

Validation of a commercial ELISA for the detection of bluetongue virus (BTV) specific antibodies in individual milk samples of Dutch dairy cows

Johannes A. Kramps, Kees van Maanen, Maria H. Mars, Johan K. Popma, Piet A. van Rijn

▶ To cite this version:

Johannes A. Kramps, Kees van Maanen, Maria H. Mars, Johan K. Popma, Piet A. van Rijn. Validation of a commercial ELISA for the detection of bluetongue virus (BTV) specific antibodies in individual milk samples of Dutch dairy cows. Veterinary Microbiology, 2008, 130 (1-2), pp.80. 10.1016/j.vetmic.2008.01.004. hal-00532386

HAL Id: hal-00532386

https://hal.science/hal-00532386

Submitted on 4 Nov 2010

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

Accepted Manuscript

Title: Validation of a commercial ELISA for the detection of bluetongue virus (BTV) specific antibodies in individual milk samples of Dutch dairy cows

Authors: Johannes A. Kramps, Kees van Maanen, Maria H. Mars, Johan K. Popma, Piet A. van Rijn


PII: S0378-1135(08)00022-9

DOI: doi:10.1016/j.vetmic.2008.01.004

Reference: VETMIC 3939

To appear in: VETMIC

Received date: 25-10-2007 Revised date: 11-1-2008 Accepted date: 17-1-2008



This is a PDF file of an unedited manuscript that has been accepted for publication. As a service to our customers we are providing this early version of the manuscript. The manuscript will undergo copyediting, typesetting, and review of the resulting proof before it is published in its final form. Please note that during the production process errors may be discovered which could affect the content, and all legal disclaimers that apply to the journal pertain.



1	validation of a commercial ELISA for the detection of bluetongue
2	virus (BTV) specific antibodies in individual milk samples of Dutch
3	dairy cows
4	
5	
6	Johannes A. Kramps ^{a,*} , Kees van Maanen ^b ,
7	Maria H. Mars ^b , Johan K. Popma ^a , and Piet A. van Rijn ^a
8	
9	
10	^a Central Veterinary Institute of Wageningen UR, P.O. box 2004,
11	8203 AA Lelystad, the Netherlands.
12	^b Animal Health Service Ltd. (GD), P.O. box 9, 7400 AA Deventer, the Netherlands.
13	
14	
15	
16	
17	* Corresponding author: Tel.: +31 320 238217; fax: +31 032 238668.
18	E-mail address: hans.kramps@wur.nl
19	
20	Corresponding address to which the proofs should be sent:
21	Central Veterinary Institute of Wageningen UR, attn. J.A. Kramps,
22	P.O. Box 2004, 8203 AA Lelystad, The Netherlands
23	e-mail address: hans.kramps@wur.nl

24	Abstract
25	
26	A recently developed indirect ELISA for the detection of bluetongue virus (BTV)-
27	specific antibodies in bovine milk samples was compared to that of the routinely used
28	competitive ELISA on serum samples.
29	During the bluetongue outbreak in the Netherlands in 2006, caused by BTV serotype
30	8, coupled serum- and milk samples were obtained from 470 individual cows from 10
31	BTV-infected farms with an average seroprevalence of 57%. In addition, bulk milk
32	samples of the same farms, and historically BT-negative samples were tested.
33	Compared to the ELISA for sera, the relative specificity and sensitivity of the ELISA
34	for milk samples is 96.5% and 98.9%, respectively when using a S/P% cut-off value
35	of 50% as advised by the manufacturer. The optimal cut-off value was found at $\ensuremath{\mathrm{S/P\%}}$
36	of 90% revealing an optimal specificity (99.0%) combined with an optimal sensitivity
37	(98.1%). Titres in positive individual milk samples ranged from 1 to 2048 with a peak
38	titre of 128. Bulk milk samples contained antibodies with titres ranging from 64 to
39	512.
40	The ELISA for milk samples was found to be a reliable and robust test. This
41	diagnostic tool is very useful, and may replace the ELISA for serum samples as first
42	choice in order to get insight into the status of lactating individual animals and
43	therewith of the entire herd with respect to BTV infection.
44	
45	
46	Keywords: Bluetongue virus; Specific antibodies; ELISA; Individual milk; Serum;
47	Sensitivity: Specificity.

Introduction

49	Bluetongue is a vector-borne viral disease of ruminants, including sheep,
50	goats and cattle, caused by bluetongue virus (BTV). Clinical signs are seen mainly in
51	sheep whereas the infection in cattle often occurs asymptomatic (Verwoerd and
52	Erasmus, 1994; Gibbs and Greiner, 1994). BTV is a member of the <i>Orbivirus</i> genus
53	of the family Reoviridae, and consists of 24 serotypes with considerable
54	immunological cross-reaction (Gorman, 1990). Animals can become infected by
55	BTV-infected midges of the genus Culicoides that have had a blood meal from
56	viraemic animals (Erasmus, 1990).
57	In 2006, an outbreak of bluetongue occurred in the Netherlands, Belgium,
58	Germany, France and Luxembourg by a BTV variant of serotype 8 (Wuijckhuise van
59	et al., 2006; Mehlhorn et al., 2007). Routine laboratory diagnosis has successfully
60	been performed in the Netherlands by an in-house real-time RT-PCR (reverse
61	transcription and polymerase chain reaction) on EDTA-blood samples (manuscript in
62	preparation) and by serology using a commercial competitive ELISA (ID.VET,
63	Montpellier, France). Positive PCR results can be found in cattle a few days post
64	infection till 200 days post infection or more. Using the ID.VET ELISA,
65	seroconversion has been detected starting 9-14 days post infection.
66	In accordance to EU legislation, control measures must be taken and disease
67	monitoring and surveillance is of fundamental importance to assess the risk posed by
68	animal movements. For dairy cattle, testing for antibodies in milk samples rather than
69	in serum samples would be very cost-effective and acceptable because of the easy
70	availability of the former sample type.
71	In this study, the performance of an indirect ELISA developed by ID.VET for
72	the detection of BTV-specific antibodies in bovine milk samples (mELISA) was

73	compared to that of the routinely used competitive ELISA on serum samples
74	(sELISA).
75	
76	Materials and Methods
77	Serum and milk BTV antibody ELISAs
78	The commercially available competitive ELISA for serum (sELISA; ID.VET,
79	Montpellier, France) was performed as described in the kit manual. Briefly, 50 μl of
80	buffer and 50µl of serum were added to the wells of a BTV-VP7 coated microtitre
81	plate. After incubation for 45 min at room temperature (rT), 100 µl anti-VP7
82	peroxidase conjugate was added and incubated for 30 min at rT. After washing, wells
83	were incubated for 15 min at rT with 100 μl TMB substrate. Colour development was
84	stopped by the addition of 100 μ I 0.5 M H ₂ SO ₄ . S/N ratios (OD _{sample} /OD _{negative control})
85	were calculated using optical density values measured at 450 nm (OD $_{450}$). S/N ratios
86	<0.4 were considered as positive.
87	
88	A recently developed, indirect ELISA for milk (mELISA) was performed according to
89	instructions of the manufacturer (ID.VET). Briefly, 50 μl skimmed milk and 50 μl
90	'wash solution' (final sample dilution 1:2) were added to the wells of a BTV-VP7
91	coated microtitre plate. After incubation for 45 min at rT, plates were washed and
92	incubated with 100 µl anti-ruminant peroxidase conjugate for 30 min at rT. After
93	washing, wells were incubated for 15 min at rT with 100 μl TMB substrate. Colour
94	development was stopped by the addition of 100 μI 0.5 M $H_2SO_4.\ S/P\%$
95	(OD _{sample} /OD _{positive control} x100%) were calculated using optical density values
96	measured at 450 nm (OD ₄₅₀). S/P% ≥50% were considered as positive.

97	Milk samples were also tested in the mELISA after a previously prepared dilution of
98	1:4 in the kit wash solution (final sample dilution 1:8).
99	Titres of individual milk samples were determined in the mELISA by testing two-fold
100	serial dilutions prepared in negative bulk milk. The titre of a milk sample was defined
101	as the highest dilution that gave a positive response (S/P% ≥50%).
102	
103	BTV-PCR
104	Viral dsRNA was automatically isolated from EDTA-blood (MagNA Pure isolation
105	robot, Roche). First strand cDNA-synthesis (RT), PCR, and real time detection was
106	performed in an one-tube, closed system (Light Cycler 2.0, Roche). The target for
107	amplification is located on genome segment 10 encoding NS3, and real time
108	detection was performed by hydrolysis of a Taqman probe. A weak positive control
109	served as cut-off for positivity for the PCR-test. This test is able to give a positive
110	signal up to 7 days before antibodies can be detected after experimental infection
111	(manuscript in preparation).
112	
113	Test Samples
114	In this study different types of samples were analysed:
115	(1) Historically negative individual milk samples (n=55) and bulk milk samples
116	(n=88) from Dutch dairy cows.
117	(2) Coupled samples (blood, EDTA blood and a milk sample) were taken
118	simultaneously of each lactating cow from 10 dairy farms selected based on a
119	high BTV seroprevalence (Table 1). In addition, from each of these 10 farms a
120	bulk milk sample was taken.

121	Milk and bulk milk samples, containing Na-azide as preservative, were defatted by
122	centrifugation and stored at -20°C in aliquots to circumvent repeated freezing and
123	thawing. Serum samples were stored at -20°C before testing. EDTA-blood samples
124	were stored at 4°C for up to 2 days before processing and testing by PCR.
125	Bulk milk taken from two BT negative farms, was used to make serial two-fold
126	dilutions of individual milk samples.
127	To test the effect of a preservative on the performance of the mELISA, either Na-
128	azide or BSM (Broad Spectrum Microtabs, D&F Control Systems Inc., San Ramon,
129	CA, USA) was added to this bulk milk.
130	
131	
132	Determination of relative sensitivity and specificity, statistical analysis
133	The relative sensitivity of the mELISA is defined as the proportion of sero-positive
134	animals giving a positive response by testing the corresponding individual milk
135	samples:
136 137 138	Seropositive animals with a positive milk sample All seropositive animals
139	The relative specificity of the mELISA is defined as the proportion of sero-negative
140	animals giving a negative response by testing the corresponding individual milk
141	samples:
142 143 144	Seronegative animals with a negative milk sample All seronegative animals
145	
146	Relative sensitivities and relative specificities of the mELISA were calculated at
147	different S/P% cut-off values. These cut-off values were used to classify milk

samples as either positive or negative. Curves were constructed by plotting the
relative specificity and the relative sensitivity against the corresponding S/P% cut-off
values ranging from 0% to 240%.
Receiver operating characteristic (ROC) curves (Kraemer, 1992) of the mELISA were
constructed by plotting the sensitivity on the ordinate as a function of the specificity at
different cut-off values. The area under the ROC curve is a quantitative measure of
the test's performance. These curves were plotted to assess the relative performance
of the mELISA when milk samples were tested undiluted or 4 times diluted.
Statistical analysis (sensitivity, specificity, confidence intervals and the agreement
index κD) was performed using Win Episcope 2.0
(http://www.clive.ed.ac.uk/winepiscope/).

Results

Specificity of the mELISA was determined by testing 55 historically negative individual milk samples and 88 historically negative bulk milk samples. All samples tested negative with a mean S/P% \pm SD of 11.0% \pm 2.8% and 10.2% \pm 0.8%, respectively. Based on the mean S/P%+3SD of these negative samples the cut-off of the mELISA should be around an S/P% of 20%.

The relative sensitivity and relative specificity of the mELISA was determined using the sELISA as the gold standard test. Coupled serum- and milk samples obtained from individual cows (n=470, see Table 1) from 10 BTV-infected farms were analysed. The overall results are shown in Table 2. The relative specificity of the mELISA is 96.5% (195/202 x 100%) with a 95% confidence interval (CI) between 94.0% and 99.1%. The 7 animals with discordant results were negative in the BTV-PCR test.

173	The relative sensitivity amounts to 98.9% (265/268 x 100%; 95% CI: 97.6% - 100%).
174	Two of the three animals with a discordant result showed a positive BTV-PCR
175	response. The overall accuracy is 97.9% (460/470 x 100%; κD=0.96).
176	A previous 4x dilution of milk samples in wash solution revealed a somewhat higher
177	specificity (97.5%; 95% CI: 95.4% - 99.7%) and a somewhat lower sensitivity (98.1%;
178	95% CI: 96.5% - 99.8%) at the prescribed cut-off SP% value of the mELISA of 50%.
179	To verify the positive/negative cut-off value of the mELISA when milk was
180	tested undiluted, a frequency distribution diagram was constructed based on the
181	results of the 470 milk samples from the 202 BT-antibody seronegative and 268 BT-
182	antibody seropositive animals (see Figure 1). In addition, relative sensitivity and
183	relative specificity at different SP% cut-off values were calculated and plotted (see
184	Figure 2). Highest relative sensitivity (98.1%; 95% CI: 96.5% – 99.8%), combined
185	with highest relative specificity (99.0%; 95% CI: 97.5% - 100%) was observed at a
186	cut-off SP% value of 90% (see also figure 3).
187	ROC curves were constructed to compare the performance of the mELISA
188	when milk samples were tested either undiluted or 4x diluted in wash solution (Figure
189	3). The test with undiluted milk samples shows a slightly greater area under the ROC
190	curve compared with the test performed with 4x diluted samples indicating a
191	somewhat better performance with undiluted milk samples.
192	To determine the BTV-specific antibody titres in the individual milk samples,
193	two-fold serial dilutions in negative bulk milk were tested in the mELISA. Typical dose
194	response curves of 8 serially diluted milk samples (1 mELISA negative and 7
195	mELISA positive samples) are depicted in Figure 4. Similar curves were found when
196	dilutions were made in bulk milk containing either azide or BSM as preservative.
197	Figure 5 depicts the frequency distribution diagram of the antibody titres of the 470

198	individual milk samples. Milk samples from seropositive animals showed a mean 2log
199	titre of 7. The 7 positive milk samples from seronegative animals had 2log titres of 0
200	(n=5) and of 1 (n=2).
201	The bulk milk samples collected from the 10 BT infected farms showed 2log
202	titres of 6 (n=5), 7 (n=4) and 9 (n=1) in the mELISA. The bulk milk sample with the
203	highest titre (29) originated from the farm with the highest seroprevalence (95%).
204	
205	Discussion
206	The commercial ID.VET mELISA is a recently developed indirect ELISA able
207	to detect BTV-VP7 specific antibodies in milk samples.
208	The performance of the mELISA was investigated using the ID.VET sELISA as the
209	reference test. In a previously performed study, the ID.VET sELISA showed both a
210	high specificity (100%) and a high sensitivity (98%) relative to BTV-PCR results on
211	field samples obtained at the early stage (Sept. and Nov. 2006) of the BT outbreak in
212	the Netherlands (manuscript in preparation). Moreover, the ID.VET sELISA showed a
213	good performance in proficiency testing (Workshop on bluetongue diagnostics and
214	epidemiology, Brussels, Belgium, November 28-29, 2006). In our laboratory (CIDC-
215	Lelystad), the sELISA is an ISO17025 accredited test and used for routine serum
216	testing of both cattle and sheep.
217	In the study described it was investigated whether the mELISA can be used as
218	a tool to identify BT-seropositive dairy cows. Analysis of historical negative individual
219	(n=55) and bulk milk samples (n=88) indicated a specificity of 100% in this small set
220	of data.
221	Testing individual milk samples obtained from seropositive (n=268) and seronegative

cattle (n=202) revealed a high sensitivity and a high specificity relative to serology

testing. Of the 470 milk samples tested in the mELISA, only 10 discordant results
(=2.1%) with the sELISA were found at the prescribed positive/negative cut-off S/P%
value of 50%. From the results as shown in Figure 2 it can be concluded that the test
shows a high sensitivity and a high specificity that change minimally between cut-off
values of 30% and 110%. Highest sensitivity combined with highest specificity
revealed 7 discordant results (5 false negatives and 2 false positives) and was found
at a cut-off S/P% value of 90%. A 4x dilution of milk samples before testing did not
result in a better performance (Figure 3): highest accuracy revealed 9 discordant
results (6 false negatives and 3 false positives) at a cut-off value of 80%.

Samples investigated were from farms with a high seroprevalence ranging from 34% to 95%. These farms were not very recently infected (samples were taken probably more than 4 months after primary infection), and most seropositive animals have reached there maximum antibody levels resulting in a minimum number of false negative test results. In case of a recently infected farm, at least part of the sera taken from infected animals may contain relatively low antibody titres giving weaker responses in the ELISA. For this reason a positive/negative cut-off value should be chosen resulting in a test with an optimal sensitivity combined with a specificity producing an acceptable low number of false positives (high sensitivity combined with an acceptable good specificity). As can be deduced from the results of this study and presented in Figures 1 and 2 using a S/P% cut-off value of 50% resulted in 3 false negatives (sensitivity = 98.9%), two of which were confirmed by a positive BTV-PCR and 7 false positives (specificity = 96.5%), all of which were confirmed by a negative BTV-PCR. In view of the impact of these false positives, one may decide to confirm (false) positive milk test results by serum antibody testing and/or by BTV-

PCR. When a higher specificity is required, the cut-off value can be cha	anged to a
S/P% of 90% with minimal effects on sensitivity.	

Antibody levels in the individual milk samples were determined by testing twofold serial dilutions, prepared in negative milk, in the mELISA. Figure 5 shows a titre
distribution of positive milk samples between 2° and 2¹¹ with a peak at 2⁷. Of the 8
milk samples with the lowest titre (2° and 2¹), 7 showed discordant results with the
serum test (positive milk samples from seronegative animals). This suggests that
animals develop BTV-specific antibody mELISA titres of 2² or higher and that by bulk
milk detection of a herd prevalence of 25% or lower seems to be feasible. BTVspecific antibody titres of the bulk milk samples of the 10 farms ranged from 2⁶ to 2⁹.
A weak positive correlation was observed between the seroprevalence and the bulk
milk titre showing the highest titre in the bulk milk sample obtained from the farm with
the highest seroprevalence. The application of the mELISA for bulk milk testing is
now being validated by testing bulk milk samples from BTV infected farms with
different percentages of seropositive lactating animals.

Due to immunological cross reaction with BTV-VP7 (Gorman, 1990) it can not be excluded that the presence of antibodies directed against BTV-related Orbiviruses like epizootic haemorrhagic disease virus (EHDV) may result in positive responses in the mELISA and/or sELISA. In case animals are tested positive in the mELISA one should consider to further investigate these animals by more specific serum tests and/or PCR. However, up to this moment EHDV never has been detected in this part of Europe.

We conclude that the mELISA using undiluted milk at cut-off S/P percentages between 30% and 110% is a very reliable diagnostic tool. Because milk samples are easy and cost-effective to obtain this could be the first choice of testing in order to get

272	insight to the status of animals and therewith of the herd concerning bluetongue virus
273	infection.
274	
275	References
276	Erasmus, B.J., 1990. Bluetongue virus. In: Dinter, Z., Morein, B., (Eds.), Virus
277	Infections of Ruminants, Elsevier, New York, USA, pp 227 – 237.
278	Gibbs, E.P., Greiner, E.C., 1994. The epidemiology of bluetongue. Comp.
279	Immunol. Microbiol. Infect. Dis. 17, 207 – 220.
280	Gorman, B.M., 1990. The bluetongue viruses. Curr. Top. Microbiol. Immunol.
281	162, 1 – 19.
282	Kraemer, H.C., 1992. Evaluating medical tests, objective and quantitative
283	guidelines. London, Sage Publications, pp 63 – 95.
284	Mehlhorn, H., Walldorf, V., Klimpel, S., Jahn, B., Jaeger, F., Eschweiler, J.,
285	Hoffmann B., Beer, M., 2007. First occurrence of Culicoides obsoletus-
286	transmitted Bluetongue virus epidemic in Central Europe. Parasitol. Res. 101,
287	219 – 228.
288	Verwoerd, D.W., Erasmus, B.J., 1994. Bluetongue. In: Coetzer, J.A.W.,
289	Thomas, G.R., Tustin, R.C. (Eds.), Infectious Disease of Livestock with
290	Special Reference to Southern Africa. Oxford University Press, Capetown
291	(Afrique du Sud), pp 443 – 459.
292	Wuijckhuise van, L., Dercksen, D., Muskens, J., Bruijn de, J., Scheepers, M.,
293	Vrouenraets, R., 2006. Bluetongue in The Netherlands; description of the first
294	clinical cases and differential diagnosis. Common symptoms just a little
295	different and in too many herds. Tijdschr Diergeneeskd. 131, 649 - 654.

296 **Table 1.**

297

Numbers of samples and anti-bluetongue seroprevalence of farms analysed in this

299 study.

300

Farm nr.	Number of lactating cows	s Seroprevalence	
1	50	94% (47/50)02	
2	25	56% (14/25) 03	
3	29	41% (12/29)04	
4	27	67% (18/273)05	
5	55	49% (27/55)06	
6	42	50% (21/42)07	
7	22	95% (21/2 2)08	
8	76	67% (51/76)09	
9	97	42% (41/973) 10	
10	47	34% (16/473)11	
All	470	57% (268/470)2	
		313	


Table 2.
 Comparison of results obtained with serum ELISA and milk ELISA of 470 coupled
 serum-milk samples.

320


		Sera tested in sELISA		321 322 323 324
		Positive (S/N<0.4)	Negative (S/N≥0.4)	325 Tot <u>326</u>
Milk samples	Positive (S/P%≥50%)	265	7	328 329 27 <u>2</u> 30 331 332
tested undiluted in mELISA	Negative (S/P%<50%)	3	195	333 334 19 8 35 336 337
	Total	268	202	338 47039 340 341

345	Figure captions
346	
347	Figure 1.
348	Frequency distribution diagram of 470 milk samples obtained from seropositive and
349	seronegative animals originating from 10 Dutch BTV infected herds. S/P% represent
350	test results obtained with the ID.VET mELISA.
351	
352	Figure 2.
353	Relative sensitivity and relative specificity of the mELISA at different S/P% cut-off
354	values.
355	
356	Figure 3.
357	ROC curves based on mELISA test results obtained with undiluted milk and milk 4x
358	diluted in buffer.
359	
360	Figure 4.
361	Typical dose response curves of serially diluted milk samples (7 positives and 1
362	negative) tested in the mELISA.
363	
364	Figure 5.
365	Frequency distribution of 2log(titres) of 470 individual milk samples collected from
366	seropositive and seronegative animals of 10 BTV-infected farms.






4/30




4/30














4/30