



**Figure S2, supplementary data. Analysis of CWPs
by SDS-PAGE after separation
by multidimensional lectin chromatography**

- A. Chromatography 1 and 2 on AIL (EA fractions) and PNA (EP fractions)
B. Chromatography 3 and 4 on ConA (EC fractions) and WGA (EW fractions)

M (kDa) Molecular mass markers; CWPs: cell wall proteins; EA(1,2): fractions eluted from the AIL column; EP(1,2): fractions eluted from the PNA column; EC(1,2): fractions eluted from the ConA column; EW (1,2): fractions eluted from the WGA column; FT4: flow-through of the WGA column. Numbering refers to the stained bands analyzed by MALDI-TOF MS (*cf* Table S2, supplementary data)