

IgA PURIFICATION FROM SOW MILK

Serge Bernard, Isabelle Lantier

▶ To cite this version:

Serge Bernard, Isabelle Lantier. IgA PURIFICATION FROM SOW MILK. Annales de Recherches Vétérinaires, 1986, 17 (1), pp.85-88. hal-00901614

HAL Id: hal-00901614 https://hal.science/hal-00901614

Submitted on 11 May 2020

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés. Ann. Rech. Vét., 1986, 17 (1), 85-88

Short note

IgA PURIFICATION FROM SOW MILK

S. BERNARD and Isabelle LANTIER

Laboratoire de Pathologie Porcine, INRA, 37380 Nouzilly, France

Résumé

PURIFICATION DES IGA DU LAIT DE TRUIE. — Les IgA du lait de truie ont été facilement purifiées par des techniques simples de chromatographie d'échange d'ions (DEAE Trisacryl M) et de filtration sur gel (Ultrogel AcA 22).

The secretory IgA is the predominant immunoglobulin in external secretions. This immunoglobulin seem very important for the protection against bacterial and viral aggression (Offit, 1985; Bohl and Saif, 1975; Tabliabue *et al.*, 1983). The purification of this class of immunoglobulin is very important to make reactif. The IgA immunoglobulin have been greatly studied, and purified with cellulose matrix ion exchange chromatography (Porter, 1969; Bourne, 1969; Debuysscher *et al.*, 1975). The use of this technique is laborious. The interaction of porcine IgA with protein A gives good yield but does not eliminate the IgG class (Bennel and Watson, 1980).

Recently, synthetic matrix (Trisacryl) have been used to easily purify the IgG class of immunoglobulin (Corthier *et al.*, 1983).

This communication will describe a very simple way to prepare IgA from milk by classical precipitation and chromatographic technique.

Materials and Methods

1) Chemicals

The ultrogel AcA 22 (LKB) and the DEAE-Trisacryl M (LKB) were purchased from reactif IBF (France); the peroxidase (grade 1) from Boehringer; the human albumin (fraction V) from Sigma.

All the other chemicals and the buffer were of an analytical grade.

2) Milk

Milk samples from various sows were manually collected after ocytocine intramuscular inoculation (2 units) and immediately used or stored at -20 °C.

3) Milk globulin extraction

We collected 45 ml of milk. The cream was separated from the skimmed milk by low-speed centrifugation (2500 g, 10 min, 5 °C, Heraeus Christ). The caseins were removed by high-speed centrifugation 45000 g on Beckman Spinco L8-55 rotor R 30. The whey proteins

| | Volume (ml) | Proteinsa | | lgAb | | Yield (%) | |
|--------------------|-------------|--------------------------|------------------|--------------------------|------------------|-----------|------|
| | | Concentration (mg/ml) | Quantity (mg) | Concentration (mg/ml) | Quantity (mg) | Protein | lgA |
| Milk | 45 | 62 | 2 790 | 14.1 | 635 | _ | _ |
| Whey | 43 | 28 | 1 200 | 7.7 | 331 | 43 | 52 |
| 50 % precipitation | 40 | 20 | 800 | 6.4 | 256 | 28.5 | 40.3 |
| DEAE start | 25 | 15 | 375 | 4.5 | 112 | 13.5 | 17.6 |
| A.C.A. 2.2. start | 10 | 15 | 150 | 6.8 | 68 | 5.4 | 10.7 |
| IgA AcA pic | | | _ | _ | 51.5° | 1.8 | 8.1 |

Table 1. - Yield of each step of IgA purification

a: The protein concentration was estimated by the Lowry method.

b: The IgA concentration was estimated by the ELISA technique described in the Materials and Method.

c: After lyophilization.



Fig. 1. — Protein elution from DEAE Trisacryl M (LKB): 25 ml of soluble material was applied to a column 1.6 × 10 cm in 0.025 M Tris-HCL buffer pH 8, 0.05 M NaCl, 90 ml/h. Two elution steps were made: 0.14 M and 1 M.



Fig. 2. — Gel filtration of the IgA DEAE eluted pic on AcA 22 Ultrogel (LKB) 10 ml of material was applied to a column (4.5 x 75 cm) in Tris-HCI 0.025 M, pH 8, 0.05 M NaCl.

(43 ml) were used immediately. The gammaglobulin fraction was obtained by ammonium sulphate precipitation (50 % of saturation V/V) at 4 °C for 2 h. After centrifugation, the pellet was resuspended in 40 ml of 0.025 M Tris-HCl buffer, pH 8.8 and reprecipitated in 33 % of ammonium sulphate concentration. The pellet was collected in 25 ml of 0.025 M Tris-HCl, pH 8.8 containing 0.05 M NaCl and extensively dialysed at 4 °C against this same buffer.

4) IgA globulin purification

The ion exchange chromatography was carried out on a column (1.6 \times 10 cm) of DEAE-Trisacryl M, equilibrated in 0.025 M Tris-HCl buffer, pH 8.8, NaCl 0.05 M. The gel filtration was performed on a column (4.5

 \times 75 cm) of AcA 22 ultrogel in the same buffer. The flow rate was 40 ml/h.

The chromatographic experiments were continuously recorded with ISCO at 280 nm (model UA 5 absorbance monitor), and the eluate affluent collected with ISCO fraction collector.

The protein concentration was determined by the Lowry technique (Lowry *et al.*, 1953) with human albumin as reference.

5) Immunoelectrophoresis (IEF)

Immunoelectrophoresis in 1 % agar was carried out according to standard procedures with total anti-swine serum proteins and anti-swine skim-milk protein antibodies (rabbits).



Fig. 3. — Immunoelectrophoresis of the different steps of IgA purification. Anti-milk antibodies; specific anti-IgA; 1: milk; 2: whey; 3: supernatant of the 50 % ammonium sulphate precipitation; 4: 50 % ammonium sulphate precipitate; 5: DEAE starting material (33 % ammonium sulphate precipitate); 6: AcA 22 starting material; 7: IgA pic on the AcA 22.

Results and Discussion

The procedure described here allows the preparation of highly purified IgA with a good yield.

The IgA purification is subordinate to the ratio between the IgA and the other components of the starting material. The swine milk contains from 10 to 15 mg/ml of IgA, and the ratio IgA/IgG is very high (about two times as much or more).

Table 1 gives the protein yield of the different steps of IgA purification. The first purification step after skimming is to separate the globulin fraction from the other components. Half of the IgA seemed trapped in the casein pellet. After the first step of ammonium sulphate precipitation (50 % V/V) we lost half of the proteins, and 60 % of IgA. The second precipitation dramatically decreased the IgA level. A better yield of IgA purification could be obtained by avoiding the second ammonium sulphate precipitation which eliminates many proteins but also eliminates a fair quantity of IgA.

The globulin preparation was passed through the DEAE gel. The elution profile (fig. 1) shows several peaks. The unadsorbed fraction is mainly constituted of IgG (Corthier *et al.*, 1984). The second pic eluted with one step gradient of 0.14 M NaCl contains IgA. After this step the purification of this material needs an exclusion chromatography (fig. 3). This fraction was concentrated to 10 ml by ammonium sulphate precipitation and layered on the top of the ultrogel column. The elution profile is shown in figure 2. The elution behaviour of the major peak (Kav = 0.41) corresponds to a molecular mass of approximately 400 000 daltons.

The proteins were concentrated by ammonium sulphate precipitation (50 % V/V). The IEP analysis gave only one precipitating arc with either the IgA specific antiserum (kindly given by J.J. Metzger) or the anti-milk antibodies (fig. 3).

The acrylamide gel electrophoresis (Dubray and Bezard, 1982) of this IgA preparation gives only three bands for the H- and L-chains and for the secretory piece (results not shown). The yield of the purification could be improved more but the purification of this IgA by this method rapidly gives IgA which only gives three bands in polyacrylamide gel.

The ELISA reactions of this immunoglobulin against specific swine anti- IgG or IgM, are very low (personal results) and the estimated contamination less than 1 %.

Received, 15th May 1985 Accepted, 13nd June 1985

Summary

IgA from swine milk may be easily purified with DEAE-Trisacryl M ion echange chromatography (upper than 75 % of purity) and AcA 22 Ultrogel gel filtration.

References

- BENNEL M.A., WATSON D.L., 1980. The interactions of porcine and ovine, serum and colostral immunoglobulins with staphylococcal protein A. *Microbiol. Immunol.*, **24**, 871-878.
- BOHL E.H., SAIF L.J., 1975. Passive immunity in TGE of swine immunoglobulin characteristics of antibodies in milk after inoculating virus by different routes. *Infect. Immun.*, **11**, 23-32.
- BOURNE F.J., 1969. IgA immunoglobulin from porcine milk. Biochem. Biophys. Acta, 181, 485-493.
- CORTHIER G., BOSCHETTI E., CHARLEY-POULAIN J., 1984. Improved method for IgG purification from various animal species by ion exchange chromatography. J. Immunol. Met., 66, 75-79.
- DE BUYSSCHER E.V., BERMAN D.T., 1975. Isolation of porcine secretory IgA by affinity chromatography and determination of its component chains. Am. J. Vet. Res., 36, 1659-1661.
- DUBRAY G., BEZARD G., 1982. A highly sensitive periodic acid silver stain for 1,2 diol groups of glycoproteins and polysaccharides in polyacrylamide gels. *Analyt. Biochem.*, **119**, 325-329.
- LOWRY O.H., ROSEBROUGH N.J., FARR A.L., RANDALL R.J., 1951. Protein measurement with the Folin phenol reagent. J. Biol. Chem., 193, 265-275.
- OFFIT P.A., CLARK F., 1985. Protection against rotavirus induced gastroenteritis in a murine model by passively acquired gastrointestinal but not circulating antibodies. J. Virol., 54, 58-64.
- PORTER P., 1969. Porcine colostral IgA and IgM antibodies to *E. coli* and their intestinal absorption by the neonatal piglet. *Immunology.*, **17**, 617-626.
- TABLIABUE E., NENCIONI L., VILLA L., KEREN D.F., LOWELL G.H., BORASHI D., 1983. Antibody-dependent cellmediated antibacterial activity of intestinal lymphocytes with secretory IgA. *Nature.*, 306, 184-186.