



HAL
open science

A substrate sensor chip to assay the enzymatic activity of Botulinum neurotoxin A

Christian Lévêque, Géraldine Ferracci, Yves Maulet, Chloé Grand-Masson, Marie-Pierre Blanchard, Michael Seagar, Oussama El Far

► **To cite this version:**

Christian Lévêque, Géraldine Ferracci, Yves Maulet, Chloé Grand-Masson, Marie-Pierre Blanchard, et al.. A substrate sensor chip to assay the enzymatic activity of Botulinum neurotoxin A. *Biosensors and Bioelectronics*, 2013, 49, pp.276-281. 10.1016/j.bios.2013.05.032 . hal-00958814

HAL Id: hal-00958814

<https://hal.science/hal-00958814v1>

Submitted on 13 Mar 2014

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

A substrate sensor chip to assay the enzymatic activity of Botulinum neurotoxin A

Christian Lévêque^{a,c,*}, Géraldine Ferracci^{b,c,&}, Yves Maulet^{a,c,&}, Chloé Grand-Masson^{b,c},
Marie-Pierre Blanchard^{b,c}, Michael Seagar^{a,c}, and Oussama El Far^{a,c}

a. INSERM, UMR_S 1072, 13015 Marseille, France

b. CNRS, UMR 7286, Plate-Forme de Recherche en Neurosciences PFRN, 13015 Marseille, France

c. Aix-Marseille Université, 13015 Marseille, France

&. The authors contributed equally to this work

*Address correspondence to this author at: INSERM 1072, Faculté de médecine secteur nord, Aix-Marseille Université CS80011 Bd Pierre Dramard 13344 Marseille Cedex 15 France. Fax +33-491090506; e-mail christian.leveque@univ-amu.fr

ABSTRACT

Botulinum neurotoxin A (BoNT/A) induces muscle paralysis by enzymatically cleaving the presynaptic SNARE protein SNAP-25, which results in lasting inhibition of acetylcholine release at the neuromuscular junction. A rapid and sensitive *in vitro* assay for BoNT/A is required to replace the mouse lethality assay (LD₅₀) in current use. We have developed a fully automated sensor to assay the endoprotease activity of BoNT/A. We produced monoclonal antibodies (mAbs) that recognize SNAP-25 neo-epitopes specifically generated by BoNT/A action. Recombinant SNAP-25 was coupled to the sensor surface of a surface plasmon resonance (SPR) system and samples containing BoNT/A were injected over the substrate sensor. On-line substrate cleavage was monitored by measuring binding of mAb10F12 to a SNAP-25 neo-epitope. The SNAP25-chip assay was toxin serotype-specific and detected 55 fM BoNT/A (1 LD₅₀/ml) in 5 minutes and 0.4 fM (0.01 LD₅₀/ml) in 5 hours. Time-course and dose-response curves were linear, yielding a limit of quantification of 0.03 LD₅₀/ml. This label-free method is 100 times more sensitive than the mouse assay, potentially providing rapid read-out of small amounts of toxin for environmental surveillance and the quality control of pharmaceutical preparations.

Nonstandard abbreviations: BoNT, botulinum neurotoxin; SPR, surface plasmon resonance; RU, resonance unit; BSA, bovine serum albumin; LD₅₀, median lethal dose; HEPES, 4-(2-hydroxyethyl)-1-piperazine ethane sulfonic acid; SNAP-25, Synaptosomal-associated protein 25; DTT, Dithiothreitol; HBS, HEPES-buffered saline; LOD, limit of detection; LOQ, limit of quantification.

Key words: botulinum neurotoxin; SNAP-25; surface plasmon resonance; toxin sensor; endoprotease; neo-epitope.

1. Introduction

Neurotoxins from *Clostridium botulinum* (BoNTs) act at the neuromuscular junction and cause the human disease botulism. BoNTs display endoprotease activity specific for the SNARE proteins that support Ca^{2+} -dependent exocytosis at the nerve terminal. Thus, cleavage of SNAREs by BoNTs inhibits acetylcholine release and compromises muscle contraction. Most cases of botulism occur through ingestion of contaminated food, with the toxin reaching peripheral nerve terminals via the circulation. Although human botulism is relatively rare, the intoxicated patient needs prolonged and extensive clinical support to counteract descending flaccid paralysis, which can lead to death by respiratory failure if untreated.

Seven different BoNT-serotypes are known, but four (BoNT/A, B, E and F) predominate in human botulism. BoNT/A is the most frequent toxin serotype detected in human botulism in the USA (Arnon et al., 2001). It is also the most potent and can persist in an active form in neurons for weeks or months. BoNT/A was the first BoNT to be introduced for clinical use in the 1980s. Since then, its indications have expanded steadily and BoNT/A is currently used therapeutically to inhibit hyper-contraction or hyper-secretion in a large range of human diseases, as well as cosmetically to reduce facial wrinkles. Moreover BoNT/A is a class A agent that constitutes a potential biological weapon (Dhaked et al., 2010).

BoNT/A occurs naturally as a high molecular mass complex of up to 900 kDa, composed of the neurotoxin (150 kDa) and a set of non-toxic components of Clostridial origin. The neurotoxin has two polypeptide chains, linked by a single disulfide bond and non-covalent interactions. Each chain has a distinct function. *In vivo* the heavy chain (100 kDa) vectors the toxin to motoneuron nerve terminals, via recognition of a specific receptor (synaptic vesicle protein 2 or SV2, associated with gangliosides). Binding results in endocytosis of the toxin followed by translocation of the light chain into the cytoplasm, concomitant with disulfide bond reduction. The light chain (50 kDa) possesses zinc-dependent endopeptidase activity, which specifically cleaves 9 amino acids from the C-terminus of the SNARE protein SNAP-25 (synaptosome-associated protein of 25 kDa). SNAP-25, which is anchored at the cytoplasmic surface of the synaptic plasma membrane, is a key component of the synaptic vesicle fusion machinery (Rossetto et al., 2013).

The standard method, used currently for serological diagnosis of botulism, surveillance and quality control of batches during pharmaceutical production, is the mouse toxicity bioassay, which has a detection limit of 1 LD_{50}/ml (10-20 pg BoNT/A). However this procedure is

increasingly considered as unethical and furthermore is relatively slow (24-48 h), as the symptoms of botulism in the mouse can take several days to appear. Thus the mouse assay is poorly adapted to emergency situations requiring rapid diagnosis and appropriate treatment.

Numerous *in vitro* assays have been proposed as substitutes for the mouse lethality assay, and these fall into two main categories. The first only detects the neurotoxin protein using specific antibodies (ELISA-based predominantly), or the nucleic acids encoding it (PCR, immuno-PCR). The second includes the methods that assay BoNT/A biological activity (for reviews, Capek and Dickerson, 2010; Hakami et al., 2010; Lindstrom and Korkeala, 2006) .

SPR is a powerful tool for studying biomolecular interactions. The versatility of SPR based-biosensors extends beyond protein-protein interactions and environmental and medical applications are expanding (Campbell et al., 2011; Helmerhorst et al., 2012). A few studies have used this label-free optical detection technique to measure protease activity using surface-captured substrate (Loew et al., 2000; Shoji et al., 2011; Steinrucke et al., 2000). However in all these reports, protease activity in the fluidic phase was monitored by a loss of substrate mass, leading to a decrease in SPR signal without revealing the site of cleavage.

In this report we describe a new simple and automated method to measure the enzymatic activity of BoNT/A entirely on-chip using SPR. On-chip cleavage of immobilized recombinant SNAP-25 was assayed by detecting neo-epitopes specific to BoNT/A endopeptidase activity. This provides a rapid and sensitive *in vitro* assay that quantifies BoNT/A activity on-line, by a gain-of-SPR signal.

2. Material and methods

2.1 Reagents: Amine coupling kits were from GE Healthcare and Biorad. ZnCl₂, Tween 20, dithiothreitol (DTT), Fc specific goat anti-mouse IgG and fatty acid-free bovine serum albumin were from Sigma. Botulinum neurotoxins: frozen aliquots of BoNT/A (Metabionics, USA) at 3.5×10^7 LD₅₀/mg (1 mg/ml from Hall strain) or BoNT/B (Metabionics, USA) at 1.1×10^7 LD₅₀/ml (1 mg/ml from Okra strain) were used immediately after thawing. BoNT/E and C (provided by Prof. S Kozaki, Osaka Prefecture University) were produced and purified as previously described (Ferracci et al., 2005). BoNT/A concentrations are indicated as molarity or mouse LD₅₀/ml. All experiments were performed in accordance with French and European guidelines for handling botulinum neurotoxin.

2.2 Monoclonal antibody production: Mice were immunized with a synthetic peptide identical to the 191-197 sequence of rat SNAP-25. This sequence is located immediately N-terminal to the peptide bond cleaved by BoNT/A (Q₁₉₇-R₁₉₈). A cysteine residue was added to the N-terminus of the peptide to allow coupling to KLH protein for immunization. Hybridomas were selected by ELISA against the coated peptide and for high affinity binding to SNAP-25 in plasma membrane vesicles treated with BoNT/A (Marconi et al., 2008). Monoclonal antibodies (mAb10F12) were purified by Protein A-Sepharose affinity chromatography and stored at -20°C.

2.3 Preparation of SNAP-25 chips: Rat 6x His-SNAP-25 was cloned in pRSFDuet (Novagen). His-SNAP-25 was expressed in BL21 and purified on Ni-NTA-agarose (Qiagen) beads under denaturing conditions according to the supplier's protocol and dialyzed against HBS (10 mM HEPES-NaOH pH 7.4, 150 mM NaCl). Purified proteins were stored in 20% glycerol at -80°C for a maximum of 6 months. SNAP-25 was amine coupled at pH 4.5 to Biacore CM5 (GE Healthcare) or pH 4 to ProteOn GLH (Biorad) sensor chips according to the manufacturer's instructions using automated programs. The final immobilization density was 9000 RU (CV = 4%, n= 15) on CM5 chips and 17000 RU (CV = 3.2%, n= 6 along a channel) on GLH chips, yielding a concentration of 10-30 µM. Chips were conditioned by injection of cleavage buffer (50 mM HEPES-NaOH pH 7.4, 0.1% BSA, 10 mM DTT, 1% Tween 20, 100 µM ZnCl₂, 10 mM DTT) and background signal was measured from five injections of mAb10F12 diluted in running buffer in the presence of anti-mouse Fc antibody, each at 10 µg/ml.

2.4 SPR analysis of BoNT/A endoprotease activity: SPR measurements were performed at 34°C using Biacore 3000, Biacore T200 (GE Healthcare) or ProteOn XPR36 (Biorad) apparatuses with TBS (10 mM Tris-HCl pH 7.4, 150 mM NaCl) or HBS as running buffer. The minimum flow rates were 1 µl/min with the Biacore and 25 µl/min with the ProteOn apparatus. BoNT/A samples were serially diluted in cleavage buffer and placed in the sample rack at 4°C before injection. In all experiments, cleavage buffer in the absence of BoNT/A was injected onto a control lane coated with SNAP-25 (control flow cell or inter-spot area for Biacore or Biorad apparatus respectively). The pre-assembled mAb10F12/anti-mouse Fc antibody complex was then injected for 2 min over all flow cells. MAb binding was measured 10 sec after the end of injection, and the chip was regenerated with an 8 sec pulse of 10 mM glycine-HCl pH 2. The data obtained from the control flow cell were automatically subtracted from experimental measurements to yield the specific signal. The pre-assembled mAb10F12/anti-mouse Fc antibody complex yielded a 3-fold increase in detection sensitivity compared to mAb10F12 alone. Data were analyzed using Biaevaluation 4.1 software (Biacore/GEHealthcare) and ProteOn Manager v. 3.1.0.6 software (Biorad). Assuming a molecular weight ratio of 6:1 between mAb and SNAP-25 and that one mAb molecule binds

to only one epitope at low rates of substrate cleavage, one RU of bound mAb10F12 corresponds to the cleavage of 0.015% of coupled substrate. The limit of detection (LOD) was measured in the absence of BoNT/A, and defined as the mean of the background plus three times the standard deviation. The limit of quantification (LOQ) was defined as the mean of the background plus ten times the standard deviation. SPR chips were used with at least 40 assay-regeneration cycles without significant loss of binding capacity. Chips were discarded when more than 10% of the substrate was cleaved.

3. Results

3.1 Production of monoclonal antibodies recognizing the reaction product of BoNT/A endopeptidase activity

The proteolytic activity of BoNT/A specifically cleaves a single peptide bond (Gln₁₉₇-Arg₁₉₈) of SNAP-25, removing a nine residue (198-206) C-terminal fragment. In order to assay BoNT/A activity using a substrate chip, we produced antibodies designed to bind specifically the C-terminal neo-epitope of SNAP-25 1-197, generated by endopeptidase activity. Mouse monoclonal antibodies were raised against a synthetic peptide corresponding to residues 191-197 of SNAP-25 and screened by ELISA, Western blotting and SPR. Six positive hybridoma clones were obtained, producing mAbs which reacted with SNAP-25 in Western blots of rat brain membrane proteins treated with BoNT/A (Fig. 1a). These mAbs displayed no cross-reactivity with full-length SNAP-25 in untreated membranes. Clone 10F12 was selected for further study on the basis of its high affinity. The specificity of mAb10F12 for SNAP-25 cleaved by BoNT/A versus full length SNAP-25 was also confirmed by SPR experiments using the native protein in membrane fractions (Fig. 1b), as described in (Marconi et al., 2008) and by immunocytochemical staining of cultured neurons incubated in the presence or absence of BoNT/A (Fig. S1).

3.2 Assay principle

Biosensors coupled to microfluidics flow cells work on the principle of continuous flow over a target coupled to the sensor surface. Samples are injected into the running buffer, at and for a given time. In this procedure, the target is the BoNT/A substrate, purified recombinant SNAP-25 immobilized on the reaction area called “flow cell” or “spot”. The method incorporates two injection steps: (i) a cleavage step in which the sample containing BoNT/A is injected over the substrate. This clips off a nine residue C-terminal peptide and generates a neo-epitope specifically recognized by mAb10F12. (ii) A detection step in which mAb10F12 is injected. MAb10F12 binding to cleaved SNAP-25 produces an SPR signal in real time

which is directly proportional to the amount of proteolytic product and thus to the concentration of BoNT/A in the initial sample (Fig. 2).

The duration of the sample injection step (ie. enzyme / substrate contact time) was varied between 5 and 300 min according to the BoNT/A concentration range to be assayed, whereas the detection step was set for 2 minutes. As the SPR signal generated is proportional to the molecular mass of bound molecules, mAb10F12 was pre-mixed with an anti-mouse Fc antibody to increase sensitivity.

3.3 Assay specificity

To confirm the specificity of BoNT/A activity detection, other BoNT serotypes were injected over SNAP-25 coated on a Biacore 3000 sensor chip at a relatively high concentrations (1 nM). No mAb10F12 binding signal was observed after BoNT/B injection (Fig. 3a), which cleaves the synaptic vesicle protein synaptobrevin / VAMP2. Similarly, BoNT/E and BoNT/C, which cleave SNAP-25 at different sites than BoNT/A (Rossetto et al., 2013), also failed to produce any signal (Fig. 3a). In contrast, injection of BoNT/A at a 1000-fold lower concentration (1 pM) generated a strong signal. These data indicate that the assay is specific for the endoprotease activity of BoNT/A, and qualifies this antibody to identify this serotype in an unknown sample.

3.4 Cumulative measurements of BoNT/A endoprotease activity

The SNAP-25 coupling procedure ensures that excess substrate ($> 5 \mu\text{M}$) is available. Experiments were systematically carried out either with a low enzyme concentration, or with limited contact time, to ensure that less $< 10\%$ of substrate was cleaved. In these conditions, substrate concentration was considered to be constant during a single or a set of toxin injections in the same flow-cell. Signal repeatability, measured by at least 3 mAb injections, separated by acid regeneration of the sensor chip, was high (CV $\sim 1\%$ from 4 to 5000 RU, and 4% for 2 RU). Hence in principle several measurements can be performed in the same flow cell by monitoring successive increments in signal. Fig. 3b illustrates 7 consecutive injections of the same toxin concentration for 5 min, followed in each case by an independent mAb10F12 injection. Product formation over each 5 min period proceeds at a linear rate indicating that the reaction is not substrate-limited (Fig. 3c). At a fixed reference time this rate will depend linearly on enzyme concentration (see below).

Endoprotease activity were found to be independent of the buffer flow rate over the sensorchip between the limits of 1 and 25 $\mu\text{l}/\text{min}$ (Fig. S2), indicating that diffusional transport of the enzyme to the solid-state substrate is non-limiting compared to the enzymatic reaction.

3.5: Rapid detection of BoNT/A activity

The rapidity of BoNT/A detection is of major importance in surveillance applications, thus the potential of SPR to achieve this aim was evaluated at 1 LD₅₀/ml, a concentration that corresponds to the limit of detection for the *in vivo* mouse assay obtained in 24 h. In order to improve measurement sensitivity, experiments were carried out on a Biacore T200 apparatus which provides very low background. As illustrated in Fig. 4, specific endoprotease activity was readily detected after a 5 min injection of reduced BoNT/A, yielding 4.5 ± 1 RU (n=3), a value largely above the LOD (dashed line). Alternatively, non-reduced BoNT/A was diluted in cleavage buffer containing DTT and immediately injected. In these conditions, although all BoNT/A was not fully reduced and thus active, a cleavage time of 15 min was sufficient to clearly detect BoNT/A activity yielding 10.7 ± 1 RU (n=4, data not shown). These data indicate that our procedure provides a sensitivity of 1 LD₅₀/ml (55 fM BoNT/A, i.e 1.3 pg of BoNT/A) within a time scale of 5-15 minutes injection.

3.6 Time-dependent sensitivity of the assay:

The dose-dependency of the method was evaluated using prolonged cleavage time with the multispot ProteOn biosensor. This device has a slightly lower sensitivity than the Biacore apparatus but has the advantage of allowing simultaneous injection of 6 samples over 6 different substrate-coated spots. This characteristic allows a direct comparison of SPR signals produced by different BoNT/A concentrations and dramatically increases the speed of sample analysis. Concentrations of BoNT/A ranging between 0.1 and 3 LD₅₀/ml were injected 4 times each for a 15 min duration. Time course and dose-response curves were linear (Fig. S3 and Fig. 5a). 1 LD₅₀ / ml was detected in 15 min (Fig. 5a) and 0.1 LD₅₀ / ml in 1 hour (Fig. 5b). A linear dose-response relationship was also measured when BoNT/A (0.01 - 0.3 LD₅₀/ml) was injected for 5 h (Fig. 5c). In this case the LOD attained 0.01 LD₅₀/ml (Fig. 5d) and the LOQ was 0.03 LD₅₀/ml. Using a Biacore T200 the LOD was 0.03 LD₅₀/ml (1.5 fM) and 0.008 LD₅₀/ml (400 attomolar; 1 aM = 10⁻¹⁸ M) for 1 and 5 hour injections respectively (data not shown).

4. Discussion

We have developed an SPR method to detect and assay the activity of botulinum neurotoxin A, one of the most potent neurotoxins in humans. This fully automated assay consists in injecting samples containing BoNT/A over a chip coated with substrate that can then be rapidly probed to measure the enzymatic product. BoNT/A cleaves a specific peptide bond, removing a six amino acid fragment and generating a new C-terminus. Production of a

custom mAb directed against the C-terminal neo-epitope was thus crucial, as the specificity of the probe/analyte interaction determines an effective assay. Only two biological reagents, substrate and mAb10F12 are required to perform the assay and samples are directly analyzed with minimal handling.

This assay is versatile as the contact time between BoNT/A and substrate can be modified depending on the purpose of the assay and the enzyme concentration to be determined. 1 LD₅₀/ml (55 fM) was detected within a few minutes, in contrast to several hours required for the *in vivo* mouse assay. For applications in which extremely low concentrations of BoNT/A were to be detected, prolonging the cleavage time achieved a sensitivity of 0.008 LD₅₀/ml (400 aM) in 5h. We predict that longer injection periods could attain even higher sensitivity in the attomolar range. A biosensor able to process simultaneous injections was used to generate linear dose response-curves for various contact times, demonstrating the potential of the method for quantitation down to 0.03 LD₅₀/ml.

Among the fastest tests for the detection of BoNT are those based on direct fluorescence measurement associated with substrate cleavage in solution. Although these assays are very appropriate for automated high-throughput screening of BoNT/A inhibitors, they lack sensitivity compared to the SPR assay: 1 nM BoNT/A detected in 2 minutes (Schmidt and Stafford, 2003), 4 pM detected in 120 min (Poras et al., 2009) versus 55 fM in 5 minutes for our assay. Although the SPR assay includes two steps: cleavage and detection, it can still compete with single-step assays in terms of rapidity, because detection is monitored on-line in the first few seconds, by a gain-of signal resulting from BoNT/A action.

The SPR procedure attains sensitivity similar to an optimized ELISA assay which constitutes to date the most sensitive enzymatic assay reported without prior concentration of BoNT/A (Ekong et al., 1997; Jones et al., 2007). However, ELISA is a multistep, labor-intensive technique that requires long incubation periods. The automated SPR method offers several advantages over ELISA including shorter cleavage time (5 h versus 18 h) and detection time (minutes versus hours).

Several factors could contribute to the high sensitivity of the assay. Substrate immobilization on negatively charged hydrogels (dextran or alginate), that maintain SNAP-25 in an accessible conformation might accelerate enzymatic cleavage. Furthermore, as BoNT/A is positively charged at pH 7.4, it is possible that the negatively charged hydrogel concentrates the toxin and promotes its enzymatic activity in a similar way to negatively charged liposomes (Caccin et al., 2003). This kind of effect has been previously described for alpha-chymotrypsin acting on immobilized substrate (Deere et al., 2008).

No loss of mass was detected during protease injection. This is due to the fact that only a small percentage of the substrate was cleaved and furthermore the released peptide (9 amino acids) accounts for only 4% of the total mass of SNAP-25. This procedure has several advantages over biosensor assays based on the measurement of mass loss upon enzymatic cleavage. Specific monitoring of the reaction product avoids false positives resulting from the action of other BoNT serotypes or non-specific proteolysis. Moreover it does not require coupling to other methods (eg. mass spectrometry) to analyze reaction products (Henn et al., 2012). Owing to the ratio of mAb/substrate mass and the fact that only few RU can be accurately measured, the method can in theory detect enzymatic cleavage of as low as 0.005% of the total immobilized substrate. This precision cannot be achieved by detecting loss of mass. Measurement of a minute amount of cleaved substrate authorizes multiple re-use of the sensor allowing us to monitor cleavage kinetics rather than a single end-point measurement or to repeat multiple assays on the same flow cell compatible with inhibitor-screening applications.

5. Conclusions

In conclusion this new biosensor method provides a simple, rapid and highly sensitive *in vitro* assay for the enzymatic activity of BoNT/A. SPR biosensors can operate in various media (Ferracci et al., 2011) and the presence of BoNT/A could be detected in liquid foods using the method proposed. Alternatively the biosensor could be directly connected to a supply of drinking water and operate automatically to continuously monitor the presence of BoNT/A by repeatedly injecting the specific mAb at fixed times. In the future the method could be extended to detect all BoNTs that affect humans, providing antibodies that specifically recognize cleaved-substrates can be generated.

ACKNOWLEDGMENTS: This work was supported by a grant from the DGA (REI N° 0634040)

LEGENDS TO FIGURES

Fig. 1: Characterization of monoclonal antibodies designed to recognize the product of SNAP-25 cleavage by botulinum neurotoxin A.

a) A rat brain membrane fraction (Marconi et al., 2008) was incubated in the presence or absence of 10 nM BoNT/A. Western blots of membranes were stained with 5 µg/ml of six mAbs raised against the 191-197 peptide of SNAP-25. In all cases the major immunoreactive bands correspond to the 1-197 fragment of SNAP-25 generated by BoNT/A endoprotease activity. Immunoreactivity was not detected in samples untreated by BoNT/A.

b) mAb10F12 was immobilized on the active flow cell of a Biacore 3000 apparatus and an irrelevant mAb on a control flow cell. Plasma membrane-derived vesicles (Marconi et al. 2008) pre-treated with BoNT/A (10 nM) were injected in the absence (upper trace) or presence (lower trace) of an excess of SNAP-25 191-197 peptide (1 µM).

Fig. 2: Assay for the endoprotease activity of botulinum neurotoxin A: general principle.

A sample containing BoNT/A was injected 300 sec over the SNAP-25 sensorchip. The deflection to approximately 600 RU and complete return to baseline reflects the change in refractive index during sample injection, irrespective of the presence (red trace) or absence (black trace) of BoNT/A (400 fM). MAb10F12 was then injected to detect epitopes generated by BoNT/A activity. The difference between red and black traces) corresponds to specific mAb binding, and indicates the number of neo-epitopes produced, which is directly related to BoNT/A activity. Each trace represents a sensorgram monitored in the single flow-cell of a Biacore sensor chip.

Fig. 3: Botulinum neurotoxin A: assay specificity and kinetics.

a) Samples containing BoNT/A (1 pM), BoNT/B, BoNT/C, BoNT/E (all at 1 nM) or without BoNT (Cont.) were diluted in cleavage buffer and injected for 10 min over a SNAP-25 sensorchip. SNAP-25 cleavage was detected by injection of mAb10F12 (\pm SD, n= 3 mAb injections). Results are representative of 3 independent experiments. b) Seven successive injections of BoNT/A (200 fM) over immobilized SNAP-25 for a duration of 5 min. After each injection of toxin, mAb10F12 was injected and the amount of cleaved SNAP-25 measured 10 sec after the end of injection plotted in c).

Fig. 4: Rapid SPR detection of botulinum neurotoxin A at the mouse bioassay threshold.

Cleavage buffer in the presence or absence of 55 fM BoNT/A (1 LD₅₀/ml) was injected over a SNAP-25 sensorchip (Biacore T200) for 5 minutes and cleavage products immediately detected using mAb10F12. Traces represent a typical experiment with superimposed signals generated by 3 injections of mAb10F12. Experiments in the absence of BoNT/A yielded a mean of -0.7 ± 0.4 RU (n = 7 independent experiments) and a LOD of 0.5 RU (dashed line).

Fig. 5: Multispot SPR analysis of botulinum neurotoxin A activity.

(a) Six samples containing the indicated concentrations of BoNT/A were injected in parallel for 15 min and then cleaved SNAP-25 was measured from mAb10F12 binding (n=3). This procedure was then repeated 3 times to reach a cleavage time of 60 min. (b) Three superimposed sensorgrams at 0.1, 0.3 and 1 LD₅₀/ml BoNT/A at the final 60 min time point in (a). (c) as in (a) with a lower BoNT/A concentration range and a single 5 h injection. (d) as in (b) with a lower BoNT/A concentration range and a single 5 h injection. Dashed lines in (a) and (c) denote the mean background level from 3 independent experiments plus 3 SD. The figure is representative of 3 independent experiments (error bars which are often smaller than data points represent SD calculated from 3 mAb injections).

Supplemental Data Legends

Fig. S1: Detection of SNAP-25 cleaved by botulinum neurotoxin A in cultured neurons.

Cultured rat hippocampal neurons were prepared as described (Ferracci et al., 2011). Cells at 15 days in vitro were either untreated (control) or exposed to BoNT/A (1 nM) for 3 hours. Cells were then fixed for immunocytochemistry. SNAP-25 cleavage was detected using mAb10F12. Asterisk indicates soma localization. Scale bar = 20 μ m.

Fig. S2: Influence of flow rate on protease activity.

BoNT/A (100 fM) was injected over surface-bound SNAP-25 (Biacore 3000) for 5 minutes at flow rates of 1, 5 or 25 μ l/ml, followed by injection of mAb10F12. Representative of 3 independent experiments.

Fig. S3: Kinetics of botulinum neurotoxin A using a multispot biosensor.

Six samples containing the indicated concentrations of BoNT/A were injected for 15 min and then cleaved SNAP-25 was measured from mAb10F12 binding (n=3). This procedure was then repeated 3 times at each concentration, ie. a total of 4 BoNT/A injections over 1 hour. Error bars are not visible. Representative of 3 independent experiments.

References

Arnon, S.S., Schechter, R., Inglesby, T.V., Henderson, D.A., Bartlett, J.G., Ascher, M.S., Eitzen, E., Fine, A.D., Hauer, J., Layton, M., Lillibridge, S., Osterholm, M.T., O'Toole, T., Parker, G., Perl, T.M., Russell, P.K., Swerdlow, D.L., Tonat, K., 2001. Botulinum toxin as a biological weapon: medical and public health management. *Jama* 285(8), 1059-1070.

Caccin, P., Rossetto, O., Rigoni, M., Johnson, E., Schiavo, G., Montecucco, C., 2003. VAMP/synaptobrevin cleavage by tetanus and botulinum neurotoxins is strongly enhanced by acidic liposomes. *FEBS letters* 542(1-3), 132-136.

Campbell, K., McGrath, T., Sjolander, S., Hanson, T., Tidare, M., Jansson, O., Moberg, A., Mooney, M., Elliott, C., Buijs, J., 2011. Use of a novel micro-fluidic device to create arrays for multiplex analysis of large and small molecular weight compounds by surface plasmon resonance. *Biosensors & bioelectronics* 26(6), 3029-3036.

Capek, P., Dickerson, T.J., 2010. Sensing the deadliest toxin: technologies for botulinum neurotoxin detection. *Toxins* 2(1), 24-53.

Deere, J., De Oliveira, R.F., Tomaszewski, B., Millar, S., Lalaoui, A., Solares, L.F., Flitsch, S.L., Halling, P.J., 2008. Kinetics of enzyme attack on substrates covalently attached to solid surfaces: influence of spacer chain length, immobilized substrate surface concentration and surface charge. *Langmuir* 24(20), 11762-11769.

Dhaked, R.K., Singh, M.K., Singh, P., Gupta, P., 2010. Botulinum toxin: bioweapon & magic drug. *The Indian journal of medical research* 132, 489-503.

Ekong, T.A., Feavers, I.M., Sesardic, D., 1997. Recombinant SNAP-25 is an effective substrate for Clostridium botulinum type A toxin endopeptidase activity in vitro. Microbiology (Reading, England) 143 (Pt 10), 3337-3347.

Ferracci, G., Marconi, S., Mazuet, C., Jover, E., Blanchard, M.P., Seagar, M., Popoff, M., Leveque, C., 2011. A label-free biosensor assay for botulinum neurotoxin B in food and human serum. Analytical biochemistry 410(2), 281-288.

Ferracci, G., Miquelis, R., Kozaki, S., Seagar, M., Leveque, C., 2005. Synaptic vesicle chips to assay botulinum neurotoxins. The Biochemical journal 391(Pt 3), 659-666.

Hakami, R.M., Ruthel, G., Stahl, A.M., Bavari, S., 2010. Gaining ground: assays for therapeutics against botulinum neurotoxin. Trends in microbiology 18(4), 164-172.

Helmerhorst, E., Chandler, D.J., Nussio, M., Mamotte, C.D., 2012. Real-time and Label-free Bio-sensing of Molecular Interactions by Surface Plasmon Resonance: A Laboratory Medicine Perspective. The Clinical biochemist 33(4), 161-173.

Henn, C., Boettcher, S., Steinbach, A., Hartmann, R.W., 2012. Catalytic enzyme activity on a biosensor chip: combination of surface plasmon resonance and mass spectrometry. Analytical biochemistry 428(1), 28-30.

Jones, O.M., 2007. The antinociceptive effects of botulinum toxin therapy for anal fissure are unproven. Diseases of the colon and rectum 50(5), 682-683; author reply 683.

Lindstrom, M., Korkeala, H., 2006. Laboratory diagnostics of botulism. Clinical microbiology reviews 19(2), 298-314.

Loew, D., Perrault, C., Morales, M., Moog, S., Ravanat, C., Schuhler, S., Arcone, R., Pietropaolo, C., Cazenave, J.P., van Dorsselaer, A., Lanza, F., 2000. Proteolysis of the exodomain of recombinant protease-activated receptors: prediction of receptor activation or inactivation by MALDI mass spectrometry. Biochemistry 39(35), 10812-10822.

Marconi, S., Ferracci, G., Berthomieu, M., Kozaki, S., Miquelis, R., Boucraut, J., Seagar, M., Leveque, C., 2008. A protein chip membrane-capture assay for botulinum neurotoxin activity. Toxicology and applied pharmacology 233(3), 439-446.

Poras, H., Ouimet, T., Orng, S.V., Fournie-Zaluski, M.C., Popoff, M.R., Roques, B.P., 2009. Detection and quantification of botulinum neurotoxin type a by a novel rapid in vitro fluorimetric assay. *Applied and environmental microbiology* 75(13), 4382-4390.

Rossetto, O., Megighian, A., Scorzeto, M., Montecucco, C., 2013. Botulinum neurotoxins. *Toxicon* 67C, 31-36.

Schmidt, J.J., Stafford, R.G., 2003. Fluorogenic substrates for the protease activities of botulinum neurotoxins, serotypes A, B, and F. *Applied and environmental microbiology* 69(1), 297-303.

Shoji, A., Kabeya, M., Sugawara, M., 2011. Real-time monitoring of matrix metalloproteinase-9 collagenolytic activity with a surface plasmon resonance biosensor. *Analytical biochemistry* 419(1), 53-60.

Steinrucke, P., Aldinger, U., Hill, O., Hillisch, A., Basch, R., Diekmann, S., 2000. Design of helical proteins for real-time endoprotease assays. *Analytical biochemistry* 286(1), 26-34.