Correlation between *in vitro* and *in vivo* data on food digestion. What can we predict with static and dynamic digestion models?



Dr Didier DUPONT

INRA, Rennes, France



Improving health properties of food by sharing our knowledge on the digestive process

International Network

Dr. Didier DUPONT, Senior Scientist, INRA, France











Industry involvement

AGRO CAMPUS 42 private companies are following INFOGEST









Guy Vergeres

Frederic Carriere

Steven Le Feunteun

In vitro gastro-intestinal digestion Consensus INFOGEST protocol

Minekus et al. 2014 Food Funct. 5, 1113-24 Highly Cited Paper (1%)



The consensus model can be learned with videos on YouTube







The International Conference on Food Digestion



The Conference has been created by Infogest and is now an event regularly followed by 200 scientists





We are pleased to announce the next 6th International Conference on Food Digestion



in Granada, Spain, 2-4 April 2019





Static *in vitro* digestion models



Static in vitro digestion models: pro's & con's



In vitro





In vivo

Main Reasons :

Ethical - Technical - Financial

Advantages:

Standardisation of the experimental conditions Good reproducibility and repeatability Easy sampling, possibility to follow kinetics

Disadvantages:

Impossible to mimic the complexity of the GI tract in a test tube!!!

What can we predict with static digestion models?







CRITICAL REVIEWS IN FOOD SCIENCE AND NUTRITION https://doi.org/10.1080/10408398.2017.1315362

REVIEW



Check for update

Correlation between in vitro and in vivo data on food digestion. What can we predict with static in vitro digestion models?

T. Bohn^a, F. Carriere^b, L. Day^c, A. Deglaire^d, L. Egger^e, D. Freitas^f, M. Golding^g, S. LeFeunteun^f, A. Macierzanka^h, O. Menard^d, B. Mirallesⁱ, A. Moscovici^j, R. Portmann^e, I. Recioⁱ, D. Rémond^k, V. Santé-Lhoutelier^l, T. J. Wooster^m, U. Lesmes^j, A. R. Mackieⁿ, and D. Dupont^d

^aLuxembourg Institute of Health, Strassen, Luxembourg; ^bCNRS UMR 7282 EIPL, Marseille, France; ^cAgresearch, Palmerston North, New Zealand; ^dINRA UMR 1253 STLO, Rennes, France; ^eAgroscope, Institute for Food Sciences, Bern, Switzerland; ^fINRA UMR GMPA 782, Grignon, France; ^gMassey University, Palmerston North, New Zealand; ^hGdansk University of Technology, Gdansk, Poland; ⁱCIAL CSIC-UAM, Madrid, Spain; ^jTechnion—Israel Institute of Technology, Haifa, Israel; ^kINRA UNH UMR 1019, Theix, France; ^IINRA UR 370 QUAPA, Theix, France; ^mNestlé Research Centre, Nestec S.A., Lausanne, Switzerland; ⁿUniversity of Leeds, Leeds, United Kingdom

CONTACT Dr. D. Dupont 🖾 didier.dupont@inra.fr 🖃 INRA UMR 1253 STLO, 65, rue de St Brieuc, 35042 Rennes Cedex, France.



Putative secondary binding site

Interaction with PC slows β-lg gastric digestion

- * Static *in vitro* digestion models can be useful to compare series of samples or understand molecular mechanisms
- * Static *in vitro* digestion models can be relevant to estimate end-point values such as:
 - Glycaemic index
 - Some micronutrient bioavailability

* Static *in vitro* digestion models are too simple to study more complex phenomena like kinetics of digestion, food structure evolution in the GI tract, bacteria survival...





Carbohydrate digestion

AGRO



In vitro model requiring fresh human saliva and human postprandial jejunal juice (no gastric phase)

Oro-gastrointestinal model with an in vivo oral phase and the combined use of amyloglucosidase and pancreatin



Comparing the peptidome obtained during Skim Milk Protein digestion



Egger et al. 2018 Food Res Int, submitted

The INFOGEST consensus protocol is relevant for protein digestion

What about for the other nutrients?

Future experiments in WG4 and WG5





Bioaccessibility/bioavailability of carotenoids



Correlation of bioaccessibility and bioavailability of various carotenoids (lycopene, beta-carotene, lutein), as presented by <u>Reboul, et al. 2006</u> (black points) and <u>Alminger, et al. 2012</u> (grey points)





Bioaccessibility/bioavailability of iron



In conclusion, for micronutrients:

Direction of bioavailability (higher, lower, equally well) appears to be generally predictable by in vitro digestion

But does mostly not predict the magnitude of bioavailability in humans

Transport, metabolism, colonic changes (especially for polyphenols), biodistribution and (renal) excretion affect bioavailability

Correlation of iron bioaccessibility (dialysis method) following *in vitro* digestion from various iron containing test meals and bioavailability (as estimated via erythrocyte incorporation) employing the same formulations. Empty black dots: Data from Walter et al. (2003). Black dots: Data from Aragon et al. (2012). Grey dots: Data from (Davidsson et al. 2002). Empty grey dots: (Sandberg 2005).









Dynamic in vitro digestion models



Dynamic Gastric Model (DGM)





Human Gastric Simulator (HGS)





Artificial Colon (ARCOL)



Digestion Dynamique Gastro-Intestinale (DIDGI)



Multi-compartmental models

Engineered Stomach and small INtestinal (ESIN)











TNO Gastro-Intestinal ModelTIM









Simulator of the Human Intestinal Microbial Ecosystem (SHIME)

Can dynamic *in vitro* digestion systems mimic the physiological reality?

Dupont D.^{a*}, Alric M.^b, Blanquet S.^b, Bornhorst G.^c, Cueva C.^d, Deglaire A.^a, Denis S.^b, Ferrua M.^e, Havenaar R.^f, Lelieveld J.^f, Mackie A.R.^g, Marzorati M.^h, Menard O.^a, Minekus M.^f, Miralles B.^d, Recio I.^d, Thuenemann E.ⁱ, Van den Abbeele P.^j







GHENT

UNIVERSITY



CONSEJO SUPERIOR DE INVESTIGACIONES CIENTÍFICAS

Submitted soon to Crit Rev Food Sci Nutr



UNIVERSITY OF LEEDS





For reprints: didier.dupont@inra.fr





.019

DGM - Simulation of grinding in the stomach

Comparison of the Mean Breaking Time of agar gel beads in low viscosity (LV) and high viscosity (HV) meals







HGS – Comparison of white rice gastric emptying and intragastric pH with growing pigs data







.021



Cumulative gastric and ileal delivery of a meal expressed as a percentage of total intake: *in vivo* (human n=7) gastric (\bullet) and ileal (O) delivery of yoghurt and gastric (\blacksquare) and ileal (\Box) delivery of blue dextran in the model simulating the slow transit of yoghurt

TIMcarbo: prediction of human glycemic response curve by combining CHO digestion in tiny-TIM, rapid glucose/fructose analysis and in silico modelling.







Glycemic response: In vivo (dotted) vs. predicted (solid) Bel

Bellmann et al., submitted for publication



In vivo \Box Ascending colon \diamond Transverse colon \triangle Descending colon

SHIME[®] – *in vitro / in vivo production* of Short-Chain Fatty Acid during fermentation of various polysaccharide sources

DIDGI®- Validation of protein digestion in infant formula



.025

Conclusion

Static in vitro digestion models are interesting for comparing large series of samples in the same conditions or understanding interaction between constituants

They could also be valuable for predicting end-point values (glycaemic index, micronutrient bioaccessibility...) but are not relevant for investigating dynamic events

Tynamic *in vitro* digestion models can perfectly mimic the physiological reality when parameters are available in the literature

But when a system has been validated for a specific food, does this means that it is relevant for all kind of foods? Can we define large categories of foods (liquids, gels, solids) and validate them?

- More validated *in vitro* dynamic digestion models are needed to simulate the digestion process of specific populations (infant, elderly...)
- The international consensus could be found on the parametering of *in vitro* dynamic systems

Next objectives of the Infogest network





The next INFOGEST Workshop will be held on the 12-13 of April 2018 in Leeds UK connected to:

17th Food Colloids Conference : Application of Soft Matter Concepts

8th - 11th April 2018 The University of Leeds, UK



ABOUT VENUE ABSTRACTS REGISTRATION ACCOMMODATION SPEAKERS ORGANISERS SPONSORS CONTACT

Food Colloids is a biennial conference in the field of physical chemistry of complex foods.

SUBMIT ABSTRACT & REGISTER

