

# Reliable Determination of Ge in Solid Environmental Samples Using a Chemical Preparation Procedure Developed for Si Isotopes and ICP-MS Analysis

Camille Delvigne, Abel Guihou, Bernard Angeletti, Isabelle Basile-Doelsch, Jean-Dominique Meunier

### ▶ To cite this version:

Camille Delvigne, Abel Guihou, Bernard Angeletti, Isabelle Basile-Doelsch, Jean-Dominique Meunier. Reliable Determination of Ge in Solid Environmental Samples Using a Chemical Preparation Procedure Developed for Si Isotopes and ICP-MS Analysis. Geostandards and Geoanalytical Research, 2018, 42 (1), pp.139 - 149. 10.1111/ggr.12197. hal-01695926

HAL Id: hal-01695926

https://hal.science/hal-01695926

Submitted on 18 May 2018

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.






# Reliable determination of Ge using the silicon isotopes chemical preparation and ICP-MS analysis for solid environmental samples.

Journal:	Geostandards and Geoanalytical Research
Manuscript ID	GGR-0514.R1
Manuscript Type:	Original Article
Date Submitted by the Author:	n/a
Complete List of Authors:	Delvigne, Camille; Centre Europeen de Recherche et d'Enseignement des Geosciences de l'Environnement, Guihou, Abel; Centre Europeen de Recherche et d'Enseignement des Geosciences de l'Environnement Angeletti, Bernard; Centre Europeen de Recherche et d'Enseignement des Geosciences de l'Environnement Basile-Doelsch, Isabelle; Centre Europeen de Recherche et d'Enseignement des Geosciences de l'Environnement Meunier, Jean-Dominique; Centre Europeen de Recherche et d'Enseignement des Geosciences de l'Environnement
Keywords:	germanium, ICP-MS, reference materials, soils, plants, sample digestion, cation exchange purification

SCHOLARONE™ Manuscripts



Germanium (Ge) exists at trace levels in the Earth's crust and is a powerful geochemical
tracer of the silicon (Si) cycle. This study proposes a simple and reliable method for
determining Ge contents in environmental samples using ICP-MS. As Si and Ge have very
similar chemical properties, we investigated the applicability of the chemical preparation
procedure developed for Si isotopes to the analysis of Ge in environmental samples.
Advantages of this procedure are: (1) efficient removal of the matrix and main interferences
affecting Ge analyses by ICP-MS; (2) a low limit of detection (6 ng l <sup>-1</sup> ); (3) relative
repeatability of approximately 3% obtained on <sup>74</sup> Ge; and (4) robustness and accuracy based
on agreement within errors with the published Ge concentrations for rock standards (BHVO-
2, AGV-2 and BCR-2). This procedure allows revising the Ge concentrations for 3 soils
standards (1.67 $\pm$ 0.09 $\mu g$ g <sup>-1</sup> ; 2.41 $\pm$ 0.18 $\mu g$ g <sup>-1</sup> ; 1.89 $\pm$ 0.10 $\mu g$ g <sup>-1</sup> for GBW 07401, GBW 07404
and GBW 07407 respectively) and proposing a value for the ERM-CD281 plant standard
$(0.06\pm0.01~\mu g~g^{-1})$ . This method provides a convenient procedure for determining Ge
concentrations in environmental samples and opens the possibility of coupling two tracers of
the Si biogeochemical cycle with a single analytical procedure.

Résumé

Le germanium (Ge), présent à des teneurs en trace dans la croute terrestre, est un traceur
géochimique performant du cycle du silicium (Si). Cette étude propose une méthode simple,
rapide et fiable pour déterminer les concentrations en Ge dans les échantillons
environnementaux par ICP-MS. Vu les propriétés géochimiques très similaires entre Ge et Si,
nous avons investigué la pertinence de la procédure de mise en solution et de purification
développée pour les isotopes du silicium à l'analyse du Ge dans des échantillons
environnementaux. Les avantages de cette technique sont : (1) une élimination efficace de la
matrice et des interférences principales affectant l'analyse du Ge par ICP-MS ; (2) une limite
de détection de 6 ng l <sup>-1</sup> ; (3) une répétabilité relative de 3% obtenue sur l'isotope <sup>74</sup> Ge; (4) une
bonne fiabilité et justesse basées sur la concordance entre les valeurs de Ge publiées pour les
standards de roche (BHVO-2, AGV-2 et BCR-2) aux incertitudes près. Cette méthode a
permis de réévaluer les concentrations en Ge de 3 standards de sol (1.67±0.09 μg g <sup>-1</sup> ;
2.41±0.18 μg g <sup>-1</sup> ; 1.89±0.10 μg g <sup>-1</sup> pour GBW 07401, GBW 07404 et GBW 07407
respectivement) et de proposer une valeur pour le standard de plante ERM-CD281 (0.06±0.01
$\mu g \; g^{\text{-1}}$ ). Cette méthode présente aussi le grand avantage de coupler deux traceurs du cycle du
silicium (à savoir les isotopes stables du silicium et le rapport Ge/Si), en une préparation
chimique unique.

1. Introduction

55	Germanium (Ge) exists at trace levels in the Earth's crust and is a powerful geochemical
56	tracer of the silicon (Si) cycle. Like Si, Ge is a Group IV element. Both elements have
57	identical outer electronic configurations and similar covalent radii (Ge: 1.22Å and Si 1.17Å;
58	Höll et al., 2007). The geochemistry of Ge is dominated by its propensity to substitute for Si
59	in silicate minerals, and thus closely follows Si through its biogeochemical cycle.
60	Fundamentals of high and low-temperature geochemistry of Ge are reviewed in Rouxel and
61	Luais (2017). To summarize, the Ge/Si ratio has been employed as a tracer of weathering
62	processes (e.g., Kurtz and Derry, 2004; Scribner et al., 2006; Opfergelt et al., 2010; Lugolobi
63	et al., 2010), hydrothermal and geothermal activity (e.g., Mortlock et al., 1993; Evans and
64	Derry, 2002), biological processes (e.g., Ellwood and Maher, 2003; Derry et al., 2005; Suttor
65	et al., 2010), and flow paths (e.g., Kurtz et al., 2011) in modern environments but also in
66	paleo-environmental studies (e.g., Filippelli et al., 2000; Hammond et al., 2004). The use of
67	Ge/Si has recently been extended to the study of the Archean Si cycle (e.g., Hamade et al.,
68	2003; Frei and Polat, 2007; Delvigne et al., 2012; 2016). Another facet of Ge environmental
69	research focuses on the potential economic resource of Ge in soils and plants as the Ge
70	demand for new technologies is exploding while Ge is rare in natural environments
71	(Rosenberg, 2009; Wiche et al., 2017). However, Ge content is often overlooked in the array
72	of analyzed elements. This arises from the analytical difficulties to obtain reliable data
73	because Ge concentrations in environmental samples are generally below the limit of
74	detection of many modern analytical methods. In addition, Ge suffers from various
75	interferences on its isotopes with ICP-MS techniques (mainly <sup>56</sup> Fe <sup>16</sup> O on <sup>72</sup> Ge, <sup>57</sup> Fe <sup>16</sup> O and
76	<sup>56</sup> FeH <sup>16</sup> O on <sup>73</sup> Ge, and <sup>74</sup> Se and <sup>58</sup> Ni <sup>16</sup> O on <sup>74</sup> Ge). Interferences due to iron oxides are
77	particularly problematic because some samples contain a few tenths of a percent of Fe while
78	Ge is present at a μg g <sup>-1</sup> level. The <sup>58</sup> Ni <sup>16</sup> O interference is also an issue as it combines the

most abundant Ni and O isotopes (58Ni=68%; 16O=99%). In addition, Ar-based molecular interferences (<sup>40</sup>Ar<sup>16</sup>O<sub>2</sub>, <sup>36</sup>Ar<sup>36</sup>Ar and <sup>38</sup>Ar<sup>36</sup>Ar on <sup>72</sup>Ge, <sup>73</sup>Ge and <sup>74</sup>Ge, respectively) should also be monitored. Germanium is then often qualified as a "tricky" element to analyze and is not routinely determined in environmental samples. The majority of studies published so far on Ge have been carried out using the isotope-dilution hydride-generation ICP-MS technique (ID-HG-ICP-MS) developed by Mortlock and Froelich (1996) (Fillipelli et al., 2000; Kurtz et al., 2002; Evans and Derry, 2002; Derry et al., 2005; Wheat and McManus, 2005; Ellwood et al., 2006; Scribner et al., 2006; Blecker et al., 2007; Makishima and Nakamura, 2009; Lugolobi et al., 2010; Sutton et al., 2010). This is probably the most reliable and precise technique for Ge determination, although it requires a hydride generation system, which is not widespread in laboratories. Alternative techniques are: (1) graphite furnace atomic absorption spectrometry (GF-AAS) (McMahon et al., 2006); (2) ICP-MS analysis with or without high resolution (Hamade et al., 2003; Delvigne et al., 2009; Cornelis et al., 2010; Tribovillard et al., 2011); (3) laser ablation ICP-MS analysis (Hamade et al., 2003; Shen et al., 2011; Belissont et al., 2014; Dong et al., 2015;); and (4) X-ray fluorescence spectrometer analysis (Frei and Polat, 2007). While valid for most geological samples, these techniques are working close to their limits of detection for environmental samples, such as plants that commonly range from 0.07 to 20 ng g<sup>-1</sup> (Derry et al., 2005; Blecker et al., 2007; Delvigne et al., 2009; Cornelis et al., 2010; Lugolobi et al., 2010). One option is to preconcentrate samples in order to reach concentrations above the detection limits. This was investigated by Soylak and Yigit (2015), who suggested a separation-preconcentration procedure using a polysulfone membrane filter combined with spectrophotometric measurements. This method has the disadvantage of requiring specific materials. For natural water samples, ID-HG-ICP-MS remains the only reliable technique thanks to its high precision and unequaled low detection limits.

The objective of this study is to propose a simple and reliable alternative method determining Ge contents in solid environmental samples, with a notable focus on soils and plants using ICP-MS, a common instrument found in geochemistry laboratories. As silicon isotopes and Ge/Si ratios are complementary tracers of the Si biogeochemical cycle, we investigate the applicability of the chemical preparation procedure developed for silicon isotopes (Georg et al., 2006) to the analysis of Ge in solid environmental samples.

2. Experimental

2.1. Material analyzed

Reference materials with the recommended values for Ge concentrations are scarce. Three international geological USGS reference materials have been chosen in this study for intercomparison with previous studies: BHVO-2 (basalt), AGV-2 (andesite), and BCR-2 (basalt). Despite certified values not being available for these standards, an increasing number of published data provide constraints on Ge concentrations (Table 1). A considerable compilation effort was undertaken by Jochum et al. (2016), which provides reference values (or informative values) determined following ISO guidelines and IAG Certification Protocol using data published between 1995 and 2015. In addition, to better fit the matrix of targeted samples of soils and plants, we selected three soil standards from the Institute of Geophysical and Geochemical Exploration (IGGE, Langfang, China) (GBW 07401, GBW 07404 and GBW 07407, also known as GSS-1, GSS-4 and GSS-7) and a plant standard from European Reference Material (ERM-CD281). Certified Ge concentrations are available for the soils standards, but with low precision (Table 1), while no Ge data is provided for ERM-CD281. Recently, additional data have been made available for GBW 07401-07408 standards (Liu et al., 2014) (Table 1).

	Sample type	Provider	Published value	rsd (2σ)	n	References
BHVO-2	basalt	USGS	1.62±0.04 μg g <sup>-1</sup>	2%	10	Jochum et al., 2016
			1.53 $\mu g g^{-1}$	0.2%	n.s.	Escoube et al., 2012
			$1.59\pm0.04~\mu g~g^{-1}$	3%	3	Scribner et al., 2006
AGV-2	andesite	USGS	$1.20\pm0.08~\mu g~g^{-1}$	7%	6	Jochum et al., 2016
BCR-2	basalt	USGS	$1.46\pm0.26~\mu g~g^{-1}$	18%	6	Jochum et al., 2016
GBW 07401	soil	IGGE	$1.50\pm0.04~\mu g~g^{-1}$	3%	6	Liu et al., 2014
			1.34±0.2 μg g <sup>-1</sup>	15% *	at least 8	Wang et al., 2013
GBW 07404	soil	IGGE	$2.10\pm0.02~\mu g~g^{-1}$	1%	6	Liu et al., 2014
			1.9±0.3 μg g <sup>-1</sup>	16% *	at least 8	Wang et al., 2013
GBW 07407	soil	IGGE	1.6±0.3 μg g <sup>-1</sup>	19% *	at least 8	Wang et al., 2013
ERM-CD281	plant	ERM	no data	-	-	-

**Table 1.** Reference materials analyzed in this study with their Ge concentrations and associated uncertainties. Values in bold are certified values while others are informative values. \* refers to an uncertainty of  $3\sigma$  instead of  $2\sigma$ ; n.s.= not specified.

2.2. Sample digestion

Samples were digested following the procedure developed for silicon isotopes using a solid sodium hydroxide (NaOH) flux (Georg et al., 2006). Between 10-20 mg (200 mg up to 600 mg for ERM-CD281, the plant standard) of powdered samples were weighed directly into silver crucibles after drying for at least 24h at 105°C. Prior to the fusion, calcination steps were carried out to ensure the combustion of organic matter into CO<sub>2</sub>. Calcinations were performed at 450°C with a 1 h stage at 250°C to ensure slow temperature increase. The 450°C calcination step lasted for at least 1 hour for soils, and 16 h for plants. In the case of plants, the total calcination procedure lasted at least 24 h. Qi et al., (2011) showed that calcination at 600°C does not induce Ge loss. Ashes were mixed with approximately 200 mg of NaOH and placed in a furnace at 720°C for 10 minutes. Fusion cakes were then allowed to dissolve overnight in ultrapure water with a brief ultrasonic bath to favor dissolution. Solutions were then quantitatively transferred into pre-cleaned HDPE bottles and acidified to pH ~2 with HNO<sub>3</sub>. Due to the volatile behavior of Ge in the presence of halogens (HCl, HClO<sub>4</sub>), only

HNO<sub>3</sub> could be used (Luais, 2007; 2012). Acidification to pH~2 was important for Fe-rich samples to ensure a complete dissolution of Fe-oxides as well as ensuring a complete recovery during the following cationic purification (Fitoussi et al., 2009). At this pH, the dominant Ge species is the Ge(OH)<sub>4</sub><sup>0</sup> form (Pokrovski and Schott, 1998).

#### 2.3. Sample purification

To overcome matrix effects and various interferences, which are complex to correct for, purification of samples appears to be an ideal approach. Perhaps more importantly, the elimination of matrix elements prior to analysis allows a substantial reduction in the dilution factor, which is of great interest for samples with low Ge contents. In addition, performance of ICP-MS can be fully exploited, as the sampling interface and the ion optics are not deteriorated due to salt depositions. A cation-free solution was obtained by passing the solution through a cation exchange column. This allowed for testing of the resin used for germanium isotopes chemical preparation, AG50W-X8 (100-200 mesh in H<sup>+</sup> form; Bio-Rad) (Luais, 2007, 2012) using the column chemistry routine of silicon isotopes (Georg et al., 2006). Sample solutions were loaded on a Bio-Rad type column filled with 1.8 ml of AG50W-X8 cation-exchange resin, pre-cleaned with several rinsing with HCl, HNO<sub>3</sub> as detailed in the study of Georg et al., (2006) (Table 2). Before loading the sample on the resin, ultrapure water was passed into the resin to remove acids (Georg et al., 2006). The matrix (sodium from fusion and other cations from the sample itself) was retained on the resin, while anionic and neutral species (e.g., Ge(OH)<sub>4</sub><sup>0</sup> and Si(OH)<sub>4</sub><sup>0</sup>) passed through. To minimize the dilution due to elution, the sample volume was maximized (here ~ 10 ml of fused solution) while the volume of eluent was kept to a minimum of one resin bed volume (~ 2 ml) (Table 2). Since the Si and Ge species do not interact with the resin, the eluent is simply ultrapure water. Ge was collected in a  $\sim$  12 ml cation-free solution. Recovered solutions were weighed to minimize uncertainties inherent to the imprecise volumes recovered. Regarding the

removal of matrix element, Na is potentially the more difficult element to eliminate as it is the most abundant and resin shows a weak relative selectivity for Na. To warrant a good purification, resin is largely in excess based on expected resin capacity and fused sample concentrations. Indeed, a 10 ml fused solution (typically 0.5 meq) saturates the resin approximately 15% of its exchange capacity. Complete removal of ambient cationic species and especially problematic cations (*i.e.*, Fe, Ni, Zn) and Na was checked with ICP-MS.

BioRad AG50W-X	8, 1.8 ml resin bed	
Separation stage	Solution matrix	Volume (ml)
Pre-cleaning	3N HCI	3 ml
Pre-cleaning	6N HCI	3 ml
Pre-cleaning	7N HNO <sub>3</sub>	3 ml
Pre-cleaning	10N HCI	3 ml
Pre-cleaning	6N HCI	3 ml
Pre-cleaning	3N HCI	3 ml
Rinse	ultrapure water	6 ml
Sample load	Acidified fused sample	10 ml
Elution	ultrapure water	2 ml

**Table 2.** Column chemistry routine of the preparation of a solid sample fused with NaOH flux adapted from Georg et al., 2006.

#### 2.4. Sample introduction and mass spectrometry

Ge concentrations were measured using an ICP-MS Nexion 300X (Perkin Elmer) using a PFA nebulizer, a quartz cyclonic spray chamber and a SC-FAST DX-4 automated introduction system. Operating conditions are detailed in Table 3. A series of five standard solutions with different concentrations (50, 100, 250, 500, 1000 ng l<sup>-1</sup>) were used to calibrate Ge concentrations. These calibration solutions were made from single element Ge standard 1000 mg l<sup>-1</sup> solution (Chem Lab, Plasma HIQU, 2-5% HNO<sub>3</sub>). To correct for instrumental drift, rhodium (Chem Lab, Plasma HIQU, 2-5% HNO<sub>3</sub>) was introduced to sample solutions through a cross flow line injection resulting in a Rh concentration of 1.58 μg l<sup>-1</sup> in the analyzed solution. Ge isotopes 72, 73 and 74 were analyzed and provided consistent values, proving the

robustness of the method. However, Ge contents were determined using the <sup>74</sup>Ge isotope, the most abundant isotopes (36.28%), given higher sensitivity and lower detection limits compared to other isotopes. In the sequence, 400 ng l<sup>-1</sup> Ge standard solutions were measured after every 15 samples to check for instrumental drift and evaluate reproducibility.


Plasma power	1600 W
Torch	Quartz glass torch
Plasma Ar gas flow rate	15 l min <sup>-1</sup>
Auxiliary Ar gas flow rate	1.2 l min <sup>-1</sup>
Nebuliser Ar gas flow rate	0.95 - 1.05 l min <sup>-1</sup>
	(optimized to minimize $^{140}\text{Ce}^{16}\text{O}^+/^{140}\text{Ce}^+<0.025$ and $^{70}\text{Ce}^{++}/^{140}\text{Ce}^+<0.03$ )
Nebuliser	PFA 0.4 ml min <sup>-1</sup>
Spray chamber	cyclonic (quartz)
Sample uptake rate	0.337 ml min <sup>-1</sup>
Interfaces cones	Nickel
Resolution	700
Typical <sup>74</sup> Ge sensitivity	~ 12000 cps μg <sup>-1</sup> l
Typical <sup>103</sup> Rh sensitivity	~ 680000 cps μg <sup>-1</sup> l

**Table 3.** Instrumental operating conditions for ICP-MS Nexion 300X

- 3. Results and discussion
- 3.1. Assessment of Ge recovery yields during chemical procedure

As the dominant Ge species after fusion is the Ge(OH)<sub>4</sub><sup>0</sup>, no affinity for cation exchange resin is expected and thus Ge should be quantitatively recovered through the chromatographic procedure. A complete Ge recovery (102±2%; n=6) during sample purification through cationic resin was obtained with single element Ge standard solution at different Ge concentrations (330 to 1100 ng l<sup>-1</sup>; Fig. 1). To take into account matrix effects, recovery yields for BHVO-2 and GBW 07401 were determined by Ge standard addition to the fused solution before purification (Ge addition from 10 ng l<sup>-1</sup> to 100 ng l<sup>-1</sup>; Fig. 1). The recovery yields for BHVO-2 (99±2%; n=4) indicated a complete recovery for Ge. In the case of GBW 07401 (Ge addition from 10 ng l<sup>-1</sup> to 1000 ng l<sup>-1</sup>; Fig. 1), recoveries of approximately 118±5% (n=6) were found when compared to the certified value (Wang et al., 2013). Considering the

Ge concentrations of Liu et al., (2014), who suggested that the Ge concentration of this standard should be revised, our recovery yield was 109±2% (Fig. 1). As the procedural blanks displayed very low levels of Ge (3±2 ng l<sup>-1</sup>; n=11) and the given large error range associated with the certified value of this standard (Table 1), we assume our recovery for this standard to be complete. A better agreement would require refining the reference value for the GBW 07401 standard (see section 3.2). Based on single element Ge standard solution and BHVO-2, we assume that Ge recovery is complete during cationic purification. The standard addition of BHVO-2 and GBW 07401 show that Si and other neutral, as well as anionic species do not create any matrix effect.



**Fig. 1** Comparison of the Ge concentrations obtained by ICP-MS after cationic purification and expected values for single element Ge standard solution at various Ge concentrations (gray circles) and reference materials (BHVO-2, black squares; GBW 07401, open diamonds) with Ge standard addition (Ge addition from 10 ng l<sup>-1</sup> to 1000 ng l<sup>-1</sup>). The dashed gray line represents a recovery yield of 100%. Expected values for reference materials are calculated

considering Ge data from the literature (Escoube et al., 2012; Scribner et al., 2006 and Liu et al., 2014, see section 3.2). Error bars represent the analytical repeatability (3%).

The separation efficiency is demonstrated by the elution curves of two samples (a single element Ge standard solution and a fused BHVO-2) that were loaded onto the pre-cleaned cation exchange resins (Fig. 2). Ge is not retained by the resin and migrates straight through the column and nearly 90% of the loaded Ge is already recovered while the sample load is still infiltrating the resin bed. The last 10% is recovered with 2 ml of ultrapure water (Fig. 2).





Fig. 2 Elution of germanium during chemical purification of BHVO-2 on AG50W-X8 resin.

The black circles corresponds to the cumulated Ge recovered whereas grey bars represent Ge recovered in each volume loaded. Error bars represent the analytical repeatability.


Furthermore, ICP-MS screening of the purified solution confirms the absence of major cations and interfering cations after purification (Fig. 3). As demonstrated for Si by Georg et al., (2006), this shows a complete separation of Ge from the ambient cation matrix. Matrix effects and interferences are then negligible and should not be corrected for, minimizing the errors in

Ge contents. In addition, Ar-based molecular interferences are not detectable in the blanks.



**Fig. 3**. Recovery yields of the main matrix elements and interfering elements of NaOH fused BHVO-2 after purification with AG50W-X8 resin.

Besides overcoming matrix effects and interferences, cationic purification has the advantage to minimize the dilution for analysis, as samples are salt-free. Indeed, the above described sample purification resulted in a dilution of the sample by a factor of 1.2. This is far less than the dilution necessary to reduce total dissolved salt to the maximum of 0.1% accepted for ICP-MS analysis (a dilution by a factor 2 being the minimum). However, some samples with low Ge content might still be close to the detection limit. To overcome this issue, a preconcentration step is helpful. A preconcentration step through low temperature evaporation was feasible, as solutions were free of major elements, limiting the risk of precipitation. For samples that might be close to detection limits, typically plants, an evaporation step at 60°C on a hotplate was performed to increase Ge concentrations. In this study, a maximal preconcentration factor of 2 was achieved without any modification of Ge recoveries. Indeed, aliquots of purified solutions of geostandards BHVO-2, AGV-2 and BCR-2, with or without evaporation, displayed similar Ge contents (Fig. 4).



**Fig. 4** Ge contents in reference materials (black; Escoube et al., 2012; Scribner et al., 2006; Jochum et al., 2016) compared to ICP-MS Ge contents obtained using our preparation protocol, with (white) or without a final evaporation step (gray). Error bars on our results represent an uncertainty of 3%, the ICP-MS analytical repeatability.

## 3.2. Detection limits, repeatability and accuracy

The contribution of Ge from the whole procedural chemistry blank was measured on 11 procedural blanks and has been found to be  $3\pm2$  ng  $\Gamma^{-1}$  and indistinguishable from the instrumental blank. According to the IUPAC definition (IUPAC, 1997), method detection limits were calculated as the mean concentration of 11 procedural blanks plus three standard deviations of these 11 procedural blanks (Table 4). The lowest detection limit was obtained for  $^{74}$ Ge (6 ng  $\Gamma^{-1}$ ), the most abundant Ge isotope. For comparison, the detection limit for  $^{72}$ Ge is 5 times higher, approximately 30 ng  $\Gamma^{-1}$  most likely due to the ArO<sub>2</sub> interference. Considering that 20 mg of powdered sample was dissolved in 100 ml, this corresponds to a detection limit of approximately 0.03  $\mu$ g g<sup>-1</sup> in the solid sample (Table 4). For the plant standard (200-600 mg dissolved in 100 ml), the detection limit is approximately 0.001  $\mu$ g g<sup>-1</sup> in the solid sample. Additionally, the linearity limit was demonstrated over the entire

calibration range. The linearity limit was above 1000 ng l<sup>-1</sup>, the highest standard in the calibration curve, as the coefficient correlation of the calibration curve is above 0.995.

Method	Dissolution	Preliminary steps	DL in solution	DL in solid samples	Reproducibility	Difference with previously published values	Samples	References
ID-HG-ICP-MS	acid	ID	0.5 ng l <sup>-1 a</sup>	-	<3%	n.s.	natural waters	Evans and Derry, 2002
HR-ICP-MS	fusion	-	10 ng l <sup>-1</sup>	n.s.	4%	n.s.	plants and soils	Delvigne et al., 2009
ICP-MS	fusion	-	2 ng l <sup>-1</sup>	0.05 μg g <sup>-1</sup>	15%	n.s.	plants and soils	Cornelis et al., 2010
LA-ICP-MS	-	-	-	< 1 μg g <sup>-1</sup>	8% <sup>b</sup>	n.s.	rocks	Belissont et al., 2011; Shen et al., 2011
X-ray fluorescence	-	-	-	1.5 μg g <sup>-1 c</sup>	1%	15%	rocks	Frei and Polat, 2007
Spectrophotometry	acid	Polysulfone membrane filter	2 ng l <sup>-1</sup>	n.s.	5%	n.s.	waters and soils	Soylak and Yigit, 2015
ICP-MS	fusion	Cationic purification	6 ng l <sup>-1</sup>	0.03 μg g <sup>-1</sup>	3%	6%	plants, soils, rocks	This study

a blank level

**Table 4.** Comparison of metrological characteristics of different techniques used for Ge determination. n.s. = not specified. DL = detection limit.

The relative repeatability of ICP-MS analysis based on a Ge single element standard solution (400 ng l<sup>-1</sup>) and analyzed six times on a single day is 3% (n=6, confidence level 95%). Similarly, four aliquots of one fused solution (BHVO-2) passed through four different columns and analyzed twice on the same day gives a relative repeatability of 2%.

Reference materials of different composition (BHVO-2; AGV-2; BCR-2; GBW 07401; GBW 07404; GBW 07407) were analyzed to evaluate the accuracy of results (Supplementary Table 1). However, the absence of properly certified Ge values complicated the exercise. As shown in Fig. 5, our BHVO-2 data agree within 10% with the published value from the compilation of Jochum et al. (2016) (Table 5). However, a careful sorting of compiled data reveals that most of the dissolution protocols are not adapted to Ge, except for that of Scribner et al. (2006). Indeed, all but one of the studies employed acid attacks using HCl, HClO<sub>4</sub> or HF (at temperature above 60°C) or a combination of these, which induce loss of Ge (Luais et al., 2012). Tests performed by Luais (2012) demonstrated that Ge loss is at 85% when HClO<sub>4</sub> is employed in the acid dissolution mixture. In addition, it is unclear if ICP-MS analysis were corrected for spectroscopic interferences and selected Ge isotopes were not specified. Thus,

b valid for Ge/Si ratio as detailed in Shen et al., 2011

<sup>&</sup>lt;sup>c</sup> minimum Ge concentration analyzed (no limit of detection provided)

Ge might be underestimated due to potential Ge volatilization during dissolution, while uncorrected interferences would overestimate Ge content. In the literature, Ge contents estimated using specific techniques devoted to Ge analysis (ID-HG-ICP-MS or double-spike MC-ICP-MS) appear to be systematically lower than Ge contents estimated with ICP-MS techniques. In the case of BHVO-2, ID-HG-ICP-MS and double spike MC-ICP-MS techniques provided a compiled value of 1.56±0.07 µg g<sup>-1</sup> (Table 5; Scribner et al., 2006; Escoube et al., 2012), while ICP-MS techniques (acid attacks and interferences corrections not specified) averaged approximately 1.62±0.04 µg g<sup>-1</sup> (Table 5; Jochum et al., 2016). Similarly, BHVO-1 data show a comparable shift between the two techniques (1.46±0.12 μg g<sup>-1</sup>, n=2 for ID-HG-ICP-MS and double spike MC-ICP-MS; 1.64±0.06 µg g<sup>-1</sup>, n=3 for ICP-MS; Jochum et al., 2016). An intermediate Ge content (1.58 µg g<sup>-1</sup>) is obtained by Luais (2012) using dedicated acid attack removing 80% of matrix analyzed by ICP-MS suggesting that discrepancy between techniques is likely induced by matrix and interferences. Despite being more limited in number, data acquired using isotope dilution techniques are the most robust (Escoube et al., 2012; Scribner et al., 2006). Considering the compiled values from these two studies, our BHVO-2 results agree within 6% (Fig. 5, Table 5).




Fig. 5 Ge content in BHVO-2 from our study (open diamonds), high quality data (gray symbols) from Escoube et al. (2012) (gray circle; double spike MC-ICP-MS), and Scribner et al. 2006 (gray triangle; ID-HG-ICP-MS), compiled data from Jochum et al. 2016 (filled circles; ICP-MS and one ID-HG-ICP-MS from Scribner et al., 2006). Bold and dashed lines represent the mean and uncertainties (95% confidence level) of each subgroup, respectively. Errors bars represent the analytical precision when provided (*i.e.*, the 3% repeatability in this study).

Regarding the other rock standards (AGV-2 and BCR-2), there is a better agreement between our data and the compiled data from Jochum et al., (2016) (Table 5). Considering these three rock standards, our data agree within ~5% with published values.

Regarding soil standards (GBW 07401, 07404, 07407), our results are systematically 20 to 30% higher than certified values provided in Wang et al. (2013). However, a better agreement is noted when compared to data from Liu et al. (2014), suggesting that certified values should

be revised. The use of HClO<sub>4</sub> in the sample decomposition procedure calls into question the accuracy of these data and may explain our systematically higher Ge contents.

	Sample type	Ge value (this study)	n	Difference with published values	Published values	n	Technique	References
				9%	1.62±0.04 μg g <sup>-1</sup>	10	ICP-MS but one	Jochum et al., 2016
BHVO-2	basalt	1.47±0.04 μg g <sup>-1</sup>	11	8%	1.56±0.04 μg g <sup>-1 a</sup>	3	ID-HG-ICP-MS and double spike MC-ICP-MS	Scribner et al., 2006 and Escoube et al., 2012
AGV-2	andesite	1.14±0.08 μg g <sup>-1</sup>	6	5%	1.20±0.08 μg g <sup>-1</sup>	6	ICP-MS	Jochum et al., 2016
BCR-2	basalt	1.43±0.14 μg g <sup>-1</sup>	5	2%	1.46±0.26 μg g <sup>-1</sup>	6	ICP-MS but one	Jochum et al., 2016
CDW 07401	i <b>BW 07401</b> soil	1.67±0.09 μg g <sup>-1</sup>	11	11%	1.50±0.04 μg g <sup>-1</sup>	6	ICP-MS	Liu et al., 2014
GBW 07401				25%	1.34±0.2 μg g <sup>-1</sup>	at least 8	n.s.	Wang et al., 2013
GBW 07404	soil	2 44 (2 42 -1	4	15%	2.10±0.02 μg g <sup>-1</sup>	6	ICP-MS	Liu et al., 2014
GBW 07404	SOII	2.41±0.18 μg g <sup>-1</sup>	4	27%	1.9±0.3 μg g <sup>-1</sup>	at least 8	n.s.	Wang et al., 2013
GBW 07407	soil	1.89±0.10 μg g <sup>-1</sup>	3	18%	1.6±0.3 μg g <sup>-1</sup>	at least 8	n.s.	Wang et al., 2013
ERM-CD281	plant	70±3 ng g <sup>-1</sup>	9				ICP-MS	This study

<sup>&</sup>lt;sup>a</sup> Compiled value of average values of each study and uncertainty values (95% confidence limit)

**Table 5.** Comparison of Ge contents in different reference materials between this study and data from the literature. n.s. = not specified. Uncertainties are provided with 95% confidence level.

To provide the first Ge data on a plant standard, the standard ERM-CD281 was analyzed. Its Ge content is estimated to be approximately 70±3 ng g<sup>-1</sup> (n=9; Supplementary Table 1).

#### 3.3. Perspectives

This method provides a convenient procedure for determining accurate and precise Ge concentrations in rocks, soils and plants. Additionally, the combination of Ge concentrations with  $\delta^{30}$ Si analyses is the great advantage of this method, as both Si cycle tracers can be analyzed from a single procedure allowing a better understanding of the continental Si cycle. As emphasized in Cornelis et al. (2011), the combination of Ge/Si ratios and silicon isotopes turned out to be a great asset for tracing the sources and fate of Si in the critical zone, where numerous co-occurring processes are difficult to decipher.

The method can likely be expanded to water samples under the condition of respecting the detection limits. The purification step of this procedure would then be suitable for hydrothermal fluids (12±13 µg l<sup>-1</sup> up to 50 µg l<sup>-1</sup>; Evans and Derry, 2002) and soil solutions

 $(7.1\pm1.4~\mu g~\Gamma^{-1}; Cornelis et al., 2010)$ . However, the analysis of silicon isotopes on the same chemical preparation may not be straightforward because of DOC and anions responsible for isotopic bias (Hughes et al., 2011). Water samples with lower Ge contents, such as freshwater and seawater, might be pre-concentrated using Mg-coprecipation (Escoube, 2008). As many other elements are scavenged by this method, the Mg-coprecipitation should be followed by the cationic purification step. As both Ge and Si recoveries are complete during Mg-coprecipitation (Escoube, 2008; Cardinal et al., 2005), the analysis of silicon isotopes and Ge concentrations could again be performed on the same preparation. However, as anionic species will not be removed completely, care should be taken to control the bias induced by anionic species on  $\delta^{30}$ Si analysis (Hughes et al., 2011).

In addition to silicon isotopes analysis, this method paves the way towards a new approach to perform Ge isotopes is soils and plants, which has never been done before. Combining silicon isotopes, germanium isotopes and Ge/Si ratios would be a powerful multi-proxy approach to tackle continental Si biogeochemical cycle.

#### 4. Conclusions

The fusion and purification technique based on the Si isotopes preparation protocol (Georg et al., 2006) allows precise determination of Ge concentrations in environmental samples down to a detection limit of approximately 6 ng  $\Gamma^1$ . Analyses performed on  $\Gamma^4$ Ge, the most abundant isotope, show the best precision of approximately 3%. Accuracy, difficult to assess given the scarcity of reliable data available for comparison, was assessed with an error of approximately 5%. This method provides data with equivalent or better metrological characteristics than other published procedures and analytical techniques except for ID-HG-ICP-MS, which shows the best precision and lowest detection limit (Table 4) but is much more time consuming. In addition to the reliability of generated data, this technique has two major

advantages: (1) it does not require peculiar instruments such a hydride generation system, not
found in most laboratories unless specific needs, and (2) Ge concentrations and Si isotopic
analysis can be performed on the same solution. This last point represents a precious gain of
time in terms of sample preparation as one procedure serves to perform both analyses. A
systematic coupling of Ge/Si ratios and $\delta^{30}$ Si compositions is a great asset in the interpretation
of the Si cycle, as both tracers are very complementary.
Acknowledgments
This study has been funded by the ANR project BIOSiSOL (ANR-14-CE01-002). The French
Agence Nationale de la Recherche supported the study through the Project EQUIPEX
ASTER-CEREGE. Authors thank Hélène Mariot for her management of the clean laboratory.
References
Blecker S.W., King S. L., Derry L.A., Chadwick O.A., Ippolito J.A. and Kelly E.F. (2007)
Blecker S.W., King S. L., Derry L.A., Chadwick O.A., Ippolito J.A. and Kelly E.F. (2007)  The ratio of germanium to silicon in plant phytoliths: quantification of biological
The ratio of germanium to silicon in plant phytoliths: quantification of biological
The ratio of germanium to silicon in plant phytoliths: quantification of biological discrimination under controlled experimental conditions. Biogeochemistry, 86,189-199.
The ratio of germanium to silicon in plant phytoliths: quantification of biological discrimination under controlled experimental conditions. Biogeochemistry, 86,189-199.  Belissont R., Boiron MC., Luais B. and Cathelineau M. (2014) LA-ICP-MS analyses of
The ratio of germanium to silicon in plant phytoliths: quantification of biological discrimination under controlled experimental conditions. Biogeochemistry, 86,189-199.  Belissont R., Boiron MC., Luais B. and Cathelineau M. (2014) LA-ICP-MS analyses of minor and trace elements and bulk Ge isotopes in zoned Ge-rich sphalerites from the
The ratio of germanium to silicon in plant phytoliths: quantification of biological discrimination under controlled experimental conditions. Biogeochemistry, 86,189-199.  Belissont R., Boiron MC., Luais B. and Cathelineau M. (2014) LA-ICP-MS analyses of minor and trace elements and bulk Ge isotopes in zoned Ge-rich sphalerites from the Noailhac – Saint-Salvy deposit (France): Insights into incorporation mechanisms and ore

Global Biogeochemical cycles, 19, GB2007, doi:10.1029/2004GB002364.

- Cornelis J.-T., Delvaux B., Cardinal D., André L., Ranger J., Opfergelt S. (2010). Tracing mechanisms controlling the release of dissolved silicon in forest soil solutions using Si isotopes and Ge/Si ratios. Geochimica et Cosmochimica Acta 74 (2010) 3913–3924. Cornelis J.-T., Delvaux B., Georg R. B., Lucas Y., Ranger J., Opfergelt S. (2011) Tracing the origin of dissolved silicon transferred from various soil-plant systems towards rivers: a review. Biogeosciences, 8, 89–112. Delvigne C., Opfergelt S., Cardinal D., Delvaux B. and André L. (2009) Distinct silicon and germanium pathways in the soil-plant system: Evidence from banana and horsetail. Journal of Geophysical Research, 114, G02013, doi:10.1029/2008JG000899. Delvigne, C., Cardinal, D., Hofmann, A., André, L., (2012) Stratigraphic changes of Ge/Si, REE+Y and silicon isotopes as insights into the deposition of a Mesoarchaean banded iron formation. Earth and Planetary Science Letters, 355-356, 109–118. Delvigne C., Opfergelt S., Cardinal D., Hofmann A. and André, L., (2016) Desilication in Archean weathering processes traced by silicon isotopes and Ge/Si ratios. Chemical Geology, 420, 139-147. Derry L.A., Kurtz A. C., Ziegler K., and Chadwick O. A. (2005) Biological control of terrestrial silica cycling and export fluxes to watersheds. Nature, 433, 728-731. Dong L., Shen B., Lee C.-T. A., Shu X.-J., Peng Y., Sun Y., Tang Z., Rong H., Lang X., Ma H., Yang F. and Guo W. (2015) Germanium/silicon of the Ediacaran-Cambrian Laobao cherts: Implications for the bedded chert formation and paleoenvironment interpretations. Geochemistry, Geophysics, Geosystems, 16, 751–763, doi:10.1002/2014GC005595.
- Ellwood M.J., Kelly M., Maher W. A. and De Deckker P. (2006) Germanium incorporation into sponge spicules: Development of a proxy for reconstructing inorganic germanium

and silicon concentrations in seawater. Earth and Planetary Science Letters, 243, 749-
759.
Ellwood M.J. and Maher W.A. (2003) Germanium cycling in the waters across a frontal zone:
the Chatham Rise. New Zealand, Marine Chemistry 80, 145–159.
Evans M. J. and Derry L.A. (2002) Quartz control of high germanium/silicon ratios in
geothermal waters. Geology, 30, 1019-1022.
Escoube R. (2008) Iron and germanium isotope geochemistry in river and hydrothermal
systems. PhD Thesis, Université de Pau et des pays de l'Adour, 275 pages.
Escoube R., Rouxel O. J., Luais B., Ponzevera E. and Donard O. F.X. (2012) An
Intercomparison Study of the Germanium Isotope Composition of Geological Reference
Materials. Geostandards and Geoanalytical Research, 36, 149-159.
Filippelli G.M., Carnahan J.W., Derry L.A., and Kurtz A. (2000) Terrestrial paleorecords of
Ge/Si cycling derived from lake diatoms. Chemical Geology, 168, 9–26.
Frei R. and Polat A. (2007) Source heterogeneity for the major components of ~3.7 Ga
Banded Iron Formations (Isua Greenstone Belt, Western Greenland): Tracing the nature
of interacting water masses in BIF formation. Earth and Planetary Science Letters, 253,
266–281.
Fitoussi C., Bourdon B., Kleine T., Oberli F. and Reynolds B. C. (2009) Si isotope
systematics of meteorites and terrestrial peridotites: implications for Mg/Si fractionation
in the solar nebula and for Si in the Earth's core. Earth and Planetary Science Letters,
287, 77–85.

431	Georg R.B., Reynolds B.C., Frank M. and Halliday A.N. (2006) New sample preparation
432	techniques for the determination of Si isotopic compositions using MC-ICPMS. Chemical
433	Geology, 235, 95–104.
434	Hamade T., Konhauser K., Raiswell R., Goldsmith S. and Morris R. (2003) Using Ge/Si
435	ratios to decouple iron and silica fluxes in Precambrian banded iron formations. Geology,
436	31, 35- 38.
437	Hammond D. E., McManus J. and Berelson W. M. (2004) Oceanic germanium/silicon ratios:
438	Evaluation of the potential overprint of temperature on weathering signals.
439	Paleoceanograpy, 19, doi/10.1029/2003PA000940.
440	Hughes H.H., Delvigne C., Korntheuer M., de Jong J., André L. and Cardinal D. (2011)
441	Controlling the mass bias introduced by anionic and organic matrices in silicon isotopic
442	measurements by MC-ICP-MS. Journal of Analytical Atomic Spectrometry, 26, 1892-
443	1896.
444	Höll R., Kling M. and Schroll E. (2007) Metallogenesis of germanium—A review. Ore
445	Geology Reviews, 30, 145–180.
446	IUPAC (1997). Compendium of Chemical Terminology, 2nd ed. (the "Gold Book").
447	Compiled by A. D. McNaught and A. Wilkinson. Blackwell Scientific Publications,
448	Oxford.
449	Jochum K.P., Weis U., Schwager B., Stoll B., Wilson S.A., Haug G.H., Andreae M.O. and
450	Enzweiler J. (2016) Reference Values Following ISO Guidelines for Frequently
451	Requested Rock Reference Materials. Geostandards and Geoanalytical Research, 40,
452	333-350.

Kurtz A.C., Lugolobi F., and Salvucci G. (2011) Germanium-silicon as a flow path tracer: Application to the Rio Icacos watershed. Water Resources Research, 47, W06516, doi:10.1029/2010WR009853, 2011 Kurtz A. C., Derry L. A., and Chadwick O. A. (2002) Germanium-silicon fractionation in the weathering environment. Geochimica et Cosmochimica Acta, 66,1525-1537. Kurtz A. C. and Derry L. A. (2004) Tracing silicate weathering and terrestrial silica cycling with Ge/Si ratios. In Wanty, R. B. and Seal, R. R., editors, Proc. 11th Int. Symp. on Water Rock Interaction, 833-836, The Netherlands. Balkema Pubs. Luais B. (2007) Isotopic fractionation of germanium in iron meteorites: Significance for nebular condensation, core formation and impact processes. Earth and Planetary Science Letters, 262, 21-36. Luais B. (2012) Germanium chemistry and MC-ICPMS isotopic measurements of Fe-Ni, Zn alloys and silicate matrices: Insights into deep Earth processes. Chemical Geology, 334, 295-311. Liu Y., Diwu C., Zhao Y., Liu X., Yuan H., and Wang J. (2014) Determination of trace and rare-earth elements in Chinese soil and clay reference materials by ICP-MS. Chinese Journal of Geochemistry, 33, 095–102. Lugolobi F., Kurtz A.C. and Derry L.A. (2010) Germanium-silicon fractionation in a tropical, granitic weathering environment. Geochimica et Cosmochimica Acta, 74, 1294–1308. Makishima A. and Nakamura E. (2009) Determination of Ge, As, Se and Te in Silicate Samples Using Isotope Dilution-Internal Standardisation Octopole Reaction Cell ICP-QMS by Normal Sample Nebulisation. Geostandards and Geoanalytical Research, 33, 369-384.

477	food samples including Chinese herbal remedies using graphite furnace atomic
478	absorption spectroscopy. Food Chemistry, 97, 411-417.
479	Mortlock R.A. and Froelich P.N. (1996) Determination of germanium by isotope dilution-
480	hydride generation inductively coupled plasma mass spectrometry. Analytica Chimica
481	Acta, 332, 277-284.
482	Mortlock R. A., Froelich P. N., Feely R. A., Massoth G. J., Butterfield D. A., and Lupton J. E.
483	(1993) Silica and germanium in Pacific-ocean hydrothermal vents and plumes. Earth and
484	Planetary Science Letters, 119, 365-378.
485	Opfergelt S., Cardinal D., André L., Delvigne C., Bremond L. and Delvaux B. (2010)
486	Variations of $\delta^{30} \text{Si}$ and Ge/Si with weathering and biogenic input in tropical basaltic ash
487	soils under monoculture. Geochimica et Cosmochimica Acta, 74, 225-240.
488	Pokrovski G.S., and Schott J. (1998) Thermodynamic properties of aqueous Ge (IV)
489	hydroxide complexes from 25 to 350°C: Implications for the behavior of germanium and
490	the Ge/Si ratio in hydrothermal fluids. Geochimica et Cosmochimica Acta, 62, 1631-
491	1642.
492	Qi HW., Rouxel O., Hu RZ., Bi XW., Wen HJ. (2011) Germanium isotopic systematics
493	in Ge-rich coal from the Lincang Ge deposit, Yunnan, Southwestern China. Chemical
494	Geology, 286, 252-265.
495	Rouxel O., and Luais B. (2017) Germanium isotope geochemistry. Reviews in Mineralogy
496	and Geochemistry, 82, 601-656.
497	Rosenberg E. (2009) Germanium: environmental occurrence, importance and speciation.
498	Reviews in Environmental Science Bio/Technology, 8, 29-57.

499	Scribner A. M., Kurtz A. C., and Chadwick O. A. (2006) Germanium sequestration by soil:
500	Targeting the roles of secondary clays and Fe-oxyhydroxides. Earth and Planetary
501	Science Letters, 243, 760-770.
502	Shen B., Lee CT. A., Xiao S. (2011) Germanium/silica ratios in diagenetic chert nodules
503	from the Ediacaran Doushantuo Formation, South China. Chemical Geology, 280, 323-
504	335.
505	Soylak M. and Yigit S. (2015) Preconcentration-separation of germanium at ultra trace levels
506	on polysulfone membrane filter and its determination by spectrophotometry. Journal of
507	Industrial and Engineering Chemistry, 24, 322-325.
508	Sutton J., Ellwood M.J., Maher W.A., and Croot P.L. (2010) Oceanic distribution of inorganic
509	germanium relative to silicon: Germanium discrimination by diatoms. Global
510	Biogeochemical Cycles, 24, GB2017, doi:10.1029/2009GB003689.
511	Tribovillard N., Bout-Roumazeilles V., Riboulleau A., Baudin F., Danelian T. and Riquier L.
512	(2011) Transfer of germanium to marine sediments: Insights from its accumulation in
513	radiolarites and authigenic capture under reducing conditions. Some examples through
514	geological ages. Chemical Geology, 282,120–130.
515	Wang Y., Gu T., Wang X., Gao YS., Jochum K.P., Müller W.E.G. (2013) Practical
516	handbook of reference materials for geoanalysis, second, Geological Publishing House,
517	Beijing.
518	Wiche O., Zertani V., Hentschel W., Achtzige R. and Midula P. (2017) Germanium and rare
519	earth elements in topsoil and soil-grown plants on different land use types in the mining
520	area of Freiberg (Germany). Journal Geochemical Exploration, 175, 120-129.

.rc. Geoch Wheat C.G. and McManus J. (2005) The potential role of ridge-flank hydrothermal systems