

Recent improvements in the differential scanning calorimetry methods applied to the study of gas hydrates.

Frédéric Plantier, Jean-Philippe Torre, Laurent Marlin, Rémi André, Pierre Le Parlouer

▶ To cite this version:

Frédéric Plantier, Jean-Philippe Torre, Laurent Marlin, Rémi André, Pierre Le Parlouer. Recent improvements in the differential scanning calorimetry methods applied to the study of gas hydrates.. 20th European Conference on Thermophysical Properties (ECTP 2014), Aug 2014, Porto, Portugal. hal-02009661

HAL Id: hal-02009661

https://hal.science/hal-02009661

Submitted on 6 Feb 2019

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Open Archive Toulouse Archive Ouverte

OATAO is an open access repo	sitory that	collects t	the work of	Toulouse
researchers and makes it freely	y available	over the	web where	possible

T : : : : : : : :			1 44 // 4		(/0 / 0 = /
This is an author's	VARCIAN	nubliched in	n' http://aata/	A LIMIV/ #ALIJALICA	ヽ tr/'ノ1()ん1
	version	DUDHSHED II	1 11110 //Oala		. 11// 195) 1
i i iio io ai i aati ioi o	* 0.0.0	pasiionoa n		orarni toaroact), _ . 00 .

To cite this version:

Plantier, Frédéric and Torré, Jean-Philippe and Marlin, Laurent and André, Rémi and Le parlouer, Pierre Recent improvements in the differential scanning calorimetry methods applied to the study of gas hydrates. (2014) In: 20th European Conference on Thermophysical Properties, 31 August 2014 - 4 September 2014 (Porto, Portugal). (Unpublished)











Aquitaine

SETARAM Instrumentation

RECENT IMPROVEMENTS IN THE HIGH PRESSURE DIFFERENTIAL SCANNING CALORIMETRY METHOD APPLIED TO THE STUDY OF GAS **HYDRATES**

Frédéric Plantier a,*, Jean-Philippe Torré a, Laurent Marlin a, André Rémib, Le Parlouer Pierreb

Outline



Thermodynamic properties of gas hydrates (e.g. phase equilibrium data, phase change enthalpies, specific heat, etc.) can be obtained by using high pressure calorimetric techniques.


However, in some cases, one of the major drawback to existing devices is the absence of in-situ agitation leading to problems such as efficient gas solubilization, long induction times, formation of an hydrate crust covering the gas/liquid interface, low hydrate to water conversion, etc.)

Development, tests and validation in our laboratory (LFCR, France) of a novel type of calorimetric cell for BT 2.15 and C80 **SETARAM calorimeters:**




industrialized and commercialized by **SETARAM Instrumentation**

Technical details of the prototype




- ◆ The mechanical agitation in the measuring cell is provided by a removable rotary shaft on which star washers have been welded (N < 200 RPM)
- ◆ The link between the wellhead and the measuring cell is provided by two concentric tubes allowing the agitator shaft to rotate freely in the middle of the tubes.
- The space let between the agitator shaft and the first tube is used for the admission of the fluids into the cell, while the annular space between the two concentric tubes allows flowing the fluids to the outlet.
- ◆ The direction of the fluid flow is completely reversible and the pressure inside the cell can be dynamic controlled during the experiment.



- (A) Picture of the prototype
- (B) High pressure cell and its agitation system
- (C) Prototype installed on a SETARAM BT 2.15 calorimeter

Dynamic control of the pressure



Fluid flow from inlet to outle (through the measuring cell)

→ Benefits:


- measurements under pressure (P < 20 Mpa)
- in-situ mechanical agitation
- dynamic control of the pressure inside the cell
- designed for BT 2.15 and C 80 SETARAM calorimeters

Experimental results

Calibration of the prototype in temperature and enthalpy


reference (NIST)	T _{fusion} _ref (°C)	T _{fusion} _exp (°C)
Water	$0,00 \pm 0,05$	0.1 ± 0.1
Cyclohexane	6.5 ± 0.3	6.7 ± 0.2
n-C ₁₂	-9.7 ± 0.3	-9.9 ± 0.2
n-C ₁₄	5.6 ± 0.9	5.3 ± 0.2
n-C ₁₆	18 ± 1	$18,1 \pm 0,2$

Precision in temperature : ± 0.2 ℃ {-10 < T (℃) < 18 }




Sensitivity constant: $K = (3.46 \pm 0.05).10^{-5} \text{ J.} \mu\text{V}^{-1}.\text{s}^{-1}$

Application to CO₂ hydrates $P = 3.05 \pm 0.05 MPa$



- the presence of the in-situ agitation triggers the crystallization (rupture of the metastability, reduction of the induction time)
- the stirring process no creates noise/perturbation of the calorimetric signal
 a good reproducibility in the results is obtained

Comparison of our results to literature data



EXP (P_{CO2} = 3.05 ± 0.05 MPa) :

number (CO₂ – N_H . .H₂0) equal to 7.3 ± 0.13 [5] s total (n = 100 %)



EXP = 62 ± 3 (kJ / mol_{cos})

References

- Conclusions
- A novel high pressure calorimetric cell equipped with a mechanical agitation and a dynamic pressure control has been developed and tested. • The temperature and enthalpy calibrations have been performed; the precision and the sensitivity constant have been determined.
- the cell prototype has been used to obtain equilibrium data and phase change enthalpy of the CO2 hydrate
- The experimental results obtained are in good agreement with literature data, and demonstrate the potentialities of this novel equipment