Development of a peptidomic workflow to apprehend dynamic digestion of a mix of proteins: application of meat product for elderly

¹Sayd T., ¹Theron L., ¹Chambon C., ¹Sicard J., ¹Mirade P.-S., ²Peyron M. A., ³Pinguet J., ³Richard D & ¹Santé-Lhoutellier V.

¹ INRA, UR370 Qualité des Produits Animaux, F-63122 Saint Genès Champanelle, France

² INRA, UMR1019, Unité de Nutrition Humaine, CRNH Auvergne, F-63000 Clermont-Ferrand, France

³ CHU Clermont-Ferrand, Service de Pharmacologie Médicale, F-63003 Clermont-Ferrand, France

Contact: thierry.sayd@inra.fr

Meat is known to be a source of protein of good nutritional quality in terms of composition in essential amino acids and digestibility. Evaluation of meat proteins digestibility is a critical issue in regard of various meat product, industrial process and according the digestive physiological status of the consumer. The peptidomic analysis is of a great interest to evaluate this digestibility. The major difficulties of these studies are reliable to the complexity of meat protein composition and kinetics of the digestion. Therefore, creating a workflow of various techniques (*in vitro* digestion model, mass spectrometry, statistic and bioinformatics) would allow to access to a multivariate proteolysis model.

The main objective of this work was to create a workflow relevant for analysing peptides from a complex meal during digestion. Meat product was chosen as an example

<u>Dynamic digestive conditions</u>: Two *in vitro* gastric physiological conditions of digestion were applied using DIDGI©: adult / elderly. Three repetitions of digestion were done. Five sampling were realized in gastric phase at 30, 60, 90, 120 and 150 min to analyze the kinetics of peptide production of meat.

<u>Mass spectrometry analysis</u>: Identification and quantification of peptides generated by gastric digestion of meat model were performed using LC-MS/MS analysis (Q Exactive plus Thermofisher) and Progenesis QI (Waters) and/or Mascot (www.matrixscience.com).

<u>Statistics/Bioinformatics:</u> Analysis of kinetic were realized using hierarchical clustering nonsupervised method and a fitting method was tested to compare proteins kinetics. Peptigram a web application was used to visualize protein hydrolysis.

This combined approach permitted to analyse senior and adult protein digestion. Specific protein hydrolysis kinetic were obtained.