

Development of a simple, low-cost and rapid thin-layer chromatography method for the determination of individual volatile fatty acids

F. Robert-Peillard, Elodie Mattio, Aïnhoa Komino, Jean-Luc Boudenne,

Bruno Coulomb

▶ To cite this version:

F. Robert-Peillard, Elodie Mattio, Aïnhoa Komino, Jean-Luc Boudenne, Bruno Coulomb. Development of a simple, low-cost and rapid thin-layer chromatography method for the determination of individual volatile fatty acids. Analytical Methods, 2019, 10.1039/c9ay00158a. hal-02090448

HAL Id: hal-02090448 https://hal.science/hal-02090448

Submitted on 4 Apr 2019

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

Development of a simple, low-cost and rapid thin-layer chromatography method for the determination of individual volatile fatty acids

- 3
- 4 Fabien Robert-Peillard^{*}, Elodie Mattio, Aïnhoa Komino, Jean-Luc Boudenne, Bruno Coulomb
- 5 Aix Marseille Univ, CNRS, LCE, Marseille, France.
- 6 ^{*}Corresponding author: fabien.robert-peillard@univ-amu.fr
- 7

8 Abstract

9 In this paper, a new thin-layer chromatography (TLC) method for the determination of individual volatile fatty acids (VFA) is presented. The experimental procedure is based on the 10 derivatization of VFA with activating agents and a naphthalene-based amine to form an 11 amide derivative, followed by fast TLC separation of these VFA derivatives and visualization 12 with phosphomolybdic acid at room temperature. Quantification is then performed by 13 densitometric measurements using a free software after TLC pictures in a home-made black 14 box for homogeneous and reproducible lighting of the plates. Optimization of the 15 experimental conditions enables fast determination (20 minutes for 3 samples) of VFA with 16 detection limits (1.5-2.5 mg.L⁻¹) and relative standard deviations (4.6-6.7 %) satisfactory 17 enough for our analytical purpose using low-cost instruments with easy to find materials and 18 19 software. The protocol was then applied to real samples (aqueous extracts of sewage sludge composts) and results were compared with gas chromatography analysis with good 20 21 agreement between the two methods.

Keywords: Volatile fatty acids, thin-layer chromatography, low-cost, environmental analysis.

22 23 24

25

1. Introduction

Anaerobic digestion of biomass leading to biogas production has raised considerable 26 27 attention over the past years, with applications such as wastewater [1], sewage sludge [2], 28 organic wastes [3] or landfill leachates [4]. The challenge for anaerobic digestion processes 29 is to maintain a stable fermentation process and consistent biogas production, and the 30 monitoring of several parameters such as pH, alkalinity or biogas production rate has been 31 proposed to this end [5]. However, it has been widely acknowledged that volatile fatty acids (VFA) concentrations (in particular acetic, propionic, isobutyric, butyric, isovaleric and valeric 32 acid) provide the most meaningful measure for assessing the effectiveness of the digestion 33 process [6]. Indeed, VFA are produced and consumed during the main stages of anaerobic 34

digestion processes (hydrolysis, acidogenesis, acetogenesis, and methanation), and this
 parameter therefore reflects best the metabolic state of the biochemical process.

VFA measurement can provide two types of information: total VFA (sum of all VFA 37 38 expressed as mg/L as acetic acid) or speciation of individual VFA. Total VFA is often used as 39 control parameter as relatively simple analytical methods such as titration [7], near-infrared spectroscopy [8] or spectrophotometry [9] are available for research or operational 40 41 laboratories. Despite this simplicity, total VFA measurements can be insufficient for precise optimization of anaerobic digesters, as understanding of the composition of VFA during all 42 stages of the digestion process is also a key point. Indeed, VFA speciation results enable 43 44 operators to monitor the current state of the anaerobic digesters and adjust the operational 45 strategies to improve performance. Accumulation of propionate is a common indicator of stress in an anaerobic digester [10-12], and traditional non-differentiating methods like 46 titration cannot give early alert for this potential stress that may lead to detrimental bacterial 47 48 community change and process failure. Moreover, analytical techniques which give individual VFA composition also show strong interests for bioplastic production field which rely on VFA 49 as carbon source, for example for polyhydroxyalkanoates production [13]. It can thus be 50 concluded that there is an urgent need for simple analytical tools for speciation of VFA in 51 operational laboratories monitoring anaerobic digesters performances. 52

Traditionally, speciation of VFA is performed by gas chromatography, liquid chromatography 53 54 or capillary electrophoresis [14-17]. These methods based on separation techniques can all 55 provide reliable results, but their relative complexity requiring skilled personnel and high purchase and operational costs (separation columns, eluents, gases...) strongly limit their 56 57 use for operational laboratories which prefer simple and low-cost analytical tools for routine 58 analysis. Among separation techniques, Thin-Layer Chromatography (TLC) is the oldest one 59 and certainly the simplest and least expensive tool in analytical chemistry. Many studies have already proposed TLC as an alternative to more modern and complex chromatographic 60 system for quality control and fast quantitative or semi-quantitative determination [18-21]. 61 Due to its simplicity, sample clean-up can be minimized, and cross-contaminations or 62 columns clogging are avoided, which are significant advantages for complex samples such 63 64 as sewage sludge or wastewater.

Our research group has developed an original procedure for the determination of total VFA by batch procedure or microplate assays [22,23]. The analytical methodology is based on a specific derivatization of VFA via a two-step mechanism: (i) activation step with two activating agents at acidic pH; (ii) amidation step with a fluorescent amine at basic pH. This methodology was also applied for the development of a multi-syringe chromatography system for the determination of individual VFA [24]. This system was based on the separation of VFA derivatives on a C_{18} monolithic column and quantification by UV or fluorescence detection, and it could serve as a relatively simple and low-cost on-line analyzer
 compared to complex on-line chromatographic systems.

74 In the present paper, a very simple and rapid TLC method for the determination of individual volatile fatty acids is presented. This method is intended as an analytical tool for operational 75 76 laboratories monitoring anaerobic digesters performances or other environmental samples who cannot afford HPLC or GC but require routine, rapid and low-cost measurement of 77 78 individual VFA (without sending samples to external laboratories). The procedure is based on previously-mentioned derivatization of VFA combined with fast TLC separation of VFA 79 80 derivatives and visualization with phosphomolybdic acid at room temperature. Development 81 and application of this procedure for quantitative determination of VFA are presented herein.

82 83

2. Experimental

84

85 **2.1 Reagents, solutions and materials**

All chemicals and solvents were of analytical or HPLC reagent grade and used without 86 further purification. Cyclohexane and ethyl acetate were purchased from Sigma-Aldrich. 87 "EDC solution" was a solution of the activating reagent N'-(3-dimethylaminopropyl)-N-88 ethylcarbodiimide (EDC) (Sigma-Aldrich) prepared at 100 g.L⁻¹ in absolute ethanol (EtOH) 89 (Sigma-Aldrich). "HOAT+EDAN solution" was a solution of the secondary activating reagent 90 7-aza-1-hydroxybenzotriazole (HOAT) (Genscript Corporation, Piscataway, USA) at 6 g.L⁻¹ 91 and the derivatizing reagent N-(1-naphthyl)ethylenediamine (EDAN) (Sigma-Aldrich) at 2 g.L⁻ 92 ¹, prepared as a mixed solution in ultrapure water (Millipore, resistivity >18 M Ω cm) with final 93 94 pH set at 3.6. "KH₂PO₄ solution" was made of potassium dihydrogen phosphate (KH₂PO₄) 95 prepared at 250 mM in NaOH 0.35 M. "PMA solution" was prepared by dissolving phosphomolybdic acid hydrate (PMA) (Sigma-Aldrich) at 3 % w/v in water. Stock solutions of 96 97 carboxylic acid standards (5000 mg L-1 for acetic acid, 1500 mg L-1 for other VFA) (Acros Organics, Geel, Belgium) were prepared in ultrapure water. 98

99 Thin layer chromatography aluminum plates (20x20 cm, 0.25 mm layer thickness) precoated 100 with silica gel 60F-254 were obtained from Merck (Darmstadt, Germany) and cut at adequate 101 dimensions for intended use (6x6 cm). Spots were applied on the plates at 1.3 cm height 102 from lower edge of the plate and were separated by a minimum distance of 0.8 cm (6 spots 103 possible on a single plate).

104

105 **2.2 Samples**

Real samples used for validation were aqueous extracts of composts (extraction protocol
described in Supplementary material) made by mixing anaerobically digested sewage sludge
(from Marseille's sewage sludge treatment plant) with green wastes (1/3-2/3, v/v). Ages of

the composts were between 1 month (AN1) and 6 months (AN 5). All samples were filtered
through 0.45 µm membrane filters (PES membrane, Millipore, France).

111

112 **2.3 Analytical protocol for VFA determination**

In a 2 mL glass vial were added 300 µL of "HOAT+EDAN solution", 250 µL of sample and 40 113 µL of "EDC solution". The vial was stirred and left to react for 5 minutes. 75 µL of "KH₂PO₄ 114 115 solution" were then added and the vial was again left to react for 1 minute after proper stirring. 3 µL of the reaction mixture was then manually spotted on the TLC plate using 3 µL 116 117 micropipettes (Brand) and the spots were dried with a simple hair-dryer set at minimum heating power until evaporation of water. The plates were developed vertically for 5 min 118 119 (solvent front at about 0.5 cm from upper edge of the plate) at room temperature (22°C) in a 100 mL polypropylene disposable tank equipped with a plastic top (Tecan) in order to 120 121 minimize organic solvent volumes used for development (1 mL is enough in these tanks; classical TLC development chambers with higher volumes can be used). Mobile phase was a 122 mixture of ethyl acetate and cyclohexane (60:40, v/v). After development, the plates were 123 dried with the hair-dryer for 1 min. Spots were visualized by dipping the plate in the PMA 124 125 solution for 10 s (in a 100 mL polypropylene disposable trough) and drying the plate on absorbent paper. A picture of the plate was taken immediately with a smartphone in a home-126 127 made black box protected from ambient light with a hole for smartphone camera. 128 Homogeneous and reproducible lighting of the plates in the box was provided by a LED strip 129 (12 V DC, 400 Lm, 24 LED/m) cut at the dimensions of the box and stuck on the upper part 130 of the box. Drawing of the box on CAD software is presented in supplementary material section. Quantitative determination of VFA was performed on a personal computer using the 131 open-source ImageJ software (National Institutes of Health) and its densitometric analysis 132 tool. Pictures were analyzed without modification using a constant rectangular selection size 133 between calibration plates and sample plates. Gel analysis tools were used to trace 134 densitometric histograms of each VFA lane and peaks areas were calculated using the wand 135 tool (results can be saved as an excel file for further analysis). 136

137

138 **2.4 Retention factors and resolution between VFA spots**

The retention factor R_f was determined directly on the TLC picture and was defined as the
ratio of the distance traveled by the center of a spot to the distance traveled by the solvent
front.

The resolutions between two neighboring spots A and B ($R_{S(A/B)}$) were calculated on the densitometric histograms given by the gel analysis tools of Image J. Resolutions were calculated as the ratio between the difference of distance of peaks maximum and the average of peak widths at base.

147 **2.5 Gas chromatography analyses**

VFA were determined under classical GC conditions using a Varian CP-3800 gas 148 chromatograph (GC) with a free fatty acid phase (FFAP) 25 m fused silica capillary column of 149 0.1 µm film thickness and an i.d. of 0.25 mm (Varian CP-WAX 58-CB). The GC oven 150 temperature was programmed to be incremented from 115 to 140 °C within 5 min, and a final 151 temperature of 140 °C was then held for 1.25 min. The injector and detector temperature 152 were set at 290 °C and 300 °C, respectively. Helium was used as carrier gas at a pressure of 153 154 1.4 bar. A flame ionisation detector (FID) was used and quantification was performed using internal standard method. Limits of detection between 3 mg L⁻¹ (valeric acid) and 10 mg L⁻¹ 155 (acetic acid) were obtained. All sample analyses were performed in duplicate 156

157

158 **3. Results**

159

Experimental conditions for the derivatization of VFA were similar as those used in a previous study with the Multi-Syringe Chromatography set-up [24]. The only differences are the concentrations of the reagents (100 g.L-1 for EDC and 2 g.L-1 for EDAN) which have been increased in order to gain sensitivity. The main experimental optimizations for the work described here are related to the TLC application which gives full interest to the analytical protocol and its intended use.

166

167 **3.1 Visualization solution of the spots**

168 Many visualization methods (or stains) can be used for TLC, which can be compound 169 specific (Dragendorff reagent, ninhydrin...) or universal (iodine, vanillin, phosphomolybdic 170 acid...), with varying preparation complexity and hazardous nature of reagents (strong acids or toxic reagents). As the method is intended for operational laboratories, visualization 171 method with adequate sensitivity, simple preparation procedure and relatively harmless 172 reagents is preferable. Phosphomolybdic acid (PMA) was thus selected based on these 173 criteria, as it is a well-known, sensitive and cheap reagent with low toxicity and wich doesn't 174 require strong acids for the preparation of the visualization solution. Usually, PMA is used as 175 176 a solution in EtOH and analytes are visualized as dark green to black spots after heating with a heat gun or a hot plate. In our first experiments, we noticed that PMA solution in water was 177 also suitable for visualization of our VFA derivatives spots and that no heating was required 178 after dipping plates in aqueous PMA, which prevents risks associated with the use of heating 179 gun and flammable organic solvents. Moreover, naked-eye comparison of spots visualized 180 with PMA in EtOH or PMA in water (Figure 1) proved that aqueous PMA was the best choice 181 for our application as the spots were much more visible and darker. 182

184 **3.2** Mobile phase selection and retention factors

Many organic solvents are used for TLC development with silica gel plates (the most 185 common and least expensive stationary phase), including highly toxic solvents such as 186 187 chloroform, benzene, toluene or hexane. Nowadays, a classical mixture of organic solvents used for TLC development is cyclohexane/ethyl acetate, due its good efficiency to separate 188 compounds with various polarities and functional groups, and also to its moderate toxicity 189 compared to previously-mentioned solvents. Mobile phases with a mixture of these two 190 191 solvents in various proportions were optimized on a VFA solution containing acetic, 192 propionic, butyric and valeric acids at 100 mg.L⁻¹ (Figure 2). Resolutions between each VFA were calculated accurately on densitometric histograms plotted from the TLC pictures and 193 are given in Table 1. As can be seen, mobile phase made with ethyl acetate/cyclohexane 194 with volume composition of 60:40 provided the best results with optimal separation between 195 196 the closest spots (propionic, butyric and valeric acid). Addition of other classical TLC mobile phase components like formic acid, diethylamine, ethanol or ammonium hydroxide didn't 197 improve results. Retention factors with this mobile phase (60:40, v/v) were also calculated for 198 comparison with real samples: acetic acid: 0.24; propionic acid: 0.43; butyric (+isobutyric) 199 200 acid: 0.55; valeric (+isovaleric) acid: 0.64.

201

202 3.3 Interferences

The optimized protocol with TLC separation was applied to a wide range of organic 203 compounds likely to be detected under the reaction conditions such as other carboxylic acids 204 205 (formic, lactic, tartaric, succinic, oxalic, malic, maleic, citric, pyruvic, benzoic acid), 206 carbohydrate (glucose) or amino acids (glycine, alanine, cysteine). Spots of potential 207 derivatives or native compounds visualized at the end of the protocol were compared to VFA 208 spots (comparison of R_f and intensity). Most of potential interfering compounds didn't lead to any spots after derivatization and TLC visualization. Only formic, lactic and succinic acid 209 were visible on the TLC pictures (Figure 3), with concentrations 2.5 times higher than VFA in 210 order to be able to visualize corresponding spots. Succinic acid gives the highest response 211 but R_f of its derivatives is much lower than VFA. Formic and lactic acid derivatives have 212 similar R_f to acetic acid, but intensity of the response is much lower. If we compare these 213 carboxylic acids with acetic acids at the same concentration, intensity response obtained with 214 densitometric histograms is <4% of that of acetic acid. Other carboxylic acids tested (tartaric, 215 oxalic, malic, maleic, citric, pyruvic, benzoic acid) led to no appreciable response even at 500 216 mg.L⁻¹. Considering the expected prevalence of VFA in real samples analyzed by this 217 method, it is reasonable to conclude that VFA measurements with this TLC protocol should 218 be free of interference from other compounds. 219

221 3.4 Analytical features

Analytical features of the method are summarized in Table 2, with TLC picture and 222 associated densitometric histograms of a VFA calibration range between 10 and 125 mg.L⁻¹ 223 represented on Figure 4 (left part). Polynomial regression equations were obtained, probably 224 due to the visualization and detection method. Calibration ranges can be easily modified 225 (increase of the lowest and highest bound of the calibration range) by reducing sample 226 volume, depending on real samples analyzed. For example, reducing sample volume to 25 227 μ L enables an extension of the calibration ranges above 1 g.L⁻¹, which will be more adapted 228 229 to VFA concentrations in sewage sludge samples. Detection limits and quantification limits have been evaluated based on repetitions (n=7) of a sample spiked at 7 mg.L⁻¹ of each VFA, 230 following U.S. EPA guidelines [25]. Detection limits and quantification limits were between 231 1.5-2.5 mg.L⁻¹ and 4.7-7.9 mg.L⁻¹, respectively. These limits were clearly low enough for our 232 analytical purposes and were in the same range as those obtained by the GC-FID method 233 234 used in our laboratory. Finally, relative standard deviation (RSD) between 4.6 and 6.7 % were obtained after application of the whole protocol on a 50 mg.L⁻¹ standard (n=7, with one 235 TLC plate for each repetition). Precision with our method is therefore lower than that usually 236 obtained with TLC scanner and dedicated data analysis software [18, 26], but is satisfactory 237 enough for our analytical purpose using low-cost instruments with easy to find materials and 238 software. 239

240

241 3.5 Analytical costs

The aim of the method described in this article is to provide a low-cost analytical tool to operational laboratories for which GC or HPLC is too expensive, both as investment and routine usage. Investment for specific equipment for our method is very low, as only standard vials, glassware and pipettes are necessary. Operational costs are summarized in Table 3 and each sample analysis costs globally only 0.12 €, which makes this method indeed very affordable

248

249 **3.6 Comparison with other TLC method for VFA determination**

Direct TLC detection of VFA (without derivatization) is of course possible but detection limits are obviously too high (g.L⁻¹ range [27]) for our analytical application. Derivatization prior to TLC separation and detection is therefore mandatory to reach mg.L⁻¹ detection range. In spite of the clear advantages that TLC techniques could provide for VFA determination (simplicity, no sample clean-up and cross contaminations, cost), only one method has been published in the literature for derivatization of VFA and ppm detection. This derivatization method is very tedious and time-consuming, as it requires neutralization of the sample, evaporation to dryness, boiling for 1 h with derivatization reagent, extraction with chloroform and evaporation to dryness of organic extracts [28]. In comparison, our method only requires reaction at room temperature for 6 min with non-toxic reagents and direct spotting of the mixture without evaporation (analytical step which has to be avoided for simplicity reasons for routine analysis as described in this study).

Therefore, no TLC method suitable for easy and practically applicable determination of 262 263 individual VFA were described so far in the literature. Regarding our method, staff effort is moderate, as it takes only about 20 minutes to perform the whole analytical protocol (from 264 265 derivatization to TLC final picture) for 3 samples (or 6 if no duplicates is performed), and 266 technicity required for the routine application of the protocol is quite low, as no complex 267 instruments are needed. This method is therefore well suited in laboratories where the global operational needs and costs for analytical tools are decisive for the routine and precise 268 269 monitoring of their process.

270

271 **3.7 Validation on real samples**

A comparative study between our TLC method and GC-FID analysis has been conducted on 272 real samples. This validation study was performed on a set of aqueous extracts of composts 273 274 (made from anaerobically digested sewage sludge and green wastes) which were analyzed in our lab for their VFA content in the frame of a research project dealing with compost 275 monitoring and early warning parameters (VFA was one of the studied parameters, as 276 277 considerable variation of VFA has already been demonstrated in composting platform [29]). TLC picture and densitometric histograms used for quantification are shown on Figure 4 278 279 (right part) for comparison with the calibration range used (scales for the left and right 280 histograms are different and peak heights can't be compared visually). As expected, it can be observed visually that acetic acid is the prevalent VFA in the compost extracts with a 281 282 noticeable evolution of VFA during the composting process (increase from AN1 to AN3 during bio-oxidation phase and decrease for older composts during maturation phase). 283 284 Quantitative results for individual C2–C5 VFA concentrations measured by our method or by GC-FID are detailed in Table 4. Results between both methods are in good agreement, with 285 only slight underestimation by our TLC method. Standard deviations are also higher for the 286 287 TLC method, but results are globally very satisfactory if we consider the simplicity of the material used for visualization and quantification. Note that semi-quantitative measurements 288 could be possible by simple naked-eye comparison with the standard TLC plates for 289 operational laboratories which just need early warning of a sudden change in VFA 290 291 concentration or composition.

292

4. Conclusions

We have demonstrated that individual VFA determination and quantification using one of the 294 simplest separation technique (thin-layer chromatography) is indeed possible for 295 296 environmental analysis. Derivatization of VFA using a simple and rapid procedure at room 297 temperature leads to the formation of derivatives which can be easily separated on silica gel 298 TLC plates and detected without heating by simple dipping in a visualization solution. Use of a home-made black box enables reproducible smartphone pictures for subsequent 299 300 quantification with a free software. As investment and operational costs are very low, this analytical method is perfectly suited for operational laboratories which can't afford on-site 301 302 HPLC or GC analysis but require routine, rapid and low-cost measurement of individual VFA. 303 The proposed method could therefore supplement classical titrimetric methods which only 304 gives access to total VFA with no details on VFA composition, although this parameter is more and more recognized for its meaningful significance for the monitoring of anaerobic 305 306 digestion processes.

307

308 Acknowledgment

This work was financially supported by the French Environment and Energy Management Agency (ADEME) through the programme "ODIAMAC" [Convention n°1506C0034].

311

312 Literature cited

313

[1] N. Narkis, S. Henefeld-Fourrier, M. Rebhun, Volatile organic acids in raw wastewater and
in physicochemical treatment. Water Res., 14 (1980), 1215-1223.

316

317 [2] B. Ruffino, G.Campo, G.Genon, E. Lorenzi, D.Novarino, G.Scibilia, M. Zanetti,

Improvement of anaerobic digestion of sewage sludge in a wastewater treatment plant by

319 means of mechanical and thermal pre-treatments: Performance, energy and economical

320 assessment. Bioresource Technol., 175 (2015), 298-308.

321

[3] B. Chatterjee, D. Mazumder, Anaerobic digestion for the stabilization of the organic
fraction of municipal solid waste: A review. Environ. Rev. 24 (2016), 426-459.

324

[4] M. Kawai, M. Kishi, M.R. Hamersley, N. Nagao, J. Hermana, T. Toda, Biodegradability
and methane productivity during anaerobic co-digestion of refractory leachate. Int. Biodeter.
Biodegrad., 72 (2012), 46-51.

328

[5] K. Boe, D.J. Batstone, J.P. Steyer, I. Angelidaki, State indicators for monitoring the
 anaerobic digestion process. Water Res.44 (2010), 5973–5980.

- 331 332
- 333 imbalance in anaerobic digestors. Appl. Microbiol. Biotechnol. 43 (1995), 559-565. 334 [7] O. Lahav, B.E. Morgan, R.E. Loewenthal, Rapid, simple, and accurate method for 335 measurement of VFA and carbonate alkalinity in anaerobic reactors. Environ. Sci. Technol. 336 337 36 (2002), 2736-2741. 338 339 [8] A.J. Ward, E. Bruni, M.K. Lykkegaard, A. Feilberg, A.P.S. Adamsen, A.P. Jensen, A.K. 340 Poulsen, Real time monitoring of a biogas digester with gas chromatography, near-infrared 341 spectroscopy, and membrane-inlet mass spectrometry. Bioresour. Technol. 102 (2011), 4098-4103. 342 343 [9] B. Chatterjee, L. Radhakrishnan, D. Mazumder, New Approach for Determination of 344 345 Volatile Fatty Acid in Anaerobic Digester Sample. Environ. Eng. Sci. 35 (2018), 1-19. 346 [10] H.B. Nielsen, H.Uellendahl, B.K.Ahring, Regulation and optimization of the biogas 347 process: Propionate as a key parameter. Biomass Bioenergy 31 (2007), 820-830. 348 349 [11] Y. Wang, Y. Zhang, J. Wang, L. Meng, Effects of volatile fatty acid concentrations on 350 351 methane yield and methanogenic bacteria. Biomass Bioenergy 33 (2009), 848-853. 352 353 [12] Y. Xiong, M. Harb, P-Y Hong, Performance and microbial community variations of 354 anaerobic digesters under increasing tetracycline concentrations. Appl. Microbiol. Biot. 101 355 (2017), 5505-5517. 356 [13] S. Munir, N. Jamil, Polyhydroxyalkanoates (PHA) production in bacterial co-culture using 357 glucose and volatile fatty acids as carbon source. J. Basic Microb. 58 (2018) 247-254. 358 359 360 [14] R. Fernandez, R.M. Dinsdale, A.J. Guwy, G.C. Premier, Critical analysis of methods for the measurement of volatile fatty acids. Crit. Rev. Env. Sci. Tec. 46 (2016), 209-234. 361 362 [15] N.T. Mkhize, T.A.M. Msagati, B.B. Mamba, M. Momba, Determination of volatile fatty 363 364 acids in wastewater by solvent extraction and gas chromatography. Phys. Chem. Earth 67 (2014), 86-92. 365 366

[6] B.K. Ahring, M. Sandberg, I. Angelidaki, Volatile fatty acids as indicators of process

- [16] J. Schiffels, M.E.M. Baumann, T. Selmer, Facile analysis of short-chain fatty acids as 4 nitrophenyl esters in complex anaerobic fermentation samples by high performance liquid
 chromatography. J. Chrom.A 1218 (2011), 5848–5851.
- 370

[17] H. Turkia, S. Holmström, T. Paasikallio, H. Sirén, M. Penttilä, J-P. Pitkänen, Online
capillary electrophoresis for monitoring carboxylic acid production by yeast during bioreactor
cultivations. Anal. Chem., 85 (2013), 9705-9712.

- 374
- [18] M. Dolowy, K. Kulpinska-Kucia, A. Pyka, Validation of a thin-layer chromatography for
 the determination of hydrocortisone acetate and lidocaine in a pharmaceutical preparation.
 Sci. World J., (2014) 1-10.
- 378

[19] K. Dynska-Kukulska, W. Ciesielski, Methods of extraction and thin-layer chromatography
determination of phospholipids in biological samples, Rev. Anal. Chem. 31 (2012), 43-56.

381

[20] P.K. Zarzycki, M.M. Ślączka, E. Włodarczyk, M.J. Baran, Micro-TLC approach for fast
screening of environmental samples derived from surface and sewage waters,
Chromatographia, 76 (2013), 1249-1259.

385

[21] M. Ajmal, A. Mohammad, N. Fatima, A.H. Khan, Determination of microquantities of
mercury(II) with preliminary thin-layer chromatographic separation from mercury(I), lead(II),
nickel(II), and copper on acid-treated silica gel layers: recovery of mercury(II) from river
waters and industrial wastewaters. Microchem. J. 39 (1989), 361-371.

390

[22] F. Robert-Peillard, E. Palacio-Barco, Y. Dudal, B. Coulomb, J.L. Boudenne, Alternative
spectrofluorimetric determination of short-chain volatile fatty acids in aqueous samples, Anal.
Chem. 81(2009), 3063–3070.

394

[23] F. Robert-Peillard, E. Palacio-Barco, B. Coulomb, J.L. Boudenne, Development of a
fluorescence-based microplate method for the determination of volatile fatty acids in
anaerobically digested and sewage sludges, Talanta 88 (2012), 230–236.

398

399 [24] F. Robert-Peillard, J.L. Boudenne, B. Coulomb, Individual volatile fatty acids
determination by chromogenic derivatization coupled to multi-syringe chromatography,
Talanta 115 (2013), 737-743.

403 [25] U.S. Environmental Protection Agency. (2016) Definition and procedure for the
404 determination of the method detection limit, Revision 2. U.S. Environmental Protection
405 Agency.

406 EPA 821-R-16-006.

407

408 [26] E. Kerr, C. West, S.K. Hartwell, Quantitative TLC-image analysis of urinary creatinine
409 using iodine staining and RGB values, J. Chrom. Sci. 54 (2014), 639-646.

410

[27] A. Pyka, K. Bober, Visualizing agents for short-chain fatty acids in TLC, J. Plan. Chrom.
18 (2005), 145-150.

413

414 [28] I.M. Seligman, F.A. Doy, Thin-layer chromatography of N,N-dimethyl-p415 aminobenzolazophenacyl esters of volatile fatty acids and hydroxy acids, Anal. Biochem. 46
416 (1972), 62-66.

417

[29] D.Plachá, H. Raclavská, M. Kucerová, J. Kucharová, Volatile fatty acid evolution in
biomass mixture composts prepared in open and closed bioreactors, Waste Manage. 33
(2013), 1104-1112.

421

422

423

424 **Figure captions**

Figure 1. Visualization of spots of VFA derivatives with PMA at 3% in EtOH (A) and PMA at
3% in water (B). All VFA at 100 mg.L⁻¹, mobile phase: ethyl acetate/cyclohexane: 60:40, v/v.
Pictures taken out of the box in preliminary experiments.

Figure 2. Influence of mobile phase composition (ethyl acetate/cyclohexane, volume composition) on the separation of VFA derivatives (acetic, propionic, butyric and valeric acid at 100 mg.L⁻¹). A = Acetic acid; P = Propionic acid; B = Butyric (+ Isobutyric) acid; V = Valeric (+ Isovaleric) acid.

432 Figure 3. Visualization of potential interfering compounds after application of the analytical

433 protocol. C = Citric acid (250 mg.L⁻¹); F = Formic acid (250 mg.L⁻¹); G = Glycine (250 mg.L⁻¹);

434 VFA = mixture of all VFA at 100 mg.L⁻¹; L = Lactic acid (250 mg.L⁻¹); S = Succinic acid (250 mg.L^{-1}).

Figure 4. TLC picture of a calibration range (10-125 mg.L⁻¹) with associated densitometric histogram (left part); TLC picture of real samples with associated densitometric histogram (right part). Scales for the two histograms are different (optimized by the software on the highest peak).

441 Tables

Table 1. Resolution between VFA peaks on densitometric histograms (with each VFA at

Mobile phase composition (Ethyl acetate/ Cyclohexane, v/v)	50:50	60:40	70:30	80:20
R _{S(V/B)} ^a	0.88	1.03	0.88	0.81
$R_{S(B/P)}^{b}$	1.21	1.32	1.15	1.23
R _{S(P/A)} ^c	1.78	2.18	2.02	2.22

100mg.L⁻¹), depending on the mobile phase composition.

444 a. Resolution between valeric and butyric acid peak.

b. Resolution between butyric and propionic acid peak.

446 c. Resolution between propionic and acetic acid peak.



448

449 Table 2. Analytical features of the proposed method

Parameter	Acetic acid	Propionic acid	Butyric acid	Valeric acid	
Rf	0.24	0.43 0.55		0.64	
Calibration range (mg.L ⁻¹)	7-125	6-125	5-125	8-125	
Polynomial regression equation	y = -1.72x ² + 638.5x + 523.7	y = -2.38x ² + 756.8x - 1812.1	y = -2.26x ² + 862.8x - 1063.8	y = -1.47x ² + 721.7x + 167.8	
Correlation coefficient	0.999	0.998	0.998	0.998	
Detection limit (mg.L ⁻¹)	2.1	2.0	1.5	2.5	
Quantification limit (mg.L ⁻¹)	6.7	6.3	4.7	7.9	
RSD (%) ^a	4.6	5.3	6.7	6.5	

450 a. Repetition of a 50 mg.L⁻¹ standard (n=7) on 7 different TLC plates.

451

452 Table 3. Operational costs of the analytical protocol described in section 2.3.

Reagent/material	Quantity	Price	Cost/analysis (€)
HOAT	1.8 mg	100€/25g	0.0072
EDAN	0.6 mg	34€/25g	0.000816
EDC	3.5 mg	50€/25g	0,007

		Total/sample	0.125
		Total for 3 samples with duplicate analysis (6 spots) for one TLC plate	0.375
TLC plates (6x6 cm)	1/9 of a full plate (20x20 cm)ª	35€/25 full plates	0,156
Phosphomolybdic acid	3 g ^b	60€/100g	0.009
Ethyl acetate	0.6 mL ^ª	8.5€/L	0.0051
Cyclohexane	0.4 mL ^a	17€/L	0.00688
KH ₂ PO ₄	2.5 mg	40€/500g	0.0002

453 a. For one TLC plate (6 spots = 3 samples with duplicate analysis).

454 b. For 100 mL of PMA solution (usable for 200 TLC plate visualization).

455

456 Table 4. Comparative results between TLC and gas chromatography (GC) analysis on

457 various compost extract samples (mg.L⁻¹, mean \pm standard deviation, n=2 replicates).

	C2		C3 C4		C5			
	GC	TLC	GC	TLC	GC	TLC	GC	TLC
AN1	68.7±2.5	63.8±3.2	9.5±1.9	7.1±2.1	16.2±1.0	14.8±4.1	15.3±0.2	12.0±1.2
AN2	80.2±1.7	88.7±10.1	11.0±0.6	8.3±1.3	19.0±1.4	17.2±0.9	18.6±0.1	14.4±0.8
AN3	131.1±5.4	121.2±5.6	14.7±2.2	10.0±2.4	23.3±0.7	18.3±3.8	14.9±0.6	12.9±1.4
AN4	41.6±4.7	36.4±2.3	12.2±0.5	11.4±1.3	13.2±1.2	9.5±2.0	3.7±0.2	<lq< th=""></lq<>
AN5	36.9±0.3	29.8±6.8	5.4±0.1	<lq< th=""><th>7.9±0.3</th><th>5.3±1.0</th><th>5.2±0.5</th><th><lq< th=""></lq<></th></lq<>	7.9±0.3	5.3±1.0	5.2±0.5	<lq< th=""></lq<>

458



469 Figure 1



- 482 Figure 2



498 Figure 3

