

Enhanced Raman Scattering of amorphous matrices for Fiber Optics Sensor

M Amraoui, M Fombonne, Marie Adier, A-M Jurdyc, F. Bessueille, Bernard Champagnon, J. Margueritat, L. Bois, Bernard Dussardier, Dominique Vouagner

▶ To cite this version:

M Amraoui, M Fombonne, Marie Adier, A-M Jurdyc, F. Bessueille, et al.. Enhanced Raman Scattering of amorphous matrices for Fiber Optics Sensor. 3rd International Workshop on Metallic Nano-Objects (MNO 2016), Nov 2016, Lyon, France. . hal-02332400

HAL Id: hal-02332400 https://hal.science/hal-02332400

Submitted on 24 Oct 2019

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.





Enhanced Raman Scattering of amorphous matrices for Fiber Optics Sensor

M.Amraoui¹, M.Fombonne¹, M.Adier¹, A-M.Jurdyc¹, F.Bessueille², B.Champagnon¹, J.Margueritat¹, L.Bois³, B,Dussardier⁴, D.Vouagner¹

¹ Institut Lumière Matière, UMR5306 Université Claude Bernard Lyon 1-CNRS, Université de Lyon, 69622 Villeurbanne Cedex, France ² Institut des Sciences Analytiques, UMR5280 UCB Lyon 1-CNRS, Université de Lyon, 69100 Villeurbanne Cedex, France ³ Laboratoire des Multimatériaux et Interfaces, UMR5615, UCB Lyon 1-CNRS, Université de Lyon, 69622 Villeurbanne Cedex, France ⁴ LPMC, CNRS, UMR 7336, Université Nice Sophia Antipolis, Parc Valrose 06108 NICE CEDEX 2, France

Introduction and objectives

- Most of Surface Enhanced Raman Scattering Sensors (SERS) is applied to detection of molecules in solution
- Development of a Fiber Optic Sensor (FOS) allowing the detection of liquid and solid analytes To Raman

Surface-Enhanced Raman Scattering (SERS)

- SERS at the interface between Raman probed materials and a noble metal with a negative real part dielectric constant $\varepsilon(\omega) < 0$ (Au, Ag...)
- Two mechanisms are responsible for SERS

 $I_{Raman} \propto \alpha^2 E_{loc}^2$



Application to sol gel amorphous matrices (TiO₂, SiO₂) covering a nanostructured gold layer deposited on glass substrate and pieces of stripped optical fibers (model samples)

Model samples elaboration

<u>Step 1</u>: Gold sputtering



Step 2 : Sol elaboration

1. Electromagnetic effect :

Surface plasmon resonance increases local electric field E_{IIOC} intensity experienced by Raman probed materials:

Localized surface plasmon resonance (LSPR)

- Propagative surface plasmon (SPP)

2. Chemical effect :

Charge transfer between metal and molecules inducing adsorbed an increase of dielectric polarizability α .





0_{3,0}

·Ē 2,5 -

û 2,0 -

In the case of solid analytes: surface (SERS) or volume (VERS*) effect? * Volume Enhanced Raman Scattering?

Experimental results

10 mA 20nm LP

80 mA 25nm HP ref1

 λ_{exc} = 780 nm, P_{inc} = 14 mW, $t_{integration}$ = 1 min, ObjX20 (NA=0,40)



10

LP

Discussion and conclusion

Raman bands exaltation of amorphous matrices (TiO₂ and SiO₂) obtained after metallization (nanostructured gold layer) of model samples (glass slide and pieces of stripped optical fibers). No Raman signal is obtained without this underlying gold layer.

- Amorphous SiO₂ identified Raman bands :

0,05

 D_0 band at 490 cm⁻¹ => O_3 Si-OH symmetric stretch

- D_2 band at 617 cm⁻¹ => ring breathing of three-membered siloxane
- Band at 793 cm⁻¹ => bond-bending vibration of Si-O-Si
- Band at 992 cm⁻¹ => symmetric stretch vibrations of silanol \equiv Si-OH
- Strong band at 264 cm⁻¹ => unidentified, might be due to phase transformations.
- Amorphous TiO₂ identified Raman bands :

Band at 639 cm⁻¹ => Ti-O stretching vibration

- Lowering the pressure during the sputtering allow to get two different nanostructures (gold nanostructures size and roughness) but both Enhanced Raman spectra are almost identical: the number of interstices between gold grains, where hot spots take place, and where sol can penetrate during the dip coating process, could be an important factor for this enhancement effect. - The next step will be to test Raman Enhancement by guided optics allowing a coupling between the optical mode of the fiber and the gold plasmon wave.

• Raman signal exaltation of amorphous matrices (TiO₂ and SiO₂) obtained after metallization of the stripped fiber pieces