

Off-line coupling of capillary isotachophoresis separation to IRMPD spectroscopy for glycosaminoglycans analysis: Application to the chondroitin sulfate disaccharides model solutes

Léna Jaravel, Baptiste Schindler, Jérôme Randon, Isabelle Compagnon, Claire

Demesmay, Vincent Dugas

▶ To cite this version:

Léna Jaravel, Baptiste Schindler, Jérôme Randon, Isabelle Compagnon, Claire Demesmay, et al.. Offline coupling of capillary isotachophoresis separation to IRMPD spectroscopy for glycosaminoglycans analysis: Application to the chondroitin sulfate disaccharides model solutes. Journal of Chromatography A, 2020, 1617, pp.460782. 10.1016/j.chroma.2019.460782 . hal-02448330

HAL Id: hal-02448330 https://hal.science/hal-02448330v1

Submitted on 20 May 2022

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Distributed under a Creative Commons Attribution - NonCommercial 4.0 International License

Version of Record: https://www.sciencedirect.com/science/article/pii/S0021967319312300 Manuscript_cc015939ddb125dab2d600f65483fb4c

1	Off-line coupling of capillary isotachophoresis separation to IRMPD spectroscopy for
2	glycosaminoglycans analysis: application to the chondroitin sulfate disaccharides model solutes
3	
4	Léna Jaravel ¹ , Baptiste Schindler ² , Jérôme Randon ¹ , Isabelle Compagnon ² , Claire Demesmay ¹ , Vincent
5	Dugas ^{1*}
6	
7	1. Université de Lyon, CNRS, Université Claude Bernard Lyon 1, Institut des Sciences Analytiques, UMR 5280, 5
8	rue de la Doua, F-69100 VILLEURBANNE, France
9	2. Univ Lyon, Université Claude Bernard Lyon 1, CNRS, Institut Lumière Matière, F-69622, LYON, France
10	
11	* Corresponding author. Tel. : +33437423552
12	E-mail address: Vincent.dugas@univ-lyon1.fr
12	
15	
14	
15	Keywords: Isotachophoresis, glycosaminoglycans, capillary electrophoresis, off-line coupling,
16	microscale separation
17	
18	

19 Abstract

20 Glycans analysis is challenging due to their immense structural diversity. Isotachophoresis was 21 investigated as separation method for the purification of isobaric sulfated disaccharides prior to their 22 characterization by Mass Spectrometry (MS) and tunable IR multiple photon dissociation (IRMPD). 23 This proof of feasibility study was applied to the separation and characterization of chondroitin 24 sulfate (CS) disaccharides. ITP separation conditions were optimized. Separation starts using a 2.5 25 mM chloride ions and 10 mM glycine at pH 3.2 solution as leading electrolyte and a terminating 26 electrolyte composed of formic acid 2.5 mM and glycine 10 mM at pH 3.5. The CS disaccharides 27 sample were prepared in the terminating electrolyte. The length of injection was also investigated in 28 order to create longer plateau-like bands of pure solutes. This strategy was helpful for collecting 29 fraction at such microseparation scale. Indeed, capillary ITP affords the injection of few tens of 30 nanoliter of sample. Fractionation of the CS disaccharides mixture in isolated ITP bands and 31 collection of solutes were successfully done using a HPC coated fused silica capillary of 1m-length 32 and 75 μ m of internal diameter. Collected fractions in a final of volume 10 μ L were analyzed by CZE, 33 tandem MS and IRMPD spectroscopy. The purity of each fraction is higher than 90% and is well-34 adapted to IRMPD characterization.

35

36 1. Introduction

37 A general trend in analytical chemistry concerns the development of original approaches opening 38 new perspectives in the context of complex molecules analysis and/or characterization. Such 39 analyses are driven, for example, by the pressing demand of biologists investigating biomolecules 40 and bio-molecular interactions as well as pharmacists qualifying new biologic drugs that enter the 41 pharma market, to mention just a few. Dealing with complex samples, oligosaccharides analysis (also 42 referred to as carbohydrates or glycans) are very exciting due to their inherent structural complexity 43 and importance in biological and physiological processes through the regulation of their various 44 protein partners. [1,2] Glycosaminoglycans (GAGs) are of particular importance in this context. Their 45 presence on the cell membranes and in the extracellular matrix plays a critical role in the intercellular communication in organisms [3,4]. The GAG family, composed of 5 members: hyaluronan (HA), 46 47 heparan sulfate (HS; including the highly sulfated form heparin), chondroitin sulfate (CS), dermatan 48 sulfate (DS) and keratin sulfate (KS), are linear, highly charged, polysaccharides, composed of 49 repeating disaccharides units. Understanding the structure of GAGs via building blocks composition, blocks connectivity, stereochemical configuration and sulfation patterns, represents a major 50 51 analytical challenge.[5] Standard approaches for GAGs characterization, rely on depolymerization

52 reactions (chemical or enzymatic reactions) followed by fractionation for their further structural 53 identification by sensitive mass spectrometry and nuclear magnetic resonance experiments [6,7]. 54 Such bottom-up approaches benefit from optimal fractionation methods. Chromatographic methods 55 such as size-exclusion chromatography (SEC) and anion exchange chromatography (AEX) as well as 56 capillary electrophoresis are widely used and even coupled to one another for the separation of 57 heterogeneous mixtures before their characterization. Each have their pros and cons and there is still 58 great way of improvement in the sample preparation [6]. Underivatized oligosaccharides are 59 hydrophilic species that can be separated in ion exchange (IEX), ion-pair chromatography (IPC) and 60 HILIC modes. With respect to mass spectrometric detection, HILIC is usually preferred as this 61 chromatographic mode, compared to IEX and IPC, does not suffer from ionization suppression effects due to mobile phase components. But this separation mode has not always been effective in 62 63 resolving isomeric oligosaccharides [6,8]. With regard to sample volume, standard chromatographic autosamplers require often few tens to hundreds microliters of sample volume even when low 64 65 volume of sample (i.e. 1 µL) is injected. Capillary electrophoresis is an alternative and powerful 66 separation technique for and ionized compounds. Its separation mechanism depends on the ratio of 67 the ionic charge to the hydrodynamic radius that determines the effective mobility of ionic species in a given background electrolyte. The high peak efficiencies in CE lead to very narrow peaks (a few 68 69 seconds) requiring high sampling rate detection systems. Isotachophoresis (ITP) is a powerful but 70 often overlooked electrokinetic separation technique because its unique plateau-like signal signature 71 may be disturbing for chromatographers accustomed to work with peak signals. This technique is 72 amenable to work with small volume of crude samples [9] (no additional sample preparation needed 73 before analysis) and produces pure solute bands (without impurities like salts or organic modifiers 74 that impede subsequent MS characterization). Capillary ITP was successfully applied to the 75 preconcentration and microseparation of heparin and heparan sulfate samples on-line coupled to NMR detection via a 1-mm microcoil probes wrapped around the fused silica capillary [10,11]. The 76 77 plateau-like signals of ITP was well-adapted to such specific microcoil NMR probes and the relatively 78 long acquisition time required for NMR (1.4 s). IRMPD spectroscopy benefits from the sensitivity of 79 Mass Spectrometry and the structural resolution of IR spectroscopy: it is a powerful tool to deduce 80 detailed structural information on the monosaccharide content and stereochemistry of the glycosidic linkages [10–12]. However, the acquisition time required to obtain an IR spectrum using this 81 82 technique is several minutes, which is very long compared to standard in-line detectors of separation 83 methods. Therefore, the online coupling of IRMPD with peak signal analytical separation methods is a tremendous challenge. [16] 84

85 In this context, we propose capillary ITP as microseparation technique prior to IRMPD 86 characterization. Capillary ITP was investigated due to its plateau-like signals that allow more time 87 for relatively long acquisition time for the detection system or for the collection of pure bands of 88 solutes for off-line coupling with sophisticated characterization methods at the microscale. The proof 89 of feasibility was demonstrated with the separation of two chondroitin sulfate disaccharide isomers, namely CSA ($\Delta U(1-3)$ GalNAc,4S) and CSC ($\Delta U(1-3)$ GalNAc,6S). Chondroitin Sulfates like other GAGs 90 91 display substantial diversity within tissues and molecular diversity across different tissues and over 92 developmental time [3] and are involved in controlling neurogenesis, neuronal plasticity and 93 regeneration after injury to the central nervous system for example. In addition, it was demonstrated 94 that the specificity of interaction of CS with proteins depends on the level and distribution of 95 sulfation [4]. The ITP conditions were optimized for the separation of these two CS disaccharides 96 isomers. Collected fractions were then characterized by CZE to determine the purity yield. Finally, the 97 collected fractions were analyzed by IRMPD by nanospray infusion, demonstrating the effectiveness 98 of such coupling at the microscale.

99

- 100 **2.** Materials and methods
- 101

102 2.1. Reagents and chemicals

103 Methanol (HPLC grade), triethylamine (TEA), dipotassium hydrogen phosphate (K₂HPO₄), O-104 phosphoric acid, sodium hydroxide (NaOH), acetic acid, hydroxypropylcellulose (HPC, Mw ~100,000), 105 hydrochloric acid, glycine, formic acid were purchased from Sigma Aldrich (France). Chondroitin 106 sulfate disaccharide standards $\Delta U(1-3)$ GalNAc,4S (CSA) and $\Delta U(1-3)$ GalNAc,6S (CSC) were purchased 107 from Iduron (ALderley Edge, UK). Stock solutions of each chondroitin sulfate disaccharides (CSA and 108 CSC) were prepared in ultra-pure water at 100 mM. Acetic acid (AcOH) was purchased from Riedel-de 109 Haënn (Riedel-de Haënn, Seelze, Germany). All aqueous solutions were prepared using >18 MΩ-cm 110 water (Millipore, Molsheim, France). Polyimide (TSP)- coated fused-silica capillaries were purchased 111 from Polymicro Technologies (Polymicro Technologies, Arizona, USA).

112 2.2. HPC Coating protocols

Capillary liquid manipulations were controlled by pressure using the microfluidic flow controller
 MFCS-100 system (Fluigent, France). This system is useful to fill, wash and dry capillaries.

The HPC coating procedure follows the protocol proposed by Koval et al [17]. Briefly, a 5 % (w/v) HPC aqueous solution was flushed through 1 m of polyimide coated fused silica capillary (75 μ m ID and 375 OD μ m) at 5 bars pressure. The capillary was placed in a gas chromatography oven (HP 5890 series II GC system) heated at 60°C, and flushed with nitrogen at 1 bar pressure, for 10 minutes. The temperature was increased up to 140°C with a 5°C/min ramp and kept at 140°C for 20 min.

120 2.3. *ITP experiments*

121 ITP experiments were carried out on an Agilent HP3DCE instrument (Agilent Technologies,
122 Waldbronn, Germany) equipped with a diode array detector (detection at 230 nm).

123 Separations were done in a 35-cm long HPC coated fused silica capillary (75-μm id / 365 μm od).

Otherwise stated in the text, ITP experiments were performed using a pH 3.2 leading electrolyte composed of HCl 2.5 mM, glycine 10mM and a pH 3.5 terminating electrolyte composed of formic acid 2.5 mM and glycine 10 mM. Mixtures of CSA and CSC at 1 mM each were prepared in the terminal electrolyte. Injections were done under hydrodynamic pressure mode.

128 2.4. *CZE experiments*

Capillary zone electrophoresis experiments were carried out on an Agilent 7100 CE instrument (Agilent Technologies, Waldbronn, Germany) equipped with a diode array detector. Separations were done in a 35-cm length HPC coated fused silica capillary (75-μm id / 365 μm od). Prior to each injection, the capillary was equilibrated by flushing with the background electrolyte (formic acid 50 mM and glycine 100 mM pH 3.24) for 30 sec. Samples were injected under hydrodynamic mode by applying a 5 mbar pressure for 1s. Separation was done at - 20 kV.

135 2.5 Mass Spectrometry and IRMPD spectroscopy

136 Off-line coupling experiments were realized by injecting 22 nl of the 1 mM chondroitin sulfate 137 disaccharides mixture (corresponding to 44 pmol of each chondroitin sulfate disaccharides). Bands collection of CSA and CSC after ITP separation were done in 10 μ L H₂O and diluted in 10 μ L MeOH. MS 138 139 analysis was carried out in a commercial mass spectrometer (Thermofinnigan LTQ XL) equipped with 140 a nanospray ionization source set in negative mode (production of deprotonated ions) and a linear 141 ion trap analyzer. Tandem MS experiments are performed via collision induced dissociation (CID) in a 142 helium collision gas. A similar setup used for IRMPD analysis was previously described by Schindler et 143 al. [18]. In brief, the mass spectrometer is modified to allow irradiation of the trapped ions by a laser 144 beam produced by a YAG-pumped tunable IR OPO/OPA. This enables in situ gas-phase IRMPD (InfraRed Multiple Photon Dissociation) spectroscopy of the trapped ions. Excitation of a vibrational 145

146 mode by the laser induces photofragmentation, which is detected by mass spectrometry. The photo-147 fragmentation yield is calculated using the formula: IRMPD yield = $-\log[I_{parent} \div (I_{parent} + I_{fragment})]$, 148 where I_{parent} is the intensity of the parent ion and I_{fragment} is the total intensity of the fragments 149 detected after laser excitation. The fragmentation yield is then monitored as a function of the 150 wavenumber to retrieve the absorption spectrum of the analyte. The spectral range 2700-3700 cm⁻¹ 151 was scanned at 1 cm⁻¹/s, which correspond to an acquisition time of 17 minutes. In these conditions 152 mass-selected ions were isolated and irradiated for 500 ms per photofragmentation MS spectrum and MS spectra were not averaged. The reference IRMPD spectra of pure CSC and CSA standards 153 154 were measured on a different setup, as described in [13].

155

156 **3. Results**

157 3.1. Determination of the optimal ITP separation conditions

158 - Choice of the electrolytes (pH and nature)

As for other differential migration methods, the condition for ITP separation is that the ionizedsolutes in the sample will separate whenever their migration rates in the mixed state are different.

161 Figure 1

The two chondroitin sulfate disaccharides will better separate if the pH in the mixing zone increases 162 163 the difference in their respective effective mobility. The electrophoretic mobility of the two 164 chondroitin sulfate disaccharides was thus first investigated in Capillary Zone Electrophoresis in acidic 165 conditions due to the weak acidic properties of the CS disaccharides (carboxylic acid and sulfate 166 groups)[19]. Figure 1 presents the relative electrophoretic mobility difference between the two CS 167 disaccharides at a pH ranging from 2.18 to 4.06. The largest relative electrophoretic mobility difference was observed between pH 3.21 to 3.45. This pH range was thus considered optimal for the 168 169 ITP separation and was fixed for the preparation of the leading electrolyte and terminal electrolyte.

A high effective mobility constituent such as chloride (strong ionic species commonly used in ITP) must be used as leader. The pH of the leading electrolyte was adjusted using glycine to pH 3.2. To cope with the objective to form long "plateau-like" bands of pure solute, low chloride leading electrolyte concentrations were investigated (ranging from 1 to 10 mM). The lowest concentration that can be used is however limited by the minimal current that can be carried by the electrolyte and controlled accurately (down to 1 μ A) by the CE instrument. The best compromise for chloride concentration versus applied current was 2.5 mM. 177 Considering the mobility of the chondroitin sulfate disaccharides in this pH range, formate anions 178 were used as constituent of the terminating electrolyte. Formic acid (pKa 3.75, formate 179 electrophoretic mobility 5.66 10^{-4} cm² V⁻¹ s⁻¹)[20] has an effective mobility below pH 3.5 lower than 180 2.18 10^{-4} cm² V⁻¹ s⁻¹, which is a lower mobility value than the mobility of CS disaccharides.

181 The terminating electrolyte was thus prepared by adjusting the pH of the formic acid solution with182 glycine as counter buffering and counter constituent.

183

184 - Sample zone optimization

ITP separation of CS disaccharides injected individually. In first approach, each chondroitin was 185 186 injected separately (Figure 2A, blue and orange plot) after conditioning the capillary with the leading 187 electrolyte (HCl 2.5 mM, Glycine (10 mM) pH 3.5). ITP separation was initiated by applying a constant 188 current between the inlet side of the capillary immersed in the terminating electrolyte and its outlet 189 sideimmersed in the leading electrolyte. While contactless conductivity was also tested, the 190 difference in conductivity of the two CS disaccharide zones is quite low. Hence it is preferred to 191 monitor the separation using UV detection. A plateau was detected at the optimal wavelength of 230 192 nm for the two CS disaccharides. The two CS disaccharides present a relatively different absorbance 193 that is useful to attribute the signal obtained to the chondroitin passing through the detector. The 194 two CS disaccharides were then injected successively (CSC was injected in first position). Two 195 plateaus were obtained and the UV signal for each chondroitin (Fig2A, purple plot) was comparable to the signal observed for the individual ITP separated CS disaccharides. UV detection is useful for 196 197 monitoring electrokinetic separation using commercial instrument as the UV detection is done in-198 capillary, the two CS disaccharides are detected at 230 nm and can be differentiated.

199 Figure 2

200 ITP separation of CS disaccharides in mixture. A mix of the two CS disaccharides solubilized in water 201 was injected according to the previous experimental conditions. Only one plateau was observed (Figure 2B, grey plot). The two CS disaccharides were not separated in these experimental conditions. 202 203 The pH of the sample zone is also an important experimental parameter for the separation of ionic 204 species in ITP [21]. In order to trigger the separation at the early stage of the ITP separation, the pH 205 of the sample zone was adjusted to 3.2-3.5. The CS were diluted in the terminating electrolyte to 206 adjust the pH. The buffering of the sample zone separation allowed to successfully separate the CS 207 disaccharides (Fig 2B, dark plot). The ITP separation leads to the formation of two plateaus of 208 increasing UV absorbance signals. The more mobile CSC solute is detected first and CSA second.

209 3.2. Fractionation/collection of the ITP separated CS disaccharides

210 - Optimization of the injection volume

For in-line or off-line coupling, preparation of long band lengths of pure solute would be obviously highly valuable. This is ever true at the miniaturized scale where collecting fraction is not easy. In order to maximize the band width, a 1-m long HPC coated capillary (75 µm ID) was used. Indeed, the length of capillary needed to separate efficiently the two CS disaccharides increases with the quantity injected. A 1-m column is the longest capillary compatible with our instrumental set-up. This is due to the low leading electrolyte concentration used. Indeed, it is difficult to maintain constant current during the separation using longer capillary.

218 Figure 3

219 The sample volume was increased linearly by adjusting the injection pressure during the injection 220 step (Figure 3). Up to 30 mbar, the separation led to the formation of two plateaus. At 40 mbar, two 221 plateaus were still detected but the boundary between the two plateaus seemed less pronounced 222 and at higher injection pressure conditions (50 mbar) 3 plateaus were observed. This intermediate 223 plateau corresponds to a mixing zone of CS A and CS S between the two pure bands in formation. For 224 the following experiments, the injection conditions were set to 30 mbar for 10 sec. The linear velocity 225 of mobile phase at 30 mbar was estimated to be 0.05 cm/sec using leading electrolyte. This 226 corresponds to an injection plug length of 5 mm. After separation, the band lengths (plateau) were 227 about 12 mm (this was estimated by measuring the duration of plateaus under pressure after their 228 formation by ITP).

In the context of miniaturized separation methods, the optimization of the ITP separation conditions led to increase the injection volume up to about 22 nL and the formation of individual bands of solutes of about 12 mm highlighting the potential interest of such methods for sample preparation at the miniaturized scale.

233

- Fractionation of the two chondroitins and analysis of each fraction by CZE

235 Figure 4

236 Under the optimized separation conditions, the two CS disaccharides were well-separated and 237 formed two plateaus of more than 1 cm length. Fractionation of the sample after separation was 238 conducted by collection of a fraction of each plateau in 10 μ L of water. From the experimental point 239 of view, the ITP separation was stopped before complete elution of the CS disaccharides from the 240 capillary. The outlet vial was then replaced by a vial containing a low volume conical insert filled with 241 10 µL of water. The fraction of pure solute selected was eluted under hydrodynamic flow mode (50 242 mbar). The two fractions were recovered using two different ITP separations. The hydrodynamic flow 243 is convenient to accurately fractionate a part of the band, however it was supposed that during this 244 step, all interfaces (CSC/CSA, CSA or CSC/electrolyte) undergo band diffusion. This is expected to mix 245 the two separated chondroitins and thus affect their purity. Purity of each collected fraction was 246 assessed by capillary zone electrophoresis (Figure 4 B). Regarding the electrophoregram obtained 247 using a chondroitins standard solution, two pairs of chondroitin peaks were detected. Each pair was 248 attributed to chondroitin sulfate anomer forms. The first peak pattern corresponds to the CSC 249 chondroitin sulfate and the second one to the CSA. Each collected fraction presents a major 250 chondroitin sulfate isomer. Fraction 1 correspond to the CSC isoform with a purity of 89.6% and 251 fraction 2 corresponds to the CSA isoform with a purity of 94.6 %. The isolation of oligosaccharide components for structural investigations must have a purity greater that 80% to 90 %. Moreover, 252 253 organic impurities or salts that can impede the characterization must also be minimized [6]. These 254 are generally observed in HPLC separations for example. ITP seems well-adapted for the separation 255 and purification of CS samples.

256

257 3.3. IRMPD analysis

258 Each collected fraction was characterized by IRMPD spectroscopy. The two disaccharides (CSA and 259 CSC) derived from chondroitin sulfate are challenging to characterize as they possess two neighboring sites of sulfation on the N-acetyl galactosamine that is O6 and O4. Infrared multiple 260 261 photon dissociation spectroscopy (IRMPD) was recently applied to the characterization of isomeric glycosaminoglycan disaccharides from heparins and chondroitin sources and has shown great 262 263 potential as a novel structural tool for carbohydrates [13]. However, this powerful technique requires 264 relatively long time of acquisition, which is not convenient for coupling with most separation 265 methods [16,22]

266 Figure 5

The MS/MS spectra of CSA and CSC are shown in the top layer of Figure 5. The singly deprotonated precursor ions are observed at m/z 458. CSA features a main fragment at m/z 300, which corresponds to the sulfated GalNAc moiety and is generally considered as its diagnostic fragment ion; and CSC features a main fragment at m/z 282, which is generally considered as its diagnostic fragment ion. [22]. [23,24]. Several authors reported that the fragmentation patterns are highly dependent on the experimental conditions and MS/MS analysis must be carried out with extreme caution. Yet, LC coupled to MS and/or tandem MS analysis is the most widespread approach for theanalysis of glycosaminoglycan composition [25].

275 The use of the IR vibrational spectrum of molecular ions(sensitive towards structural details) enables 276 the investigation of functional groups, intra- or intermolecular interactions and molecular 277 conformations bringing additional information compared to the MS/MS analysis. The IR spectrum 278 presents intrinsic features, which do not depend on the MS conditions, that benefits to the 279 saccharides analysis by recording IR fingerprints of the CS disaccharides for unambiguous 280 characterization. The IRMPD spectra of the two collected fractions show an apex at 3410 cm⁻¹ and 281 3375 cm⁻¹ attributed to the CSA and CSC disaccharides respectively, which is in good agreement with 282 reference spectra of CS disaccharides standards [13]. In the specific case of CSA and CSC, the MS/MS 283 fingerprint is sufficient to distinguish between the isomers, and the differences in the IRMPD 284 signatures are relatively subtle. However, IRMPD has proven useful to distinguish between a variety 285 of isomers when tandem MS is ambiguous. [14] Therefore, the presented example supports that the 286 fractionation of tens of nL of isobaric saccharides mixtures using ITP can be applied to the off-line 287 collection of high purity compounds for further IR spectroscopic characterization. This demonstrates the pertinence of such coupling strategy at the microscale separation. 288

289

4. Conclusion

291 Electrokinetic separation of isomeric chondroitin sulfate disaccharides (CSA and CSC) at the microscale was achieved by isotachophoresis and fractions of high purity were collected for further 292 293 structural analysis. ITP generates plateau-like bands of purified compounds which is particularly 294 advantageous for collecting sample at the microscale or for coupling with long acquisition time 295 detection methods. In contrast with ion exchange chromatography, ITP does not require the use of 296 high salt concentration that can interfere severely with Mass Spectrometry. Such proof of feasibility 297 aims at bringing an original coupling strategy with sophisticated spectroscopic characterization 298 method of oligosaccharides. The molecular characterization of such complex molecules requires 299 combined approaches of separation and characterization. Optimization of the ITP conditions allows 300 the separation of two isobaric CS disaccharides from few tens of nL sample solution and their 301 collection in 10 µL. The purity of the collected fraction, more than 90 %, is well-adapted to off-line 302 coupling with gas phase IR action spectroscopy. This case study validates this proof of feasibility 303 study using two similar CS disaccharides GAGs precursor. In order to fully respond to biological 304 demand, the in-line coupling ITP-IRMPD spectroscopy would be further investigated allowing 305 miniaturized and automated analysis of complex samples.

306

307 References

- N.S. Gandhi, R.L. Mancera, The Structure of Glycosaminoglycans and their Interactions with
 Proteins, Chem Biol Drug Des. 72 (2008) 455–482. https://doi.org/10.1111/j.1747 0285.2008.00741.x.
- A. Varki, Biological roles of oligosaccharides: all of the theories are correct, Glycobiology. 3
 (1993) 97–130. https://doi.org/10.1093/glycob/3.2.97.
- R.A. Townley, H.E. Bülow, Deciphering functional glycosaminoglycan motifs in development,
 Curr Opin Struct Bio. 50 (2018) 144–154. https://doi.org/10.1016/j.sbi.2018.03.011.
- R.P. Richter, N.S. Baranova, A.J. Day, J.C. Kwok, Glycosaminoglycans in extracellular matrix
 organisation: are concepts from soft matter physics key to understanding the formation of
 perineuronal nets?, Curr Opin Struct Bio. 50 (2018) 65–74.
 https://doi.org/10.1016/j.sbi.2017.12.002.
- 319 [5] I. Compagnon, B. Schindler, G. Renois-Predelus, R. Daniel, Lasers and ion mobility: new
 320 additions to the glycosaminoglycanomics toolkit, Current Opinion in Structural Biology. 50
 321 (2018) 171–180. https://doi.org/10.1016/j.sbi.2018.06.006.
- 322 [6] C.J. Jones, S. Beni, J.F.K. Limtiaco, D.J. Langeslay, C.K. Larive, Heparin Characterization:
 323 Challenges and Solutions, Annual Rev. Anal. Chem. 4 (2011) 439–465.
 324 https://doi.org/10.1146/annurev-anchem-061010-113911.
- 325 [7] E. Sisu, C. Flangea, A. Serb, A.D. Zamfir, Modern developments in mass spectrometry of
 326 chondroitin and dermatan sulfate glycosaminoglycans, Amino Acids. 41 (2011) 235–256.
 327 https://doi.org/10.1007/s00726-010-0682-4.
- V.L. Gill, U. Aich, S. Rao, C. Pohl, J. Zaia, Disaccharide Analysis of Glycosaminoglycans Using
 Hydrophilic Interaction Chromatography and Mass Spectrometry, Anal. Chem. 85 (2013) 1138–
 1145. https://doi.org/10.1021/ac3030448.
- S.L. Eldridge, V.K. Almeida, A.K. Korir, C.K. Larive, Separation and Analysis of Trace Degradants
 in a Pharmaceutical Formulation Using On-Line Capillary Isotachophoresis-NMR, Anal. Chem. 79
 (2007) 8446–8453. https://doi.org/10.1021/ac7016629.
- A.K. Korir, V.K. Almeida, D.S. Malkin, C.K. Larive, Separation and Analysis of Nanomole
 Quantities of Heparin Oligosaccharides Using On-Line Capillary Isotachophoresis Coupled with
 NMR Detection, Anal. Chem. 77 (2005) 5998–6003. https://doi.org/10.1021/ac050669u.
- A.K. Korir, C.K. Larive, On-line NMR detection of microgram quantities of heparin-derived
 oligosaccharides and their structure elucidation by microcoil NMR, Anal Bioanal Chem. 388
 (2007) 1707–1716. https://doi.org/10.1007/s00216-007-1400-2.
- E. Mucha, A. Stuckmann, M. Marianski, W.B. Struwe, G. Meijer, K. Pagel, In-depth structural
 analysis of glycans in the gas phase, Chem. Sci. 10 (2019) 1272–1284.
 https://doi.org/10.1039/C8SC05426F.
- 343 [13] B. Schindler, L. Barnes, C.J. Gray, S. Chambert, S.L. Flitsch, J. Oomens, R. Daniel, A.R. Allouche, I.
 344 Compagnon, IRMPD Spectroscopy Sheds New (Infrared) Light on the Sulfate Pattern of
 345 Carbohydrates, J. Phys. Chem. A. 121 (2017) 2114–2120.
- 346 https://doi.org/10.1021/acs.jpca.6b11642.
- [14] B. Schindler, L. Barnes, G. Renois, C. Gray, S. Chambert, S. Fort, S. Flitsch, C. Loison, A.-R.
 Allouche, I. Compagnon, Anomeric memory of the glycosidic bond upon fragmentation and its
 consequences for carbohydrate sequencing, Nature Com. 8 (2017) 973.
 https://doi.org/10.1038/s41467-017-01179-y.
- [15] B. Schindler, G. Renois-Predelus, N. Bagdadi, S. Melizi, L. Barnes, S. Chambert, A.-R. Allouche, I.
 Compagnon, MS/IR, a new MS-based hyphenated method for analysis of hexuronic acid
- epimers in glycosaminoglycans, Glycoconjugate J. 34 (2017) 421–425.
- 354 https://doi.org/10.1007/s10719-016-9741-8.

- B. Schindler, G. Laloy-Borgna, L. Barnes, A.-R. Allouche, E. Bouju, V. Dugas, C. Demesmay, I.
 Compagnon, Online Separation and Identification of Isomers Using Infrared Multiple Photon
 Dissociation Ion Spectroscopy Coupled to Liquid Chromatography: Application to the Analysis of
 Disaccharides Regio-Isomers and Monosaccharide Anomers, Anal. Chem. 90 (2018) 11741–
 11745. https://doi.org/10.1021/acs.analchem.8b02801.
- 360 [17] D. Koval, V. Kašička, H. Cottet, Analysis of glycated hemoglobin A1c by capillary electrophoresis
 361 and capillary isoelectric focusing, Analytical Biochemistry. 413 (2011) 8–15.
 362 https://doi.org/10.1016/j.ab.2011.01.048.
- 363 [18] B. Schindler, A.D. Depland, G. Renois-Predelus, G. Karras, B. Concina, G. Celep, J. Maurelli, V.
 364 Loriot, E. Constant, R. Bredy, C. Bordas, F. Lépine, I. Compagnon, FAIMS-MS-IR spectroscopy
 365 workflow: a multidimensional platform for the analysis of molecular isoforms, International
 366 Journal for Ion Mobility Spectrometry. 20 (2017) 119–124. https://doi.org/10.1007/s12127367 017-0225-8.
- 368 [19] B. Larsson, M. Nilsson, H. Tjälve, The binding of inorganic and organic cations and H+ to
 369 cartilage in vitro, Biochem Pharmacol. 30 (1981) 2963–2970. https://doi.org/10.1016/0006370 2952(81)90260-4.
- 371 [20] J.L. Beckers, Steady-state models in electrophoresis: From isotachophoresis to capillary zone
 372 electrophoresis, Electrophoresis. 16 (1995) 1987–1998.
 373 https://doi.org/10.1002/elps.11501601327.
- F.E.P. Mikkers, E.M. Everaerts, J.A.F. Peek, Isotachophoresis: The concepts of resolution, load
 capacity and separation efficiency I. Theory, J Chromatogr A. 168 (1979) 293–315.
 https://doi.org/10.1016/0021-9673(79)80001-1.
- 377 [22] G. Renois-Predelus, B. Schindler, I. Compagnon, Analysis of Sulfate Patterns in
 378 Glycosaminoglycan Oligosaccharides by MSn Coupled to Infrared Ion Spectroscopy: the Case of
 379 GalNAc4S and GalNAc6S, J Am Soc Mass Spectrom. 29 (2018) 1242–1249.
 380 https://doi.org/10.1007/s13361-018-1955-5.
- 381 [23] H. Desaire, J.A. Leary, Detection and quantification of the sulfated disaccharides in chondroitin
 382 sulfate by electrospray tandem mass spectrometry, J Am Soc Mass Spectrom. 11 (2000) 916–
 383 920. https://doi.org/10.1016/S1044-0305(00)00168-9.
- Y. Zhang, A.H. Conrad, E.S. Tasheva, K. An, L.M. Corpuz, Y. Kariya, K. Suzuki, G.W. Conrad,
 Detection and Quantification of Sulfated Disaccharides from Keratan Sulfate and
 Chondroitin/Dermatan Sulfate during Chick Corneal Development by ESI-MS/MS, Invest
- 387 Ophthalmol Vis Science. 46 (2005) 1604–1614. https://doi.org/10.1167/iovs.04-1453.
- 388 [25] G.O. Staples, J. Zaia, Analysis of Glycosaminoglycans Using Mass Spectrometry, Curr
 389 Proteomics. 8 (2011) 325–336. https://doi.org/10.2174/157016411798220871.

390

392 Figure Captions:

Figure 1. Plot of the relative electrophoretic mobility difference between CS A and CS C as a function
of the pH of the glycine/HCl electrolyte. Electrophoretic mobilities were determined by capillary zone
electrophoresis using the procedure described by Vigh et al. [18].

396 Figure 2. Isotachophoregrams of each individual chondroitin (CSA and CSC at 1 mM) injected separately and injected in stacked configuration (CSC + CSA) (A) and ITP separation of CS 397 398 disaccharides in mixture in water or in the terminating electrolyte (B). Conditions for (A): HPC-coated 399 capillary length: 50 cm, Injection 50 mbar for 3 sec. Applied Current = $-4 \mu A$ and in mixture (CSA and 400 CSC at 1 mM each in water (grey) and solubilized in the terminating electrolyte (black)). Conditions 401 for (B) : ITP of CS disaccharides solubilized in water was realized in a 45.9 cm-length HPC-coated 402 capillary, injection 5 mbar for 8 sec. Applied Current = $-5\mu A$. ITP of CS disaccharides solubilized in 403 terminating electrolyte was realized in a 58.7 cm length HPC-coated capillary, injection 20 mbar for 404 4sec). Applied Current = $-5\mu A$.

Figure 3. Influence of the injection pressure of CS disaccharides mixture at 1 mM solubilized in the terminating electrolyte, separation at -3 μ A, HPC-coated fused silica capillary (75 μ m id, total length 1m), UV detection at 230 nm, injection time (10 sec), corresponding injection volumes 7.8 to 39 nL.

408 Figure 4. ITP separation of chondroitin sulfate disaccharides (CSA and CSC) in a 1-m long HPC-coated 409 capillary and fractionation of the separated solutes. Injection of chondroitin sulfate disaccharides at 1 410 mM each in terminating electrolyte, separation under -3 µA, detection at 230 nm, injection 30 mbar/ 10 sec (A). CZE separation of the two collected fractions and of the mixture of the two chondroitin 411 412 sulfate disaccharides. HPC Fused silica capillary (75 µm id, 30-cm total length, 21.5 cm-length to 413 detection), UV detection at 230 nm, BGE: Formate 50 mM, glycine 10 mM, pH 3.5, injection CS 414 standard mixture (20 µM in terminating electrolyte) (-10kV, 5 sec) injection of the CS fraction (-10kV, 415 10 sec), response time 0.063 sec (40 Hz).

Figure 5. IRMPD analysis of fraction 2 (A) and fraction 1 (B). Top panels: CID spectra of singly deprotonated disaccharides (m/z 458). Center panels: IRMPD signatures of the m/z 458 ions. Lower panels: reference IRMPD spectra of pure CSA (A) and CSC (B) disaccharide standards. Strips corresponding to the maximum of absorption of CSA (light red) and CSC (light blue) standards are overlaid for facile structural identification.









