

A simple GC–MS method for the determination of diphenylamine, tolylfluanid propargite and phosalone in liver fractions

Ali Kadar, Ludovic Peyre, Henri Wortham, Pierre Doumenq

To cite this version:

Ali Kadar, Ludovic Peyre, Henri Wortham, Pierre Doumenq. A simple GC–MS method for the determination of diphenylamine, tolylfluanid propargite and phosalone in liver fractions. Journal of Chromatography B Biomedical Sciences and Applications, 2019, 1113, pp.69-76. 10.1016/j.jchromb.2019.03.005. hal-02469623

HAL Id: hal-02469623 <https://hal.science/hal-02469623>

Submitted on 6 Feb 2020

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

A simple GC-MS method for the determination of diphenylamine, tolylfluanid propargite and phosalone in liver fractions

3 Ali Kadar^{a*}, Ludovic Peyre^{b*}, Henri Wortham^a, Pierre Doumenq^a

 ^a Aix Marseille Univ, CNRS, LCE, Marseille, France, ^b Cabinet conseil Santé-Environnement, 43 Rue Antoine Brun, 06150 Cannes, France * Corresponding authors. *E-mail addresses:* ali.kadar@etu.univ-amu.fr (A. Kadar) and ludovic.peyre-teisseire@hotmail.fr (L. Peyre)

ABSTRACT

 In this study, an accurate and robust gas chromatography/mass spectrometry method was developed for quantitative analysis of diphenylamine, tolylfluanid, propargite and phosalone in liver fractions. Different injector parameters were optimized by an experimental design of experiment technique (central composite design). An optimal combination of injector temperature (°C), splitless time (min) and overpressure (KPa) values enabled to maximize the chromatographic responses. Sample preparation was based on protein precipitation using trichloroacetic acid followed by liquid-liquid extraction (LLE) of the pesticides with hexane. All compounds were quantified without interference in selected ion monitoring mode using endrin as internal standard. The calibration curves for diphenylamine, tolylfluanid, propargite and phosalone compounds were linear over the concentration range of 19 0.1 to $25 \mu M$ (*i.e.*0.1, 0.5, 2.0, 5.0, 15.0 and $25.0 \mu M$) with determination coefficients (R²) higher than 0.999. A lower limit of quantification of 0.1 µM was obtained for all analytes, *i.e.* 422.5, 868.0, 876.2 and 919.5 µg/kg of liver fraction (hepatocytes) for diphenylamine, tolylfluanid, propargite and phosalone, respectively. All compounds showed extraction recoveries higher than 93%, with a maximum RSD of 3.4%. Intra- and inter-day accuracies varied from 88.4 to 102.9% and, imprecision varied from 1.1 to 6.7%. Stability tests demonstrated that all pesticides were stable in liver extracts 25 during instrumental analysis (20° C in the autosampler tray for 72 h) following three successive freeze- thaw cycles and, at -20°C for up to 12 months. This simple and efficient analytical procedure is thus suitable for assessing mammals liver contamination or metabolism studies.

 Keywords: Bioanalytical method, Pesticides mixture, Experimental design, GC-MS, Liver fractions, Protein precipitation.

1. Introduction

 For more than half a century, agricultural practices have involved the use of a large amounts of pesticides. These pesticides provide better agricultural yields and allow extending the shelf life of perishable fruits and vegetables to fulfill the need of the growing worldwide population [1]. Unfortunately, these helpful active compounds may at the same time constitute a significant risk to animal and human health [2]. Indeed, pesticide residues are widespread in the whole environment and are likely to be present in water resources and agricultural products [3-5]. Water and food crop consumption is the predominant pathway of exposure for the general population [6-12].

 On this population scale, several authors evaluated the dietary exposure to pesticides residues and pointed out that the consumers may be simultaneously exposed to different residues [13-22]. In their work [23], Crepet et al. concluded that depending on the foodstuff consumed the French general population was mainly and most heavily exposed to 7 different pesticide mixtures composed of 2 to 6 compounds. Diphenylamine, tolylfluanid, propargite and phosalone belong to a mixture that was significantly correlated to common fruits such as apples and pears. Once these potentially contaminated fruits have been consumed and the compounds have passed into the body, these latter are transported by the blood flow to the liver for metabolization [24]. Therefore, in the scope of metabolism [24] or biomonitoring [25-27] studies, a sensitive and reliable analytical method is needed to assess the 48 pesticides level in liver fractions. Until now, as far as we know, bibliographic research shows that no work has been published on the simultaneous determination of diphenylamine, phosalone, propargite and tolylfluanid in human biological samples.

 However, Oliveira et al. and Kaczynski et al. have recently published multiresidues analytical protocols including the determination of propargite and phosalone [24] or the quantification of propargite and tolylfluanid [25] in fish liver. These methods were based on the use of liquid chromatography coupled to tandem mass spectrometry (LC-MS/MS). The study reported by Russo et al. [26] is, to the best of our knowledge, the only one dealing with real human liver samples. Authors carried out the development of a gas chromatography coupled to a negative chemical ionization mass spectrometry (GC-NCI-MS) method for the analysis of several organophosphorus pesticides including phosalone. Valente et al. [27] quantified this latter in human blood along with other organophosphates thanks to a GC-MS method. In the same matrix, Sharma et al. [28] successfully quantitated phosalone using gas chromatography coupled to a flame thermionic detector (GC-FTD).

 For sample purification, works published by Valente et al. [27] and Sharma et al. [28] depicted a one- step phosalone extraction/clean-up from blood liquid samples using either silica based C 18 reversed- phase [27] or florisil [28]. Robles-Molina et al. [x] extracted diphenylamine and phosalone from water with Oasis HLB sorbent. Lehotay et al. [y] validated automated mini-SPE with MgSO4/primary secondary amine (PSA)/C18/CarbonX sorbent for the analysis of diphenylamine and propargite in food

samples.

 On the other hand, solvent extractions proved to be suitable for the analysis of phosalone, propargite and tolyfluanid in liver samples [24-26].

 An analytical protocol using solvent extraction followed by a GC separation and MS detection seemed convenient for assessing liver contamination by the present pesticide mixture.

Thus, this work aimed at developing and validating a simple, accurate and robust analytical method for

 the analysis of diphenylamine, tolylfluanid, propargite and phosalone in human liver fractions (hepatocytes).

2. Experimental

- *2.1. Chemicals, materials and biological samples*
- Certified standards of diphenylamine, tolylfluanid, propargite, phosalone and endrin with purities higher
- than 99.5 %; ammonium sulfate salts and trichloroacetic acid (TCA) of research grade purity were
- purchased from Sigma-Aldrich (St Quentin Fallavier, France). *n*-Hexane, ethyl acetate and acetonitrile
- of ultra-trace analysis grade were purchased from Carlo Erba (Val de Reuil, France). Standard stock
- solutions were prepared from pure compounds and appropriately diluted in acetonitrile.
- A high-throughput tissue grinder (MM 300) produced by Retsch (Haan, Germany) was used as a
- powerful transversal shaker for the LLE experiments. Phase separation was achieved using a Thermo
- IEC Micromax RF benchtop centrifuge from Thermo Fisher Scientific (Illkirch, France).

 All experiments on human tissues were carried out according to the ethical standards of the responsible committee on human experimentation and the Helsinki Declaration. Cellular (hepatocytes) fractions were obtained after mechanical decomposition of liver tissues. Here we decided to carry out the 87 experiments with thermally (100°C for 3 min) inactivated hepatocytes previously isolated following 88 Berry and Friend procedure [29].

2.2 Sample treatment

 A 400 µL volume of thermally inactivated liver cells (16 mg) in 100 mM phosphate potassium buffer (pH 7.4) was added to 1.8 mL Eppendorf® tubes. Prior to a brief vortex mix of the samples, the required amounts of diphenylamine, tolylfluanid, propargite, phosalone and endrin were added as internal 93 standard (IS). Then, 100 μ L of ice-cold trichloroacetic acid (20%) and 400 μ L of Hexane were added to the tubes. The samples were vigorously shaken for 3 minutes and centrifuged at 16000g and 4°C for 4 min, enabling the precipitation of the denatured proteins and full separation of the liquid phases. The 96 supernatant was then transferred into a 200 µL GC vial insert for GC analysis.

2.3. GC-MS Analysis

Separation and quantification of the pesticides mixture were performed using a TRACETM ultra gas chromatograph coupled to a DSQ II single quadrupole mass spectrometer (Thermo Fisher Scientific, Courtaboeuf, France). The gas chromatograph was equipped with a programmable split/splitless 101 injector, a capillary column and a programmable oven. A sample volume of $2 \mu L$ was injected at $271^{\circ}C$, 102 in splitless mode, in a baffle Siltek-deactivated liner $(2 \text{ mm} \times 2.75 \text{ mm} \times 120 \text{ mm})$ provided by Thermo Fisher Scientific . A surge pressure of 490 kpa was applied for a 1.25 min period right after injection. The separation of the four compounds was achieved on a TRACE TR5-MS (Thermo Fisher Scientific) 105 capillary column (30 m \times 0.25 mm, *i.d.*, 0.25 µm film thickness). The electron impact ion source and transfer line temperatures were set at 250 and 280°C, respectively. The elution was carried out using a 1.0 ml/min constant flow rate of helium carrier gas and a temperature gradient. After optimization of the standards mixture separation, the oven temperature was programmed as follows: the initial 109 temperature set at 50 °C for 2.0 min followed by a first ramp of 35 °C .min⁻¹ up to 210 °C and held for 110 4.0 min; then a second ramp of 35 $^{\circ}$ C min⁻¹ up to 280 $^{\circ}$ C was set and maintained for 4.1 min. The total optimized run time was 17 min. The criteria for positivity used to identify the compounds were both retention times and characteristic m/z signals. Quantification was accomplished in the selected ion monitoring (SIM) mode and the monitored ions are displayed in **Table 1**. 114

115 Table 1. Ions monitored under the [S](https://www.sciencedirect.com/science/article/pii/S157002321831208X#tf0005)IM mode by GC–MS^a and their relative

116 intensities (%).

117 a The compounds were quantified with the sum of both base peak and fragment peak signals.

118 **Ionized in the positive mode with an electron energy of 70 eV.**

119

120 *2.4. Optimization of splitless injection parameters*

130 tolylfluanid, propargite and phosalone.

121 To achieve the lowest instrumental detection limits, in addition to preliminary optimization of 122 chromatographic and ionization source settings, the splitless injection was optimized. For that purpose, 123 a central composite design (CCD) with three independent variables $(X_1,$ Temperature; X_2 , Surge 124 pressure; X_3 , Splitless Time) was performed. Preliminary experiments with one factor at a time allowed 125 to determine the appropriated ranges of these independent variables: X_1 : 133–318 °C; X_2 : 0–490 kPa; 126 X3: 0–2.50 min. As shown in **Table 2A**, five levels were used for optimization of the three factors. 127 Design and analysis of the central composite experiment were carried out using Statistica 8.0 software 128 (Statsoft, Maison Alfort, France). A total of 48 assays described in **Table 2B** were carried out by 129 automated injections of 2μ L of a 10 μ M multi-compounds standard solution containing diphenylamine,

Table 2. Values of the factors at the five levels examined (A) and experiments undertaken for the central composite design (B).

2.5. Method validation

Validation was achieved according to the Food and Drug Administration (FDA) guidance for

225 2.51 245

bioanalytical methods [30].

The LLOQ, established as the lowest concentration that can be measured with acceptable imprecision

141 $({\langle}20\%)$ and accuracy $({\pm}20\%)$, was assessed with four serial dilutions of spiked sample extract

Commenté [ZZ1]: Si plusieurs échantillons spikés mettre sample au pluriel

 containing 0.400 mM of each compound (six replicates). To validate a specific LLOQ, the chromatographic peak area of the corresponding analyte should be at least five times higher than the noise background of blank samples (*i.e.* peak to peak signal-to-noise ratio, S/N=5). The LLOQ were established as 422.5, 868.0, 876.2 and 919.5 µg/kg of liver for diphenylamine, tolylfluanid, propargite and phosalone, respectively.

 Selectivity was evaluated on ten different blank human cells extracts by checking that chromatographic peaks of pesticides and endrin (IS) did not co-elute with any endogenous compound. The evaluation consisted in verifying that specific retention times and SIM responses of pesticides were well discriminated from potential interfering signals.

 In addition, the matrix effect criterion was studied: extracting ten different blank matrices were extracted and spiked with the pesticides mixture at the predetermined LLOQ concentrations and quantified. Figure 2. presents both a typical chromatogram at the LLOQ and a chromatogram of a blank sample extract. The results were compared with those obtained for fortified aqueous extract of the same concentration levels. In the absence of acceptable limits defined in the FDA guidance and, as previously reported by Kadar et al. [31], to be acceptable, the deviation between the calculated and the nominal values should 157 be less than \pm 5%.

 To study the linearity for diphenylamine, tolylfluanid, propargite and phosalone, 6 replicates of 159 calibration standard samples from the validated LLOQ (0.1 μ M) to 25 μ M were prepared in hexane and analyzed by GC-MS.

161 Recovery rates of analytes were evaluated at low (0.1 μ M), medium (5.0 μ M) and high (50.0 μ M) 162 concentration levels. Despite the detector linearity being validated up to $25 \mu M$ only, a validation of the last concentration allowed expanding the concentration range of applicability of the method. Three replicates were prepared by fortifying thermally inactivated human hepatocytes samples before extraction and analysis. The ratio between mean peak areas from these samples and the one from post-extracted spiked samples allowed determining each pesticide percentage recovery.

 The imprecision (intra- and inter-day) and accuracy of the method were studied at four different 168 concentration levels: low (0.1 μ M), medium (5.0 μ M) and high (50.0 μ M). For intra-day imprecision and accuracy, five replicate samples per concentration level were prepared and analyzed on the same

 day. For inter-day imprecision, six different days on a 15 days period were chosen to repeat the preparation and analysis of sample duplicates at the same spiking levels. Intra- and inter-day imprecision were considered acceptable if the relative standard deviations (RSD%) were below 15% [31]. On the other hand, accuracy, expressed as the mean percentage deviation (Intra-day: Ar% and Inter-day: Br%) from the spiked value, was in accordance with the guidance when the yields ranged between 85 and 115% of the nominal concentrations.

 To check the stability of the compounds, degradation tests were carried out in triplicate using processed samples previously spiked at 5 µM with the pesticides mixture. Various storage conditions were experimented: 96 h in the autosampler tray at +25°C, after three freeze – thawing cycles from -20°C to 179 room temperature during 15h or, long term storage for either 1 month, 3 months or 9 months at -20°C. For this criterion, the FDA did not set acceptable limits. Then, after samples' analysis using freshly prepared calibration curves, imprecision and accuracy were respectively considered acceptable if below 15% and if between 85% and 115% of their nominal values as depicted by Gonzalez et al. [3]. In addition, after consecutive analysis of 8 samples, a blank sample and a standard at 0.1 µM were analyzed to check for any pollution or significant drift (> 10%) of the instrument signal.

3. Results and discussions

3.1. Method development

3.1.1. Sample treatment

 On the basis of previously published works [32-35], hexane and ethyl acetate seemed to be good solvents candidates for liquid/liquid extraction of the studied pesticides mixture from a 400 µL sample of thermally inactivated human liver cells. As presented in **Table 3**, the assay was optimized through a sequence of experiments aiming at establishing the best sample treatment prior to GC-MS analysis. First, the extraction efficiencies of the two hydrophobic solvents were compared on 400 µL spiked phosphate potassium buffer samples. As displayed in **Table 3**, in the absence of matrix**,** hexane proved to be the best compromise for the whole pesticides mixture. Then, the experiment was repeated on a sample spiked liver cells sample. However, due to the formation of a very compact emulsion, no clear layers appeared, thus preventing GC-MS analysis. Partial solubilization of proteins, highly present, could

 explain the outbreak of this emulsion. Thanks to their combined hydrophobic and hydrophilic properties, these proteins are attracted to both the non-polar solvent and the aqueous medium. Afterwards, comparison of simultaneous protein precipitation and LLE versus sequential protein precipitation and LLE was carried out. For this purpose, two efficient protein precipitators were compared: TCA and ammonium sulfate [36]. Globally, the results in **Table 3** show that regardless of the assay procedure (either simultaneous or sequential), recoveries obtained were higher with TCA compared to ammonium sulfate. In fact, contrarily to ammonium sulfate [37], TCA is known to help with the release of bound analytes into solution [38]. Indeed, at pH lower than the protein isoelectric point, TCA interacts with the positively charged amine group of proteins to form an insoluble salt [39]. In addition, as displayed in **Table 3**, simultaneous protein precipitation and LLE gave the best recoveries with TCA. This is probably because the vigorous shaking enabled a rapid solvent exchange between aqueous and organic solution while precipitation of protein was taking place simultaneously.

209

210 Table 3. Pesticides recoveries of a 10 μM concentration spiked samples 211 regarding the chemical treatment applied.

212 ^a LLE.

213 b Proteins precipitation followed by LLE.

214 C Simultaneous proteins precipitation and LLE.

215 Highest recoveries were obtained using simultaneous proteins precipitation and LLE after TCA and

216 hexane were added. Recoveries varied from 86 % for tolylfluanide to 98 % for endrin, with a maximum

- 217 RSD value of 1.8 %.
- 218 The final protocol was as follows: 100 µL of 20 % TCA and 400 µL of hexane were added to the human
- 219 hepatic preparation before 3 minutes vigorous shaking. Then, the sample was centrifuged for 4 minutes
- 220 at 16000 g and 4°C. Finally, the supernatant was transferred into a vial prior to GC-MS analysis.
- 221 *3.1.3 Optimization of splitless time, temperature and surge pressure using response surface* 222 *methodology*
- 223 The experimental results were analyzed by multiple linear regression to fit to the postulated model (1)

224 where Y is the instrumental response, X_1, X_2 and X_3 are the three independent variables described above 225 and β_i , β_{ii} and β_{ij} are the fitting coefficients.

226
$$
Y = \beta_0 + \beta_1 X_1 + \beta_2 X_2 + \beta_3 X_3 + \beta_{12} X_1 X_2 + \beta_{13} X_1 X_3 + \beta_{23} X_2 X_3 + \beta_{11} X_1^2 + \beta_{22} X_2^2 + \beta_{33} X_3^2
$$
 (1)

227

228 Table 4. Second polynomial equations obtained for the different pesticides.

230 Calculated coefficients of the pesticides response model, corresponding determination coefficients $(R²)$ 231 and adjusted determination coefficients (Adj. R²) are reported in **Table 4**. ANOVA results presented in **Table 5** demonstrated that the model was highly significant for each compound. The values of $R²$ and 233 Adj. R² (0.9180 and 0.8985 for diphenylamine; 0.8531 and 0.8183 for tolylfluanid; 0.8922 and 0.8667

234 for propargite; 0.9363 and 0.9212 for phosalone:) indicated a satisfying degree of correlation between

235 the observed and the predicted values.

Phosalone

239 SS = sum of squares, ddl = degree of freedom, MS = mean square, F = ratio and p = probability 240 level.

241

 Even if it was possible to find the optimal injector settings for a single response using surface response design, here the study aimed at finding a compromise for the simultaneous optimization of responses from diphenylamine, tolylfluanid, propargite and phosalone. Consequently, the multicriteria methodology established by Derringer and Suich [40] was implemented for the ensuing work. This methodology first involves the construction of an individual desirability function (di) for each compound. Then, as depicted by Eq. (2), the overall desirability function (D) is defined as the weighted geometric average of the individual desirability (di).

249
$$
D = \sqrt[4]{d_1 d_2 d_3 d_4}
$$
 (2)

250 A value of D close to 1 means that each individual pesticide's response is maximized and that the 251 corresponding levels of factors are globally optimum. As can be observed in **Fig. 1**, the desirability

- surfaces obtained by Eq. (2) show the effect of operating conditions on the overall pesticides
- chromatographic responses in the range of the investigated variables.
- Therefore, applying the numerical optimization function from the Statistica software to the Derringer's
- 255 desirability function allowed to find the following injection experimental conditions: X1 (271.7°C), X2
- (1.25 min) and X3 (490 KPa); which simultaneously maximized the signals of all the compounds.
-
- *3.2 Performance of the analytical method*
- No instrumental carry-over or significant drift were noticed during all the validation process.
- The absence of interfering peaks on chromatograms from ten different blank human hepatocyte samples,
- compared chromatograms acquired from reference standard in neat solvents led to the full validation of
- the method selectivity. Indeed, chromatographic signals of pesticides were always satisfactorily
- discriminated on the basis of their specific retention times and SIM responses. In addition, no matrix
- effects waere observed at the LLOQ level for all pesticides.
- **Table 6** lists the key elements of the method validation.
- 266 The determination coefficient values (R^2) were always greater than 0.999 indicating that there was a
- good correlation of linearity through the concentrations range used and a homoscedastic distribution of
- replicates at all levels.On the basis of recovery results ranging from 96.5% to 100.0%, with a maximum
- RSD of 3.4%, the implemented sample treatment was considered satisfactory.
- In addition, intra- and inter-day imprecision and accuracy exhibited acceptable values, ranging from 271 1.1% to 6.7%.
- Moreover, stability tests demonstrated that all the analytes did not degrade in the final extract while
- stored at -20°C or manipulated at room temperature. Finally, all the data met the validation criteria set
- by the FDA guidelines, rending this method suitable for future researches.
- 275 The evaluated LLOQ for diphenylamine (422.5 µg/kg), tolylfluanid (868.0 µg/kg), propargite (976.2
- 276 μ g/kg) and phosalone (919.5 μ g/kg) are higher than the values published by other authors in the context
- of biomonitoring studies. Indeed, for diphenylamine and propargite, Marschner et al. [41] obtained an
- LLOQ of 166.5 µg/kg on animal liver. Moreover, the phosalone LLOQ published by Russo et al. [26]
- was 15 µg/kg of human liver.

 Fig. 1. Response surface and contour plots of (a) injection time and injection temperature, (b) surge pressure and injection temperature and (c) surge pressure and injection time on the

283 overall desirability of *diphenylamine*, tolylfluanid, propargite and phosalone mixture

284 response.

287 Fig. 2. Hepatocytes liver extract chromatograms, A: blank sample, B: sample at the LLOQ.

288

289	Table 6. Results of the analytical method validation: linearity $(n = 6)$, recovery	

291 R%: percent recovery; RSD%: percent relative standard deviation.

292 Ar%: intra-day percent accuracy rate; Br%: inter-day percent accuracy rate.

293 SCt%: cold-war percent stability; SLt%: long term percent stability; SA%: autosampler 294 percent stability.

295

296 **4. Conclusion**

 In this paper, a simple and selective GC-MS method was developed to simultaneously determine the amounts of diphenylamine, tolylfluanid, propargite and phosalone in human liver samples. The instrumental settings were optimized to obtain the highest chromatographic responses for the above- mentioned compounds. The optimum splitless injection conditions obtained through RSM and global Derringer's desirability function were 271.7°C, 1.25 min and 490 KPa for injection temperature, splitless time and surge pressure, respectively. Convenient optimization of simultaneous protein precipitation and LLE extraction conditions allowed to achieve very good recoveries and efficient sample purification using a user-friendly sample treatment. This one-step LLE and cleanup method followed by GC-MS analysis complied with current FDA

306 requirements and showed satisfying selectivity, linearity, recovery, precision and accuracy. The limits

- of quantification obtained in this work are satisfying, considering previously published works related to
- human and animal liver samples.
- This procedure is thus appropriate for the monitoring of the parent compound loss during *in vitro* human
- liver metabolism studies. Finally, upgrading it to a more sensitive technology (MS/MS) it may also be
- useful for conducting human or mammalian biomonitoring studies.
-

Acknowledgements

- This work was funded by the Agence Nationale de la Recherche (under reference ANR-2008-CESA-
- 016-01) by the Office National de l'Eau et des Milieux Aquatiques (ONEMA) and by the Agence
- Nationale de SÉcurité Sanitaire de l'alimentation, de l'environnement et du travail (ANSES) in the
- frame of the MEPIMEX Project.

 33. Huang, Z., Li, Y., Chen, B., Yao, S., 2007. J. Chromatogr. B Analyt. Technol. Biomed. Life Sci. Simultaneous determination of 102 pesticide residues in Chinese teas by gas chromatography-mass spectrometry. 853, 154-62.

 34. Nakajima, T., Tsuruoka, Y., Kanda, M., Hayashi, H., Hashimoto, T., Matsushima, Y., Yoshikawa, S., Nagano, C., Okutomi, Y., Takano, I., 2015. Determination and surveillance of nine acaricides and one metabolite in honey by liquid chromatography-tandem mass spectrometry. Food Addit. Contam. Part A. Chem. Anal. Control Expo. Risk Assess. 32, 1099-104.

 35. Picó, Y., la Farré, M., Segarra, R., Barceló, D. 2010. Profiling of compounds and degradation products from the postharvest treatment of pears and apples by ultra-high pressure liquid chromatography quadrupole-time-of-flight mass spectrometry. Talanta, 81, 281-93.

 36. Polson, C., Sarkar, P., Incledon, B., Raguvaran, V., Grant, R., 2003. Optimization of protein precipitation based upon effectiveness of protein removal and ionization effect in liquid chromatography–-tandem mass spectrometry. J. of Chromatogr. B. 785, 263-275.

 37. Kaufmann, A., Butcher, P., Maden, K., Walker, S., Widmer, M., 2011. Development of an improved high-resolution mass spectrometry based multi-residue method for veterinary drugs in various food matrices. Anal Chim Acta., 700, 86-94.

- 38. Chiaochan, C., Koesukwiwat, U., Yudthavorasit, S., Leepipatpiboon, N., 2010. Efficient hydrophilic interaction liquid chromatography-tandem mass spectrometry for the multiclass analysis of veterinary drugs in chicken muscle. Anal. Chim. Acta. 682, 117-29.
- 39. Souverain, S., Rudaz, S., Veuthey, J.L., 2004. Protein precipitation for the analysis of a drug
- cocktail in plasma by LC-ESI-MS. J. Pharm. Biomed. Anal.35, 913-20.

