

Development of a simple multiresidue extraction method for the quantification of a wide polarity range list of pesticides and transformation products in eggs by liquid chromatography and tandem mass spectrometry

Vincent Dufour, Laure Wiest, Sylvain Slaby, François Le Cor, Lucile Auger, Olivier Cardoso, Laurence Curtet, Laure Pasquini, Xavier Dauchy, Emmanuelle Vulliet, et al.

▶ To cite this version:

Vincent Dufour, Laure Wiest, Sylvain Slaby, François Le Cor, Lucile Auger, et al.. Development of a simple multiresidue extraction method for the quantification of a wide polarity range list of pesticides and transformation products in eggs by liquid chromatography and tandem mass spectrometry. Journal of Chromatography A, 2020, 1628, pp.1-9. 10.1016/j.chroma.2020.461447 . hal-03011261

HAL Id: hal-03011261 https://hal.science/hal-03011261

Submitted on 20 Nov 2020

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

1	Development of a simple multiresidue extraction method for the quantification of
2	a wide polarity range list of pesticides and transformation products in eggs by
3	liquid chromatography and tandem mass spectrometry
4	Vincent DUFOUR ^{a b} , Laure WIEST ^b , Sylvain SLABY ^a , François LE COR ^{a c} , Lucile AUGER ^{a d} ,
5	Olivier CARDOSO ^e , Laurence CURTET ^d , Laure PASQUINI ^c , Xavier DAUCHY ^c , Emmanuelle
6	VULLIET ^b , Damien BANAS ^a
7	
8	^a Université de Lorraine, INRAE, URAFPA, F-54000 Nancy, France
9	^b Univ Lyon, CNRS, Université Claude Bernard Lyon 1, Institut des Sciences Analytiques, UMR
10	5280, 5 Rue de la Doua, F-69100, Villeurbanne, France
11	^c LHN, Laboratoire d'Hydrologie de Nancy, ANSES, 40 Rue Lionnois, 54000 Nancy, France
12	^d Office Français de la Biodiversité - Montfort, 01330 Birieux, France
13	^e Office Français de la Biodiversité – Unité Sanitaire de la Faune, 9 avenue Buffon, 45071 Orléans,
14	France
15	
16	VD: <u>vincent.dufour@univ-lorraine.fr</u>
17	LW: <u>laure.wiest@isa-lyon.fr</u>
18	SS: <u>sylvain.slaby@univ-lorraine.fr</u>
19	FLC: <u>francois.le-cor@univ-lorraine.fr</u>
20	LA: <u>lucile.auger@univ-lorraine.fr</u>
21	OC: <u>olivier.cardoso@ofb.gouv.fr</u>
22	LC: <u>laurence.curtet@ofb.gouv.fr</u>
23	LP: <u>laure.pasquini@anses.fr</u>

- 24 XD: <u>xavier.dauchy@anses.fr</u>
- 25 EV: <u>emmanuelle.vulliet@isa-lyon.fr</u>

- 26 DB: <u>damien.banas@univ-lorraine.fr</u>
- 27
- 28 Corresponding author: Damien BANAS, URAFPA, Unité de Recherche Animal et Fonctionnalités
- des Produits Animaux, Université de Lorraine, INRA, 2 Avenue de la Forêt de Haye, 54500
- 30 Vandœuvre-lès-Nancy, Nancy, France; damien.banas@univ-lorraine.fr; +33 372 744 151;

31

33 Abstract

Waterfowl populations have been decreasing in Europe for the last years and pollution appears to be 34 one of the main factors. This study was conducted to develop a single sensitive and robust analytical 35 method for the monitoring of 2 fungicides, 15 herbicides, 3 insecticides and 24 transformation 36 products in wild bird eggs. One of the major challenges addressed was the characterization of 37 chemicals with large logP range (from -1.9 to 4.8). A total of 11 different extraction parameters were 38 39 tested in triplicate to optimize the extraction protocol, on generic parameters, buffer addition and use of clean-up steps. Quantification was based on matrix-match approach with hen eggs as reference 40 41 matrix (34 analytes with r²>0.99). Particular attention was payed to matrix effects (72 % on average), quantification limits (0.5 to 25 ng.g⁻¹ dry mass / 0.2 to 7.5 ng.g⁻¹ fresh mass) and extraction yields 42 (46 to 87 % with 25 analytes up to 70 %) to ensure the relevance of the method and its compatibility 43 with ultra-trace analysis. It led to a simple solid/liquid low temperature partitioning extraction 44 method followed by LC-MS/MS. Analysis of 29 field samples from 3 waterfowl species revealed 45 that eggs were slightly contaminated with pesticides as only one egg presented a contamination 46 (terbutryn, herbicide, 0.7 ng.g⁻¹) and confirmed the relevance of the method. 47

48

49 **Keywords**: eggs; pesticides; multiresidue; SLE-LTP; LC-MS/MS; polar transformation products

51

52 **1. Introduction**

Nowadays we assist to the decline of bird population and diversity across the world [1,2], and plant 53 protection products (PPP) are identified as one of the main factor of this decline [3–5], affecting 54 pollination, seed dispersal, pest control and, thus, economy [6]. The analysis of wild bird eggs could 55 allow to acquire more impregnation data by PPP while preserving reproductive capacities of sampled 56 57 species and considering a highly vulnerable stage of development. Egg represents an important analytical challenge: it contains ≥ 10 % of both proteins and lipids [7] which are well known to cause 58 59 matrix effects in extractions and analysis of organic chemicals [8-11], affecting sensitivity, robustness and reliability of methods. Moreover, as their contamination depend on the ability of 60 pesticides to cross different biological membranes [12], contamination is expected to be low and the 61 final protocol should adress trace levels (ng.g⁻¹). 62

Several works already focus on the caracterisation of PPP in eggs. Most [7,12–17] are using a 63 64 QuEChERS (quick, easy, cheap, effective, rugged, and safe) approach, firstly developed by Anastassiades and coworkers [18]. They monitored between 10 and >100 different pesticides based 65 on LC-MS/MS or GC-MS/MS, with limits of quantification (LOQ) from 0.5 to 10 ng.g⁻¹. However, 66 nature of solvent, buffer salt and clean-up step can vary a lot between studies, according to the nature 67 of monitored PPP. On the other side, Luo et al. [19] successfully developped a solid-liquid 68 extraction with low temperature partitionning (SLE-LTP) that have to be mentioned as it allow the 69 quantification of 325 pesticides and some transformation products (TP) with most of LOQ below 1 70 ng.g⁻¹ without using the salting-out effect. 71

Despite the very good performances of the mentioned studies, 3 major limits can be put forward. First, most of characterized analytes are relatively hydrophobic as they are often organochlorines, organophosphorines, permethrines, triazines or chloracetanilides, with logP close or up to 3. The only exception is the caracterization of neonicotinoids that are much more polar ($2 > \log P > -1$

according to [20]). Indeed, if less-polar pesticides are more incline to accumulate in organisms, 76 especially via trophic chain, they tend to be ban and replace with more polar ones which can impact 77 78 aquatic ecosystems and wildlife in general and, thus, have to be monitored. Secondly, recent studies showed the importance of not only considering PPP in environmental investigations, but also their 79 TP as they may have lowest predicted no effect concentrations than active substances [21]. 80 Moreover, they often are more polar than the active substance (weak acids or basis), which reinforce 81 82 the first argument of disposing from a method for the caracterization of large polarity range analytes. In studies mentioned above, no TP were monitored. Zhang et al. [22] integrated TP in their analyse 83 84 but they represent less than 10 % of the analytes. Finally, all previous prototocols were based on the analysis of mixed fresh eggs stored at 4 °C or -20 °C. In some cases this is not the easiest way to 85 store them over the long term as they deteriorate at room temperature and necessitate both 86 temperature control and energy source. This makes it difficult to store large series of samples over 87 the long term until they are analysed, while at the same time ensuring their proper conservation. 88

As a consequence, it is necessary to dispose of a quick and reliable method for the monitoring of 89 traces of large polarity range PPP and TP in complex matrices, and with an appropriate long term 90 storage. This study focuses on the development of a sole extraction method for 44 chemicals in 91 freeze-dried eggs: 20 active substances (herbicides, fungicides, and insecticides) but also 24 highly 92 polar and ionic TP. The method, based on QuEChERS extraction and a matrix-matched 93 quantification, should allow to quantify these chemicals at traces levels (<10 ng.g⁻¹ dry mass), with 94 good extraction recoveries (ER, comprised between 80 % and 120 %) and least possible matrix 95 effects, taking in account the samples complexity. It was tested on a cohort of 29 eggs of 3 different 96 waterfowl species obtained from 4 different ponds in Sologne region (France) to be sure to provide a 97 functional tool for evaluating the impregnation of bird eggs by PPP and TP. 98

100 **2. Materials and methods**

101 2.1. Standards and reagents

Ultrapure water (Fisher Chemical, Geel, Belgium), acetonitrile (Honeywell, Seelze, Germany) and 102 heptane (Merck, Darmstadt, Germany) were of LC-MS quality. QuEChERS citrate buffer (4 g 103 magnesium sulfate, 1 g sodium chloride, 1 g sodium citrate dihydrate, 0.5 g sodium hydrogen 104 citrate), acetate buffer (6 g magnesium sulfate, 1.5 g sodium acetate) and NaCl buffer (6 g 105 106 magnesium sulfate, 1.5 g sodium chloride), and dispersive solid phase extraction salts EMR-lipids (Enhanced Matrix Removal-Lipid) were purchased from Agilent (Santa Clara, USA). Dispersive 107 108 solid phase extraction reagent PSA (primary-secondary amine) was obtained from Chromebond (Düren, Germany) and magnesium sulfate (MgSO₄) from Sigma-Aldrich (Saint Louis, USA). Formic 109 acid (FA) and ammonium sulfate were purchased from Biosolve-chemicals (Dieuze, France). 110 Zirconium oxide grinding balls (diameter 3 mm) were obtained from Fisher Scientific (Illkirch, 111 France). 112

Analytical standards were purchased from A2S (Martignas-sur-Jalle, France), Honeywell (Seelze, Germany), HPC standards (Cunnersdorf, Germany), LGC group (Teddington, UK), Sigma-Aldrich (Saint-Louis, USA) and TCI Europe (Anvers, Belgium). All were high purity level (> 96 %), with most >99 %. Individual solutions made in acetonitrile (ACN) (100 mg.L⁻¹) were stored at -18 °C. They were mixed to constitute a 5 μ g.L⁻¹ solution used to develop and validate the extraction protocol. Mix solutions were regularly characterized to check their stability over time and data were compiled in a control chart. Stored at 18 °C, they were stabled at least for 6 months.

120

121 2.2. Sampling strategy

Egg sampling was carried out in accordance with the country's legislation (Interministerial decree NOR: DEVP0650375A, Environmental Code L424-11, [23]). It was performed during spring 2018 in the Sologne region (France). The prospection was done in four ponds to increase the probability to find nests and to improve representativeness of studied sites. The sampled species were typical water birds found in some French regions with different ecological behaviours: *Anas platyrhynchos* (Mallard duck, Anatidae), *Aythya ferina* (Common pochard, Anatidae) and *Fulica atra* (Eurasian coot, Rallidae). Only eggs incubated for less than 7 days were sampled with gloves. In the laboratory, each egg was weighed before and after being broken. The content was roughly mixed and frozen at -20 °C in aluminum tray. Then, it was freeze-dried (Reacti-Vap from PIERCE, Bellafonte, USA), grinded and stocked at room temperature.

Eggs from hens kept in battery were used as reference matrix for protocol development, protocol blanks, or positive controls in extraction batches when spiked with analytical standards. They were pretreated as field samples.

135

136 2.3. QuEChERS-based extraction

Different conditions were tested to optimize the extraction procedure. The first applied protocol 137 (Ref#1) (Fig. 1) was based on a simple SLE-LTP, aiming to increase its complexity until reaching 138 OuEChERS approaches. Two hundred milligrams of reference matrix were weighted in a 50 mL 139 polypropylene centrifuge tube (VWR, Radnor, USA) and spiked with 50 µL of a mix solution 140 containing all analytes and 75 µL of a mix solution containing all isotope labeled standards, both at 141 nominal concentration 500 µg.L⁻¹. Solution solvent was evaporated under gentle nitrogen flux for 10 142 min. A volume of 5 mL of ultrapure water was added to the sample, mixed together for 5 min at 16 143 Hz with Geno/Grinder (SamplePrep 2010 Geno/Grinder from SPEX, Costa Mesa, USA) to re-144 hydrate the sample. Then 10 mL of ACN and 5 mL of heptane were added in the centrifuge tube and 145 mixed together for 5 min at 16 Hz with Geno/Grinder to extract analytes and remove part of non-146 polar interferents. Samples were then centrifuged for 5 min at 9500 RCF at 20 °C (SIGMA 3K3OH, 147 Fisher, Osterode am Harz, Germany). Six milliliters of the ACN/water mix were transferred in new 148 centrifuge tube containing PSA, handshaked for 30 s and centrifuged 2 min 9500 RCF to separate 149

salts from extract. Final extract was aliquoted: 1 mL was transferred in a 2 mL glass injection vial 150 (from Supelco, USA) and stored at -20 °C until analysis. Generic parameter modifications were 151 tested to improve the protocol: replacing ACN with methanol (MeOH), increasing the mass of 152 sample to 0.5 g, increasing the extraction time to 10 min, acidification of ACN (0.2 %FA), 153 introducing QuEChERS buffers (NaCl, acetate or citrate) and replacing PSA d-SPE salts with EMR-154 lipid salts or MgSO₄. The comparison of samples spiked before or after extraction (in triplicate and 155 156 by a single operator) allowed to calculate extraction yield and estimate limits of quantification (LOQ) of each tested condition. Twenty microliters of the final extract were diluted in 180 µL of 157 158 ultrapure water + 0.1%FA (v:v) before injection.

The optimized protocol was used to evaluate if sample preparation (e.g. freeze drying) affected the representativeness of samples. Fresh hen eggs were broken, weighted and mixed and 57 g of it was spiked with 1 mL of the 500 μ g.g⁻¹ standards mix solution while 54 g were supplemented with 1 mL of ACN as reference. Both were strongly mixed for 2 min to ensure good samples homogeneity, aliquoted in 3 distinct samples and then freeze dried (Water represented 80 % ± 2 % of the total egg mass; n = 6). Reference samples were spiked with 50 μ L of mix solution (taking into account dehydration) and all samples were extracted following final protocol (section 3.2. and Fig. 1).

166

167 2.4. LC-MS/MS

The following analytical method was initially developed in the Laboratory of Hydrology of Nancy (ANSES-LHN, Nancy, France) for the quantification of 88 pesticides and transformation products in water bodies, among which the 44 herein targeted chemicals. Analysis was performed with a 1290 Infinity UPLC (Agilent Technologies, Massy, France) coupled with 5500 QTRAP (Sciex, Villebon sur Yvette, France) in both positive and negative modes, using electrospray ionization. Separation was performed at 0.4 mL.min⁻¹ at 50 °C over a Xselect® HSS T3 column 2.1x5 mm, 2.5 μm (Waters, Saint Quentin-en-Yvelines, France). The column was protected with a KrudKatcher Ultra

(Phenomenex, Le Pecq, France) and a Xselect® HSS T3 precolumn 2.1x5 mm, 2.5 µm (Waters, 175 Saint Quentin-en-Yvelines, France). Ultrapure water + 0.1 %FA (A) and ACN + 0.1 %FA (B) were 176 177 used as mobile phases in positive mode while 1 mM ammonium acetate in ultrapure water (A) and ACN (B) were used in negative mode. The gradient was the same for both methods, started with 10 178 % B for 1 min, increased to 100 % B in 7 min, held 6 min, back to 10 % B in 1 min and held 4 min. 179 Temperature of the mass spectrometer source, nebulizer gas, heating gas, were respectively settled at 180 181 550 °C, 40 psi and 50 psi. Ion spray voltages were optimal at 5500 eV for positive ionization and -4500 eV, for negative ionization. 182

183 Analysis parameters of each compound, including multiple reaction monitoring (MRM) transitions, are summarized in SI Table S1. Analyte infusions were performed in MeOH:H₂O mix (50:50), with 184 formic acid in positive mode and with ammonium acetate in negative mode, to ensure good signal 185 intensity. Transitions were automatically optimized with Analyst 1.6.2 and the 2 most intense and 186 sensitive ones were kept for each analyte. Nature of phases, flow, injection volume and column 187 temperature were tested to ensure stable retention times, good separation, and good intensity on the 188 chromatographic peaks in both modes. The most sensitive transition of each analyte was kept for 189 quantification and the other as confirmation. 190

191

192 2.5. Method validation

The matrix-matched approach was validated over six days with two different operators (three days each), following the experimental organization described by Giroud *et al.* [24] and the SANTE Guidance [25], and changing extraction materials every day. Hen eggs samples spiked with analytical standards from LOQ/2 to 7xLOQ (6 concentration levels) and a matrix blank (reference matrix extracted without standards spiking) were extracted each day for three days by each of the operators to determine regression law for each of the analytes. Each day one of the concentration levels (LOQ, 3xLOQ, 7xLOQ) was extracted in triplicate to document the intermediate precision of all parameters. Matrix effect (ME, in %), extraction recoveries (ER, in %) were calculated as follows
(equation 1 and 2):

$$202 \qquad ME = \frac{Area \ of \ post-extraction \ spike}{Area \ of \ solution \ reference} \times 100 - 100 \tag{1}$$

203
$$ER = \frac{Area \ of \ pre-extraction \ spike}{Area \ of \ post-extraction \ spike} \times 100$$
(2)

Herein ER were determined by comparison of artificial samples spiked before and after the extraction (LOQ, 3xLOQ, 7xLOQ) on two triplicates for each concentration level. Reference matrix samples spiked after extraction were also compared with the analytical response of a solution of equivalent constitution to calculate ME for each analyte.

208

209 2.6. *Method application*

The developed method was applied on 29 waterfowl eggs (section 2.2.) for pesticides quantification, with a six points matrix match calibration curve for each sequence (same concentrations as described in section 2.3.) and following the European Commission recommendations [25]. A blank (to track contamination) and an artificial sample made from reference matrix spiked with analytes at LOQ level (to calculate recoveries and LOQ) were also extracted. To preserve analytical relevance of the results, protocol blanks were taken into account by increasing LOQ at three times the concentration in blanks when a contamination was identified.

217

218 **3. Results and discussion**

219 *3.1. Optimization of extraction conditions*

Sequential modifications of the protocol were tested compared to the protocol described in section2.3. of this article. Performances were evaluated with the aim of centering ER between 80 and 120

222 %, and keeping LOQ below 10 ng.g⁻¹ to analyze contaminant traces. All these performances are 223 described in sections below and are summarized in Fig. 2.

The nature of the solvent, the mass of matrix extracted, contact times between sample and solvents can all play role in the determination of a performant extraction protocol. A total of 4 different versions of the reference protocol based on generic parameter changes were tested in triplicates to evaluate their impact on LOQ and ER: (1) ACN was replaced with MeOH, (2) the mass of matrix was increased from 0.2g to 0.5g, (3) the mixing time was increased to 10 min, and (4) ACN was acidified with 0.2 % AF.

230

231

3.1.1. <u>Preliminary tests</u>

ACN is a commonly used organic solvent for the extraction of polar pesticides from solid matrices 232 and more particularly in QuEChERS approaches [12–14,17,18,24,26] as it is easily separable from 233 water with the adjunction of salts [18]. MeOH was tested in parallel as it is also a polar solvent that 234 may be useful for the extraction of most polar molecules. The use of MeOH instead of ACN 235 decreased the ER by 40 %, and only 2 analytes presented good ER (120 % > ER > 70 %) while 8 236 were concerned for Ref#1. Otherwise, the solvent change did not affect the global LOQ as 26 and 22 237 analytes presented LOQ below 10 ng.g⁻¹ for Ref#1 and MeOH condition respectively. The difference 238 in ER was probably due to matrix effects: (1) MeOH can extract a more important part of the matrix, 239 leading to strong signal extinction and low ER determination, and (2) MeOH is well known to 240 increase the ionization of chemicals in the mass spectrometer source which may have led to more 241 important noise on chromatograms even in MRM acquisition. This solvent may also induce an 242 esterification of some of the acidic molecules as proposed by Dong et al. [27] leading to some of the 243 low ER observed. 244

The mass of analyzed sample has to be considered: this parameter increases the quantity of available analytes for extraction and can therefore decrease the LOQ, but it may also increase the quantity of interferents extracted, leading to more matrix effects. In our case, increasing sample mass from 0.2 g
to 0.5 g strongly reduced the LOQ with respectively 22 and 35 analytes quantifiable below 10 ng.g⁻¹.
Moreover, the number of analytes quantified above 50 ng.g⁻¹ decreased from 4 to 2. At the same
time, 5 more molecules presented good ER when matrix mass was increased.

Mixing time is also an important parameter: an increase leads to better interaction between solvent and matrix and can contribute to a better ER. Extend the extraction time from 5 to 10 min increased the LOQ by a factor 2 on average with still 21 analytes quantifiable under 10 ng.g⁻¹. However, ER were improved as 14 analytes presented good recoveries and only 4 of them were up to 120 %, with no outlier >200 % which ensure a better quantification. Moreover, compared to the Ref#1, longer mixing time limited the variability of the extraction process as the %RSD on extraction recoveries decreased from 40 % on average to only 11 %.

Some of the analyzed chemicals (especially TP) are weak acids that may bonded to the matrix in the 258 Ref#1 conditions. Acidification at pH close or lower than the pKa of acidic molecules may increase 259 extraction performances by fostering their neutral form, less susceptible to interact with matrix or 260 water and ease the transfer to organic solvent [28]. The acidification of the ACN with 0.2 % of FA 261 (v:v) increased LOQ by a factor 2 on average but 21 chemicals presented LOQ below 10 ng.g⁻¹, but 262 ER were close to Ref#1 with only 7 analytes with ER comprised between 70 % and 120 %. 263 Nevertheless, acidification reduced the outliers as no analyte presented ER up to 200 %. It is also 264 important to mention that acidification induced better separation between the different fractions and 265 a more compact solid deposit after centrifugation, leading to quicker and easier-to-perform clean-up 266 step. 267

268

269

3.1.2. Determination of a new reference protocol

The reference protocol presented higher LOQ and lower ER performances than any of the alternative extraction tested. As a consequence, Ref#1 was modified to integrate key strengths of previously evaluated parameters before testing new buffers and clean-up steps: mixing times were increased to 10 min and the ACN was acidified with FA (+0.2 % v:v). As extracts were significantly more pigmented with 0.5 g of sample than with 0.2 g, and as pigmentation could be linked with matrix effects and presence of unwanted molecules, this parameter was temporarily set at 0.2 g for Ref#2 to preserve the analytical system. New protocol presented good quantification limits as 30 chemicals had LOQ inferior to 10 ng.g⁻¹. In parallel ER comprised between 80 % and 120 % concerned 19 chemicals and only 5 ones are up to 120 %.

279

280

3.1.3. Impact of buffer addition and selection of d-SPE cleanup

All parameters tested in this section were built from the new reference Ref#2 (section 3.1.5.) as it 281 represented a significant improvement of initial protocol (Fig. 1). The addition of salts as in 282 QuEChERS approaches can lead to improve ER of the protocol, especially for the most polar 283 chemicals [30,33,36,40,41,49] because of the salting-out effect that force the transfer of chemicals 284 from the water layer to the organic solvent layer. The buffer used in the original QuEChERS method, 285 NaCl [18], the acetate buffer recommended by the American standard (AOAC) [29] and the citrate 286 buffer used in the European standard EN 15662 [30] were tested. Compared to Ref#2, quantification 287 limits were higher when a buffer was used: the number of molecules with LOQ inferior to 10 ng.g⁻¹ 288 for NaCl, acetate and citrate condition was respectively 23, 18 and 19, and ones with LOQ up to 100 289 ng.g⁻¹ concerned 2, 6 and 5 analytes respectively. Concerning ER, protocols had 18, 22 and 17 290 molecules with good rates. Acetate buffer which is recommended in standardized American 291 OuEChERS protocols [29] presented the best performances in terms of ER while NaCl allowed to 292 get the lowest LOQ. However, compared to Ref#2, adjunction of any buffer almost systematically 293 decreased the average extraction recovery of acetochlor ESA /alachlor ESA, acetochlor OXA, CGA 294 357704, CGA 368208, dimetachlor ESA, dimethenamid ESA, flufenacet ESA, flufenacet OXA, 295 MCPA, metazachlor ESA, metolachlor ESA and metolachlor OXA to 15 % on average. The addition 296

of important quantities of salts may increased the pH up to the pKa of these molecules, leading to the dominance of the ionic forms, more strongly bounded to the water layer and more hardly extractible with ACN [28], which may explain the poor observed recoveries. Buffer was removed from final protocol because transformation products were strongly impacted by this phenomenon and because the priority of this study was to monitor these TP in biota as very few data are available.

D-SPE step can increase method performance by purifying the extract and removing interferences, 302 303 leading to less matrix effects, but the use of unsuitable salts can also remove analytes. Lipids are often considered as responsible for matrix effect in biotics matrices and can be partially removed 304 305 with d-SPE [10]. PSA is one of the most widely used d-SPE cleanup mix for the removal of fatty acids [18,31,32] and was *de facto* integrated into the protocols Ref#1 and Ref#2. EMR-lipids mix 306 was tested as it is specifically designed to remove lipids without affecting extraction recoveries of 307 PPP [32] and presented good results for extraction of PPP in eggs [19]. A condition with only 308 MgSO₄ was also tested as this salt is often used in association with the PSA and EMR-lipids to 309 remove water from samples [18,24,32] and force the transfer of analytes from water to organic 310 solvents. The two new conditions tested respectively had 23 and 27 molecules with LOO lower than 311 10 ng.g⁻¹. ER were acceptable (80-120 %) respectively for 17 and 35 chemicals and none of the 312 experimental conditions presented ER up to 120 %. If all conditions were comparable in terms of 313 LOQ, the addition of MgSO₄ clearly stands out with both good ER for large range of chemicals, 314 including TP, while those of EMR-lipid were close to Ref#2 with no great improvement, nor 315 degradation of performances apart for benzamide (0 ER, no LOQ definable). In opposition to Luo et 316 al. [19], it seemed that EMR-Lipid causes the removal of many chemicals causing the diminution of 317 ER by average 24 % compared to Ref#2, probably because we studied a wider range of polarity. On 318 319 the contrary, the use of MgSO₄ only, forced the transfer of chemicals from water to ACN, even for the most polar ones, and part of interferences may have precipitated as chromatograms showed less 320 background noise and as no perturbation of LOQ was highlighted. 321

322 *3.2. Impact of freeze drying on extraction yields*

As specified in section 2.2., due to logistical limitations, sampled eggs were broken, mixed and freeze-dried in laboratory for appropriate storage. As the analytical method was not developed yet when samples were collected (from April 2018 to July 2018 and analyzed in May 2019), they could not be analyzed fresh as in other studies [13,14,16,17] or after being stored at -20 °C [12,15]. The impact of freeze drying was evaluated as stated in section 2.3. and results are presented in Fig. 3.

328 Freeze drying did not affect the representativeness of the sample for most contaminants as they presented extraction yields between 80 % and 120 %. However, it led to a loss for flufenacet (53 %), 329 330 dimetachlor (66 %), fenthion (73 %) and s-metolachlor (78 %). CGA 357704 was also affected as it presented extraction yields of 80 % but with an RSD of 39 % while it was below 20 % for other 331 chemicals. These chemicals may have bounded more strongly with the matrix during the time of 332 freeze-drying (48 h) that was much longer than the usual contact time between the native solution 333 and the matrix for the constitution or artificial samples (10 min). In this case, the developed protocol 334 was not efficient for the characterization of these molecules. Studies on the effect of an increasing 335 contact time between matrix and analytes would confirm these observations to evaluate the role of 336 this parameter in the representativeness of artificial samples toward field samples, as it was identified 337 as an important limit in matrix match approaches [33]. 338

339 3.3. Analytical performances and validation of the method

340 Method performance parameters are presented in Table 1.

Linear regression with r² up to 0.99 was successfully determined for 34 chemicals. Among the 44 analytes, 33 analytes presented LOQ $\leq 10 \text{ ng.g}^{-1}$ and all were $\leq 25 \text{ ng.g}^{-1}$. Moreover, taking into account the hydration rate of samples, these LOQ are 72 % lower on average when considering the fresh mass concentration which is in the high LOQ trends of Luo *et al.* [19] (0.03 to 23 ng.g⁻¹). Considering matrix effects, only metazachlor ESA (225 %) presented a signal enhancement and only 3 chemicals presented signal extinction: prosulfocarb (34 %), dimethenamid ESA (40 %) and metolachlor ESA (48 %). Indeed, matrix effects were low, with 72 % average, which is consistent with literature as signal extinction seemed to affected most part of chemicals in complex matrices [26,34]. However, matrix effect did not affect the relevance of the analysis as matrix match quantification should take them into account if the reference matrix is well chosen (discussed section 3.4.).

ER at LOQ were comprised between 46 % (metolachlor OXA) and 87 % (imidacloprid), with 25 analytes with ER comprised between 70 % and 120 % as recommended by European Commission [25].However, the inter-day precision was below or close to 20 % (upper values were 21 % for CGA 357 357704 at LOQ, 22 % for metsulfuron-methyl at LOQ and 22 % for metolachlor OXA at 3xLOQ), assuming that the quantification was precise for the tested chemicals. Definitive LOQ according to SANTE guidelines are referenced Table 1.

No signal was detected in blanks for most of the chemicals monitored (Table 1). Contaminations ranged from 0.1 ng to 8.9 ng, corresponding to 0.2 and 17.8 ng.g⁻¹ dry mass respectively. No chemical was systematically quantified in blanks but 7 pesticides and TP presented unique contamination over the 6 extraction days, 2 were quantified twice in blanks, and 4 were quantified in more than half of blanks (fipronil sulfone, prosulfocarb, s-metolachlor and terbutryn).

363

364 3.4. Field sample analysis

Spiking with isotope labeled standards was used to check if the reference matrix (hen eggs) used for quantification was close enough to field samples, by comparing their respective areas. Average ratio was comprised between 80 % and 120 % for all internal standards except for D5-acetochlor ESA (124 ± 19 %, ESI neg, 4.9 min) and D6-chlorpyriphos methyl (123 ± 45 %, ESI pos, 8.6 min). It showed a good correspondence between the reference matrix and field samples, and so proved the good reliability of the method. PCA (Fig. S1) was performed on these ratios to determine if clusters of samples can be identified taking into account egg species (XLstats by Addinsoft Inc, New York, USA). The PCA was based on Pearson regression ($\alpha = 0.005$), axis 1 (47.2 %) and 2 (29.0 %) explain 76.1 % of the total variability but no clear clusters could be identified. It attested of a relatively homogeneous matrix complexity between different species, but also between all samples. It confirmed the good reliability of the developed protocol with field samples, and so the suitability of this method for the monitoring of wide polarity range of pesticides in eggs.

Over the 29 collected and analyzed environmental samples, only 1 coot egg presented a contamination superior to LOQ for terbutryn (3.1 ng.g⁻¹ dry mass) for LOQ at 2.4 ng.g⁻¹ taking in account contamination as mentioned sections 2.6. and 3.3. Considering the hydration ratio of the sample (78.5 %), it corresponded to 0.7 ng.g⁻¹ fresh mass contamination or 22 ng for the whole egg.

This molecule presented a contamination induced by the extraction protocol (estimated to 0.13 \pm 381 0.02 ng.g⁻¹ dry mass in the extraction batch of the considered sample). It was maybe due to the actual 382 use of terbutryn as a biocide in some paintings, for the protection of building materials against 383 moisture or microalgae [35-37] which can create contamination in laboratories. However, the 384 quantity found in blanks was quite constant and their area was systematically removed for 385 quantification. Moreover, the LOQ was increased at 3 times the concentration of terbutryn found in 386 blanks, so the communicated results were reliable. The quantification of terbutryn in coot egg 387 confirmed that the protocol developed is adapted to quantification of pesticides in waterfowl eggs 388 with good sensitivity and ER, and with acceptable matrix effects even if the contamination was 389 below the quantification limit in most cases. 390

391

392 **4. Conclusion**

We successfully developed a sensitive, robust, accurate and reliable extraction protocol for the characterization of 20 PPP and 24 of their TP in eggs, based on matrix matched quantification. About 80 % of LOQ were below 10 $ng.g^{-1}$ dry mass and more than 50 % of extraction yields were comprised between 70 % and 100 %, which makes it compatible with trace monitoring in complex matrices. Special attention was paid to matrix effect as they were evaluated on artificial samples for
method validation, but also directly in samples with the survey of isotope labeled standards signals.
The application of the developed method to 29 real waterfowl eggs sampled in the French region
Sologne, revealed very low contamination for the monitored PPP and TP in this biotic matrix. Only
one egg presented a concentration higher than the LOQ of terbutryn.

402

403 Acknowledgement

This project was funded by the "Office Français de la Biodiversité" (Vincennes, France). The authors want to thank Serge Bourdais from the "Fédération Départementale des Chasseurs du Loir et Cher" (Vineuil, France), the "Fédération des Chasseurs de l'Ain" (Bourg-en-Bresse, France), the owners of ponds who have granted us access to their properties, as well as our collaborators who helped us search and collect eggs. Thanks also to the "Laboratoire d'hydrologie de Nancy" of the ANSES (Nancy, France) for the provision of their analytical methods. The authors finally want to thank Alain Iuretig (Université de Lorraine, Vandoeuvre-lés-Nancy, France) for his technical support.

411

412 **References**

- 413 [1] E. Pennisi, Billions of North American birds have vanished, Science. 365 (2019) 1228–1229.
 414 https://doi.org/10.1126/science.365.6459.1228.
- 415 [2] Y. Si, Y. Xu, F. Xu, X. Li, W. Zhang, B. Wielstra, J. Wei, G. Liu, H. Luo, J. Takekawa, S.

416 Balachandran, T. Zhang, W.F. de Boer, H.H.T. Prins, P. Gong, Spring migration patterns,

- 417 habitat use, and stopover site protection status for two declining waterfowl species wintering
- 418 in China as revealed by satellite tracking, Ecol. Evol. 8 (2018) 6280–6289.
- 419 https://doi.org/10.1002/ece3.4174.
- 420 [3] A.M. Calvert, C.A. Bishop, R.D. Elliot, E.A. Krebs, T.M. Kydd, C.S. Machtans, G.J.
- 421 Robertson, A Synthesis of Human-related Avian Mortality in Canada, Avian Conserv. Ecol. 8

422 (2013) art11. https://doi.org/10.5751/ACE-00581-080211.

- P. Mineau, M. Whiteside, Pesticide Acute Toxicity Is a Better Correlate of U.S. Grassland
 Bird Declines than Agricultural Intensification, PLoS One. 8 (2013) e57457.
- 425 https://doi.org/10.1371/journal.pone.0057457.
- 426 [5] R.L. Stanton, C.A. Morrissey, R.G. Clark, Analysis of trends and agricultural drivers of
- farmland bird declines in North America: A review, Agric. Ecosyst. Environ. 254 (2018) 244–
 254. https://doi.org/10.1016/j.agee.2017.11.028.
- 429 [6] K. V. Rosenberg, A.M. Dokter, P.J. Blancher, J.R. Sauer, A.C. Smith, P.A. Smith, J.C.
- 430 Stanton, A. Panjabi, L. Helft, M. Parr, P.P. Marra, Decline of the North American avifauna,

431 Science. 366 (2019) 120–124. https://doi.org/10.1126/science.aaw1313.

- F. Hildmann, C. Gottert, T. Frenzel, G. Kempe, K. Speer, Pesticide residues in chicken eggs –
 A sample preparation methodology for analysis by gas and liquid chromatography/tandem
 mass spectrometry, J. Chromatogr. A. 1403 (2015) 1–20.
- 435 https://doi.org/10.1016/j.chroma.2015.05.024.
- B. Holmes, A. Dunkin, R. Schoen, C. Wiseman, Single-Laboratory Ruggedness Testing and
 Validation of a Modified QuEChERS Approach To Quantify 185 Pesticide Residues in
 Salmon by Liquid Chromatography– and Gas Chromatography–Tandem Mass Spectrometry,

439 J. Agric. Food Chem. 63 (2015) 5100–5106. https://doi.org/10.1021/jf5055276.

- P. Kaczyński, B. Łozowicka, M. Perkowski, J. Szabuńko, Multiclass pesticide residue analysis
 in fish muscle and liver on one-step extraction-cleanup strategy coupled with liquid
 chromatography tandem mass spectrometry, Ecotoxicol. Environ. Saf. 138 (2017) 179–189.
 https://doi.org/10.1016/j.ecoenv.2016.12.040.
- 444 [10] P. Kaczyński, Clean-up and matrix effect in LC-MS/MS analysis of food of plant origin for

- 445 high polar herbicides, Food Chem. 230 (2017) 524–531.
- 446 https://doi.org/10.1016/j.foodchem.2017.03.091.
- [11] R. Romero-González, A. Garrido Frenich, Application of HRMS in Pesticide Residue
 Analysis in Food From Animal Origin, in: Appl. High Resolut. Mass Spectrom., Elsevier,
- 449 2017: pp. 203–232. https://doi.org/10.1016/B978-0-12-809464-8.00007-5.
- [12] E. Bro, J. Devillers, F. Millot, A. Decors, Residues of plant protection products in grey
 partridge eggs in French cereal ecosystems, Environ. Sci. Pollut. Res. 23 (2016) 9559–9573.
 https://doi.org/10.1007/s11356-016-6093-7.
- [13] K. Lichtmannegger, R. Fischer, F.X. ave. Steemann, H. Unterluggauer, S. Masselter,
 Alternative QuEChERS-based modular approach for pesticide residue analysis in food of
 animal origin, Anal. Bioanal. Chem. 407 (2015) 3727–3742. https://doi.org/10.1007/s00216015-8594-5.
- [14] S. Choi, S. Kim, J.Y. Shin, M. Kim, J.-H. Kim, Development and verification for analysis of
 pesticides in eggs and egg products using QuEChERS and LC–MS/MS, Food Chem. 173
 (2015) 1236–1242. https://doi.org/10.1016/j.foodchem.2014.10.143.
- 460 [15] Y. Li, Z. Chen, R. Zhang, P. Luo, Y. Zhou, S. Wen, M. Ma, Simultaneous Determination of
- 461 42 Pesticides and Herbicides in Chicken Eggs by UHPLC–MS/MS and GC–MS Using a
- 462 QuEChERS-Based Procedure, Chromatographia. 79 (2016) 1165–1175.
- 463 https://doi.org/10.1007/s10337-016-3132-y.
- [16] X. Xu, X. Xu, M. Han, S. Qiu, X. Hou, Development of a modified QuEChERS method based
 on magnetic multiwalled carbon nanotubes for the simultaneous determination of veterinary
 drugs, pesticides and mycotoxins in eggs by UPLC-MS/MS, Food Chem. 276 (2019) 419–
- 467 426. https://doi.org/10.1016/j.foodchem.2018.10.051.

468	[17]	NE. Song, J.Y. Lee, A.R. Mansur, H.W. Jang, MC. Lim, Y. Lee, M. Yoo, T.G. Nam,
469		Determination of 60 pesticides in hen eggs using the QuEChERS procedure followed by LC-
470		MS/MS and GC-MS/MS, Food Chem. 298 (2019) 125050.
471		https://doi.org/10.1016/j.foodchem.2019.125050.
472	[18]	M. Anastassiades, S.J. Lehotay, D. Štajnbaher, F.J. Schenck, Fast and Easy Multiresidue
473		Method Employing Acetonitrile Extraction/Partitioning and "Dispersive Solid-Phase
474		Extraction" for the Determination of Pesticide Residues in Produce, J. AOAC Int. 86 (2003)
475		412-431. https://doi.org/10.1093/jaoac/86.2.412.
476	[19]	P. Luo, X. Liu, F. Kong, L. Tang, Q. Wang, W. Li, W. Xu, S. Wen, L. Chen, Y. Li, Multi-
477		residue determination of 325 pesticides in chicken eggs with EMR-Lipid clean-up by
478		UHPLC-MS/MS and GC-MS/MS, Chromatographia. 83 (2020) 593-599.
479		https://doi.org/10.1007/s10337-020-03876-1.

- [20] K.A. Lewis, J. Tzilivakis, D.J. Warner, A. Green, An international database for pesticide risk
 assessments and management, Hum. Ecol. Risk Assess. An Int. J. 22 (2016) 1050–1064.
 https://doi.org/10.1080/10807039.2015.1133242.
- 483 [21] Z. Tousova, P. Oswald, J. Slobodnik, L. Blaha, M. Muz, M. Hu, W. Brack, M. Krauss, C. Di
- 484 Paolo, Z. Tarcai, T.-B. Seiler, H. Hollert, S. Koprivica, M. Ahel, J.E. Schollée, J. Hollender,
- 485 M.J.F. Suter, A.O. Hidasi, K. Schirmer, M. Sonavane, S. Ait-Aissa, N. Creusot, F. Brion, J.
- 486 Froment, A.C. Almeida, K. Thomas, K.E. Tollefsen, S. Tufi, X. Ouyang, P. Leonards, M.
- 487 Lamoree, V.O. Torrens, A. Kolkman, M. Schriks, P. Spirhanzlova, A. Tindall, T. Schulze,
- 488 European demonstration program on the effect-based and chemical identification and
- 489 monitoring of organic pollutants in European surface waters, Sci. Total Environ. 601–602
- 490 (2017) 1849–1868. https://doi.org/10.1016/j.scitotenv.2017.06.032.
- 491 [22] X. Zhang, Y. Song, Q. Jia, L. Zhang, W. Zhang, P. Mu, Y. Jia, Y. Qian, J. Qiu, Simultaneous

determination of 58 pesticides and relevant metabolites in eggs with a multi-functional filter
by ultra-high performance liquid chromatography-tandem mass spectrometry, J. Chromatogr.
A. 1593 (2019) 81–90. https://doi.org/10.1016/j.chroma.2019.01.074.

- 495 [23] Arrêté du 7 juillet 2006 portant sur l'introduction dans le milieu naturel de grand gibier ou de
 496 lapins et sur le prélèvement dans le milieu naturel d'animaux vivants d'espèces dont la chasse
 497 est autorisée (NOR: DEVP0650375A, Art. L424-11), Code l'Environnement. July 7, 2006.
- 498 [24] B. Giroud, A. Vauchez, E. Vulliet, L. Wiest, A. Buleté, Trace level determination of

499 pyrethroid and neonicotinoid insecticides in beebread using acetonitrile-based extraction

500 followed by analysis with ultra-high-performance liquid chromatography-tandem mass

- 501 spectrometry, J. Chromatogr. A. 1316 (2013) 53–61.
- 502 https://doi.org/10.1016/j.chroma.2013.09.088.
- European Commission DG-SANTE, Analytical Quality Control and Method Validation
 Procedures for Pesticide Residues Analysis in Food and Feed, N° SANTE/12682/2019, (2019)
 49.
- A. Lazartigues, L. Wiest, R. Baudot, M. Thomas, C. Feidt, C. Cren-Olivé, Multiresidue
 method to quantify pesticides in fish muscle by QuEChERS-based extraction and LC-MS/MS,
 Anal. Bioanal. Chem. 400 (2011) 2185–2193. https://doi.org/10.1007/s00216-011-4945-z.
- 509 [27] X. Dong, S. Liang, Z. Shi, H. Sun, Development of multi-residue analysis of herbicides in
 510 cereal grain by ultra-performance liquid chromatography–electrospray ionization-mass
 511 spectrometry, Food Chem. 192 (2016) 432–440.
- 512 https://doi.org/10.1016/j.foodchem.2015.07.025.

513 [28] M. Anastassiades, E. Scherbaum, B. Tadelen, D. Štajnbaher, Recent Developments in

514 QuEChERS Methodology for Pesticide Multiresidue Analysis, in: H. Ohkawa, H. Miyagawa,

515 P.W. Lee (Eds.), Pestic. Chem. Crop Prot. Public Heal. Environ. Saf., WILEY-VCH,

516		Cambridge, 2007: pp. 1–30. https://doi.org/10.1017/CBO9781107415324.004.
517	[29]	S.J. Lehotay, M. O'Neil, J. Tully, A.V. García, M. Contreras, H. Mol, V. Heinke, T. Anspach,
518		G. Lach, R. Fussell, K. Mastovska, M.E. Poulsen, A. Brown, W. Hammack, J.M. Cook, L.
519		Alder, K. Lindtner, M.G. Vila, M. Hopper, A. de Kok, M. Hiemstra, F. Schenck, A. Williams,
520		A. Parker, Determination of Pesticide Residues in Foods by Acetonitrile Extraction and
521		Partitioning with Magnesium Sulfate: Collaborative Study, J. AOAC Int. 90 (2007) 485–520.
522		https://doi.org/10.1093/jaoac/90.2.485.
523	[30]	NF EN 15662:2018, Aliments d'origine végétale - Multiméthode de détermination des résidus
524		de pesticides par analyse CG et CL après extraction/partition avec de l'acétonitrile et
525		purification par SPE dispersive - Méthode modulaire QuEChERS, France, 2018.
526	[31]	O. Shimelis, Y. Yang, K. Stenerson, T. Kaneko, M. Ye, Evaluation of a solid-phase extraction
527		dual-layer carbon/primary secondary amine for clean-up of fatty acid matrix components from
528		food extracts in multiresidue pesticide analysis, J. Chromatogr. A. 1165 (2007) 18-25.
529		https://doi.org/10.1016/j.chroma.2007.07.037.
530	[32]	J. V. Dias, V. Cutillas, A. Lozano, I.R. Pizzutti, A.R. Fernández-Alba, Determination of
531		pesticides in edible oils by liquid chromatography-tandem mass spectrometry employing new
532		generation materials for dispersive solid phase extraction clean-up, J. Chromatogr. A. 1462
533		(2016) 8–18. https://doi.org/10.1016/j.chroma.2016.07.072.
534	[33]	M. Chys, K. Demeestere, A.S. Ingabire, J. Dries, H. Van Langenhove, S.W.H. Van Hulle,
535		Enhanced treatment of secondary municipal wastewater effluent: comparing (biological)
536		filtration and ozonation in view of micropollutant removal, unselective effluent toxicity, and
537		the potential for real-time control, Water Sci. Technol. (2017) wst2017207.
538		https://doi.org/10.2166/wst.2017.207.


539 [34] M.V. Barbieri, C. Postigo, N. Guillem-Argiles, L.S. Monllor-Alcaraz, J.I. Simionato, E. Stella,

540	D. Barceló, M. López de Alda, Analysis of 52 pesticides in fresh fish muscle by QuEChERS
541	extraction followed by LC-MS/MS determination, Sci. Total Environ. 653 (2019) 958–967.
542	https://doi.org/10.1016/j.scitotenv.2018.10.289.

- [35] M. Burkhardt, T. Kupper, S. Hean, R. Haag, P. Schmid, M. Kohler, M. Boller, Biocides used
 in building materials and their leaching behavior to sewer systems, Water Sci. Technol. 56
- 545 (2007) 63–67. https://doi.org/10.2166/wst.2007.807.
- 546 [36] U. Schoknecht, J. Gruycheva, H. Mathies, H. Bergmann, M. Burkhardt, Leaching of Biocides
 547 Used in Façade Coatings under Laboratory Test Conditions, Environ. Sci. Technol. 43 (2009)
 548 9321–9328. https://doi.org/10.1021/es9019832.
- [37] U.E. Bollmann, C. Tang, E. Eriksson, K. Jönsson, J. Vollertsen, K. Bester, Biocides in urban
 wastewater treatment plant influent at dry and wet weather: Concentrations, mass flows and
 possible sources, Water Res. 60 (2014) 64–74. https://doi.org/10.1016/j.watres.2014.04.014.

552

554 Figure and table legends



556 Fig. 1. Detailed steps of the extractions Ref#1, Ref#2 and the final protocol, with each of the

⁵⁵⁷ modified steps in grey.






Fig. 2. Evaluation of global LOQ in ng.g⁻¹ (A) and Extraction Recoveries (ER) in % (B) for each tested parameter (n=3), with Ref#1: the first protocol tested, and Ref#:2 the protocol built from general parameters' optimization.







		quantinoution			ER	\pm %RSD ((%)	Inter-da	y precision	(%RSD)	Blanks	
Chemicals	Range (ng.g ⁻¹)	r ² ± ‰RSD	LOQ (ng.g ⁻¹)	ME ± SD (%)	LOQ	(n=3) 3xLOQ	7xLOQ	LOQ	(n=3) 3xLOQ	7xLOQ	Max - Min (ng)	QF
Fungicides												(%)
Boscalid	5 - 100	0.995 (5)	5.0	-30 ± 42	82 ± 22	78 ± 19	85 ± 10	6	6	3	_	-
Tebuconazole	0.5 - 10	0.994 (4)	0.5	-10 ± 17	75 ± 15	75 ± 18	83 ± 12	8	5	3	0.1	17
Herbicides												
Acetochlor	10 - 200	0.994 (7)	25.0	-35 ± 30	67 ± 28	66 ± 30	76 ± 9	13	12	2	-	-
Atrazine	0.5 - 10	0.995 (4)	0.5	-39 ± 30	76 ± 25	71 ± 26	84 ± 7	15	5	2	-	-
Bentazone	0.5 - 10	0.992 (3)	0.5	-26 ± 49	74 ± 38	68 ± 20	78 ± 13	16	5	6	-	-
Benzamide	25 - 500	0.995 (5)	25.0	-35 ± 41	79 ± 20	78 ± 23	87 ± 9	7	7	3	5.6	17
Chlorotoluron	1 - 20	0.996 (3)	5.0	-32 ± 37	82 ± 17	77 ± 21	86 ± 10	8	6	4	-	-
Dimethachlor	5 - 100	0.994 (6)	5.0	-24 ± 10	74 ± 19	74 ± 20	69 ± 27	12	8	2	-	-
Dimethenamid	5 - 100	0.993 (9)	5.0	-35 ± 8	79 ± 17	74 ± 21	71 ± 32	8	6	2	-	-
Flufenacet	0.5 - 10	0.992 (4)	2.5	-26 ± 25	68 ± 22	69 ± 21	69 ± 32	14	9	3	-	-
Isoproturon	0.5 - 10	0.994 (8)	0.5	-42 ± 46	80 ± 14	77 ± 18	91 ± 10	7	7	1	-	-
Metazachlor	0.5 - 10	0.993 (4)	2.5	-39 ± 31	87 ± 32	81 ± 18	75 ± 30	10	7	5	0.4	17
Metsulfuron-methyl	5 - 100	0.99 (5)	25.0	-31 ± 33	68 ± 29	80 ± 24	73 ± 37	22	11	6	-	-
Prosulfocarb	5 - 100	0.991 (11)	5.0	-66 ± 27	39 ± 23	38 ± 24	64 ± 34	9	12	4	4 - 1	50
S-metolachlor	1 - 20	0.993 (7)	1.0	-29 ± 26	80 ± 27	69 ± 27	75 ± 15	14	10	4	0.2	50
Terbuthylazine	0.5 - 10	0.996 (4)	0.5	-37 ± 34	76 ± 20	72 ± 22	80 ± 10	10	8	2	-	-
Terbutryn	0.5 - 10	0.995 (4)	0.5	-43 ± 27	72 ± 20	69 ± 27	78 ± 8	8	9	3	0.4 - 0.1	67
Insecticides												
Fenthion	5 - 100	0.986 (9)	25.0	-36 ± 29	71 ± 25	65 ± 30	76 ± 9	7	11	8	-	-
Fipronil	0.5 - 10	0.992 (8)	0.5	-6 ± 13	75 ± 28	79 ± 23	83 ± 12	11	9	5	0.4 - 0.1	33
Imidacloprid	5 - 100	0.991 (7)	25.0	-29 ± 38	87 ± 27	77 ± 28	92 ± 14	17	10	4	-	-

Table 1. Analytical performances of the validated method with n=6 unless specified (Range: range of tested concentrations (ng; LOQ: calculated limit of quantification; ME:
 matrix effect; ER: extraction recovery; QF: quantification frequency)

Transformation products												
Acetochlor OXA	5 - 100	0.985 (10)	25.0	-43 ± 33	60 ± 33	51 ± 33	62 ± 14	13	10	7	-	-
Alachlor-Acetochlor ESA	2 - 40	0.985 (7)	10.0	-32 ± 13	51 ± 33	53 ± 23	57 ± 7	14	7	7	-	-
Atrazine-2-hydroxy	1 - 20	0.989 (9)	1.0	-9 ± 27	93 ± 26	78 ± 30	77 ± 7	7	16	6	-	-
CGA 357704	5 - 100	0.987 (7)	5.0	-4 ± 25	49 ± 34	50 ± 25	57 ± 9	21	11	7	1.2	17
CGA 368208	1 - 20	0.978 (15)	1.0	-13 ± 70	61 ± 36	54 ± 22	62 ± 14	24	9	13	-	-
CGA 50267	0.5 - 10	0.995 (5)	0.5	-17 ± 21	84 ± 22	73 ± 20	78 ± 12	9	9	4	-	-
DEA-OH	5 - 100	0.995 (3)	5.0	-50 ± 26	77 ± 20	72 ± 23	69 ± 33	3	5	3	3.5 - 1.8	33
Desmethyl-isoproturon	1 - 20	0.991 (8)	5.0	-31 ± 14	80 ± 23	79 ± 22	80 ± 23	20	8	3	-	-
Dimethachlor ESA	5 - 100	0.992 (5)	25.0	-30 ± 14	56 ± 30	51 ± 24	65 ± 11	13	10	5	-	-
Dimethenamid ESA	5 - 100	0.986 (5)	25.0	-60 ± 32	54 ± 35	47 ± 23	72 ± 15	10	10	7	-	-
Fipronil sulfone	0.5 - 10	0.995 (5)	0.5	-20 ± 38	82 ± 24	77 ± 20	78 ± 23	13	8	1	0.3 - 0.1	83
Flufenacet ESA	5 - 100	0.99 (6)	5.0	-43 ± 29	53 ± 19	48 ± 22	70 ± 17	15	9	5	5.1	17
Flufenacet OXA	10 - 200	0.988 (8)	10	-40 ± 17	52 ± 35	48 ± 25	60 ± 16	18	13	6	-	-
Imidacloprid desnitro	0.5 - 10	0.994 (6)	0.5	-26 ± 49	74 ± 28	78 ± 24	87 ± 11	3	7	3	-	-
MCPA	5 - 100	0.991 (5)	25.0	-36 ± 35	65 ± 43	66 ± 20	79 ± 8	15	8	6	-	-
Me-DPC	5 - 100	0.995 (6)	5.0	-3 ± 32	84 ± 22	74 ± 23	83 ± 18	9	6	2	-	-
Metazachlor ESA	5 - 100	0.989 (7)	5.0	125 ± 359	55 ± 29	52 ± 23	73 ± 17	20	10	4	-	-
Metazachlor OXA	25 - 500	0.994 (7)	25.0	-26 ± 32	55 ± 22	53 ± 24	78 ± 12	5	10	3	5.3	17
Metolachlor ESA	5 - 100	0.99 (9)	5.0	-52 ± 23	53 ± 21	49 ± 21	61 ± 7	6	8	6	-	-
Metolachlor OXA	5 - 100	0.971 (28)	25.0	-46 ± 28	41 ± 26	47 ± 28	49 ± 15	20	22	10	8.9	17
Saccharin	5 - 100	0.99 (5)	25.0	-48 ± 24	55 ± 34	61 ± 23	57 ± 39	17	10	6	-	-
ТСР	10 - 200	0.995 (4)	10.0	0 ± 31	68 ± 24	70 ± 21	50 ± 82	9	9	2	-	-
Terbuthylazine desethyl	0.5 - 10	0.993 (7)	0.5	-29 ± 31	86 ± 24	74 ± 20	86 ± 9	9	6	3	-	-
Terbuthylazine-desethyl-OH	5 - 100	0.991 (7)	5.0	-28 ± 30	76 ± 21	71 ± 24	83 ± 8	7	11	4	1.6	17

Compound	CAS	Uses	Log P		Transition	t _R	DP	EP	CE	СХР	Ion rat
Acetochlor ESA + Alachlor ESA*	947601-84-5	Acetochlor and alachlor TP	2.2	MRM 1	314.0→80.0	4.9	-135	-10	-52	-15	2.1
				MRM 2	314.0→120.9	4.9	-135	-10	-30	-13	
Acetochlor OXA*	194992-44-4	Acetochlor TP	2.9	MRM 1	264.1→144.0	4.7	-55	-10	-38	-13	0.4
				MRM 2	264.1→146.0	4.7	-65	-10	-16	-19	
Bentazon*	25057-89-0	Herbicide	-0.5	MRM 1	239.0→132.0	4	-100	-10	-38	-4	2.0
				MRM 2	239.0→175.0	4	-100	-10	-28	-6	
CGA 357704*	1217465-10-5	S-metolachlor TP	1.9	MRM 1	278.0→206.0	2.7	-45	-10	-18	-15	22.8
				MRM 2	278.0→162.0	2.7	-45	-10	-26	-15	
CGA 368208*	1173021-76-5	S-metolachlor TP	1.9	MRM 1	256.0→79.8	4.2	-145	-10	-58	-11	2.0
				MRM 2	256.0→120.8	4.2	-145	-10	-28	-11	
Dimethachlor ESA*	1231710-75-0	Dimethachlor TP	-0.2	MRM 1	300.0→80.0	4.4	-25	-10	-66	-13	1.9
				MRM 2	300.0→120.9	4.4	-25	-10	-30	-15	
Dimethenamid ESA*	205939-58-8	Dimethenamid TP	-0.6	MRM 1	319.9→120.7	4.6	-50	-10	-30	-17	1.5
				MRM 2	319.9→80.0	4.6	-50	-10	-66	-9	
Fipronil sulfone*	120068-36-2	Fipronil TP	4.6	MRM 1	450.8→414.9	8.5	-65	-10	-24	-21	1.2
				MRM 2	450.8→281.9	8.5	-65	-10	-38	-19	
Fipronil*	120068-37-3	Insecticide	3.8	MRM 1	434.9→329.9	8.2	-110	-10	-22	-21	3.1
				MRM 2	434.9→249.9	8.2	-110	-10	-38	-17	
Flufenacet ESA*	947601-87-8	Flufenacet TP	1.2	MRM 1	274.1→121.0	4.5	-25	-10	-28	-11	1.1
				MRM 2	274.1→79.9	4.5	-25	-10	-54	-29	
Flufenacet OXA*	201668-31-7	Flufenacet TP	1.6	MRM 1	224.0→152.0	4.2	-20	-10	-15	-10	1.2
				MRM 2	152.0→136.0	4.2	-90	-8	-17	-10	
MCPA*	94-74-6	Herbicide	-0.8	MRM 1	199.0→141.0	4.5	-65	-10	-20	-12	2.5
				MRM 2	201.0→143.0	4.5	-65	-10	-20	-12	
Metazachlor ESA*	172960-62-2	Metazachlor TP	-0.1	MRM 1	322.0→120.9	4.4	-110	-10	-30	-11	2.2
				MRM 2	322.0→79.9	4.4	-110	-10	-80	-11	
Metolachlor ESA*	947601-85-6	S-metolachlor TP	2.1	MRM 1	328.1→79.8	4.9	-135	-10	-72	-11	2.3
				MRM 2	328.1→121.0	4.9	-135	-10	-30	-13	

569 **Table S1.** Analytes and analytical parameters associated, (t_R: retention time; DP: declustering potential; EP: entrance potential; CE: collision energy; CXP: collision cell exit 570 potential; TP: transformation product; IS: isotope labeled standard, MRM1: used for quantification; MRM2: used for confirmation)

Compound	CAS	Uses	Log P	Transition		t _R	DP	EP	CE	СХР	lon ratio
Metolachlor OXA*	152019-73-3	S-metolachlor TP	2.3	MRM 1	278.1→174.1	4.8	-65	-10	-22	-15	7.6
				MRM 2	278.1→206.1	4.8	-65	-10	-14	-13	
Saccharin*	81-07-2	Metsulfuron-methyl TP	0.9	MRM 1	181.9→41.9	2	-150	-10	-88	-19	9.7
				MRM 2	181.9→61.9	2	-115	-10	-26	-9	
13C4-Fipronil*	2140327-54-2	IS	-	MRM 1	438.9→334.0	8.2	-50	-10	-22	-17	-
D5-acetochlor ESA*	-	IS	-	MRM 1	319.0→79.9	4.9	-135	-10	-66	-11	-
D6-Bentazone*	-	IS	-	MRM 1	245.0→131.9	4	-40	-10	-36	-15	-
D6-Metazachlor ESA*	-	IS	-	MRM 1	328.2→120.9	4.4	-115	-10	-32	-13	-
Acetochlor	34256-82-1	Herbicide	4.1	MRM 1	270.1→223.9	8.1	61	10	13	16	1.6
				MRM 2	270.1→148.1	8.1	61	10	27	10	
Atrazine	1912-24-9	Herbicide	2.7	MRM 1	216.1→174.0	6.8	86	4.5	23	4	3.1
				MRM 2	216.1→104.1	6.8	86	4.5	39	4	
Atrazine-2-hydroxy	2163-68-0	Atrazine TP	2.1	MRM 1	198.1→156.0	4.1	71	10	25	12	2.6
				MRM 2	198.1→114.0	4.1	71	10	31	20	
Benzamide	55-21-0	Herbicide	0.6	MRM 1	122.0→77.0	3.8	95	10	33	13	1.6
				MRM 2	122.0→79.0	3.8	95	10	19	13	
Boscalid	188425-85-6	Fungicide	3	MRM 1	343.0→307.0	7.7	111	10	29	20	2.3
				MRM 2	343.0→140.0	7.7	111	10	29	12	
CGA 50267	82508-02-9	S-metolachlor TP	2.4	MRM 1	208.0→135.0	6.4	61	10	22	13	1.1
				MRM 2	208.0→120.0	6.4	61	10	41	13	
Chlorotoluron	15545-48-9	Herbicide	2.5	MRM 1	213.1→72.0	6.6	74	6	33	4	2.6
				MRM 2	213.1→46.2	6.6	74	6	27	4	
DEA-OH	19988-24-0	Atrazine TP	-1	MRM 1	170.0→128.0	1.5	96	10	21	12	1.3
				MRM 2	170.0→86.0	1.5	96	10	29	14	
desmethyl isoproturon	34123-57-4	Isoproturon TP	2.8	MRM 1	193.0→94.0	6.4	40	10	32	13	2.6
				MRM 2	193.0→136.0	6.4	40	10	22	13	
Dimethachlor	50563-36-5	Herbicide	2.2	MRM 1	256.0→224.2	7.2	64	10.5	19	6.5	2
				MRM 2	256.0→148.0	7.2	64	10	33	6	
Dimethenamid	87674-68-8	Herbicide	2.2	MRM 1	275.9→244.1	7.6	16	10	19	16	3.3
				MRM 2	275.9→168.0	7.6	16	10	35	20	

Compound	CAS	Uses	Log P	Log P Transition		t _R	DP	EP	CE	СХР	lon ratio
Fenthion	55-38-9	Insecticide	4.8	MRM 1	278.9→169.1	8.3	71	10	23	12	1
				MRM 2	278.9→246.9	8.3	71	10	17	18	
Flufenacet	142459-58-3	Herbicide	3.5	MRM 1	364.0→152.0	8	91	10	27	14	1.1
				MRM 2	364.0→194.1	8	91	10	15	16	
Imidacloprid desnitro	127202-53-3	Imidacloprid TP	-0.8	MRM 1	211.0→126.0	2.6	96	10	33	12	3.7
				MRM 2	211.0→90.0	2.6	96	10	49	8	
Imidacloprid	138261-41-3	Insecticide	0.6	MRM 1	255.9→175.1	5	66	10	29	14	1.0
				MRM 2	255.9→209.1	5	66	10	23	14	
Isoproturon	34123-59-6	Herbicide	2.5	MRM 1	207.2→72.1	6.8	86	4	29	4	5.5
				MRM 2	207.2→165.2	6.8	86	10.5	19	6.5	
Me-DPC	64-55-1	Chloridazon TP	-1.4	MRM 1	160.0→87.9	2.2	61	10	39	10	2.6
				MRM 2	160.0→130.0	2.2	61	10	33	12	
Metazachlor	67129-08-2	Herbicide	2.5	MRM 1	278.1→210.1	7.1	58	3.5	15	4	0.9
				MRM 2	278.1→134.2	7.1	58	3.5	27	4	
Metazachlor OXA	1231244-60-2	Herbicide	1.5	MRM 1	274.0→69.0	4.9	56	10	31	10	1.3
				MRM 2	274.0→134.1	4.9	56	10	29	18	
Metsulfuron-methyl	74223-64-6	Herbicide	-1.9	MRM 1	382.0→167.0	6.2	86	10	23	14	1.4
				MRM 2	382.0→141.0	6.2	86	10	21	10	
Prosulfocarb	52888-80-9	Fungicide	4.5	MRM 1	252.3→91.1	8.9	71	10	35	14	3.5
				MRM 2	252.3→128.1	8.9	71	10	17	14	
S-metolachlor	87392-12-9	Herbicide	3	MRM 1	284.1→252.1	8.1	101	10	21	18	2.8
				MRM 2	284.1→176.2	8.1	101	10	35	14	
ТСР	6515-38-4	Chlorpyrifos-ethyl TP	2.2	MRM 1	197.9→107.0	6.8	56	10	43	8	1.9
				MRM 2	197.9→134.0	6.8	56	10	31	10	
Tebuconazole	107534-96-3	Fungicide	3.7	MRM 1	308.1→70.1	7.8	61	10	39	4	8.7
				MRM 2	308.1→125.0	7.8	61	10	51	4	
Terbuthylazine desethyl	30125-63-4	Terbuthylazine TP	2.3	MRM 1	202.0→146.0	6.2	46	10	23	16	5.9
				MRM 2	202.0→103.9	6.2	46	10	37	12	
Terbuthylazine -desethyl-2-OH	66753-06-8	Terbuthylazine TP	-1.2	MRM 1	184.0→128.1	2.2	81	10	21	10	2.3
				MRM 2	184.0→86.1	2.2	81	10	35	10	

Compound	CAS	Uses	Log P		Transition	t _R	DP	EP	CE	СХР	lon ratio
Terbuthylazine	5915-41-3	Herbicide	3.4	MRM 1	230.1→174.1	7.4	81	11.5	23	6.5	6.6
				MRM 2	230.1→104.1	7.4	86	7.5	41	4	
Terbutryn	886-50-0	Herbicide and terbuthylazine TP	3.7	MRM 1	242.0→186.1	6.8	41	10	27	14	5.2
				MRM 2	242.0→68.0	6.8	41	10	61	8	
13C-DEDIA	1216850-33-7	IS	-	MRM 1	149.0→105.8	1.7	16	10	27	14	-
D10-Simazine	220621-39-6	IS	-	MRM 1	212.0→137.1	6.1	66	10	29	10	-
D11-Acetochlor	1189897-44-6	IS	-	MRM 1	281.1→235.0	8.1	71	10	15	14	-
D3-Isoproturon	352438-80-3	IS	-	MRM 1	210.1→75.1	6.8	86	10	25	10	-
D4-carbendazim	291765-95-2	IS	-	MRM 1	196.1→164.2	3.2	60	10	25	12	-
D4-Imidacloprid	1015855-75-0	IS	-	MRM 1	260.0→213.1	5	61	10	23	14	-
D5-Atrazine	163165-75-1	IS	-	MRM 1	221.1→179.1	6.8	121	10	25	16	-
D5-Atrazine-2-hydroxy	2163-68-0	IS	-	MRM 1	203.1→161.4	4.1	126	10	23	10	-
D5-DIA	1189961-78-1	IS	-	MRM 1	179.0→100.9	4.3	101	10	27	10	-
D5-Metazachlor	1246816-51-2	IS	-	MRM 1	284.1→216.1	7.1	58	3.5	15	4	-
D5-Terbutryn	1219804-47-3	IS	-	MRM 1	247.4→190.9	6.8	96	10	23	14	-
D6-chlorpyrifos Me	2083629-84-7	IS	-	MRM 1	327.9→130.9	8.6	76	10	31	10	-
D9-Tebuconazole	1246818-83-6	IS	-	MRM 1	317.1→70.0	7.8	86	10	57	12	-

*: chemicals analyzed in negative electrospray mode, otherwise analyzes were performed in positive mode




Fig. S1. Representation of isotope labeled standards area ratio between reference and field samples, PCA based on Pearson regression ($\alpha = 0.005$).