Centre de Recherche sur les lons, les Matériaux et la Photonique





Proteomic analysis of the Radiation-Induced Bystander Effects

Lepleux C.^{1,2} and Chevalier F.^{1,2}

edf RÉGION NORMANDIE

Funded by :

Experimental evidence shows that, in addition to this direct and targeted effect of radiation, a side or secondary effect is observed within the un-irradiated area. Irradiated cells (1 and 2) relay a stress signal in the close proximity. Neighboring un-irradiated cells (3) to 7) react to this signal, i.e., BE with a specific "bystander response" (BR), characteristic of a biological context, but with a close relationship to the biological response typically associated with direct radiation exposure. Bystander responses and bystander factors (BF) secreted by irradiated cells were observed and studied in many physical and biological conditions, in vitro and in vivo.



The cells were irradiated (X-rays / Carbon ions) and the Bystander supernatant was transferred on non-irradiated cells. To study the effects of these treatments on the survival of cell lines, we carried out clonogenic assays and in order to study the impact of these treatments (X-rays) on the cellular proteome, we carried out proteomic analysis by two-dimensional electrophoresis (2DE), starting

1 LARIA, iRCM, François JACOB Institute, DRF-CEA, Caen, France; 2 Normandie Univ, ENSICAEN, UNICAEN, CEA, CNRS, CIMAP UMR6252, 14000 Caen, France

In the present study, we investigated the capacity of chondrocytes (T/C28-A2) in responding to bystander factors released by irradiated chondrosarcoma cells (SW1353) using a medium transfer protocol.

from the cellular protein lysates.

The 2DE analysis showed the effect of irradiation on the proteome of irradiated chondrosarcoma cells as well as the Bystander effect on the proteome of bystander chondrocytes. We then selected variant spots that were analyzed by mass spectrometry.

An in-depth analysis of the identified proteins will provide a better understanding of the mechanisms involved in cell irradiation as well as the radio-induced Bystander effect.



Secreted proteins : potential bystander factors

Conditioned medium (SW1353 0,1Gy) vs conditioned medium (SW1353 0Gy)



Q9NZN4-01	EH domain-containing protein 2	1,50	
P54577-01	TyrosinetRNA ligase, cytoplasmic	-1,64	
Q9UBM7-01	7-dehydrocholesterol reductase	-1,67	
Q14152-03	Eukaryotic translation initiation factor 3 subunit A	-1,80	
P13861	cAMP-dependent protein kinase type II-alpha regulatory subur	it -2,00	MAPK pathway regulation
P49189-01	4-trimethylaminobutyraldehyde dehydrogenase	-2,67	Mitochondrial permeability
Ductoin outvoction and concretion	green : increased in bystander ; red : decreased in by	stander	
Proteins were TCA precipitated and the Mass spectrometry Each bands from the SDS-PAGE gels w	en solubilized in laemmli buffer, and separated by SDS-PAGE.		_

Modulated proteins in bystander cells : potential bystander responses

Bystander T/C28-A2 (0,1Gy) vs control T/C28-A2 (0Gy)





Modulate protein folding and traffic within the cell, can be secreted, and is involved in several human pathologies Enzyme with several functions, including control growth and hypoxia (interact

with MBP-1)



Oxidative stress response

Protein extraction and solubilisation

Proteins were extracted from TC-28/Ac cells in sample buffer containing 7 M urea, 2M thiourea, 4% CHAPS, 0.05% Triton X100, 65 mM DTT, 40mM spermine, protease and phosphatase inhibitor cocktails. Proteins were then precipitated using the 2D clean-up kit (GE Healthcare) and the pellet was solubilized with TUC solution (7M urea, 2M thiourea, 4% CHAPS) and quantified with the 2D quant kit (GE Healthcare).

Protein load : 250 μg and 2D electrophoresis (6 technical / biological replicates)

IPG strips : pI 4-7 ; 1st dimension : Protean IEF Cell (Bio-Rad) ; 2nde dimension : 12% acrylamide gel (dodeca Cell electrophoresis unit, Bio-Rad) Gel staining and picture acquisition : Silver stain, stained gels scanned on a GS 800 densitometer (Bio-Rad)

2D gel analysis and Mass spectrometry

Stained gels were analysed with the Progenesis Samespots software and variant spots were picked and analysed by LC-MS/MS (Q-EX) mass spectrometry

123	1,3127831	6,56	136	P49368	T-complex protein 1 subunit gamma OS=Homo sapiens GN=CCT3 PE=1 SV=4 - [TCPG_HUMAN]	
126	1,6276638	4,67	56	P09493	Tropomyosin alpha-1 chain OS=Homo sapiens GN=TPM1 PE=1 SV=2 - [TPM1_HUMAN]	
				P31946	14-3-3 protein beta/alpha OS=Homo sapiens GN=YWHAB PE=1 SV=3 - [1433B_HUMAN]	
127	1,3070507	4,91	87	P60709	Actin, cytoplasmic 1 OS=Homo sapiens GN=ACTB PE=1 SV=1 - [ACTB_HUMAN]	
				Q16186	Proteasomal ubiquitin receptor ADRM1 OS=Homo sapiens GN=ADRM1 PE=1 SV=2 - [ADRM1_HUMAN]	
128	1,3553405	5,03	87	Q16186	Proteasomal ubiquitin receptor ADRM1 OS=Homo sapiens GN=ADRM1 PE=1 SV=2 - [ADRM1_HUMAN]	
134	1,2603369	5,68	73	Q13347	Eukaryotic translation initiation factor 3 subunit I OS=Homo sapiens GN=EIF3I PE=1 SV=1 - [EIF3I_HUMAN]	
138	1,9817232	4,18	145	P04264	Keratin, type II cytoskeletal 1 OS=Homo sapiens GN=KRT1 PE=1 SV=6 - [K2C1_HUMAN]	
145	1,4712852	4,67	98	Q9BTY7	Protein HGH1 homolog OS=Homo sapiens GN=HGH1 PE=1 SV=1 - [HGH1_HUMAN]	
147	1,4495543	6,56	144	x		
149	1,3028428	6,09	51	х		
152	1,3145383	6,72	51	x		
155	1 3942451	5.43	82	×		



Several protein were observed increased or decreased and could be related to a bystander response in cells receiving conditioned medium from irradiated cells. These preliminary results need now to be validated with functional assays and the factors will be tested on other cells.



ENLIGHT Annual Meeting and Training 2019, Caen



DRF / iRCM – CIMAP





Université de Caen Normandie