

Sonoporation: from the lab to human clinical trials

Spiros Kotopoulis, Georg Dimcevski, Bjørn Tore Gjertsen, Odd Helge Gilja,

Emmet Mc Cormack, Michiel Postema

▶ To cite this version:

Spiros Kotopoulis, Georg Dimcevski, Bjørn Tore Gjertsen, Odd Helge Gilja, Emmet Mc Cormack, et al.. Sonoporation: from the lab to human clinical trials. 2014 IEEE International Ultrasonics Symposium (IUS), Sep 2014, Chicago, United States. pp.846-849, 10.1109/ULTSYM.2014.0208 . hal-03193064

HAL Id: hal-03193064 https://hal.science/hal-03193064v1

Submitted on 11 Apr 2021

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

Sonoporation: From the lab to human clinical trials

Spiros Kotopoulis^{*†¶}, Georg Dimcevski^{*}, Bjørn Tore Gjertsen[‡], Odd Helge Gilja^{*§},

Emmet Mc Cormack[‡] and Michiel Postema[†]

*National Centre for Ultrasound in Gastroenterology, Haukeland University Hospital,

Jonas Lies vei 65, 5021 Bergen, Norway

[†]Department of Physics and Technology, University of Bergen, Allégaten 55, 5007 Bergen, Norway

[‡]Department of Clinical Science, University of Bergen, Jonas Lies vei 65, 5021 Bergen, Norway

[§]Department of Clinical Medicine, University of Bergen, Jonas Lies vei 65, 5021 Bergen, Norway

[¶]Corresponding author. Email: Spiros.Kotopoulis@uib.no

Abstract—Therapeutic ultrasound has been in use for over 70 vears but has primarily been a thermal modality. Sonoporation, 2 the use of ultrasound and stable gas microbubbles in the size 3 range of 2-10 µm to form transient pores in cell membranes, has been of great interest in the past 15 years. This technique could 5 be used to improve the delivery of current drugs in very localised 6 regions. There are several phenomena behind sonoporation 7 that all occur non-exclusively: push, pull, jetting, inertial cavitation, shear and, translation. Pre-clinical work has shown that sonoporation can be used to reduce primary tumour burden 10 and inhibit metastatic development. Our clinical trial showed that 11 ultrasound in combination with microbubbles and chemotherapy 12 can effectively double the number of chemotherapy cycles 13 patients can undergo, meaning that the patients were healthier 14 for a longer period of time. Nevertheless, sonoporation is still in 15 its infancy and there is vast room for improvement in both the 16 areas of microbubbles and ultrasound. 17

I. INTRODUCTION

18

Therapeutic ultrasound has been in use since the early 1940s 19 [1], primarily as method to heat deep tissue to promote healing. 20 Nevertheless, in the past 15 years, due to advancements 21 in technology, new ultrasound-based therapeutic modalities 22 have surfaced; specifically High Intensity Focused Ultrasound 23 (HIFU) and sonoporation. HIFU is used to heat targeted tissue 24 in very specific regions to such a high temperature that thermal 25 necrosis occurs. This is often done with the aid of Magnetic 26 Resonance Imaging to visualise the target tissue and measure 27 the thermal deposition [2]. Sonoporation on the other hand 28 aims to solve the problem of time and space controlled drug 29 delivery [3]. One of the major barriers in effective drug 30 delivery is the cell membrane; it acts as wall selectively 31 choosing what to allow in and out of the cell [4]. Sonoporation 32 can be defined as the *the transient formation of pores in a cell* 33 membrane due to ultrasound or a combination of ultrasound 34 and microbubbles. By forming pores in the cell membranes we 35 can effectively bypass this barrier that is the cell membrane, 36 resulting in highly efficient targeted drug delivery. In our work 37 we explore the phenomenon of sonoporation, from the lab 38 bench all the way to the World's first clinical trials. 39

40 II. MICROBUBBLES IN A SAFE CLINICAL SETTING

The first clinical use of microbubbles was reported in 1969 [5]. Microbubbles in the blood pool are used to improve the signal-to-noise ratio or contrast of blood when imaged

 TABLE I

 CLINICALLY APPROVED ULTRASOUND CONTRAST AGENTS.

Contrast agent	Shell	Gas	Size range Ø(µm)	Manufacturer
Optison TM	Albumin	C_3F_8	2.2–4.5	GE Healthcare
SonoVue®	Lipid	SF_6	2.5	Bracco Int.
Definity®	Lipid	C_3F_8	1.1–3.3	Lantheus Medical Imaging
Sonazoid®	Lipid	C_4F_{10}	1.9–2.4	Daiichi Pharmaceutical Co.

using B-mode sonography. Hence, these microbubbles are 44 also known as ultrasound contrast agents. Their sizes range 45 between 2-5 µm, typically with over 95% below 10 µm so 46 that the microbubbles can pass through the lung capillaries. 47 Microbubbles, like most physical objects, have a resonant 48 frequency; a frequency where its oscillation amplitude is 49 maximal. Similar to a bell, when a force strikes it, the 50 microbubble will volumetrically oscillate and generate sound 51 at its resonance frequency. For a free gas bubble this frequency 52 can be approximated by [6]: 53

$$f_r \approx \frac{6.5}{D},\tag{1}$$

64

65

66

67

where f_r is the resonant frequency and D is the diameter 54 of the bubble. From this we can see that a 2.5 µm gas 55 bubble will resonate at approximately 2.6 MHz, a typical 56 diagnostic ultrasound frequency. Free gas microbubbles are 57 inherently unstable, and dissolve almost instantly. For this 58 reason clinical diagnostic ultrasound contrast agents typically 59 have a lipid or albumin shell with a slowly diffusing gas core, 60 increasing stability, allowing for efficient clinical imaging. 61 Table I summarises clinically approved ultrasound contrast 62 agents and their contents. 63

In clinical diagnostic ultrasound imaging, the Mechanical Index (MI) gives an indication of the mechanical damage possible due to inertial cavitation when an ultrasound contrast agent is present. It is defined by:

101

105

106

107

108



Fig. 1. Photo micrographs depicting clusters of Definity $^{(\! R)}$ microbubbles attracting and merging in under 2 seconds.

$$MI = \frac{p^-}{\sqrt{f_c}},$$
(2)

where p^- is the peak-negative acoustic pressure in-situ 68 measured in MPa, and f_c is the centre frequency in MHz. 69 An MI<0.3 is considered safe, an MI between 0.3 and 0.7 70 considered moderately safe as there is a risk of damage is 71 to neonatal lung tissue or intestine, and an MI>0.7 there is 72 a high risk of inertial cavitation of the ultrasound contrast 73 agents, and a theoretical risk of cavitation formation without 74 ultrasound contrast agents present [7]. For this reason, in 75 all our work, we attempt to work at MI < 0.2, ensuring no 76 destructive and currently uncontrollable inertial cavitation is 77 present and allowing for faster transition from lab to clinic. 78

The first step in being able to induce localised sonoporation 79 was to attempt to control the location of the microbubbles [8]. 80 A custom made experimental setup incorporating a 2.2 MHz 81 ultrasound transducer, a synthetic Ø200-µm capillary and a 82 high-speed camera was built on top of an inverted microscope. 83 Continuous-wave ultrasound with centre frequencies of 2-MHz 84 and 7-MHz with peak-to-peak acoustic pressures of 20 kPa 85 were used in conjunction with Definity[®] ultrasound contrast 86 agent. Upon insonation, the microbubbles were seen to attract 87 each other forming small spherical clusters of microbubbles. 88 These clusters of 1-2 thousand microbubbles would form 89 within several seconds, and were spaced $\frac{1}{4}\lambda$ apart. As 90 the attracted microbubbles were oscillating in phase, when 91 continuing sonication, the clusters were attracted to each other, 92 forming even larger clusters (c.f. Fig. 1). When increasing 93 the acoustic pressure the clusters could be radiated towards 94 the membrane wall. These results indicated that if necessary, 95 microbubbles could be accumulated to specific regions and 96 radiated toward a vascular wall. 97

98

III. In-vitro WORK

In-vitro work has been the hallmark for evaluating cellular
 effects for years. In our work, we attempted to understand why

specific acoustic conditions were optimal for increased drug delivery.

As can be expected, injecting a gas bubble into the blood stream, that generates a detectable acoustic force, may have some inherent side-effects. In 1997, upon studying such side-effects [9] and [10] showed that ultrasound in combination with microbubbles could increase the uptake of fluorescent genes.

Following this exciting discovery, several research groups, both experimental and theoretical rushed to understand the physical behaviours behind this increased uptake. It accepted that the oscillatory behaviour of the microbubbles under sonication forms small transient pores in the cell membrane, yet the question of *how* still remained.

Hence, the next step was to attempt to understand the 115 physical mechanisms of microbubble-cell interaction using 116 optimised acoustic conditions. Whilst a vast majority of 117 high-speed imaging is performed at several million frames 118 per second, we attempted to record at slower frame rates 119 giving us longer imaging duration but with sufficient temporal 120 resolution to track microbubble motion. For this work we 121 used immortalised cervical cancer cells (HeLa cells). $1.6 \times$ 122 10^6 cells were seeded into an OptiCell $^{\rm \tiny IM}$ and incubated for 123 24 hours. Following incubation, low concentrations of DiD 124 lipophilic fluorescent probe (Vybrant^{1M} Molecular Probes) 125 tagged lipid shelled microbubbles, Definity®, SonoVue® or 126 MicroMarker^{1M} (Bracco) were injected into the OptiCell^{1M}. 127 Concentrations of 1 bubble to 1 cell were aimed for. The 128 OptiCell^{1M} was then placed in a water bath cell-side up, 129 ensuring the floating bubbles were in contact with the cells. 130 Two ultrasound frequencies were evaluated: 1.0 MHz and 131 6.6 MHz with duty cycles of 40% and 50%, respectively. 132 The MI for both frequencies was ≤ 0.2 and acoustic 133 intensity (I_{SPTA}) was ≪1 W/cm². Imaging was performed 134 using a high-magnification, high-numerical aperture (NA) 135 setup using a $60 \times / 0.90$ -NA objective or using a confocal 136 fluorescence microscope with a $40 \times /1.30$ -NA objective. 137 Following sonication, a range of phenomena was observed, 138 including bubbles hammering back and forth into cells [11], 139 to complete bubble translation into a cell. To validate if the 140 microbubble had truly translated into the cell, 3D confocal 141 fluorescence imaging was performed. Results showed that the 142 lipid bonded fluorescent probe was then within the cell. To 143 exclude the chance of endocytosis, experiments were re-done 144 with endocytosis inhibitors where no change in localised 145 microbubble translation was seen. 146

In conclusion, to date there are several non-exclusive 147 phenomena that can be directly related to this increased up 148 take [12]: Push and pull, where the microbubble near, or in 149 contact with the cell, stretches the bilayer membrane opening 150 pores. Jetting, when a bubble collapses asymmetrically, 151 making a needle-like point towards the cell. Streaming, 152 where the fluid flow around the bubble near a cell is 153 strong enough to shear the membrane apart. Translation, 154 where a microbubble with a therapeutic load can be forced 155 inside a cell. Inertial cavitation, where the violent formation 156



Fig. 2. Illustration depicting five of the sonopration mechanisms.

and destruction of bubbles generates shock-waves damaging 157 everything in the local area. It is important to understand 158 that all these phenomena are non-exclusive, and each and 159 all have advantages and disadvantages. In addition, when 160 sonoporation will occur most effectively when exciting the 161 microbubbles at their resonance frequency. Figure 2 illustrates 162 five of the sonoporation mechanisms that induce transient cell 163 membrane porosity. In addition to these phenomena there are 164 several controversial yet thought-provoking theories regarding 165 the acoustic activity of the cell membranes themselves [13]. 166

IV. PRE-CLINICAL WORK

167

Pre-clinical mouse work was performed to evaluate the effect of sonoporation in combination with a chemotherapeutic on an orthotopic pancreatic ductal carcinoma (PDAC) model [71] [14].

The orthotopic xenograft model was developed by injecting 172 1×10^{6} bioluminescent human PDAC cells (MIA PaCa-2^{luc}) 173 into the head of the pancreas of immunodeficient (NOD-scid 174 IL2r γ^{null}) mice. The bioluminescent cells allowed for 2D full 175 body imaging helping visualise the spread of the tumour cells. 176 The 3D development of the primary tumour was measured 177 using high-resolution 3D ultrasound (MS250 probe and Vevo 178 2100, VisualSonics Inc, Ontario, Canda). 179

Mice were imaged and treated once weekly using Gemzar (Eli Lilly and Company, Indianapolis, IN) the current best chemotherapeutic for treating PDAC, or Gemzar in combination with ultrasound and SonoVue[®] microbubbles. Acoustic conditions remained the same as in the *in-vitro* work, where bubbles were forced into the cells, *i.e.*, 1-MHz, 40



Fig. 3. Bioluminesence imaging and tumour volume of mice with orthotopic pancreatic cancer treated with Gemzar alone or Gemzar with microbubbles and ultrasound.

cycles, MI=0.2, with a 40% duty cycle. Treatment was started 3 weeks after xenografting.

After only two treatments, *i.e.*, 5 weeks after xenografting, 188 a statistical significant difference in tumour volume was seen 189 between the sonoporated group and the chemotherapy alone 190 group. Figure 3 shows the bioluminescent images and tumour 191 volumes 10 weeks after xenografting (after six treatments). 192 By week 10, the primary tumour in the sonoporation group 193 was nearly 4 times smaller than the control and chemotherapy 194 alone tumour. In addition the metastatic spread was limited 195 to the liver in the sonoporation group whereas the Gemzar 196 and control group showed whole body metastatic spread. 197 In conclusion, translating the identical acoustic conditions 198 seen to force microbubbles into cells in vitro resulted in 199 primary tumour inhibition and delayed metastatic spread when 200 combining a chemotherapeutic. 201

V. CLINICAL TRIALS

202

186

187

In many cases pre-clinical work does not translate well into the clinic. In our work here we perform, to our knowledge, the world first clinical trial on sonoporation [15].

A total of five patients with PDAC were treated using 206 Gemzar, SonoVue[®] and ultrasound. The chemotherapeutic 207 was administered in accordance to the manufacturer's 208 guidelines; infusion over 30 minutes. At the completion 209 of the infusion, the concentration of the chemotherapeutic 210 was maximal in the blood. At this point sonoporation was 211 initiated. Ultrasound was generated using a clinical diagnostic 212 scanner (GE LOGIQ 9 + 4C probe, GE Healthcare, Waukesha, 213 WI), with optimised acoustic output conditions. The clinical 214 diagnostic machine allowed for simultaneous imaging and 215 treatment of the primary tumour. The output frequency was 216 1.8 MHz with a duty cycle of 1%. To compensate for the 217

reduced duty cycle the ultrasound treatment lasted 31.5 min. 218 SonoVue[®] were injected intravenously every 3.5 minutes 219 resulting in 9 doses of 0.5 mL each over 31.5 minutes. 220

To evaluate the effect of the treatment, we compared the 221 number of chemotherapy cycles our patients could undergo 222 versus 80 historical patients treated at the same hospital. The 223 healthier the patients are, the more treatment cycles they 224 can undergo. The historical group was able to receive 9 ± 6 225 treatment cycles. The patients treated with sonoporation were 226 able to undergo almost double the amount of treatment with 227 16 ± 7 cycles of chemotherapy. It is important to note that 228 several patients were removed from the study due to the 229 efficacy of the treatment. These patients would be able to go 230 substantially more cycles, hence the treated cycles number is 231 artificially low. 232

Maximum tumour diameters were measured using CT 233 imaging from independent radiologists. The primary tumour 234 growth was inhibited and in some cases, even reduced. 235

In conclusion, the clinical trials showed that there are 236 no immediate adverse side effects to sonoporation, and that 237 sonoporation can help increase the efficacy of the treatment. 238

239

VI. FUTURE WORK

Whilst sonoporation has shown success in the lab and in 240 the clinic, there is still a lot of room for improvement. To 241 date, sonoporation has primarily been performed by combining 242 existing materials and technologies aimed for other tasks, for 243 example the ultrasound contrast agents and clinical diagnostic 244 machines are designed for imaging not treatment. Both these 245 aspects need to be drastically improved. 246

A. Microbubbles 247

In the area of microbubbles, new agents need to be 248 developed that can carry a therapeutic load to the target region. 249 These bubbles need to be highly acoustically sensitive, and be 250 able to release their load on demand. Current concepts allow 251 loading inside the shell, or on the surface of the shell, hence 252 only small loads can be carried. New methods and concepts are 253 needed that allow for larger volumes of therapeutic loads and 254 specific release mechanisms. Such a method would only treat a 255 localised region, and would have minimal affect on metastasis, 256 for this reason, it may be ideal to combine sonoporation with 257 traditional chemotherapy, ensuring maximum affect on both 258 primary tumour and metastasis. 259

B. Ultrasound 260

Clinical diagnostic machines are highly specialised and have 261 inverse acoustic properties than those needed for sonoporation. 262 Sonoporation requires a large area to be treated, so a large 263 footprint probe, that can treat a whole tumour is necessary. In 264 addition, longer duty cycles are required to continuously excite 265 the microbubbles in the target region increasing the cellular 266 permeability. 267

VII. CONCLUSION

Although sonoporation was discovered only in 1997 and 269 clinical trials were first done from 2013, we have shown 270 that there is great potential in sonoporation. In-vitro results 271 have shown that microbubbles can be forced into cells, 272 whereas pre-clinical mice studies under identical conditions 273 have shown to reduce primary tumour burden and metastatic 274 development. Our clinical trial showed that combining 275 sonoporation and chemotherapy can help the patients to 276 undergo more chemotherapy cycles, presumably prolonging 277 survival. 278

ACKNOWLEDGEMENTS

This study has been supported by MedViz 280 (http://medviz.uib.no/), an interdisciplinary research cluster 281 from Haukeland University Hospital, University of Bergen, 282 and Christian Michelsen Research AS. The authors would like 283 to thank The Michelsen Centre for Industrial Measurement 284 Science and Technology for their support throughout this 285 project.

REFERENCES

- [1] V. J. Robertson and K. G. Baker, "A review of therapeutic ultrasound:
- effectiveness studies," Phys Ther, vol. 81, no. 7, pp. 1339-1350, 2001. [2] K. Hynynen, "MRI-guided focused ultrasound treatements," Ultrasonics,
- vol. 50, no. 2, pp. 221-229, 2010. [3] B. Geers, H. Dewitte, S. C. De Smedt, and I. Lentacker, "Crucial factors and emerging concepts in ultrasound-triggered drug delivery," J Control Release, vol. 164, no. 3, pp. 248-255, 2012.
- [4] C. Peetla, S. Vijayaraghavalu, and V. Labhasetwar, "Biophysics of cell membrane lipids in cancer drug resistance: Implications for drug transport and drug delivery with nanoparticles." Adv. Drug Deliv. Rev., vol. 65, no. 13-14, pp. 1686-98, Nov. 2013.
- [5] R. Gramiak and P. M. Shah, "Echocardiography of the aortic root," Invest. Radiol., vol. 3, no. 5, pp. 356-366, 1968.
- M. Postema, A. Bouakaz, M. Versluis, and N. de Jong, [6] "Ultrasound-induced gas release from contrast agent microbubbles," IEEE Trans. Ultrason. Ferroelectr. Freq. Control, vol. 52, no. 6, pp. 1035-1041, 2005.
- [7] G. ter Haar, "Safety and bio-effects of ultrasound contrast agents," Med. Biol. Eng. Comput., vol. 47, no. 8, pp. 893-900, 2009.
- S. Kotopoulis and M. Postema, "Microfoam formation in a capillary," [8] Ultrasonics, vol. 50, no. 2, pp. 260-268, 2010.
- [9] S. Bao, B. D. Thrall, and D. L. Miller, "Transfection of a reporter plasmid into cultured cells by sonoporation in vitro," Ultrasound Med. Biol., vol. 23, no. 6, pp. 953-959, 1997.
- [10] Unger and E, "SonoPoration and Gene Delivery with Acoustically Active Carriers," p. S2, 1997.
- P. Marmottant and S. Hilgenfeldt, "Controlled vesicle deformation and [11] lysis by single oscillating bubbles," Nature, vol. 423, pp. 153-156, 2003.
- [12] M. Postema, S. Kotopoulis, A. Delalande, and O. H. Gilja, "Sonoporation: Why microbubbles create pores," Ultraschall Med., vol. 33, no. 1, pp. 97-98, 2012.
- [13] B. Krasovitski, V. Frenkel, S. Shoham, and E. Kimmel, "Subcellular sonophores: ultrasound induced intramembrane cavitation," Proc. Nat. Acad. Sci., vol. 108, no. 8, pp. 3258-3263, 2011.
- [14] S. Kotopoulis, A. Delalande, M. Popa, V. Mamaeva, G. Dimcevski, O. H. Gilja, M. Postema, B. T. Gjertsen, and E. Mc Cormack, "Sonoporation-enhanced chemotherapy significantly reduces primary tumour burden in an orthotopic pancreatic cancer xenograft," Mol. Imaging Biol., vol. 16, no. 1, pp. 53-62, 2014.
- [15] S. Kotopoulis, G. Dimcevski, O. H. Gilja, D. Hoem, and M. Postema, "Treatment of human pancreatic cancer using combined ultrasound, microbubbles, and gemcitabine: A clinical case study," Med. Phys., vol. 40, no. 7, pp. 72902-72909, 2013.

268

279

286

287

288

289

290

291

292

293

294

295

296

297

298

299

300

301

302

303

304

305

306

307

308

309

310

311

312

313

314

315

316

317

318

319

320

321

322

323

324

325

326

327

328

329

330