

Synthetic cannabinoids in e-liquids: A proton and fluorine NMR analysis from a conventional spectrometer to a compact one

Nao Wu, Saïda Danoun, Stéphane Balayssac, Myriam Malet-Martino, Véronique Gilard, Catherine Lamoureux

▶ To cite this version:

Nao Wu, Saïda Danoun, Stéphane Balayssac, Myriam Malet-Martino, Véronique Gilard, et al.. Synthetic cannabinoids in e-liquids: A proton and fluorine NMR analysis from a conventional spectrometer to a compact one. Forensic Science International, 2021, 324, pp.110813. 10.1016/j.forsciint.2021.110813. hal-03245107

HAL Id: hal-03245107

https://hal.science/hal-03245107

Submitted on 24 May 2023

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



- 1 Synthetic cannabinoids in e-liquids: a proton and fluorine NMR analysis from a
- 2 conventional spectrometer to a compact one
- 3 Nao Wu^{1#}, Saida Danoun¹, Stéphane Balayssac^{1§}, Myriam Malet-Martino¹, Catherine
- 4 Lamoureux², Véronique Gilard^{1§}
- ¹Groupe de RMN Biomédicale, Laboratoire SPCMIB (UMR CNRS 5068), Université Paul
- 6 Sabatier, Université de Toulouse, 118 route de Narbonne, 31062 Toulouse Cedex, France
- ² SCL, Laboratoire de Paris, 25 Avenue de la République, 91300 Massy
- 8 **Present address: National R&D Center for Se-rich Agricultural Products Processing, School
- 9 of Modern Industry for Selenium Science and Engineering, Wuhan Polytechnic University,
- 10 Wuhan 430023, China.
- 11 §Present address: Laboratoire des IMRCP, Université de Toulouse, CNRS UMR 5623,
- Université Toulouse III Paul Sabatier 118 Route de Narbonne, 31062, Toulouse Cedex 9,
- 13 France
- 16 Corresponding author: Pr Véronique GILARD, Laboratoire des IMRCP, Université de
- 17 Toulouse, CNRS UMR 5623, Université Toulouse III Paul Sabatier 118 Route de Narbonne,
- 18 31062. Toulouse Cedex 9. France
- 19 Tel: 33561558281; e-mail: gilard@chimie.ups-tlse.fr

20

14

15

Abstract

- 2 The ¹H NMR profiles of 13 samples of e-liquids supplied by French customs were obtained
- with high-field and low-field NMR. The high-field ¹H NMR spectra allowed the detection of
- 4 matrix signals, synthetic cannabinoids, and flavouring compounds. Quantitative results were
- 5 obtained for the five synthetic cannabinoids detected: JWH-210, 5F-MDMB-PICA, 5F-ADB,
- 6 5F-AKB48, and ADB-FUBINACA. Conventional GC-MS analysis was used to confirm
- 7 compound identification. Fluorine-19 NMR was proposed for the quantification of fluorinated
- 8 synthetic cannabinoids and was successfully implemented on both 400 MHz and 60 MHz
- 9 NMR spectrometers. This study based on few examples explored the potentiality of low-field
- 10 NMR for quantitative and quantitative analysis of synthetic cannabinoids in e-liquids.

Keywords: e-liquid, e-cigarette, qNMR, compact NMR, 19F NMR, cannabinoid

Introduction

Electronic cigarettes or e-cigarettes have spread rapidly in recent years as an alternative method for nicotine delivery. Compared to traditional tobacco, they are inexpensive, easy to use, and described as safer for public health [1, 2] although they have been associated with a variety of health effects, particularly the exposure to some constituents of aerosols that might result in respiratory complications [3, 4]. Electronic liquids or e-liquids are formulations used in e-cigarettes which can vaporize aerosol for inhalation when the e-cigarette is activated. E-liquids usually contain a mixture of propylene glycol (PG) and glycerol as the main matrix but also flavours and in some case a pharmaceutical and/or herbal agent inside [5, 6].

Synthetic cannabinoids (SCs) are the most rapidly growing class of recreational designer drugs used in e-liquids [6-8]. Their chemical structures differ from naturally occurring cannabinoids, but they have similar pharmacological activities by binding with CB1 and/or CB2 receptors [9]. SCs could enhance the excitatory effect by promoting the ability to pass the blood-brain barrier because of chemical structure modifications reducing the molecular polarity [10]. SCs were originally used as therapeutic drugs for their anti-convulsant, anti-nociceptive, and anti-psychotic properties, but long time inhalation or high dose intake can induce addiction, cardiovascular disease, mental confusion, high blood pressure, seizures, acute kidney injury, cancer, and even death [11]. SCs usually exhibit more potency than natural cannabis as they are designed to bind more strongly to CB1/CB2

receptors in the neuronal system [12], and they also may interact with other receptors in the brain that natural cannabis does not [13]. Due to the fast generation of SCs, their structures are quite diverse and a large number of new psychoactive SCs have appeared on sale, which makes it even more difficult to detect and distinguish them [14]. Over the past decade, various studies have described the detection of SCs, most of them focusing on the identification and structural characterization of new cannabinoids [15-21]. In addition to these researches, the use of different analytical methods for the detection or/and quantification of SCs in commercial herbal mixture blends has become a vibrant area of research [22-29]. However, only a few studies on screening for SCs or natural cannabinoids in e-liquids have been reported. Peace and co-workers have identified SCs and other components in commercially available e-liquid formulations using direct analysis in real time mass spectrometry (DART-MS) and gas chromatography-mass spectrometry (GC-MS), as well as naturally occurring cannabinoids quantified by high-performance liquid chromatography-triple quadrupole mass spectrometry (HPLC-MS/MS) [5, 7, 30]. Poklis et al. evaluated 5F-ADB in 9 cannabidiol e-liquids from the same manufacturer by DART-MS and GC-MS in 2019 [8].

Nuclear Magnetic Resonance (NMR) spectroscopy has already been described as a powerful method for the analysis of complex mixtures and recent review articles summarized the promising capabilities of NMR in the fields of forensic and fraud analysis for both qualitative and quantitative purposes [31-33]. Concerning SCs, NMR was used for their structural identification and for the analysis of herbal blends [18, 22, 25-27, 34]. Although high-field (HF) NMR is not yet considered as a routine tool in analytical chemistry mainly because of its high cost, benchtop cryogen-free low-field (LF) NMR which uses much cheaper instruments can partially solve this cost issue and be implemented in forensic or customs laboratories. Indeed, compact NMR spectroscopy has been applied recently for the quality control of various matrices, including meat [35], edible oils [36], drugs [37], dietary supplements [38, 39] and herbal blends containing cannabinoids [24].

The aim of the present study was to assess the capacity of HF and LF NMR for the characterization of e-liquid contents. Proton (1 H) and fluorine (19 F) NMR analyses were carried out with both HF (400 MHz for 1 H) and LF (60 MHz for 1 H) spectrometers and GC-MS was used as a confirmation and reference method. To our knowledge, this study is the first to use compact LF NMR for the characterization of SCs and flavouring agents in e-liquids.

2. Experimental

2.1. Materials

67 13 e-liquid samples from customs seizures were provided by the SCL laboratory
68 (French Customs). The standard SCs JWH-210, 5F-AKB48, and 5F-MDMB-PICA were
69 obtained from Cayman Chemical (Ann Arbor, MI, USA). 3-(trimethylsilyl)propionic-2,2,3,370 d4 acid sodium salt (TSP, purity≥99%, isotopic purity 98 atom % D), α,α,α-trifluorotoluene
71 (TFT, purity≥99%), vanillin, ethyl maltol, limonene, anisyl alcohol, raspberry ketone (4-(472 hydroxyphenyl)butan-2-one), menthol, linalool and chromium(III) acetylacetonate (Cr(acac)₃)
73 were purchased from Sigma-Aldrich (St. Louis, MO, USA).

74 2.2. NMR analysis

2.2.1. Sample preparation for NMR

200 μl of each e-liquid withdrawn with a glass precision syringe were diluted with 400 μl of methanol-d4 (MeOD) and mixed with 30 μl of a TSP solution (10 mM in MeOD) as reference for ¹H NMR experiments. ¹⁹F NMR experiments were run for the 10 e-liquid samples **4-13** containing the fluorinated SCs 5F-MDMB-PICA, 5F-ADB, 5F-AKB48, and ADB-FUBINACA. For these ¹⁹F NMR experiments, 200 μl of sample were diluted with 400 μl of MeOD and then mixed with 9 μl of a TFT solution (1M in MeOD; final concentration 15 mM). Two mg of Cr(acac)₃ as a relaxation agent were added into the NMR tube. The HF and LF ¹⁹F NMR experiments were recorded on the same NMR tube. The exact concentrations of the TSP and TFT solutions used were determined relative to that of a solution of yohimbine hydrochloride (European Pharmacopoeia Reference Standard). All quantitative experiments were done in triplicate from sampling the e-liquid to processing the NMR spectrum.

2.2.2. HF NMR analysis

HF ¹H NMR experiments were performed on a Bruker Avance 400 spectrometer (Bruker Biospin AG, Fallanden, Switzerland) equipped with a 5 mm triple resonance probe (TXO) with ¹⁹F direct detection. Typical acquisition parameters optimized for quantitative experiments were as follows: number of scans (NS) 128, pulse width 4.52 μs (flip angle 30°), acquisition time 3.40 s, spectral width 12 ppm, 32K data points, and an additional relaxation delay of 10 s; the recording time was thus ca. 29 min. The signal of TSP set at 0 ppm was used as an internal reference for chemical shift measurement and quantitative analysis. Three samples with lower concentrations of SCs (1-3) were recorded by increasing the NS to 1024

for obtaining a correct signal-to-noise (S/N) ratio. HF ¹H NMR data were processed using the Bruker TopSpin 4.0.8 software with one level of zero-filling and Fourier transformation after multiplying FIDs by an exponential line-broadening function (LB) of 0.3 Hz.

HF ¹⁹F NMR experiments were recorded with the same spectrometer and probe at 376 MHz. First, a basic zg sequence allowed to acquire ¹H-coupled ¹⁹F NMR spectra, then quantitative experiments were carried out using an inverse gated decoupling pulse sequence with the following parameters: NS 32, pulse width 4.52 μs (flip angle 30°), acquisition time 1.57 s, spectral width 221 ppm, 256K data points, and an additional relaxation delay of 3 s for complete longitudinal relaxation; the recording time was thus less than 3 min. The signal of TFT set at -63.72 ppm with respect to CFCl₃ (0 ppm) was used as an internal reference for chemical shift measurement and quantitative analysis.

2.2.3. LF NMR analysis

Spectra were acquired on a PulsarTM benchtop NMR spectrometer (Oxford Instruments, Abingdon, UK) operating at a frequency of 59.7 MHz for ¹H and 56.17 MHz for ¹⁹F. The temperature inside the spectrometer was 310 K. Acquisition was performed with the SpinFlow 1.2.0.1 software (Oxford Instruments) and data processing with the MNova 11.0 software (Mestrelab Research, Santiago de Compostela, Spain).

LF 1 H NMR spectra were recorded with a flip angle of 90° (12.5 µs), a spectral width of 5000 Hz, and 8K data points (acquisition time of 1.64 s). The relaxation delay was set at 2 s, and 1024 transients were recorded leading to a total acquisition time of \approx 1 h. For data processing, the FIDs with one level of zero-filling were apodized with an exponential (0.3 Hz) filter and an automatic Whittaker smoother baseline correction was applied.

LF ¹⁹F NMR spectra were recorded with a flip angle of 90°, a spectral width of 20833 Hz, and 32K data points (acquisition time of 1.57 s). The relaxation delay was set at 1 s, and 2048 transients were recorded leading to a total acquisition time of 1 h 27 min. For data processing, the FIDs were apodized with an exponential (3.0 Hz) filter and an automatic baseline correction was applied.

2.2.4. T_1 measurements

T₁s of protons targeted for HF quantification (written in bold in Table 1) and of ¹⁹F nuclei (HF and LF) were measured by the inversion-recovery pulse sequence with 20 inversion delays from 0.001 to 25 s. For these T₁ measurements, the acquisition parameters were as follows: 64K data points, acquisition time 6.81 s, spectral width 12 ppm, number of

scans 32, and relaxation delay 40 s for HF ¹H NMR; 128K data points, acquisition time 1.75

s, spectral width 100 ppm, number of scans 16, and relaxation delay 5 s for HF ¹⁹F NMR;

131 16K data points, acquisition time 2.05 s, spectral width 133 ppm (8000 Hz), number of scans

96, and relaxation delay 5 s for LF ¹⁹F NMR.

133

134

138

139

140

145

148

149

150

151

152

153

154

155

156

157

158

159

132

2.2.5. NMR quantification

135 Concentrations were measured by comparing the areas of convenient signals (written 136 in bold in Table 1) of targeted compounds with that of the reference signal, the area of each 137 NMR peak being directly proportional to the number of nuclei giving rise to it.

After phasing and baseline correction of the spectra, the targeted signals were integrated. The amount of each compound expressed in mg/ml of e-liquid sample was determined using the following general equation:

141
$$Amount(mg/ml) = \frac{A_{SCs}}{N_{SCs}} \times [REF] \times \frac{N_{REF}}{A_{REF}} \times V_{total} \times \frac{Mw_{SCs}}{V_{e-liquid}}$$

142 [*REF*]: concentration of the reference (TSP or TFT),

143 A_{SCs} and A_{REF} : signal areas of SCs and reference signals, respectively,

 N_{SCs} and N_{REF} : number of nuclei giving rise to the signal considered,

 V_{total} : total volume of solution in the NMR tube,

146 Mw_{SCs} : molecular weight of SCs,

 $V_{e-liquid}$: volume of e-liquid used for the NMR assay.

Linear regression and the Wilcoxon signed-rank test were used to compare the amounts of fluorinated SCs measured by HF ¹⁹F NMR versus HF ¹H NMR and versus LF ¹⁹F NMR. The significance of the slope and the fit of the regression line were tested by t-test and F-test respectively.

2.3. GC-MS analysis

Each sample was vortexed for 30 s before aliquoting. 250 μl of each e-liquid sample were then mixed with 500 μl of CHCl₃. The analyses were performed using an Agilent 6890 N Gas Chromatograph coupled with a GCT Premier CAB109 TOF mass spectrometer (Waters Corporation, Milford, MA, USA). An Agilent HP-5MS capillary column (30 m × 0.25 mm, 0.25 μm film thickness) was used to perform the chromatographic separation. Helium as the carrier gas was delivered at a constant flow of 1.0 ml/min. The temperature program was as follows: initial temperature 70 °C for 2 min, increase to 200 °C at 5 °C/min

- then to 280 °C at 10 °C/min, 24 min holding time at 280 °C for a total run time of 60 min.
- 161 Injections were performed with a 1:20 split ratio and injection volumes of 1.0 μL. Mass
- spectra were recorded in the range m/z 30-650. Compounds were identified using the
- SWGDRUG 3.3 (www.swgdrug.org) and NIST 08 databases.

3. Results and discussion

164

165

166

167

168

169

170

171

172

173

174

175

176

177

178

179

180

181

182

183

184

185

186

187

188

189

190

191

3.1. HF ¹H NMR qualitative analysis of e-liquids

In the first step of this study, the 13 e-liquid samples were screened by HF ¹H NMR in order to characterize their contents in SCs and other chemicals. GC-MS was used as the gold standard method to confirm the identification of SCs as well as flavouring components.

The HF ¹H NMR spectra of all 13 samples allowed the detection of major signals of propylene glycol (PG) and glycerol which constitute the matrix of e-liquids. As an illustration, Figure 1A shows the complete HF ¹H NMR spectrum of one e-liquid sample (3). Matrix signals of PG and glycerol have the most prominent signals in the upfield area. For PG, a doublet (d; CH₃) was detected at 1.11 ppm (J = 6.4 Hz), a multiplet at 3.77 ppm (<CH), an ABd system centered at 3.41 ppm (J = 11.1, 6.5, 4.8 Hz) (-CH₂), and broad singlets for hydroxyl protons at 4.60 and 5.24 ppm. The methine and methylene groups of glycerol led to a multiplet at 3.66 ppm, an ABd system at 3.58 and 3.51 ppm (J = 11.3, 6.1, 4.8 Hz) whereas hydroxyl groups resonated as doublet and triplet at 5.29 and 5.19 ppm respectively (depending on the rate of hydroxyl proton exchange, these signals were more or less broad). Moreover, ethylene glycol and/or polyethylene glycol were identified in samples 9-12 from their characteristic intense broad signals around 3.6 ppm. At last, ethanol was detected in samples 6, 8, 10, and 11 by its methyl signal at 1.2 ppm whereas the quadruplet of its methylene group overlapped with the signal of glycerol at 3.59 ppm. The presence of ethanol in e-liquid increases the solubility of SCs and flavors and can significantly lower the viscosity of glycerol, which improves the formation of the aerosol for inhalation. Peace and co-workers previously evaluated alcohol content in e-liquid by headspace gas chromatography with flame ionization detector (HS-GC-FID) [7].

In addition to the signals of the matrix, the resonances of psychoactive SCs were detected in the NMR spectra. Five SCs were identified: JWH-210 was detected in samples **1-3**, 5F-MDMB-PICA in samples **4-6**, 5F-ADB in samples **7-9**, ADB-FUBINACA in samples **10-12**, and 5F-AKB48 in sample **13**. Structures of these SCs (Table 1) are based on indole or indazole rings with differently substituted carboxamide and alkyl side chains, four SCs being

fluorinated. The HF ¹H NMR profiles (5.6-10.0 ppm) of samples **3**, **4**, **8**, **11**, and **13** are reported in Figure 1B and the assignments of SCs signals are given in Table 1. Characteristic signals with close chemical shifts and multiplicity are observed for the indazole protons (H4, H5, H6, and H7) of 5F-ADB, ADB-FUBINACA and 5F-AKB48 whereas signals of the indole ring (H2, H4, H5, H6, and H7) in JWH-210 and 5F-MDMB-PICA have a different spectral signature. For SCs bearing fluorine on the alkyl chain (5F-MDMB-PICA, 5F-ADB, and 5F-AKB48), a doublet of triplet is observed for H5" due to both ²J coupling with fluorine nucleus and ³J_{HH} with neighboring methylene. Protons from the methylene groups H1" adjacent to the nitrogen atom of indazole or indole ring give rise to characteristic deshielded triplets (≈4.3-4.5 ppm). Some signals could not be detected due to their overlap with matrix resonances; for instance, the singlet of methyl ester (5F-ADB or 5F-MDMB-PICA) overlaps with the signal of propylene glycol around 3.7-3.8 ppm. NMR characteristics for SCs are in accordance with previously published data [15, 16, 19, 24].

In addition to matrix and SCs signals, resonances of other chemicals were detected in e-liquid HF ¹H NMR spectra. Indeed, flavouring ingredients are widely added in e-liquids to give fruit aromas or other pleasant scents. Flavouring compounds like limonene, vanillin, ethyl maltol, raspberry ketone, anisyl alcohol, linalool, and menthol were successfully identified by spiking with standard compounds and their HF ¹H NMR assignments are reported in Table 2. Vanillin, the most common flavor agent, was detected in 8 out of 13 samples. Its more significant signal is that of the aldehyde proton which leads to a deshielded singlet at 9.74 ppm (Figure 1B); protons on the aromatic ring resonate in the range 6.9-7.5 ppm whereas the methoxy group singlet is detected at 3.92 ppm. Ethyl maltol which has sweet caramelized sugar smell was detected in six samples. It has two characteristic doublets at 8.00 and 6.42 ppm corresponding to aromatic protons, whereas protons from the ethyl group are detected at 2.74 ppm (q, H2) and 1.22 ppm (t, H1). Characteristic signals of raspberry ketone were observed in three samples with doublets of the para-substituted aromatic ring at 7.00 and 6.69 ppm and singlet of the methyl group at 2.11 ppm. Linalool was found in samples 10 and 11 from its characteristic doublet of doublet at 5.89 ppm as shown in Figure 1B. The characteristic upfield signals of menthol allowed its detection in sample 13. Anisyl alcohol and limonene were also each detected in only one sample.

Besides NMR analysis, GC-MS was performed for verification of SCs and flavouring compounds identification. Retention times, exact masses and characteristic fragments are reported in Table S1 for SCs and Table 2 for flavours. These data are also supported by previous GC-MS characterization of e-liquids [5, 30].

3.2. LF ¹H NMR profiling

All samples were then analysed on the benchtop 60 MHz ¹H NMR spectrometer using the same solutions in the same NMR tubes than for HF recording. As an illustration, the HF and LF ¹H NMR spectra of sample 8 are shown in Figure 2A. Compared to the HF spectrum, signals are much broader with more overlap due to the lower resolution at LF. The resonances are especially crowded in the low frequency region (0.5-4.5 ppm) where targeted signals are erased by intense matrix signals, but some clues can nevertheless be drawn from the ethylenic/aromatic region i.e. above 5.6 ppm. Figure 2B shows the five typical LF profiles of samples 3, 4, 8, 11, and 13 with the same zoom than for their HF spectra illustrated in Figure 1B. The SC JWH-210 being at a low concentration in sample 3 could not be detected under the chosen experimental conditions. For the four other samples, even if signals are broad and similarities appear, the LF profile of each sample differs from the others and gives specific spectral signatures. For instance, singlets detected at 8.0 and 5.7 ppm in samples 4 and 11 are characteristic of 5F-MDMB-PICA (H2) and ADB-FUBINACA (H1") respectively. The profiles of 5F-ADB and 5F-AKB48 in samples 8 and 13 are close even if the multiplet at \approx 7.2 ppm is more spread out in sample 13. Concerning flavouring compounds, the singlet of vanillin is observed at 9.74 ppm in samples 4 and 8, likewise for the two characteristic doublets of ethyl maltol at 8.0 and 6.4 ppm in sample 8.

From LF ¹H NMR analyses of all samples, it appears that despite intense matrix signals, some characteristic signals of SCs can be detected if their concentration is not too low. These LF ¹H NMR experiments constitute a proof of concept that direct profiling of e-liquids containing SCs could be relevant and should be applied on a larger number of samples for considering a chemometrics treatment as recently proposed in agrifood and pharmaceutical applications [39-41]. However, we did not further investigate the ¹H NMR experiments on the compact spectrometer because the method was not suitable for quantitative purposes in this study.

3.3. ¹⁹F NMR analysis

Because of their enhanced potency on CB1 receptors, an increasing number and rate of emergence of fluorinated SCs have been reported over the last decades [42, 43]. In this context, it is therefore of interest to propose a specific method for detecting fluorinated molecules. Insofar the ¹⁹F nucleus has favourable NMR characteristics (nuclear spin of 1/2, relatively narrow lines, 100% natural abundance, good sensitivity (83% that of proton) and

large chemical shift range), the fluorine NMR analysis can be used to study these fluorinated cannabinoids as reported recently by Naqi et al. who applied HF ¹⁹F NMR spectroscopy to quantify the amounts of two fluorinated SCs, AM-694 and 5F-ADB, in herbal incense packages [43].

¹⁹F NMR spectra were recorded for all samples at both HF and LF, i.e. 376 and 56.2 MHz respectively. Typical spectra without proton decoupling are presented in Figure 3, with HF and LF spectra stacked. TFT was added as a reference for chemical shift (-63.72 ppm with respect to CFCl₃) and quantification. The fluorine atom of ADB-FUBINACA linked to the aromatic ring leads to a singlet at -116.0 ppm. Fluorine nuclei at the end position of the alkyl chains of 5F-AKB48, 5F-MDMB-PICA, and 5F-ADB give rise to signals at ≈-220 ppm. The chemical shifts of 5F-AKB48 and 5F-ADB are almost identical (Δv < 1 Hz at HF), while that of 5F-MDMB-PICA slightly differs with Δv ≈37 Hz at HF and ≈6 Hz at LF. They appear as triplet of triplet at HF due to the coupling of the ¹⁹F nucleus with the neighbouring methylene protons H5" (2 J_{HF} = 47.5 ± 0.1 Hz) and H4" (3 J_{HF} = 25.5 ± 0.2 Hz) while, at LF, an apparent septuplet (3 J_{HF} = 23.7 ± 0.4 Hz) is observed because signals are much more spread out as coupling constants are independent of the applied magnetic field.

Unlike LF proton NMR analyses of e-liquids that have limitations for quantitative purposes and cannot be implemented due to strong signal overlap, ¹⁹F NMR opens up the possibility of setting up quantitative experiments.

3.4. Quantitative analysis

3.4.1. Setting up of quantitative conditions

The aim of the quantitative study was to consider HF proton and fluorine NMR as reference quantitative methods to evaluate the ability of quantification of LF ¹⁹F NMR. Indeed, HF NMR is recognized as a powerful method for quantification of compounds in various mixtures for both ¹H [24, 25, 43-45] and ¹⁹F [43, 44, 46, 47] nuclei.

To accurately quantify SCs in e-liquids, MeOD was chosen as the solvent due to its good miscibility with e-liquids, and T_1 relaxation times were determined by the classical inversion-recovery pulse sequence method in order to define the relevant acquisition parameters for quantitative 1H and ^{19}F NMR experiments. HF T_1 values of protons selected for the quantification of SCs and internal reference were measured in samples 3, 4, 7, 11 and 13 (Table 3). All 1H resonances of SCs were considered as fully relaxed since 99.9% of the signal intensity of the proton with the longest T_1 (2.7 s) was recovered in the recording

conditions used (flip angle of 30° and repetition time of 13.4 s). For ¹⁹F NMR quantifications, the relaxation agent Cr(acac)₃ was added into the NMR tube to shorten the relaxation times of the ¹⁹F nuclei. The T₁s of the ¹⁹F of the four fluorinated SCs in e-liquids are reported in Table 3. The longest T₁ for SCs was found for 5F-ADB in sample **7** (0.14 s and 0.18 s at HF and LF respectively) and the T₁ of the internal reference TFT was measured at 0.27 s at LF. For both HF and LF, spectra recording conditions were chosen to gain full relaxation of the ¹⁹F nuclei. Repetition times of 4.57 s and 2.57 s were used for HF and LF NMR with flip angles of 30° and 90°, respectively. All ¹⁹F nuclei were thus completely relaxed since more than 99.9% of the ¹⁹F magnetization was recovered.

Recording times were set at 2.5 min (32 scans) and 90 min (2048 scans) for HF and LF ¹⁹F NMR experiments leading to S/N ratios of 200 and 50, respectively. The limits of detection (LOD) and quantitation (LOQ) of the LF analyses were estimated by recording diluted solutions of standard SCs (5F-AKB48 and 5F-MDMB-PICA) of known concentrations. The LODs were evaluated at 0.33 and 0.25 mg/ml and the LOQs at 0.92 and 1.0 mg/ml for 5F-MDMB-PICA and 5F-AKB48, respectively, i.e. significantly lower than the concentrations measured in e-liquids (Table 4). We did not measure the LOQ of 5F-ADB which bears a fluorinated chain similar to 5F-AKB48, nor that of ADB-FUBINACA as the absence of multiplicity of the signal (singlet) enhances the S/N ratio of the resonance.

3.4.2. Quantitative analyses of SCs in e-liquids

The amounts of SCs in each e-liquid sample analysed in triplicate are reported in Table 4. For the HF ¹H NMR quantification, a characteristic signal could be integrated for each SC despite some overlaps. The ¹⁹F NMR quantification was straightforward as the NMR signals of the fluorinated SCs and the reference TFT are easily integrated at both HF and LF.

HF proton and fluorine NMR are recognized as quantitative methods in forensic [25, 31, 43] but have never been applied to e-liquids to our knowledge. The results reported in Table 4 show that, even if a matrix effect could bias the 1 H NMR quantification, data obtained by both techniques are in good agreement. Indeed, the linear regression equation of the HF 1 H versus 19 F values display a slope of 0.99 ± 0.03 (p-value = 0), a y-intercept of -0.34 ± 0.46 and a correlation coefficient of 0.98 (p-value 1.1×10^{-22}). The p-value of 0.55 given by the Wilcoxon signed-rank test demonstrates the absence of significant differences between these two data sets. The same observation can be done when comparing LF 19 F versus HF 19 F quantifications: a slope of 0.90 ± 0.03 (p-value = 0), a y-intercept of 1.22 ± 0.47 , a correlation coefficient of 0.97 (p-value = 1.6×10^{-21}) and a p-value of 0.74 for the Wilcoxon signed-rank

test. These results demonstrate that SCs in e-liquids can be successfully quantified by NMR and that an accurate ¹⁹F NMR quantification of fluorinated SCs can be implemented on a compact LF NMR spectrometer without requiring any specific standard compound.

As shown in Table 4, the 13 e-liquids contained variable concentrations of SCs: samples 1-3 contained the lowest amount of SCs with JWH-210 at around 1 mg/mL, whereas sample 13 contained more than 20 mg/mL of 5F-AKB48. The contents in 5F-MDMB-PICA, 5F-ADB, and ADB-FUBINACA ranged from \approx 7 to \approx 18 mg/ml in samples 4-12. To our knowledge, there are no published quantitative data on SCs in e-liquids for discussion of these results. However, even if the toxicological incidence is not the direct purpose of this study, we must highlight the high CB1/CB2 binding potency of the third generation of SCs which include the fluorinated SCs detected in this study [10, 15, 48, 49]. For instance, ADB-FUBINACA is a highly potent agonist of CB1 (EC50 = 1.2 nM) [15] by comparison to JWH 210 (EC50 = 25.3 nM) [49] or THC (EC50 = 172 nM) [15]. This high affinity of SCs for cannabinoids receptors along with their significant concentration in e-liquids obviously poses a threat to the health of consumers.

Conclusion

In this study, 13 commercial e-liquid formulations were screened by ¹H and ¹⁹F NMR. HF ¹H NMR allowed a holistic analysis of SCs and flavouring compounds in the e-liquids. In addition, the ability of benchtop LF ¹H and ¹⁹F NMR spectroscopy was evaluated. To the best of our knowledge, this is the first study dealing with the analysis of e-liquids containing SCs by LF NMR. Although LF ¹H NMR has shown limitations for quantitative applications, in particular due to (i) the huge difference in concentrations between the matrix and the active compounds, and (ii) signal overlaps, some characteristic signals could be distinguished. The interesting novelty of this study is that NMR and in particular LF ¹⁹F NMR can easily allow the quantification of fluorinated SCs in e-liquids without complex sample preparation and with no need for specific standard compounds. Thanks to the detection specificity of LF fluorine NMR, this proof-of-concept study opens the way to routine screening for the detection of illicit fluorinated SCs in various matrices such as e-liquids, herbal blends, or bath salts in forensic or customs laboratories.

Captions for Figures

- Figure 1. Proton NMR spectra of e-liquids (400 MHz, in MeOD): (A) Whole NMR spectrum
- of sample 3, (B) Overlap of 5 spectra of e-liquids containing different SCs (5.6-10.0 ppm).
- EM: ethyl maltol, V: vanillin, RK: raspberry ketone, Lin: linalool, AA: anisyl alcohol.
- Figure 2. Low-Field (60 MHz) proton NMR spectra of e-liquids: (A) Comparison of the ¹H
- NMR spectra of sample 8 recorded in MeOD on the high-field (400 MHz) and low-field (60
- 358 MHz) spectrometers, (B) Overlap of 5 spectra (60 MHz) with different SCs (5.6-10.0 ppm).
- 359 EM: ethyl maltol, V: vanillin, RK: Raspberry Ketone, Lin: Linalool.
- 360 **Figure 3.** ¹⁹F NMR spectra of e-liquid samples recorded in MeOD at 376 MHz (HF) and
- 361 56.17 MHz (LF).

362 **References**

- 363 [1] M.R. Peace, T.R. Baird, N. Smith, C.E. Wolf, J.L. Poklis, A. Poklis, Concentration of
- Nicotine and Glycols in 27 Electronic Cigarette Formulations, J Anal Toxicol 40(6) (2016)
- 365 403-7. 10.1093/jat/bkw037.
- 366 [2] A. Breland, E. Soule, A. Lopez, C. Ramoa, A. El-Hellani, T. Eissenberg, Electronic
- cigarettes: what are they and what do they do?, Ann Ny Acad Sci 1394(1) (2017) 5-30.
- 368 https://doi.org/10.1111/nyas.12977.
- 369 [3] G. Kaur, R. Pinkston, B. McLemore, W.C. Dorsey, S. Batra, Immunological and
- toxicological risk assessment of e-cigarettes, Eur Respir Rev 27(147) (2018).
- 371 10.1183/16000617.0119-2017.
- 372 [4] S.V. Cherian, A. Kumar, Y.M.R.M. Estrada, E-Cigarette or Vaping Product-Associated
- 373 Lung Injury: A Review, Am J Med 133(6) (2020) 657-663. 10.1016/j.amjmed.2020.02.004.
- 374 [5] M.R. Peace, J.W. Stone, J.L. Poklis, J.B. Turner, A. Poklis, Analysis of a Commercial
- Marijuana e-Cigarette Formulation, J Anal Toxicol 40(5) (2016) 374-8, 10.1093/jat/bkw021.
- 376 [6] J. Pourchez, V. Forest, E-cigarettes: from nicotine to cannabinoids, the French situation,
- 377 Lancet Resp Med 6(5) (2018) e16. 10.1016/S2213-2600(18)30069-9.
- 378 [7] M.R. Peace, R.I. Krakowiak, C.E. Wolf, A. Poklis, J.L. Poklis, Identification of MDMB-
- 379 FUBINACA in commercially available e-liquid formulations sold for use in electronic
- 380 cigarettes, Forensic Sci Int 271 (2017) 92-97. 10.1016/j.forsciint.2016.12.031.
- [8] J.L. Poklis, H.A. Mulder, M.R. Peace, The unexpected identification of the
- cannabimimetic, 5F-ADB, and dextromethorphan in commercially available cannabidiol e-
- 383 liquids, Forensic Sci Int 294 (2019) e25-e27. 10.1016/j.forsciint.2018.10.019.
- 384 [9] L. Hondebrink, A. Zwartsen, R.H.S. Westerink, Effect fingerprinting of new psychoactive
- substances (NPS): What can we learn from in vitro data?, Pharmacol Ther 182 (2018) 193-
- 386 224. 10.1016/j.pharmthera.2017.10.022.
- [10] C. Hess, C.T. Schoeder, T. Pillaiyar, B. Madea, C.E. Muller, Pharmacological evaluation
- of synthetic cannabinoids identified as constituents of spice, Forensic Toxicology 34(2)
- 389 (2016) 329-343. 10.1007/s11419-016-0320-2.
- 390 [11] K. Cohen, A.M. Weinstein, Synthetic and Non-synthetic Cannabinoid Drugs and Their
- 391 Adverse Effects-A Review From Public Health Prospective, Front Public Health 6 (2018)
- 392 162. 10.3389/fpubh.2018.00162.
- 393 [12] S.M. Gurney, K. Scott, S. Kacinko, B. Presley, B. Logan, Pharmacology, toxicology, and
- adverse effects of synthetic cannabinoid drugs, Forensic Sci Rev 26(1) (2014) 53-78. (no
- 395 DOI)
- 396 [13] A.M. Kemp, M.S. Clark, T. Dobbs, R. Galli, J. Sherman, R. Cox, Top 10 facts you need
- to know about synthetic cannabinoids: not so nice spice, The American journal of medicine
- 398 129(3) (2016) 240-244. 10.1016/j.amjmed.2015.10.008.

- 399 [14] D. Debruyne, R. Le Boisselier, Emerging drugs of abuse: current perspectives on
- 400 synthetic cannabinoids, Subst Abuse Rehabil 6 (2015) 113-29. 10.2147/SAR.S73586.
- 401 [15] S.D. Banister, M. Moir, J. Stuart, R.C. Kevin, K.E. Wood, M. Longworth, S.M.
- Wilkinson, C. Beinat, A.S. Buchanan, M. Glass, M. Connor, I.S. McGregor, M. Kassiou,
- 403 Pharmacology of Indole and Indazole Synthetic Cannabinoid Designer Drugs AB-
- 404 FUBINACA, ADB-FUBINACA, AB-PINACA, ADB-PINACA, 5F-AB-PINACA, 5F-ADB-
- 405 PINACA, ADBICA, and 5F-ADBICA, ACS Chem Neurosci 6(9) (2015) 1546-59.
- 406 10.1021/acschemneuro.5b00112.
- 407 [16] V. Shevyrin, V. Melkozerov, A. Nevero, O. Eltsov, Y. Shafran, Y. Morzherin, A.T.
- 408 Lebedev, Identification and analytical characteristics of synthetic cannabinoids with an
- indazole-3-carboxamide structure bearing a N-1-methoxycarbonylalkyl group, Anal Bioanal
- 410 Chem 407(21) (2015) 6301-15. 10.1007/s00216-015-8612-7.
- 411 [17] J. Wallgren, S. Vikingsson, A. Astrand, M. Josefsson, H. Green, J. Dahlen, X.Y. Wu, P.
- 412 Konradsson, Synthesis and identifications of potential metabolites as biomarkers of the
- 413 synthetic cannabinoid AKB-48, Tetrahedron 74(24) (2018) 2905-2913.
- 414 10.1016/j.tet.2018.04.026.
- 415 [18] N. Uchiyama, S. Matsuda, M. Kawamura, R. Kikura-Hanajiri, Y. Goda, Two new-type
- cannabimimetic quinolinyl carboxylates, QUPIC and QUCHIC, two new cannabimimetic
- carboxamide derivatives, ADB-FUBINACA and ADBICA, and five synthetic cannabinoids
- 418 detected with a thiophene derivative α-PVT and an opioid receptor agonist AH-7921
- 419 identified in illegal products, Forensic Toxicology 31(2) (2013) 223-240. 10.1007/s11419-
- 420 013-0182-9.
- 421 [19] R. Denooz, J.-C. Vanheugen, M. Frederich, P. de Tullio, C. Charlier, Identification and
- 422 Structural Elucidation of Four Cannabimimetic Compounds (RCS-4, AM-2201, JWH-203
- and JWH-210) in Seized Products, Journal of Analytical Toxicology 37(2) (2013) 56-63.
- 424 10.1093/jat/bks095.
- 425 [20] N. Uchiyama, M. Kawamura, R. Kikura-Hanajiri, Y. Goda, URB-754: a new class of
- designer drug and 12 synthetic cannabinoids detected in illegal products, Forensic science
- 427 international 227(1-3) (2013) 21-32. 10.1016/j.forsciint.2012.08.047.
- 428 [21] M. Bovens, C. Bissig, S.N. Staeheli, M. Poetzsch, B. Pfeiffer, T. Kraemer, Structural
- characterization of the new synthetic cannabinoids CUMYL-PINACA, 5F-CUMYL-
- 430 PINACA, CUMYL-4CN-BINACA, 5F-CUMYL-P7AICA and CUMYL-4CN-B7AICA,
- 431 Forensic Sci Int 281 (2017) 98-105. 10.1016/j.forsciint.2017.10.020.
- 432 [22] M.D.P. Risseeuw, P. Blanckaert, V. Coopman, S. Van Ouekelberghe, S. Van
- 433 Calenbergh, J. Cordonnier, Identification of a new tert-leucinate class synthetic cannabinoid
- in powder and "spice-like" herbal incenses: Methyl 2-[[1-(5-fluoropentyl)indole-3-
- carbonyl]amino]-3,3-dimethyl-butanoate (5F-MDMB-PICA), Forensic Sci Int 273 (2017) 45-
- 436 52. 10.1016/j.forsciint.2017.01.023.
- 437 [23] A.D. Lesiak, R.A. Musah, M.A. Domin, J.R. Shepard, DART MS as a Preliminary
- 438 Screening Method for "Herbal Incense": Chemical Analysis of Synthetic Cannabinoids,
- 439 Journal of forensic sciences 59(2) (2014) 337-343. 10.1111/1556-4029.12354.
- 440 [24] G. Assemat, F. Dubocq, S. Balayssac, C. Lamoureux, M. Malet-Martino, V. Gilard,
- Screening of "spice" herbal mixtures: From high-field to low-field proton NMR, Forensic Sci
- 442 Int 279 (2017) 88-95. 10.1016/j.forsciint.2017.08.006.
- 443 [25] S.J. Dunne, J.P. Rosengren-Holmberg, Quantification of synthetic cannabinoids in herbal
- smoking blends using NMR, Drug Test Anal 9(5) (2017) 734-743. 10.1002/dta.2032.
- [26] F. Fowler, B. Voyer, M. Marino, J. Finzel, M. Veltri, N.M. Wachter, L. Huang, Rapid
- screening and quantification of synthetic cannabinoids in herbal products with NMR
- spectroscopic methods, Anal Methods-Uk 7(18) (2015) 7907-7916. 10.1039/c5ay01754h.

- 448 [27] M.A. Marino, B. Voyer, R.B. Cody, A.J. Dane, M. Veltri, L. Huang, Rapid Identification
- of Synthetic Cannabinoids in Herbal Incenses with DART-MS and NMR, J Forensic Sci 61
- 450 Suppl 1 (2016) S82-91. 10.1111/1556-4029.12932.
- 451 [28] K.G. Shanks, G.S. Behonick, T. Dahn, A. Terrell, Identification of Novel Third-
- 452 Generation Synthetic Cannabinoids in Products by Ultra-Performance Liquid
- 453 Chromatography and Time-of-Flight Mass Spectrometry, Journal of Analytical Toxicology
- 454 37(8) (2013) 517-525. 10.1093/jat/bkt062.
- 455 [29] N.K. Burns, A.G. Theakstone, H. Zhu, L.A. O'Dell, J.R. Pearson, T.D. Ashton, F.M.
- 456 Pfeffer, X.A. Conlan, The identification of synthetic cannabinoids surface coated on herbal
- substrates using solid-state nuclear magnetic resonance spectroscopy, Anal Chim Acta 1104
- 458 (2020) 105-109. 10.1016/j.aca.2019.12.051.
- 459 [30] M.R. Peace, K.E. Butler, C.E. Wolf, J.L. Poklis, A. Poklis, Evaluation of Two
- 460 Commercially Available Cannabidiol Formulations for Use in Electronic Cigarettes, Frontiers
- 461 in Pharmacology 7(279) (2016). 10.3389/fphar.2016.00279.
- 462 [31] A.D.C. Santos, L.M. Dutra, L.R.A. Menezes, M.F.C. Santos, A. Barison, Forensic NMR
- spectroscopy: Just a beginning of a promising partnership, TrAC Trends in Analytical
- 464 Chemistry 107 (2018) 31-42. 10.1016/j.trac.2018.07.015.
- 465 [32] A.P. Sobolev, F. Thomas, J. Donarski, C. Ingallina, S. Circi, F.C. Marincola, D. Capitani,
- 466 L. Mannina, Use of NMR applications to tackle future food fraud issues, Trends Food Sci
- 467 Tech 91 (2019) 347-353. 10.1016/j.tifs.2019.07.035.
- 468 [33] T. Castaing-Cordier, V. Ladroue, F. Besacier, A. Bulete, D. Jacquemin, P. Giraudeau, J.
- 469 Farjon, High-field and benchtop NMR spectroscopy for the characterization of new
- psychoactive substances, Forensic Sci Int 321 (2021) 110718.
- 471 10.1016/j.forsciint.2021.110718.
- 472 [34] O. Apirakkan, A. Frinculescu, H. Denton, T. Shine, D. Cowan, V. Abbate, N. Frascione,
- 473 Isolation, detection and identification of synthetic cannabinoids in alternative formulations or
- dosage forms, Forensic Chemistry 18 (2020) 100227. 10.1016/j.forc.2020.100227.
- 475 [35] W. Jakes, A. Gerdova, M. Defernez, A. Watson, C. McCallum, E. Limer, I. Colquhoun,
- D. Williamson, E. Kemsley, Authentication of beef versus horse meat using 60 MHz 1H
- NMR spectroscopy, Food Chemistry 175 (2015) 1-9. 10.1016/j.foodchem.2014.11.110.
- 478 [36] T. Parker, E. Limer, A. Watson, M. Defernez, D. Williamson, E.K. Kemsley, 60 MHz
- 479 1H NMR spectroscopy for the analysis of edible oils, TrAC Trends in Analytical Chemistry
- 480 57 (2014) 147-158. 10.1016/j.trac.2014.02.006.
- 481 [37] G. Assemat, S. Balayssac, A. Gerdova, V. Gilard, C. Caillet, D. Williamson, M. Malet-
- 482 Martino, Benchtop low-field (1)H Nuclear Magnetic Resonance for detecting falsified
- 483 medicines, Talanta 196 (2019) 163-173. 10.1016/j.talanta.2018.12.005.
- 484 [38] G. Pagès, A. Gerdova, D. Williamson, V. Gilard, R. Martino, M. Malet-Martino,
- Evaluation of a benchtop cryogen-free low-field 1H NMR spectrometer for the analysis of
- 486 sexual enhancement and weight loss dietary supplements adulterated with pharmaceutical
- 487 substances, Analytical chemistry 86(23) (2014) 11897-11904. 10.1021/ac503699u.
- 489 [39] N. Wu, S. Balayssac, S. Danoun, M. Malet-Martino, V. Gilard, Chemometric Analysis of
- 490 Low-field 1H NMR Spectra for Unveiling Adulteration of Slimming Dietary Supplements by
- 491 Pharmaceutical Compounds, Molecules 25(5) (2020) 1193. 10.3390/molecules 25051193.
- 492 [40] W. Zhu, X. Wang, L. Chen, Rapid detection of peanut oil adulteration using low-field
- 493 nuclear magnetic resonance and chemometrics, Food Chem 216 (2017) 268-74.
- 494 10.1016/j.foodchem.2016.08.051.

488

- 495 [41] P.H.J. Keizers, F. Bakker, J. Ferreira, P.F.K. Wackers, D. van Kollenburg, E. van der Aa,
- 496 A. van Beers, Benchtop NMR spectroscopy in the analysis of substandard and falsified

- 497 medicines as well as illegal drugs, J Pharm Biomed Anal (2019) 112939.
- 498 10.1016/j.jpba.2019.112939.
- 499 [42] S.D. Banister, J. Stuart, R.C. Kevin, A. Edington, M. Longworth, S.M. Wilkinson, C.
- Beinat, A.S. Buchanan, D.E. Hibbs, M. Glass, Effects of bioisosteric fluorine in synthetic
- 501 cannabinoid designer drugs JWH-018, AM-2201, UR-144, XLR-11, PB-22, 5F-PB-22,
- 502 APICA, and STS-135, ACS chemical neuroscience 6(8) (2015) 1445-1458.
- 503 10.1021/acschemneuro.5b00107.
- 504 [43] H.A. Naqi, T.J. Woodman, S.M. Husbands, I.S. Blagbrough, 19F and 1H quantitative-
- NMR spectroscopic analysis of fluorinated third-generation synthetic cannabinoids, Anal
- 506 Methods-Uk 11(24) (2019) 3090-3100. 10.1039/c9ay00814d.
- 507 [44] G. Fardella, P. Barbetti, I. Chiappini, G. Grandolini, Quantitative analysis of
- 508 fluoroquinolones by 1H-and 19F-NMR spectroscopy, International journal of pharmaceutics
- 509 121(1) (1995) 123-127. 10.1016/0378-5173(95)00038-K.
- 510 [45] C. Simmler, J.G. Napolitano, J.B. McAlpine, S.-N. Chen, G.F. Pauli, Universal
- quantitative NMR analysis of complex natural samples, Current Opinion in Biotechnology 25
- 512 (2014) 51-59. 10.1016/j.copbio.2013.08.004.
- 513 [46] R. Martino, V. Gilard, F. Desmoulin, M. Malet-Martino, Fluorine-19 or phosphorus-31
- NMR spectroscopy: a suitable analytical technique for quantitative in vitro metabolic studies
- of fluorinated or phosphorylated drugs, J Pharm Biomed Anal 38(5) (2005) 871-91.
- 516 10.1016/j.jpba.2005.01.047.
- 517 [47] J.J. Nasr, S. Shalan, Validated 1H and 19F nuclear magnetic resonance for the
- quantitative determination of the hepatitis C antiviral drugs sofosbuvir, ledipasvir, and
- daclatasvir in tablet dosage forms, Microchemical Journal 152 (2020) 104437.
- 520 10.1016/j.microc.2019.104437.
- 521 [48] S.D. Banister, M. Longworth, R. Kevin, S. Sachdev, M. Santiago, J. Stuart, J.B. Mack,
- 522 M. Glass, I.S. McGregor, M. Connor, M. Kassiou, Pharmacology of Valinate and tert-
- Leucinate Synthetic Cannabinoids 5F-AMBICA, 5F-AMB, 5F-ADB, AMB-FUBINACA,
- MDMB-FUBINACA, MDMB-CHMICA, and Their Analogues, ACS Chem Neurosci 7(9)
- 525 (2016) 1241-54, 10.1021/acschemneuro.6b00137.
- 526 [49] A. Cannaert, J. Storme, F. Franz, V. Auwarter, C.P. Stove, Detection and Activity
- 527 Profiling of Synthetic Cannabinoids and Their Metabolites with a Newly Developed
- 528 Bioassay, Anal Chem 88(23) (2016) 11476-11485. 10.1021/acs.analchem.6b02600.

Figure 1




Figure 1




Figure 2




Figure 2




Figure 3




Table 1. High-field ¹H NMR (400 MHz) characteristic signals of synthetic cannabinoids detected in e-liquid samples recorded in MeOD. Signals in bold were used for quantification.

	5 6 4 7 7 2" 4" 5" 6' 5' 3'		22.7	5 6 4 7 2" 3" 5" 2' 1" 2' 1" 5"		5 6 4 7 2" 3" 5"		5 6 3" F 4 7 2" F NH 1"		5 6 7 2" 3" 5" 4" 5" 4" 5" 4" 5" 4" 5" 4" 5" 4" 5"					
	JWH-210				IDMB-P	ICA		5F-ADB			FUBINA	CA	5	5F-AKB48	
¹ H label ^a	δ (ppm)	M^b	J (Hz)	δ (ppm)	M^{b}	J (Hz)	δ (ppm)	M^{b}	J (Hz)	δ (ppm)	M^{b}	J (Hz)	δ (ppm)	M^{b}	J (Hz)
2	7.56	S		8.06	S										
4	8.31	dd	7.0/1.3	8.02	d	7.8	8.20	d	8.2	8.22	d	8.2	8.20	d	8.2
5	7.30-7.60	m		7.21	m		7.31	m		7.31	m		7.26	m	
6	7.30-7.60	m		7.27	m		7.48	m		7.44	ddd	8.6 7.0 1.1	7.44	ddd	8.6 6.9 1.1
7	7.30-7.60	m		7.49	d	8.2	7.68	d	8.6	7.61	d	8.5	7.62	d	8.5
2'	7.44-7.60	m			Overlap		Overlap			Overlap		2.20	m		
3'	7.44-7.60	m		/			/			/			2.13	br s	
4'	/			/			/			/			1.78	br s	
5'	8.04	d	8.4	/			/			/			/		
6'	7.30-7.60	m		/			/			/			/		
7'	7.30-7.60	m		/			/			/			/		
8'	8.18	d	8.6	/			/			/			/		
1"		verlap		4.26	t	7.1	4.54	t	7.0	5.72	S		4.48	t	7.0
2"	1.78	quin	7.5	1.92	quin	7.5	2.01	quin	7.5	7.31	m		1.98	quin	7.5
3''	Overlap		1.42	m		1.43	m		7.06	app t	8.8	1.42	m		
4''		verlap		1.69	m		1.72	m		/			1.72	m	
5"	0.82	t	7.0	4.40	dt	47.5/6.0	4.39	dt	47.5/6.0	/			4.39	dt	47.5/6.0
CH ₃ -CH ₂	1.42	t	7.5	/			/			/			/		
$(CH_3)_3$	/			1.09	S		1.08	S		1.09	S		/		

^a Proton signals mentioned as "overlap" or missing (CH₃-OCO-, CH₃-CH₂) could not be readily detected due to their overlap with other signals.

^b M: multiplicity; s: singlet; d: doublet; dd: doublet of doublet; t: triplet; dt: doublet of triplet; quin: quintuplet; m: multiplet; app t: apparent triplet; br s: broad signal.

Table 2. High-field ¹H NMR (400 MHz, MeOD) and GC-MS data of flavouring compounds detected in the e-liquids analyzed. Only detected ¹H NMR signals are reported.

Compound	Structure	¹H NMR		GC-MS		Detected in samples
		δ (ppm), multiplicity ^a (J (Hz)), attribution	Retention time (min)	Formula Experimental accurate mass m/z (mass error in ppm)	m/z of major characteristic fragment ions observed	
Ethyl maltol	5 OH CH ₃	8.00, d (5.5), H6 6.42, d (5.5), H5 2.74, q (7.6), CH ₂ 1.22, t (7.6), CH ₃	9.89	C ₇ H ₈ O ₃ 140.0467 (0.7)	55, 69, 71, 97, 139	7-12
Vanillin	O H 6 2 5 OH CH ₃	9.74, s, CHO 7.46, m, H6 7.44, s, H2 6.96, d (8.5), H5 3.92, s, -O-CH ₃	15.00	C ₈ H ₈ O ₃ 152.0462 (3.9)	81, 109, 123, 151	4-9, 11, 12
Raspberry ketone	HO 3' 2' CH ₃	7.00, d (8.5), H2' 6.69, d (8.5), H3' 2.11, s, CH ₃	19.18	$C_{10}H_{12}O_2$ 164.0834 (1.2)	107, 131, 149	8, 10, 11
Anisyl alcohol	OH 2 3 OCH3	6.88, d (8.7), H3 4.52, s, CH ₂	12.04	C ₈ H ₁₀ O ₂ 138.0676 (0.7)	77, 109, 121	4

Menthol	9 CH ₃ 10 H ₃ C 4 HO 6 CH ₃	1.61, m, H3 (1H) 0.92, d (7.1), CH ₃ 9 0.91, d (6.6), CH ₃ 7 0.84, m, H3 (1H) 0.79, d (6.9), CH ₃ 10	9.65	C ₁₀ H ₂₀ O 156.1519 (3)	71, 81, 95, 123, 138	13
Limonene	H ₃ C 2' CH ₃	4.68, m, H1' 1.62, broad s, CH ₃ 8	5.50	C ₁₀ H ₁₆ 136.1245 (1.5)	68, 79, 93	6
Linalool	8' H ₃ C 6 4 2	5.89, dd (17.4, 10.8), H2 5.01, dd (10.8, 1.6), H1 (1H) 1.65;1.58, two broad s, CH ₃ 8;8' 1.23, broad s, CH ₃ 9	7.19	C ₁₀ H ₁₈ O 154.1360 (2)	71, 80, 93, 121, 136	10, 11

^a s: singlet; d: doublet; dd: doublet of doublet; t: triplet; q: quadruplet; m: multiplet.

Table 3. T₁ values measured in e-liquid samples.

Cannabinoids	HF	LF NMR		
(sample number)	¹ H NMR	¹⁹ F NMR ^c	¹⁹ F NMR ^c	
	$T_1(s)$	$T_1(s)$	$T_1(s)$	
JWH-210 (3)	1.5 (H5')	-	-	
5F-MDMB-PICA (4)	0.5 (H2'')	0.13	0.17	
5F-ADB (7)	1.1 (H5)	0.14	0.18	
ADB-FUBINACA (11)	0.4 (H1'')	0.09	0.14	
5F-AKB48 (13)	2.0 (H5'')	0.14	0.16	
TSP^a	2.7	-	-	
$\mathrm{TFT}^{\mathrm{b}}$	-	0.20	0.27	

^a Mean values from samples 3, 4, 7, 11 and 13.
^b Mean values from samples 4, 7, 11 and 13.
^c ¹⁹F NMR measurements were done on samples containing the relaxation agent Cr(acac)₃.

Table 4. Amounts of synthetic cannabinoids (SC) in 13 e-liquids measured with high-field ¹H and ¹⁹F NMR and low-field ¹⁹F NMR. All experiments were done in triplicate from sampling the e-liquid to processing the NMR spectrum.

		Aı	SD		
		HF N	HF NMR		
Sample	Cannabinoid	¹ H NMR	¹⁹ F NMR	¹⁹ F NMR	
1	JWH-210	0.9 ± 0.1	_a	-	
2		0.9 ± 0.1	-	-	
3		1.1 ± 0.1	-	-	
4	5F-MDMB-PICA	18.4 ± 0.3	19.3 ± 0.3	18.4 ± 0.3	
5		17.2 ± 0.3	18.0 ± 0.2	17.3 ± 0.3	
6		15.2 ± 0.6	14.9 ± 0.4	14.3 ± 1.3	
7	5F-ADB	15.6 ± 0.3	15.6 ± 0.3	16.0 ± 0.7	
8		7.3 ± 0.1	7.7 ± 0.3	8.0 ± 0.2	
9		10.9 ± 0.4	11.7 ± 0.5	11.8 ± 0.6	
10	ADB- FUBINACA	15.1 ± 0.3	15.4 ± 0.1	15.2 ± 0.5	
11		16.6 ± 0.4	17.8 ± 0.5	17.5 ± 0.5	
12		15.7 ± 0.3	16.3 ± 0.3	15.0 ± 0.2	
13	5F-AKB48	24.1 ± 0.1	24.1 ± 0.1	22.7 ± 0.5	

^a JWH-210 is not fluorinated.

