

Quantification of curium isotopes in environmental samples: drawbacks, speciation and specific tracer

Azza Habibi, Roselyne Le Corre, Christele Wampach Aubert, Sandrine Fleury, Patricia Hennequet, Celine Augeray, Kevin Galliez, Michelle Agarande

▶ To cite this version:

Azza Habibi, Roselyne Le Corre, Christele Wampach Aubert, Sandrine Fleury, Patricia Hennequet, et al.. Quantification of curium isotopes in environmental samples: drawbacks, speciation and specific tracer. Journal of Radioanalytical and Nuclear Chemistry, 2021, 329 (2), pp.545-554. 10.1007/s10967-021-07751-7. hal-03368652

HAL Id: hal-03368652 https://hal.science/hal-03368652v1

Submitted on 6 Oct 2021

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

11 12	Quantification of curium isotopes in environmental samples: drawbacks, speciation and specific tracer
13 14	A. HABIBI ¹ , R. LE CORRE ¹ , C. WAMPACH AUBERT ¹ , S. FLEURY ² , P. HENNEQUET ² , C. AUGERAY ¹ , K. GALLIEZ ² , M. AGARANDE ³
15	¹ IRSN/PSE-ENV/SAME/LERCA, 31 rue de l'écluse, 78116, Le Vésinet (France)
16	² IRSN/PSE-ENV/SAME/LMN, 31 rue de l'écluse, 78116, Le Vésinet (France)
17	³ IRSN/PSE-ENV/SAME, 31 rue de l'écluse, 78116, Le Vésinet (France)

18 Keywords

19 Curium, americium, environmental samples, speciation, ²⁴⁸Cm

20 Abstract

In this study, drawbacks of the use of ²⁴³Am as a tracer in order to quantify curium isotopes during a complex chemical procedure were highlighted. These drawbacks were investigated through a speciation computation simulation based on the compilation of over 50 stability constants. A new curium standard was then characterized and used to quantify a wide ²⁴⁴Cm activities range in environmental samples. Obtained results were in excellent agreement with spiked activities and more consistent than those obtained when using ²⁴³Am as tracer.

27 Introduction

Curium isotopes are mostly produced by irradiation of uranium and/or plutonium oxides with neutrons in nuclear power reactors [1]. In fact, during the nuclear fuel cycle, the production of transuranic elements with mass numbers exceeding 240 by multiple neutron capture of uranium and plutonium can be very substantial (e.g. ~ 20 g.tonne⁻¹ for curium isotopes) [2]. Over twenty isotopes of curium are known with mass numbers ranging from 232 to 252, all of these isotopes are radioactive and primarily α -particle emitters, with ²⁴²Cm and ²⁴⁴Cm (respective half-lives of 163 days and 18.1 years) accounting for 90 % of those produced in the nuclear fuel cycle [2]. Indeed, in the fallout following the Chernobyl accident, 242 Cm was the dominating alphaparticle emitter (14 times higher than $^{239+240}$ Pu) [2]. Furthermore, the high specific activity of 242 Cm and 244 Cm (> 10¹² Bq/g) has mainly been considered as an issue because it results in significant contribution to the radiotoxicity of nuclear waste. Regulatory and non-regulatory releases of curium in the environment, especially near reprocessing power plants such as La Hague (France), presents a considerable public health concern [1]. Curium isotopes quantification in environmental samples is then fundamental to monitor this releases.

42 Reliable detection and quantification of curium isotopes in complex environmental matrices require extensive removal of chemical and radiometric interferences. The procedure yield must 43 be determined using an elemental or isotopic analog. The common yield monitor of curium 44 separation is ²⁴³Am since considered as a chemical analog. The use of this isotope is therefore 45 known for its drawbacks. A bias could indeed be observed between americium and curium 46 behaviors [3-7]. This difference generates bias between chemical yields of americium and 47 curium which leads to inexact results. To understand the observed bias and behaviors of 48 americium and curium, their speciation during a complex purification procedure was 49 investigated in this study by computer simulation. This simulation was executed using the 50 MEDUSA software and the associated HYDRA equilibrium database, compiled with published 51 data[8].52

To avoid the risk of a bias between the chemical yields and therefore an incorrect curium isotopes activities determination in environmental samples, an alternative approach using ²⁴⁸Cm as tracer to enable a determination based on the isotopic dilution technique was studied. The concentration and isotopic composition of a commercial curium solution with ²⁴⁸Cm enrichment of 95.82 % (w/w) was standardized by alpha spectrometry. Various environmental samples, each spiked with known amounts of ²⁴⁸Cm and ²⁴⁴Cm were processed following our radiochemical procedure involving alpha spectrometry counting. Experimental results and method performances observed when analyzing various types of environmental samples (water, vegetables and sediment) are presented. A detailed methodology including curium and americium speciation, ²⁴⁸Cm solution qualification, interferences corrections, results and uncertainties calculation are presented and discussed.

64 **Reagents and instrumentation**

65 **Reagents**

- 66 All used reagents were prepared from analytical-grade powders and analytical reagent-grade
- 67 for acids. Extraction resins were from Triskem, Int. (FR) (TEVA resin and TRU resin). Anionic
- resin was purchased from BIO-RAD (AG^{\mathbb{R}} 1X4 50-100 mesh).
- 69 Tracer solutions were conserved in acidic media $(1 4 \text{ M of HNO}_3)$ and purchased from the
- 70 NIST, NPL, AEA Technology and ORNL for ²⁴²Pu (SRM 4334j), ²⁴³Am (A14063), ²⁴⁴Cm
- 71 (CLP 10010) and ²⁴⁸Cm (8R103CMPROD), respectively.

72 Instrumentation

Alpha source filters are counted by alpha spectrometry for approximately 69 h.

74 Experimental procedure

75 The major steps of the experimental procedure are summarized in figure 1.



Figure 1: experimental procedure for americium and curium quantification in environmental
 samples

The first step of the experimental procedure consists on the addition of tracers (²⁴²Pu and ²⁴³Am)
to 5 g of solid samples (e.g. soil, sediment, vegetables...) or 500 mL of liquid samples (e.g.
drinking water, surface water...).

When organic solid samples are analyzed (e.g. vegetables, fish, seaweed...), an acidic leaching
involving a mixture of concentrated nitric acid and hydrochloric acid is performed.

Furthermore, when soil or sediment samples are analyzed, the leaching is performed using 9 M
hydrochloric acid.

After the elimination, by centrifugation, of the undissolved residue, plutonium and americium/curium phosphates are precipitated by adding 10 mL of calcium phosphates and ammonium hydroxide (20 %) to adjust the pH at 8-9. Once the precipitate is centrifuged, it's dissolved by adding 100 mL of 7 M HNO₃. Valence adjustment is then performed with $\sim 200 \text{ mg of NaNO}_2$.

Plutonium is extracted from the sample using a column containing ~ 5 g of AG[®] 1-X4 50 – 100 91 92 mesh (Anion Exchange Resin) previously washed with H₂O and pre-conditioned with 30 mL of 7 M HNO₃. The sample is loaded through the column. The resin is then rinsed with 2 x 30 93 mL and 40 mL of 7 M HNO₃. The loaded sample and the rinsing fractions are collected to be 94 95 treated afterward since containing americium and curium (fraction 2). The column is then rinsed with 2 x 30 mL of 9 M HCl before plutonium elution with 2 x 30 mL of a solution of 0.2 M of 96 hydroxylamine Hydrochloride. The elution fraction is evaporated to dryness and the residue is 97 dissolved with 100 mL of 1 M HNO₃. Plutonium is precipitated by adding 5 mL of iron (III) 98 nitrate, previously prepared by dissolving 30 g of Fe(NO₃)₃ in 1 L of 0.01 M of HNO₃, and 99 ammonium hydroxide (20 %) to adjust the pH at least at 9. The precipitate is centrifuged and 100 dissolved with 30 mL of 7 M HNO₃. 10 mL of 750 g/L of Al(NO₃)₃. 9H₂O and ~ 100 mg of 101 NaNO₂ are then added to the previously dissolved precipitate. The sample is purified from the 102 remaining uranium and thorium with a TEVATM resin. After the conditioning with 20 mL of 3 103 M HNO₃ and the sample loading, the TEVATM resin is rinsed with 20 mL of 3 M HNO₃, 20 mL 104 of 0.5 M HNO₃ to eliminate uranium, 50 mL of 6 M HCl to eliminate thorium and again with 105 106 20 mL of 3 M HNO₃. Plutonium is finally eluted with 2 x 30 mL of a solution of 0.2 M of hydroxylamine Hydrochloride. The elution fraction is evaporated to dryness. 107

To the fraction containing americium and curium (fraction 2), ~ 100 mg of ascorbic acid, 3 mL 108 of the previously used solution of iron (III) nitrate and 12 g of oxalic acid are added. The 109 solution volume is then extended to 200 mL and the pH is adjusted to 1.48 - 1.5 with 20 % 110 ammonium hydroxide. The precipitate is then collected by centrifugation, dissolved with 20 111 mL of 7 M HNO₃ and heated to dryness in order to convert the oxalate complexes on carbonate 112 complexes. Carbonate complexes are then dissolved with 20 mL of 1 M HNO₃. Americium and 113 curium are extracted using a TRUTM resin. The sample is loaded through a prepacked column 114 purchased from Triskem international. The resin is rinsed with 2 x 10 mL of 1 M HNO₃ and 2 115 mL of 9 M HCl. The loading and rinsing fractions are discarded. Americium and curium are 116 eluted from the TRU resin with 10 mL of 4 M HCl. The elution fraction is heated to dryness 117 and the residue is dissolved with 20 mL of a freshly prepared solution containing 0.1 M formic 118 acid and 2 M ammonium thiocyanate. The sample is loaded and the TEVATM resin is rinsed 119 120 with 20 mL of 0.1 M formic acid and 2 M ammonium thiocyanate and 30 mL of 0.1 M formic acid and 1 M ammonium thiocyanate in order to remove the remaining lanthanides. Americium 121 and curium are eluted from the TEVATM resin with 30 mL of 0.25 M HCl. The elution fraction 122 is evaporated to dryness. 123

Both residues containing plutonium and americium/curium are dissolved with 50 mL of 0.2 M HCl and lanthanum fluoride microprecipitation is performed by adding 1 mL of 0.29 g/L lanthanum oxide, prepared in 0.2 M HCl, and 1 mL of 40 % hydrofluoric acid. pH is then adjusted to 1.89 - 1.91.

Quantification is based on isotopic dilution technique where ²⁴²Pu tracer is used for ²³⁸Pu and ²³⁹⁺²⁴⁰Pu quantification and ²⁴³Am tracer used for ²⁴¹Am, ²⁴²Cm and ²⁴³⁺²⁴⁴Cm quantification. Indeed, americium and curium have similar chemical properties in solution as a result of predominantly trivalent oxidation state with comparable ionic radii [9] and their extraction behavior is often perceived as indistinguishable [10]. Some differences were nonetheless observed between americium and curium behaviors. These differences were then exploited in
chromatographic systems to achieve separations of Cm from Am [3–7]. Due to these
differences, the ISO 13167:2016 [11] recommended, if curium isotopes are needed in drinking
water and ²⁴³Am is used as tracer, the determination of a correction factor (CF) using reference
materials or spiked water samples.

$$CF = \frac{Cm \text{ isotope theoretical activity}}{Cm \text{ isotope calculated activity with 243}_{Am} \text{ tracer}} \quad Eq. (1)$$

139 The calculated correction factor (Eq. (1)) has then to be applied to curium isotopes measured140 activities.

An attempt to experimentally determine the correction factor was done by our team. This determination was done using 11 spiked water samples and samples from proficiency tests with ²⁴⁴Cm activities ranging from 2 mBq.L⁻¹ to 10 Bq.mL⁻¹. The calculated correction factors are presented in figure 2.







147 The presented results were used to calculate the average correction factor which is 148 1.74 ± 15 %, the standard deviation is however relatively high (~ 0.68).

149 In order to determine the accuracy of the results when the determined correction factor (CF) is applied, 11 water samples, containing known activities of ²⁴⁴Cm, were analyzed using the 150 procedure described in figure 1 and ²⁴³Am as a tracer. The standardized deviations (SD) (Eq. 151 (2)) and relative deviations (RD) (Eq. (3)) between theoretical activities and experimentally 152 determined ones when using the correction factor (e.g. 1.74) are then calculated. The SD and 153 RD have to be, respectively, lower than 1 and 15 % in order to consider the experimental activity 154 of ²⁴⁴Cm acceptable. The obtained calculated SD and RD [6], when the determined correction 155 factor is applied, were compared to those obtained without the application of the correction 156 157 factor (table 1).

158
$$SD = \frac{|\text{calculated } 244_{Cm} \text{ activity} - \text{theoretical } 244_{Cm} \text{ activity}|}{\sqrt{(\text{uncertainty of the calculated activity}^2 + \text{uncertainty of the theoretical activity}^2)}} Eq. (2)$$

$$RD = \frac{|\text{calculated } 244_{Cm} \text{ activity-theoretical } 244_{Cm} \text{ activity}|}{\text{theoretical } 244_{Cm} \text{ activity}} \times 100 \quad \text{Eq. (3)}$$

161

Sample	Theoretical ²⁴⁴ Cm activities, Bq/L (k=2)	Calculated ²⁴⁴ Cm activities, Bq/L (k=2)	SD of ²⁴⁴ Cm	SD of ²⁴⁴ Cm (With correction factor)	RD of ²⁴⁴ Cm	RD of ²⁴⁴ Cm (With correction factor)
1	14.6 ± 2.2	25.4 ± 5.4	7.96	2.22	62 %	34 %
2	213.0 ± 19.2	369.8 ± 65.8	6.33	0.34	46 %	6 %
3	2.2 ± 0.3	3.9 ± 0.8	5.30	0.32	46 %	6 %
4	6.8 ± 0.7	11.9 ± 2.2	5.75	0.50	37 %	10 %
5	20.2 ± 4.4	35.0 ± 9.4	4.55	0.82	53 %	18 %
6	106. 7 ± 13.9	185.3 ± 37.3	4.39	0.14	41 %	3 %
7	2.4 ± 0.3	4.2 ± 0.9	2.93	1.02	28 %	26 %
8	7.8 ± 1.0	13.5 ± 2.7	6.38	0.81	51 %	14 %
9	197.1 ± 23.7	342.3 ± 66.7	1.46	1.57	16 %	46 %
10	9.1 ± 1.2	15.7 ± 3.2	4.53	0.07	42 %	2 %
11	2.4 ± 0.3	4.2 ± 0.8	0.84	1.97	8%	59 %

correction factor in water samples

Although a downward trend of the standardized deviations and the relative deviations is observed when the correction factor is applied, some calculated ²⁴⁴Cm activities are not acceptable since the associated SD is higher than 1 and/or the associated RD is higher than 15 %. For the sample n°11 the SD and the RD even increased after the correction. The use of ²⁴³Am when curium isotopes are analyzed could generate a bias even when a correction factor is experimentally determined and applied.

In order to understand the physico-chemical causes of this bias during the analysis of spiked aqueous solution, Cm and Am speciation, when using the previously described protocol, was computed. This speciation was performed using the HYDRA database associated with the MEDUSA software. Missing species and stability constants were added to complete the database. Activity coefficients were calculated by using Davies equations [14]. This speciation study was performed using a total of 50 stability constants. During the speciation computation, two hypotheses were made. The first hypothesis was that reactions between and with not studied ions were neglected. The second hypothesis was that 243 Am and 244 Cm concentration is 10^{-12} M (respectively, 2.43 x 10^{-10} g/L and 2.47 x 10^{-10} g/L) through the entire protocol.

The first studied step was the calcium phosphates co-precipitation (figures 3 and 4) [15–17]. This co-precipitation is performed at $pH \ge 8$. According to the previously described procedure, during this step, PO_4^{3-} concentration is ~ 30 mM and Ca^{2+} concentration is 45 mM. At pH higher than 8, the present calcium phosphate complex is $Ca_5(PO_4)_3OH_{(s)}$. Figures 3 and 4 show the fractions of Am and Cm complexes during the co-precipitation step.





184 *Figure 3: speciation of americium during the calcium phosphates co-precipitation step*



185

186

Figure 4: speciation of curium during the calcium phosphates co-precipitation step

Figures 3 and 4 show that at a pH between 8 and 10, Am and Cm phosphate complexes are AmPO₄. $x H_2O_{(s)}$ and CmPO₄. $x H_2O_{(s)}$. At pH higher than 10, Cm(OH)₂⁺ becomes a dominate species and its fraction reach 100 % at pH 11.2. In this range of pH, almost all Am is precipitated under the form of Am(OH)₃ while Cm(OH)₃ remains a minor species until pH 13.5. When a calcium phosphates co-precipitation of Am and Cm is performed, the pH has a significant influence on Am and Cm forms which are not necessarily the same and so can generate a bias.

The second studied step is the precipitated complexes dissolution using 7 M HNO₃ [18–20]. 193 During this step we suppose that all the present Am³⁺ and Cm³⁺ are under the forms AmPO₄. x 194 $H_2O_{(s)}$ and CmPO₄. x $H_2O_{(s)}$. However, since during the previous step, 100 % of Ca is 195 precipitated under the form $Ca_5(PO_4)_3OH_{(s)}$, Ca^{2+} and PO_4^{3-} concentrations are respectively at 196 45 mM and 30 mM during this step. NO₃⁻ concentration is 7 M during this step. The aim of 197 198 this step is to convert the solid AmPO₄. $x H_2O_{(s)}$ and CmPO₄. $x H_2O_{(s)}$ to aqueous complexes. Figures 5 and 6 present the speciation of Am³⁺ and Cm³⁺ during this step as a function of the 199 200 pH.



201

202

Figure 5: speciation of americium during the phosphates dissolution step



203

204

Figure 6: speciation of curium during the phosphates dissolution step

Figures 5 and 6 shows that during the phosphate complexes dissolution with 100 mL 7 M HNO₃ (pH < 1), 100 % of curium is under the form of aqueous nitrate complexes. On the other hand, ~ 6 % of Am is not under the form of aqueous nitrate complexes but remains under the form

- AmH₂PO₄²⁺. This aqueous americium complex fraction will not have the same behavior during the separation step on the AG[®] 1-X4 resin.
- The third studied step was the calcium oxalates co-precipitation [20, 21]. Figures 7 and 8 show the speciation of americium and curium in function of the pH.
- During this step, NO_3^- , Ca^{2+} and PO_4^{3-} concentrations are still at respectively 7 M, 45 mM and 30 mM. Moreover, according to the previously described procedure, oxalic acid (OxH₂) concentration is at 0.7 M.



216 *Figure 7: speciation of americium during the calcium oxalates co-precipitation step*



217

218 *Figure 8: speciation of americium during the calcium oxalates co-precipitation step*

For pH \geq 1, 100 % of the present calcium is precipitated as Ca(ox).H₂O_(s). According to figures 7 and 8, at pH 1, 97 % of Cm is precipitated as Cm₂(ox)₃. 10H₂O_(s) while only 83.2 % of Am is precipitated as Am₂(ox)₃. 7H₂O_(s). Once again, Am and Cm behaviors present a slight difference. Cm and Am remaining in solution are under the form of aqueous nitrate complexes. Furthermore at pH 1.5, 99.9 % of both Cm and Am precipitate as oxalate complexes. The pH adjustment step is then very important and a slight variation could influence Am and Cm forms.

After the recovery of the carrier, oxalate complexes were converted to carbonate complexes with 7 M HNO₃ under heating [20, 22, 23]. Speciation of americium and curium during this step are presented in figures 9 and 10. Since during the previous step, 100 % of Ca²⁺ is precipitated under the form of CaC₂O₄(s), Ca²⁺ concentration is still at 45 mM during this step. At [NO₃⁻] = 7 M, ~ 99 % of curium and americium are present as aqueous nitrate complexes. In contrary to the previously steps, the dissolution of the carbonate complexes carrier does not imply any differences between americium and curium behaviors.





Figure 9: speciation of americium during the carbonates dissolution step



Figure 10: speciation of curium during the carbonates dissolution step

Speciation of americium and curium during the separation using the TRUTM resin is not possible 236 without a perfect knowledge of the resin properties (e.g. quantity of coated CMPO, packing 237 procedure...). However, behaviors of americium and curium on the TRUTM resin were studied 238 by Gharibyan et al.[10]. In this study, minor variations between Cm and Am behaviors at 1 M 239 HNO₃ were described. This medium is the most common one used for the sample loading when 240 analyzing Am and Cm. The authors described batch studies and deduced a separation factor of 241 Am from Cm of 1.45. This type of behavior, where a heavier neighboring element has a lower 242 retention factor than its lighter neighbor, is also observed in the lanthanide series between Eu 243 (III) and Gd (III) which are isoelectronic to Am (III) and Cm (III), respectively. This separation 244 factor was then confirmed by column studies and other publications [24]. 245

246 In order to verify if this behavior difference is impacting in our case, several water samples spiked with Am and Cm were analyzed using TRUTM resin. Conditions of the separation were 247 the same as our routine protocol (e.g. 2 mL pre-packed columns, a flow of 1.5 mL/min 248 249 controlled by a peristaltic pump). Elution fractions (4 M HCl) were analyzed by ICP-MS in order to prevent any behavior difference due to the lanthanum fluoride co-precipitation step. 250 Figure 11 presents the obtained results. These results highlight that, within the uncertainties, 251 extraction results obtained for Am and Cm are similar. A slight trend of higher Am recoveries 252 is nonetheless observed. Since both actinides are not retained over the entire HCl concentration 253 range [24, 25], this difference is obviously due to a difference between retention factors during 254 the sample loading in 1 M HNO₃. 255





Figure 11: Am and Cm extracted fraction with the TRU^{TM} resin

258 This speciation computation study points up a noticeable difference in the hydrolysis constants of Am and Cm. Moreover, $AmH_2PO_4^{2+}$ seems to be more stable than $CmH_2PO_4^{2+}$ and a small 259 difference between formation's constants of Cm₂(Ox)₃.10 H₂O and Am₂(Ox)₃.7 H₂O was 260 observed. The literature described a higher affinity of americium to the TRUTM resin [10]. Cm 261 and Am behaviors are indistinguishable during all the other steps of our protocol. In fact, in 262 spite of a slight difference between extractions of Am and Cm by ammonium 263 dinonylnaphthalene sulphonate in moderate thiocyanate concentration solutions [26], behaviors 264 of Am and Cm during the separation using TEVATM resin were similar. Finally, since 265 americium and curium complexes are not retained on the AG® 1-X4 resin over the entire HNO3 266 concentration range, any difference between their behaviors during this step can be dismissed 267 268 [27].

When analyzing environmental samples, difference between Am and Cm behaviors could be more noticeable and not predictable. Due to the importance of curium analysis in some of the environmental samples (e.g. samples from the surrounding of reprocessing plants), the quantification has to be reliable. The only known standard containing ²⁴⁸Cm was then bought from the Oak Ridge National Laboratory. Activity concentrations of this standard are not 274 certified and associated uncertainties are not communicated. The total ²⁴⁸Cm content was 275 however measured to be 6.6 μ g (~ 1012 Bq) and the total volume of the solution was 5 mL at 276 2.0 M HCl. Table 2 presents the announced isotopic composition of the material.

277

Table 2: certified isotopic composition of the material

Isotope	Atom (%)	Weight (%)
²⁴⁴ Cm	< 0.02	< 0.02
²⁴⁵ Cm	0.06	0.06
²⁴⁶ Cm	4.12	4.09
²⁴⁷ Cm	0.02	0.02
²⁴⁸ Cm	95.78	95.82

278

Figure 12 presents the alpha spectra of the material. This spectra was obtained by counting a

lanthanum fluoride precipitate of ~ 0.02 Bq of the material during 69h.





Figure 12: alpha spectra of the used standard

Alpha spectrometry highlighted that despite the fact ²⁴⁴Cm presence is low, it represents 5.4 % of ²⁴⁸Cm net count. ²⁴⁴Cm presence in the radioactive material has then to be taken into account. However, ²⁴³Am emits at 5275.3 keV with an intensity of 93.2 % while ²⁴⁵Cm and ²⁴⁶Cm emit at respectively 5361.8 keV (93.2 %) and 5387.5 keV (86.74 %). The addition of both tracers (²⁴³Am and ²⁴⁸Cm) is then not possible and two replicates of the sample are needed in order to analyze americium and curium.

In order to verify activity concentrations of curium isotopes and determine uncertainties associated to each concentration, a diluted fraction of the purchased radioactive material was doped with a known amount of a certified ²⁴⁴Cm tracer solution and analyzed. This ²⁴⁴Cm tracer solution was previously analyzed and no significant ²⁴⁸Cm activity concentration was detected above our detection limit. Table 3 shows expected and measured activities concentrations of ²⁴⁸Cm and ²⁴⁵⁺²⁴⁶Cm.

295

Table 3: expected and measured activity concentrations

	expected activit	y concentrations	Measured activity concentrations		
	Activity (Bq.g ⁻¹)	Uncertainty (%)	Activity (Bq.g ⁻¹)	Uncertainty (%)	
²⁴⁸ Cm	6.29×10^{-2}	-	6.26×10^{-2}	7.52	
^{245 + 246} Cm	2.29×10^{-1}	-	2.27×10^{-1}	6.81	

296

Results presented in table 3 show that the measured activity concentrations are consistent with the calculated ones. Furthermore, this characterization step allowed the determinations of associated uncertainties.

In order to validate the use of the newly purchased ²⁴⁸Cm tracer solution for curium isotopes
quantification, water, soil and tea samples were spiked to be analyzed. In fact, although a few

environmental reference materials are containing curium, none has been reliably certified. Tea
 and soil samples were from previous proficiency tests intended to americium isotopes analysis.
 Water, soil and tea samples were analyzed before the spike and no significant ²⁴⁴Cm activity
 concentration was detected above our detection limit.

The quantification was based on the isotopic dilution technique, where the activity of each isotope is calculated using equation 4 [11].

Where A_T is the added activity of the tracer (²⁴⁸Cm), N_{net} are counts of the studied isotope (²⁴⁴Cm) corrected from the background and from the impurity in the tracer solution, N_{netT} are the counts of the tracer (²⁴⁸Cm), m is the sample quantity, E and E_T are respectively emitting rates of the studied isotope and the tracer.

Table 4 presents spiked and measured activities of ²⁴⁴Cm in tea, soil and water samples.

Spiked activity, ²⁴⁴ Cm (k=2)			Experimental activity, ²⁴⁴ Cm (k=2)			Standard deviation	Relative deviation
Tea samples (Bq/kg, ashes)							
0,29	±	0,01	0,32	±	0,09	0,34	11 %
0,28	±	0,01	0,28	±	0,08	0,05	1 %
0,91	±	0,05	0,90	±	0,18	0,04	1 %
0,91	±	0,05	0,92	±	0,17	0,06	1 %
4,28	±	0,21	4,32	±	0,78	0,05	1 %
4,11	±	0,21	4,07	±	0,81	0,05	1 %
	S	oil samples	(Bq/kg, ash	es)			
0,29	±	0,01	0,31	±	0,09	0,31	10 %
0,91	±	0,05	0,82	±	0,17	0,47	9 %
4,18	±	0,21	3,68	±	0,74	0,66	12 %
5,00	±	0,25	4,88	±	0,49	0,23	2 %
1,20	±	0,06	1,15	±	0,10	0,42	4 %
1,20	±	0,06	1,09	±	0,10	0,96	9 %
1,20	±	0,06	1,12	±	0,10	0,68	7 %
0,50	±	0,03	0,48	±	0,05	0,34	4 %
0,50	±	0,03	0,45	±	0,05	0,90	10 %
0,50	±	0,03	0,52	±	0,05	0,35	4 %
0,20	±	0,01	0,18	±	0,03	0,69	10 %
0,20	±	0,01	0,23	±	0,03	0,89	14 %

Table 4: spiked and measured activity concentrations of ²⁴⁴Cm in tea, soil and water samples

315

According to results obtained when analyzing tea, soil and water samples, the use of ²⁴⁸Cm as a tracer allows having excellent compatibility between spiked and measured activities even when analyzing complex matrices such as soil samples. These performances are characterized by standard deviations lower than 1 and relative deviations lower than 15 %. Results presented in table 4 are more reliable and accurate since activity concentrations determination is based on an isotopic tracer.

322 Conclusion

Curium origins in environmental samples are mostly due to regulatory releases from fuel 323 reprocessing plants (e.g. La Hague, Sellafield...) and fallout following major nuclear accidents. 324 Quantification of curium isotopes in environmental samples is frequently based on isotopic 325 dilution using ²⁴³Am as a tracer. This procedure shows a few drawbacks since it generates a 326 327 bias between americium and curium behaviors. The speciation computation studied in this work permitted to conclude that this bias is mostly due to differences in the hydrolysis constants of 328 americium and curium, in the formation constants of Cm₂(Ox)₃. 10 H₂O and Am₂(Ox)₃. 7 H₂O, 329 in the stability constants of AmH₂PO₄²⁺ and CmH₂PO₄²⁺ and in affinities of americium and 330 curium for the TRUTM resin in 1 M HNO₃ solution. In order to prevent the observed bias in 331 quantified curium activities, a new standard was reliably characterized and used to quantify 332 curium isotopes in spiked environmental samples (e.g. water, tea and soil samples). The new 333 procedure is based on isotopic dilution using ²⁴⁸Cm as tracer. Excellent results were obtained 334 335 for the analyzed samples. The characterization and the use of a specific curium tracer, even when it's not reliably certified, to quantify curium isotopes is then recommended since it allows 336 having more consistent results. 337

338 **Bibliography**

- Schneider DL, Livingston HD (1984) Measurement of curium in marine samples. Nucl Instrum Methods Phys Res 223:510–516. https://doi.org/10.1016/0167-5087(84)90700-2
- Holm E, Roos P, Aarkrog A, et al (2002) Curium isotopes in Chernobyl fallout. J
 Radioanal Nucl Chem 252:211–214. https://doi.org/10.1023/A:1015786431984
- 343 3. Moore FL (1971) New method for separation of americium from curium and associated
 add elements in the zirconium phosphate-nitric acid system. Anal Chem 43:487–489
- Myasoedov BF (1994) New methods of transplutonium element separation and determination. J Alloys Compd 213–214:290–299. https://doi.org/10.1016/0925-8388(94)90918-0
- Buijs K, Muller W, Reul J, Toussaint JC (1976) The separation and purification of
 americium and curium on the multigram scale. Proc Mosc Symp Chem Transuranium
 Elem 209–213. https://doi.org/10.1016/B978-0-08-020638-7.50043-9
- Nakayama S, Nelson D (1988) Comparison of distribution coefficients for americium
 and curium: Effects of pH and naturally occurring colloids. J Environ Radioact 8:173–
 181. https://doi.org/10.1016/0265-931X(88)90024-0
- Nash KL, Madic C, Mathur JN, Lacquement J (2010) Actinide Separation Science and Technology. In: The Chemistry of the Actinide and Transactinide Elements. Springer, Dordrecht, pp 2622–2798
- Puigdomenech I (2000) Windows software for the graphical presentation of chemical speciation. https://doi.org/10.13140/RG.2.2.33538.91842
- 9. Shannon RD (1976) Revised effective ionic radii and systematic studies of interatomic distances in halides and chalcogenides. Acta Crystallogr A 32:751–767.
 https://doi.org/10.1107/S0567739476001551
- Gharibyan N, Dailey A, McLain DR, et al (2014) Extraction Behavior of Americium and Curium on Selected Extraction Chromatography Resins from Pure Acidic Matrices.
 Solvent Extr Ion Exch 32:391–407. https://doi.org/10.1080/07366299.2014.884888
- 365 11. NF ISO 13167 January 2016
- Modolo G, Kluxen P, Geist A (2010) Demonstration of the LUCA process for the
 separation of americium(III) from curium(III), californium(III), and lanthanides(III) in
 acidic solution using a synergistic mixture of bis(chlorophenyl)dithiophosphinic acid and
 tris(2-ethylhexyl)phosphate. Radiochim Acta Int J Chem Asp Nucl Sci Technol 98:193–
 201. https://doi.org/10.1524/ract.2010.1708
- 371 13. NF T90-210 Novembre 2018
- 14. Davies CW (1938) The extent of dissociation of salts in water. Part VIII. An equation for
 the mean ionic activity coefficient of an electrolyte in water, and a revision of the

dissociation constants of some sulphates. J Chem Soc Resumed 2093–2098. 374 https://doi.org/10.1039/JR9380002093 375 15. Grivé M, García D, Campos I, et al (2013) Simple Functions Spreadsheet tool: 376 Phosphates update and Temperature assessment. Final Rep Amphos21 377 378 16. Moll H, Brendler V, Bernhard G (2011) Aqueous curium(III) phosphate species characterized by time-resolved laser-induced fluorescence spectroscopy. Radiochim 379 Acta Int J Chem Asp Nucl Sci Technol 99:775-782. 380 https://doi.org/10.1524/ract.2011.1878 381 17. Fanghänel Th, Kim JI, Paviet P, et al (1994) Thermodynamics of Radioactive Trace 382 Elements in Concentrated Electrolyte Solutions: Hydrolysis of Cm3+ in NaCl-Solutions. 383 Radiochim Acta 66-67:81-88. https://doi.org/10.1524/ract.1994.6667.special-issue.81 384 385 18. Skerencak A, Panak PJ, Hauser W, et al (2009) TRLFS study on the complexation of Cm(III) with nitrate in the temperature range from 5 to 200 °C. Radiochim Acta Int J 386 Chem Asp Nucl Sci Technol 97:385-393. https://doi.org/10.1524/ract.2009.1631 387 19. Bansal BML, Patil SK, Sharma HD (1964) Chloride, nitrate and sulphate complexes of 388 europium (III) and americium (III). J Inorg Nucl Chem 26:993-1000. 389 https://doi.org/10.1016/0022-1902(64)80261-X 390 20. Kitamura A, Doi R, Yoshida Y (2014) Update of JAEA-TDB. Update of thermodynamic 391 data for palladium and tin, refinement of thermodynamic data for protactinium, and 392 preparation of PHREEQC database for use of the Brønsted-Guggenheim-Scatchard 393 model. Japan Atomic Energy Agency 394 21. Burney GA, Porter JA (1967) Solubilities of Pu(III), Am(III), and Cm(III) oxalates. 395 Inorg Nucl Chem Lett 3:79-85. https://doi.org/10.1016/0020-1650(67)80128-4 396 397 22. Fanghänel Th, Weger HT, Schubert G, Kim JI (1998) Bicarbonate Complexes of Trivalent Actinides - Stable or Unstable? Radiochim Acta 82:55-58. 398 https://doi.org/10.1524/ract.1998.82.special-issue.55 399 23. Fouillac C, Criaud A (1984) Carbonate and bicarbonate trace metal complexes: critical 400 reevaluation of stability constants. Geochem J 18:297-303 401 24. Janda J, Sas D, Tokárová A (2015) Characterization of the extraction properties of the 402 selected extraction sorbents TRU, TEVA, UTEVA and DIPHONIX. J Radioanal Nucl 403 Chem 304:301-312. https://doi.org/10.1007/s10967-014-3538-7 404 25. Könnecke Th, Fanghänel Th, Kim JI (1997) Thermodynamics of Trivalent Actinides in 405 Concentrated Electrolyte Solutions: Modelling the Chloride Complexation of Cm(III). 406 Radiochim Acta 76:131-136. https://doi.org/10.1524/ract.1997.76.3.131 407 408 26. Khopkar PK, Mathur JN (1974) Thiocyanate complexing of some trivalent actinides and 409 lanthanides. J Inorg Nucl Chem 36:3819-3825. https://doi.org/10.1016/0022-1902(74)80171-5 410 27. NAVRATIL JD (1989) Ion Exchange Technology in Spent Fuel Reprocessing. J Nucl 411 Sci Technol 26:735-743. https://doi.org/10.1080/18811248.1989.9734377 412