

Can real-time polymerase chain reaction allow a faster recovery of hospital activity in cases of an incidental discovery of carbapenemase-producing Enterobacteriaceae and vancomycin-resistant Enterococci carriers?

R. Saliba, C. Neulier, D. Seytre, A. Fiacre, F. Faibis, P. Leduc, M. Amara, F. Jauréguy, E. Carbonnelle, J.-R. Zahar, et al.

▶ To cite this version:

R. Saliba, C. Neulier, D. Seytre, A. Fiacre, F. Faibis, et al.. Can real-time polymerase chain reaction allow a faster recovery of hospital activity in cases of an incidental discovery of carbapenemase-producing Enterobacteriaceae and vancomycin-resistant Enterococci carriers?. Journal of Hospital Infection, 2019, 103, pp.115 - 120. 10.1016/j.jhin.2019.07.004. hal-03487531

HAL Id: hal-03487531

https://hal.science/hal-03487531

Submitted on 20 Jul 2022

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



- 1 Can Real-Time PCR Allow a Faster Recovery of Hospital Activity in Cases of an
- 2 Incidental Discovery of Carbapenemase-Producing Enterobacteriaceae and
- 3 Vancomycin-Resistant Enterococci Carriers?
- 4 Rindala Saliba^{1,2}, Caroline Neulier³, Delphine Seytre², Alain Fiacre⁴, Frédéric Faibis⁴,
- 5 Pierre Leduc³, Marlène Amara⁵, Françoise Jauréguy^{1,2}, Etienne Carbonnelle^{1,2}, Jean-
- 6 Ralph Zahar^{1,2*}, Laurence Marty⁶,
- 7 ¹ IAME, UMR 1137, Université Paris 13, Sorbonne Paris Cité, France
- 8 ² Service de Microbiologie Clinique et Unité de Contrôle et de Prévention du risque
- 9 Infectieux, Groupe Hospitalier Paris Seine Saint-Denis, AP-HP, 125 rue de Stalingrad
- 10 93000, Bobigny, France
- ³ Service de Prévention du Risque Infectieux, CH André Mignot, Versailles, France
- ⁴ Département de Microbiologie, Grand Hôpital de l'Est Francilien site Marne la Vallée,
- 13 France

- ⁵ Département de Microbiologie et Hygiène, CH André Mignot, Versailles, France
- ⁶ Unité d'Hygiène Inter Hospitalière Nord Seine et Marne Grand Hôpital de l'Est Francilien
- site Marne la Vallée, France
- * Corresponding author:
- 19 Jean-Ralph Zahar
- 20 125 rue de Stalingrad 93000, Bobigny, France
- 21 +33 6 03 07 90 88
- 22 jrzahar@gmail.com
- 24 Running title: RT-PCR for Infection Control Management

Summary
Background. Detection of faecal carriers of carbapenemase-producing Enterobacteriaceae
(CPE) and vancomycin-resistant enterococci (VRE) has become a routine medical practice
many countries. In an outbreak setting, several public health organizations recommend
threenweekly rectal screenings to rule-out acquisition in contact patients. This strategy,
associated with bed closures and reduction of medical activity for a relatively long time,
seems costly.
Aim. The objective of this study was to test the positive and negative predictive values of I
PCR(GeneXpert®) carried-out at Day-0, compared to conventional three weekly culture-
based rectal screenings, in identifying, among contact patients, those who acquired
CPE/VRE.
<i>Methods.</i> We conducted a multicentre retrospective study from January2015 to October20
We included all contact patients (CP) identified from index patients (IP) colonized or infection
with CPE/VRE, incidentally discovered. Each CP was investigated at Day-0 by PCR
(GeneXpert®), and by the recommended three weekly screenings.
Findings. Twenty-two IP and 159 CP were included. An average of 0.77 secondary case p
patient was noted, with a mean duration of contact of 10 days [range 1-64]. Among the 15
CP, 16 (10%) had a CPE/VRE-positive culture during the monitoring period. Rectal
screenings were positive at Day-0 (10 patients), Day-7 (2 patients), Day-14 (4 patients).
Thirteen out of 16 patients with positive culture had at Day-0 a positive PCR. Overall, a
concordance of 97.5% (155/159) was observed between the three-weekly screenings and
Day-0 PCR results. When performed on CP at Day-0 of the identification of an IP,
PCR(GeneXpert®) allowed to shorten turnaround time by 5 to 27days, compared to three
weekly screenings. PPV and NPV were respectively 100% and 98%.
Conclusions. RT-PCR (GeneXpert®) can avoid three weekly rectal samplings needed to
rule-out acquisition of CPE/VRE.

Keywords. Real-time PCR, Carbapenemase-Producing Enterobacteriaceae, Vancomycin-

Resistant Enterococci, infection control, outbreak, rectal screening.

53

55	Text:
56	
57	
58	Introduction:
59	Carbapenemase-producing Enterobacteriaceae (CPE) and vancomycin-resistant enterococci
60	(VRE) emerged globally over the last decade, resulting in antibiotic-resistant infections that
61	are difficult to treat and costly to healthcare facilities. ¹ The search and isolate policy is
62	considered the main infection prevention and control strategy. Rapid, accurate detection of
63	intestinal CPE and VRE colonization is critical to minimize transmission, and hence reduce
64	costly, difficult to treat infections caused by these emerging extensively drug-resistant (XDR)
65	bacteria. ²
66	The 2013 French National Guidelines recommend implementation of a strict isolation policy,
67	including cohorting patients with CPE in a dedicated ward with dedicated healthcare workers
68	(HCWs) and an extensive screening policy of contact patients. In case of identifying a non-
69	cohorted index patient, recommendations are to close the ward concerned and apply a
70	screening policy to all contact patients. Weekly rectal swabbing is recommended to
71	investigate possible acquisition among contact patients. Contact patients are considered non-
72	carriers if three consecutive weekly negative rectal swabs are obtained; they are then no
73	longer monitored. This strategy seems to be costly as it is associated with a reduction of
74	medical activity. Moreover, it could expose patients to higher risk of complications as they
75	will usually remain isolated until screening is complete. ³
76	The current process for detecting CPE in clinical samples is usually performed using
77	selective cultures and may require up to 72h for obtaining the results. ⁴
78	Different molecular assays have been developed and have shown their effectiveness in
79	identifying the resistance genes from bacterial colonies. Real-time PCR assays have been
80	reported to have lower limits of detection than conventional agar-based culture
81	methodologies. 4-7 When used directly on rectal samples, they offer laboratories the ability to
82	reduce turnaround times. However, there is no data on their effectiveness in this specific
83	situation. Moreover, no previous study has demonstrated if they have a potential added value
84	for infection control management in an epidemiological setting. Therefore, evaluating the

- 85 clinical utility of PCR-based assays as a screening tool, is needed.
- 86 The main objective of this study was to test in a multicentre retrospective survey the
- 87 effectiveness of real time PCR (GeneXpert®) carried out at Day 0 in identifying among
- 88 contact patients those who acquired XDR. An additional objective was to define the negative
- 89 predictive value (NPV) of this technique; a high NPV would support using this strategy in
- 90 place of the standard three week sampling and testing regimen.
- 91 Part of the results was presented at the 29th ECCMID Congress.

Materials/methods:

93

94 We conducted a multicentre retrospective study in three large French hospitals (Avicenne Hospital AP-HP, Mignot Hospital-Hospital Center of Versailles, Grand Hospital Eastern IIe 95 Site de Marne-la-Vallée) from January 2015 to October 2018. In these hospitals, the 2013 96 French National Guidelines recommendations ⁸ for the prevention of spread of XDR are 97 systematically applied. Thus, in case of identifying an index patient, patients are cohorted and 98 an extensive screening policy for all isolated contact patients is applied. Contact patients are 99 only considered non-carriers, and isolation discontinued, if three consecutive weekly rectal 100 swabs screenings are negative by recommended culture techniques. 101 An index patient was defined as a non-cohorted patient colonized or infected with CPE or 102 VRE incidentally discovered post hospitalization. We defined contact patients as all the 103 104 adults (> 18 years old) staying in the same unit during the same period as an index patient, sharing the same HCWs and for whom, at least three-weekly rectal screenings were 105 106 performed during the same hospitalization and after the end of their exposure using 107 conventional culture techniques, plus a real-time PCR (GeneXpert®) assay done on the first 108 screening sample. A secondary case patient was defined as a contact patient for whom CPE or VRE was later identified by the weekly rectal screenings. A non-colonized patient was 109 defined as a contact patient who remained CPE- or VRE-negative during his stay in the 110 hospital. Exposure time was defined as the duration of contact between the index patient and 111 a secondary case patient before being identified as CPE or VRE positive or a control patient 112 throughout his stay in the same ward. ⁹ 113 Data were collected retrospectively in case and control patients by reviewing medical and 114 microbiological records. Collected data included demographic characteristics, including age, 115 gender, ward of contact (medical, surgical, or intensive care unit), exposure time, duration of 116 hospitalization in the ward, and antibiotic administration during the exposure period. 117 CPE and VRE detections were performed using respectively Xpert® Carba-R and Xpert® 118 vanA/vanB (Cepheid, GeneXpert) following manufacturer recommendations, on the first 119 screening sample. CPE detection by culture consisted of direct plating of the rectal swab on 120 ChromID® CARBA SMART medium (BioMérieux, Marcy-l'Étoile, France) for 24h, 121 likewise, screening of VRE was carried out on ChromID® VRE medium (BioMérieux, 122 Marcy-l'Étoile, France). Antibiotic susceptibility testing according to Comité de 123

124	l'antibiogramme de la Société Française de Microbiologie - European Committee on
125	Antimicrobial Susceptibility Testing (CA-SFM - EUCAST) guidelines 10 was performed on
126	each type of suspected colony grown on selective media. Confirmation of resistance was
127	accomplished using Carba NP test and $\beta\text{-}Carba$ test (Bio-Rad, Marnes-la-Coquette, France)
128	^{11,12} and Xpert® Carba-R on suspected CPE isolates and Xpert® vanA/vanB on suspected
129	VRE isolates.
130	Statistical method:
131	The sensitivity, specificity, positive predictive value (PPV) and negative predictive value
132	(NPV) were calculated for all tested isolates. Statistical analysis was performed using 95%
133	confidence intervals (CIs) calculated using Student's t-test.

Results:

134

During the study period 22 non-cohorted patients, hospitalized in different wards, carrying or 135 infected with CPE (59%) or VRE (41%) were incidentally discovered post hospitalization. 136 The mean age of these index patients was 68 years-old and their average length of hospital 137 stay, 12 days (range between 1 and 50 days). 16 index patients over 22 (72.7%) were exposed 138 to antibiotics during their hospital stay. CPE accounted for 13 (59%) cases, with OXA-48 139 140 being the most prevalent mechanism of resistance (Table I). Among the 159 contact patients (as defined previously) only 17 (10.7%) acquired an XDR. In 141 fact, 8 index patients of 22 generated at least one secondary case patient, the average being 142 0.77 secondary case per patient (range between 0 and 3). The mean duration of contact in the 143 whole population was 10 days (range between 1 and 64 days). Among the 159 patients 96 144 (60.4%) patients were receiving antibiotic therapy during the contact period (Table II). 145 According to French recommendations, all contact patients were followed-up weekly and 146 were screened by rectal swabbing. Furthermore, real-time PCR (GeneXpert®) carried out on 147 rectal swabs of the 159 contact patients at Day 0 and 14 patients were detected as CPE 148 149 carriers. Conventional agar-based methodology used for weekly rectal screenings, was positive in 16 patients. Ten had a positive culture result at first day of screening (Day-0), 2 at 150 the second screening week and 4 at the third screening week. In one case, PCR was positive 151 152 at Day 0, whereas, conventional culture remained negative for the three-weekly rectal screenings. In 3 other cases, all related to VRE acquisition, PCR was negative at Day 0, while 153 154 cultures were positive at the 2nd or 3rd round screening. Overall, a concordance of 97.5% (155/159) was observed between the three weekly rectal 155 screens using culture techniques and Day 0 real-time PCR (GeneXpert®) results. When 156 performed on all contact patients at Day 0 after identifying of an index patient, real-time PCR 157 (GeneXpert®) allowed to shorten turnaround time by 5 to 27 days (6 days in average) 158 compared to the conventional three weekly culture-based rectal screening method. Among 159 160 the 4 discordances, one case had a positive PCR at Day 0 for CPE OXA-48, whereas, conventional culture methodology stayed negative for the threenweekly rectal screens. In the 161 162 3 remaining discordant cases, all VRE acquisition, PCR was negative at Day 0, while cultures were positive during the 2nd screening in one case, and the 3rd screening in the two other 163 cases (Table III). 164

We did not find any statistically significant differences in age, duration of contact and antibiotic consumption between secondary case patients diagnosed during the 1^{st} screening and those found negative during the 1^{st} screening but identified as carriers later at the 2^{nd} or 3^{rd} screening.

Discussion and Conclusion:

170

In our multicentre retrospective study a concordance of 97.5% (155/159) was observed 171 between the three-weekly rectal screenings using culture techniques and Day 0 real-time PCR 172 (GeneXpert®) results. All, secondary cases among contact patients were considered having 173 174 acquired the XDR bacteria within the first day of contact with the index patient and before the identification of carriage. Indeed, after identification, all case index were systematically 175 cohorted according to the French recommendations and contact patients were not in contact 176 with any other index patient or known carrier during their hospital stay. Among the 4 177 178 discordances, one case had a positive PCR at Day 0, whereas, conventional culture methodology stayed negative for the three-weekly rectal screenings. This may be explained 179 180 by the fact that the patient had been receiving antibiotics, affecting the expression of the resistance genes in culture, or, since it was an OXA-48 gene detected, the possibility of a 181 transiently carriage of *Shewanella spp.* which naturally produce OXA-48 variants. ⁵ It needs 182 to be considered also that molecular methods may be more sensitive than the culture method 183 that was the reference standard in our study. Knowing that, heavily colonized patients are 184 more likely to yield opportunities for horizontal transmission, it raises the issue of whether 185 the excellent sensitivity of PCR and detection of very low levels of colonization is 186 significant, regarding the consequent costly infection control measures. ² 187 In the 3 remaining discordant cases, all VRE acquisition, PCR was negative at Day 0, while 188 cultures were positive during the 2nd screening in one case, and the 3rd screening in the two 189 other cases. The absence of detection in these 3 cases at Day 0 may be associated with low 190 undetectable bacterial load in the rectal swab specimens, becoming detectable at the 2nd or 3rd 191 screenings or even to sample techniques. ¹³ However, looking closer to those cases, all 3 are 192 193 parts of complicated outbreaks generating several secondary cases. Therefore, we can also suggest that those 3 patients were "true-negative" at Day 0 and have acquired the VRE later 194 195 during their hospital stay. Detection of faecal carriers of XDR bacteria has become a routine clinical practice in many 196 197 parts of the world and is recommended by public health organizations for the containment of 198 their spread. Among different approaches used, more recently, molecular approaches have been developed to increase detection sensitivity and decrease reporting time. 14 199 When performed on all contact patients at Day 0 of the discovery of an index patient, real-200 time PCR (GeneXpert®) allowed to shorten turnaround time by 5 to 27 days (6 days in 201

202 average) compared to the conventional three weekly culture-based rectal screenings method. The result of a diagnostic assay has a clinical implication, and as a screening tool for XDR 203 bacteria, that implication is the implementation of infection control measures. The rapid 204 cohorting of carriers is crucial to prevent outbreaks and reduce the social impact and the cost 205 206 of infection control measures empirically implemented for high risk patients, who end up being negative for XDR bacteria colonization. ⁵ Moreover, the knowledge of the PPV and 207 NPV is essential for decision making by the infection control team in respect to the value of a 208 positive and a negative test result. Some authors have studied the performance of real time 209 210 PCR (GeneXpert®), but most have studied its effectiveness in detecting resistance mechanisms on bacterial strains. ^{15–19} In addition, the limited number of authors who have 211 been interested in evaluating the performance of GeneXpert® directly on rectal swabs, 212 compared to the routine culture techniques, have studied its effectiveness only at the 213 admission of patients considered as 'high-risk' for CPE carriage. ^{2,4–7,14,20–22} Hence, there are 214 few data on the effectiveness of real-time PCR in this specific situation and no previous study 215 216 has demonstrated if they have a potential added value for infection control management in an epidemiological setting. To our best knowledge, this is the first report evaluating the clinical 217 utility of PCR-based assays as a screening tool. Our study suggested that real time PCR 218 (GeneXpert®), carried out at Day 0 is efficient in identifying among contact patients those 219 who acquired XDR, with an NPV of 98%. 220 According to the 2013 French National Guidelines recommendations ⁸ and the international 221 recommendations, ²³ in cases of identifying an index patient, an extensive screening policy 222 for all isolated contact patients must be applied. Contact patients are only considered non-223 224 carriers and no more isolated, if three weekly rectal swabs screenings came back with a negative result, by recommended culture techniques. This strategy seems to be costly as it is 225 226 associated with bed closures and reduction of medical activity, for a relatively long duration of time. In addition to that, it could expose the isolated patients to a higher risk of 227 complications as they are possibly receiving less optimal management for their medical 228 condition, compared to non-isolated patients with the same medical condition. 229 230 The major strength of our study is the finding that, in an outbreak context, if the screening real-time PCR (GeneXpert®) test for a patient at Day 0 is negative, the risk that it will be 231 discovered later colonized remains weak. Thereby, contact patient's isolation will only be 232

limited to those who have their screening real-time PCR positive, means that they have acquired the XDR bacteria. We acknowledge that this study suffers from several limitations. First, its retrospective nature and the lack of retrospective clinical data allowing for better targeting patients at risk of acquisition of XDR bacteria, knowing that Hilliquin et al demonstrated recently in their large multicentre case-control study, that at least three risk factors are associated with KP OXA- 48 acquisition in contact patients of a non-cohorted index patient: (i) geographical proximity (ii) antibiotic therapy during exposure time and (iii) at least one invasive procedure. ⁹ Second, the sample size with the limited number of secondary cases. Third, we could not exclude that the negative results on culture and/or on PCR were possibly the consequence of a poor quality of sampling (no macroscopic staining was visualized), information not available retrospectively. In conclusion, PCR (GeneXpert®) can provide valuable information for infection control programs, and avoid three-weekly rectal samplings needed to rule out acquisition of CPE or VRE. When included in a clinical algorithm, it can reduce the cost behind the infection control measures implementation. Further studies need to be conducted comparing different real-time PCR techniques to the three-weekly rectal screenings using culture techniques, on a larger number of patients, taking into account different risk factors of XDR acquisition and transmission.

233

234

235

236

237

238

239

240

241

242

243

244

245

246

247

248

249

250

252	Funding
253	
254	
255	This research did not receive any specific grant from funding agencies in the public,
256	commercial, or not-for-profit sectors.
257	
258	Conflicts of Interest
259	
260	
261	None to report.

2	62	References	;

263		
203		

- 1. Wilson HJ, Khokhar F, Enoch DA, Brown NM, Ahluwalia J, Dougan G, et al. Point-prevalence survey of carbapenemase-producing Enterobacteriaceae and vancomycin-resistant enterococci in adult inpatients in a university teaching hospital in the UK. J Hosp Infect. 2018 Sep;100(1):35–9.
- Lowman W, Marais M, Ahmed K, Marcus L. Routine active surveillance for carbapenemase-producing Enterobacteriaceae from rectal swabs: diagnostic implications of multiplex polymerase chain reaction. J Hosp Infect. 2014 Oct;88(2):66–71.
- Otter JA, Burgess P, Davies F, Mookerjee S, Singleton J, Gilchrist M, et al. Counting the cost of an outbreak of carbapenemase-producing Enterobacteriaceae: an economic evaluation from a hospital perspective. Clin Microbiol Infect Off Publ Eur Soc Clin Microbiol Infect Dis. 2017 Mar;23(3):188–96.
- 4. Huang T-D, Bogaerts P, Ghilani E, Heinrichs A, Gavage P, Roisin S, et al. Multicentre evaluation of the Check-Direct CPE® assay for direct screening of carbapenemase-producing Enterobacteriaceae from rectal swabs. J Antimicrob Chemother.
 2015;70(6):1669–73.
- 5. Hoyos-Mallecot Y, Ouzani S, Dortet L, Fortineau N, Naas T. Performance of the Xpert® Carba-R v2 in the daily workflow of a hygiene unit in a country with a low prevalence of carbapenemase-producing Enterobacteriaceae. Int J Antimicrob Agents. 2017 Jun;49(6):774–7.
- Tenover FC, Canton R, Kop J, Chan R, Ryan J, Weir F, et al. Detection of colonization
 by carbapenemase-producing Gram-negative Bacilli in patients by use of the Xpert
 MDRO assay. J Clin Microbiol. 2013 Nov;51(11):3780–7.
- 7. Otter JA, Dyakova E, Bisnauthsing KN, Querol-Rubiera A, Patel A, Ahanonu C, et al. Universal hospital admission screening for carbapenemase-producing organisms in a low-prevalence setting. J Antimicrob Chemother. 2016;71(12):3556–61.
- 8. Lepelletier D, Berthelot P, Lucet J-C, Fournier S, Jarlier V, Grandbastien B, et al. French recommendations for the prevention of "emerging extensively drug-resistant bacteria" (eXDR) cross-transmission. J Hosp Infect. 2015 Jul;90(3):186–95.
- Hilliquin D, Le Guern R, Thepot Seegers V, Neulier C, Lomont A, Marie V, et al. Risk factors for acquisition of OXA-48-producing Klebsiella pneumonia among contact patients: a multicentre study. J Hosp Infect. 2018 Mar;98(3):253–9.
- 10. European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for
 interpretation of MICs and zone diameters. Version 9.0, valid from 2019-01-01
 [Internet]. Available from: https://www.sfm-microbiologie.org/casfm-eucast

- 300 11. Dortet L, Cuzon G, Plésiat P, Naas T. Prospective evaluation of an algorithm for the 301 phenotypic screening of carbapenemase-producing Enterobacteriaceae. J Antimicrob 302 Chemother. 2016 Jan;71(1):135–40.
- 12. Compain F, Gallah S, Eckert C, Arlet G, Ramahefasolo A, Decré D, et al. Assessment of
 Carbapenem Resistance in Enterobacteriaceae with the Rapid and Easy-to-Use
 Chromogenic β Carba Test. J Clin Microbiol. 2016;54(12):3065–8.
- 13. Dyakova E, Bisnauthsing KN, Querol-Rubiera A, Patel A, Ahanonu C, Tosas Auguet O,
 et al. Efficacy and acceptability of rectal and perineal sampling for identifying
 gastrointestinal colonization with extended spectrum β-lactamase Enterobacteriaceae.
 Clin Microbiol Infect Off Publ Eur Soc Clin Microbiol Infect Dis. 2017
 Aug;23(8):577.e1-577.e3.
- 14. Tato M, Ruiz-Garbajosa P, Traczewski M, Dodgson A, McEwan A, Humphries R, et al.
 Multisite Evaluation of Cepheid Xpert Carba-R Assay for Detection of Carbapenemase Producing Organisms in Rectal Swabs. J Clin Microbiol. 2016;54(7):1814–9.
- 15. Nijhuis R, Samuelsen O, Savelkoul P, van Zwet A. Evaluation of a new real-time PCR
 assay (Check-Direct CPE) for rapid detection of KPC, OXA-48, VIM, and NDM
 carbapenemases using spiked rectal swabs. Diagn Microbiol Infect Dis. 2013
 Dec;77(4):316–20.
- 16. Traczewski MM, Carretto E, Canton R, Moore NM, Carba-R Study Team. Multicenter
 Evaluation of the Xpert Carba-R Assay for Detection of Carbapenemase Genes in
 Gram-Negative Isolates. J Clin Microbiol. 2018 Aug;56(8).
- 17. Smith M, Diederen B, Scharringa J, Leversteijn-van Hall M, Fluit AC, Cohen Stuart J.
 Rapid and accurate detection of carbapenemase genes in Enterobacteriaceae with the
 Cepheid Xpert Carba-R assay. J Med Microbiol. 2016 Sep;65(9):951–3.
- 18. Tanner H. Verification of the Cepheid Xpert Carba-R assay for the detection of
 carbapenemase genes in bacterial isolates cultured on alternative solid culture media. J
 Hosp Infect. 2017 Nov;97(3):254–7.
- 19. Findlay J, Hopkins KL, Meunier D, Woodford N. Evaluation of three commercial assays
 for rapid detection of genes encoding clinically relevant carbapenemases in cultured
 bacteria. J Antimicrob Chemother. 2015 May;70(5):1338–42.
- 20. Moore NM, Cantón R, Carretto E, Peterson LR, Sautter RL, Traczewski MM, et al. Rapid
 Identification of Five Classes of Carbapenem Resistance Genes Directly from Rectal
 Swabs by Use of the Xpert Carba-R Assay. J Clin Microbiol. 2017;55(7):2268–75.
- 21. Kim D-K, Kim HS, Pinto N, Jeon J, D'Souza R, Kim MS, et al. Xpert CARBA-R Assay
 for the Detection of Carbapenemase-Producing Organisms in Intensive Care Unit
 Patients of a Korean Tertiary Care Hospital. Ann Lab Med. 2016 Mar;36(2):162–5.
- 22. Rezende TFT, Doi AM, Quiles MG, Pignatari ACC, Manfrendi S, Grothe C, et al.
 Detection of colonization by carbapenem-resistant organisms by real-time polymerase chain reaction from rectal swabs in patients with chronic renal disease. J Hosp Infect.
 2017 Jun;96(2):123–8.

23. Schwaber MJ, Carmeli Y. An ongoing national intervention to contain the spread of
 carbapenem-resistant enterobacteriaceae. Clin Infect Dis Off Publ Infect Dis Soc Am.
 2014 Mar;58(5):697–703.

Table I. Description of the characteristics of the index patients

	Total Index Patients		
	N=22 (%)		
Female sex	11 (50)		
Age, mean	68		
Hospitalization duration (days), mean	12		
Ward of admission			
Medical	15 (68.2)		
Surgical	4 (18.2)		
Intensive Care Unit	3 (13.6)		
Antibiotic therapy during hospitalization	16 (72.7)		
Secondary cases generated	8 (36.4)		
Mechanisms of resistance			
VRE ¹ (vanA)	9 (40.9)		
CPE ² (OXA 48)	11 (50)		
CPE ² (NDM)	1 (4.5)		
CPE ² (KPC)	1 (4.5)		

¹ VRE: Vancomycin-resistant enterococci

350 ² CPE: Carbapenemase-producing Enterobacteriaceae

353 Table II. Description of the characteristics of the contact patients

	Total Contact Patients N=159 (%)	Secondary Cases N=17 (%)	Negative Cases N=142 (%)	p
Female sex	78 (49)	10 (58.8)	68 (47.9)	0.44
Age, median CI [min25%-75%,max]	70	80 [61-85]	72 [60,2-83,7]	0.54
Duration of contact with IP				
(days), median CI [25%-	10.4 [1,64]	5 [3,7]	10 [4,15]	0.03
75%]				
Ward of admission				
Medical	120 (75.5)	12 (70.6)	108 (76)	0.56
Surgical	27 (17)	5 (29.4)	22 (15.5)	0.17
Intensive Care Unit	12 (7.5)	0 (0)	12 (8.5)	0.36
Antibiotic therapy during hospitalization	96 ³ (60.4)	11 (64.7)	85 (59.9)	0.79
Mechanisms of resistance	17 (10.7)			
VRE ¹ (vanA)	10 (6.3)	10 (58.8)	0	
CPE ² (OXA 48)	6 (3.8)	6 (35.3)	0	
CPE ² (NDM)	1 (0.6)	1 (5.9)	0	

^{355 &}lt;sup>1</sup> VRE: Vancomycin-Resistant-Enterococci

^{356 &}lt;sup>2</sup> CPE: Carbapenemase-Producing-Enterobacteriaceae

³ Penicillin + β-lactamase inhibitor 26; cephalosporins 32; carbapenems 5; fluoroquinolones

^{358 13;} aminoglycosides 7; metronidazole 13; other antibiotics 27

Table III. Description of the characteristics of discordant cases between the three weekly rectal screening using culture techniques and Day 0 real-time PCR (GeneXpert®)

	Case 1	Case 2	Case 3	Case 4
Day 0 real-time PCR	Positive	Negative	Negative	Negative
1st week screening	Negative	Negative	Negative	Negative
2 nd week screening	Negative	Positive	Negative	Negative
3 rd week screening	Negative	Negative	Positive	Positive
Sex	M	M	M	F
Age	55	81	69	91
Duration of contact with index patient (days)	2	4	14	7
Ward of admission	Surgical	Medical	Medical	Medical
Antibiotic therapy during hospitalization	Yes	No	No	Yes
Mechanisms of resistance	CPE ¹ (OXA 48)	VRE ² (vanA)	VRE ² (vanA)	VRE ² (vanA)

362 ¹ CPE: Carbapenemase-producing Enterobacteriaceae

363 ² VRE: Vancomycin-resistant enterococci

359

360