



**HAL**  
open science

## Automatic quantification of uracil and dihydrouracil in plasma

Tiphaine Robin, Franck Saint-Marcoux, Doriane Toinon, Naïma Tafzi, Pierre Marquet, Souleiman El Balkhi

► **To cite this version:**

Tiphaine Robin, Franck Saint-Marcoux, Doriane Toinon, Naïma Tafzi, Pierre Marquet, et al.. Automatic quantification of uracil and dihydrouracil in plasma. *Journal of Chromatography B - Analytical Technologies in the Biomedical and Life Sciences*, 2020, 1142, pp.122038 -. 10.1016/j.jchromb.2020.122038 . hal-03489995

**HAL Id: hal-03489995**

**<https://hal.science/hal-03489995>**

Submitted on 22 Aug 2022

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Distributed under a Creative Commons Attribution - NonCommercial 4.0 International License

## **Automatic quantification of uracil and dihydrouracil in plasma.**

Tiphaine Robin<sup>1</sup>; Franck Saint-Marcoux<sup>1,2</sup>; Doriane Toinon<sup>3</sup>; Naïma Tafzi<sup>1</sup>; Pierre Marquet<sup>1,2</sup>; Souleiman El Balkhi<sup>1,2</sup>

- (1) Department of pharmacology, toxicology and pharmacovigilance, Limoges University Hospital, France
- (2) INSERM UMR 1248
- (3) Shimadzu corporation, Kyoto, Japan

**Corresponding author :**

Franck Saint-Marcoux  
Department of pharmacology, toxicology and pharmacovigilance  
Limoges University Hospital  
France

Tel: +33555056143  
Fax: +33555056162  
[franck.saint-marcoux@unilim.fr](mailto:franck.saint-marcoux@unilim.fr)

**Abstract:**

Fluoropyrimidines-based chemotherapies are the backbone in the treatment of many cancers. However, the use of 5-fluorouracil and its oral pre-prodrug, capecitabine, is associated with an important risk of toxicity. This toxicity is mainly due to a deficiency of dihydropyrimidine dehydrogenase (DPD). This deficiency may be detected by using a phenotypic approach that consists in the measurement of uracilemia or the calculation of dihydrouracil (UH<sub>2</sub>)/uracil (U) ratio. For uracilemia, a threshold value of 16 ng/ml has been proposed for partial deficiency, while a value of 150 ng/ml has been proposed for complete deficiency

We have developed a rapid, accurate and fully-automated procedure for the quantification of U and UH<sub>2</sub> in plasma. Sample extraction was carried out by a programmable liquid handler directly coupled to a liquid chromatography - tandem mass spectrometry (LC-MS/MS) system. The method was validated according to the EMA guidelines and ISO 15189 requirements and was applied to real patient samples (n=64).

The limit of quantification was 5 and 10 ng/ml for U and UH<sub>2</sub> respectively. Imprecision and inaccuracy were less than 15% for inter and intra-assay tests. Comparison with dedicated routine method showed excellent correlation.

An automated procedure perfectly fulfills the need of low inaccuracy and CVs at the threshold values (less than 5% at 16 ng/ml) and is highly suitable for the characterization of DPD deficiency.

Automatization should guaranty reliable and robust performances by minimizing the sources of variation such as volume inaccuracies, filtration or manual extraction related errors.

**Keyword** : LC-MS/MS, automated sample preparation, 5-FU, U, UH<sub>2</sub>, DPD deficiency

## 1. Introduction

Fluoropyrimidines-based chemotherapies remain important antineoplastic agents since first approval for testing in humans about 60 years ago [1, 2]. Indeed, 5-fluorouracil and its oral prodrug capecitabine are the backbone in the treatment of colorectal, pancreatic, gastric, breast, head and neck cancers [3, 4]. However, the use of fluoropyrimidines is associated with an important risk of toxicity, including a painful skin condition known as hand-foot syndrome, diarrhea, mucosal ulceration and immune suppression [5-7]. As a result, this toxicity limits or delays the administration of optimal or successive courses, which impacts the clinical outcome of patients with cancer. The mechanisms of fluoropyrimidines-based chemotherapies toxicity and the different strategies to prevent this toxicity are of great importance and numerous works have been published [5-10]. The dihydropyrimidine dehydrogenase (DPD), a rate-limited enzyme, is responsible of 80% of the catabolism of fluoropyrimidines. Most of the toxicities under fluoropyrimidines treatment can be explained by an inherited deficiency in DPD activity. Up to 15% of patients could exhibit a partial deficit exposing them to moderate or even severe adverse effects. However, in 0.1 to 0.5% of patient, the complete DPD deficiency could lead to lethal toxicity. This deficiency may be detected by using genotyping that detect the DPYD gene inactivating polymorphisms or by using a phenotyping by measuring uracilemia. The calculation of the dihydrouracil (UH<sub>2</sub>)/uracil (U) ratio to determine the phenotypes has been also used but it has raised some controversies.

Currently, there is no clear European regulatory obligation for DPD deficiency screening in patients given a fluoropyrimidine-based chemotherapy. However, based on the levels of evidence from the literature, the Group of Clinical Pharmacology in Oncology (GPCO)-UNICANCER and the French Network of Pharmacogenetics (RNPGx) has recently established the following recommendations [5] : (i) to explore DPD deficiency before initiating any chemotherapy containing fluoropyrimidines-based chemotherapies by measuring plasma U concentrations (possibly associated with UH<sub>2</sub>/U ratio); (ii) to genotype DPYD by exploring different variants; (iii) to reduce the initial dose, when partial deficiency is diagnosed or to proscribe fluoropyrimidines-based chemotherapies when the deficiency is complete. For uracilemia, a threshold value of 16 ng/ml has been proposed for partial deficiency, while a value of 150 ng/ml has been proposed for complete deficiency [11-13]. No consensus has been defined for the interpretation of UH<sub>2</sub>/U ratio, due to a great interlaboratory variability in the measurements [13-16].

In the past years, many labs have developed methods with UV detection [17-26] and more recently some have reported the use of mass spectrometry (MS) detection [26-32]. All these methods were based on classical manual sample preparations with protein precipitation, solid-phase or liquid-liquid extractions. With regards to the rapidity needed to report results and the increase in the number of patients to be monitored, the time dedicated to sample extraction may be a limiting step and its automation will offer a great added value. To our knowledge, no such automated extraction procedure coupled to LC-MS has been proposed yet for the measurement of U and UH<sub>2</sub>. Most recent extraction procedures involving any automated step of sample preparation (for the measurements of drugs or toxic compounds) were extensively reviewed by Vogeser and Kirchhoff [33], Zheng et al [34], Prabhu and Urban [35] or Wei [36]. Their main drawback is that they were not fully automated as the extracts needed at least to be transferred into the autosampler. In this context, we recently published a very high throughput method for the measurement of illicit drugs where extraction procedure was fully synchronized to LC-MS and where no human intervention was necessary after the primary tube was loaded on the system [37]. The main objective of this study was to develop a LC-MS method for the measurement of U and UH<sub>2</sub> in plasma with a fully-automated sample preparation.

## **2. Material and methods**

### *2.1 Chemicals and reagents:*

Dihydrouracil, dihydrouracil-<sup>13</sup>C<sup>15</sup>N<sub>2</sub>, uracil and uracil-<sup>13</sup>C<sup>15</sup>N<sub>2</sub> were supplied by LGC Standards (Molsheim, France). These analytes were stable in stock solution for 6 months at -20°C [29]. Acetonitrile and isopropanol were purchased from Carlo Erba Réactifs (Val de Reuil, France). Bovin serum albumin was purchased from Sigma Aldrich (ref. A7030) (Saint-Quentin Fallavier, France). Formic acid and acetic acid were supplied from Fisher Chemical (Illkirch, France). Pure water was obtained using a Millipore Direct Q purification system (Saint-Quentin-en-Yvelines, France).

### *2.2 Sample preparation:*

Extraction was performed by an automated sample preparation system, the CLAM-2030 (Shimadzu Corporation, Marne-la-Vallée, France), coupled to a LC-MS/MS system (see below). Briefly, once the primary (or secondary) tube was loaded onto the automated system, no further human intervention was required.

Blood samples were collected in heparinized tubes and plasma was separated within 1h30 by centrifugation at +4°C and stored at -20°C according to the recommendations [6, 38].

Once on board, 20 µL plasma (store at +4°C) were automatically pipetted in a pre-conditioned tube containing a filter, to which reagents were then added, mixed and filtered. Precisely: 40 µL of acetonitrile/acid formic (1%) was added to a PTFE (polytetrafluoroethylene) filter vial (0.45 µm pore size) previously conditioned with 20 µL of isopropanol. Then, the 20 µL of plasma and 10 µL of isotopically labelled internal standards (at 1500 ng/ml: namely, 167 ng/mL in the tube) were added. The mixture was agitated for 120 seconds (1900 rpm) then filtered by application of vacuum pressure (-60 to -65 kPa) for 120 seconds into a collection vial. Finally, the extract was diluted with 550µL of water and 50 µL were injected in the LC-MS/MS system.

### *2.3 LC-MS/MS conditions:*

The chromatographic system consisted of two Shimadzu LC-30 AD pumps (Nexera X2), a CTO 20AC oven and a SIL-30 AC-MP autosampler (Shimadzu Corporation, Marne-la-Vallée, France). Chromatographic separation was performed using a Hypercarb 3 µm, 150\*2.1 mm column (ThermoFisher, Illkirch, France), and a gradient of (A) water containing 0.5 % acetic acid, and (B) acetonitrile containing 0.5 % acetic acid as mobile phase at a flow-rate of 250 µL/min, as follows: 0.00-3.00 min, 5% (B); 3.00-5.00 min, 5 to 25% (B); 5.00-9.1 min, 25% (B); 9.1-9.2 min, 25 to 100% (B); 9.2-10.5 min, 100% (B); 10.5-10.6, 100 to 5% (B); 10.6-15 min, 5% (B). Oven temperature was set at 25°C.

LC-MS/MS analysis was performed on an LCMS-8060 triple quadrupole mass spectrometer (Shimadzu Corporation, Marne-la-Vallée, France) used in the positive electrospray ionization modes. The interface parameters and common settings were as follows: interface voltage: 1 kV; nebulizing gas flow: 3 L/min; heating gas flow: 14 L/min; drying gas flow: 10 L/min; interface temperature: 380°C; DL (desolvation line) temperature: 300 °C; heat block temperature: 500°C; collision gas pressure 170 kPa.

All compounds were measured by scheduled-MRM (Multiple Reaction Monitoring) using 1 ms pause time and 30 ms dwell time.

### *2.4 Validation procedure:*

The calibration range was 5 to 320 ng/mL for U concentration and 10 to 640 ng/ml for UH<sub>2</sub>. The calibration standards and the 3 quality controls (QC) prepared in BSA at 80 g/L were used:

(i) QC1 13 and 23 ng/ml, (ii) QC2 16 and 32 ng/ml and (iii) QC3 150 and 300 ng/ml for U and UH<sub>2</sub>, respectively.

The laboratory of Pharmacology-toxicology of the Limoges University Hospital activities has been accredited following the International Standards Organization (ISO) 15189 standard (accreditation number: 8-2607). The methods validation protocol respects also the EMA guidelines and examines the following criteria:

- The lower limit of quantitation (LLOQ) was defined as the lowest concentration that could be measured with intra-assay and inter-assay imprecision (CV%), and inaccuracy (bias) all < 20%.
- The intra-assay precision and accuracy (n=6) and the inter-assay precision and accuracy (n=6) were assessed by analysis of the 3 different QC. To assess the intra-assay precision and accuracy 6 replicates of each QC were analyzed in one day. To assess the inter-assay precision and accuracy, each QC were analyzed each day for 6 days. The CV% values and bias values had to be less than 20% at the LLOQ and less than 15% for the other levels.
- Calibration standards were prepared by adding appropriate working standard solutions to 100 µL of BSA prior to extraction in order to obtain the following 6 concentration levels for uracil: 5, 10, 15, 20, 80 or 320 ng/mL; and the following for dihydrouracil: 10, 20, 30, 40, 160 and 640 ng/mL. Calibration curves of the compounds-to-internal standard peak-area ratios of the quantification transition versus expected compound concentrations were constructed using a linear with 1/x weighting regression analysis.
- Recoveries were determined for the 3 QC by comparing analyte / internal standard peak area ratios for sample spiked either before or after the automated sample preparation.
- The carry-over effect was evaluated by analyzing the blank samples that were systematically analyzed after the injection of calibration standards at the upper limit of quantification (respectively 320 and 640 ng/ml for U and UH<sub>2</sub>). According to EMA guidelines [39] carry-over should not be greater than 20% of the LLOQ for the 2 analytes and 5% of their internal standards.
- Matrix effects were evaluated using deuterated internal standards. In fact, U and UH<sub>2</sub> are endogenous compounds. Six different human plasma and 6 different water samples were extracted and secondary spiked with internal standards at 100 ng/ml. Potential ion suppression or enhancement was evaluated by comparing the average peak intensity obtained in the two different conditions. The difference had to be less than 15%.
- Selectivity examination is not applicable for U and UH<sub>2</sub> since the analytes are naturally present in human plasma at concentrations of approximately 10 ng/mL and 100 ng/mL, respectively. However, we examined the selectivity of IS by evaluating the interference in 6

different blank plasmas. Selectivity of IS was proven by the absence of any peaks at the same retention time.

- The stability of U and UH<sub>2</sub> were studied by others [25-31]. According to Jacobs et al [31] the U and UH<sub>2</sub> are stable in whole blood for 4 hours at +4°C at maximum. Noteworthy, 13% increases in the U concentration after 4h has been observed. Both U and UH<sub>2</sub> were stable in plasma at -20°C after long term storage. This is why all samples have to be centrifuged within 1h30 and stored at -20°C.

#### *Application of the whole procedure to real patient samples:*

The automated sample preparation and LC-MS/MS analysis procedure was tested using 64 patient samples by comparing its quantitative results to those of a pre-existing and validated quantitative LC-MS/MS method routinely used in our lab and considered as a reference here after. Briefly, this method was based on a liquid-liquid extraction of 500µL of plasma with ethyl acetate/2-propanol preceded by a protein precipitation with ammonium sulfate. The chromatographic separation was performed using a XSelect®CSH™ C18, 2.5 µm (100×2.1 I.D.) column and mass detection and quantification were performed on an AB-Sciex QTRAP®4500MD triple quadrupole mass spectrometer (Sciex, Villebon sur Yvette) operated in the positive electrospray ionization (ESI+) with selected reaction monitoring mode.

### **3. Results and discussion**

As DPD deficiency screening in patients receiving a fluoropyrimidine-based chemotherapy is an obligation in France and highly recommended in other countries, most labs in charge of the measurement of U (and UH<sub>2</sub>) will or are already facing a great increase of this activity. We propose here a fully-automated solution able to ensure an accurate and robust measurement without requiring precious manpower and allowing high throughput by running samples continuously when necessary. The method was validated according to EMA guidelines. A summary of the results of validation and the MRM conditions are reported in Table 1 and in Table 2.

- Selectivity

Among validation criteria, the specificity examination is challenging since U and UH<sub>2</sub> are both endogenous compounds in human plasma and that most determination methods use BSA to build matrix matched calibrations. We examined the selectivity of the IS by evaluating



interferences in 6 different plasma. No significant signal (<5%) were observed at the RT of the IS.

An elegant approach has been proposed by Jacobs et al [31] where control human plasma was spiked with U-<sup>15</sup>N<sub>2</sub> and UH<sub>2</sub>-<sup>13</sup>C<sub>4</sub>-<sup>15</sup>N<sub>2</sub> and quantified using U and UH<sub>2</sub> as internal standards for examination of the selectivity, the recovery and the matrix effects. However, for selectivity, only the selectivity of the IS was examined but any interference from possible co-administered medications could not be verified. As patients receiving fluoropyrimidines-based chemotherapies are generally poly-medicated, this verification is necessary, though. Therefore, we examined, in our in-house built library, a list of 900 drugs and metabolites of interest to identify a potential interfering compound. None of these drugs presented a mass of 114.0 or 112.0. Collision induced dissociation ions, produced in the ESI source, with these m/z could be another source of interference. But it is still unlikely that such a compound would have the same retention time and transitions as U and HU<sub>2</sub>. In addition, the use of a column allowing a significant retention of very polar molecules, such Hypercarb, contributes to improve the selectivity of the method. Table 3 resumes the published methods assessing U and UH<sub>2</sub>.

- Linearity and limits of quantification:

Correlation coefficients values were greater than 0.99 for each standard calibration curves in the concentration range from 5 to 320 ng/ml for U and from 10 to 640 ng/ml for UH<sub>2</sub>. The LLOQs were 5 ng/ml for U and 10 ng/ml for UH<sub>2</sub>. All LLOQs were below the physiological concentrations and were quantifiable with an accuracy and a precision within 20%. ULOQs were above the concentrations observed in DPD deficient patients, namely: 320 and 640 ng/ml for U and UH<sub>2</sub> respectively.

- Carry over, recovery and matrix effect:

By inspecting blank samples used during the validation procedure, no carry over was observed after the injection of calibration standards prepared at the ULOQ. A typical example of chromatogram is presented in Figure 1 for U and UH<sub>2</sub>. The mean recoveries were very stable for all the tested concentrations: 65 to 70% for UH<sub>2</sub> and 75 to 80% for U, with low CV% values for the 2 molecules. The mean matrix effects were 10.5 and 10.2% for U <sup>13</sup>C<sup>15</sup>N<sub>2</sub> and UH<sub>2</sub> <sup>13</sup>C<sup>15</sup>N<sub>2</sub>, respectively.

- Accuracy and precision

Intra and inter-assay precisions (%RSD) varied from 2.8 to 8.8 complying with EMA validation guidelines. Intra and inter-assay accuracies varied from 85.3 to 104.6 and were also within these acceptance criteria. The inaccuracy and CV were less than 5% for intra and inter-assay tests for the threshold value 16 ng/ml.

Recent recommendations proposed threshold values of 16 and 150 ng/ml for uracilemia to characterize a partial or a complete DPD deficiency, respectively. This implies that analytical methods have to be the most accurate, reliable and robust as possible. Indeed, a suboptimal quantification, in particular around threshold values may change or delay patient's care and medical decisions. Generally, inaccuracy and CV less than 15% are the common acceptance criteria to validate each parameter of a method [39]. Unfortunately, the French recommendations do not mention any acceptable inaccuracies for the 2 thresholds. In this case, a range of quantification from 13.6 to 18.4 ng/ml could be expected for the threshold value set at 16 ng/ml. When the recommendations are strictly respected, for patient's care, such range may have different consequences: (i) Wrongly measured less than 16 ng/mL, a patient could be considered as having no DPD deficiency. The clinician may not reduce the dose of 5-FU and the patient could be overexposed with a risk of toxicity;(ii) Wrongly measured more than 16 ng/mL, a patient could be considered as having a partial DPD deficiency. The 5-FU dose may be reduced and the patient could be underexposed with a suboptimal treatment. Practically, to prevent these risks it may be recommended to re-analyze the sample or to explore the deficiency by genotyping. Both solutions delay the treatment whereas it's highly recommended that patients benefit from it as soon as they are diagnosed.

Obviously, a method with very low inaccuracy is needed to explore DPD deficiency. Even if they can provide such performances, the analytical methods reported up to now for the measurement of U and UH<sub>2</sub> are particularly long and tedious (Table 3). Most of these methods required 200 to 500 µL of plasma (except for [33]) and a preparation of the sample with a protein precipitation was always necessary before the extraction. Each also needed a liquid/liquid extraction (LLE) or solid phase extraction (SPE) with multiple manual steps including centrifugation, filtration and evaporation. In this study, we have automated the sample preparation in order to insure and maintain low accuracy. The automated sample preparation required 8 minutes. Then, after the automated transfer of the extract, the chromatographic separation was performed in 15 minutes. About 23 minutes were needed to obtain the first

result. Then, as extraction and separation were performed in parallel, the system produced a result every 15 minutes.

Once the sample is on board of the system, this approach drastically decreases the manpower dedicated to sample preparation. In addition, automation is an ideal solution allowing the analysis of patients' samples rapidly/immediately in labs willing to perform this analysis every day/at all times.

- Clinical practice

The relevance of UH<sub>2</sub>/U ratio was already reviewed by French group (ref HAS). The UH<sub>2</sub>/U ratio as a biomarker for severe toxicities under fluoropyrimidines treatment is still controversial. In addition, no consensual value could be set between the different laboratories. The main reason is the presence of interfering compounds of UH<sub>2</sub> when analyzed by LC-UV. The genotyping approach for the diagnosis of DPD deficiency focuses on only 4 variants leading to a lack of sensitivity (false negatives). Therefore, the use of U concentrations as phenotypic approach of DPD deficiency is recommended currently by the French health authorities. The validated method here is now used for routine analysis and allowed the comparison of results obtained for 64 patients with a reference method. Our results show reliable quantification. The concentrations of U in 63 patients ranged between 4.9 and 36.6 ng/mL. Seventeen patients have concentrations between 14 and 18 ng/mL. Among these, the 2 methods would have proposed different diagnosis for 5 patients (based on the threshold of 16 ng/ml). The difference varied from 1.8 to 5.6 ng/ml. In this study, 100% and 93.7% of the Bland Altman plots were within the  $\pm 2SD$  interval for U and UH<sub>2</sub> concentrations, respectively. A summary of this study is presented in Figure 2.

#### **4. Conclusion**

To our knowledge, a fully-automated LC-MS/MS method for accurate and robust quantification of U and UH<sub>2</sub> has not been published yet. This method could be used in routine labs facing a great increase in this activity since it allows high throughput by running samples continuously when necessary. Its application in routine clinical for the characterization of DPD deficiency is an ideal solution to reduce incidence of errors and thus insure best reproducibility, robustness and reliability [33, 35].

## References

- [1] R.B. Diasio, B.E. Harris, Clinical pharmacology of 5-fluorouracil, *Clin Pharmacokinet*, 16 (1989) 215-237.
- [2] C. Heidelberger, N.K. Chaudhuri, P. Danneberg, D. Mooren, L. Griesbach, R. Duschinsky, R.J. Schnitzer, E. Plevin, J. Scheiner, Fluorinated pyrimidines, a new class of tumour-inhibitory compounds, *Nature*, 179 (1957) 663-666.
- [3] F. Chionh, D. Lau, Y. Yeung, T. Price, N. Tebbutt, Oral versus intravenous fluoropyrimidines for colorectal cancer, *Cochrane Database Syst Rev*, 7 (2017) CD008398.
- [4] V.C. Jordan, A Retrospective: On Clinical Studies with 5-Fluorouracil, *Cancer Res*, 76 (2016) 767-768.
- [5] F. Lemaitre, F. Goirand, M. Launay, E. Chatelut, J.C. Boyer, A. Evrard, M.N. Paludetto, R. Guilhaumou, J. Ciccolini, A. Schmitt, [5-fluorouracil therapeutic drug monitoring: Update and recommendations of the STP-PT group of the SFPT and the GPCO-Umicancer], *Bull Cancer*, 105 (2018) 790-803.
- [6] M.A. Lorient, J. Ciccolini, F. Thomas, C. Barin-Le-Guellec, B. Royer, G. Milano, N. Picard, L. Becquemont, C. Verstuyft, C. Narjoz, A. Schmitt, C. Bobin-Dubigeon, A. Harle, A. Paci, V. Poinson, S. Quaranta, A. Evrard, B. Hennart, F. Broly, X. Fonrose, C. Lafay-Chebassier, A.S. Wozny, F. Masskouri, J.C. Boyer, M.C. Etienne-Grimaldi, [Dihydropyrimidine dehydrogenase (DPD) deficiency screening and securing of fluoropyrimidine-based chemotherapies: Update and recommendations of the French GPCO-Umicancer and RNPx networks], *Bull Cancer*, 105 (2018) 397-407.
- [7] A.B. van Kuilenburg, J.G. Maring, Evaluation of 5-fluorouracil pharmacokinetic models and therapeutic drug monitoring in cancer patients, *Pharmacogenomics*, 14 (2013) 799-811.
- [8] F. Goirand, F. Lemaitre, M. Launay, C. Tron, E. Chatelut, J.C. Boyer, M. Bardou, A. Schmitt, How can we best monitor 5-FU administration to maximize benefit to risk ratio?, *Expert Opin Drug Metab Toxicol*, 14 (2018) 1303-1313.
- [9] D. Meulendijks, A. Cats, J.H. Beijnen, J.H. Schellens, Improving safety of fluoropyrimidine chemotherapy by individualizing treatment based on dihydropyrimidine dehydrogenase activity - Ready for clinical practice?, *Cancer Treat Rev*, 50 (2016) 23-34.
- [10] J. Ciccolini, E. Gross, L. Dahan, B. Lacarelle, C. Mercier, Routine dihydropyrimidine dehydrogenase testing for anticipating 5-fluorouracil-related severe toxicities: hype or hope?, *Clin Colorectal Cancer*, 9 (2010) 224-228.
- [11] D. Meulendijks, L.M. Henricks, B.A.W. Jacobs, A. Aliev, M.J. Deenen, N. de Vries, H. Rosing, E. van Werkhoven, A. de Boer, J.H. Beijnen, C. Mandigers, M. Soesan, A. Cats, J.H.M. Schellens, Pretreatment serum uracil concentration as a predictor of severe and fatal fluoropyrimidine-associated toxicity, *Br J Cancer*, 116 (2017) 1415-1424.
- [12] M.C. Etienne-Grimaldi, J.C. Boyer, C. Beroud, L. Mbatchi, A. van Kuilenburg, C. Bobin-Dubigeon, F. Thomas, E. Chatelut, J.L. Merlin, F. Pinguet, C. Ferrand, J. Meijer, A. Evrard, L. Llorca, G. Romieu, P. Follana, T. Bachelot, L. Chaigneau, X. Pivot, V. Dieras, R. Largillier, M. Mousseau, A. Goncalves, H. Roche, J. Bonnetterre, V. Servent, N. Dohollou, Y. Chateau, E. Chamorey, J.P. Desvignes, D. Salgado, J.M. Ferrero, G. Milano, New advances in DPYD genotype and risk of severe toxicity under capecitabine, *PLoS One*, 12 (2017) e0175998.
- [13] M. Boisdron-Celle, G. Remaud, S. Traore, A.L. Poirier, L. Gamelin, A. Morel, E. Gamelin, 5-Fluorouracil-related severe toxicity: a comparison of different methods for the pretherapeutic detection of dihydropyrimidine dehydrogenase deficiency, *Cancer Lett*, 249 (2007) 271-282.
- [14] E. Gamelin, M. Boisdron-Celle, V. Guerin-Meyer, R. Delva, A. Lortholary, F. Genevieve, F. Larra, N. Ifrah, J. Robert, Correlation between uracil and dihydropyrimidine plasma ratio, fluorouracil (5-FU) pharmacokinetic parameters, and tolerance in patients with advanced colorectal cancer: A potential interest for predicting 5-FU toxicity and determining optimal 5-FU dosage, *J Clin Oncol*, 17 (1999) 1105.
- [15] M. Launay, J. Ciccolini, C. Fournel, C. Blanquicett, C. Dupuis, N. Fakhry, F. Duffaud, S. Salas, B. Lacarelle, Upfront Dpd Deficiency Detection to Secure 5-Fu Administration: Part 2- Application to Head-and-Neck Cancer Patients, *Clin Cancer Drugs*, 4 (2017) 122-128.

- [16] M. Launay, L. Dahan, M. Duval, A. Rodallec, G. Milano, M. Duluc, B. Lacarelle, J. Ciccolini, J.F. Seitz, Beating the odds: efficacy and toxicity of dihydropyrimidine dehydrogenase-driven adaptive dosing of 5-FU in patients with digestive cancer, *Br J Clin Pharmacol*, 81 (2016) 124-130.
- [17] E. Gamelin, M. Boisdron-Celle, A. Turcant, F. Larra, P. Allain, J. Robert, Rapid and sensitive high-performance liquid chromatographic analysis of halogenopyrimidines in plasma, *J Chromatogr B Biomed Sci Appl*, 695 (1997) 409-416.
- [18] M.B. Garg, J.C. Sevester, J.A. Sakoff, S.P. Ackland, Simple liquid chromatographic method for the determination of uracil and dihydrouracil plasma levels: a potential pretreatment predictor of 5-fluorouracil toxicity, *J Chromatogr B Analyt Technol Biomed Life Sci*, 774 (2002) 223-230.
- [19] J. Ciccolini, C. Mercier, M.F. Blachon, R. Favre, A. Durand, B. Lacarelle, A simple and rapid high-performance liquid chromatographic (HPLC) method for 5-fluorouracil (5-FU) assay in plasma and possible detection of patients with impaired dihydropyrimidine dehydrogenase (DPD) activity, *J Clin Pharm Ther*, 29 (2004) 307-315.
- [20] G. Remaud, M. Boisdron-Celle, C. Hameline, A. Morel, E. Gamelin, An accurate dihydrouracil/uracil determination using improved high performance liquid chromatography method for preventing fluoropyrimidines-related toxicity in clinical practice, *J Chromatogr B Analyt Technol Biomed Life Sci*, 823 (2005) 98-107.
- [21] R. Deporte, M. Amiand, A. Moreau, C. Charbonnel, L. Champion, High-performance liquid chromatographic assay with UV detection for measurement of dihydrouracil/uracil ratio in plasma, *J Chromatogr B Analyt Technol Biomed Life Sci*, 834 (2006) 170-177.
- [22] M.H. Kristensen, P. Pedersen, J. Mejer, The value of dihydrouracil/uracil plasma ratios in predicting 5-fluorouracil-related toxicity in colorectal cancer patients, *J Int Med Res*, 38 (2010) 1313-1323.
- [23] R.Z. Hahn, A.F. Galarza, A. Schneider, M.V. Antunes, G. Schwartzmann, R. Linden, Improved determination of uracil and dihydrouracil in plasma after a loading oral dose of uracil using high-performance liquid chromatography with photodiode array detection and porous graphitic carbon stationary phase, *Clin Biochem*, 48 (2015) 915-918.
- [24] F. Thomas, I. Hennebelle, C. Delmas, I. Lochon, C. Dhelens, C. Garnier Tixidre, A. Bonadona, N. Penel, A. Goncalves, J.P. Delord, C. Toulas, E. Chatelut, Genotyping of a family with a novel deleterious DPYD mutation supports the pretherapeutic screening of DPD deficiency with dihydrouracil/uracil ratio, *Clin Pharmacol Ther*, 99 (2016) 235-242.
- [25] H. Jiang, J. Jiang, P. Hu, Y. Hu, Measurement of endogenous uracil and dihydrouracil in plasma and urine of normal subjects by liquid chromatography-tandem mass spectrometry, *J Chromatogr B Analyt Technol Biomed Life Sci*, 769 (2002) 169-176.
- [26] R. Svobaite, I. Solassol, F. Pinguet, L. Ivanauskas, J. Bres, F.M. Bressolle, HPLC with UV or mass spectrometric detection for quantifying endogenous uracil and dihydrouracil in human plasma, *Clin Chem*, 54 (2008) 1463-1472.
- [27] B. Buchel, P. Rhyn, S. Schurch, C. Buhr, U. Amstutz, C.R. Largiader, LC-MS/MS method for simultaneous analysis of uracil, 5,6-dihydrouracil, 5-fluorouracil and 5-fluoro-5,6-dihydrouracil in human plasma for therapeutic drug monitoring and toxicity prediction in cancer patients, *Biomed Chromatogr*, 27 (2013) 7-16.
- [28] I.C. Cesar, G.F. Cunha-Junior, R.M. Duarte Byrro, L.G. Vaz Coelho, G.A. Pianetti, A rapid HPLC-ESI-MS/MS method for determination of dihydrouracil/uracil ratio in plasma: evaluation of toxicity to 5-fluorouracil in patients with gastrointestinal cancer, *Ther Drug Monit*, 34 (2012) 59-66.
- [29] O. Chavani, B.P. Jensen, R.M. Strother, C.M. Florkowski, P.M. George, Development, validation and application of a novel liquid chromatography tandem mass spectrometry assay measuring uracil, 5,6-dihydrouracil, 5-fluorouracil, 5,6-dihydro-5-fluorouracil, alpha-fluoro-beta-ureidopropionic acid and alpha-fluoro-beta-alanine in human plasma, *J Pharm Biomed Anal*, 142 (2017) 125-135.
- [30] F. Coudore, D. Roche, S. Lefeuvre, D. Faussot, E.M. Billaud, M.A. Lorient, P. Beaune, Validation of an ultra-high performance liquid chromatography tandem mass spectrometric method for quantifying uracil and 5,6-dihydrouracil in human plasma, *J Chromatogr Sci*, 50 (2012) 877-884.

- [31] B.A. Jacobs, H. Rosing, N. de Vries, D. Meulendijks, L.M. Henricks, J.H. Schellens, J.H. Beijnen, Development and validation of a rapid and sensitive UPLC-MS/MS method for determination of uracil and dihydrouracil in human plasma, *J Pharm Biomed Anal*, 126 (2016) 75-82.
- [32] C. Stentoft, M. Vestergaard, P. Lovendahl, N.B. Kristensen, J.M. Moorby, S.K. Jensen, Simultaneous quantification of purine and pyrimidine bases, nucleosides and their degradation products in bovine blood plasma by high performance liquid chromatography tandem mass spectrometry, *J Chromatogr A*, 1356 (2014) 197-210.
- [33] M. Vogeser, F. Kirchhoff, Progress in automation of LC-MS in laboratory medicine, *Clin Biochem*, 44 (2011) 4-13.
- [34] N. Zheng, H. Jiang, J. Zeng, Current advances and strategies towards fully automated sample preparation for regulated LC-MS/MS bioanalysis, *Bioanalysis*, 6 (2014) 2441-2459.
- [35] G.R.D. Prabhu, P.L. Urban, The dawn of unmanned analytical laboratories, *TRAC-TRENDS IN ANALYTICAL CHEMISTRY*, 88 (2017) 41\_52.
- [36] D. Wei, M. Li, K.W. King, L. Yang, Online and automated sample extraction, *Bioanalysis*, 7 (2015) 2227-2233.
- [37] T. Robin, A. Barnes, S. Dulaurent, N. Loftus, S. Baumgarten, S. Moreau, P. Marquet, S. El Balkhi, F. Saint-Marcoux, Fully automated sample preparation procedure to measure drugs of abuse in plasma by liquid chromatography tandem mass spectrometry, *Anal Bioanal Chem*, 410 (2018) 5071-5083.
- [38] R.F. Murphy, F.M. Balis, D.G. Poplack, Stability of 5-fluorouracil in whole blood and plasma, *Clin Chem*, 33 (1987) 2299-2300.
- [39] T.C.O.T.E. COMMUNITIES, implementing Council Directive 96/23/EC concerning the performance of analytical methods and the interpretation of results, in: E.u. commission (Ed.), 2002.

**Table 1:** Results of the validation study

Compound	Associated IS	Recovery (%) (%RSD) (n=3)			Intra-assay precision (%) (%RSD) (n=6)			Inter-assay precision (%) (%RSD) (n=6)			Matrix effect
		QC 1	QC 2	QC 3	QC 1	QC 2	QC 3	QC 1	QC 2	QC 3	
Dihydrouracil	Dihydrouracil <sup>13</sup> C <sup>15</sup> N <sub>2</sub>	64.8 (7.7)	68.9 (6.7)	69.6 (4.2)	85.3 (7.7)	86.5 (5.7)	101.3 (8.8)	89.5 (3.0)	93.8 (2.8)	100.6 (4.8)	10.2 (4.9)
Uracil	Uracil <sup>13</sup> C <sup>15</sup> N <sub>2</sub>	76.5 (14.2)	72.0 (5.6)	78.8 (8.9)	104.6 (4.7)	96.6 (4.4)	98.4 (3.1)	100.8 (6.6)	97.1 (3.8)	98.3 (6.0)	10.5 (2.2)

RSD: relative standard deviation; QC: quality control

**Table 2:** Optimised mass spectrometry parameters

Compounds	Precursor ion		Product ion						Retention time (min)
	m/z	Q1 pre-bias (V)	Quantitation		Reference				
			m/z	Collision energy (V)	Q3 (V)	m/z	Collision energy (V)	Q3 (V)	
Dihydrouracil	114.9	-12.0	30.1	-18.0	-13.0	55.0	-23.0	-21.0	4.26
Dihydrouracile <sup>13</sup> C <sup>15</sup> N <sub>2</sub>	118.3	-10.0	55.0	-25.0	-20.0	76.1	-14.0	-12.0	4.28
Uracil	113.1	-14.0	70.1	-19.0	-27.0	96.0	-22.0	-17.0	8.71
Uracil <sup>13</sup> C <sup>15</sup> N <sub>2</sub>	116.2	-10.0	71.1	-19.0	-28.0	98.0	-20.0	-18.0	8.71



**Table 3:** State of the art for methods assaying U and UH<sub>2</sub> by LC-MS/MS or LC-UV

Compound	Matrix	ISTD	Sample preparation	Extraction	Column	LC run (min)	Mobile phase	Detection mode	Limit of quantification (ng/ml)	Reference
U, UH <sub>2</sub> , 5FU, 5FUH <sub>2</sub> , FUPA, Fβala	Plasma 50 μL	U- <sup>13</sup> C <sup>15</sup> N <sub>2</sub> UH <sub>2</sub> - <sup>13</sup> C <sup>15</sup> N <sub>2</sub> 5FU- <sup>13</sup> C <sup>15</sup> N <sub>2</sub> FUH <sub>2</sub> - <sup>13</sup> C <sup>15</sup> N <sub>2</sub> FUPA- <sup>13</sup> C <sub>3</sub> Fβala- <sup>13</sup> C <sub>5</sub>	Prot prec: 100 μL ZnSO <sub>4</sub> 1M	L/L (*2): 600 μL 95% Acetonitrile/5% Methanol 50 μL Acetic acid 0.25 M 400 μL 15% Isopropanol/85% Ethyl acetate	Luna-PFP	8	0.25% Acetic acid/0.05% Formic acid/99.7% Water Acetonitrile	LC-MS/MS ESI+	U = 5 UH <sub>2</sub> = 10 5FU = 50 FUH <sub>2</sub> = 50 FUPA = 50 Fβala = 50	32
U, UH <sub>2</sub>	Plasma 300 μL	U- <sup>15</sup> N <sub>2</sub> UH <sub>2</sub> - <sup>13</sup> C <sub>4</sub> <sup>15</sup> N <sub>2</sub>	Prot prec: 900 μL 50% Acetonitrile/50% Methanol	Dryness, recovery and centrifugation	Acquity UHPLC HSS T3	5	0.1% Formic acid Acetonitrile/0.1% Formic acid	LC-MS/MS ESI : U : - UH <sub>2</sub> : +	U = 1 UH <sub>2</sub> = 10	31
U, UPA, βala	Plasma 300 μL	U- <sup>13</sup> C <sub>4</sub> <sup>15</sup> N <sub>2</sub> UPA βala U- <sup>13</sup> C <sub>3</sub> <sup>15</sup> N	Dilution: 4 times Prot prec: 1.8 mL Ethanol	Centrifugation and ultrafiltration	Synergi Hydro-RP	18	50 mM Ammonium acetate (pH 2.8)/10% Methanol 50 mM ammonium acetate (pH 2.8)/50% Methanol	LC-MS/MS ESI +	U = 73 UPA = 616 Bala = 1157	30
U, UH <sub>2</sub> , 5FU, 5FUH <sub>2</sub>	Plasma 200 μL	U- <sup>13</sup> C <sub>2</sub> UH <sub>2</sub> - <sup>13</sup> C <sub>4</sub> <sup>15</sup> N <sub>2</sub> 5FU- <sup>15</sup> N <sub>2</sub>	Prot prec : 200 μL H <sub>2</sub> SO <sub>4</sub> NH <sub>3</sub>	L/L (*2) : 10% Isopropanol/90% Ethyl acetate	C18 Atlantis	10	1 mM Ammonium acetate/0.5 mM formic acid/3.3% Methanol	LC-MS/MS ESI+	U = 4.1 UH <sub>2</sub> = 12.1 5FU = 13 FUH <sub>2</sub> = 100	29
U, UH <sub>2</sub>	Plasma 500 μL	U-Cl	Prot prec: 50 μL phosphonic acid	L/L: 40% Ethyl acetate/30% Dichloromethane/30% Methyl-ter-butyl-ether	C18X-Terra	5	Methanol/0.1% Sodium hydroxide (15:85)	LC-MS/MS ESI-	U=5 UH <sub>2</sub> = 10	28
U, UH <sub>2</sub>	Plasma 500 μL	U-Br	Prot prec : 300 mg H <sub>2</sub> SO <sub>4</sub> NH <sub>3</sub>	L/L : 15% Isopropanol/85% Ethyl acetate and derivation	Acquity UPLC BEH	4.5	0.5% Acetic acid 99.5% Acetonitrile/0.5% acetic acid	LC-MS/MS ESI +	U = 0.625 UH <sub>2</sub> = 0.625	27
U, UH <sub>2</sub>	Plasma 500 μL	5FC 5B <sub>3</sub> C	Prot prec: 500 μL KH <sub>2</sub> PO <sub>4</sub> 10 mM	SPE: Atoll Xtrem 500 μL Methanol for elution, dryness, 200 μL mobile phase A	C18 Atlantis	12	Water/0.1% Acetic acid Acetonitrile	LC-MS/MS ESI+ LC-UV	U = 4.1 UH <sub>2</sub> = 12.1	26
U, UH <sub>2</sub>	Plasma 200 μL	U-Br	Prot prec: 150 mg (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>	L/L : 5 mL 85% Ethyl acetate/15% Isopropanol.	Discovery Amide C16		3% Methanol	LC-MS/MS ESI -	U = 0.5 UH <sub>2</sub> = 5	25
U, UH <sub>2</sub>	Plasma 500 μL	5FU	Prot prec: 600 mg (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>	Centrifugation, dryness, 120 μL water, filtration.	Hypercarb	42	Water Acetonitrile	LC-UV	U = 1.25 UH <sub>2</sub> = 0.675	24
U, UH <sub>2</sub>	Plasma 500 μL	5FU	Prot prec: 500 mg (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>	L/L : 4 mL 85% Ethyl acetate/15% Isopropanol	Hypercarb	30	Water Acetonitrile	LC-UV	U = 0.2 UH <sub>2</sub> = 0.2	23
U, UH <sub>2</sub>	Plasma 1 mL	5-bromo-2'-deoxyuridine	Prot prec: 1 mL perchloric acid 1 M	SPE: non polar high capacity polymeric sorbent	RP-18 Supelcosil		50% Acetonitrile 0.1% trifluoroacetic acid	LC-UV		22

				One column volume Methanol for elution. Dryness. 150 µL water						
U, UH <sub>2</sub>	Plasma 500 µL	5FC	Prot prec: 500 µL KH <sub>2</sub> PO <sub>4</sub> 10 mM	SPE: Atoll Xtrem 1000 µL ammonium formate buffer (10 mM, pH 5.0) for elution, dryness, 200 µL mobile phase A	dC18 Atlantis	15	10 mM potassium phosphate buffer (pH 3.0)	LC-UV	U = 2.5 UH <sub>2</sub> = 6.75	21
U, UH <sub>2</sub>	Plasma 500 µL	FUH <sub>2</sub>	Prot prec: 600 mg (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>	L/L : 4 mL 85% Ethyl acetate/15% Isopropanol	Hypercarb	77	Water Acetonitrile	LC-UV	U = 1.25 UH <sub>2</sub> = 0.675	20
U, UH <sub>2</sub> , 5FU	Plasma 500 µL	U-Br	Prot prec: 20 µL orthophosphorique acid (5%)	L/L : 6 mL 16% Propanol/84% Ether	RP-18X Terra	35	50 mM potassium phosphate, 0.1% triethylamine	LC-UV	5FU = 5	19
U, UH <sub>2</sub>	Plasma 500 µL	U-Cl	Prot prec: 200 µL (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>	L/L : 6 mL 85% Ethyl acetate/15% Isopropanol	Spherisorb ODS2	45	10 mM potassium phosphate buffer (pH 3.0)	LC-UV	U = 2 UH <sub>2</sub> = 10	18
U, 5FU	Plasma 500 µL	U-Br U-Cl	Prot prec: 600 mg (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>	L/L : 1200 µL 85% Ethyl acetate/15% Isopropanol	Spherisorb ODS2 Phase Sep	25	10 mM potassium dihydrogen phosphate (pH 3.0)	LC-UV	5FU = 6	17

Prot prec: protein precipitation; L/L: liquid/liquid; SPE : solid phase extraction ; LC : liquid chromatography ; MS/MS : tandem mass spectrometry; UV: ultraviolet; ESI: electrospray ionization

Figure 1: Typical MRM chromatograms of uracil (bottom) and dihydrouracil (top) at the lower and upper limits of quantification (LLOQ and ULOQ). Blank samples were obtained after ULOQ.

Figure 2: Bland Altman analysis for uracil and dihydrouracil quantified with the newly procedure and a validated LC-MS method with manual sample preparation.

d: mean difference

ssd: standard deviation



