

# Performance of a multiplex polymerase chain reaction panel for identifying bacterial pathogens causing pneumonia in critically ill patients with COVID-19

François Caméléna, Anne-Clotilde Moy, Emmanuel Dudoignon, Thibaut Poncin, Benjamin Deniau, Lucie Guillemet, Jérôme Le Goff, Mélissa Budoo, Mourad Benyamina, Maïté Chaussard, et al.

## ▶ To cite this version:

François Caméléna, Anne-Clotilde Moy, Emmanuel Dudoignon, Thibaut Poncin, Benjamin Deniau, et al.. Performance of a multiplex polymerase chain reaction panel for identifying bacterial pathogens causing pneumonia in critically ill patients with COVID-19. Diagnostic Microbiology and Infectious Disease, 2021, 99, pp.115183 -. 10.1016/j.diagmicrobio.2020.115183 . hal-03493840

## HAL Id: hal-03493840 https://hal.science/hal-03493840v1

Submitted on 17 Oct 2022

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Distributed under a Creative Commons Attribution - NonCommercial 4.0 International License

1	Performance of a multiplex polymerase chain reaction panel for identifying
2	bacterial pathogens causing pneumonia in critically ill patients with COVID-
3	19
4	François Caméléna <sup>1,2</sup> *, Anne-Clotilde Moy <sup>3*</sup> , Emmanuel Dudoignon <sup>3,4</sup> , Thibaut Poncin <sup>1,2</sup> , Benjamin
5	Deniau <sup>3,4</sup> , Lucie Guillemet <sup>3</sup> , Jérôme Le Goff <sup>1</sup> , Mélissa Budoo <sup>1</sup> , Mourad Benyamina <sup>3</sup> , Maïté
6	Chaussard <sup>3</sup> , Maxime Coutrot <sup>3</sup> , Matthieu Lafaurie <sup>5</sup> , Benoît Plaud <sup>3,4</sup> , Alexandre Mebazaa <sup>3,4</sup> , François
7	Depret <sup>3,4</sup> and Béatrice Berçot <sup>1,2</sup>
8	
9	*these authors contributed equally to this work.
10	
11	<sup>1</sup> Department of Microbiology, Saint-Louis-Lariboisière Hospital Group, Assistance Publique - Hôpitaux
12	de Paris, Paris, France;
13	<sup>2</sup> University of Paris, INSERM 1137, IAME, Paris, France;
14	<sup>3</sup> Department of Anesthesiology and Critical Care and Burn Unit, Saint-Louis-Lariboisière Hospital Group,
15	Assistance Publique - Hôpitaux de Paris, Paris, France.
16	<sup>4</sup> University of Paris, INSERM 942, INI-CRCT network, Paris, France.
17	<sup>5</sup> Department of Infectious Disease, Saint-Louis-Lariboisière Hospital Group, Assistance Publique -
18	Hôpitaux de Paris, Paris, France.
19	
20	Corresponding author: Prof. Béatrice Berçot
21	Department of Microbiology, Saint-Louis-Lariboisière Hospital Group, Assistance Publique - Hôpitaux de
22	Paris, Paris, France
23	Tel +33 1 42 38 50 96; Fax +33 1 42 49 92 00; E-mail: beatrice.bercot@aphp.fr

- 24 ABSTRACT
- 25

The FilmArray® Pneumonia Plus (FA-PP) panel can provide rapid identifications and semi-quantitative 26 27 results for many pathogens. We performed a prospective single-center study in 43 critically ill patients with COVID-19, in which we performed 96 FA-PP tests and cultures of blind bronchoalveolar lavage 28 29 (BBAL). FA-PP detected one or more pathogens in 32% (31/96 of samples), whereas culture methods 30 detected at least one pathogen in 35% (34/96 of samples). The most prevalent bacteria detected were Pseudomonas aeruginosa (n=14) and Staphylococcus aureus (n=11) on both FA-PP and culture. The FA-31 PP results from BBAL in critically ill patients with COVID-19 were consistent with bacterial culture 32 findings for bacteria present in the FA-PP panel; showing sensitivity, specificity, positive and negative 33 34 predictive value of 95%, 99%, 82% and 100%, respectively. Median turnaround time for FA-PP was 5.5 35 hours, which was significantly shorter than for standard culture (26 hours) and AST results (57 hours).

36

37 Keywords: FilmArray Pneumonia plus, SARS-CoV-2, Bacterial pneumonia, Rapid diagnostics,

#### 39 1. Introduction

40 Severe acute respiratory syndrome coronavirus-2 (SARS-CoV-2), which causes coronavirus disease 2019 (COVID-19), has rapidly spread worldwide. Several studies have reported complications of COVID-19, 41 42 such as bacterial pneumonia, acute respiratory distress syndrome (ARDS) and multiple organ failure 43 syndromes [1–3]. Recent guidelines for the management of adults critically ill with COVID-19 have 44 suggested the empiric use of antimicrobial agents in patients with respiratory failure [4]. The accurate and timely diagnosis of bacterial pneumonia, particularly in cases of hospital-acquired pneumonia (HAP) and 45 46 ventilator-associated pneumonia (VAP), is particularly challenging, and this condition remains a major 47 cause of morbidity and mortality [5,6]. Molecular tests provide a rapid turnaround time (TAT), together with identifications and semi-quantitative results for many pathogens responsive to antimicrobial therapy. 48 Multiplex testing may provide additional information concerning the presence of antibiotic resistance 49 50 genes, thereby improving antibiotic management [7,8]. We performed a prospective single-center study on 51 critically ill patients with COVID-19, in which we conducted parallel tests of blind bronchoalveolar lavage (BBAL) by conventional culture and FilmArray® Pneumonia Plus (FA-PP) panel (BioFire, Salt 52 53 Lake City, US). The aim of this study was to evaluate the performance of FA-PP and to compare its TAT 54 with that of conventional cultures.

55

#### 56 2. Materials and Methods

#### 57 2.1. Study population

We conducted a prospective single-center study of consecutive patients admitted to the ICU at the Saint-Louis Hospital (Assistance Publique - Hôpitaux de Paris, Paris, France) between March 22, 2020 and April 15, 2020, for COVID-19 complicated by respiratory failure. We included all patients on invasive mechanical ventilation for whom BBAL was performed. The diagnosis of COVID-19 was based on the detection of SARS-CoV-2 from nasopharyngeal samples. This detection relies on viral RNA amplification by using RT-PCR (RealStar Altona®) to amplify SARSCoV2 E gene. Patients with the following criteria were excluded from the analysis: aged under 18 years, pregnant or dying within 48 hours of admission.
For each patient, BBAL was performed immediately after intubation (or on admission to the ICU if the patient was already intubated) and in cases in which VAP was suspected. We recorded clinical, laboratory and antimicrobial therapy data. The study was approved by the local ethics committee (the ethics committee of *Société Française d'Anesthésie Réanimation* IRB 00010254-2019-203).

69

70

#### 2.2. Conventional microbiological analysis

71 BBAL was performed with a sterile catheter and a specimen trap kit. The suction tube was introduced 72 blindly through the endotracheal tube and wedged into the tracheobronchial tree before suction; 50 mL of saline were successively injected and aspirated twice, and the contents of the second aspiration were 73 analyzed, the contents of the first being discarded [9]. Samples were transported immediately to the 74 laboratory. Samples received from 6 pm to 8 am were conserved at +4°C and tested the following 75 76 morning. We tested the BBAL samples in parallel by FA-PP and conventional culture, and the results were released as soon as they become available. The decision to use FA-PP results for the management of 77 78 antimicrobial therapy was left to the discretion of the treating physician.

79

We plated 10 µL of BBAL on agar media (BioMérieux, Marcy l'Etoile, France) including blood, colistin-80 81 nalidixic acid blood and chocolate agars, and incubated the plates in ambient air or under an atmosphere enriched with 5% CO<sub>2</sub> at 35°C. No Gram staining was performed on the BBAL, to protect laboratory 82 staff, as this procedure requires centrifugation, which may generate aerosols. Bacteria were semi-83 quantified by calculating colony-forming units (CFU) per ml (CFU/ml). Culture results were interpreted 84 85 according to the recommendations of the European Society for Clinical Microbiology and Infectious Diseases [10], after observation at 18-24 hours, and negative plates were incubated for a further two days, 86 to ensure an absence of bacterial growth, before being classified definitively as negative. A cut-off of  $10^4$ 87 CFU/ml for significant growth of potential pathogens for BBAL was used. Bacteria were identified by 88

89	using MALDI-TOF MS (VITEK-MS, BioMérieux, Marcy-l'Étoile, France), according to the					
90	manufacturer's instructions. Antimicrobial susceptibility testing (AST) was performed by the disk					
91	diffusion method or with the VITEK 2 (BioMérieux, Marcy-l'Étoile, France), in accordance with the					
92	recommendations of the European Committee on Antimicrobial Susceptibility Testing.					
93						
94	2.3. Molecular investigation					
95	As recommended by the manufacturer, a dry swab was dipped into the BBAL and then immersed in					
96	dilution buffer before injection of the fluid into the FA-PP cartridge. The analysis was carried out in 67					
97	minutes on the Filmarray® apparatus (version 2.1.336.0), with the Pneumoplus v2.0 pouch system					
98	(version 2.1.0.5).					
99						
100	2.4. Statistical analysis					
101	Sensitivity (Se), specificity (Sp), positive and negative predictive values (PPV and NPV), likelihood ratio					
102	and $\kappa$ index were calculated by comparing the results for conventional culture with those of FA-PP only					
103	for bacterial pathogens present in the FA-PP panel. Performance was measured only for bacterial analysis,					
104	with bacterial culture used as the gold standard reference method. Two-tailed Mann-Whitney tests were					
105	used to compare TAT between the two methods.					
106						

#### 107 **3. Results**

In total, 96 samples from 43 patients were tested in parallel by FA-PP and conventional BBAL culture.
The median age of the patients was 64 years, 86% were male, 27% were obese and 72% had hypertension
(Table 1). On admission to the ICU, 37 (86%) patients were already on invasive mechanical ventilation.
BBAL were performed under antibiotic treatment in 55% (53/96) of samples.

FA-PP detected one or more pathogens in 32% (31/96) of samples (corresponding to 18 patients), whereas
culture methods detected at least one pathogen in 35% (34/96) of samples (corresponding to 20 patients).

114 At first BBAL, 25% (10/43) of FA-PP tests and 21% (9/43) of cultures gave positive results. The most prevalent bacteria detected were Staphylococcus aureus (n=5 on both FA-PP and culture) followed by 115 Haemophilus influenzae (n=4 and n=2 on FA-PP and culture, respectively) (Table 2). In cases in which 116 others BBAL were performed, 40% (21/53) of samples tested positive by FA-PP whereas 47% (25/53) 117 118 tested positive on culture. The most prevalent bacteria detected was *Pseudomonas aeruginosa* (n=11 on 119 both FA-PP and culture), followed by S. aureus (n=6 on both FA-PP and culture) (Table 2). The performance of FA-PP for each pathogen is shown in Table 2. Of the 43 patients, 11/43 (25.5%) had only 120 121 one sampling, and in each case, a negative result was found for both FA-PP tests and cultures. For the 122 remaining 32/43 (74.5%) patients, iterative sampling was performed; all FA-PP test and culture results for 123 these patients are presented in Supplementary Table S1. Neither atypical bacteria nor viruses were 124 detected by FA-PP. The pathogen most frequently detected by FA-PP in the absence of positive cultures was *H. influenzae* (*n*=3). Positive cultures were obtained for eight samples that tested negative by FA-PP. 125 In six of these eight samples, the pathogens present were not included in the FA-PP panel (Citrobacter 126 127 koseri in 4 cases, Stenotrophomonas maltophilia in one case and Enterococcus faecium in one case). S. aureus and Klebsiella aerogenes were recovered from single culture from two samples that tested also 128 negative by FA-PP. Bacterial load was 10<sup>3</sup> colony-forming units (CFU)/mL of culture for both these 129 130 pathogens, corresponding to a low organism burden. We found that 28 of the bacteria of the FA-PP panel grew to significant levels (*i.e.*  $\ge 10^4$  CFU/mL) in culture, with 100% (28/28) presenting  $\ge 10^4$  copies of 131 bacterial nucleic acid/mL on FA-PP. Another 10 bacteria grew to non-significant levels (*i.e.*  $< 10^4$ 132 CFU/mL), and 80% (8/10) of these bacteria were found to be present at a level  $\geq 10^4$  copies/mL by FA-133 134 PP. FA-PP detected two *bla*<sub>CTX-M</sub> extended spectrum  $\beta$ -lactamase (ESBL), one *bla*<sub>VIM</sub> carbapenemase and 135 two S. aureus methicillin resistance genes. AST detected the presence of an ESBL encoded by a bla<sub>CTX-M</sub> 136 gene in Escherichia coli and methicillin-resistant S. aureus (MRSA) in two isolates corresponding to two cases of mecA/C detection by FA-PP. The detection of the other resistance genes, one bla<sub>CTX-M</sub> and one 137  $bla_{\rm VIM}$ , remained unexplained because no other viable bacteria with a phenotype corresponding to the 138 139 presence of these genes were obtained in culture.

Median TAT was significantly shorter for FA-PP tests than for standard culture and AST results (5.5 hours [IQR: 3.4-12.4] vs. 25.9 hours [IQR: 22.5-35.7] and 57 hours [IQR: 47.7-70.8], p<0.001 and p<0.001, respectively).

143

#### 144 4. Discussion

145 We studied 43 critically ill patients with COVID-19 who underwent BBAL at least once for suspected bacterial pneumonia. A few cases of bacterial pneumonia have been reported in critically ill patients with 146 COVID-19. The diagnosis of bacterial pneumonia in such patients remains challenging, as obtaining distal 147 respiratory samples by bronchoscopy exposes intensivists to a high risk of contamination [1,2,11]. We, 148 149 therefore, performed BBAL, which is simple, suitable for bacterial culture and does not require bronchoscopy [9]. Multiplex syndromic panels for the detection of bacterial pathogens responsible for 150 151 pneumonia have recently been developed and shown to yield more sensitive results more rapidly than 152 conventional culture, particularly in cases in which antibiotics were administered before sampling [6–8]. However, there are currently no recommendations for routine molecular testing and multiplex PCR is 153 rarely used for detection of the bacteria typically responsible for pneumonia. Questions about the clinical 154 usefulness of this approach remain, particularly in the context of critically ill patients with COVID-19. 155

156 This is, to our knowledge, the first study to evaluate, in real-time, FA-PP relative to culture, using exclusively distal respiratory samples for the diagnosis of bacterial pneumonia in critically ill patients. It 157 158 shows that the results of FA-PP are consistent (Se 95%, Sp 99%, PPV 82% and NPV 100%) with those of conventional culture for the bacterial pathogens present in the FA-PP panel. FA-PP significantly 159 decreased TAT and increased the number of bacterial detections in critically ill patients with COVID-19. 160 However, attention must be paid to bacteria not present in the FA-PP panel (e.g. Citrobacter spp.) that 161 might affect antibiotic therapy and to calculations of FA-PP performance that were performed 162 163 independently of the consideration of bacterial thresholds and antibiotic administration prior to sampling.

Bacterial cultures of BBAL grew with significant amount of bacteria (ie,  $\geq 10^4$  CFU/mL) in 37% (n =16) of patients. Among the 16 bacterial pneumonia, 4 (25%) were classified as community-acquired pneumonia, 1 (6%) as HAP, and 11 (69%) as VAP. In our study, 6% (6/96) of the samples were found to contain a bacterium not present in the FA-PP panel, and this bacterium was present at a high abundance (*i.e.*  $\geq 10^4$  CFU/mL) in four of these samples. These bacteria were involved in 1 HAP and 3 VAP.

169 Furthermore, FA-PP may provide a major argument for very early de-escalation and narrowing of 170 antimicrobial therapy in specific situations (such as to withdraw anti-MRSA), if negative FA-PP results 171 are obtained, or for the faster initiation of appropriate antimicrobial therapy, if positive FA-PP results are 172 obtained. We also found that most of the positive FA-PP results from culture-negative samples corresponded to H. influenzae; the presence of this bacterium may reflect contamination by the 173 174 oropharyngeal flora. Antimicrobial therapy is also known to impact bacterial growth, as 88% (7/8) of the false positive FA-PP results were obtained from BBAL performed under antibiotic treatment. Similarly, 175 63% (5/8) of bacteria, present at a level  $\geq 10^4$  copies/mL by FA-PP and grew at non-significant level in 176 culture, were from BBAL performed under antimicrobial therapy. Furthermore, others studies have 177 178 described that the bacterial burden could be overestimate by FA-PP compared to culture [7,8]. In addition, 179 resistance genes detected by FA-PP were not confirmed by AST results in 2 samples suggesting 180 limitations to predict phenotypic susceptibility from molecular tests. These results raise important point 181 that PCR can detect both dead and viable bacteria [7,8,12]. It remains unclear for how long bacterial loads 182 remain detectable after the initiation of appropriate antibiotic treatment in critically ill patients with COVID-19 and whether the monitoring of bacterial load by FA-PP would be useful in these patients and, 183 more broadly, in all patients with suspected VAP. Further studies are required to investigate this issue. 184 185 Additionally, FA-PP has a cost per test much greater than the costs of routine culture; future studies should focus on how best to use this test in a cost-effective manner. 186

187 Finally, the strengths of this study are the well-characterized population of patients (COVID-19 in ICU),188 the use of a single type of specimen (BBAL) and prospective testing by FA-PP and culture in parallel. Its

189 limitations are that it was a single-center study and that the COVID-19 population may not be190 representative of the ventilated patients generally treated in the ICU

### 191 5. Conclusion

In conclusion, in this prospective study of critically ill patient with COVID-19, FA-PP results were in agreement with standard culture for identifying bacterial pathogens in BBAL specimens, with a significant decrease of TAT. The use of FA-PP may enable a more complete and faster approach for the diagnosis of bacterial pneumonia in patient hospitalized in ICU. Further studies are required to address the clinical impact of FA-PP on the management of VAP in critically ill patients.

197

#### 198 Acknowledgments

199 The authors thank the intensive care teams and the clinical microbiology staff of Saint-Louis Hospital for

their contributions. They also thank Samuel Amar for his help with the collection of microbiological data.

201

#### 202 **Declarations**

- **Funding:** This work was supported by internal funding.
- 204 **Competing Interests**: F.C. received conference invitations from BioMérieux.
- 205 Ethical Approval: The study was approved by the local ethics committee (the ethics committee
- of Société Française d'Anesthésie Réanimation IRB 00010254-2019-203).

#### 207 **REFERENCES**

- 208 [1] Lescure F-X, Bouadma L, Nguyen D, Parisey M, Wicky P-H, Behillil S, et al. Clinical and
- 209 virological data of the first cases of COVID-19 in Europe: a case series. Lancet Infect Dis
- 210 2020;2:1–10. doi:10.1016/S1473-3099(20)30200-0.
- 211 [2] Rodriguez-Morales AJ, Cardona-Ospina JA, Gutiérrez-Ocampo E, Villamizar-Peña R, Holguin-
- 212 Rivera Y, Escalera-Antezana JP, et al. Clinical, laboratory and imaging features of COVID-19: A
- systematic review and meta-analysis. Travel Med Infect Dis 2020;101623.
- doi:10.1016/j.tmaid.2020.101623.
- Shen C, Wang Z, Zhao F, Yang Y, Li J, Yuan J, et al. Treatment of 5 Critically Ill Patients With
  COVID-19 With Convalescent Plasma. JAMA 2020:1–8. doi:10.1001/jama.2020.4783.
- 217 [4] Alhazzani W, Møller MH, Arabi YM, Loeb M, Gong MN, Fan E, et al. Surviving Sepsis
  218 Campaign: guidelines on the management of critically ill adults with Coronavirus Disease 2019

219 (COVID-19). Intensive Care Med 2020:1–41. doi:10.1007/s00134-020-06022-5.

- [5] Nair GB, Niederman MS. Ventilator-associated pneumonia: present understanding and ongoing
  debates. Intensive Care Med 2015;41:34–48. doi:10.1007/s00134-014-3564-5.
- Papazian L, Klompas M, Luyt C-E. Ventilator-associated pneumonia in adults: a narrative review.
  Intensive Care Med 2020;s00134-20. doi:10.1007/s00134-020-05980-0.
- 224 [7] Lee SH, Ruan S-Y, Pan S-C, Lee T-F, Chien J-Y, Hsueh P-R. Performance of a multiplex PCR
- pneumonia panel for the identification of respiratory pathogens and the main determinants of
- 226 resistance from the lower respiratory tract specimens of adult patients in intensive care units. J
- 227 Microbiol Immunol Infect 2019;52:920–8. doi:10.1016/j.jmii.2019.10.009.
- 228 [8] Yoo IY, Huh K, Shim HJ, Yun SA, Chung YN, Kang OK, et al. Evaluation of the BioFire®
- 229 FilmArray® Pneumonia Panel for rapid detection of respiratory bacterial pathogens and antibiotic

230		resistance genes in sputum and endotracheal aspirate specimens. Int J Infect Dis 2020;S1201-9712.
231		doi:10.1016/j.ijid.2020.03.024.
232	[9]	Mentec H, May-Michelangeli L, Rabbat A, Varon E, Le Turdu F, Bleichner G. Blind and
233		bronchoscopic sampling methods in suspected ventilator-associated pneumonia. A multicentre
234		prospective study. Intensive Care Med 2004;30:1319-26. doi:10.1007/s00134-004-2284-7.
235	[10]	Cornaglia G, Courcol R, Herrmann J. European Manual of Clinical Microbiology. European
236		Society for Clinical Microbiology and Infectious Diseases, Basel; 2012.
237	[11]	Bhatraju PK, Ghassemieh BJ, Nichols M, Kim R, Jerome KR, Nalla AK, et al. Covid-19 in
238		Critically Ill Patients in the Seattle Region - Case Series. N Engl J Med 2020:1-11.
239		doi:10.1056/NEJMoa2004500.
240	[12]	Clavel M, Barraud O, Moucadel V, Meynier F, Karam E, Ploy M-C, et al. Molecular quantification
241		of bacteria from respiratory samples in patients with suspected ventilator-associated pneumonia.
242		Clin Microbiol Infect 2016;22:812.e1-812.e7. doi:10.1016/j.cmi.2016.06.013.

## 244 TABLE AND FIGURE

245 Table 1. Main characteristics of patients on admission to the ICU.

Characteristics	Patients (n=43)
Median age, years (IQR)	64 [56-70]
Sex, male, no. (%)	37 (86)
Median BMI, kg per m <sup>2</sup> (IQR)	27 [26-30]
Comorbidities, no. (%)	
Hypertension	31 (72)
Heart failure	1 (2.3)
ARBS	17 (39.5)
Coronary disease	5 (11.6)
Diabetes mellitus	16 (37.0)
Pulmonary disease	7 (16.0)
Immunosuppression	3 (7.0)
Severity of illness	
Duration from onset of symptoms to ICU admission, days (IQR)	8.5 [5.3-11.6]
SAPS II (IQR)	37 [28-49]
SOFA score on day 1 (IQR)	5 [4-7]
Acute kidney injury, no. (%)	28 (65.1)
KDIGO 1	6 (13.9)
KDIGO 2	6 (13.9)
KDIGO 3	16 (34.9)
Renal replacement therapy, no. (%)	8 (18.6)
Mechanical ventilation on ICU admission, no. (%)	37 (86)
ARDS, no. (%)	43 (100)
Ratio of $PaO_2$ to $FiO_2$ , mmHg ( $PaO_2/FiO^2$ ) at admission (IQR)	136 [82.5-174.5]
Ratio of PaO <sub>2</sub> to FiO <sub>2</sub> nadir, mmHg (IQR)	82 [70-103]
Median number of prone position ventilation sessions (IQR)	1 [0-2]
Duration of mechanical ventilation, days (IQR)	10 [7-12]
Use of noradrenaline during the first 48 hours, no. (%)	36 (83.7)
Length of stay in ICU, days (IQR)	11 [8-13]
Death, no. (%)	12 (27.9)
Admission biological data	
WBC count, x $10^9$ cells per L (IQR)	9.52 [5.2-12.4]
Lymphocyte count, x $10^9$ cells per L (IQR)	0.79 [0.5-1.1]
Serum creatinine concentration, µmol per L (IQR)	75 [57-100.5]

D-dimers, ng per mL (IQR)	2610 [1490-4820]		
Radiological data	logical data       4 [2-4]         X-ray, number of dials affected at admission       4 [2-4]         CT during hospitalization, no. (%)       29 (67.4)         (10-25]*       2 (6.9)         (25-50]*       9 (31.0)         (50-75]*       10 (34.5)         >75*       8 (27.6)		
X-ray, number of dials affected at admission	4 [2-4]		
CT during hospitalization, no. (%)	29 (67.4)		
(10-25]*	2 (6.9)		
(25-50]*	9 (31.0)		
(50-75]*	10 (34.5)		
>75*	8 (27.6)		
Antimicrobial treatment, no. (%)			
Antibiotic therapy before admission	29 (67.4)		
3 <sup>rd</sup> generation cephalosporin	22 (51.1)		
Macrolide	21 (48.3)		
Other antibiotics	11 (25.6)		

246 **BMI**: Body mass index; **ICU**: Intensive care unit; **SAPS II**: Simplified Acute Physiology Score II;

247 SOFA: Sepsis-Related Organ Failure Assessment; KDIGO: Kidney Disease Improving Global

248 Outcomes; ARDS: Acute respiratory distress syndrome; ARBS: Angiotensin-renin-aldosterone

blocker system; WBC. White blood cells; CT: computerized tomography; \*Specific COVID-19
 pulmonary lesions

Bacterial pathogen	FA-PP+ & Culture+	FA-PP+ & Culture-	FA-PP- & Culture+	FA-PP- & Culture-	Se (95% CI)	Sp (95% CI)	PPV (95% CI)	NPV (95% CI)	LR+ (95% C
Acinetobacter baumannii	0	0	0	96	-	1.00 (-)	-	1.00 (-)	-
Staphylococcus aureus	10	1	1	84	0.91 (0.58-1.0)	0.99 (0.94 - 1.0)	0.91 (0.59 - 0.98)	0.99 (0.93 - 1.0)	77 (11 - 54
Streptococcus pneumoniae	3	1	0	92	1.00 (0.29 - 1.0)	0.99 (0.94 - 1.0)	0.75 (0.30-0.95)	1.00	93 (13 - 65
Escherichia coli	3	0	0	93	1.00 (0.29 - 1.0)	1.00 (0.96 - 1.0)	1.00	1.00	-
Enterobacter cloacae	0	1	0	95	-	0.99 (0.97 - 1.0)	0.00	1.00	-
Klebsiella aerogenes	3	0	1	92	0.75 (0.19 - 0.99)	1.00 (0.96 - 1.0)	1.00 (-)	0.99 (0.94 - 1.0)	-
Klebsiella pneumoniae	0	1	0	95	-	0.99 (0.97 - 1.0)	0.00 (-)	1.00	-
Klebsiella oxytoca	0	0	0	96	-	1.00 (-)	-	1.00	-
Pseudomonas aeruginosa	14	0	0	82	1.00 (0.77 - 1.0)	1.00 (0.96 - 1.0)	1.00 (-)	1.00	-
Haemophilus influenzae	2	3	0	91	1.00 (0.16 - 1.0)	0.97 (0.91 - 0.99)	0.40 (0.18 - 0.67)	1.00 (-)	31 (10 - 95
Streptococcus agalactiae	1	0	0	95	1.00 (0.03 - 1.0)	1.00 (0.96 - 1.0)	1.00 (-)	1.00 (-)	-
Streptococcus pyogenes	0	0	0	96	-	1.00 (-)	-	1.00 (-)	-
Serratia marcescens	0	0	0	96	-	1.00 (-)	-	1.00 (-)	-
Proteus spp.	0	1	0	95	-	0.99 (0.97 - 1.0)	0.00 (-)	1.00 (-)	-
Moraxella catarrhalis	0	0	0	96	-	1.00 (-)	-	1.00	-
Total (per analysis)	36	8	2	1394	0.95 (0.82 - 0.99)	0.99 (0.99 - 1.0)	0.82 (0.69 - 0.90)	1.00 (0.99 - 1.0)	166 (83 - 33

### 252 Table 2 Performance of FA-PP relative to conventional culture

253 FA-PP: FilmArray Pneumonia Plus; Se: Sensitivity; Sp: Specificity; PPV: Positive predictive value; NPV:

254 Negative predictive value; LR: Likelihood ratio