



HAL
open science

Evaluation of 4 quantification methods for monitoring 16 antibiotics and 1 beta-lactamase inhibitor in human serum by high-performance liquid chromatography with tandem mass spectrometry detection

Patrick Seraissol, Thomas Lanot, Sarah Baklouti, Camille Mané, Stéphanie Ruiz, Michel Lavit, Pascale de Riols, Jean-Christophe Garrigues, Peggy Gandia

► To cite this version:

Patrick Seraissol, Thomas Lanot, Sarah Baklouti, Camille Mané, Stéphanie Ruiz, et al.. Evaluation of 4 quantification methods for monitoring 16 antibiotics and 1 beta-lactamase inhibitor in human serum by high-performance liquid chromatography with tandem mass spectrometry detection. *Journal of Pharmaceutical and Biomedical Analysis*, 2022, 219, pp.114900. 10.1016/j.jpba.2022.114900 . hal-03814985

HAL Id: hal-03814985

<https://hal.science/hal-03814985v1>

Submitted on 22 Jul 2024

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Distributed under a Creative Commons Attribution - NonCommercial 4.0 International License

1 **Evaluation of 4 quantification methods for monitoring 16 antibiotics and 1 beta-**
2 **lactamase inhibitor in human serum by high-performance liquid chromatography with**
3 **tandem mass spectrometry detection**

4

5 Patrick Seraisol^{a,b}, Thomas Lanot^a, Sarah Baklouti^a, Camille Mané^a, Stéphanie Ruiz^c, Michel
6 Lavit^a, Pascale De Riols^a, Jean-Christophe Garrigues^{d*}, Peggy Gandia^{a,e}

7

8 ^a *Laboratoire de Pharmacocinétique et Toxicologie, IFB, Hôpital Purpan, 330 Avenue de*
9 *Grande-Bretagne, 31059 Toulouse Cedex 9, France.*

10 ^b *Laboratoire Départemental 31, Eau – Vétérinaire – Air, 76 chemin de Boudou, CS 50013,*
11 *31140 Launaguet, France.*

12 ^c *Service de Réanimation Polyvalente Adulte, Hôpital Rangueil, 1 avenue du Professeur Jean*
13 *Poulhès, 31059 Toulouse Cedex 9, France*

14 ^d *Laboratoire des IMRCP, Université de Toulouse, UMR 5623, 118 Route de Narbonne,*
15 *31062 Toulouse Cedex 9, France.*

16 ^e *INTHERES, Université de Toulouse, INRA, ENVT, 23 Chemin des Capelles, BP 87614, 31*
17 *076, Toulouse Cedex 3, France.*

18

19 * Corresponding author : Jean-Christophe Garrigues

20 Laboratoire des IMRCP (CNRS UMR 5623)

21 118 route de Narbonne,

22 31062 Toulouse cédex , France

23 Email : garrigues@chimie.ups-tlse.fr

24 Phone : +33 (0) 5 61 55 62 69

25

26 **Abstract**

27 Antibiotic (ATB) prescription in an intensive care unit (ICU) requires continuous monitoring
28 of serum dosages due to the patient's pathophysiological condition. Dosing adjustment is
29 necessary to achieve effective targeted concentrations. Since ICUs routinely use a large
30 number of ATBs, global monitoring needs to be developed. In the present study, we
31 developed a global analytical method for extracting, separating and quantifying the most
32 widely used ATBs in ICUs: amoxicillin, piperacillin, cefazolin, cefepime, cefotaxime,
33 ceftazidime, ceftolozane, ceftriaxone, ertapenem, meropenem, ciprofloxacin, moxifloxacin,
34 levofloxacin, daptomycin, dalbavancin, linezolid and a beta-lactamase inhibitor: tazobactam.
35 To guarantee the robustness of the quantification, we differentiated the 16 ATBs and the beta
36 lactamase inhibitor into 4 pools (ATB1 to ATB4), taking into account prescription frequency
37 in the ICU, the physicochemical properties and the calibration ranges of the ATBs selected.
38 The whole ATB was then separated with two LC columns in reversed phase: Kinetex Polar-
39 C18 100 Å and Polar-RP-80 synergy, in less than 6.5 min. Detection was carried out by
40 electrospray in positive ion mode, by tandem mass spectrometry (LC-MS/MS). The four
41 quantification methods were validated according to the European guidelines on bioanalytical
42 method validation (EMA guide), after determining the extraction yields, matrix effects,
43 recovery, precision, accuracy, within-run precision and between-run precision. For all
44 analyses, bias is <15% and is comparable to the literature and LOQs vary from 0.05 mg.L⁻¹
45 for ciprofloxacin to 1.00 mg.L⁻¹ for ceftriaxone and dalbavancin. The stability time of
46 cefepime and piperacillin is 3 hrs and for the other ATBs 6 hrs in serum at room temperature.
47 For long-term stability, freezing at -80°C guarantees 3 months of stability for ceftriaxone and
48 dalbavancin and more than 6 months for the other ATBs.

49 **Keywords**

50 LC–MS/MS, Antibiotics, Therapeutic drug monitoring, Critical care medicine, Beta-
51 lactamase, Sample preparation

52

53 **Highlights**

- 54 • Development of a monitoring method for 16 ATBs and beta lactamase inhibitor
- 55 • The 16 ATBs were grouped into 4 pools and quantified by specific methods
- 56 • Application of the method to 51 plasma samples from ICU patients
- 57 • Serum stability of Cefepime and Piperacillin is 3 hrs at room temperature
- 58 • The long-term stability of all ATBs is greater than 3 months at -80°C

59

60 **1. Introduction**

61 In intensive care units (ICUs), infection is a major problem due mainly to nosocomial
62 infections. Many antibiotics from different families are then administered to treat patients.
63 The various drug dosage regimens applied to hospitalized patients require precise monitoring
64 of the antibiotic concentrations in plasma [1]. Thus, accurate and precise quantification of
65 antibiotics in plasma is needed in order to avoid underexposure or overexposure of patients,
66 which could lead to inefficacy and bacterial resistance or toxicity, respectively. This is further
67 compounded by the frequent pathophysiological changes observed in critically ill patients,
68 which lead to modifications of antibiotic pharmacokinetics (PK) [2–4]. Nowadays, the
69 majority of analytical methods used to perform therapeutic drug monitoring (TDM) are
70 dedicated to the homogeneous chemical family of ATBs (e.g. β -lactams [5–8],
71 fluoroquinolones [9,10], oxazolidinones [11]). These analytical methods are essentially liquid
72 chromatography techniques with UV detection. The methods involve specific sample
73 preparation for each class of ATBs and do not allow rapid monitoring. Moreover, large-scale
74 studies have shown that 16% of patients in ICUs have effective plasma concentrations below

75 the pharmacokinetic/pharmacodynamic (PK/PD) targets [12]. Thus, there is a need for precise
76 therapeutic drug monitoring (TDM) of antibiotics in plasma in order to individualize the
77 treatment regimen [1,13,14].

78 The limit of these analytical methods is found in the specific wavelengths used for each
79 therapeutic class, with UV quantification and dedicated sample treatment for specific
80 chemical families of ATBs. These methods are developed for only one chemical class, so they
81 are time-consuming and not adapted to ICU specifications. Developments in mass
82 spectrometry and tandem mass spectrometry have been published for the analysis of β -
83 lactams [15][16], peptides [17] and multi-class components [18], showing the possibility of
84 developing a comprehensive method, from extraction to separation, and optimizing the
85 detection parameters.

86 A few methods have been developed to ensure monitoring of a large panel of ATBs
87 commonly used in the ICU and are based on detection by high resolution mass spectrometry
88 (HRMS), a system that is rarely used in hospitals [19]. To take of all these constraints into
89 account, the present study focuses on developing an analytical methodology applied to
90 quantifying of 16 antibiotics and 1 beta-lactamase inhibitor. This methodology uses four
91 quantification methods based on reversed phase separation with tandem mass spectrometry
92 detection for monitoring the most commonly used ATBs in our ICU: two penicillins:
93 amoxicillin and piperacillin; six cephalosporins: cefazolin, cefepime, cefotaxime, ceftazidime,
94 ceftolozane and ceftriaxone; two carbapenems: ertapenem and meropenem; a beta-lactamase
95 inhibitor: tazobactam; three fluoroquinolones: ciprofloxacin, moxifloxacin and levofloxacin; a
96 glycopeptide: daptomycin; a lipoglycopeptide: dalbavancin; and an oxazolidinone: linezolid.
97 In order to take into account the different physicochemical properties of the analytes, the
98 specific concentration ranges for each ATB and the priorities for application in relation to the

99 frequency of administration for treated patients, four differentiated extraction and
100 quantification methods were developed.

101

102 **2. Material and methods**

103 2.1 Chemicals

104 Amoxicillin (CAS-N. 61336-70-7 – purity 99.6%), piperacillin sodium salt (CAS-N.
105 59703-84-3 – purity 93.8%), cefazolin sodium salt (CAS-N. 27164-46-1– purity 97.7%),
106 cefotaxime sodium (CAS-N. 64485-93-4 – purity 100%), ceftazidime hydrate (CAS-N.
107 120618-65-7– purity 98.4%), ceftriaxone disodium salt hemi (heptahydrate) (CAS-N.
108 104376-79-6 – purity 99.6%), meropenem trihydrate (CAS-N. 119478-56-7 – purity 99.1%),
109 tazobactam sodium salt (CAS-N. 89785-84-2 – purity 92.5%), ciprofloxacin (CAS-N. 85721-
110 33-1 – purity 99.5%), moxifloxacin hydrochloride (CAS-N. 186826-86-8 – purity 100%),
111 levofloxacin (CAS-N. 100986-85-4 – purity 99.6%), dalbavancin (CAS-N. 171500-79-1 –
112 purity 92.6%) and linezolid (CAS-N. 165800-03-3 – purity 99.5%) were purchased from
113 Sigma-Aldrich (Saint Quentin Fallavier, France). $^2\text{H}_4$ -amoxicillin (CAS- 61336-70-7 – purity
114 96%), $^2\text{H}_5$ -piperacillin (CAS-59703-84-3 – purity 95%), $^{13}\text{C}_2$ ^{15}N -cefazolin sodium salt (CAS-
115 27164-46-1 – purity 98%), cefepime dihydrochloride monohydrate (CAS-N. 123171-59-5 –
116 purity 95%), $^2\text{H}_3$ -Cefepime sulfate (CAS- 123171-59-5 – purity 92%), $^2\text{H}_3$ -cefotaxime (CAS-
117 64485-93-4 – purity 95%), $^2\text{H}_5$ -ceftazidime (CAS- 120618-65-7 – purity 90%), $^2\text{H}_3$ -
118 ceftriaxone disodium salt (CAS-N.1132650-38-4 – purity 96%), $^2\text{H}_6$ -meropenem (CAS-
119 119478-56-7 – purity 97%), $^{15}\text{N}_3$ -tazobactam sodium salt (CAS- 89785-84-2 – purity 98%),
120 $^2\text{H}_8$ -ciprofloxacin (CAS-N.1130050-35-9 – purity 96%), daptomycin (CAS-N. 103060-53-3 –
121 purity 98%) and $^2\text{H}_3$ -linezolid (CAS-N.1127120-38-0 – purity 99.5%) were purchased from
122 LGC (Molsheim, France). Ceftolozane trifluoroacetate salt (CAS-N.1628046-32-1 – purity

123 98.7%), ¹⁵N₂ ²H₄-ceftolozane trifluoroacetate salt (CAS-N.1628046-32-1 – purity 98.1%),
124 ertapenem sodium salt (CAS-N.153773-82-1 – purity 92.7%), ²H₄-ertapenem sodium salt
125 (CAS-N.153773-82-1 – purity 90.6%), ²H₅-moxifloxacin trifluoroacetate salt (CAS-
126 N.1092356-43-8 – purity 99.7%), ¹³C ²H₃-levofloxacin (CAS-N.1261398-33-7 – purity 99%),
127 ²H₆-dalbavancin trifluoroacetate salt (CAS-N.1588823-86-2 – purity 95.3%) and ²H₅-
128 daptomycin trifluoroacetate salt (CAS-N.13565796-56-0 – purity 97.3%) were purchased
129 from ALSACHIM (Illkirch Graffenstaden, France). LC-MS-grade formic acid, acetonitrile
130 and methanol were supplied by Honeywell-Fluka (Illkirch, France), isopropanol was
131 purchased from VWR (Fontenay-sous-Bois, France) and water (UHPLC-MS quality) came
132 from Scharlau (Barcelona, Spain).

133 2.2 Biological samples

134 Blank drug-free serum was provided by the French blood establishment. Clinical samples
135 were collected as part of standard care from patients hospitalized in an ICU at Toulouse
136 University Hospital between 01/01/2020 and 30/06/2021. This retrospective cohort study was
137 approved by the Toulouse University Hospital review board (registration number: RnIPH
138 2021-78; CNIL number: 2206723 v 0).

139

140 2.3 Sample preparation

141 In order to overcome stability and weighing problems, two types of standard stock
142 solutions and quality control for each molecule were weighed by two different manipulators
143 and frozen at -80°C. One was used to prepare calibration ranges and the other for the
144 preparation of the controls. Before use, they were tested by several injections at the same
145 concentration in our analytical system and compared to give a difference of less than 5%.
146 Each marked antibiotic used as an internal standard was prepared under the same

147 solubilization conditions and stored at -80°C (Table S4). The antibiotics used had to reach
148 different areas of therapeutic concentrations and the range amplitudes were adapted for each
149 mixture from ATB1 to ATB4, while three levels of quality control were prepared to cover the
150 full amplitudes of the different calibration ranges (Table S5). In order to reproduce a
151 representative patient sample, we were required not to exceed a ratio of 1/10 or 1 volume of
152 overload of stock solution and 9 volumes of matrix (serum, plasma) [19]. Low quality control
153 was prepared at 3 or 4 times the low quantification limit, with average control between 30%
154 and 60% of calibration ranges, and high control between 75% and 85% of calibration ranges
155 (Table S5).

156 Given the characteristics of solubility, the low and high quantification limits to be
157 reached, the large number of these molecules, and our choice not to exceed a ratio of 1/10 or 1
158 volume of stock solution overload and 9 matrix volumes (serum and plasma), it was very
159 difficult to monitor all the molecules in a single extraction procedure. To take into account the
160 differences in physicochemical properties related in particular to differences in the solubilities
161 of analytes, the ranges of concentrations sought and sensitivity, the 16 ATBs and the beta-
162 lactamase inhibitor were extracted differentially into 4 groups ATB1, ATB2, ATB3 and
163 ATB4 following 4 specific protocols. Each group is organized according to the frequency of
164 requests for monitoring ATBs in the Toulouse ICU and the physicochemical properties of the
165 16 ATBs and beta lactamase inhibitor. The ATB1 pool was built with ciprofloxacin,
166 tazobactam, amoxicillin, ceftolozane, cefepime, ceftazidime, piperacillin and cefazolin; the
167 ATB2 group contained cefotaxime, daptomycin, levofloxacin, linezolid and moxifloxacin;
168 ATB3 was built with ertapenem and meropenem; ATB4 contained the last 2 molecules
169 ceftriaxone and dalbavancin.

170 Sample pretreatment of all biological samples was carried out in 2 steps. In the first step,
171 the proteins were precipitated. Depending on the pool of ATBs, precipitation was carried out

172 with a mixture of Methanol with 0.1% formic acid for ATB1 to ATB3. For ATB4
173 precipitation was carried out with a mixture of acetonitrile and water (50/50) (vol./vol.) with
174 0.1% formic acid. The samples were then centrifuged. In a second step dilution with water
175 with 0.1% formic acid was applied at differentiated ratios for ATB1 to ATB3. For ATB4 no
176 dilution is carried out. All internal standards were contained in the precipitation solvent added
177 in the first step.

178 For ATB1 to ATB3 protocols, the molecules were extracted and diluted with the same
179 solvent mixture composed of Methanol / 0.1% formic acid and water / 0.1% formic acid,
180 respectively. A simple deproteinization and adapted dilution of 1/39 for ATB1, 1/65 for
181 ATB2 and 1/91 for ATB3 were used. For the ATB 4 group, we opted for a much milder
182 precipitation.

183 ATB1 group: ATB1 internal standard stock solutions: $^{15}\text{N}_2$, $^2\text{H}_4$ -Ceftolozane, $^2\text{H}_4$ -
184 Amoxicillin, $^2\text{H}_8$ -Ciprofloxacin, $^{13}\text{C}_2$, $^{15}\text{N}_1$ -Cefazolin, $^2\text{H}_3$ -Cefepime, $^2\text{H}_5$ -Ceftazidime, $^2\text{H}_5$ -
185 Piperacillin and $^{13}\text{C}_2$, $^{15}\text{N}_1$ Tazobactam were mixed extemporaneously and diluted in methanol
186 with 0.1% (vol.) formic acid to 417, 208, 56, 278, 417, 556, 556 and 208 $\mu\text{g/L}$ respectively.
187 300 μL of the stock solution of ATB1 internal standards was mixed with 25 μL of serum
188 samples, standards or controls, vortexed 1 min and centrifuged at 10,000 g for 10 min. 200 μL
189 of supernatant was recovered, diluted with 400 μL of 0.1% (vol.) formic acid in water,
190 vortexed, transferred into autosampler vials and 2 μL was used for injection.

191 ATB2 group: ATB2 internal standard stock solutions: $^2\text{H}_3$ -Cefotaxime, ^{13}C , $^2\text{H}_3$ -
192 Levofloxacin, $^2\text{H}_3$ -Linezolid, $^{13}\text{C}_2$, $^2\text{H}_5$ -Moxifloxacin and $^2\text{H}_5$ -Daptomycin were mixed
193 extemporaneously and diluted in methanol with 0.1% (vol.) formic acid to 1000, 325, 250,
194 500 and 2500 $\mu\text{g/L}$ respectively. 240 μL of the stock solution of ATB2 internal standards
195 were mixed with 20 μL of serum samples, standards or controls, vortexed 1 min and

196 centrifuged at 10,000 g for 10 min. 150 μ L of supernatant were recovered, diluted with 600
197 μ L of 0.1% (vol.) formic acid in water, vortexed, transferred into autosampler vials and 2 μ L
198 were used for injection.

199 ATB3 group: ATB3 internal standard stock solutions: $^2\text{H}_6$ -Meropenem and $^2\text{H}_4$ -Ertapenem
200 were mixed extemporaneously and diluted in methanol with 0.1% (vol.) formic acid to 501
201 μ g/L. 240 μ L of the stock solution of ATB3 internal standards were mixed with 20 μ L of
202 serum samples, standards or controls, vortexed 1 min and centrifuged at 10,000 g for 10 min.
203 100 μ L of supernatant were recovered, diluted with 600 μ L of 0.1% formic acid (vol.) in
204 water, vortexed, transferred into autosampler vials and 2 μ L were used for injection.

205 ATB4 group: ATB4 internal standard stock solutions: $^2\text{H}_3$ -Ceftriaxone and $^2\text{H}_6$ -Dalbavancin
206 were mixed extemporaneously and diluted in 0.1% (vol.) formic acid in water / 0.1% (vol.)
207 formic acid in acetonitrile (50/50) (vol./vol.) to 1250 and 2400 μ g/L respectively. 250 μ L of
208 the stock solution of ATB4 internal standards were mixed with 20 μ L of serum samples,
209 standards or controls, vortexed 1 min and centrifuged at 10,000 g for 10 min. 100 μ L of
210 supernatant were recovered and transferred into autosampler vials and 2 μ L were used for
211 injection.

212

213 2.4 Instrumentation and analytical conditions

214 The analyses were carried out by HPLC consisting of 2 high-pressure binary pumps
215 (LC20ADXR) that can support up to 660 bars, an injector (SIL-20AXR) and a multi-channel
216 column oven (6-way) (CTO-20AC) from Shimadzu (Marne la Vallée, France), coupled with a
217 triple quadrupole mass spectrometer (Sciex API 4500) from ABSCIEX (Villebon-sur-Yvette,
218 France). Given the concentrations to be achieved, the physical and chemical characteristics,
219 and stability, 3 separation methods were developed to measure 16 molecules on a Kinetex

220 Polar-C18 100 Å 2.1 x 100 mm, 2.6 µm column from Phenomenex (Le Pecq, France) and a
221 method using a Synergi Polar-RP-80 Å 2 x 50 mm, 5 µm column from Phenomenex for the
222 two remaining molecules.

223 The mobile phase was composed of ultrapure water added with 0.1% formic acid (solvent A)
224 and acetonitrile added with 0.1% formic acid (solvent B), in order to obtain a pH of the
225 mobile phase below the lowest pKa of all the ATBs studied (Table S1). Solvents A and B
226 were delivered using dedicated gradients (Table S2), followed by an equilibration step for the
227 column before the next injection. The analysis times were respectively 6 min, 6.5 min, 5 min
228 and 4 min for each ATB pool. ATB1, ATB2 and ATB3 were separated on Kinetex 2.6 µm
229 Polar-C18 100 Å 100 x 2.1mm and ATB4 on Synergi 4 µm Polar-RP 80 Å 50 x 2 mm.

230 Given strong demand for analyses at our laboratory and thanks to our instrumentation
231 consisting of a column selector, our strategy was to minimize mobile phase changes for a
232 sequence overnight. Four different ATB groups (ATB1 to ATB4) were created to respond to
233 optimized separation, as well as to obtain concentration levels suitable for therapeutic
234 monitoring while having the chromatographic conditions best suited in terms of separation
235 and peak qualities (shape, asymmetry). Our choice was therefore based on 2 columns meeting
236 these criteria.

237 The LC/MS/MS conditions were set as follows: electrospray ionization (ESI) in positive
238 mode, spray voltage: 4500 V, gas pressure (nitrogen): curtain gas: 35 psi, ion source gas 1: 10
239 psi, ions source gas 2: 60 psi, collision gas: 10 psi, capillary temperature: 450 °C. All specific
240 parameters were set as described in Table S3.

241

242 2.5 Analytical method validation

243 Analytical validation was carried out according to the guidelines of the European
244 Medicines Agency (EMA) [20], indicating the set parameters that we had to test with the
245 acceptability criteria that we had to meet to prove the robustness of our method before it goes
246 into production. The selectivity of the method was studied by verifying its ability to quantify a
247 reference element (here ATBs) without interference from other reference elements or
248 endogenous compounds present in the matrix. It was tested using the appropriate blank matrix
249 from at least six different sources. The signal of any interfering peak in each blank matrix
250 should be less than 20.0% of the signal obtained at the LOQ for the reference element and
251 5.0% of the signal obtained for the internal standard.

252 The isotopic contribution was studied by verifying that no signal from the labeled internal
253 standard (IS) interfered with the ATB to be quantified and vice versa. A first blank matrix was
254 overloaded with the reference element to the concentration of the highest point of the range
255 without IS, a second with the IS without the reference element (P0), and finally a third with
256 the reference element to the LOQ without IS. The peak area obtained at the concentration of
257 the highest point of the range and at the retention time of the IS should be less than 5.0% of
258 the area of the peak of the IS obtained in the P0. The peak area observed in the P0 and the
259 retention time of the reference element should be less than 20.0% of the peak area of the
260 reference element obtained at the LOQ.

261 Carry-over was studied by verifying that no signal from the highly concentrated reference
262 element (or its EI) from a previous injection would contaminate a blank matrix. The same
263 blank was injected 3 times in a row after the highest standard of a calibration range. The peak
264 signal interfering in each blank matrix should be less than 20.0% of the signal obtained at the
265 LOQ for the reference element and 5.0% of the signal obtained for the IS without the
266 reference element (P0).

267 The influence of the anticoagulant agent was measured by comparing the quantification
268 of a serum sample with a heparinized plasma with the same concentration. Sample spikes
269 were performed extemporaneously in triplicate with the low internal quality control (IQC) and
270 high IQC level on our validation matrix (serum) as well as on lithium heparin plasma. The
271 accuracy criteria were calculated relative to the concentration determination of IQC serum in
272 triplicate. Accuracy and precision were to be less than 15%. A dilution study had to be
273 validated beforehand in order to be able to measure samples that were higher than the last
274 point of range. The dilution was done with the same matrix used for the overloads and the
275 volume of the final test take had to correspond exactly to that of the method. A sample was
276 prepared so that six extractions could be made. Its concentration is twice the highest point of
277 the range. During analysis, it was diluted to 1/10. Accuracy and precision were to be less than
278 15%.

279 In initial validation, the global measurement uncertainty (GMU) was calculated as a
280 percentage corresponding to the dispersion of values assigned to the measurements of our
281 IQC (low, medium, high) during the between-run precision study. It is estimated in initial
282 validation from eq. 1:

$$283 \quad GMU = \frac{(|\bar{X} - X_{ref}| + 2 * \sigma)}{X_{ref}} * 100$$

284 Where x_{ref} : target value, \bar{x} : average value, σ : standard deviation, GMU: Global Measurement
285 Uncertainty

286 Then, this value will be regularly updated by incorporating the results of external quality
287 control (EQC) obtained within the year by inter-laboratory controls.

288 A calibration curve was used to establish a relationship between the response of the
289 instrument and the concentration of the analyte. The calibration curves included six levels of

290 concentrations, a blank (matrix not overloaded with analyte or internal standard), and a P0
291 (matrix not overloaded with analyte but overloaded with internal standard). The P0 result was
292 not included in the calibration curve. During validation, this was redone every day for five
293 days and had to include the same points from one day to the next. The regression model was a
294 linear model with a 1/X weighting and the regression coefficient had to be >0.995. The values
295 of the first point and the last point of the range were not supposed to be excluded from the
296 calculation of the calibration curve equation. Recalculated concentrations were expected to
297 show a maximum accuracy of $\pm 20.0\%$ for the first point of range and $\pm 15\%$ for the other
298 points.

299 The limit of quantification (LOQ) was the smallest amount of ATB that could be
300 quantified with defined accuracy and precision. It corresponded to the first point of the
301 calibration range with accuracy and precision of less than 20%. The upper limit of
302 quantification (ULOQ), was the highest calibration standard.

303

304 The mathematical calculation compared with the measurements made showed us that the
305 estimate gave a good reflection of reality.

306 The within-run precision concerned appreciating intra-series variability. For all three
307 levels of controls, the samples were measured within a series with an extemporaneous
308 calibration.

309 Precision is expressed as the coefficient of variation (CV) of peak areas from eq. 4.

310
$$(CV)\% = \frac{\text{Standard Deviation}}{\text{Average}} * 100$$

311 Accuracy describes the closeness of the determined value obtained by the method to the
312 nominal concentration of the analyte (expressed in percentage or 'bias') from eq. 5.

313
$$Bias = \frac{m - v}{v} * 100$$

314

315 Between-run precision consisted of appreciating inter-series variability (CV, eq. 4) with
316 the notion of time or even analyst or equipment (estimation of robustness). Five sets of
317 analyses were performed over five different days, each with an extemporaneous calibration,
318 and each level of IQC was extracted in duplicate. The mean values for each range point were
319 calculated over all five series. Each point had to show a precision and accuracy within the
320 limit of 15% (20% for the LOQ). The precision of the slopes of the calibration curves should
321 not exceed 15% with a regression coefficient >0.995. In each series, the accuracy of the CIQs
322 had to be within the 15% limit. For each concentration level, the average accuracy obtained
323 had to be less than 15% as well as the precision calculated from the ten IQCs. All the
324 quantitative methods developed on the 2 LC columns were validated according to these
325 criteria.

326 2.6 Extraction efficiency and ionization loss (matrix effect in LC-MS-MS)

327 The matrix effect, determined as previously described for quantitative bioanalytical
328 methods [21], corresponds to the potential interference in the ionization of molecules in the
329 source and may increase or decrease ionization efficiency, thus altering the device's signal.
330 Three parameters were studied:

331 The matrix effect, by comparing the signal (X) obtained from a blank extracted and then
332 overloaded with the reference element, with the signal of a pure solution (100).

333
$$Matrix\ effect = \frac{X}{100}$$

334 The extraction efficiency, by comparing the signal (Y) obtained from a sample overloaded
335 with the reference element and then extracted, with the signal (X) obtained from a blank
336 extracted and then overloaded with the reference element.

337
$$\text{Extraction efficiency} = \frac{X}{Y}$$

338 The recovery rate, by comparing the signal (Y) obtained from a sample overloaded with the
339 reference element and then extracted, with the signal of a pure solution (100).

340
$$\text{Recovery Rate} = (\text{Extraction efficiency} * \text{Matrix effect}) = \frac{Y}{100} = \frac{Y}{X} * \frac{X}{100}$$

341 The three parameters were calculated in triplicate for each of the three levels of control and
342 for the internal standard. The extraction efficiency and matrix effect were expected to be less
343 than 15% (CV) for each level of concentration on the analyses. A matrix effect greater than
344 100% reflects an increase in ionization, while an effect less than 100% is the expression of an
345 ion suppression phenomenon.

346

347 2.7 Stability

348 The stability of an ATB in a biological matrix depends on its chemical properties, matrix
349 and storage conditions. Different stability conditions were tested to cover the different cases
350 of pre-dosing storage.

351 2.7.1 Stability in whole blood

352 Overloads were carried out extemporaneously in triplicate at the level of low IQC and
353 high IQC on our validation matrix (serum) as well as on whole blood taken on lithium heparin
354 and extract (H0). The overloaded samples were stored for 24 hours at 4°C and then extracted
355 in triplicate. The accuracy criteria were calculated based on the concentration of serum IQC
356 (H0). Accuracy and precision were to be less than 15%.

357 2.7.2 Short-term stability

358 Overloads were performed extemporaneously in triplicate at the level of low IQC and
359 high IQC on our validation matrix (serum) and extract (H0). Overloaded samples were stored
360 for 3 hours and 6 hours at room temperature, for 24 hours, 48 hours and 72 hours at 4°C and
361 then extracted in triplicate. The accuracy criteria were calculated based on the concentration
362 of serum IQC (H0). Accuracy and precision were to be less than 15%.

363 2.7.3 Long-term stability

364 Overloads were performed extemporaneously in triplicate at the level of low IQC and
365 high IQC on our validation matrix (serum) and extract (H0). Overloaded samples were stored
366 for 1 month, 3 months, 6 months and 12 months at -20°C and -80°C and then extracted in
367 triplicate. The accuracy criteria were calculated based on the concentration of serum IQC
368 (H0). Accuracy and precision were to be less than 15%.

369 2.7.4 Freeze-thaw stability

370 Overloads were performed extemporaneously in triplicate at the level of low IQC and
371 high IQC on our validation matrix (serum) and extract (H0). The overloaded samples were
372 frozen for at least 24 hours and then completely thawed at room temperature and finally
373 refrozen for a period of 24 hours (1st cycle). The second and third cycles were carried out in
374 the same way. At the end of the third cycle, they were extracted in triplicate. The accuracy
375 criteria were calculated based on the concentration of serum IQC (H0). Accuracy and
376 precision were to be less than 15%.

377 2.7.5 Injector stability

378 Overloads were performed extemporaneously in triplicate at the level of low IQC and
379 high IQC on our validation matrix (serum) and extract (H0). The extracts were stored for 24
380 hours at room temperature in the vials on the injector and then re-injected in triplicate. The

381 accuracy criteria were calculated based on the concentration of serum IQC (H0). Accuracy
382 and precision were to be less than 15%.

383 2.7.6 Analysis of clinical samples

384 These analytical methods are now validated and used daily to determine antibiotic
385 concentrations in plasma from patients at Toulouse University Hospital and nearby care
386 facilities. As a proof of concept, we thus analyzed antibiotic concentrations determined for
387 ICU patients. We chose to focus on cefepime for two reasons. On the one hand, cefepime has
388 a narrow therapeutic index and a well-documented link between drug exposure and
389 efficacy/neurotoxicity, and on the other hand it is by far the most requested antibiotic. We
390 selected patients treated by continuous cefepime perfusion. Blood was collected on
391 heparinized collection tubes when concentration reached the steady state. After arrival at the
392 laboratory, samples were centrifuged at 3,500 rpm for 10 min and stored at -20°C before
393 analysis. Measured concentrations were compared with guidelines published by two French
394 working groups (i.e. the French Society of Anesthesia and Intensive Care Medicine (SFAR)
395 and the French Society of Pharmacology and Therapeutics (SFPT)) [22].

396

397 **3. Results and discussion**

398 3.1 Chromatographic optimization

399 All the ATBs chosen for this study were in perfect harmony with the daily administration
400 practiced in our intensive care units. In this period of COVID 19, many of these molecules
401 had been administered and required rapid dosing to ensure optimum efficacy. The goal was to
402 obtain selective methods that were quick and easy to implement with all specific classes of
403 antibiotics (penicillins, cephalosporins, carbapenems, fluoroquinolones, etc.). All of the ATBs

404 studied contained specific chemical structures with very different physicochemical properties,
405 as shown in Fig. S1 (supplementary data). The values of the solubility product constant (K_{sp}),
406 as well as the octanol/water partition coefficients of each molecule, made it possible to choose
407 the solvents best suited to optimal dissolution (Table S1). The solubilization properties of the
408 analytes, described in Table S1, do not allow the use of a single generic method for the 17
409 molecules. We decided to develop four different quantification methods, grouping all 16
410 ATBs and 1 beta lactamase inhibitor into four pools, taking into account the frequency of
411 requests for analysis of each molecule and their physicochemical properties. The
412 compositions of the four ATB pools were mainly based on clinical priorities and practices:
413 thus, the ATB1 pool was built to include the most prescribed antibiotic drugs in our hospital
414 (ciprofloxacin, tazobactam, amoxicillin, ceftolozane, cefepime, ceftazidime, piperacillin and
415 cefazolin). Overall, this pool represents more than 80% of antibiotic dosing demands from
416 clinical departments. While most of the remaining antibiotics were then included in the ATB2
417 pool, some of them were isolated due to logistical and analytical constraints. The ATB3 pool
418 included carbapenems (meropenem and ertapenem), which display poor stability in biological
419 fluids [23] and therefore must be prioritized for analysis. The ATB4 pool, including
420 daptomycin and ceftriaxone, had to be developed with different chromatographic conditions
421 due to the poor performances of the initial conditions in terms of separation and peak aspect.

422 To mix a maximum of compounds, the concentrations to be reached for stock solutions
423 had to be high. We therefore adapted the concentrations in relation to the high limits of
424 linearity to be reached and established the antibiotic pools according to these different
425 constraints. Typically, the upper limits of quantification for piperacillin and cefazolin were set
426 at 40 mg.L^{-1} and serum samples were diluted to quantitate the analytes in the linear domain
427 (Table S4, supplementary data). During method development, we controlled the absence of

428 interference peaks with the ATBs studied in order to promote selectivity. Moreover, all ATBs
429 had an isotopic contribution <20%.

430 Due to the complex structure of all the antibiotics studied, with negative or very high log
431 P, we favored columns capable of performing hydrophilic and hydrophobic interactions and
432 with the ability to work in 100% aqueous phase with maximum latitude in optimizing the
433 elution gradient. The Kinetex Polar C18 column gave the best separation and peak shape
434 results for ATB1 to ATB3 pools. However, it was not possible to obtain correct peaks with
435 this column for ceftriaxone, probably due to the dioxotriazide group with this polar graft [19].
436 Better results were observed on a Synergi Polar-RP column, however, with an ether bonded
437 phenyl phase and with a hydrophilic endcapping giving a satisfactory peak shape quality. We
438 therefore developed the ceftriaxone method with this column, associated with dalbavancin,
439 which does not pose any particular problems. The chromatographic separations of each pool
440 was optimized using the iterative software Osiris[®]. The latter provided us with good
441 predictions with the best separations ($RS > 1.5$ and $2 < K < 10$), while having a short analysis
442 time. However, it was not always possible for certain molecules (ceftolozane, cefepime,
443 amoxicillin, tazobactam) (Table S1, supplementary data) with low values of octanol/water
444 partition coefficients, to obtain a correct separation with resolution >1.5 (Fig S2,
445 supplementary data), in less than 3 minutes of elution time, but in no way interfering with
446 quantification. On the other hand, the ATB2 to ATB4 pools gave good separations (Fig S2,
447 supplementary data). Similar problems of co-elution and separation with low resolution were
448 obtained for these same classes of ATBs and beta lactamase inhibitor on semi-porous
449 columns, with particle sizes equivalent to that used in our study, using buffered eluents.
450 Amoxicillin, tazobactam and cefepime were also the polar compounds that show the least
451 selectivity [24].

452 For the point linked to the MS transitions, ATBs and beta lactamase inhibitor are small
453 molecules with molecular masses of less than 700 g/mol, except for daptomycin (1620.7
454 g/mol) and dalbavancin (1816.7 g/mol). Only one transition was followed, due to the high
455 dose given to patients in ICU and the comparison made between the analytical result and the
456 dose administered. For low masses, the mass spectrometer followed simply charged
457 transitions protonated into the ESI source, and for the high masses, it followed doubly-
458 charged protonated transitions. MRM transition tracking is listed on Table S3 (supplementary
459 data).

460 3.2 Sample preparation

461 The results of the extraction yields, matrix effects and recovery rates for the sets of
462 molecules are presented in Table 1. The very wide range of hydrophobicity in the ATBs,
463 determined by log P at between -6.17 and 3.8 (Table S1) for the ATBs studied, led us to
464 develop a method of extraction as simple as possible in terms of the number of steps and types
465 of solvents used. The difficulty was to best adapt the matrix volumes collected for each group,
466 on the one hand, and the composition and volume of the dilution reagent of the supernatant
467 after precipitation, on the other hand, in order to respect the low and high quantification
468 limits. As the amplitudes of the ranges were sometimes very large, we had to be able to
469 correctly quantify very low levels of antibiotics, but also very high levels without risking
470 saturating the detector. The matrix proteins of the first three ATB groups were precipitated
471 with a mixture of methanol and formic acid at 0.1% and then the supernatant was diluted with
472 the same solvent. For the ATB 4 group, we opted for a much gentler precipitation to avoid
473 problems of a low recovery rate due to the adsorption of dalbavancin on proteins [24].

474 Ceftolozane, ceftazidime, levofloxacin and cefepime had ionization gains of 225, 174,
475 158 and 154%, respectively. The matrix effect is important for these 4 ATBs, which are polar,

476 with log P values of -6.17, -1.60, -0.40, -0.37, respectively (Table S1, supplementary data).
477 Lefeuvre et al. found values of >100% for ceftazidime, levofloxacin and cefepime, but of
478 <110%, using high resolution mass spectrometry and XIC mode for quantification, but no
479 information was given on the extraction efficiency and extraction yield [19]. Rigo-Bonnin et
480 al. developed a method showing a lower matrix effect (108%) for the quantification of
481 ceftolozane, in tandem mass spectrometry, but separating only this ATB and tazobactam [25].
482 For levofloxacin, a lower matrix effect was determined by Fang et al. in mouse plasma (94%),
483 with strong precipitation conditions, unusable with all ATBs [26]. The percentages of other
484 ATBs ranged from 90% to 130%, except for dalbavancin, which has a very significant
485 ionization loss at 20%.

486 Extraction yields were >65% except for ceftolozane (41%) and levofloxacin (65%). The
487 extraction yield was lower for ceftolozane than for the other ATBs, but for the method
488 developed in this study, it was higher than that found by Putnam et al. (20%), which used a
489 precipitation of serum proteins using 90% acetonitrile [27]. In the method developed by
490 Putnam et al. the very low extraction yield is observed as a very low matrix effect (96%).
491 Finally, recovery for all ATBs had always been between 64% and 138%, except for
492 dalbavancin at 20%, despite mild deproteinization, with precision on the three levels of
493 concentration (7%) <15%. However, the internal standards have always corrected this
494 characteristic phenomenon of the sometimes consequential LC-MS-MS. Precision measured
495 for extraction efficiency, matrix effect and recovery is <15% for all ATBs, except for
496 ceftolozane, cefepim, levofloxacin, linezolid, and moxifloxacin linked to the higher values for
497 matrix effects (Table 1).

498 Table 1: Extraction yields, matrix effects and recovery obtained with the 16 ATBs and 1
499 beta-lactamase with the differentiated quantification methods ATB1 to ATB4.

Antibiotics	Extraction efficiency		Matrix effect		Recovery	
	Rate	Precision	Rate	Precision	Rate	Precision
	%	%	%	%	%	%
ATB 1						
Ceftolozane	41	11	225	22	92	19
Cefepime	89	13	154	16	137	19
Amoxicillin	92	8	105	8	96	8
Tazobactam	91	7	91	7	92	6
Ceftazidime	70	6	174	9	122	8
Ciprofloxacin	93	15	122	15	110	12
Cefazolin	88	5	105	5	92	5
Piperacillin	86	15	106	15	90	7
ATB 2						
Cefotaxime	68	9	90	3	64	10
Daptomycin	88	10	94	7	82	11
Levofloxacin	65	20	158	25	107	40
Linezolid	80	9	99	13	80	18
Moxifloxacin	76	15	130	20	100	27
ATB 3						
Ertapenem	68	4	127	7	87	3
Meropenem	81	5	130	2	105	5
ATB 4						
Ceftriaxone	102	3	120	4	122	5
Dalbavancin	99	4	20	7	20	7

500

501 3.3 Method validation

502 3.3.1 Carry-over

503 Ertapenem and daptomycin had a significant memory effect of 169% and 30% requiring,
504 after each CQI and depending on the patient, 3 injections of 10 µL of phase rinse
505 (H2O/methanol/acetonitrile/isopropanol 25-25-25-25). The other ATBs had no memory effect
506 phenomena.

507 3.3.2 Precision and accuracy

508 The analytical validation results are summarized in Table 2. For all analyses, bias is
509 <15%. Ciprofloxacin has the largest bias observed in inter-day between-run precision (-
510 11.7%) and piperacillin has the largest bias observed for intra-day within-run precision (-
511 13.3%). For these same molecules, Lefeuvre and al. [19] observed lower biases, but
512 ceftriaxone was analyzed with a high intra-day bias (14.1%), while for our study we observe a
513 low value of bias (<2.8%) for this ATB. The relative standard deviation was acceptable for all
514 ATBs. The highest values were observed for piperacillin (11.8%) for between-run precision
515 and for daptomycin (7.8%) for within-run precision, which is consistent with the literature.
516 The specifications to be met were reached in terms of within-run precision, between-run
517 precision and accuracy.

518

519 Table 2: Accuracy, Within-run and between-run precision results.

520

Accuracy and within-run precision (Intra-day)			Accuracy and between-run precision (Inter-day)		
Low QC *	Medium QC *	High QC *	Low QC*	Medium QC *	High QC *

	Molecules	Bias %	CV %	Bias %	CV %	Bias %	CV %	Bias %	CV %	Bias %	CV %	Bias %	CV %
ATB1	Ciprofloxacin	-11.1	7.1	4.8	5.1	11.0	3.3	-11.7	5.5	-2.2	7.6	3.3	7.1
	Tazobactam	-6.3	4.9	-5.3	6.7	-0.5	7.1	-10.0	4.9	-9.9	5.1	-6.8	5.7
	Amoxicillin	-8.8	4.8	-10.0	3.8	-6.7	2.7	-3.3	6.8	-0.6	9.3	-1.4	10.0
	Ceftolozane	-8.2	5.5	-11.0	4.0	-7.8	3.1	-8.6	5.9	-3.7	10.3	-0.6	8.7
	Cefepime	-11.8	6.0	-8.4	3.7	-1.6	4.8	-7.6	6.5	-3.9	8.8	-2.3	7.9
	Ceftazidime	-1.8	2.6	-3.5	5.1	-1.1	6.8	0.2	6.7	-3.9	5.7	-2.0	8.9
	Piperacillin	-11.3	2.4	-13.3	2.6	-12.5	2.7	-3.5	8.6	-2.7	9.1	0.9	11.8
	Cefazolin	-9.0	3.8	-4.2	6.2	-3.9	4.7	-9.9	3.8	-4.3	4.0	-5.7	4.4
ATB2	Levofloxacin	4.7	3.6	7.4	2.8	7.4	2.8	-1.5	5.2	2.0	4.8	2.9	4.1
	Ofloxacin	-11.6	1.9	-2.1	7.6	-13.0	0.9	-8.6	8.0	-9.5	5.8	-10.4	4.8
	Moxifloxacin	-11.6	4.2	-5.7	1.6	-5.1	1.1	-7.5	7.7	-2.4	5.6	-2.1	7.0
	Linezolid	-10.5	3.9	-7.4	2.0	-4.8	2.4	-10.1	5.5	11.4	2.8	-8.7	3.0
	Cefotaxime	-6.8	5.4	0.5	5.5	-3.8	2.5	-7.9	4.0	-0.4	4.1	-6.4	3.5
	Daptomycin	8.6	7.8	3.4	6.1	6.1	4.9	6.3	10.4	0.1	5.5	3.1	6.6
ATB3	Meropenem	4.3	2.2	3.1	3.7	2.3	2.1	4.2	2.7	4.3	3.2	4.0	2.8
	Ertapenem	-5.7	5.8	-6.8	2.8	-7.7	1.6	-2.9	6.3	-4.9	5.7	-5.3	7.1
ATB4	Ceftriaxone	-0.7	1.4	0.8	1.6	-2.8	1.7	2.5	4.1	2.7	4.4	-1.9	4.1
	Dalbavancin	1.9	2.5	5.6	0.6	3.7	1.7	3.6	5.2	6.0	1.7	3.2	1.9

521 * Low QC, Medium QC and High QC are given for each ATBs in Table S5.

522

523 3.3.3 Limit of quantification (LOQ)

524 The LOQs were studied and all corresponded to the first point of the ATB calibrations,
525 for which accuracy and precision were measured with a standard deviation of <20% (Table 3).
526 These LOQs were adapted to each low value of the ATBs sought in patients' serums. They
527 ranged from 0.05 mg.L⁻¹ for ciprofloxacin to 1 mg.L⁻¹ for ceftriaxone and dalbavancin. Higher
528 values were expected in serums from ICU patients with pathophysiological features, which
529 are integrated into the final dosage by the systematic dilution applied to all samples (1/10), in
530 order to meet the guidelines of specific ICU therapeutic practices [22]. For most literature
531 studies, LOQs are set at 0.5 mg.L⁻¹ or 1 mg.L⁻¹ for ATBs. [16,19].

532 The ULOQs were also studied and all corresponded to the last point of the ATB
533 calibrations, for which accuracy and precision were measured with a standard deviation of

534 <15% (Table 3). They ranged from 2.5 mg.L⁻¹ for ciprofloxacin to 100 mg.L⁻¹ for
535 Cefotaxime, Daptomycin, Meropenem and Ertapenem.

536 By differentiating the quantification into 4 pools, it is possible to monitor of 16 ATBs and
537 the beta lactamase inhibitor using an MS/MS detection based on a single transition, with
538 within-run precision, between-run precision and accuracy close to those using an HRMS
539 detection [19]. If dilution of the serum is an additional step, it is however mandatory, because
540 of the high doses administered. The robustness of this step is correctly evaluated for High QC
541 values, with CVs<15% (Table 2). The validation of these 4 quantification methods shows
542 performances equivalent to those obtained for the monitoring of 15 ATBs in plasma
543 developed by Lefeuvre et al. [19], in which all compounds were detected by HRMS and
544 quantified with a linear regression between 0.5 and 32 mg/L ($R^2 > 0.998$). The 4 pools ATB1,
545 ATB2, ATB3 and ATB4 show a quantification with a linear regression adapted to each ATB
546 as presented in Table 3. The LOQ and ULOQ are modulated according to the prescriptions
547 and target values [12]. In all the ranges sought, the linearity of the calibration is validated
548 ($R^2>0.998$) (Table 3). The slopes of the calibration lines are in the same order of magnitude as
549 those obtained with an HRMS detection mode, although comparisons are limited, since the
550 internal standards used are not all identical or used at the same concentration [19]. The
551 advances obtained with this new methodology are the modulation of the calibration ranges
552 and the validation of the robustness of the 4 extraction and quantification pools.

553 Table 3: Measured LOQ, ULOQ and linear regression parameters for the 16 ATBs and 1
554 beta-lactamase inhibitor.

555

Measured LOQ (mg.L ⁻¹)	Measured ULOQ (mg.L ⁻¹)	Linear regression
--	---	-------------------

Molecules		Slope	Intercept	R ²
ATB1	Ciprofloxacin	0.05	2.5	y = 1.418x - 0.0853 (R ² = 0.999)
	Tazobactam	0.10	5	y = 1.044x + 0.0935 (R ² = 0.999)
	Amoxicillin	0.20	10	y = 0.540x + 0.0239 (R ² = 0.999)
	Ceftolozane	0.40	20	y = 0.202x - 0.0009 (R ² = 0.998)
	Cefepime	0.40	20	y = 0.189x + 0.0089 (R ² = 0.999)
	Ceftazidime	0.40	20	y = 0.146x - 0.0010 (R ² = 0.999)
	Piperacillin	0.80	40	y = 0.097x - 0.0019 (R ² = 0.999)
	Cefazolin	0.80	40	y = 0.130x + 0.0168 (R ² = 0.999)
ATB2	Levofloxacin	0.13	25	y = 0.290x + 0.0014 (R ² = 0.999)
	Moxifloxacin	0.13	25	y = 0.267x - 0.0039 (R ² = 0.999)
	Linezolid	0.13	25	y = 0.235x + 0.0150 (R ² = 0.999)
	Cefotaxime	0.50	100	y = 0.076x + 0.0109 (R ² = 0.998)
	Daptomycin	0.50	100	y = 0.144x + 0.0421 (R ² = 0.999)
ATB3	Meropenem	0.20	100	y = 0.139x + 0.0023 (R ² = 0.999)
	Ertapenem	0.20	100	y = 0.129x + 0.0041 (R ² = 0.999)
ATB4	Ceftriaxone	1.00	50	y = 0.064x + 0.0109 (R ² = 0.999)
	Dalbavancin	1.00	50	y = 0.044x + 0.0033 (R ² = 0.999)

556

557 3.3.4 Anticoagulant, dilution, global measurement uncertainty and stability.

558 This accuracy and precision study first showed that heparinized tubes can be accepted for
559 monitoring 16 ATBs and the beta lactamase inhibitor, bias and CV being <15% for the high
560 and low QC values (Table 4). Global measurement uncertainty was small <30%, without
561 reference in the literature.

562 Table 4: Lithium Heparin Validation Results, Dilution, and Global Uncertainty Calculation.
563

Molecules	Comparison of doped serum versus doped heparinized plasma				Validation of dilution		Global measurement uncertainty			
	Low QC		High QC		Dilution (1/10 ULOQ)		Low QC	Medium QC	High QC	
	Bias %	CV %	Bias %	CV %	Bias %	CV %	CV %	CV %	CV %	
ATB1	Ciprofloxacin	0.5	4.9	-1.5	1.9	12.6	2.0	21.3	17.1	17.9
	Tazobactam	-13.4	1.2	-9.5	4.7	1.9	7.0	18.9	19.0	17.3
	Amoxicillin	-0.3	7.6	-7.7	0.3	6.0	4.1	16.4	19.1	20.0
	Ceftolozane	-8.3	9.6	-6.1	0.3	4.7	5.9	19.4	23.6	17.9
	Cefepime	-11.5	4.0	-9.3	0.3	9.8	3.2	19.6	20.8	17.7
	Ceftazidime	-9.3	4.8	-5.2	3.2	12.9	3.1	13.5	14.8	19.5
	Piperacillin	-8.9	3.7	-3.4	0.3	6.1	2.6	20.2	20.3	24.7
	Cefazolin	-6.1	5.0	-6.9	4.4	8.7	2.8	16.7	12.1	13.9

ATB2	Levofloxacin	10.4	0.8	4.4	1.8	12.0	2.5	11.8	11.7	11.4
	Ofloxacin	1.1	0.8	-0.6	1.2	2.5	-0.6	23.2	20.0	19.0
	Moxifloxacin	-5.5	6.2	1.0	3.0	8.9	5.4	21.8	13.4	16.9
	Linezolid	-1.6	3.7	4.9	3.2	-7.8	3.3	20.0	16.3	14.2
	Cefotaxime	-8.1	3.2	3.4	2.3	8.7	3.6	15.2	8.6	13.0
	Daptomycin	-6.8	2.8	0.0	4.8	7.8	4.8	28.3	11.2	16.7
ATB3	Meropenem	0.6	2.2	1.9	0.2	1.3	5.3	9.7	10.9	9.9
	Ertapenem	7.9	1.1	1.3	1.9	3.5	4.6	15.0	15.8	18.8
ATB4	Ceftriaxone	3.1	0.9	-3.9	2.0	11.5	1.9	11.0	11.7	10.0
	Dalbavancin	4.8	7.9	-2.6	3.9	3.4	5.4	14.4	9.5	7.1

564

565 A study of the stability of all the antibiotics was carried out and is shown in Table 5.

566 Table 5: Short-term, long-term, autosampler and freeze-thaw stability results in serum and
567 whole blood, for the 4 differentiated extraction and quantification ATBs pools (ATB1 to
568 ATB4).

569

		Stability						
		Short-term		Long-term		Autosampler	Freeze-thaw	
		Room temperature	4°C	-20°C	-80°C	Room temperature	Room temperature / -80°C	
Spiked sample	Serum (hours)	Whole blood (hours)	Serum (days)	Serum (months)	Serum (months)	Extract (hours)	Serum (number)	
ATB 1	Ciprofloxacin	6	24	2	12	12	24	3
	Tazobactam	6	24	2	1	6	24	3
	Amoxicillin	6	24	5	1	6	24	3
	Ceftolozane	6	24	2	6	6	24	3
	Cefepime	3	24	2	1	6	24	3
	Ceftazidime	6	24	5	3	6	24	3
	Piperacillin	3	24	1	1	6	24	3
	Cefazolin	6	24	2	12	12	24	3
ATB 2	Levofloxacin	6	24	3	6 (test 12 months in progress)	6 (test 12 months in progress)	24	3
	Ofloxacin	6	24	3	6 (test 12 months in progress)	6 (test 12 months in progress)	24	3

	Moxifloxacin	6	24	3	6 (test 12 months in progress)	6 (test 12 months in progress)	24	3
	Linezolid	6	24	3	6 (test 12 months in progress)	6 (test 12 months in progress)	24	3
	Cefotaxime	6	QC Low not stable 24 for QC High	3	3	6 (test 12 months in progress)	24	3
	Daptomycin	6	24	3	6 (test 12 months in progress)	6 (test 12 months in progress)	24	3
ATB 3	Meropenem	6	not stable 24 H	1	not stable 1 M	6 (test 12 months in progress)	not stable 24 H	3
	Ertapenem	6	24	2	not stable 1 M	6 (test 12 months in progress)	not stable 24 H	3
ATB 4	Ceftriaxone	6	24	3	3 (test 6 and 12 months in progress)	3 (test 6 and 12 months in progress)	24	3
	Dalbavancin	6	24	3	3 (test 6 and 12 months in progress)	3 (test 6 and 12 months in progress)	24	3


570

571 The stability study shows that some ATBs and the beta lactamase inhibitor cannot be
572 quantified after 3 hours of stacking of samples at room temperature, such as cefepime and
573 piperacillin. For piperacillin, Mortensen et al. showed 6-hr stability at room temperature for
574 human plasma treated with EDTA [23]. The stability of these sensitive ATBs and the beta
575 lactamase inhibitor is always low (1 month), after freezing at -20°C, but above 6 months, if
576 the stacking is carried out at -80°C. The remaining stability durations are consistent with
577 those given in the literature [15,17–19]. On total blood, we find 2 ATBs that are not stable
578 over a duration of 24 hrs at 4°C: cefotaxime for the low QC and meropenem.

579 3.4 Application to clinical samples.

580 Cefepime was the most prescribed antibiotic for patients in ICUs, representing around 80% of
581 total antibiotic requests. As such, we chose to focus our study on this antibiotic, found in the
582 ATB1 pool. Cefepime levels were measured at steady state in patients treated with continuous
583 infusion of the antibiotic. Daily dose of cefepime varied between 2000 mg and 8000 mg

584 according to identified bacteria, infection site and the patients renal function. In total, 50
585 cefepime plasma concentrations were measured with the ATB1 quantification method. Their
586 distribution is represented in Fig.1.



587
588 Fig.1. Distribution of measured concentrations of cefepime in ICU patients (in green:
589 recommended concentration window from the SFAR-SFPT consensus; in red: median value
590 of measured concentrations).

591 The recommended concentrations window is 5 to 35 mg.L⁻¹ for unknown bacterial infections
592 in the case of a continuous infusion [22]. While no underexposure was detected, indicating
593 sufficient antibiotic efficacy, around 75% of measured concentrations were above the
594 recommended concentrations range. This point is particularly important as our results
595 highlight an increased risk of toxicity, and more precisely neurotoxicity [28]. This can be
596 easily explained by the pharmacokinetic behavior of the drugs administered to ICU patients.
597 Indeed, this population tends to have damaged renal functions, leading to a decrease in
598 cefepime elimination. Monitoring cefepime concentrations is thus essential to optimize patient
599 care.

600 4. Conclusion

601 In this work, we developed the analysis of 16 ATBs from different classes and one beta-
602 lactamase inhibitor, using a simple protein precipitation method and grouping the ATBs into 4
603 pools, to rapidly respond to requests for assays, prioritizing the most requested ATBs in the
604 ATB1 pool. For the 4 quantification methods evaluated in this study, sampling was limited to
605 20 μ L of serum, whereas conventional methods use 50 μ L to 100 μ L. These quantification

606 methods allow wide spectrum monitoring from the most commonly administered ATBs
607 (ATB1) to the last molecules grouped in ATB3 and ATB4. The separation methods use the
608 same solvents, to simplify the column change needed to quantify all 17 molecules, while LC-
609 MS-MS technology enabled us to optimize the analysis time as much as possible (extraction,
610 injection, integration and validation). Performance was evaluated and validated according to
611 the EMEA guidelines and demonstrated a robust, easy-to-implement method. As we had also
612 seen widely varying stabilities for the different ATBs, we improved sample flow and storage
613 upon receipt in the laboratory. In addition, as this technology allowed for extremely low
614 volume (a few μL), the amount of sample was reduced. This was particularly interesting for
615 pediatric samples and pharmacokinetic studies.

616 Developing this new method was motivated by an optimization of patient care. In this context,
617 therapeutic drug monitoring for antibiotics has to be performed 7 days a week to adjust drug
618 dosage as quickly as possible. Our goal has been achieved, since our strategy was to extend
619 this methodology to transfer our old methods and optimize them as best we could with our
620 more efficient instruments.

621 **Acknowledgements**

622 We especially thank Mr. Clément De Almeida for his technical support for this work.

623

624 **References**

- 625 [1] A. Charmillon, E. Novy, N. Agrinier, M. Leone, A. Kimmoun, B. Levy, B. Demoré, J.
626 Dellamonica, C. Pulcini, The ANTIBIOPERF study: a nationwide cross-sectional
627 survey about practices for β -lactam administration and therapeutic drug monitoring
628 among critically ill patients in France, *Clin. Microbiol. Infect.* 22 (2016) 625–631.

- 629 <https://doi.org/10.1016/j.cmi.2016.04.019>.
- 630 [2] J.A. Roberts, M.H. Abdul-Aziz, J. Lipman, J.W. Mouton, A.A. Vinks, T.W. Felton,
631 W.W. Hope, A. Farkas, M.N. Neely, J.J. Schentag, G. Drusano, O.R. Frey, U.
632 Theuretzbacher, J.L. Kuti, Individualised antibiotic dosing for patients who are
633 critically ill: Challenges and potential solutions, *Lancet Infect. Dis.* 14 (2014) 498–509.
634 [https://doi.org/10.1016/S1473-3099\(14\)70036-2](https://doi.org/10.1016/S1473-3099(14)70036-2).
- 635 [3] A. Borsuk-De Moor, E. Rypulak, B. Potręć, P. Piwowarczyk, M. Borys, J. Sysiak, D.
636 Onichimowski, G. Raszewski, M. Czuczwar, P. Wiczling, Population pharmacokinetics
637 of high-dose tigecycline in patients with sepsis or septic shock, *Antimicrob. Agents*
638 *Chemother.* 62 (2018) 1–11. <https://doi.org/10.1128/AAC.02273-17>.
- 639 [4] J.L. Vincent, M. Bassetti, B. François, G. Karam, J. Chastre, A. Torres, J.A. Roberts,
640 F.S. Taccone, J. Rello, T. Calandra, D. De Backer, T. Welte, M. Antonelli, Advances in
641 antibiotic therapy in the critically ill, *Crit. Care.* 20 (2016) 1–13.
642 <https://doi.org/10.1186/s13054-016-1285-6>.
- 643 [5] F. Wolff, G. Deprez, L. Seyler, F. Taccone, M. Hites, B. Gulbis, J.L. Vincent, F.
644 Jacobs, F. Cotton, Rapid quantification of six β -lactams to optimize dosage regimens in
645 severely septic patients, *Talanta.* 103 (2013) 153–160.
646 <https://doi.org/10.1016/j.talanta.2012.10.024>.
- 647 [6] M.C. Verdier, O. Tribut, P. Tattevin, Y. Le Tulzo, C. Michelet, D. Bentué-Ferrer,
648 Simultaneous determination of 12 β -lactam antibiotics in human plasma by high-
649 performance liquid chromatography with UV detection: Application to therapeutic drug
650 monitoring, *Antimicrob. Agents Chemother.* 55 (2011) 4873–4879.
651 <https://doi.org/10.1128/AAC.00533-11>.

- 652 [7] S.E. Briscoe, B.C. McWhinney, J. Lipman, J.A. Roberts, J.P.J. Ungerer, A method for
653 determining the free (unbound) concentration of ten beta-lactam antibiotics in human
654 plasma using high performance liquid chromatography with ultraviolet detection, *J.*
655 *Chromatogr. B Anal. Technol. Biomed. Life Sci.* 907 (2012) 178–184.
656 <https://doi.org/10.1016/j.jchromb.2012.09.016>.
- 657 [8] A. Kratzer, S. Schießer, P. Matzneller, B. Wulkersdorfer, M. Zeitlinger, J.
658 Schlossmann, F. Kees, C. Dorn, Determination of total and free ceftolozane and
659 tazobactam in human plasma and interstitial fluid by HPLC-UV, *J. Pharm. Biomed.*
660 *Anal.* 163 (2019) 34–38. <https://doi.org/10.1016/j.jpba.2018.09.044>.
- 661 [9] F.U. Khan, F. Nasir, Z. Iqbal, I. Khan, N. Shahbaz, M. Hassan, F. Ullah, Simultaneous
662 determination of moxifloxacin and ofloxacin in physiological fluids using high
663 performance liquid chromatography with ultraviolet detection, *J. Chromatogr. B Anal.*
664 *Technol. Biomed. Life Sci.* 1017–1018 (2016) 120–128.
665 <https://doi.org/10.1016/j.jchromb.2016.03.002>.
- 666 [10] Y. Zheng, Z. Wang, G. Lui, D. Hirt, J.M. Treluyer, S. Benaboud, R. Aboura, I. Gana,
667 Simultaneous quantification of levofloxacin, pefloxacin, ciprofloxacin and
668 moxifloxacin in microvolumes of human plasma using high-performance liquid
669 chromatography with ultraviolet detection, *Biomed. Chromatogr.* 33 (2019).
670 <https://doi.org/10.1002/bmc.4506>.
- 671 [11] D. Cattaneo, S. Baldelli, F. Conti, V. Cozzi, E. Clementi, Determination of linezolid in
672 human plasma by high-performance liquid chromatography with ultraviolet detection,
673 *Ther. Drug Monit.* 32 (2010) 520–524.
674 <https://doi.org/10.1097/FTD.0b013e3181d5eeee>.
- 675 [12] J.A. Roberts, S.K. Paul, M. Akova, M. Bassetti, J.J. De Waele, G. Dimopoulos, K.M.

676 Kaukonen, D. Koulenti, C. Martin, P. Montravers, J. Rello, A. Rhodes, T. Starr, S.C.
677 Wallis, J. Lipman, A. Margarit Ribas, L. De Crop, H. Spapen, J. Wauters, T.
678 Dugernier, P. Jorens, I. Dapper, D. De Backer, F.S. Taccone, L. Ruano, E. Afonso, F.
679 Alvarez-Lerma, M.P. Gracia-Arnillas, F. Fernández, N. Feijoo, N. Bardolet, A. Rovira,
680 P. Garro, D. Colon, C. Castillo, J. Fernado, M.J. Lopez, J.L. Fernandez, A.M. Arribas,
681 J.L. Teja, E. Ots, J. Carlos Montejo, M. Catalan, I. Prieto, G. Gonzalo, B. Galvan, M.A.
682 Blasco, E. Meyer, F. Del Nogal, L. Vidaur, R. Sebastian, P.M. Garde, M.D.M. Martin
683 Velasco, R. Zaragoza Crespo, M. Esperatti, A. Torres, O. Baldesi, H. Dupont, Y.
684 Mahjoub, S. Lasocki, J.M. Constantin, J.F. Payen, J. Albanese, Y. Malledant, J.
685 Pottecher, J.Y. Lefrant, S. Jaber, O. Joannes-Boyau, C. Orban, M. Ostermann, C.
686 McKenzie, W. Berry, J. Smith, K. Lei, F. Rubulotta, A. Gordon, S. Brett, M. Stotz, M.
687 Templeton, C. Ebm, C. Moran, V. Pettilä, A. Xristodoulou, V. Theodorou, G.
688 Kouliatsis, E. Sertaridou, G. Anthopoulos, G. Choutas, T. Rantis, S. Karatzas, M.
689 Balla, M. Papanikolaou, P. Myrianthefs, A. Gavala, G. Fildisis, A. Koutsoukou, M.
690 Kyriakopoulou, K. Petrochilou, M. Kompoti, M. Michalia, F.M. Clouva-Molyvdas, G.
691 Gkiokas, F. Nikolakopoulos, V. Psychogiou, P. Malliotakis, E. Akoumianaki, E.
692 Lilitsis, V. Koulouras, G. Nakos, M. Kalogirou, A. Komnos, T. Zafeiridis, C.
693 Chaintoutis, K. Arvaniti, D. Matamis, C. Kydona, N. Gritsi-Gerogianni, T.
694 Giasnetsova, M. Giannakou, I. Soultati, I. Chytas, E. Antoniadou, E. Antipa, D.
695 Lathyris, T. Koukoubani, T. Paraforou, K. Spiropoulou, V. Bekos, A. Spring, T.
696 Kalatzi, H. Nikolaou, M. Laskou, I. Strouvalis, S. Aloizos, S. Kapogiannis, O.
697 Soldatou, C. Adembri, G. Villa, A. Giarratano, S. Maurizio Raineri, A. Cortegiani, F.
698 Montalto, M.T. Strano, V.M. Ranieri, C. Sandroni, G. De Pascale, A. Molin, P. Pelosi,
699 L. Montagnani, R. Urbino, I. Mastromauro, F.G. De Rosa, T. Cardoso, S. Afonso, J.
700 Gonçalves-Pereira, J.P. Baptista, A. Özveren, DALI: Defining antibiotic levels in

- 701 intensive care unit patients: Are current β -lactam antibiotic doses sufficient for
702 critically ill patients?, *Clin. Infect. Dis.* 58 (2014) 1072–1083.
703 <https://doi.org/10.1093/cid/ciu027>.
- 704 [13] J.A. Roberts, M. Ulldemolins, M.S. Roberts, B. McWhinney, J. Ungerer, D.L.
705 Paterson, J. Lipman, Therapeutic drug monitoring of β -lactams in critically ill patients:
706 Proof of concept, *Int. J. Antimicrob. Agents.* 36 (2010) 332–339.
707 <https://doi.org/10.1016/j.ijantimicag.2010.06.008>.
- 708 [14] S.L. Parker, F.B. Sime, J.A. Roberts, Optimizing dosing of antibiotics in critically ill
709 patients, *Curr. Opin. Infect. Dis.* 28 (2015) 497–504.
710 <https://doi.org/10.1097/QCO.0000000000000206>.
- 711 [15] F.B. Sime, M.S. Roberts, J.A. Roberts, T.A. Robertson, Simultaneous determination of
712 seven β -lactam antibiotics in human plasma for therapeutic drug monitoring and
713 pharmacokinetic studies, *J. Chromatogr. B Anal. Technol. Biomed. Life Sci.* 960
714 (2014) 134–144. <https://doi.org/10.1016/j.jchromb.2014.04.029>.
- 715 [16] R. Cazorla-Reyes, R. Romero-González, A.G. Frenich, M.A. Rodríguez Maresca, J.L.
716 Martínez Vidal, Simultaneous analysis of antibiotics in biological samples by ultra high
717 performance liquid chromatography-tandem mass spectrometry, *J. Pharm. Biomed.*
718 *Anal.* 89 (2014) 203–212. <https://doi.org/10.1016/j.jpba.2013.11.004>.
- 719 [17] I.L. Tsai, H.Y. Sun, G.Y. Chen, S.W. Lin, C.H. Kuo, Simultaneous quantification of
720 antimicrobial agents for multidrug-resistant bacterial infections in human plasma by
721 ultra-high-pressure liquid chromatography-tandem mass spectrometry, *Talanta.* 116
722 (2013) 593–603. <https://doi.org/10.1016/j.talanta.2013.07.043>.
- 723 [18] M. Paal, M. Zoller, C. Schuster, M. Vogeser, G. Schütze, Simultaneous quantification

- 724 of cefepime, meropenem, ciprofloxacin, moxifloxacin, linezolid and piperacillin in
725 human serum using an isotope-dilution HPLC–MS/MS method, *J. Pharm. Biomed.*
726 *Anal.* 152 (2018) 102–110. <https://doi.org/10.1016/j.jpba.2018.01.031>.
- 727 [19] S. Lefeuvre, J. Bois-Maublanc, L. Hocqueloux, L. Bret, T. Francia, C. Eleout-Da
728 violante, E.M. Billaud, F. Barbier, L. Got, A simple ultra-high-performance liquid
729 chromatography-high resolution mass spectrometry assay for the simultaneous
730 quantification of 15 antibiotics in plasma, *J. Chromatogr. B Anal. Technol. Biomed.*
731 *Life Sci.* 1065–1066 (2017) 50–58. <https://doi.org/10.1016/j.jchromb.2017.09.014>.
- 732 [20] C. Wharf, U. Kingdom, Guideline on bioanalytical method validation Guideline on
733 bioanalytical method validation Table of contents, 44 (2012) 1–23.
- 734 [21] B.K. Matuszewski, M.L. Constanzer, C.M. Chavez-Eng, Strategies for the assessment
735 of matrix effect in quantitative bioanalytical methods based on HPLC-MS/MS, *Anal.*
736 *Chem.* 75 (2003) 3019–3030. <https://doi.org/10.1021/ac020361s>.
- 737 [22] R. Guilhaumou, S. Benaboud, Y. Bennis, C. Dahyot-Fizelier, E. Dailly, P. Gandia, S.
738 Goutelle, S. Lefeuvre, N. Mongardon, C. Roger, J. Scala-Bertola, F. Lemaitre, M.
739 Garnier, Optimization of the treatment with beta-lactam antibiotics in critically ill
740 patients, *Crit. Care.* 23 (2019) 104.
741 <https://ccforum.biomedcentral.com/articles/10.1186/s13054-019-2378-9>.
- 742 [23] J.S. Mortensen, B.P. Jensen, M. Zhang, M. Doogue, Preanalytical Stability of
743 Piperacillin, Tazobactam, Meropenem, and Ceftazidime in Plasma and Whole Blood
744 Using Liquid Chromatography-Tandem Mass Spectrometry, *Ther. Drug Monit.* 41
745 (2019) 538–543. <https://doi.org/10.1097/FTD.0000000000000650>.
- 746 [24] T. Alebic-Kolbah, R. Demers, L. Cojocar, Dalbavancin: Quantification in human

- 747 plasma and urine by a new improved high performance liquid chromatography-tandem
748 mass spectrometry method, *J. Chromatogr. B Anal. Technol. Biomed. Life Sci.* 879
749 (2011) 2632–2641. <https://doi.org/10.1016/j.jchromb.2011.07.027>.
- 750 [25] R. Rigo-Bonnin, J. Gomez-Junyent, L. García-Tejada, E. Benavent, L. Soldevila, F.
751 Tubau, O. Murillo, Measurement of ceftolozane and tazobactam concentrations in
752 plasma by UHPLC-MS/MS. Clinical application in the management of difficult-to-treat
753 osteoarticular infections, *Clin. Chim. Acta.* 488 (2019) 50–60.
754 <https://doi.org/10.1016/j.cca.2018.10.034>.
- 755 [26] P.F. Fang, H.L. Cai, H. De Li, R.H. Zhu, Q.Y. Tan, W. Gao, P. Xu, Y.P. Liu, W.Y.
756 Zhang, Y.C. Chen, F. Zhang, Simultaneous determination of isoniazid, rifampicin,
757 levofloxacin in mouse tissues and plasma by high performance liquid chromatography-
758 tandem mass spectrometry, *J. Chromatogr. B Anal. Technol. Biomed. Life Sci.* 878
759 (2010) 2286–2291. <https://doi.org/10.1016/j.jchromb.2010.06.038>.
- 760 [27] W.C. Putnam, R.R. Kallem, V. Edpuganti, I. Subramaniyan, R.G. Hall, Development
761 and validation of a quantitative LC-MS/MS method for the simultaneous determination
762 of ceftolozane and tazobactam in human plasma and urine, *J. Chromatogr. B Anal.*
763 *Technol. Biomed. Life Sci.* 1159 (2020).
764 <https://doi.org/10.1016/j.jchromb.2020.122354>.
- 765 [28] L.E. Payne, D.J. Gagnon, R.R. Riker, D.B. Seder, E.K. Glisic, J.G. Morris, G.L. Fraser,
766 Cefepime-induced neurotoxicity: A systematic review, *Crit. Care.* 21 (2017) 1–8.
767 <https://doi.org/10.1186/s13054-017-1856-1>.

768