Laboratoire de Chimie Physique -Matière et Rayonnement **Philippe JONNARD** Activités de l'équipe 'Interfaces, Multimatériaux, Sources et **Optiques** X Sorbonne Universités, UPMC Univ Paris 06, LCPMR, Paris, France CNRS UMR 7614, LCPMR, Paris, France





Equipe IMSOX

« Interfaces, Multimatériaux, Sources et Optiques X »

https://lcpmr.cnrs.fr/content/interfaces-mutimatériaux-sources-et-optique-x

Laboratoire de Chimie Physique - Matière et Rayonnement

Département/UFR de Chimie

Faculté des Sciences et Ingénierie

Sorbonne Université

(fusion Paris VI – Paris IV au 1er janvier 2018)










Topics IMSOX team

Characterization of nanometric multimaterials by x-ray spectroscopy and x-ray optics


Development of spectroscopic and optical devices for materials characterization



Applications of multilayers

- Solar images from spatial telescopes
- Photolithography @ 13.5 nm
- Optical components for synchrotron beamlines
- Optical components for new x-ray sources (x-ray laser, high harmonic generation, ...)
- Optics for diffraction apparatus


X-ray emission spectroscopy



Imaging at a given wavelength

Photolithography @ 13.5 nm, 6.x nm

Motivation and challenges



	λ= 500nm	λ= 30.4nm	F
d	~125nm	Co_2.6nm/Mg_13.6nm	
σ	~10 nm	~0.5-0.8nm	

R multiplied by DW factor exp(-2[2π σ sinθ / λ]²) exp(-2[π σ / d]²)

Interface is the largest challenge for EUV multilayer!

X-ray reflectivity $R_{-}\theta = I_{r}(\lambda)/I_{0}(\lambda)$ $I_{r}(\lambda)$ $R_{-}\lambda = I_{r}(\theta)/I_{0}(\theta)$

 $\begin{array}{l} \text{2d sin}\theta = p \; \lambda \\ \text{Bragg law} \end{array}$

Ni (50 nm) 1E-8 épaisseur Reflectivity (a.u.) 1E-9 densité 1E-10 1E-11 rugosité 1E-12 2,0 0,0 0,5 1.5 2,5 1.0 Angle (°)


multilayer

 $I_0(\lambda)$

Géométrie

Determination of parameters of the stack

- thickness
- roughness
- density / optical index of the various layers



Observation of the energy distribution of the emitted radiation Information about the chemical state of the emitting element


Incoming beam

- Electrons
- Particles
- X-rays

Detection / dispersion

- WDS : crystal spectrometer
- EDS : Si(Li), SDD, bolometer

NMR spectroscopy




HF: hyperfine field

In ferromagnets, the magnetic field H₀ is the field inside the material on the nuclei sites (20 T for Co nuclei in bulk Co)

The HF strength depends on the local symmetry and local chemical environment of the probed nuclei: nearest neighbour local structure and/or local chemical environment.

Mg/Co multilayers (1)



• Design

Si / [Co (2.67 nm)/Mg (13.61 nm)]x30

Fabrication
 Deposition by DC-MS
 Ultra-smooth Si substrates (σ~0.3nm)
 High purity sputtering targets
 High purity sputtering gas
 High precision deposition rates


• Measurement High reflectiviy: 40.3% Narrow band-width: 1.32 nm

Co/Mg is a promising multilayer for EUV applications

Mg/Co results (2)

Design

Si substrate / [Co (2.55 nm)/Mg (14.45 nm)]x30/ B_4C (3.5 nm) B_4C : capping layer to prevent oxidation




R_sim = 56.5% @ 25.2 nm (49.2 eV) 45° s-pol

for perfect multilayer: without roughness or interdiffusion

X-ray Reflectivity

- High structural quality
- Thickness close to the design
- Interfacial roughness 0.6 nm

Mg/Co results (3)




EUV reflectivity

- High reflectance
- 75% of the value simulated for perfect multilayer
- Assymetrical peak owing to Mg L edge


Mg/Co results (4)

X-ray Emission Spectroscopy



- Mg atoms as in pure Mg
- No interaction between Mg and Co layers

Mg/Co results (5)




Nuclear Magnetic Resonance

Spectroscopy

- Co atoms as in pure Co
- No interaction between
 - Co and Mg layers
 - Faint peak due to Co


atoms in contact with the Mg layer

Mg/Co results (6)




TEM 100 keV

Sample prepared by focused ion beam Thinning down to ≈40 nm High-resolution images Co and Mg layers crystallized EELS Abrupt interfaces




Mg/Co results (annealing)




EUV reflectivity

- Stable up to 280° C
- R divided by 2 after 305° C
- No more reflectivity after 400° C

Mg/Co results (annealing)






Nuclear Magnetic Resonance Spectroscopy


- Co atoms as in pure Co
- No more interfacial Co atoms

XES

• Mg atoms as in pure Mg


Mg/Co results (annealing)





Mechanical deficiency due to Co clustering

X-ray standing waves



Berkeley Lab.

In (or close to) **Bragg conditions**

- a strong standing wave develops inside and outside the multilayer
- having the period of the multilayer

Generally used

- In hard x-ray range
 - In glancing incidence
 - With crystal or
 - With multilayer (set of **bilayers**) To probe
 - thin layer on top of the multilayer
 - **Interfaces of the multilayer** itself By using
 - Fluorescence or photoelectrons

What we do

- In soft x-ray range
- In glancing incidence and exit
- With multilayer


(set of **bilayers or trilayers**)

To probe


- Interfaces of the multilayer itself By using
- X-ray fluorescence

X-ray standing waves

Depth distribution of the electric field



Two XSW modes



Glancing incidence mode "Normal" or GI

Glancing exit mode "Kossel" or GE

A bit of history



M. von Laue

M. Von Laue Nobel prize 1914 for his discovery of the x-ray diffraction by crystals


First x-ray diffraction pattern


Hydrated copper sulfate P1 triclinic 1912: von Laue, Friedrich, Knipping



W. Kossel 1935 Diffraction of the f


Diffraction of the fluorescence produced within the crystal by the crystal itself






Fluorescence of Mg/Co multilayer

Si substrate / [Mg (5.1nm) / Co (2.6nm)]x30 / B4C (3.5nm)



GE experiment-simulations Mg/Co



abrupt interfaces

• Effect of the refractive indices

Introduction of Zr into Mg/Co stack


- Makes the stack stable upon annealing
- Improves the reflectivity (optimization on the optical path within the stack)

Si substrate / [**Mg (5.45nm) / Co (2.55nm) / Zr (1.5nm)**]x30 / B4C (3.5nm) Si substrate / [**Mg (5.45nm)) / Zr (1.5nm) / Co (2.55nm)**]x30 / B4C (3.5nm)

Two different orders of the layers

NMR spectroscopy (1)


Probe of the nearest neighbour local structure around the Co atoms



Mg/Co well-defined layers and interfaces

NMR spectroscopy (2)


Probe of the nearest neighbour local structure around the Co atoms



Mg/Co and Mg/Co/Zr well-defined layers and interfaces


NMR spectroscopy (3)

Probe of the nearest neighbour local structure around the Co atoms



GE measurements – Mg-Co-Zr


Different contrasts of the features for Co L emission depending on the sample



Same behavior in GI mode

Simulations – Mg-Co-Zr

Mg/Zr/Co – Co L emission



- Mg/Zr/Co described as Mg/Co₄Zr
- Mg/Co/Zr trilayer with nominal parameters

Mg-Co-Zr interfaces

From XRR, NMR, XES, XSW-XRF

- Sharp Zr-on-Co interfaces in Mg/Co/Zr
- Mixing at Co-on-Zr interfaces in Mg/Zr/Co

Surface free energy

- Co 2.0 J.m⁻²
- Zr 1.6 J.m⁻²

Asymmetric behavior of Zr-on-Co and Co-on-Zr interfaces

Diffraction de Kossel Caméra CCD

Pd/Y based periodic multilayers

Problem: the original design has bad optial performance

• Grazing incident X-ray reflectometry (Cu K radiation 8048 eV)

Low reflection, bad periodicity


• Severe interdiffusion between Pd and Y layers.

Solution: Derivative systems


• N₂ in the sputtering gas [1].

0% 2% 6%


- B₄C barrier layers [2].
- Pd-on-Y interfaces
- Y-on Pd interfaces
- all interfaces




[2] Prasciolu et al., "Thermal stability studies of short period Sc/Cr and Sc/B4C/Cr multilayers," Appl. Opt. 53(10), 2126 (2014).



Combine PIXE with Kossel diffraction



Proton X-ray emission beam Detection angle θ Sample



Energy dispersive CCD


Advantage: spatially resolved.

- Measurement of X-ray emission intensity in a squared area of about 13 x 13 mm.
- Native: 1024 x 1024 sensor array with 13 x 13 µm pixels.
- Our selection: 256 x 256 of 52 x 52 μm.

>>> compromise of angular/energy resolution

- No need to scan the angle.
- We place a 200 µm beryllium film in front of the camera to filter the scattered protons. Protection of the camera.
- Acquisition time significantly reduced compared to SDD.

Days to 2 hours >>>potential of in-situ measurements (annealing test, oxidation test)




X-ray emission spectrum

- B₄C(2.5 nm)/[Pd(2 nm)/B₄C(2nm)/Y(2nm)]_{x40}/Si
- Energy resolution is about 150 eV at 5.9 keV (Mn Kα).
- Pd Lα Lβ, Y Lα Kα, Si Kα can be observed
- Pile up effect due to multiple photons detection
- Energy calibration of the spectrum is carried out with Pd Lα peak (2.838 keV) and Y Kα peak (14.958 keV)



- Y Lα, Si Kα are in the low efficiency part of the X-ray camera.
- Pd Lβ does not bring additional infomation
- We measure Pd Lα emission





X-ray spectrum of Pd(2 nm)/ $B_4C(2 nm)/Y(2 nm)$ multilayer:


- a) overview (log scale)
- b) zoomed range of interest (linear scale).

Kossel curves

- The original image obtained by X-ray color camera contains the intensity of emitted photons of all energies.
- Kossel curve: the angular distribution of the characteristic X-ray emission.
- Too get the Kossel curve of Pd Lα emission, we need 3 steps:
- 1. Select region of interest (ROI) on the spectrum.
- 2. "Filter" the original image by applying ROI.
- 3. Integrate the "filtered" image along the vertical pixels.









Calibration of detection angle

- Pixel of the detector grazing exit angle
- 4-layer sample: B₄C(2.5 nm)/[B₄C(1 nm)/Pd(2 nm)/B₄C(1 nm)/Y(2 nm)]_{x40}/Si
- Ist and 2nd order of Kossel ocillation
- a) experimental curve in channel (pixel)
- b) simulated curve in angle (degree)
- Precise values obtained from the derivatives of the curves




Non-destructive test




Kossel curves of Pd La emission

- 2-hour acquisition time. Curve statistics can be further improved by increasing the acquisition time.
- Angles are adjusted to the Bragg angle, which is in most cases located in the center of the Kossel oscillation.
- The periodicity of the originally designed Pd/Y multilayer is totally compromised.
- Nitrogen reduces the interdiffusion.
- Kossel features can be distinguished for samples with different B₄C barrier layers.
 Structural sensitivity.



Kossel curves of the series of samples:


a) Deposited with nitrogen in the sputtering gas.




Kossel curves of the series of samples:

b) With B₄C barrier layers. The notation gives multilayer structure and layer thickness respectively.

Simulation




- Code for XRF-Kossel curve.
- Approach.
- Uniform ionization.
- Still not perfect due to background.
- Improvement needed.



B4C/Pd/B4C/Y Pd L alpha emission

Multicouches gravées


Multilayer mirror (MM)



Bad spectral resolution Owing to large diffraction pattern

1D periodical structure used at the 1st (2nd) diffraction order

Bragg law $2d \sin\theta = n\lambda$



d: multilayer periodD: grating period ΓD : grating line width1 - ΓD : grating line spacingCalculations: $D = 1 \ \mu m$ $\Gamma = 0.35$


2D periodical structure used:

- at the 0th order of the grating
- at the 1st (2nd) diffraction order

Bragg law $2d \sin\theta = n\lambda$



Reflectivity curves: comparison MM/EM




Etching process


Scanning electronic microscopy images


Aspect ratio:

- aimed 300:1200
- experimental ~ 310:1210




Acc.V Spot Magn Det WD |------- 2 μm 15.0 kV 2.0 8000x TLD 6.2 ech 30






Example of an etched multilayer in a curve-crystal spectrometer: C Kα emission band of graphite





et le plan (0001) du graphite

anisotropic graphite crystal



Réseau de diffraction en transmission






Tracé de rayons Spectromètre Johann

RAY-UI: A powerful and extensible user interface for RAY

P. Baumgärtel', M. Witt, J. Baensch, M. Fabarius, A. Erko, F. Schäfers, and H. Schirmacher

Citation: AIP Conference Proceedings 1741, 040016 (2016); doi: 10.1063/1.4952888



Johann R = 500 mm InSb (111) Si Ka = 1740 eV

