

Motor imagery-based Brain Computer Interface: theory and practice with OpenVibe

C. Zielinski

Let's start with a little demo... EEG in real time





OPENBCI DIY NEUROTECHNOLOGIST' S STARTER KIT

\$ 1,499.99 USD

Split your purchase into monthly installments with **shop** Pay Learn more

of Channels

Up to 8 channels

Quantity

- 1 +

Add to cart

Product Overview

This bundle includes a Cyton, Ultracortex, and additional electrodes & accessories for EEG, EMG, and ECG sensing. Available in 8 or 16 channels.



OPENBOLI CYTON BIOSENSING BOARD (8-CHANNELS)

\$ 999.00 USD

Pay in 4 interest-free installments of \$249.75 with shop@@_Learn more Quantity - 1 + Add to cart

▲ Share

Product Overview

The board that started it all! The Cyton allows you to gather 8-channels of scientifically-validated physiological data.

Cyton board 8-channels Dry electrodes 250 Hz Ultracortex Markup IV

OpenBCI



- Time series widget
- Frequency spectra
- Band power
- Networking: protocols to
 - output data (Serial, LSL, OSC)





Brain Computer Interface





by thought



lindAffect

Succeed in having a common language between the brain (seen through the EEG) and the computer

Co-adaptation: Human learning / BCI processing adjustment

Various technics based on different forms of electrical activities, **spontaneous** (wave rhythms linked to mental imagery) or **evoked** (by visual or acoustic stimuli) (ERP, P300)

Active BCI: Mental imagery

Mental tasks: motor-imagery of the limb, mental calculation or navigation

Active BCI that employs the user's endogenous brain activity in the absence of any external stimulus (*Pfürtscheller & Da Silva, 1999*)

Initially: communication for paralyzed people (yes/no)

Today: neurorehabilitation of stroke patients – MI neurofeedback effective in reestablishing motor movements immediately after a stroke (Mattia 2016, Pichiorri 2015)

Active BCI: Mental imagery



(Pfürtscheller et al 2011)



Active BCI: Mental imagery

Mental task will induce in specific cortical regions and in specific frequency bands:

- 1°) Event-Related Desynchronisation (ERD while being performed)
- 2°) Event-Related Synchronisation (ERS once the user has stopped)



(adapted from Lemm et al., 2009)

Event-Related Desync/Synchronization

Decrease or increase in synchrony of the underlying neuronal populations reflecting cortical activation and deactivation (*Pfurstcheller & Da Silva, 1999*)





Clinical Neurophysiology Volume 110, Issue 11, 1 November 1999, Pages 1842-1857

Invited review

Event-related EEG/MEG synchronization and desynchronization: basic principles

G. Pfurtscheller ^a [∧] [∞], F.H. Lopes da Silva ^b

Neuronal networks become asynchronous during mental activity; **ERD** – correlated with mental **activity**; **ERS** – correlated with mental **inactivity**





Figure 2.2: The Wilder Penfield Homunculus: the drawing of the limbs along the M1 region shows the regions responsible for generating a motor control signal to different parts of the body.[Adapted from [19]]

The Primary Motor Cortex (PMC) is responsible for planning and executing movements

Correspondence between the PMC areas and the various muscle groups

Movements tasks induce changes to brain activity (ERS/ERD)



Fig. 10. Movement-specific location of the beta ERS in one subject displaying a somatotopic organization of the beta oscillations after finger, arm and foot movement. Note the different subject-specific frequency bands lowest with finger and highest with arm and foot movement, respectively. 'Black' indicates location of maximal ERS.

Pfurstcheller & Da Silva, 1999



Synchro and Desynchronisation of brain waves at motor cortex areas, whether by preparing the movement of the body or by imagining it

Movements tasks induce changes to brain activity visible in the EEG





Hoffmann et al., 2006 : "A good feature extraction method should reduce the dimension of the input data as much as possible while keeping all the information necessary for classification."

Motor imagery: Graz protocol

I- training the system

Series of MI-task (ex. Left/Right) [feature extraction to train a classifier]





II- training the user

MI-task with feedback.

User goal: to develope effective strategies to increase system performance

Pfurtscheller & Neuper, "Motor Imagery and Direct Brain-Computer Communication", 2001

Graz protocol





Multiple sessions, each of which multiple sequences (4-6 runs - fatigue>6) One run:

10-20 trials per class (Left/Right) ~ 7 min Indication of which MI task is to perform





Feedback



- (1) an offline calibration stage, during which the system's settings are determined
- (2) an online operational stage, during which the system recognizes the user's brain activity patterns and translates them into application commands

BCI research community is currently searching for solutions to help avoid the costly offline calibration stage

OpenVibe



openvibe.inria.fr

An open-source software platform

- Design, test and use BCIs
- Real-time EEG acquisition, processing and visualization
- Offline or batch analysis of large datasets

Wide EEG device compatibility (EGI, Brainamp, Biosemi,... C++ driver)

Compatibility with 🔁 & 📣

BCI paradigm demos available + tutorials

Visual programmation (GUI boxes) + own plugin development in C++





scenarios > bci-examples	
Nom	
	erp-recording
	filter-bank-CSP
	handball
1	motor-imagery
1	motor-imagery-CSP
	neurofeedback
	p300-magic-card
	p300-speller
	p300-speller-xDAWN
	riemannian
	spaceship
	ssvep
	ssvep-mind-shooter
	ssvep-mind-shooter-classic

OpenVibe working flow



Dotted gray: communication protocol, solid blue: data flow, and colored boxes represent processes.

Acquisition server: Configuration to receive raw data depending on the device system

Designer: Define and Launch stimulation, online and offline data processing

Acquisition server: Configuration to receive raw data depending on the device system

Designer: Define and Launch stimulation, online and offline data processing

OPENVIBE – BOX TUTORIALS

power-spectrum.xml

You can browse each box' documentation by selecting the box and pressing F1

OpenViBE Designer

OPENVIBE – BOX

connectivity_measure.xml

OPENVIBE – BOX TUTORIALS

wavelet-decomposition.xml

Artifact removal, PCA, Machine learning...

Wavelet decomposition example

The scenario produces uniform random noise signal which rides on a low-frequency sine-wave carrier.

This signal is decomposed and then reconstructed by *Discrete Wavelet Transform* and *Inverse DWT* boxes.

Observing the outputs A and D shows that the low-frequency part A is now more clean of the random noise, some of which is carried now in the high-frequency part D. The reconstructed signal should closely resemble the original one.

Openvibe: Graz protocol + OpenBCI

- mi_1_stim-acq: stimulation (L/R) while recording the training data set
- mi_2_classifier-trainer: classifier parameters estimate
- mi_3_online_osc: classification in real time with feedback
- + send command by OSC to a web app

Ultracortex Mark IV Node Locations (35 total)

Based on the internationally accepted **10-20 System** for electrode placement in the context of EEG research

Motor imagery: data processing pipeline

Linear discriminant Analysis

- Find the projector line that gives:
- maximum separation between means of projected classes
- minimum variance within each projected class
- Solution: eigenvectors based on between-class and within-class covariance matrices

Example - Projector 2: best discriminant

Openvibe: Graz protocol + OpenBCI

- Multiple sessions with Franziska
- Low classification accuracy by imagining the movement (60% with the Left class being better discriminated) Cross-validation test accuracy is 55.836040% (st

Cross-validation test accuracy is 55.836040% (sigma = 5.424622%) Cls vs cls 1 2 Target 1: 58.5 41.5 %, 1020 examples Target 2: 46.9 53.1 %, 1020 examples Training set accuracy is 60.441176% (optimistic) Cls vs cls 1 2 Target 1: 63.4 36.6 %, 1020 examples Target 2: 42.5 57.5 %, 1020 examples

 By actually doing the movement: good accuracy (92%!) – classification based on the mu-beta features is very reliable

> Classifier trainer> Cross-validation test accuracy is 92.354288% (sigma = 7.365529%) Classifier trainer> Cls vs cls 1 2 Classifier trainer> Target 1: 97.1 2.9 %, 1020 examples Classifier trainer> Target 2: 12.4 87.6 %, 1020 examples Classifier trainer> Training set accuracy is 92.156863% (optimistic) Classifier trainer> Cls vs cls 1 2 Classifier trainer> Target 1: 96.9 3.1 %, 1020 examples Classifier trainer> Target 2: 12.5 87.5 %, 1020 examples

Co-adaptation user/machine - substantial training to find the correct imagination strategy that leads to the most distinct ERD/ERS

Currently it can only be used with a maximum of 3 or 4 different MI tasks to ensure maximal performances (classification accuracy)

BCIs improvement

« Illiterate users »

Most current systems do not properly establish this mutual understanding

- 10 to 30 % of users unable to control a BCI: "BCI illiteracy" or "BCI deficiency" phenomenon
- MI-BCI : even for not "illeterate", performance rather low
 - 75% of classification accuracy for 2 class MI-BCIs
 - 20% reach 80 to 100%
- 2 main factors

(1) Signal processing (Stage I) – many studies over the last years
(2) The user-training role: mostly neglected – inefficient training protocol (instructions, training tasks, feedback and training environment)

- Controlling an BCI requires acquisition of specific skills = to generate stable and distinct brain activity patterns while performing the tasks
- Understand how the humain learn and how to adapt training process

BCIs improvement

Summary guidelines for designing more effective training protocols (C. Jeunet)

 \bigcirc

Training task

Progression (increasing difficulty), adaptative user-specific Self-paced and asynchronous sessions Preparatory training tasks (meditation)

Feedback

Visual with emotional connotations (smileys) – tactile is promising (channel much less saturated in interactive situations) / increase the quantity and quality of information (topography of cerebral activity...)

Training environment

Improvement of user motivation & experience: gamified protocols (ex. ball-basket paradigm, spaceship that must avoid asteroids) Integration of virtual reality. Ex "Use the force" – allows users to levitate a spaceship by imagining moving their feet

MI-based BCIs improvement

Continuous Tactile Feedback for Motor-Imagery based Brain-Computer Interaction in a Multitasking Context Camille Jeunet, Chi Vi, Daniel Spelmezan, Bernard N'Kaoua, Fabien Lotte, Sriram Subramanian

interest; blue - referent electrodes; green - ground

electrode.

motor imagery EEG mental task visual cue online vibrotactile reedback

 b) Experimental setup of BCI training.
 Red dots - electrodes of interest; blue dot - referent electrode; yellow dots - vibrotactile actuators.

"Vibrotactile feedback during MI training induced significant enhancement of ERD activity only for non-dominant, left hand over contralateral motor cortex area measured in C4 electrode"

New Results

A Follow this preprint

A BCI-based vibrotactile neurofeedback training improves motor cortical excitability during motor imagery

N. Grigorev, A. Savosenkov, M. Lukoyanov, A. Udoratina, N. Shusharina, A. Kaplan, A. Hramov, V. Kazantsev, S. Gordleeva doi: https://doi.org/10.1101/2021.02.28.433220

Now published in IEEE Transactions on Neural Systems and Rehabilitation Engineering doi: 10.1109/tnsre.2021.3102304

