

Coupled dielectric spectroscopy and DSC for ferroelectric materials

Gwenn Morvézen, N Bréfuel, Daniel Bourgault, Hervé Guillou, A Sylvestre

► To cite this version:

Gwenn Morvézen, N Bréfuel, Daniel Bourgault, Hervé Guillou, A Sylvestre. Coupled dielectric spectroscopy and DSC for ferroelectric materials. 7th Central and Eastern European Conference on Thermal Analysis and Calorimetry (CEEC-TAC7), Aug 2023, Brno, Czech Republic. , 2023. hal-04197731

HAL Id: hal-04197731 https://hal.science/hal-04197731

Submitted on 12 Sep 2023

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Coupled dielectric spectroscopy and DSC for ferroelectric materials

<u>Gwenn Morvézen^{1,2}</u>, N. Brefuel¹, D. Bourgault², H. Guillou², A. Sylvestre¹ ¹ Univ. Grenoble Alpes, CNRS, Grenoble INP, G2Elab, 38000 Grenoble, France ² Univ. Grenoble Alpes, CNRS, Grenoble INP, Institut Néel, Grenoble, France


ABSTRACT


We present here a new planar calorimeter coupling dielectric spectroscopy and DSC measurements permitting to link dielectric and thermodynamic behaviors. It is especially suited for sintered ferroelectric pellets or electroactive polymers with transition temperatures from 10 to 140°C and could be applied to thin films with high sensitivity. To show these possibilities, a hybrid molecular ferroelectric material named dabcoHBF₄ is studied and new features are measured in its ferroelectric-paraelectric transition.

Studied material


Developped devices

Calibration







Calibration process : 1) Voltage induced by power pulses @ constant Temperature



2) Sensitivity S $@T_i = slope of U_i vs P$







- Dielectric + DSC measurement
- DC voltage + DSC measurement
- Controlled gas athmosphere


(M/Lm)

Dielectric-DSC Coupled measurements


3) Sensitivity S vs T







New DSC-dielectric correlated structures shown at low heating/cooling rate



DabcoHBF₄

Lead free ferroelectrics Aqueous synthesis High T_c for organics $Pr = 4.9 \ \mu C/cm^2$ Large dielectric constant range

Ref: 1) P.P. Shi et. al . Am. Chem. Soc. 139, 1319–1324 (2018) 2) A. Olejniczak et. al Crystal Growth & Design, 18, 11, 6488-6496 (2018) 3) N. Bréfuel et.al J. Mat. Science, 56, 18582-18591. (2021)



Conclusion and perspectives

We successfully developed a new calorimeter based on thermoelectric sensors made for planar samples (such as pellets, thick films or polymer sheets) coupled with dielectric spectroscopy or voltage applications. First, we calibrated it and present its specifications. We used that coupled DSC-dielectric spectroscopy measurement to study a hybrid molecular ferroelectric material named DabcoHBF₄ at low scanning speeds and showed new corelated substructures in the ferroelectric-paraelectric transition.

Now, some of the following steps planned to be made soon to develop that new tool :

Device optimizations to be made

- Lower thermal noise
- Minimize temperature ramps drifts
- Upgrade cooling power

New measurements

- On thin films and polymers
- With voltage, for electrocaloric effects
- On liquids and solutions



