

A vascularized glioblastoma multiforme within a 3D perfused microphysiological system: combining a self-organized microvasculature and a central venule in a hydrogel

Agathe Figarol, Roumaïssa Mosbah, Muhammad Hamidullah, Léo Durand, Marion Pouit, Sylvain Chamouton, Franck Chollet, Thérèse Leblois, Vincent Humblot

▶ To cite this version:

Agathe Figarol, Roumaïssa Mosbah, Muhammad Hamidullah, Léo Durand, Marion Pouit, et al.. A vascularized glioblastoma multiforme within a 3D perfused microphysiological system: combining a self-organized microvasculature and a central venule in a hydrogel. microphysiological systems society (MPSS), Jun 2023, BERLIN, Germany. hal-04224584

HAL Id: hal-04224584 https://hal.science/hal-04224584

Submitted on 12 Oct 2023

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

femto-st SCIENCES & **TECHNOLOGIES**

A vascularized glioblastoma multiforme within a 3D perfused microphysiological system: combining a self-organized microvasculature and a central venule in a hydrogel.



Figarol Agathe ¹, Mosbah Roumaïssa ^{1,2}, Hamidullah Muhammad ¹, Durand Léo ^{1,3}, Pouit Marion ^{1,3}, Chamouton Sylvain ^{1,3}, Chollet Franck¹, Leblois Thérèse¹, & Humblot Vincent¹.

Introduction

The most frequent and aggressive brain tumor is the glioblastoma multiforme (GBM). After diagnostic, the median patient survival is only of 15 months after diagnosis. There is a crucial need in unlocking the development of innovative drug treatments.

The blood-brain barrier (BBB):

- inhibits > 98% drugs to reach brain tissues
- hampers the development of new chemotherapy against GBM.



The clinical failure rate could be lowered with:

- drug nano-vectorization
- screening on pertinent preclinical models.



The 3D-Glimpse project

- Development of a BBB and GBM-on-chip
- Dynamic culture under perfusion to mimic the blood flow
- Instrumentation with
- biosensors for an integrated detection of nanocarriers transport.
- Validation as an alternative to animal testing for drug screening

Development of the BBB-on-chip.

3D coculture of human brain cells (HBEC:AC:PC at 1:0.5:2 ratio) in an ECM hydrogel (collagen I + fibrin) enabled their self-organization as a vascularized tissue. Optimized as free drops, the hydrogel was then poured in the polydimethylsiloxane (PDMS) microchip prototype over a needle, to create a central venule. The first perfusion assays of this venule confirmed the microchip design but induced as yet some degradation of the hydrogel.



After software-assisted design, molds were 3D printed in polylactic acid to cast PDMS microchips.









The central venule carved into the collagenfibrin hydrogel retracted in the first 24 h then remained stable over time.

Similarly, a guide was designed for the reproductive and controlled needle insertion.

Hydrogel after 1 week of static culture under bright field and fluorescent microscope. Tight junctions (CD31) of endothelial cells were immunostained and showed the self-organized microarchitecture as a capillary network.





Preliminary results have been obtained for endothelial cells seeding and perfusion.

BBB-on-chip in a healthy versus tumoral context.

The microchip will be adapted to fit the microenvironment of a vascularized brain tumor, with adjunction of GBM cells (U87) to create a GBM-on-chip.



Instrumentation of the BBB and GBM-on-chip.

Functionalized piezo-electrical substrates will be used as biosensors to detect the transport of nano-vectorized drugs through the BBB in healthy and pathological conditions.



A micro-processed porous membrane was developed for a controlled separation of the microchip chambers.





Activation of PC and AC due to the presence of GBM cells? Enhanced permeability of the capillaries?

Perspectives: transport and therapeutic efficiency screening.

Besides the transport, nanocarriers will be screened for their delivery specificity, and their toxicity towards targeted and off-target cells. The goal is to guide their dosage adjustment to prevent the strong side-effects of those chemotherapeutic treatments. ondatior



