Development of high resolution and radiation hard electromagnetic calorimeter for the DVCS experiment

Ho San KO

Institut de Physique Nucléaire d'Orsay

2019 IEEE NSS MIC

Co-authors : J. Bettane, G. Hull, C. Munoz Camacho, T. Nguyen Trung, and E. Rindel











In collaboration with

- Catholic University of America, USA
 - PbWO₄ light yield & transmittance measurement cross-check
- Justus Liebig University Giessen, Germany (many works have been done)
 - PbWO₄ light yield & transmittance measurement cross-check
 - PbWO₄ radiation hardness measurement cross-check
 - Optical bleaching method



- Deeply Virtual Compton Scattering(DVCS) for GPDs measurements
- Neutral Particle Spectrometer(NPS) for DVCS experiments in Hall C
- NPS simulations and crystal optical properties measurement
 - Energy resolution simulation
 - Radiation hardness measurements
 - Crystal optical properties recovery tests

Generalized Parton Distributions



Form Factors: <u>via elastic scattering</u> -charge & magnetization spatial distribution



Parton distribution: <u>via deep inelastic scattering</u> -Longitudinal momentum & helicity distribution of partons Generalized Parton Distributions: <u>via exclusive reactions</u> [eg. DVCS] -Transverse position distribution and longitudinal momentum of partons

 $f(x,r_1)$

r

 $\delta z_{\perp} \sim 1/Q$

Neutron Particle Spectrometer(NPS) in Hall C



 $e \ p \to e'p'\gamma$



(Exclusive process)

To detect neutral particles (γ , π_0), NPS will be installed on to the SHMS SHMS will be used as a cantilever

HMS detects scattered electrons NPS detects neutral particles

HMS : High Momentum Spectrometer SHMS : Super High Momentum Spectrometer

Continuous 11 GeV electron beam

Neutron Particle Spectrometer(NPS) in Hall C

<<pre><<picture : Emmanuel Rindel>>



Highest luminosity($\sim 10^{38} cm^{-2} s^{-1}$) of DVCS ever before

• 1080 (30X36) PbWO₄ blocks

each crystals : 2X2X20 cm3, ~670 gram

• PbWO₄ as a scintillator :

-Radiation hard
-140~200 photons/MeV
-Widely used in high energy physics as a calorimeter

Electronics : PMT, LED, and optical fiber

NPS Energy Resolution Simulation

 σ



$$\frac{\sigma}{E} = A \oplus \frac{B}{\sqrt{E}} \oplus \frac{C}{E}$$



<<pre><<picture : Julien Bettane>>

Purpose :

-NPS is made of 1080 crystals (30X36, stacked)

-Need structures for uniform arrays \rightarrow Better energy resolution

-More structure between crystals \rightarrow Worse energy resolution

 \rightarrow Check the dependency of the energy resolution of the NPS on the distance between the crystals(gap)



Recent NPS prototype : No material in the middle part

<< PrimEx also used PbWO₄ crystals for the calorimeter>>

PbWO₄ light yield's temperature dependence



<<Light yield measurement in temperature controlled box>>



The light yield has temperature dependence. \rightarrow The temperature of the calorimeter will be fixed at 18°C

31/10/2019

PbWO₄ Irradiation Test



Irradiation facility at Laboratoire de Chimie Physique (LCP), Orsay, France

- Radiation hardness measurement
 - Check crystal performance in high radiation environment
- Optical bleaching
 - Method to recover the optical properties of the crystals after the irradiation damage
 - 222TBq ⁶⁰Co source
 - Vary the distance from the source to change the dose rate
 - Can simultaneously irradiate 9 crystals

PbWO₄ Irradiation Test

Light transmittance measurement

J23

The radiation damage differs for every crystal

J43





Picture : \sim 20Gy/min \sim 4hrs with 60 Co at Laboratoire de Chimie Physique

31/10/2019

900

PbWO₄ Radiation Damage Recovery

- Ionizing radiation creates color centers
 - By electron traps or point structure defects
- Relaxation of color centers
 - Thermo-activation
 - Injection of specific energy
 - -Blue light : PMT sensitive

-Infrared light : PMT insensitive \rightarrow possibility of curing during data taking



12

PbWO₄ Optical Bleaching

PMTs, LEDs, optical fibers on the back side of the calorimeter

<<Crystal curing using blue light with optical fiber attached to the crystal>>







PbWO₄ Optical Bleaching



-The crystal was irradiated with 30Gy of dose -Applied two types of light illumination recovery

-Infrared(IR) light with optical fiber for 2hrs-Blue light with optical fiber for 2hrs

 \rightarrow Blue light works better in recovering the transmittance.

→Blue light for the crystal optical property damage recovery system of NPS Optical Bleaching

Quartz (\emptyset : 1000µm & 800µm) and plastic (\emptyset : 1000µm) optical fiber have no difference in the amount of recovery for 2hrs of curing after 30Gy of irradiation.

LED and Optical Fiber irradiation test

- LED and quartz optical fiber will also be used as a calibration tool for the calorimeter
 - Tests these equipments' radiation hardness



- Quartz optical fiber is radiation hard up to $\sim 1.5 \times 10^5 \ Gy$
- LED is not radiation hard
 - The LED's casing color changes after $\sim 1.0 \times 10^3 Gy$

Expected dose in the experiment : < 300 Gy

LED and Optical Fiber irradiation test



- The LED's plastic casing changes its color due to irradiation.
- Its light output changes.
 However, increases first and then decreases.

Expected dose in the experiment : < 300 Gy

➔ Further LED irradiation tests are planned with less radiation dose



- DVCS access to the GPDs
- NPS is needed in Hall C in order to perform the DVCS experiment
- NPS energy resolution calculation with Geant4
 - Detector structure : 0.5mm carbon material only in the front and the back side, 2cm each of the crystals
- Optical bleaching method for maintaining calorimeter's energy and position resolution
- Blue light curing (optical bleaching) via quartz optical fiber to be used for the NPS

Time Line of the NPS

- Detector construction at the end of 2019
- Delivery of NPS to Jefferson Lab at the beginning of 2020
- NPS tests at Jefferson Lab in 2020
- DVCS experiment expected in 2021