

## Studying Muscle Transcriptional Dynamics at Single-molecule Scales in Drosophila

Emma Leroux, Nourhene Ammar, Savannah Moinet, Thierry Pecot, Hadi

Boukhatmi

#### ► To cite this version:

Emma Leroux, Nourhene Ammar, Savannah Moinet, Thierry Pecot, Hadi Boukhatmi. Studying Muscle Transcriptional Dynamics at Single-molecule Scales in Drosophila. Journal of visualized experiments: JoVE, 2023, 199, 10.3791/65713. hal-04260328

### HAL Id: hal-04260328 https://hal.science/hal-04260328

Submitted on 26 Oct 2023

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

# Journal of Visualized Experiments An optimized method to study muscle transcriptional dynamics at single-molecule scales in Drosophila.

Manager and Duraft	oouroc	100	oprina.
N/IODUOORIDT LIROTT	Ma	int D	roft

Article Type:	Invited Methods Article - Author Provided Footage (APF)
Manuscript Number:	JoVE65713R2
Full Title:	An optimized method to study muscle transcriptional dynamics at single-molecule scales in Drosophila.
Corresponding Author:	Hadi Boukhatmi, Ph.D IGDR: Institut de Genetique et Developpement de Rennes Rennes, FRANCE
Corresponding Author's Institution:	IGDR: Institut de Genetique et Developpement de Rennes
Corresponding Author E-Mail:	hadi.boukhatmi@univ-rennes1.fr
Order of Authors:	Emma Leroux
	Nourhene Ammar
	Savannah Moinet
	Thierry Pecot
	Hadi Boukhatmi
Additional Information:	
Question	Response
Please indicate whether this article will be Standard Access or Open Access.	Open Access (\$3900)
Please confirm that you have read and agree to the terms and conditions of the author license agreement that applies below:	I agree to the <u>Author License Agreement</u>
Please confirm that you have read and agree to the terms and conditions of the video release that applies below:	I agree to the <u>Video Release</u>
Please provide any comments to the journal here.	

1	TITLE:			
2	Studying Muscle Transcriptional Dynamics at Single-molecule Scales in <i>Drosophila</i>			
3				
4	AUTHORS AND AFFILIATIONS:			
5	Emma Leroux <sup>1</sup> , Nourhene Ammar <sup>1</sup> ,	Savannah Moinet <sup>1</sup> , Thierry Pecot <sup>2</sup> , Hadi Boukhatmi <sup>1, *</sup>		
6				
7	<sup>1</sup> Institut de Génétique et Développ	ement de Rennes (IGDR), Université de Rennes, CNRS, 2 rue		
8	Prof. Leon Bernard, Rennes, France			
9	<sup>2</sup> BIOSIT, UAR 3480 US 018, Université de Rennes, 2 rue Prof. Leon Bernard, Rennes, France			
10				
11	Email addresses of the co-authors:			
12	Emma Leroux	emma.leroux@univ-rennes.fr		
13	Nourhene Ammar	nourhene.ammar@univ-rennes.fr		
14	Savannah Moinet	savannah.moinet@etudiant.univ-rennes.fr		
15	Thierry Pecot	thierry.pecot@univ-rennes.fr		
16	* Common and the south one			
1/	*Corresponding author:	hadi havulu atasi Quain asa ata fa		
18	Hadi Boukhatmi	nadi.bouknatmi@univ-rennes1.fr		
19				
20	Droconhila is a well established m	adal for studying key malasulas that regulate myogenesis		
21	However current methods are ins	ufficient to determine mPNA transcriptional dynamics and		
22	spatial distribution within syncytia	To address this limitation, we entimized an PNA fluerosconce		
23	spatial distribution within syncytia. To address this limitation, we optimized an RNA fluorescence			
24	moloculo scalo			
26				
27	ABSTRACT:			
28	Skeletal muscles are large syncytia r	nade up of many bundled myofibers that produce forces and		
29	enable body motion. <i>Drosophila</i> is a classical model to study muscle biology. The combination of			
30	both <i>Drosophila</i> genetics and advanced omics approaches led to the identification of key			
31	conserved molecules that regulate muscle morphogenesis and regeneration. However, the			
32	transcriptional dynamics of these molecules and the spatial distribution of their messenger RNA			
33	within the syncytia cannot be assessed by conventional methods. Here we optimized an existing			
34	single-molecule RNA fluorescence <i>in situ</i> hybridization (smFISH) method to enable the detection			
35	and quantification of individual mRNA molecules within adult flight muscles and their muscle			
36	stem cells. As a proof of concept, we have analyzed the mRNA expression and distribution of two			
37	evolutionary conserved transcription factors, Mef2 and Zfh1/Zeb. We show that this method can			
38	efficiently detect and quantify single mRNA molecules for both transcripts in the muscle			
39	precursor cells, adult muscles, and muscle stem cells.			
40				
41	INTRODUCTION:			
42	Adult skeletal muscles are made of	f differentiated multinucleated myofibers whose contractile		
43	properties generate movements. Muscle growth, maintenance, and regeneration rely upon			
44	muscle progenitors and muscle stem cells (MuSCs) that are specified during embryonic			

- development <sup>1</sup>. The myogenic program is finely controlled by a set of core myogenic transcription
   factors (TFs) (e.g., Pax3/7, MYOD, Mef2, and ZEB) <sup>2-4</sup>. Deciphering the molecular mechanisms
   regulating muscle biology is important for elucidating fundamental myology-related questions
- 48 and for possible therapeutic use in treating muscle-degenerative diseases.
- 49

50 Drosophila has a long history as a genetic model to study myogenesis <sup>5</sup>. It has recently emerged as a new model to investigate muscle regeneration <sup>6-8</sup>. The muscle structure and core myogenic 51 programs are highly conserved between flies and mammals. For example, the TFs Mef2 and 52 53 Zfh1/ZEB have a conserved function in regulating muscle development and regeneration <sup>3,9-11</sup>. 54 Adult Drosophila muscles, such as the Indirect Flight Muscles (IFMs), are formed from a specific 55 population of MuSCs, known as Adult Muscle Progenitors (AMPs) <sup>12</sup>. These AMPs are specified during embryogenesis and associate with epithelial structures such as wing and leg discs. 56 57 Throughout embryonic and larval stages, AMPs remain undifferentiated until metamorphosis when they engage in differentiation and fusion to form the IFMs <sup>13,14</sup>. The TF Spalt major (*spalt*, 58 salm) is expressed during flight muscle development and is necessary to determine their 59 structural identity <sup>15</sup>. Mef2 is another important myogenic TF that is essential for adult muscle 60 formation <sup>16,17</sup>. It is expressed in the AMPs and maintained in the adult IFMs <sup>10,18</sup>. While the 61 majority of AMPs differentiate into functional muscles, a subset escapes differentiation and 62 forms the adult MuSCs<sup>11</sup>. Similar to vertebrates, the TF Zfh1/ZEB is required to prevent 63 64 premature differentiation of the adult MuSCs and to maintain their stemness <sup>9,10</sup>.

65

Gene expression dynamics have been proven to regulate various muscle biological processes <sup>19,20</sup>. 66 67 The advent of high-throughput single-cell and single-nuclei RNA sequencing techniques has 68 enabled a comprehensive exploration of these transcriptional dynamics <sup>21,22</sup>. One notable 69 limitation of these approaches is their inability to provide the spatial distribution of mRNA 70 molecules in the multinucleated muscle fibers. These features can be investigated by singlemolecule fluorescence in situ hybridization (smFISH) that reveals two types of mRNA bodies: 1. 71 72 Individual mRNA molecules spread out in the cytoplasm and representing the mature RNA and 73 2. A maximum of two bright nuclear foci correspond to the nascent transcript and reveal the 74 transcriptionally active alleles <sup>23</sup>. Therefore, smFISH is a method of choice for individual mRNA 75 molecule quantification, investigating their spatial distribution and providing snapshots of the 76 gene transcriptional dynamics.

77

78 The smFISH method relies on a set of short fluorescent oligonucleotide probes specifically 79 designed to be complementary to the target transcript. Upon annealing, it generates high-80 intensity point sources permitting the detection of mRNA at a single-molecule level using confocal microscopy <sup>24</sup>. This method is increasingly applied for a wide variety of cell types, 81 including mammalian muscle tissues <sup>19,20</sup>. However, in *Drosophila* like other animal models, most 82 83 of the knowledge about adult muscle gene expression is derived from bulk molecular assays, 84 which lack quantitative information on the spatial location of mRNA molecules. Here we optimized a method for performing smFISH on muscle precursor cells and adult Drosophila 85 muscles <sup>23,25</sup>. This protocol includes a fully automated analysis pipeline for nuclei segmentation 86 and mRNA counting and localization. 87

88

89	
90	PROTOCOL:
91	
92	1. Wing disc and adult muscle dissection and preparation
93	
94	NOTE: Clean up the dissection bench and all the dissection instruments with RNAse inhibitor
95	solution ( <b>Table of Materials</b> ). Laboratory gloves should be worn during the whole experimental
96	procedure.
97	
98	1.1. Wing disc
99	
100	1.1.1. Place the L3 staged larvae in cold 1x phosphate-buffered saline (PBS) ( <b>Figure 1A</b> ).
101	
102	1.1.2. Rinse the larvae 3x with 1x PBS and keep them in this solution for 5 min.
103	
104	1.1.3. Using two pairs of sharp forceps, hold the larva from the anterior part and pull and discard
105	the rest of the body tissues (around 2/3).
106	
107	1.1.4. Hold both the first $1/3$ of the larva with one pair of forceps and the mouth books with the
108	other pair. Retract the mouth books inside the larva until it is fully inverted.
109	
110	NOTE: Dissect approximately 30 larvae per experiment, to ensure adequate sample size for
111	analysis.
112	
113	1.1.5. Observe that the wing discs are closely associated with the two primary branches of the
114	tracheae, which run along each side of the larva.
115	
116	1.1.6. Remove all the other larval components (leg discs and brain) to obtain a clean inverted
117	larval carcass with a pair of wing discs attached to the trachea ( <b>Figure 1B</b> ).
118	
119	1.1.7. Place the samples in a 1 mL centrifuge tube and fix them in 4% paraformaldehyde (PFA)
120	diluted in 1x PBS for 45 min at room temperature under agitation.
121	
122	NOTE: The sample agitation is achieved by an oscillating stirrer ( <b>Table of Materials</b> ).
123	
124	1.1.8. Wash the samples 3 x 5 min with 70% ethanol (EtOH, diluted in RNAse-free water) and
125	store them in the same washing solution at 4 °C for at least 2 days and up to 1 week.
126	
127	1.2. Indirect flight muscles
128	
129	1.2.1. Anesthetize the adult flies on a $CO_2$ fly pad and remove the heads, abdomens, legs, and
130	wings.
131	

132	1.2.2. Place the dissected thoraxes in cold 1x PBS and continue doing so until all samples are
133	dissected (Figure 1D).
134	
135	1.2.3. Prefix the thoraxes in 4% PFA, (1% Triton) for 20 min at room temperature under agitation.
136	
137	1.2.4. Rinse the samples 3 x 5 min with 1x PBT washing solution (1% Triton in 1x PBS).
138	
139	1.2.5. Position the thoraxes on a double-sided tape (Table of Materials) on a glass slide and bisect
140	them with a sharp microtome blade to produce two hemi-thoraces (Figure 1E,F).
141	
142	1.2.6. Fix the samples in 4% PFA (1% Triton) for 45 min at room temperature under agitation.
143	
144	1.2.7. Wash 2 x 20 min with 1x PBT washing solution (1% Triton, 1x PBS) at room temperature
145	under agitation.
146	
147	1.2.8. Place the hemi-thoraces in cold 1x PBS and use the forceps to gently isolate the IFMs from
148	the hemi-thorace (Figure 1G).
149	
150	NOTE: To achieve this, it is important to remove the cuticle as thoroughly as possible.
151	
152	1.2.9. To permeabilize the samples, store the isolated IFMs in 70% EtOH solution for at least 2
153	days and up to 1 week at 4 °C.
154	
155	NOTE: Isolate the IFMs from approximately 10 thoraxes per experiment to ensure adequate
156	sample size for analysis.
157	
158	2. Hybridization, immunostaining, and mounting
159	
160	NOTE: Similar hybridization procedure is applied for both wing discs and IFMs samples. It is critical
161	to isolate the wing discs from the larval carcasses and put them in 200 $\mu$ L of Buffer A prior to
162	starting the hybridization.
163	
164	2.1. Resuspend the probes in 400 µL of TE buffer (10 mM Tris-HCl. 1 mM EDTA, pH 8.0) for a final
165	stock solution of 12.5 µM.
166	
167	2.2. Transfer the samples to a 0.5 mL centrifuge tube containing Buffer A ( <b>Table of Materials</b> ).
168	
169	NOTE: By decreasing the volume of the container, it becomes easier to visualize the samples and
170	avoid losing them during the washing process.
171	
172	2.3. Wash the samples 2 x 20 min with Buffer A.
173	
174	2.4. Aspirate Buffer A, replace it with 400 $\mu$ L of the hybridization buffer ( <b>Table of Materials</b> ), and
175	prehybridize the samples for 30 min at 37 °C in a thermal mixer (Table of Materials).

176	
177	2.5. Dilute the probes to a final concentration of 125 nM for wing discs and 200 nM for adult
178	muscles, in the hybridization buffer and add the primary antibody.
179	
180	NOTE: The antibody dilution varies depending on the specific antibody being used (Table of
181	Materials). Adding antibodies is optional. Users may opt to forego protein staining and instead
182	use a similar smFISH protocol that does not include the addition of antibodies.
183	
184	2.6. Aspirate the hybridization buffer and replace it with 100 $\mu$ L of hybridization buffer containing
185	the probe and the antibody.
186	
187	2.7. Incubate the samples in the dark at 37 °C for at least 16 h under 300 rpm agitation in a
188	<mark>thermal mixer.</mark>
189	
190	<mark>2.8. Warm up an aliquot of Buffer A at 37 °C.</mark>
191	
192	2.9. Wash the samples for 3 x 10 min in the thermal mixer at 37 °C under 300 rpm agitation.
193	
194	2.10. During the last wash, dilute the secondary antibody (1/200) and DAPI (1/10,000) in Buffer
195	A and incubate the samples in this solution at 37 °C for at least 1 h.
196	
197	2.11. Wash the samples for 3 x 20 min with Buffer B ( <b>Table of Materials</b> ).
198	
199	2.12. Carefully transfer the samples on a microscope slide, wipe off the residual Buffer B, and add
200	30 µL of mounting medium, and cover with an 18 x 18 mm coverslip. Use a nail polish to seal the
201	preparation.
202	
203	2.13. Store the samples for up to 1 week at 4 °C.
204	
205	NOTE: However, it is recommended to perform the imaging as soon as possible to avoid probe
206	signal degradation.
207	
208	3. Imaging
209	
210	3.1. Acquire images with an xyz acquisition mode using a confocal microscope equipped with 40x
211	and 63x oil immersion objectives ( <b>Table of Materials</b> ). The smFISH signal is detected by the HyD
212	detector via a photon counting mode.
213	2.2. Leasts the complex using the DADI signal and UV least
214	3.2. Locate the samples using the DAPI signal and UV lamp.
215	2.2. To image the wing disc associated muscle progenitors and IFMs (Figure 2 and Figure 2) apply
210	the following setup for lasor lines and omission filters, select <b>DADL</b> (excitation (ex.) <b>405</b> pm
21/ 210	omission (om) 450 pm). Alora Elucr 488 (or 496 pm, om 510 pm), and Ouerer 570 for the
210 210	smEISH probes (av 647 nm om 670 nm)
213	sini isri probes (ex. <b>047 mii</b> , em. <b>070 mii</b> ).

220	
221	3.4. To image the adult MuSCs (Figure 2D), apply the following setup for laser lines and emission
222	filters: select DAPI (ex. 405 nm, em. 450 nm), Alexa Fluor 488 (ex. 496 nm, em. 519 nm), Alexa
223	Fluor 555 (ex. 555 nm, em. 565 nm), and Quasar 670 for the smFISH probes (ex. 647 nm, em. 670
224	nm).
225	
226	NOTE. For the imaging, the wavelengths can be adapted according to the dyes associated with
220	secondary antibodies and probes
227	secondary antibodies and probes.
220	2 E. Save the confecal images as TIEE files
229	S.S. Save the comocal images as this mes.
230	A Dest imaging analysis
231	4. Post iniaging analysis
252	4.1. To get started lownsh imagel and pavigate to the <b>Divging</b> many. From there, colort <b>Magne</b>
233	4.1. To get started, launch images and havigate to the <b>Plugins</b> menu. From there, select <b>Macros</b>
234	<b>Edit</b> to open the macros source code.
235	
236	NOTE: This will enable adjustment of the parameters as needed, and if required, modify the
237	folder directory for each channel within the macro's source code.
238	
239	4.2. To follow this protocol and quantify <i>Mef2</i> smFISH spots in larvae, use the following settings:
240	Log_radius = 3.0 and Log_quality = 20.0. Segement Mef2-positive nuclei with blur = 2;
241	nucleus_scale_parameter = 30; nucleus_threshold = -8.0; nucleus_size = 300.
242	
243	4.3. In muscle fibers, use the following parameters: Log_radius = 2.5; Log_quality = 60.0; blur =
244	2; nucleus_scale_parameter = 100; nucleus_threshold = 0.0; nucleus_size = 300.
245	
246	4.4. To initiate the macro, simply execute the <b>run</b> command and wait for the macro to
247	automatically load all images from the designated folder and quantify them in a sequential
248	manner.
249	
250	NOTE: Please note that the duration of this process may vary from a few seconds to several hours,
251	depending on the volume of data being analyzed.
252	
253	4.5. Look for the results displayed in two new .csy files: file FISH results and file nuclei results.
254	The first shows the intensities of FISH spots and distances with their nearest nucleus. The second
255	one gives the number of EISH-detected spots inside and outside the nucleus. It takes into account
255	one gives the number of hor detected spots inside and outside the nucleus. It takes into decount
250	the distances between the detected snots and the segmented nucleus
15/	the distances between the detected spots and the segmented nucleus.
257	the distances between the detected spots and the segmented nucleus.
257 258 250	the distances between the detected spots and the segmented nucleus. <b>REPRESENTATIVE RESULTS:</b>
257 258 259	the distances between the detected spots and the segmented nucleus. <b>REPRESENTATIVE RESULTS:</b> In this protocol, we used commercially produced probes targeting <i>Mef2</i> and <i>zfh1</i> mRNA. We designed the probes with the manufacturer's FIGH maskes designer (Table of Patiential). The
257 258 259 260	the distances between the detected spots and the segmented nucleus. <b>REPRESENTATIVE RESULTS:</b> In this protocol, we used commercially produced probes targeting <i>Mef2</i> and <i>zfh1</i> mRNA. We designed the probes with the manufacturer's FISH probes designer ( <b>Table of Materials</b> ). The
257 258 259 260 261	the distances between the detected spots and the segmented nucleus. <b>REPRESENTATIVE RESULTS:</b> In this protocol, we used commercially produced probes targeting <i>Mef2</i> and <i>zfh1</i> mRNA. We designed the probes with the manufacturer's FISH probes designer ( <b>Table of Materials</b> ). The <i>Mef2</i> set targets the exons common to all <i>Mef2</i> RNA isoforms (from exon 3 to exon 10). The <i>zfh1</i>
257 258 259 260 261 262	the distances between the detected spots and the segmented nucleus. <b>REPRESENTATIVE RESULTS:</b> In this protocol, we used commercially produced probes targeting <i>Mef2</i> and <i>zfh1</i> mRNA. We designed the probes with the manufacturer's FISH probes designer ( <b>Table of Materials</b> ). The <i>Mef2</i> set targets the exons common to all <i>Mef2</i> RNA isoforms (from exon 3 to exon 10). The <i>zfh1</i> set targets the third exon common to both <i>zfh1-RB</i> and <i>zfh1-RE</i> isoforms <sup>10</sup> . Both probes are
257 258 259 260 261 262 263	the distances between the detected spots and the segmented nucleus. <b>REPRESENTATIVE RESULTS:</b> In this protocol, we used commercially produced probes targeting <i>Mef2</i> and <i>zfh1</i> mRNA. We designed the probes with the manufacturer's FISH probes designer ( <b>Table of Materials</b> ). The <i>Mef2</i> set targets the exons common to all <i>Mef2</i> RNA isoforms (from exon 3 to exon 10). The <i>zfh1</i> set targets the third exon common to both <i>zfh1-RB</i> and <i>zfh1-RE</i> isoforms <sup>10</sup> . Both probes are conjugated to the Quasar 670 fluorescent dye.

#### 264

We generated a macro in-house that was compatible with ImageJ software to automatically analyze the raw .tif data. The macro allows the 3D segmentation of the nuclei by Cellpose <sup>26</sup> with the Fiji plugin BIOP and the transcriptional spot quantification by a Laplacian of Gaussian as implemented in the Trackmate plugin<sup>27</sup>. MorpholibJ plugin<sup>28</sup> is used for postprocessing (**Supplemental File 1**).

270

271 As a first step, we performed smFISH on wing imaginal discs, where *Mef2* and *zfh1* are known to 272 be expressed in the AMPs. These data show that both transcripts are uniformly detected in the 273 AMPs population co-stained with either Mef2 or Zfh1 antibodies (Figure 2A,B). Active 274 transcription sites (TS) can be revealed and distinguished from mature mRNA by smFISH as they 275 tend to be larger and have more intense signals than individual cytoplasmic transcripts or mature nuclear transcripts. Consistently, higher magnifications of AMPs distinguish between TS foci and 276 277 mature mRNA scattered in the cytoplasm (Figure 2A',B'), validating the sensitivity of this smFISH 278 protocol. It should, however, be mentioned that we observed one active TS per nucleus for both 279 *zfh1* and *Mef2* (Figure 2A',B'). These data further validate the accuracy of the probe design since 280 Mef2 and zfh1 transcriptional patterns in the AMPs colocalize with the detection of their 281 respective proteins (Figure 2).

282

In a second step, we examined the transcription site and distribution of *Mef2* and *zfh1* mRNAs in differentiated adult IFMs and associated stem cells (**Figure 2C,D**). These data clearly illustrate *Mef2* TSs in the syncytial muscle nuclei and *Mef2* mRNA distribution throughout the cytoplasm (**Figure 2C,C'**). Zfh1 specifically marks the MuSCs population <sup>10,11</sup>. Using this smFISH protocol, we successfully detected *zfh1* transcription in MuSCs marked by Zfh1-Gal4 > GFP expression. Higher magnification illustrates the detection of both *zfh1* TS and cytoplasmic single mRNAs (**Figure 2D'**).

Finally, we used our in-house-built ImageJ macro to quantitatively analyze *Mef2* transcriptional dynamics and spatial distribution. The computational pipeline can efficiently detect and segment *Mef2* spots and muscle nuclei (**Figure 3A–C**). However, due to signal-to-noise ratio issues, some spots may not be detected in certain cases (**Figure 3A',B'**). We used this automated method to investigate whether *Mef2* mRNA abundance varies between muscle precursor cells (AMPs) and differentiated adult IFMs. Our analysis shows that the number of *Mef2* mRNA per nuclei in adult IFMs and the AMPs are not significantly different (**Figure 3D**).

297

To gain insight into the spatial distribution of *Mef2* transcripts in adult IFM muscles, we quantified the fraction of cytoplasmic versus nuclear *Mef2* mRNA spots (**Figure 3E,F**). Using our program, we counted both the total number of *Mef2* mRNA spots and the number of *Mef2* mRNA spots overlapping with nuclei staining (Mef2). Subtracting the nuclei-associated mRNA from the total mRNA number provides the number of cytoplasmic mRNA spots. Our findings reveal that approximately 92% of *Mef2* mRNA is present in the cytoplasm, while the remaining 8% is associated with muscle nuclei, as shown in **Figure 3E,F**.

305

**FIGURE AND TABLE LEGENDS:** 

307

**Figure 1: Procedure to dissect and prepare wing discs and IFM samples for smFISH. (A)** L3 staged larvae. (**B**) Dissected and inverted anterior end of the larva. Arrows indicate the wing discs. (**C**) Isolated wing disc. (**D**) Dissected adult thoraxes. (**E**) Adult thorax oriented on a doublesided tape prior to bisection (dotted line). (**F**) Adult hemithorax. (**G**) Isolated IFMs. (**H**) Workflow of the smFISH protocol. Abbreviations: IFMs = indirect flight muscles; smFISH = single-molecule RNA fluorescence *in situ* hybridization; EtOH = ethanol.

314

315 Figure 2: Representative images of *Mef2* and *zfh1* smFISH in wing discs and adult IFMs. (A,B) 316 Mef2 and zfh1 mRNAs (purple) are detected uniformly in AMPs and co-localize with (A,A') Mef2 317 and (B,B') Zfh1 proteins (green), respectively. (A",B") Higher magnification of the AMPs. (C) Mef2 318 mRNAs and Mef2 protein are detected in adult IFMs. (C') Higher magnification of the boxed 319 region in C. (D) Zfh1-Gal4 >UAS-mCD8GFP expression (green) labels adult MuSCs. zfh1 RNA is 320 detected in adult MuSCs. (D') Higher magnification of the boxed region in D. Zfh1 transcription 321 start sites and single mRNAs are indicated by arrows and arrowheads, respectively. DAPI staining (blue). Scale bars = 50 μm (A,B,A',B'), 10 μm (A",B",C,D), 5 μm (C',D'). Abbreviations: IFMs = 322 323 indirect flight muscles; smFISH = single-molecule RNA fluorescence in situ hybridization; AMPs = 324 adult muscle progenitors; MuSCs = muscle stem cells; DAPI = 4',6-diamidino-2-phenylindole; GFP 325 = green fluorescent protein.

326

327 Figure 3: Quantification of *Mef2* transcription and cytoplasmic mRNA accumulation by smFISH.

328 (A) Distribution of *Mef2* transcripts (purple) and Mef2 protein (green) in wildtype adult IFMs. (A') 329 Higher magnification of boxed region in A. Transcription start sites and single mRNAs are 330 indicated by arrows and arrowheads, respectively. (B,C) Automatic Mef2 spots finding and 331 nuclear segmentation. Cytoplasmic and nuclear Mef2 mRNAs are indicated in green and 332 magenta, respectively. (B',C') Higher magnifications of boxed regions in B and C. Asterisks 333 indicate spots that are not counted by the macro. (D) Example of Mef2 spot quantification in 334 AMPs and adult IFMs, relative to the total number of nuclei. (p = 0,0785, Student's t-test, wing 335 discs (n = 4), IFMs (n = 11)). (E) Example of quantification of *Mef2* spot distribution in adult IFMs. (\*\*\* p< 0.0007, Student's t-test, n = 5). (F) Percentage of Mef2 mRNA spots detected inside and 336 337 outside the nuclei (n = 5). Scale bars =  $10 \mu m$  (A,A'). Abbreviations: IFMs = indirect flight muscles; 338 smFISH = single-molecule RNA fluorescence in situ hybridization; AMPs = adult muscle 339 progenitors.

340

#### 341 Supplemental File 1: MorpholibJ plugin used for postprocessing.

342

#### 343 **DISCUSSION:**

344

The application of the smFISH method has gained popularity in recent times, with its use extending to various cell types and model organisms. However, in the case of *Drosophila*, the majority of knowledge about adult muscle gene expression is derived from bulk molecular assays, which fail to provide quantitative information on the precise spatial location of mRNA molecules. To address this gap, we have optimized a method for performing smFISH on muscle precursor cells and adult *Drosophila* muscles. This approach has been adapted from a previously published protocol <sup>23</sup> and optimized for *Drosophila* muscle tissues. 352

The major obstacle in obtaining high-quality smFISH images is the thickness of adult muscles, which hinders the optimal penetration of probes. Therefore, it is crucial to separate the adult muscles from the other tissues of the animal and carry out the hybridization process with mild agitation. This particular step ensures effective permeabilization of the tissues.

357

Another critical aspect is that the signal-to-noise ratio can cause the computational pipeline to fail in detecting certain mRNA spots. We found that increasing the signal-to-noise ratio can be achieved by finding the optimal concentration of probes. The optimal dilution may differ depending on the composition of each probe set, including the conjugated dye and oligonucleotide composition. We recommend trying different dilutions; a final dilution of 200 nM yielded the best signal-to-noise ratio for adult tissues in these experiments.

364

365 With the smFISH method, it is possible to quantify the number of newly synthesized RNAs at the 366 TS foci. When a gene is actively transcribed, multiple nascent RNAs are produced concurrently at 367 the TS. As a result, the TS's intensity will surpass that of mature cytoplasmic RNA, and this feature, 368 combined with its nuclear localization, can differentiate the TS from individual cytoplasmic RNAs. 369 In the context of muscle biology, TS detection and quantification are particularly important for 370 determining the synchronization of transcription of a specific gene among muscle nuclei within 371 the same syncytium. However, this computational pipeline is not designed to differentiate 372 between cytoplasmic mRNA and TS signals. As an alternative, we suggest combining this smFISH 373 method with other well-established smFISH analysis tools, such as BayFISH or FISH-quant <sup>29,30</sup>. 374 These tools have been proven to facilitate automated segmentation and fluorescence intensity 375 calculations of RNA aggregates with remarkable precision.

376

Finally, while smFISH detects mRNA molecules with high spatial resolution, it is limited to analyzing a small number of mRNA at the same time. Multiscale methods like merFISH (multiplexed error-robust fluorescence *in situ* hybridization) enable the simultaneous analysis of a large number of different mRNA <sup>31</sup>. By combining oligonucleotide probes and error-correcting codes, this method facilitates the detection of hundreds of RNA species in single cells.

382

#### 383 ACKNOWLEDGMENTS:

We thank Alain Vincent for his critical reading of the manuscript. We thank Ruth Lehmann for providing the Zfh1 antibody. We acknowledge the help of Stephanie Dutertre and Xavier Pinson for confocal microscopy at the MRic Imaging Platform. EL is supported by PhD fellowship from ministère de l'enseignement supérieur et de la recherche scientifique. NA is supported by the AFM-Telethon Ph.D. fellowship 23846. TP is funded by a Chan Zuckerberg Initiative DAF grant to T.P. (2019-198009). HB is supported by CNRS, the Atip Avenir program, and the AFM-Telethon trampoline grant 23108.

391

#### 392 **DISCLOSURES:**

- 393 The authors have no conflicts of interest to disclose.
- 394
- 395 **REFERENCES** :

396		
397	1	Relaix, F. et al. Perspectives on skeletal muscle stem cells. Nature Communications. 12 (1),
398		692 (2021).
399	2	Chal, J. & Pourquie, O. Making muscle: skeletal myogenesis in vivo and in vitro.
400		Development. <b>144</b> (12), 2104–2122 (2017).
401	3	Taylor, M. V. & Hughes, S. M. Mef2 and the skeletal muscle differentiation program.
402		Seminars in Cell and Developmental Biology. <b>72</b> , 33–44 (2017).
403	4	Siles, L. et al. ZEB1 imposes a temporary stage-dependent inhibition of muscle gene
404		expression and differentiation via CtBP-mediated transcriptional repression. Molecular
405		and Cellular Biology. <b>33</b> (7), 1368–1382 (2013).
406	5	Weitkunat, M. & Schnorrer, F. A guide to study Drosophila muscle biology. Methods. 68
407		(1), 2–14 (2014).
408	6	Boukhatmi, H. Drosophila, an integrative model to study the features of musclestem cells
409		in development and regeneration. Cells. 10 (8), 2112 (2021).
410	7	Bothe, I. & Baylies, M. K. Drosophila myogenesis. Current Biology. 26 (17), R786-791
411		(2016).
412	8	Catalani, E. et al. RACK1 is evolutionary conserved in satellite stem cell activation and
413		adult skeletal muscle regeneration. Cell Death Discovery. 8 (1), 459 (2022).
414	9	Siles, L., Ninfali, C., Cortes, M., Darling, D. S., Postigo, A. ZEB1 protects skeletal muscle
415		from damage and is required for its regeneration. Nature Communications. 10 (1), 1364
416		(2019).
417	10	Boukhatmi, H., Bray, S. A population of adult satellite-like cells in Drosophila is maintained
418		through a switch in RNA-isoforms. <i>Elife.</i> <b>7</b> , e35954 (2018).
419	11	Chaturvedi, D., Reichert, H., Gunage, R. D., VijayRaghavan, K. Identification and functional
420		characterization of muscle satellite cells in Drosophila. Elife. 6, e30107 (2017).
421	12	Gunage, R. D., Reichert, H., VijayRaghavan, K. Identification of a new stem cell population
422		that generates <i>Drosophila</i> flight muscles. <i>Elife</i> . <b>3</b> , e03126 (2014).
423	13	Figeac, N., Daczewska, M., Marcelle, C., Jagla, K. Muscle stem cells and model systems for
424		their investigation. Developmental Dynamics. 236 (12), 3332–3342 (2007).
425	14	Aradhya, R., Zmojdzian, M., Da Ponte, J. P., Jagla, K. Muscle niche-driven insulin-Notch-
426		Myc cascade reactivates dormant adult muscle precursors in Drosophila. Elife. 4, e08497
427		(2015).
428	15	Schonbauer, C. et al. Spalt mediates an evolutionarily conserved switch to fibrillar muscle
429		fate in insects. <i>Nature.</i> <b>479</b> (7373), 406–409 (2011).
430	16	Bryantsev, A. L., Baker, P. W., Lovato, T. L., Jaramillo, M. S., Cripps, R. M. Differential
431		requirements for Myocyte Enhancer Factor-2 during adult myogenesis in Drosophila.
432		Developmental Biology. <b>361</b> (2), 191–207 (2012).
433	17	Soler, C., Han, J., Taylor, M. V. The conserved transcription factor Mef2 has multiple roles
434		in adult Drosophila musculature formation. Development. 139 (7), 1270–1275 (2012).
435	18	Soler, C., Taylor, M. V. The Him gene inhibits the development of Drosophila flight muscles
436		during metamorphosis. Mechanisms of Development. <b>126</b> (7), 595–603 (2009).
437	19	Pinheiro, H. et al. mRNA distribution in skeletal muscle is associated with mRNA size.
438		Journal of Cell Science. <b>134</b> (14), jcs256388 (2021).
439	20	Denes, L. T., Kelley, C. P., Wang, E. T. Microtubule-based transport is essential to distribute

- 440 RNA and nascent protein in skeletal muscle. *Nature Communications*. **12** (1), 6079 (2021).
- 441 21 Li, H. et al. Fly Cell Atlas: A single-nucleus transcriptomic atlas of the adult fruit fly. *Science*.
  442 **375** (6584), eabk2432 (2022).
- 443 22 Santos, M. D. et al. Extraction and sequencing of single nuclei from murine skeletal 444 muscles. *STAR Protocols.* **2** (3), 100694 (2021).
- Raj, A., van den Bogaard, P., Rifkin, S. A., van Oudenaarden, A., Tyagi, S. Imaging individual
  mRNA molecules using multiple singly labeled probes. *Nature Methods.* 5 (10), 877–879
  (2008).
- Femino, A. M., Fay, F. S., Fogarty, K., Singer, R. H. Visualization of single RNA transcripts
  in situ. *Science.* 280 (5363), 585–590 (1998).
- Raj, A., Tyagi, S. Detection of individual endogenous RNA transcripts in situ using multiple
  singly labeled probes. *Methods in Enzymology.* 472, 365–386 (2010).
- 452 26 Stringer, C., Wang, T., Michaelos, M., Pachitariu, M. Cellpose: a generalist algorithm for 453 cellular segmentation. *Nature Methods.* **18** (1), 100–106 (2021).
- 45427Tinevez, J. Y. et al. TrackMate: An open and extensible platform for single-particle455tracking. *Methods.* **115**, 80–90 (2017).
- Legland, D., Arganda-Carreras, I., Andrey, P. MorphoLibJ: integrated library and plugins
  for mathematical morphology with ImageJ. *Bioinformatics.* 32 (22), 3532–3534 (2016).
- 458 29 Mueller, F. et al. FISH-quant: automatic counting of transcripts in 3D FISH images. *Nature*459 *Methods.* 10 (4), 277–278 (2013).
- Gomez-Schiavon, M., Chen, L. F., West, A. E., Buchler, N. E. BayFish: Bayesian inference of
  transcription dynamics from population snapshots of single-molecule RNA FISH in single
  cells. *Genome Biology.* 18 (1), 164 (2017).
- 463 31 Chen, K. H., Boettiger, A. N., Moffitt, J. R., Wang, S., Zhuang, X. RNA imaging. Spatially
  464 resolved, highly multiplexed RNA profiling in single cells. *Science*. **348** (6233), aaa6090
  465 (2015).
- 466



## Figure 1





Figure 3

Name of Material/ Equipment	Company	Catalog Number
Confocal microscope SP8	Leica	SP8 DMI 6000
DAPI	ThermoFisher	62248
Double-sided tape	Tesa	57910-00000
Dumont #55 Forceps	Fine Science Tools	11255-20
Formaldehyde 16%	Euromedex	EM-15710
Mef2 antibody	DHSB	NA
PBS 10x	Euromedex	ET330
RNaseZAP	Sigma-Aldrich	R2020
Stellaris probes	LGC	NA
Stellaris RNA FISH Hybridization Buffer	LGC	SMF-HB1-10
Stellaris RNA FISH Wash Buffer A	LGC	SMF-WA1-60
Stellaris RNA FISH Wash Buffer B	LGC	SMF-WB1-20
Swann-Morton Blades	Fisher scientific	0210
ThermoMixer	Eppendorf	5382000015
Triton X-100	Euromedex	2000
UAS-mCD8::GFP line	Bloomington	RRID:BDSC_5137
Ultrapure water	Sigma-Aldrich	95284
Watch Glass, Square, 1 5/8 in	Carolina	742300
Zfh1 antibody	Gifted by Ruth Leahman	
Zfh1-Gal4	Bloomington	BDSC: 49924, FLYB: FBtp0059625

Supplemental File 1

Click here to access/download Supplemental Coding Files Le Roux et al-Macro.txt