

Late Clinical and Echocardiographic Results with the Magna Ease© Pericardial Aortic Bioprosthesis

Amedeo Anselmi, Marie Aymami, Jacques Tomasi, Gemma D'alessandro, Thierry Langanay, Hervé Corbineau, Julien Mancini, Erwan Flécher, Jean-Philippe Verhoye

▶ To cite this version:

Amedeo Anselmi, Marie Aymami, Jacques Tomasi, Gemma D'alessandro, Thierry Langanay, et al.. Late Clinical and Echocardiographic Results with the Magna Ease© Pericardial Aortic Bioprosthesis. European Journal of Cardio-Thoracic Surgery, 2023, European Journal of Cardio-Thoracic Surgery: Official Journal of the European, 10.1093/ejcts/ezad351. hal-04338096

HAL Id: hal-04338096 https://hal.science/hal-04338096

Submitted on 29 Jan2024

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

1

Anselmi A.

1 2	Late Clinical and Echocardiographic Results with the Magna Ease© Pericardial Aortic Bioprosthesis
3	
4 5 6 7	Amedeo ANSELMI ^{1,2} , MD PhD, Marie AYMAMI, MD PhD, Jacques TOMASI, MD PhD, Gemma D'ALESSANDRO ¹ , MD, Thierry LANGANAY, MD, Hervé CORBINEAU, MD, Julien MANCINI ³ , MD PhD, Erwan FLECHER, MD PhD, Jean-Philippe VERHOYE ^{1,2} , MD PhD
9 10	¹ Division of Thoracic and Cardiovascular Surgery, Pontchaillou University Hospital, Rennes, France
11	² Univ Rennes, CHU Rennes, Inserm, LTSI - UMR 1099, F-35000 Rennes
12 13	³ Aix-Marseille Univ, Inserm, IRD, ISSPAM, APHM, Biostatistics Dept., UMR1252 SESSTIM Research Unit, Marseille, France
14 15 16 17 18 19 20 21 22 23 24	Address for correspondence: Amedeo Anselmi, MD, PhD Division of Thoracic and Cardiovascular Surgery Pontchaillou University Hospital – Rennes, France E-mail <u>amedeo.anselmi@chu-rennes.fr</u> Telephone: 00 33 2 99 28 24 16
25	Meeting presentation: this content will be presented at the 2022 EACTS Annual
26	Meeting to be held in Milan, Italy, from October 5 th to October 8 th , 2022.
27	WORD COUNT:
28	KEY WORDS: Surgical Aortic Valve Replacement; Bioprostheses; Outcomes
29	Visual Abstract.
30 31	Key Question: Which are the long-term outcomes with the Perimount Magna $\mbox{Ease}\xspace$ bioprosthesis for SAVR?
~ ~	Kan Finding and A 040 OAV/Da after medicar 0.0 second aligned extension and anti-

- Key Findings: In 1,016 SAVRs after median 9.8 years, clinical outcomes were optimal and freedom from SVD (10-yrs) was 96.5% and 97.4% (Kaplan-Meier and competing
- 34 risks).
- Take-Home message: The Magna Ease© bioprosthesis can be safely employed for
- 36 SAVR in patients with life expectancy of at least 10 years. Continued investigation is
- 37 required for further follow-up timepoints.

38

39

40 **Abstract.**

Objectives. The population of candidates to surgical aortic valve replacement (SAVR)
 is evolving. The Perimount Magna Ease© bioprosthesis has been introduced relatively
 recently in the practice. We aimed at evaluating its long-term results.

Methods. Single-center cohort of 1,016 consecutive SAVRs with the Magna Ease©
valve (2008-2014), all-comers population. Prospective collection of in-hospital data,
systematic clinical and echocardiographic follow-up. Evaluation of valve-related
events: SVD, (structural valve deterioration, according to modified definition criteria),
NSVD (nonstructural valve dysfunction), PPM (patient-prosthesis mismatch).

Results. Age at SAVR was 73.4±9.5 years; calcified aortic stenosis was the indication 49 to surgery in 59.6%. 974 patients entered the follow-up, 564 were alive at last follow; 50 up (median duration: 9.8 years) (up to 13.4 years). NYHA class was I or II in 92.1%. 51 Overall survival at 10 years was 56.8%±1.8. Freedom from SVD at 10 was 96.5%±0.8 52 53 (Kaplan-Meier), and 97.4%±0.6 (competing risks) (28 SVD events after 6.9 years ± 3.3). There were 15 reinterventions for SVD (redo-SAVR and TAVI); 10-year freedom 54 from reintervention was 97.8% ± 0.6. Moderate and severe PPM occurred in 26.8% 55 and 5.4%, respectively; without association with late mortality (p=0.12 for moderate 56 and p=0.70 for severe PPM). Freedom from valve-related mortality was 97.8%±0.5 at 57 10 years. 58

Conclusions. In this follow-up of the Magna Ease bioprosthesis for SAVR, data indicate good late outcomes (30-days outcomes are excluded). Continued follow-up is required to further support its use in patients with life expectancy greater than 10-12 years.

63

- 64
- 65
- 66

67

68 Introduction.

Surgical Aortic Valve Replacement (SAVR) remains the treatment of choice for 69 severe aortic stenosis in patients younger than 75 and presenting acceptable operative 70 risk, and for operable patients affected by severe aortic regurgitation [1]. While older 71 72 patients are preferentially treated through transcatheter aortic valve replacement (TAVR), the population of candidates to SAVR is being characterized by younger age. 73 Additionally, societal evolutions and the importance attributed to lifestyle trigger 74 frequent demands from informed patients younger than 60 to receive a bioprosthesis 75 when SAVR is required [2], in order to avoid the limitations associated with 76 anticoagulant therapy after implantation of a mechanical valve. Therefore, the 77 appreciation of long-term durability of bioprostheses acquires a growing importance. 78 The Perimount Magna Ease[©] bioprosthesis is widespread in the clinical practice. It 79 80 has been the object of investigations focusing on its mid-term durability and clinical outcomes [3] even in comparison to other pericardial bioprostheses. Nonetheless, its 81 introduction following its predecessor Perimount Magna© valve is relatively recent. In 82 83 this context, the aim of the current study is to evaluate the long-term outcomes of SAVR with this valve device under the profiles of overall survival, freedom from SVD and from 84 other valve-related outcomes. 85

86

87

88

89

90

91 **Patients and Methods.**

92 <u>Ethical Statement.</u>

The institutional database, ensuring consecutiveness and completeness of the cohort, is declared in the CNIL (*Commission Nationale des Informatiques et Libertés*) under the number 1207754, according to compliance regulations and to the French law. Patients signed an informed consent for the employment of personal data for research purposes, and all data were analyzed anonymously.

98

99 <u>Criteria of inclusion and management of data.</u>

We performed a retrospective recall of data about 1,017 SAVR procedures 100 performed with the Perimount Magna Ease[©] bovine pericardial bioprosthesis in 1,016 101 patients. We considered only the first procedure for the single patient with twice 102 procedures in the database. Although this device is implanted in our department even 103 during current years, we enrolled patients undergoing SAVR between January 2nd, 104 2008 and December 31st, 2014 in order to focus on long-term follow-up. Combined or 105 reoperative procedures did not represent an exclusion criterion. The electronic 106 database of Rennes University Hospital, entailing prospective collection of in-hospital 107 data of all patients undergoing cardiac surgery at our department since 1978, has been 108 109 used as the source of all pre, intra- and early postoperative variables. Surgical methods for SAVR were homogeneous across the study period; bioprostheses were implanted 110 using a supra-annular technique with a series of non-everting U-shaped pledgetted or 111 112 non-pledgetted sutures.

Systematic follow-up of the included patients was performed in 2020 and 2021. At 113 first, we interrogated the INSEE (Institut National de la Statistique et des Etudes 114 Economiques – National Institute of Statistics and Economical Studies) database to 115 obtain the patients' vital status and the date of death. Following collection of the 116 available in-house follow-up data, questionnaires were sent by surface mail to general 117 practitioners and/or cardiologists. Such questionnaires enquired about the results of 118 the last available echocardiography, NYHA class and causes of death, if applicable. 119 Questionnaires inquired about any adverse event occurring during the follow-up; both 120 death and adverse events were categorized as valve- or non-valve related. 121 122 Controversies were adjudicated by consensus. Any alteration related to the bioprosthesis as assessed by echocardiography, clinical examination or reoperation 123 was recorded. Echocardiography date were interpreted according to standard criteria. 124 When the questionnaire was not informative or when we had no response, the 125 cardiologists or referring physicians were directly contacted, or even the 126 patients/families. For individual patients to be considered not lost at follow-up, 127 echocardiography results had to be available; the date of the last available 128 echocardiography was considered as the date of last follow-up. 129

130

131 Outcomes and definitions.

Definitions of valve-related events were formulated with the basis of the 2008 Akins guidelines and of more recent consensus documents as well [4, 5, 6]. SVD was defined as a persistent dysfunction or deterioration of the bioprosthesis evident at echocardiography and/or reoperation, but exclusive of infection or thrombosis. SVD could be under the form of stenosis, with significant and persistent increase in average gradient (>10 mmHg increase in average gradient vs. reference echocardiography,

unless attributable to other transient causes such as anemia or abnormally increased 138 cardiac output), and/or decrease in effective orifice area (EOA) compared to reference 139 echocardiography. SVD could otherwise be under the form of intraprosthetic leak ≥ 140 grade 2/4. Morphologic valve deterioration was therefore a component of the definition 141 of SVD ("Bioprosthetic Valve Dysfunction" in a previous consensus paper [7]). Recent 142 granular definitions have been proposed for SVD. These comprise the evidence of 143 average gradient at any timepoint >20 mmHg to define SVD [6, 8, 9]. In the view of the 144 Authors, and consistently with a large body of previous literature, elevated 145 transvalvular gradient evident at hospital discharge (reference echocardiography) in 146 147 absence of bioprosthetic dysfunction and corresponding to PPM (patient-prosthesis mismatch), persisting during the follow-up with no significant change in the average 148 transvalvular gradients the EOA, and with no evidence 149 or at echocardiography/reoperation of evolutive structural dysfunction, does not constitute 150 an SVD event. These events were computed as PPM. Additionally, we distinguished 151 between SVD and reintervention for SVD: such outcome measures are presented 152 separately. Reintervention for SVD included both redo aortic valve surgery and 153 transcatheter Valve-in-Valve procedure. NSVD (Non-Structural Valve Deterioration) is 154 any dysfunction in valve function not linked to a structural failure, including 155 periprosthetic leak, thrombosis, hemolysis, pannus. PPM was categorized as severe 156 when the indexed Effective Orifice Area (iEOA) was ≤0.65 cm²/m², moderate when 157 iEOA was ≤0.85 cm²/m², or absent when iEOA was >0.85 cm²/m² [4, 10]. Infectious 158 endocarditis (IE) was any infection involving the bioprosthesis. Early post-operative 159 mortality was death within 30 days following surgery or later if during the index 160 hospitalization. We defined valve-related mortality as due to SVD, NSVD, IE or 161 reintervention on the index bioprosthesis leading to death, or as due to embolism, 162

hemorrhage and cardiac insufficiency as well, unless the cause-effect relationship with
 any failure of the index bioprosthesis could be excluded.

165

166 <u>Endpoints</u>

167 The primary endpoint was the freedom from SVD at follow-up, using both actuarial 168 (Kaplan-Meier) and actual (competing risks) methods.

The secondary endpoints were the overall follow-up survival, the freedom for reintervention for SVD, the rate of PPM and its impact on late survival, the freedom from NSVD, the freedom from IE.

172

173 Statistical analysis

Continuous variables were described as mean ± standard deviation (SD) in case 174 of normal distribution, or as median in case of non-normal distribution. The 175 Kolmogoroff-Smirnoff test was used to assess normality. Categorical variables were 176 177 expressed as percentages. Time-to-event analyses for non-lethal valve-related events were built employing both the Kaplan-Meier ("actuarial") and competing risks ("actual") 178 methodologies in complementary manner [5, 11]. Non-lethal events at follow-up (such 179 as SVD) were also analyzed with the 'actual' method (competing risks). In fact, the 180 description of SVD using the "actuarial" (Kaplan-Meier) method alone would provide 181 excessively pessimistic estimations for the purposes of counseling individual patients 182 concerning the expected risk of SVD, since death due to other causes censors the 183 occurrence of SVD itself. On the contrary, the Kaplan-Meier method is appropriate to 184 compare outcomes among different models of bioprostheses. In the graphical 185

representation of Kaplan-Meier curves, confidence intervals were shown as dotted 186 lines corresponding to each curve. For the purposes of the analysis of late survival, 187 early mortality was excluded from Kaplan-Meier curves. The median follow-up duration 188 was also estimated using reversed Kaplan-Meier method. To study the impact of PPM 189 on late survival, a Fine and Gray model was used. We performed stratified analysis by 190 size of implanted prosthesis. Intergroup comparison was conducted with ANOVA for 191 continuous data and chi-square test for categorical data. Opposed Kaplan-Meier 192 curves were compared with the log-rank statistic. The alpha level was set at 0.05. The 193 analyses were conducted using IBM SPSS Statistics 27.0 (IBM Inc., New York, USA) 194 195 and survminer and cmprsk packages for R 4.0 (The R Foundation for Statistical Computing Platform, Vienna, Austria), within an external, independent and academic 196 Biostatistics Department (JM). 197

198

199

200

201

202

203

204

205 **Results.**

- The dataset comprised 1,016 patients receiving SAVR with the Magna Ease© valve;
- 207 Figure 1 displays the study workflow. The baseline characteristics are summarized in

Table 1; the population is mainly represented by patients affected by calcified aortic 208 stenosis, while isolated aortic regurgitation or infectious endocarditis were markedly 209 less frequent. Table 2 reports the early in-hospital findings. There was a considerable 210 rate of concomitant procedures to SAVR (34.2%) (concomitant coronary bypass in 211 23.1%, other valve surgery in 5.7%, ascending aortic replacement in 47 pts 4.9%), 212 while redo surgery occurred in 3.4% of instances. Operative mortality was 4.1% with 213 214 42 deaths. No early death was valve-related. The rate of implantation of a new permanent pacemaker was 3% in this cohort including patients affected by infectious 215 endocarditis and undergoing reoperative valvular operations. 216

Stratified analysis by prosthesis size suggested that the recipients of smaller sized valves tended to be older, to be female, to have smaller BSA and higher baseline LVEF, to present more frequently calcified aortic stenosis and less frequently bicuspid aortic valve (Supplementary Table 1). At discharge, recipients of smaller-sized valves presented higher average and peak aortic gradients, and lower EOA. They also displayed a higher rate of new permanent pacemaker implantation (Supplementary Table 2).

224

225 Follow-up Results

Of the 974 patients discharged alive from the hospital, none was lost at follow-up (100% completeness); median follow-up duration was 9.8 years according to reversed Kaplan-Meier method (up to 13.4 years). During the follow-up period, we noted 410 deaths (42.1%); therefore, 564 patients were alive at the end of the follow-up. The functional result was NYHA class I or II in 92.1% of patients; in the entire population. Overall survival (Table 4 and Figure 2B) was $83.2\% \pm 1.2$ at 5 years and $56.8\% \pm 1.8$

at 10 years. At last echocardiography, mean average transvalvular gradient in the 232 entire population was 14.4 mmHg ± 7.5, and average EOA was 1.46 cm² ± 0.42; aortic 233 regurgitation was noted to be of grade 3 or 4 in 0.4% of cases. Recipients of smaller-234 sized prostheses presented higher average and peak gradients, and lower EOA than 235 recipients of larger-sized valves (Supplementary Table 3). At follow-up, we observed 236 no statistically significant differences in valve-related events by size of implanted 237 prosthesis (Supplementary Table 4). We observed a total of 28 SVD events during the 238 follow-up, occurring after 6.9 years ± 3.3 (longest interval: 12.7 years); freedom from 239 SVD was 96.3% \pm 0.8 (Kaplan-Meier) and 97.2% \pm 0.6 (competing risks) at 10 years 240 (Table 4 and Figure 2A). Of the 28 SVD events, 15 were treated with reintervention 241 (10 redo-SAVR and 5 Valve-in-Valve TAVR); the remainders were not treated due 242 either to moderate SVD with limited clinical impact, or to severely compromised clinical 243 conditions in very aged patients. The average age of redo-SAVR patients was 244 50.1±18.2, and 66.8±8.7 for valve-in-valve TAVR patients (p=0.077). Freedom from 245 reintervention (any) for SVD was 97.8% ± 0.6 at 10 years (Kaplan-Meier) (Table 4). 246 Freedom rates from other valve-related events (NSVD, IE) are reported in Table 4; of 247 note, there were a total of 20 late valve-related late deaths (10-years freedom: 97.8% 248 ± 0.58) (Figure 2C). There were no significant intergroup differences at stratified 249 Kaplan-Meier analysis by size of implanted prosthesis, in terms of freedom from late 250 valve-related events (Supplementary Figure 1). At last echocardiography, data about 251 EOA allowing estimation of PPM were available in 88& patients (90.5% of patients 252 entering the follow-up). Among these, no PPM was noted in 567 (85.2%), moderate 253 PPM in 261 (26.8%), and severe PPM in 53 (5.4%). In Fine & Gray model, PPM was 254 not identified as a significant predictor of late mortality (p=0.12 for moderate PPM; 255 p=0.70 for severe PPM). From discharge echocardiography to last follow-up, we noted 256

10

257	no statistically significant difference in average aortic gradient, but a statistically
258	meaningful decrease in peak aortic gradient and in EOA (Supplementary Figure 2).
259	
260	
261	
262	
263	
264	
265	
266	
267	
207	
268	
269	
270	
271	
272	Discussion.
273	Surgical aortic valve replacement (SAVR) is still the treatment of choice for severe
274	aortic valve disease, in younger patients (aged <75 according to the recent
275	ESC/EACTS recommendations) [1] with low and intermediate operative risk.

a diffusing tendency towards the employment of bioprostheses for SAVR over 277 mechanical valves in patients younger than 60 years after appropriate patients' 278 information. Given the employment of bioprostheses in patients with important life 279 expectancy, often beyond 15 or 20 years, the appreciation of very late reliability and 280 valve-related outcomes acquires greater significance in the current era. The device 281 investigated herein is one of the last-generation stented pericardial bioprostheses with 282 inner-mounted pericardial leaflets. It can be distinguished from its predecessor 283 (Magna© valve) by a lower profile, sleek commissure posts and a low stent base, 284 aimed at facilitating implantation namely in smaller root anatomies (Figure 3). It entails 285 286 a flexible cobalt-chromium alloy stent, designed to absorb energy changes during the cardiac cycle and reduce leaflet stress. It mounts three independent bovine pericardial 287 leaflets. Given its relatively recent introduction in clinical practice (2007) and its 288 widespread diffusion, it is important to evaluate its late performance in the current 289 context. In previous investigations, we have showed the optimal hemodynamic 290 performance of the Magna Ease bioprostheses in small aortic annuli (19 to 23 mm 291 nominal size valves) [2], which is a distinguished feature over porcine valve devices. 292 Nonetheless, the aforementioned investigation was limited by a 3.7 ± 2 years average 293 follow-up. Previous studies with large sample size have been so far performed to 294 evaluated mid-term durability (ex. Thorp and associates, reporting a 0.3% rate of 295 reoperation for SVD in 1,126 recipients after a mean follow-up time of 5 ± 2 years) [12]. 296 Such investigation was further limited, in the Authors' view, by the use of VARC-2 297 criteria for definition of SVD (including average transprosthetic gradient >20 mmHg 298 irrespective of the cause), leading to an excessively pessimistic 28.7% SVD rate in this 299 mid-term study [9]. An important series [3] comprised 1,365 recipients of Magna Ease© 300 valve, but presented a 4.3 years only average follow-up time. Additionally, this paper 301

presented only freedom from SVD requiring reintervention, which likely underestimates 302 the true SVD burden. In the view of the Authors, elevated transvalvular gradient should 303 be adjudicated between PPM, valve thrombosis (or other potential causes such as 304 subvalvular pannus) vs. SVD, based on transthoracic and transesophageal 305 echocardiography evidence, morphology and movement of leaflets, evolution of 306 gradients / EOA over time. Thereof, we herein considered as SVD only cases with 307 valve dysfunction due to inherent structural failure (such as, but not limited to, fracture, 308 poppet escape, calcification, leaflet tear, stent creep, and suture line disruption inside 309 the prosthetic valve). The SVD definition herein corresponds to the "Structural 310 311 Bioprosthetic Valve Dysfunction" in the 2020 consensus paper [7], which confirms the importance to distinguish between SVD and PPM, and between echocardiographically 312 evident SVD and its clinical consequences (with or without reintervention). Our 313 approach is, therefore, essentially consistent with the consensus document. We 314 described event-freedom rates which actually correspond to the outcome required for 315 appropriate patient information at the time of preoperative counseling. More recently, 316 other studies have addressed the 10-year durability of the Magna Ease© valve (mean 317 6.6 ± 2.6 years follow-up) [13], but was based on a limited cohort of 338 patients. With 318 1,016 included SAVR cases, our study features one of the largest cohorts and one of 319 the longest follow-up (average 8,9 years among survivors). We report good late 320 durability and reliability in terms of freedom from valve-related adverse events with the 321 investigated device, supporting its reasonable employment in patients whose life 322 expectancy is 10 to 12 years. The current study, yet one of the largest in terms of 323 cohort and follow-up duration, cannot provide robust indications about longer follow-324 up (i.e. with respect to patients with longer life expectancy), since after the 10th follow-325 up year less than one fourth of the initial cohort is still at risk. Continued investigation 326

is required to this purpose. To this respect, we underline the originality of the current
series which provides a longer follow-up than previous studies with a comprehensive
assessment of valve-related events. This rationally supports clinical decision-making
on the basis of demonstrated durability at longer timepoints.

Structural valve deterioration is a relatively infrequent event often occurring after 331 several years after implantation. Therefore, sufficient sample size and follow-up 332 duration are required for appropriate appreciation [14]. Also, the features, timing of 333 appearance and modalities of clinical presentation of SVD can be often associated with 334 individual models of bioprostheses. In example, noncalcified leaflet tear potentially 335 occurring 3 to 5 years after implantation with abrupt clinical manifestation has been 336 associated to some devices with externally mounted pericardium [4, 15], and 337 accelerated calcification to others [16]. In the current study, the form of SVD was more 338 frequently represented by progressive leaflet calcification leading to stenosis. Porcine 339 bioprostheses feature lesser hemodynamic performance than pericardial valves 340 (notably in smaller annuli) but SVD seems to occur more frequently under the form of 341 progressive, indolent calcification and to show longer delay until occurrence (ex. 9.3 ± 342 4.9 years after implantation in a third-generation device followed up during an average 343 8.5 ± 3.9 years period [17]). Similar considerations can be done for other porcine 344 xenografts [18]. Such average time to appearance of SVD is greater that observed in 345 the current study (6.9 years \pm 3.3); this point deserves particular attention and possibly 346 comparative investigations for the very long-term follow-up (15 to 20 years). Therefore, 347 appropriate evaluation of the burden of SVD and of its clinical consequences need to 348 be considered when discussing the late effectiveness of SAVR as a therapy vs. 349 transcatheter techniques (TAVR). Similarly, such comparative evaluation cannot be 350

14

based only on the rate of freedom from valve-related reintervention, which excludes
SVD cases not undergoing reoperation or valve-in-valve TAVI [19].

The clinical significance of PPM remains object of debate. Some Authors have 353 354 associated PPM with increased overall and cardiac mortality, mainly in patients with baseline impaired left ventricular ejection fraction [20, 21]. The appreciation of PPM 355 depends on individual model of bioprostheses [22]. In the present large cohort, we 356 failed to identify a significant association between PPM and overall late mortality (Fine-357 Gray model), although no subgroup-based analysis was performed. Given the known 358 very good hemodynamic characteristics of the investigated bioprosthesis [2], in our 359 center we tended to its frequent employment in patients with small annuli (19-21 mm), 360 thus leading to a potential overestimation of the rate of PPM in the overall population 361 herein (28.4% and 5.9%, respectively for moderate and severe PPM). Stratified 362 analysis described that the distribution of other valve-related events than PPM across 363 subgroups by prosthesis size. Our series is homogeneous in technical aspects 364 potentially impacting the hemodynamic results (all valves implanted supra-annular). 365

Age at surgery influences the lifetime likelihood of incurring in SVD due to longer life 366 expectancy. Herein, we observed a significantly higher risk of SVD or reoperation for 367 SVD among patients younger than 65 at surgery (Supplementary Table 5). This further 368 underlines the usefulness of the competing risks method when addressing non-lethal 369 valve-related events. Additionally, we observed a significant association at Kaplan-370 Meier analysis between mean aortic gradient >20 mmHg at discharge and occurrence 371 of SVD (Supplementary Figure 3). This can be interpreted as a consequence of smaller 372 'reserve' in EOA predisposing to evident SVD in case of evolution towards leaflet 373 calcification. 374

The current investigation is limited by its single-center nature. Although the included 375 patients were identified retrospectively, their in-hospital clinical data had been collected 376 prospectively. Although no patient was lost at follow-up and control echocardiography 377 was available for all of them, causes of death at follow-up remained unknown for 17.8% 378 of cases, with a potential underestimation of the rate of valve-related death. We missed 379 EOA information for 9.5% of followed-up patients, which should be considered in the 380 interpretation of PPM results. The average age at SAVR in the current population (73.4 381 ± 9.5 years) is potentially higher and therefore less representative of the population of 382 patients undergoing SAVR in the current era. Nonetheless, good life expectancy in a 383 384 large proportion of our patients facilitated the performance of long-term investigations. In the interpretation of long-term results, the exclusion of immediate postoperative 385 events from Kaplan-Meier curves (as we focused herein on late valve-related events) 386 should be considered. Therefore, the current conclusions can only be applied to late 387 outcomes. 388

In conclusion, we herein report the long-term valve-related outcomes with the Perimount Magna Ease© bioprosthesis for SAVR. Durability of this device is good at 10 to 12 years, with a 2.9% rate of SVD. The hemodynamic characteristics are optimal, with low rates or PPM; rates of other valve-related adverse events (NSVD, IE) are limited. It is to be considered that early (30-days) events are excluded from this analysis. Continued evaluation (up to 15 and 20 years) is required for appropriate decision-making in the currently evolving population of candidates to SAVR.

396

397

398

399 Acknowledgments.

- 400 We thank patients, cardiologists and referring physicians who helped us to collect data
- 401 of interest. The Authors sincerely thank Mrs. Sylvie Marie for her work.

402

- 403 **Data Availability Statement.**
- 404 The current database's data are available from the corresponding author upon 405 reasonable request.

406

- 407 Authors' contribution statement.
- 408 Conceptualization: AA, JPV
- 409 Data curation: AA, MA, JT, GDA
- 410 Formal analysis: AA, GDA, JM
- 411 Funding acquisition: AA, JPV
- 412 Investigation: AA, MA, JT, GDA, JM, JPV
- 413 Methodology: AA, JM, JPV
- 414 Project administration: AA, JPV
- 415 Resources: GDA
- 416 Software: JM
- 417 Supervision: JPV, HC, TL, EF
- 418 Validation: MA, JT, EF
- 419 Visualization:
- 420 Writing Original draft: AA
- 421 Writing Review & Editing: AA, MA, JT, GDA, JPV, TL, HC, EF

422

- 423
- 424

Funding statement.

This work was partially supported by an unrestricted research grant from EdwardsLifesciences inc.

Conflict of Interest statement.

430 This work was partially supported by an unrestricted research grant from Edwards431 Lifesciences inc.

447 **References.**

448	1.	Vahanian A, Beyersdorf F, Praz F, Milojevic M, Baldus S, Bauersachs J et al.
449		ESC/EACTS Scientific Document Group. 2021 ESC/EACTS Guidelines for the
450		management of valvular heart disease. Eur J Cardio-Thorac Surg 2021;60:727-
451		800.
452	2.	Anselmi A, Ruggieri VG, Belhaj Soulami R, Flécher E, Langanay T, Corbineau
453		H et al. Hemodynamic Results and Mid-term Follow-up of 850 19 to 23 mm
454		Perimount Magna Ease Valves. Thorac Cardiovasc Surg. 2019;67:274-281.
455	3.	Biancari F, Valtola A, Juvonen T, Husso A, Dahlbacka S, Laakso T et al. Trifecta
456		Versus Perimount Magna Ease Aortic Valve Prostheses. Ann Thorac Surg.
457		2020;110:879-888.
458	4.	Anselmi A, Flecher E, Chabanne C, Ruggieri VG, Langanay T, Corbineau H et
459		al. Long-term follow-up of bioprosthetic aortic valve replacement in patients
460		aged ≤60 years. J Thorac Cardiovasc Surg. 2017;154:1534-1541.e4.
461	5.	Akins CW, Miller DC, Turina MI, Kouchoukos NT, Blackstone EH, Grunkemeier
462		GL et al. Guidelines for reporting mortality and morbidity after cardiac valve
463		interventions. J Thorac Cardiovasc Surg. 2008;135:732-8.
464	6.	Généreux P, Piazza N, Alu MC, Nazif T, Hahn RT, Pibarot P et al. Valve
465		Academic Research Consortium 3: Updated Endpoint Definitions for Aortic
466		Valve Clinical Research. J Am Coll Cardiol. 2021;77:2717-2746.
467	7.	Pibarot P, Herrmann HC, Wu C, Hahn RT, Otto CM, Abbas AE, et al.
468		Standardized Definitions for Bioprosthetic Valve Dysfunction Following Aortic or
469		Mitral Valve Replacement: JACC State-of-the-Art Review. J Am Coll Cardiol.
470		2022 Aug 2;80(5):545-561.

19

- 8. Dvir D, Bourguignon T, Otto CM, Hahn RT, Rosenhek R, Webb JG et al.
 Standardized Definition of Structural Valve Degeneration for Surgical and
 Transcatheter Bioprosthetic Aortic Valves. *Circulation. 2018;137:388-399.*
- 474 9. Kappetein AP, Head SJ, Généreux P, Piazza N, van Mieghem NM, Blackstone
 475 EH et al. Updated standardized endpoint definitions for transcatheter aortic
 476 valve implantation: the Valve Academic Research Consortium-2 consensus
 477 document. *Eur Heart J. 2012 ;33:2403-18.*
- 478 10. Blais C, Dumesnil JG, Baillot R, Simard S, Doyle D, Pibarot P. Impact of valve
 479 prosthesis-patient mismatch on short-term mortality after aortic valve
 480 replacement. *Circulation. 2003;108:983-8*.
- 11. Bourguignon T, Lhommet P, El Khoury R et al. Very long-term outcomes of the
 Carpentier-Edwards Perimount aortic valve in patients aged 50-65 years. *Eur J Cardiothorac Surg. 2016;49:1462-8.*
- 12. Thorp SD, Khazaal J, Yu G, Parker JL, Timek TA. Magna ease bioprosthetic
 aortic valve: mid-term haemodynamic outcomes in 1126 patients. *Interact Cardiovasc Thorac Surg. 2021;32:839-845*.
- 487 13. Kermen S, Strella J, Aupart A, Espitalier F, Aupart M, Bernard A et al. Durability
 488 of a bovine pericardial aortic bioprosthesis based on Valve Academic Research
 489 Consortium-3 echocardiographic criteria. *JTCVS Open. 2022;11:72-80*.
- 490 14. Rajab TK, Ali JM, Hernández-Sánchez J, Mackie J, Grimaudo V, Sinichino S et
 491 al. Mid-term follow-up after aortic valve replacement with the Carpentier
 492 Edwards Magna Ease prosthesis. *J Cardiothorac Surg. 2020;15:209*.
- 493 15. Werner P, Gritsch J, Scherzer S, Gross C, Russo M, Coti I et al. Structural valve
 494 deterioration after aortic valve replacement with the Trifecta valve. *Interact* 495 *Cardiovasc Thorac Surg.* 2021;32:39-46.

- 496 16. Sénage T, Le Tourneau T, Foucher Y, Pattier S, Cueff C, Michel M et al. Early
 497 structural valve deterioration of Mitroflow aortic bioprosthesis: mode, incidence,
 498 and impact on outcome in a large cohort of patients. *Circulation.*499 2014;130:2012-20.
- 17. Anselmi A, Flécher E, Ruggieri VG, Harmouche M, Langanay T, Corbineau H
 et al. Long-term results of the Medtronic Mosaic porcine bioprosthesis in the
 aortic position. *J Thorac Cardiovasc Surg. 2014;147:1884-91*.
- 18. Lehmann S, Merk DR, Etz CD, Oberbach A, Uhlemann M, Emrich F et al.
 Porcine xenograft for aortic, mitral and double valve replacement: long-term
 results of 2544 consecutive patients. *Eur J Cardiothorac Surg. 2016;49:1150-6*.
- 19. Persson M, Glaser N, Nilsson J, Friberg Ö, Franco-Cereceda A, Sartipy U.
 Comparison of Long-term Performance of Bioprosthetic Aortic Valves in
 Sweden From 2003 to 2018. *JAMA Netw Open. 2022;5:e220962*.
- 20. Chen J, Lin Y, Kang B, Wang Z. Indexed effective orifice area is a significant
 predictor of higher mid- and long-term mortality rates following aortic valve
 replacement in patients with prosthesis-patient mismatch. *Eur J Cardiothorac Surg. 2014;45:234-40.*
- 21. Mohty D, Dumesnil JG, Echahidi N, Mathieu P, Dagenais F, Voisine P et al.
 Impact of prosthesis-patient mismatch on long-term survival after aortic valve
 replacement: influence of age, obesity, and left ventricular dysfunction. *J Am Coll Cardiol. 2009;53:39-47*.
- 22. Durko AP, Pibarot P, Atluri P, Bapat V, Cameron DE, Casselman FPA et al.
 Essential information on surgical heart valve characteristics for optimal valve
 prosthesis selection: expert consensus document from the European
 Association for Cardio-Thoracic Surgery (EACTS)-The Society of Thoracic

521	Surgeons (STS)-American Association for Thoracic Surgery (AATS) Valve
522	Labelling Task Force Eur J Cardiothorac Surg. 2021;59:54-64.
523	
524	
525	
526	
527	
528	
529	
530	
531	
532	
533	
534	
535	
536	
537	
538	
539	
540	

Table 1. Baseline characteristics in the entire study population (N=1,016).

Characteristic	N (%) or mean ±
	SD
Age (years)	73.4±9.5
$BSA(m^2)$	1.80 ± 0.21
Creatinine (µmol/L)	97.2 ± 67.4
Male gender	577 (56.8)
NYHA class	
NYHA I	58 (5.7)
NYHA II	595 (58.7)
NYHA III	328 (32.4)
NYHA IV	32 (3.2)
Syncope	72 (7.1)
Previous stroke	35 (3.4)
Systemic hypertension	623 (61.3)
Smoker: past or present	203 (20)
Dyslipidemia	562 (55.3)
Diabetes mellitus	156 (15.4)
Renal failure	
- not requiring hemodialysis	45 (4.4)
- requiring hemodialysis	7 (0.7)
COPD	109 (10.7)
Cancer : past or present	86 (8.5)
LVEF(%)	60.5±11.2
LV end-Diastolic diameter (mm)	49.7±7.9
SPAP (mmHg)	36.9±11.2
Average Aortic Gradient (mmHg)	49.3±16.2
Aortic surface (cm ²)	0.74 ± 0.24
Aortic insufficiency	
AI grade III	40 (3.9)
AI grade IV	34 (3.3)
Etiology of aortic valve disease	
Calcified aortic stenosis	606 (59.6)
Bicuspid morphology	120 (11.8)
Prolapse/annular dilation	1 (0.1)
Endocarditis	40 (3.9)
Other	28 (2.9)

542

NYHA: New York Heart Association. BSA: Body Surface Area. COPD: Chronic
Obstructive Pulmonary Disease. LVEF: Left Ventricular Ejection Fraction. SPAP:
Systolic Pulmonary Artery Pressure. AI: Aortic Insufficiency.

546

547

23

548

549

Table 2. In-hospital results in the entire study population (N=1,016).

Characteristic	N (%) or
	mean ± SD
Timing	
Elective	966 (95.1)
Non-elective	50 (4.9)
CPB time (min)	75±36.7
Cross-clamp time (min)	59.9±30.3
Concomitant procedures	
Any	347 (34.2)
CABG	235 (23.1)
Other valve	58 (5.7)
Thoracic aorta	50 (4.9)
Other	25 (2.5)
Redo surgery	35 (3.4)
Prosthesis diameter (mm)	22.2±1.82
Early mortality	42 (4.1 <u>)</u>
Discharge echocardiography	
Mean Gradient (mmHg)	14±5
Max Gradient (mmHg)	25.7 ± 8.8
$EOA (mm^2)$	1.67±0.4
Aortic insufficiency grade I	111 (10.9)
Aortic insufficiency grade II	19 (1.9)
Aortic insufficiency grade III	2 (0.2)
Revision for bleeding	4 (0.4)
Gastrointestinal bleeding	5 (0.5)
Coma	
-transient	2 (0.2)
-permanent	4 (0.4)
<u>Stroke</u>	
-transient	5 (0.5)
-permanent	8 (0.8)
New Permanent Pacemaker	31 (3.0)
Renal failure	
- Not requiring hemodialysis	91 (9.0)
- Requiring hemodialysis	12 (1.2)

⁵⁵¹

552 CPB: Cardiopulmonary Bypass. CABG: Coronary Artery Bypass Grafting. EOA. 553 Effective Orifice Area.

554

Α.

555

556

Characteristic N (%) or mean ± SVD **NYHA class** 370 (46.7) NYHA I 361 (45.4) NYHA II NYHA III 60 (7.5) NYHA IV 3 (0.4) LVEF(%) 61.4 ± 9.7 LV end-Diastolic diameter (mm) 49.2±12.4 SPAP (mmHg) 34.3±9.9 Average Aortic Gradient (mmHg) 14.4 ± 7.5 Max Aortic Gradient (mmHg) 24.8±10.6 $EOA (mm^2)$ 1.46 ± 0.42 Aortic insufficiency AI grade I 95 (21.0) AI grade II 17 (3.8) AI grade III 1 (0.2) AI grade IV 1 (0.2)

557 **Table 3.** Follow-up results (N=974).

558

- NYHA: New York Heart Association. LVEF: Left Ventricular Ejection Fraction. SPAP:
 Systolic Pulmonary Artery Pressure. EOA: Effective Orifice Area. AI: Aortic
 Insufficiency.
- 562

563

564

565

566

567

Table 4. Follow-up results in time-to-event analyses (overall survival, freedom from
SVD and from any reintervention for SVD, freedom from NSVD and from IE) (N=974).
Actuarial (Kaplan-Meier) rates are provided; if applicable, freedom rates calculated
according to the competing risks method are provided in *italics*.

N (%)		
-	5 y	10 y
Deaths: 410	83±1.2	56.8 ± 1.8
(42.1)		
SVD :	99.3±0.3	96.5 ± 0.8
28 (2.9)	99.4±0.3	97.4±0.6
Reintervention for	99.4±0.3	97.8 ± 0.6
SVD:		
15 (1.5)		
NSVD :	99.2±0.3	99±0.4
10 (1)		
IE :	97.8±0.5	96.2±0.7
29 (3)		
Valve-related	98.6±0.4	97.8±0.5
death : 20 (2.1)		
N at risk	703	256
	195	230

573 SVD: Structural Valve Deterioration. NSVD: Nonstructural Valve Dysfunction. IE:

- 574 Infectious Endocarditis.

- **Table 5.** Follow-up valve-related events described as cumulative incidence (competing
- 584 risks).

Characteristic	5 y	10 y
SVD	0.6±0.3	2.6±0.6
Reintervention for SVD	0.5±0.2	1.6±0.4
NSVD	0.6±0.3	0.8±0.3
IE	2.1±0.5	3.2±0.6
Any reintervention	1.4±0.4	3.5±0.6
Valve-related death	1.3±0.4	1.9±0.5

596 **Figure Legends.**

597 **Central Image Legend.** Summary of SVD results: number of events, time to 598 occurrence, freedom rates according to both the Kaplan-Meier and competing risks 599 methods.

600 **Graphical Abstract legend**. Summary of the main study methods and findings 601 concerning valve durability.

Figure 1. Study workflow. SAVR: Surgical Aortic Valve Replacement. FU: Follow-up.

Figure 2. Incidence and Kaplan-Meier survival curves (operative mortality cases
excluded). A: Incidence of non-lethal events. B: Overall survival. C: Survival free from
any reintervention for SVD.

Figure 3. The Carpentier-Edwards Perimount Magna Ease[©] aortic bioprosthesis compared to its predecessor the Carpentier-Edwards Perimount Magna[©] aortic bioprosthesis.

609











Magna Ease valve (23 mm) (23 mm)

Click here to access/download Supplementary material Supplementary Table 1.docx

Click here to access/download Supplementary material Supplementary Table 2.docx

Click here to access/download Supplementary material Supplementary Table 3.docx

Click here to access/download Supplementary material Supplementary Table 4.docx

Click here to access/download Supplementary material Supplementary Table 5.docx

Supplementary Figure 1

Click here to access/download **Supplementary material** Supplementary Figure 1.tif

Supplementary Figure 2

Click here to access/download Supplementary material Supplementary Figure 2.tif

Supplementary Figure 3

Click here to access/download Supplementary material Supplementary Figure 3.tif