

232.038

# **Projet BePAT :**

BeaQuant, a digital autoradiograph for the production of radionuclides for nuclear medicine applications

### **Arthur Bongrand**

Melissa Azzoune, Samuel Duval, Arnaud Guertin Julie Champion

### Boulogne-Billancourt, 21 – 25 november







Jubatech



Uranium 238.029

Astatine

85

**Physics of Radiation InteractionS with Matter and Applications** 







# My career path



# **BeaQuant : real-time autoradiography system for "small" samples**

### Autoradiography

 $\rightarrow$  An **emission imaging technique** carried out from a radioactive source placed in contact with an emulsion, a photographic film or an adapted detector

 $\rightarrow$  Area of interest located at the surface ( $\approx$  10 µm deep)



**Glass slide** 





Lefeuvre et al., under review NIMA

Metal foil



# Samples

Applications

 $\rightarrow$  "Thin" (<1 cm thick)  $\approx$  size of a small sugar square

 $\rightarrow$  Pre-clinical (oncology + drug discovery)

Visualization of molecules or fragments of molecules

 $\rightarrow$  Can be electroplated on a glass slide, a "sugar" or a **metal foil (new)** 

 $\rightarrow$  Geology (study of uraniferous rocks, environmental pollution...)

Specificities of the BeaQuant

labeled with radioactive elements

- $\rightarrow$  Real-time direct counting
  - ✓ Quantification
  - ✓ Tracking



Technology	Direct gaseous counter	
Performance	<ul> <li>Spatial: Alpha 20 μm – Beta 20 μm – High energy β or β<sup>+</sup> 50 μm</li> <li>Linearity over 5 orders of magnitude</li> <li>Sensitivity: 0.0005 cpm/mm<sup>2</sup> (≈ mBq)</li> </ul>	
Markers	All radioisotopes	







# Radionuclide production scheme and expertise of PRISMA@Subatech



### **BePAT**

### **Properties of Alpha- and Beta-Emitting Radionuclides**



#### α therapy benefit

- Short range in human tissues : ≈ 100µm (≈ cell diameter)
- High Linear Energy Transfer (LET) : 60-200 keV/µm
- Great DNA damages  $\rightarrow$  no cell repair

#### <sup>230</sup>U/<sup>226</sup>Th from <sup>230</sup>Pa

considering Arronax's characteristics 17 and 30-70 MeV Proton beam



### **BePAT**



#### • α therapy benefit

- Short range in human tissues : ≈ 100µm (≈ cell diameter)
- High Linear Energy Transfer (LET) : 60-200 keV/µm
- Great DNA damages → no cell repair

#### <sup>230</sup>U/<sup>226</sup>Th from <sup>230</sup>Pa

considering Arronax's characteristics 17 and 30-70 MeV Proton beam



### **Application of BePAT**







Irradiation station and beam line







### Verification of production cross section



# Some rules for the preparation of the BePAT samples



Capsule and metal foils



Metal foil can be electrodeposited after digestion and chemical separation



### Samples

- $\rightarrow$  Metal foil
- → Metal foil electrodeposited after digestion and chemical separation (in the future)
- Sample activity
- → Insensitive to γ and X →  $A(\alpha) + A(\beta+) + A(\beta-) < 20 \text{ kBq}$  - Specific calculations required !
- Challenge
- $\rightarrow$  Short-lived isotope !



# **Experimental study conditions**

Cible	Source	
Th 232 (11,72g/cm <sup>3</sup> )	Protons	
6,25 cm <sup>2</sup> (2,5.2,5 cm <sup>2</sup> ) epaisseur 43,2 μm	2,5 cm d'extension latérale	

$$\boldsymbol{\sigma}(\mathbf{E}) = \frac{M. Y(E)}{Na \ e \ x \ \rho}$$
  
ou  
$$\boldsymbol{\sigma}(\mathbf{E}) = \frac{M. Act(tirr)}{Na \ e \ x \ \rho \ (1 - exp(-\lambda tirr)))}$$



Modèle physique  $\rightarrow$  INCL Cascade Intranucléaire de Liège tdchain.lst  $\rightarrow$  rendement [atm/sources] si >10<sup>-8</sup> sinon  $\rightarrow$ tdchain.act  $\rightarrow$  activité [Bq]

nucl	lide	yield	nuclide	yield
	 2	1 0005 09	V:: 91	1 0005 07
Ho	2	1.000E-08	Kr 85	2 100E-07
Co	70	1.000E-08	Kr 86	1.210E-06
Ni	71	1.000E-08	Kr 87	1.310E-06
Cu	68	1.000E-08	Kr 88	2.450E-06
Cu	70	2.000E-08	Kr 89	2.130E-06
Cu	71	2.000E-08	Kr 90	1.880E-06

nuclide		atoms	radioac	tivity
		[atoms/cc]	[Bq/cc]	[Bq]
Н	3	3.2542E+03	5.8015E-06	1.5664E-07
He	6	2.9263E+00	2.5144E+00	6.7890E-02
He	8	8.0768E-03	4.7006E-02	1.2692E-03
Li	8	5.2072E-02	4.2973E-02	1.1603E-03
Li	9	1.0599E-03	4.1203E-03	1.1125E-04
Be	7	1.9679E+01	2.9666E-06	8.0097E-08
Be	8	4.9126E-18	4.1621E-02	1.1238E-03
Be	10	1.3998E+01	2.0341E-13	5.4922E-15
Be	11	2.6077E-19	1.3089E-20	3.5339E-22





### **Experimental procedure**



### **Irradiation scenarios**



**AZZOUNE 2022** 

Contraintes  $\rightarrow A(\alpha)+A(\beta) \le 20kBq$  $\rightarrow I = 1 nA$  $\rightarrow tirr = 15 min$ 





### Localized $\alpha$ -spectrometry



MATERIAUX







# **Energy spectrum**

gamn	ia-ray spect	trum weig	nted by energ	y ( 42-groups)	
grp.	energy	[MeV]	flux	energy flux	relative
	upper	lower	[n/s/cc]	[MeV/s/cc]	error
1	50.000	30.000	0.0000E+00	0.0000E+00	0.0000E+00
2	30.000	20.000	0.0000E+00	0.0000E+00	0.0000E+00
3	20.000	14.000	0.0000E+00	0.0000E+00	0.0000E+00
4	14.000	12.000	0.0000E+00	0.0000E+00	0.0000E+00
5	12.000	10.000	3.4274E-08	3.7702E-07	3.9448E-02
6	10.000	8.000	2.5845E-08	2.3260E-07	3.9432E-02
7	8.000	7.500	2.6211E-08	2.0313E-07	3.9404E-02
8	7.500	7.000	2.4590E-08	1.7828E-07	3.9308E-02
9	7.000	6.500	9.7197E+04	6.5608E+05	6.6655E-02
10	6.500	6.000	1.0106E+06	6.3162E+06	4.0591E-02
11	6.000	5.500	3.3388E+06	1.9198E+07	3.4109E-02
12	5.500	5.000	6.5112E+06	3.4184E+07	2.5902E-02
13	5.000	4.500	7.2631E+06	3.4500E+07	1.9996E-02
14	4.500	4.000	1.3590E+07	5.7756E+07	2.0798E-02
15	4.000	3.500	2.2669E+07	8.5009E+07	1.8877E-02
16	3.500	3.000	3.6049E+07	1.1716E+08	2.1157E-02
17	3.000	2.500	6.8483E+07	1.8833E+08	1.4998E-02
18	2.500	2.000	8.6090E+07	1.9370E+08	1.2366E-02
19	2.000	1.660	8.0379E+07	1.4709E+08	1.0404E-02
20	1.660	1.500	7.6818E+07	1.2137E+08	2.7601E-02
21	1.500	1.340	8.0274E+07	1.1399E+08	1.9271E-02
22	1.340	1.330	7.7190E+06	1.0305E+07	1.0502E-01
23	1.330	1.000	2.1836E+08	2.5439E+08	1.4559E-02
24	1.000	0.800	1.8044E+08	1.6239E+08	1.4900E-02







# Conclusion

- Monte Carlo simulation software PHITS
- $\rightarrow$  Irradiation scenarios + Energy spectrum
- $\rightarrow$  Tally: [T-DCHAIN] + DCHAIN <3



# Perspectives

- 1. Comparison with experimental measurements
- 2. Study of the isotopic purity
- 3. Influence of other parameters : Target thickness, incident particles, other production channels (Ac-226/Th-226)
- 4. Spectrometry





# Backup



