

Simple rheological tests to assess the swallowing capacity of foods.

Rémi Lecanu, Cassandre Leverrier, Guy Della Valle, Marco Ramaioli

▶ To cite this version:

Rémi Lecanu, Cassandre Leverrier, Guy Della Valle, Marco Ramaioli. Simple rheological tests to assess the swallowing capacity of foods.. Annual European Rheology Conference 2022,, Apr 2022, Sevilla (Spain), Spain. hal-04522986

HAL Id: hal-04522986

https://hal.science/hal-04522986

Submitted on 27 Mar 2024

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



INRAO

Simple Rheological tests to assess the swallowing capacity of foods







Everything for the Rheologist



10 mL -

8 mL -

Level 3

Level 2

4 mL -

Level 1

1 mL -

Level 0

Rémi Lecanu¹, Cassandre Leverrier¹, Guy Della Valle², Marco Ramaioli¹

¹ Université Paris-Saclay, INRAE, AgroParisTech, UMR 782 SayFood, 91300 Massy & 78850 Thiverval-Grignon, France. ² INRAE UR 1268, Biopolymers, Interactions & Assemblies (BIA), 44316 Nantes, France

INTRODUCTION

Empirical rheological tests are commonly performed in the industry, in order to predict the processability of materials. For practical reasons, they have to be easily handled, to require short measurement time, and, ideally, their relevance should be ascertained by a comparison with scientific rheological measurements.

Regarding food, these constraints are not only requested for processability, but also to predict sensory properties and the capacity of the food to be swallowed safely. For instance, this last property is essential for individuals suffering from dysphagia.

Objectives of this work:

- To propose a classification of empirical rheological tests
- To present two different methods used for the characterization of food products, here thickened drinks.
- To delineate their main limitations and the necessity of carefully choosing the rheological test according to the needs and the type of product to assess.

COLLAPSING TESTS: The Line-Spread Test

Classification of simple rheological tests

We propose a non-exhaustive classification of methods that are not only related to foods. Methods developed in other fields are of great interest particularly in the oil, paint and concrete industries where the **control of texture** or **consistence** is mandatory for a steady processing. The main purpose of covering a wide variety of tests is to be able to select the most relevant one according to the needs of applications.


1.Imposed Stress			2.Imposed Strain Rate
Collapsing Tests	Gravity Flow Tests	Induced Flow Tests	
Slump test	Capillary Viscometer (U-shaped, Ostwald)	Capillary Piston Viscometer (MFI,)	Truck Mix
			Rotational Viscometer (Brookfield, Couette)
Line-Spread Test	Funnels (Marsh, Posthumus, Flow Cups, IDDSI Flow Test)	Falling Ball Viscometer	Vane Test
Bostwick Consistometer			Vibrational Viscometer

TESTS OVERVIEW


FLOW TESTS: Funnels and the Syringe Flow Test

The syringe flow test (SFT) is used to characterize the texture of drinks in the frame established by the International Dysphagia Standardization Diet Initiative (IDDSI)¹.

By allowing 10mL of fluid to flow during 10 seconds under hydrostatic pressure, the users can report the volume left in the syringe and link it to a texture level, from thin drinks (Level 0) to moderately thick (Level 3).



It is recommended to respect a precise length of syringe¹ (BD Syringe). But, it has been shown² that the flow in the syringe is highly sensitive to change in the nozzle diameter because the flow rate depend on the nozzle diameter power 3.



However these tests

Are dependent on geometry

Do not consider the effect of saliva

Lack objectivity

Both the IDDSI Flow test and the Line-Spread test are used to characterize the same type of food product, but, depending on the thickeners used, one might be more relevant than the others, for instance the IDDSI tests proved its use for xanthan gum-based thickeners^{5,6} while the line-spread test has shown to be more relevant for **starch**-based thickeners⁵.

NDDTF classification, but cannot provide an accurate assessment of

viscosity for fluids of similar texture³. Although close to the slump

test where the final height is also measured, and linked to a yield


stress τ_0^4 , the LST cannot be predictive of rheological parameters.

비스코업

Kim et al. (2018)

The line-spread test (LST) consists of a cylindrical cookie cutter of 3.5cm height and 5cm diameter which is filled with the thickened beverage, lifted, and allowed to flow on a plane surface marked with concentric circles for 60s.

This method characterize texture of thickened drinks in the classification of the National Dysphagia Diet Task Force (NDDTF), thin, nectar-like, honey-like.



Most other funnel test link the total efflux **time** t_f to a viscosity or apparent viscosity η . For Newtonian fluid a simple calibration with water allows to obtain α , a coefficient accounting for the **geometry**.

CONCLUSION

These tests are highly valuable because

- Not time consuming
- Do not need extensive trainings
- Are affordable

REFERENCES

- 1. IDDSI framework, https://iddsi.org/framework/
- 2. Matsuyama et al., International Journal of Food Science and Technology (2020)
- 3. Nicosia & Robbins. Dysphagia (2007)
- 4. Roussel & Coussot. Journal of Rheology (2005)
- 5. Kim et al., Journal of Texture Studies (2018)

Those tests have greatly improved the care of people with dysphagia disorder and thus their quality of life, but it is important to carefully choose them.

6. Jeong et al. MDPI Foods (2021)