



LC-HRMS-based metabolomics as a tool to develop analytical methods: How to choose the best extraction protocol when it comes to untargeted analysis?

Hikmat Ghosson, Delphine Raviglione, Cédric Bertrand, Marie-Virginie Salvia

► To cite this version:

Hikmat Ghosson, Delphine Raviglione, Cédric Bertrand, Marie-Virginie Salvia. LC-HRMS-based metabolomics as a tool to develop analytical methods: How to choose the best extraction protocol when it comes to untargeted analysis?. 16èmes Journées Scientifiques du Réseau Francophone de Métabolomique et Fluxomique, Jun 2024, Saint-Malo, France. , pp.S5-P22. hal-04606057

HAL Id: hal-04606057

<https://hal.science/hal-04606057>

Submitted on 9 Jun 2024

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Distributed under a Creative Commons Attribution - NonCommercial - NoDerivatives 4.0 International License

LC-HRMS-based metabolomics as a tool to develop analytical methods: How to choose the best extraction protocol when it comes to untargeted analysis?

Hikmat Ghosson^{1,2,†,*}, Delphine Raviglione^{1,2}, Cédric Bertrand^{1,2,3}, Marie-Virginie Salvia^{1,2}

*: hikmat.ghosson@univ-perp.fr; hikmatghosson@gmail.com

1: Centre de Recherches Insulaires et Observatoire de l'Environnement (CRIODE), Perpignan, France

2: UFR Sciences Exactes et Expérimentales (SEE), Université de Perpignan, Perpignan, France

3: S.A.S. AkiNaO, Perpignan, France

†: Actual affiliation: Chromatographies et Spectrométrie de Masse (CxM), Expertise Physico-Chimique (EPC), Michelin, Clermont-Ferrand, France

Context

In the current poster, different computational and graphical tools applied to LC-HRMS-based untargeted analysis are presented. These tools are used to compare 5 different extraction protocols developed to perform exhaustive untargeted analysis of formulated pesticides residues and microbial metabolites in soil. The comparisons are based on 4 analytical and contextual criteria: (1) the coverage of different types/families of compounds, (2) the quantitative performance, (3) the repeatability, and (4) the capacity to discriminate between contaminated and non-contaminated soils. Each of these tools needs specific usage precautions and presents advantages, limitations and complementarities with the other suggested tools.

The toolbox


Performance



Heatmap + HCA Dendrogram

- ❖ **Adjustable parameters:**
- ✓ Distance measuring method and clustering algorithm
 - ✓ Metadata organization (based on experimental design)
 - ✓ Intensity scale and color contrasts
 - ✓ Configuration architecture (holistic, polarity- or origins-segmented, etc.)


Advantages	Limitations
✓ Exhaustive 3D cartography of data	✗ Needs removal of ionization redundancies
✓ Compounds annotation is not mandatory	✗ No information on chemical families



Van Krevelen Diagram

- ❖ **Mandatory precautions:**
- ✓ Molecular formula annotation (level 4) is needed
 - ✓ Needs removal of ionization redundancies (advantage to have LC coupling to HRMS)
 - ✓ Intensity scaling needs to be carefully optimized

Advantages	Limitations
✓ Chemical families coverage assessment	✗ Non-exhaustive cartography
✓ Revealing heteroatoms-specific families	✗ Results quality depends on HRMS power



References

- Ghosson, H. et al. 2024. *Anal. Chem.*, 96(7):2810-2821
 Ghosson, H. 2020. Thesis, Université de Perpignan
 Gu, Z. 2022. *iMeta.*, 1(3):e43
 Pang, Z. et al. 2024. *Nucleic Acids Res.*, gkae253

Repeatability



All ions-RSDs vs. Euclidean Distances


❖ **Mandatory precautions:**

➢ All ions-RSDs

✓ Adjusting noise to zero to avoid insignificant amplification of RSDs

➢ Euclidean Distances

✓ Standardizing ion intensities to avoid “size effects” (e.g. logarithmic transformations)



Discrimination Power


R²X: Percentage of the explained variation that is causing the discrimination

R²Y: Correlation between the discrimination and the tested factor (here: contamination)

R²Y - Q²: Significance of the discrimination (overfitting assessment)

Q²: Discrimination predictivity: The chance (in %) to reliably determine the sample type

O: Systemic variations that are orthogonal to the tested factor



For more information

analytical chemistry

pubs.acs.org/ac


Article

LC-HRMS-Driven Computational Toolbox to Assess Extraction Protocols Dedicated to Untargeted Analysis: How to Ease Analyzing Pesticide-Contaminated Soils?

Hikmat Ghosson,^{*} Delphine Raviglione, Cédric Bertrand, and Marie-Virginie Salvia

Cite This: [Anal. Chem.](https://doi.org/10.1021/acs.analchem.4c00000) 2024, 96, 2810–2821

Read Online



Acknowledgments

Interreg POCTFE PALVIP



UNION EUROPÉENNE



AKINAO



INSERM

bioMarie

Web Networks

