

### Consequences of a telomerase-related fitness defect and chromosome substitution technology in yeast synIX strains

Laura Mcculloch, Vijayan Sambasivam, Amanda Hughes, Narayana Annaluru, Sivaprakash Ramalingam, Viola Fanfani, Evgenii Lobzaev, Leslie Mitchell, Jitong Cai, Hua Jiang, et al.

### ▶ To cite this version:

Laura Mcculloch, Vijayan Sambasivam, Amanda Hughes, Narayana Annaluru, Sivaprakash Ramalingam, et al.. Consequences of a telomerase-related fitness defect and chromosome substitution technology in yeast synIX strains. Cell Genomics, 2023, 3 (11), pp.100419. 10.1016/j.xgen.2023.100419 . hal-04658941

### HAL Id: hal-04658941 https://hal.science/hal-04658941

Submitted on 23 Jul 2024

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

# Consequences of a telomerase-related fitness defect and chromosome substitution technology in yeast *synIX* strains

Laura H. McCulloch, <sup>1,22</sup> Vijayan Sambasivam, <sup>2,22</sup> Amanda L. Hughes, <sup>3</sup> Narayana Annaluru, <sup>2,17</sup> Sivaprakash Ramalingam, <sup>2,21</sup> Viola Fanfani, <sup>4,19</sup> Evgenii Lobzaev, <sup>4,5</sup> Leslie A. Mitchell, <sup>1,20</sup> Jitong Cai, <sup>6</sup> The Build-A-Genome Class <sup>6,7,8</sup> Hua Jiang, <sup>9</sup> John LaCava, <sup>9,10</sup> Martin S. Taylor,<sup>11</sup> William R. Bishai,<sup>12</sup> Giovanni Stracquadanio,<sup>4</sup> Lars M. Steinmetz,<sup>3,14,15</sup> Joel S. Bader,<sup>6,13</sup> Weimin Zhang, <sup>1,\*</sup> Jef D. Boeke, <sup>1,13,16,23,\*</sup> and Srinivasan Chandrasegaran <sup>2,18,\*</sup>

The Build-A-Genome Class authors: Breeana G. Anderson, Abena Apaw, Pavlo Bohutskyi, Erin Buchanan, Daniel Chang, Melinda Chen, Eric Cooper, Amanda Deliere, Kallie Drakos, Justin Dubin, Christopher Fernandez, Zheyuan Guo, Thomas Harrelson, Dongwon Lee, Jessica McDade, Scott Melamed, Héloise Muller, Adithya Murali, José U. Niño Rivera, Mira Patel, Mary Rodley, Jenna Schwarz, Nirav Shelat, Josh S. Sims, Barrett Steinberg, James Steinhardt, Rishi K. Trivedi, Christopher Von Dollen, Tianyi Wang, Remus Wong, Yijie Xu, Noah Young, Karen Zeller, Allen Zhan.

1 Institute for Systems Genetics and Department of Biochemistry and Molecular Pharmacology, NYU Langone Health, New York, NY 10016, USA

2 Department of Environmental Health and Engineering, Bloomberg School of Public Health, Johns Hopkins University, Baltimore, MD 21205, USA

3 European Molecular Biology Laboratory (EMBL), Genome Biology Unit, 69117 Heidelberg, Germany

4 School of Biological Sciences, The University of Edinburgh, Edinburgh EH9 3BF, UK

5 School of Informatics, The University of Edinburgh, Edinburgh EH8 9AB, UK

6 Department of Biomedical Engineering, Whiting School of Engineering, Johns Hopkins University, Baltimore, MD 21218, USA

7 Department of Biology, Krieger School of Arts and Sciences, Johns Hopkins University, Baltimore, MD 21218, USA 8 Department of Chemical and Biomolecular Engineering, Whiting School of Engineering, Johns Hopkins University, Baltimore, MD 21218, USA

9 Laboratory of Cellular and Structural Biology, Rockefeller University, New York, NY 10065, USA

10 European Research Institute for the Biology of Ageing, University Medical Center Groningen, Groningen, the Netherlands

11 Department of Pathology, Massachusetts General Hospital, Harvard Medical School, Boston, MA 02114, USA 12 Department of Medicine/Division of Infectious Diseases, Johns Hopkins University School of Medicine, Baltimore, MD

21205, USA

13 High Throughput Biology Center, Johns Hopkins University School of Medicine, Baltimore, MD 21205, USA

14 Stanford Genome Technology Center, Stanford University, Palo Alto, CA 94304, USA

15 Department of Genetics, School of Medicine, Stanford University, Stanford, CA 94305, USA

16 Department of Biomedical Engineering, NYU Tandon School of Engineering, Brooklyn, NY 11201, USA

17 Present address: L'Oréal Research and Innovation, Clark, NJ 07066, USA

18 Present address: Pondicherry Biotech Private Limited, Pondicherry Engineering College Campus, East Coast Road, Pillaichavady, Puducherry 605014, India

19 Present address: Department of Biostatistics, Harvard T.H. Chan School of Public Health, Boston, MA 02115, USA 20 Present address: Neochromosome, Inc., Long Island City, NY 11101, USA

21 Present address: Center for Scientific and Industrial Research, Institute of Genomics & Integrative Biology, Sukhdev Vihar, Mathura Road, New Delhi 110025, India

22 These authors contributed equally

23 Lead contact

\*Correspondence: weimin.zhang@nyulangone.org (W.Z.), jef.boeke@nyulangone.org (J.D.B.), chandra@jhmi.edu (S.C.)

### 1 Summary

- 2
- We describe the complete synthesis, assembly, debugging, and characterization of a synthetic 3 4 404,963 bp yeast chromosome, synIX. Combined chromosome construction methods were used to 5 synthesize and integrate the left arm of synIX (synIXL) into a strain containing previously described synIXR. During synIX strain characterization, we identified and resolved a bug related to EST3, a 6 gene involved in telomerase function, producing a *synIX* strain with near wild-type fitness. To 7 8 facilitate future synthetic chromosome consolidation and increase flexibility of chromosome transfer 9 between distinct strains, we combined chromoduction, a method to transfer a whole chromosome 10 between two strains, with conditional centromere destabilization to substitute a chromosome of 11 interest for its native counterpart. We found that both synIX transfer via chromoduction and wild-12 type IX destabilization were efficient methods. We observed that wild-type II tended to co-transfer with synIX and was co-destabilized with wild-type IX, suggesting a potential gene dosage 13 compensation relationship between these two chromosomes. 14 15
- 16 **Keywords:** *Saccharomyces cerevisiae*, *synIX*, megachunk assembly, chromosome debugging,
- 17 EST3, chromoduction, chromosome substitution, centromere destabilization, transcriptomics
- 18

### 19 Introduction

#### 20

The Synthetic Yeast Genome Project (Sc2.0) marks a key milestone in the development of 21 22 "designer" eukaryotic genomes. This global effort seeks to produce a modified version of the ~12 23 Mb Saccharomyces cerevisiae genome from the bottom up that retains native strain fitness while 24 eliminating repetitive elements to improve genome stability and adding custom features to endow 25 new genome functionalities. These changes include deletion of tRNA genes and introduction of 26 loxPsym sites (nondirectional loxP sites that, being symmetric, can recombine with each other in 27 two distinct orientations) downstream of non-essential genes, enabling genome rearrangement via 28 the SCRaMbLE (Synthetic Chromosome Rearrangement and Modification by LoxP-mediated 29 Evolution) system and subsequent study of how different DNA sequence changes affect organismal fitness and function.<sup>1-4</sup> More broadly, by establishing a set of common design principles 30 for all sixteen yeast chromosomes, Sc2.0 aims to provide a blueprint for future eukaryotic genome 31 engineering endeavors. As part of this project, each chromosome in the yeast genome has been 32 33 redesigned and built from synthesized blocks of DNA. Following the previous construction of 34 several synthetic chromosomes, we report the assembly and characterization of synthetic 35 chromosome IX (synIX), coinciding with completion of the assembly phase of the Sc2.0 project as a whole.5-17 36

37

38 Our group initially assembled a circular version of the right arm of the chromosome, synIXR, as a 39 proof-of-principle for feasibility of the Sc2.0 project,<sup>6</sup> and subsequently incorporated that 40 chromosome arm into an otherwise wild-type and linear version of chromosome IX.<sup>2</sup> Here, we describe the construction and characterization of the left arm of chromosome IX, synIXL, and 41 consequently, the completion of synIX in its entirety. As with many of the other completed Sc2.0 42 chromosomes and in accordance with the Sc2.0 project objectives, synIX-containing strains, 43 44 following a "debugging" process to resolve any deficiencies in fitness across multiple growth 45 conditions, behave similarly to the ancestral BY4741 yeast strain with respect to growth fitness and 46 transcriptional profiling.

47

48 Over the past decade, fitness defects have emerged during completion of many of the Sc2.0 49 chromosome strains. We attribute these growth issues, or "bugs," to our incomplete understanding 50 of Saccharomyces cerevisiae and its genome during the project's initial design phase. Identifying 51 and resolving such problems may improve future genome engineering efforts and help us 52 anticipate and avoid potential complications. More specifically, bugs in previous Sc2.0 strains have 53 commonly arisen from design choices (sequence additions, deletions, or alterations) that 54 unexpectedly affect protein or RNA levels. For example, PCRTag recoding (intended to alter DNA 55 sequence without changing the resulting amino acid sequence, producing a watermark to 56 distinguish synthetic and wild-type sequences) has lowered protein expression or interfered with transcription factor binding sites in other synthetic chromosome strains.<sup>7,9,11,12</sup> Addition of loxPsym 57 58 sites and removal of introns as part of synthetic chromosome design have also affected RNA expression of genes in Sc2.0 strains, altering strain fitness.<sup>7,11</sup> In the case of synIX, we identify and 59 60 characterize a fitness defect related to a sequence adjacent to the EST3 gene, one of the 61 components of yeast telomerase holoenzyme that is involved in yeast telomere replication.<sup>18-20</sup> We 62 show that this was the consequence of a specific intentional design step, namely removal of a

- 63 tRNA gene and associated repetitive DNA sequences. Deletion of a tRNA gene, Ty1 retroelement,
- and the DNA between them led to reduced average telomere length and impaired the *synIX*
- 65 strain's ability to grow at higher temperature (37°C), a likely consequence of significant Est3p
- 66 reduction. Reintroduction of part or all of this sequence dramatically improved strain fitness, with
- associated partial or complete restoration of telomere length and Est3p levels depending on
- 68 sequence composition and variable effects on the strain's transcriptional profile.
- 69

70 The ultimate goal of Sc2.0 is to create a high-fitness yeast strain with a fully synthetic genome. A 71 previously described chromosome consolidation method that relies on chromosome 72 endoreduplication followed by sporulation becomes laborious and time-consuming as the synthetic chromosome number increases in the strain.<sup>2,21</sup> Here we present a strategy which combines 73 74 chromoduction of the chromosome of interest and subsequent loss of the native counterpart to 75 accelerate consolidation. This approach takes advantage of a karyogamy defect that arises when either parent in a yeast mating has a kar1-1 mutation.<sup>22</sup> Such crosses yield heterokaryons 76 77 consisting of nuclear material from one yeast strain plus a mixed cytoplasm of both parental 78 strains. In rare cases, whole chromosomes can transfer between strains following a kar1-1 x KAR1 79 cross.<sup>23</sup> While chromoduction occurs more frequently with smaller chromosomes than with larger 80 ones,<sup>23</sup> prior work has often employed a selection strategy to identify chromoductants in which a 81 particular chromosome of interest has been transferred.<sup>24,25</sup> Conversely, whole chromosome 82 destabilization can be achieved by activating a GAL promoter immediately upstream of yeast centromeres.<sup>7,26</sup> Here, we demonstrated the feasibility of synthetic chromosome transfer using 83 svnIX, and produced a final haploid veast strain in which svnIX has replaced wild-type 84 85 chromosome IX. This represents a promising proof-of-principle for targeted transfer of 86 chromosomes among yeast strains; indeed, this method is in active use to perform the final 87 consolidation of Sc2.0 chromosomes into a single veast strain.<sup>21</sup> 88

### 89 **Results**

90

### 91 **Design and assembly of** *synIX*

92 The design of *synIX* adheres to the same principles used for the Sc2.0 project.<sup>2</sup> The designed 93 sequence for synIX contains a variety of modifications, including deletion of ten tRNA genes 94 (relocated to a separate tRNA neochromosome) and 11,632 bp of repetitive DNA, recoding of 54 95 stop codons from TAG to TAA, 436 bp of restriction enzyme sites, and 7,943 bp of PCRTags 96 (recoded DNA sequences designed to distinguish between synthetic and native DNA without 97 altering the resulting amino acid sequence), and addition of 142 loxPsym sites downstream of 98 nonessential genes that allow for the Cre-mediated SCRaMbLE evolution system (Figure 1A).<sup>1-4</sup> 99 These changes reduced the total size of chromosome IX from a wild-type length of 439,885 bp to a 100 synthetic length of 404,963 bp.

101

102 Following the synthesis of *synIXR*,<sup>6</sup> we began construction of *synIXL* using a strain containing a

- 103 linear, partially synthetic chromosome *synIXR* in which the left arm was in a native state
- 104 (yLM461).<sup>2</sup> Assembly of *synIXL* proceeded in several stages (**Figure 1B**). Initially, we assembled
- 105 750 bp building blocks from synthetic 60-70 bp oligonucleotides through polymerase cycling
- assembly (PCA) in the context of the Build A Genome undergraduate class (**Figure S1A**).<sup>5,6</sup> We

107 then combined these building blocks to form 2-4 kb minichunks. In later stages of the project, as

108 DNA synthesis technology improved and costs dropped, we obtained the remaining minichunks 109 directly from a DNA synthesis vendor.

110

We divided *synIXL* into nine 30-60 kb megachunks, each made up of a series of minichunks that 111 overlapped adjacent ones. To assemble the first six megachunks (megachunks A-F), we 112 transformed component minichunks into our entry strain one megachunk at a time, including an 113 114 auxotrophic marker added to the final minichunk of each region (Figure 1C). Through yeast 115 homologous recombination, each assembly step overwrote the wild-type chromosome segment 116 and auxotrophic marker attached to it with the synthetic segment and a new auxotrophic marker in 117 accordance with the Switching Auxotrophies Progressively for Integration (SwAP-In) method,<sup>2,6</sup> 118 enabling selection of strains containing newly integrated DNA (Figure 1C, Table S1). After each 119 round of integration, the PCRTag watermarking system allowed identification of colonies 120 containing desired synthetic DNA and lacking corresponding wild-type chrIX DNA (Supplementary 121 File 1).

122

123 The minichunk integration approach used in the assembly of megachunks A-F of synIXL left some 124 unwanted patches of wild-type *chrIX* sequence (Figure 1B, top). To improve the efficiency of synthetic DNA incorporation and to avoid such wild-type DNA patches, later portions of synIXL 125 used the "megachunk cloning" integration strategy.<sup>12</sup> For these sections of the chromosome, we 126 127 first assembled minichunks into 40-60 kb megachunks (megachunks G-I) as extrachromosomal plasmids in yeast, with an assembly success rate of approximately 20-45% (Figure 1D, Table S2), 128 129 and subsequently transformed them into *E. coli* for plasmid extraction. Following sequence verification of the megachunks, these larger assembled pieces were released from their plasmid 130 131 backbones using restriction enzyme digestion, and delivered directly into our semisynthetic (A-F) 132 synIX strain for SwAP-In. We verified successful integration of the three megachunks via PCRTagging analysis and whole genome sequencing. All three megachunks constructed in this 133 134 manner were integrated into the synIX strain in their entirety and lacked novel mutations (Figure 135 1B, middle).

136

To convert the wild-type DNA patches seen in megachunks A-F of *synIXL* to synthetic, we used a CRISPR-Cas9 editing strategy to selectively target *synIX* at residual *chrIX* PCRTags within each patch (**Figure S1B**) and repair these sequences with appropriate synthetic donor DNAs.<sup>9,10</sup>
Following replacement of ten minichunks over seven rounds of CRISPR-mediated editing, we obtained a version of *synIX* that contained all the synthetic minichunks included in the

- 142 chromosome's design (**Figure 1B, bottom, Table S3**).
- 143

We continued CRISPR-mediated editing of this draft strain to correct point mutations, small

deletions, and duplications in *synIX*, as detailed in **Table S4**. Additionally, we discovered a
 discrepancy that stemmed from inaccuracies in the wild-type yeast reference genome sequence

146 discrepancy that stemmed from inaccuracies in the wild-type yeast reference genome sequence 147 that had originally been used for synthetic chromosome design. In *synIX*, the designer sequence

147 that had originally been used for synthetic chromosome design. In *synIX*, the designer sequence 148 for the non-essential *FAA3* gene contained an extra seven base pairs missing from the most

recent version of the wild-type reference sequence, resulting in a frameshift (**Figure 1E**).<sup>27,28</sup> A

149 recent version of the whit-type reference sequence, resulting in a frameshift (**Figure 1E**).<sup>-1</sup>

150 previous genetic screen linked the  $faa3\Delta$  null allele to decreased vegetative growth;<sup>29</sup>

consequently, we restored proper *FAA3* function in our *synIX* strain by deleting the unwanted
 sequence via CRISPR-mediated editing, reestablishing the normal reading frame of *FAA3* (Figure
 S2).

154

155 During this editing process, we unexpectedly identified whole-chromosome disomy of synlX (Figure 1F, top). To restore the normal karyotype of the synIX strain, we integrated a URA3-pGAL 156 cassette upstream of CEN9 into one synIX copy,<sup>26</sup> induced chromosome destabilization and loss 157 158 via growth in media containing galactose, and then selected against the URA3 marker to identify 159 strains with a single copy of synIX (Figure S3). Whole genome sequencing confirmed that the resulting synIX strain contained approximately equal read depth for each chromosome, as 160 expected for a monosomic haploid yeast strain (Figure 1F, bottom). There were no overt 161 162 phenotypic differences between the disomic and monosomic strains, suggesting that the disomy was an inadvertent consequence of CRISPR-mediated editing rather than selection for a second 163 copy of *synIX* due to haploinsufficiency. This monosomic strain was used for further debugging. 164 165



<sup>166</sup> 167

Figure 1. Design and assembly of *synIX*. (A) Diagram of the differences between wild-type (WT) *chrIX* and *synIX*.
 Designer features specific to *synIX* include the addition of PCRTags and loxPsym sites, the recoding of TAG stop

170 codons to TAA, the replacement of wild-type telomeres with universal telomere caps (UTCs), and the removal of tRNA

171 genes (relocated to a tRNA neochromosome), introns, retrotransposons, and subtelomeric repeats. (B) Schematic of 172 the synIXL construction process, reflecting the synthetic and (unwanted) WT composition of chromosome IX in the 173 synIX strain after each stage of assembly. Stage 1: SwAP-In with minichunks (A-F). Stage 2: SwAP-In with pre-174 assembled megachunks (G-I, opague with blue lines). Stage 3: Error correction with CRISPR-Cas9 (primarily in 175 highlighted megachunks A-F). Blue: syn minichunks, red: WT minichunks, gray: megachunks. Red patches in stages 1 176 and 2 reflect segments of DNA that were not replaced by their expected synthetic MC counterparts during SwAP-In 177 with minichunks in stage 1 (A-F); appropriate synthetic sequences were integrated at these sites during stage 3. (C) 178 Minichunk integration strategy for assembling megachunks A-F of svn/X in yeast (used for stage 1 in (B)). A kanMX 179 marker (black) was integrated into the left end of wild-type chrIX (red). Individual minichunks (blue, synthetic) 180 comprising one megachunk were then co-transformed into the in-progress synIX strain and used to overwrite the wild-181 type chrIX sequence (red) via homologous recombination. Alternating auxotrophic markers (LEU2, brown, and URA3. 182 orange) were used and overwritten at each step, allowing for selection of the new marker and counterselection against 183 the previous round's marker. (D) Megachunk plasmid assembly and integration approach for building and integrating 184 megachunks G-I in synIX (used for stage 2 in (B)). Minichunks (blue) were assembled in yeast, and an auxotrophic 185 marker (here, URA3, orange) was added to the end of each assembly, followed by a segment of chrIX homology. The 186 megachunk assemblon overwrote the prior round auxotrophic marker (brown, LEU2) and wild-type DNA (red) by 187 homologous recombination. (E) Schematic depicting a frameshift mutation in the FAA3 gene in the synIX design. The 188 synthetic sequence included seven base pairs not present in the S288C wild-type yeast reference genome. The added 189 bases, shown in red, cause a frameshift, with the resulting synIX amino acid sequence (red) varying from the expected 190 chrlX amino acid sequence (black, bottom). (F) Coverage plots showing read depth along each yeast chromosome in 191 disomic synIX strain yLHM0588 (top) and monosomic synIX strain yLHM0721 (bottom). X axis: Yeast chromosome 192 (not to scale). Y axis: Relative depth based on number of reads at each position divided by average read depth across 193 the sixteen yeast chromosomes.

194

### 195 **Debugging of** *synlX*

196 After correcting the above *synIX* sequence issues, we found that this strain, coined

197 yeast chr09 9 1 (also referred to as yLHM1192), grew more poorly than the parental BY4741 strain under high-temperature (37°C) conditions on both dextrose-containing (YPD) and glycerol-198 containing (YPG) media (Figure 2A). A similar phenotype was observed in additional strain 199 isolates equivalent to yLHM1192 (Figure S5). Following a pooled sequencing bug mapping 200 201 method,<sup>9</sup> we backcrossed our synIX strain to a wild-type strain to generate a heterozygous diploid veast strain for sporulation, and subsequently dissected tetrads. The goal of this process was to 202 203 obtain spores containing a randomized mixture of recombinant synIX and chrIX DNA as a 204 consequence of meiotic recombination (Figure S4A). By conducting whole genome sequencing of 205 healthy and sick spores, we identified a region spanning from approximately 15-20 kb upstream of 206 CEN9 to 5 kb downstream of CEN9 that appeared to be associated with our synIX fitness defect 207 (Figures 2B, S4B). Consistent with this region's involvement in synIX's growth deficiency, our synIX strain initially exhibited a marked decline in fitness on YPD at 37°C following the integration 208 209 of megachunk I, which contained the synthetic sequence corresponding to this portion of the chromosome (Figure S6). To determine the gene responsible for the bug within this ~20-25 kb 210 211 region, we assembled the wild-type version of each of the genes in this region into a CEN plasmid 212 for complementation (Figure S7). The strains containing a plasmid with EST3 showed a 213 pronounced improvement in fitness on both YPD and YPG at 37°C compared with the buggy synIX 214 strain (Figures 2B, S7).

215

To investigate which modifications of the synthetic *EST3* gene might contribute to this fitness

217 defect, we replaced sections of the synthetic region with corresponding wild-type sequences in our

218 synIX strain and monitored strain fitness. Replacing the coding sequence of EST3 was insufficient

219 to restore fitness on YPD or YPG at 37°C (Figure 2C, yLHM1429); however, restoring both a 220 deleted segment containing an aspartyl tRNA gene and a nearby Ty1 long terminal repeat (LTR) 221 element dramatically improved *svnIX* strain fitness on both YPD and YPG at 37°C (Figure 2C. **vLHM1430** and **vLHM1431**). To better assess the contributions of these *EST3*-adjacent elements 222 to strain fitness, we next introduced smaller modifications to EST3's upstream region in the synIX 223 224 strain, including different combinations of the aspartyl tRNA gene, nearby Tv1 LTR, and an 225 intervening single-copy "gap" DNA sequence with no known function that separates the tRNA gene 226 and Ty1 LTR in the native *chrIX* sequence (Figure 3A, Table S5). We found that all of these edits 227 improved strain growth fitness (Figure 3A). Similar phenotypes were observed in additional strains 228 generated for each modification, and did not noticeably change following extended strain 229 passaging (Figure S8).

230

As EST3 is involved in telomere maintenance, we hypothesized that EST3-associated changes in 231 232 fitness might correlate with changes in average strain telomere length. Consequently, we next assessed telomere length of our synIX and variant strains via Southern blotting with a probe 233 specific for telomeric repeats (Figure 3B). Following enzymatic digestion of yeast gDNA with Xhol, 234 a smear of telomeric fragments from chromosome ends containing Y' elements can be seen at the 235 bottom or each lane; these fragments were used to estimate average telomere length, with other 236 237 telomeric fragments appearing at larger sizes in each lane. As expected, vLHM1192's telomeres were markedly shorter than those seen in the wild-type BY4741 strain. Interestingly, however, the 238 239 telomeres in the modified synIX strains split into two length groups; three groups of strains, yLHM1504 (containing the entire tRNA, intervening "gap" sequence between the tRNA and Ty1 240 241 with no known function, and Ty1 sequence), yLHM1506 (containing the tRNA and gap sequences, but not the Ty1 sequence), and yLHM1601 (containing a tRNA with a transcription-inactivating 242 point mutation<sup>30-33</sup> and the aforementioned gap sequence) had telomeres similar in length to or 243 slightly longer than wild-type yeast, while two other groups of strains, yLHM1505 (containing just 244 245 the gap sequence and Tv1 LTR) and vLHM1591 (containing just the gap sequence), had 246 telomeres that were shorter than BY4741, but slightly longer than yLHM1192. 247

248 Next, we wanted to assess the role of these upstream EST3 modifications on strain fitness and 249 telomere length independent of any other synIX-associated strain changes. To do this, we 250 introduced each of our synthetic modifications upstream of EST3 (from yLHM1192, yLHM1504-251 1506, yLHM1591, and yLHM1601) to wild-type yeast, while otherwise maintaining the native chrIX 252 sequence (Figure 4A). As in our synIX strains, the yLHM1192-equivalent wild-type edited strains showed a decline in fitness on YPG plates at 37°C, while the other modified strains showed fitness 253 254 similar to wild-type BY4741 (Figure S9). Fitness did not change appreciably following extended strain passaging (**Figure S9**). We again examined average telomere length of these strains via 255 Southern blot (**Figure 4B**). As with the synthetic strains, we saw three groups of telomere lengths: 256 257 wild-type (BY4741, as well as the strains with the vLHM1504, vLHM1506, and vLHM1601-258 equivalent modifications), shortened (strains with the yLHM1192-equivalent modifications), and 259 intermediate (strains with the yLHM1505-equivalent and yLHM1591-equivalent modifications). 260

We next sought to understand why all of our modifications improved strain fitness, but only some fully restored telomere length to wild-type levels. We hypothesized that differences in *EST3* 

translation and/or transcription might result in variable telomere length across strains. To examine 263 Est3p expression in our modified synIX strains, we tagged the 3' end of EST3 with a 3X FLAG tag 264 and monitored protein levels after strain growth at 37°C (Figure 5A, Table S6). FLAG-tagged 265 synthetic strains were crossed to an untagged wild-type strain of opposite mating type (BY4742) to 266 complement any other recessive svn/X-related effects on strain growth while allowing for 267 monitoring of Est3p production from the synthetic locus. Consistent with the telomere length data, 268 269 the vLHM1504-, vLHM1506-, and vLHM1601-derived strains expressed higher levels of FLAG-270 tagged Est3p than did the vLHM1192-, vLHM1505-, and vLHM1591-derived strains (Figure 5A). A 271 follow-up comparison of the vLHM1192-, vLHM1505-, and vLHM1591-derived strains revealed 272 higher Est3p levels in the yLHM1505-derived and yLHM1591-derived strains than in the 273 yLHM1192-derived strains, although all three strains had markedly lower Est3p levels than did wild-type yeast (Figure 5B). We saw neither longer nor shorter isoforms of FLAG-tagged Est3p 274 associated with the low-abundance strains (Figures 5B, S10). Similar results were observed in 275 276 FLAG-tagged heterozygous diploid strains generated by crossing FLAG-tagged versions of the 277 wild-type EST3 variant strains from Figure 4A and Figure S9 to untagged BY4742 and monitoring 278 Est3p levels (Figures S11A. S11B).

279

280 To further elucidate what might contribute to the variations in Est3p levels in our engineered synIX 281 strains, we examined RNA sequencing data to look for changes in RNA expression patterns. Of 282 note, strand-specific RNA sequencing data revealed differences in the pattern of EST3 transcripts in the wild-type and synIX strains at both 30°C and 37°C (Figures 5C, 5D, S12A, S12B). In 283 particular, while both strains showed sequencing read depth throughout the EST3 coding region, 284 285 the synIX strain contained transcripts upstream of the expected EST3 transcription start site that were absent from the wild-type strain. Importantly, given that these transcripts were transcribed 286 from the same strand as native EST3, some EST3 transcripts might have had an extended 5' UTR. 287 predicted to encode upstream AUG codons that would not produce functional Est3p. In our 288 289 modified *svnIX* strains, we saw two classes of upstream transcript profiles (Figures 5E, 5F). 290 yLHM1504, yLHM1506, and yLHM1601, the three strains with wild-type-length or slightly longer 291 telomeres, showed a relative reduction in RNA sequencing reads directly upstream of EST3 292 compared with strand-specific total reads aligned to the upstream plus coding sequence of EST3 293 (Figures 5E, 5F, S13A. S13B). By contrast, yLHM1505 and yLHM1591, the two strains with 294 intermediate-length telomeres, still showed substantial transcript read depth upstream of EST3. The wild-type strains with modifications upstream of EST3 yielded a similar transcriptional pattern 295 296 (Figures S14A, S14B).

297

298 Nanopore direct RNA sequencing, which evaluates full-length transcript patterns, revealed a 299 distinct population of EST3-spanning transcripts in vLHM1192 containing at least one out-of-frame 300 AUG upstream of the expected EST3 transcription start site (Figure 5G). Many of these transcripts started 250-500 bp upstream of the AUG that normally gives rise to the first methionine of 301 302 functional Est3p (Figures S15). We hypothesize that such "non-translatable" upstream transcripts fail to yield properly translated Est3p, contributing to reduced telomerase function and decreased 303 strain fitness (Figure S16). By contrast, while vLHM1505 and vLHM1591 contained a few non-304 translatable transcripts spanning EST3, most of the transcripts upstream of EST3 that were 305 observed in these strains terminated prior to the expected EST3 transcription start site (gray in 306

9

Figure 5G). yLHM1192, yLHM1505, and yLHM1591 all had relatively fewer full-length, presumably translatable *EST3* transcripts than did the tRNA-containing strains (yLHM1504, yLHM1506, and yLHM1601).

310

Taken together, these results point to a model in which upstream initiated transcripts are unable to 311 produce Est3p, and only reads initiating at or near the native transcription start site produce 312 functional EST3 transcripts and resulting protein. The original synIX strain (yLHM1192) appears to 313 314 yield almost no functional Est3p, while variant strains produce either low (yLHM1505 and yLHM1591) or near-wild-type (yLHM1504, yLHM1506, and yLHM1601) Estp3 levels. However, 315 316 strains with even a low level of functional Est3p and a moderate reduction in telomere length 317 appear to surpass a threshold for maintaining normal strain fitness. Thus, our three classes of Est3p levels and telomere length yield only two fitness profiles; healthy strains and sick ones. 318 319



#### 320 321

322 Figure 2. Identifying and mapping an EST3-related bug. (A) Spot assays comparing the growth of BY4741 (wild-type) 323 to yLHM1192 (synthetic) across a variety of growth conditions. Each column represents a tenfold dilution. RT: Room 324 temperature, approximately 22°C. (B) (Top) Schematic of region near CEN9 identified as containing the gene 325 responsible for a synIX fitness defect. The unhealthy synIX strain (yLHM1192) was transformed with plasmids 326 containing transcription units for each of the genes shown. (Bottom) Spot assays comparing the growth of BY4741 and 327 yLHM1192 to yLHM1192 transformed with a plasmid containing the EST3 coding sequence plus 500 bp of upstream 328 sequence and 200 bp of downstream sequence. Each column represents a tenfold dilution. (C) Spot assays 329 comparing growth of strains with different combinations of synthetic and wild-type features in and upstream of EST3. 330 Each column represents a tenfold dilution. Schematics on left illustrate wild-type (red) and synthetic (blue) features 331 present in each strain. The asterisk (\*) marks the site of a programmed +1 ribosomal frameshift during EST3 332 translation.

333

	YPD, 30°C, 2d			YP	YPD, 37°C, 2d				PG, 3	7°C,	7d		
BY4741		۲	-	:		۲	皺				橡	-	PCRTag + PCRTag tRNA Gap seq Ty1 LTR
yLHM1192			*	22	•	۲	-	:					Buggy synIX strain
yLHM1504		۲	鐮	- ,2			*	ġ.	$\bigcirc$		575	1.	-K L EST3 L
yLHM1505			物	•		۲	*	-					Gap + Ty1 LTR restored
yLHM1506		۲	*	• •		۲	÷.		$\bigcirc$			~~	tRNA + gap restored
yLHM1591		۲	-	: %			*	**	0		1		Gap restored
yLHM1601		۲	Ser.	7:		۲		1	$\bigcirc$			2.0	Mutated tRNA + gap restored

В



#### 334 335

336 Figure 3. Modifying synIX strains upstream of EST3 alters fitness and telomere length. (A) Spot assays and 337 accompanying schematics of the region upstream of EST3 for strains BY4741 (wild-type chrIX), yLHM1192 (synthetic), 338 and variant strains derived from yLHM1192. Wild-type PCRTags, tRNA gene (Asp), and Ty1 long terminal repeat 339 (LTR) are depicted in red, synthetic PCRTags and loxPsym sites in blue, mutation in tRNA<sup>Asp</sup> gene as a vertical black 340 line. Each spot assay column represents a tenfold dilution. (B) Southern blot of Xhol-digested yeast gDNA derived 341 from BY4741 (wild-type) and replicate strains with modifications equivalent to yLHM1192 (synthetic) and variant strains 342 (yLHM1504, yLHM1505, yLHM1506, yLHM1591, and yLHM1601, as depicted in (A)). DNA was probed with a 343 digoxigenin-labeled fragment specific for telomeric repeats. Black arrow (right): Average size of Y'-containing telomeric 344 fragments for BY4741. Gray arrow (right): Average size of Y'-containing telomeric fragments for yLHM1192. Left and 345 middle panels: Left and right halves of one Southern blot, cropped to remove one lane in the middle. Right panel: Right 346 side of a second Southern blot, cropped to remove left-hand ladder, BY4741, and yLHM1192 lanes.



#### 347 348

349 Figure 4. Modifying wild-type yeast strains upstream of EST3 alters telomere length. (A) Schematics of the region 350 upstream of EST3 for strains BY4741 (wild-type chrIX), yLHM1192-equivalent (synthetic modification matching 351 yLHM1192 in an otherwise wild-type yeast strain), and variant strains derived from the yLHM1192-equivalent strain. 352 Wild-type PCRTags, tRNA gene (Asp), and Ty1 long terminal repeat (LTR) are depicted in red, loxPsym sites in blue, 353 mutation in tRNA<sup>Asp</sup> gene as a vertical black line. (B) Southern blot of Xhol-digested yeast gDNA derived from BY4741 354 (wild-type), yLHM1192 (synthetic), and replicate strains with synthetic-equivalent modifications made upstream of 355 EST3, as depicted in (A)). DNA was probed with a digoxigenin-labeled fragment specific for telomeric repeats. Black 356 arrow (right): Average size of Y'-containing telomeric fragments for BY4741. Gray arrow (right): Average size of Y'-357 containing telomeric fragments for BY4741 + yLHM1192-equivalent deletion strain.

358



359 360

361 Figure 5. Dissecting the effects of EST3-adjacent features on EST3 protein and RNA expression. (A) Immunoblotting 362 for immunoprecipitated (IP) Est3p (top) and H3 histone pre-IP loading control (bottom). Est3p was C-terminally tagged 363 with 3X FLAG in each synthetic strain, and diploids were generated via mating with wild-type BY4742. Two 364 independent colonies were grown in YPD at 37°C for each strain. Right-most lane (second ladder) was cropped out of 365 image. (B) Immunoblotting for immunoprecipitated (IP) Est3p (top) and H3 histone pre-IP loading control (bottom) for 366 low-expressing Est3p strains plus controls. Est3p was C-terminally tagged with 3X FLAG in each synthetic strain, and 367 diploids were generated via mating with wild-type BY4742. Three independent colonies were grown at YPD at 37°C for 368 each yLHM1192-, yLHM1505-, and yLHM1591-derived diploid strain. Right-most lane (second ladder) was cropped 369 out of image. (C) Strand-specific RNA sequencing alignments for the reverse strand of BY4741 and yLHM1192 370 (strand from which EST3 transcription is expected, i.e., "reverse strand reads") grown at 37°C in YPD to a composite 371 reference containing both synthetic and native IX features. Schematic on top illustrates wild-type chrIX (red) and 372 synthetic (blue) features. Black box on the RNA sequencing alignment plot highlights region between the expected 373 EST3 start codon and the upstream tRNA. The asterisk (\*) marks the site of a programmed +1 ribosomal frameshift 374 during EST3 translation. (D) Quantification of reverse strand reads mapping directly upstream of EST3 (between tRNA 375 and start codon, corresponding to black box in C) as a fraction of the total reverse strand EST3 coding and upstream 376 reads for each strain shown in C. Read counts based on alignment to native (BY4741) or synthetic (yLHM1192) 377 reference sequence. Error bars represent standard deviation of three biological replicates. (E) Spot assays and strand-378 specific RNA sequencing alignments for the reverse strand reads aligning to the EST3 region for samples grown at 379 37°C on YPD plates or in YPD liquid media. Strains: BY4741 (wild-type *chrlX*), yLHM1192 (synthetic), and variant 380 strains derived from yLHM1192. Wild-type PCRTags, tRNA gene (Asp), and Ty1 long terminal repeat (LTR) are

381 depicted in red, synthetic PCRTags and loxPsym sites in blue, mutation in tRNA<sup>Asp</sup> gene as a vertical black line. Black 382 box on the RNA sequencing alignment plot highlights region between the expected EST3 start codon and the 383 upstream tRNA. The asterisk (\*) marks the site of a programmed +1 ribosomal frameshift during EST3 translation. (F) 384 Quantification of reverse strand reads mapping directly upstream of EST3 (between tRNA and start codon, 385 corresponding to black box in E) as a fraction of the total reverse strand EST3 coding and upstream reads for each 386 strain shown in E. Read counts based on alignment to native (BY4741) or synthetic (yLHM1192) reference sequence. 387 Error bars represent standard deviation of three biological replicates. (G) Nanopore direct RNA sequencing reads 388 aligned to the EST3 region (±3kb). Strains: BY4741, yLHM1192, and variant strains derived from yLHM1192, all grown 389 at 37°C in YPD liquid media. X axis: Chromosome coordinate. Y axis: Number of reads. Forward strand reads are 390 above the Y axis, reverse strand reads below the Y axis. Dotted lines indicate the boundaries of the EST3 coding 391 region. EST3 transcripts are expected to lie on the reverse strand. Yellow reads: "Non-translatable" reads that contain 392 at least one AUG upstream of the expected EST3 start codon plus a minimum of 20 bp of 5' UTR sequence. Blue 393 reads: Reads mapping to EST3 that start no more than 20 bp upstream of the first AUG upstream of the expected 394 EST3 start codon, including "translatable/functional" reads spanning the entirety of EST3 and reads starting 395 downstream of the EST3-initiating AUG.

396

#### 397 Characterization of the healthy *synIX* strain

398 We ultimately proceeded with the strain in which a mutated version of the deleted tRNA gene and 399 the adjacent gap sequence were reintroduced upstream of EST3 (yLHM1601 in Figure 5E). The tRNA sequence contained a point mutation previously shown to block tRNA transcription,<sup>30-33</sup> 400 allowing functional tRNA relocation to a tRNA neochromosome<sup>15</sup> while maintaining Est3p levels 401 similar to those of wild-type yeast (Figure 5A). This strain is most consistent with Sc2.0's design 402 403 principles of maximizing strain fitness while eliminating repetitive elements and tRNA genes in particular. This final synIX version, coined yeast chr09 9 2, grew similarly to wild-type yeast 404 405 across a variety of media types and temperature conditions, as did additional sequence-matched replicate strains (Figures 6A, S17, S18, S19). The final strain largely matched the expected 406 407 synthetic reference sequence, with slight deviations from the original design reflecting changes 408 made during bug fixing and other minor alterations during the assembly process (Table S7). The 409 final strain also showed even coverage across *svnIX* and normal chromosomal copy number 410 (Figure 6B). SynIX migrated slightly faster on pulsed-field gel electrophoresis (PFGE) than did 411 native IX (Figure 6C). We also observed an unusual banding pattern among a couple of the higher 412 bands in yLHM1601, possibly related to an expanded repetitive sequence near the telomere of chromosome XVI. This difference did not appear to affect strain fitness, and disappeared following 413 synIX transfer to an alternative matched-background wild-type yeast strain via the chromosome 414 415 substitution method discussed in the next section (Figure 6C, lanes 3 and 4).

416

417 Transcriptional analysis of the completed synIX strain (yLHM1601, or yeast chr09 9 2) identified a small number of differentially expressed genes compared to the wild-type strain (Figures 6D, 6E, 418 419 S20A, S20B). Several hits found on synIX were located in or near telomeres, with some of these 420 genes showing upregulation and other showing downregulation. On the left arm of synIX, SOA1 (YIL166C) showed decreased expression in the synIX strain at 37°C (~5x lower than in BY4741). 421 Several of the other genes originally found between SOA1 and TEL09L in native chromosome IX 422 423 were deleted in the synIX design. Consequently, we hypothesize that SOA1 expression may be decreased due to increased silencing associated with its new position near TEL09L, as has been 424 425 observed on other synthetic chromosomes where genes moved closer to the telomeres due to subtelomeric DNA deletion.<sup>12</sup> Conversely, right arm subtelomeric genes YIR042C and YIR043C 426

427 showed expression increases at 37°C (~12x increase for YIR042C) and both 30°C and 37°C (~15x 428 and ~34x increases for YIR043C), respectively. Unlike SOA1, YIR042C and YIR043C were already located very close to TEL09R in native yeast chromosome IX. As a result, it is possible that 429 430 the universal telomere caps introduced to synIX per Sc2.0 design principles failed to silence nearby genes to the same degree as the native telomeres, since the earliest version of the UTC 431 used in this case lacked the X sequence used at the left telomere and elsewhere in the project. 432 433 Since the X sequence, a known subtelomeric anti-silencing region.<sup>34</sup> is present in the native TEL09R, its absence could explain the selectively elevated expression of these TEL09R-adjacent 434 435 genes.

436



437 438

439 Figure 6. SynIX characterization. (A) Spot assays comparing the growth of BY4741 (wild-type) to yLHM1601 440 (synthetic) across a variety of temperature conditions. Each column represents a tenfold dilution. RT: Room 441 temperature, approximately 22°C. (B) DNA sequencing coverage plot for synIX strain (yLHM1601). X axis: Yeast 442 chromosome (not to scale), or coordinate along synIX. Relative depth (Y axis) based on reads mapped to the 443 yeast chr09 3 55 reference sequence divided by the average read depth across the sixteen yeast chromosomes 444 (top) or svn/X (bottom), (C) Pulsed-field gel showing chromosomes from BY4741 (chrlX), vLHM1601 (svn/X). 445 yLHM2337 (BY4742-background strain with copies of chrIX and synIX), and yLHM2401 (yLHM2337 after loss of chrIX) 446 strains. Vertical black line indicates the location of a potential discrepancy in band intensity for larger chromosomes 447 between BY4741 and vLHM1601; synIX was moved to an alternative strain (vLHM2337) with larger bands resembling 448 BY4741 in intensity, and *chrIX* was subsequently lost from this strain (producing yLHM2401) via the strategy described 449 in the next section. (D-E) Volcano plots of differentially expressed genes obtained from RNA sequencing data for synIX 450 strain (yLHM1601) vs. BY4741 measured at 30°C (D) and 37°C (E) in YPD. Upregulated genes in yLHM1601 are 451 depicted in red (chromosome IX), downregulated genes in medium blue (chromosome IX) or light blue (other 452 chromosomes). Transcript counts are based on alignment to the S288C reference transcriptome. Auxotrophic gene 453 LYS2 was present in BY4741, but not yLHM1601. Fold change cutoff is 4, adjusted p-value cutoff is 0.01. Three

biological replicates were used for each strain. For corresponding unfiltered plots of RNA sequencing data depicting all
 measured genes, see Figures S21A and S21B.

456

### 457 Chromosome substitution of synIX

Here, we sought to combine chromoduction with subsequent loss of the corresponding native 458 459 chromosome from the recipient strain to seamlessly move svn/X between yeast strains while 460 maintaining normal karyotype in the final strain of interest (Figure 7A). As a proof-of-concept, a disomic synIX strain was used as the donor, as an extra copy of synIX was deemed likely to 461 produce higher chromosome transfer efficiency.<sup>24</sup> A kar1-1 mutant strain was used as the recipient 462 463 strain for blocking karyogamy during mating. Notably, the recipient strain also contained two mutations, can1 and cyh2, as recessive selectable markers to ensure robust exclusion of donor 464 cells or potential diploid cells on canavanine- and cycloheximide-containing medium. To facilitate 465 selection of recipient strains containing synIX transferred from the donor strain following strain 466 467 crossing, we deleted the LYS12 gene on the recipient strain's wild-type chromosome IX (Figure **S21A**). This allowed facile selection for lysine prototrophy to ensure *synIX* was present in the 468 chromoductants. Additionally, to enable eventual loss of the wild-type chrIX from the recipient 469 470 strain, we subsequently knocked in a URA3-pGAL cassette upstream of CEN9 in the recipient

- 471 strain for *chrIX* destabilization and counterselection (**Figure S21B**).<sup>26</sup>
- 472

After modifying the recipient strain, we mated the donor and recipient strains (yLHM0387 and 473 yWZ601, respectively) and selected for recipient cells that had acquired at least one copy of synIX 474 via chromoduction on SC–Lys–Arg+canavanine+cycloheximide plates. A PCRTagging assay was 475 476 performed using selected PCRTags across *chrIX* and *synIX*. Four chromoductants (yWZ610-613) containing both chrIX and synIX PCRTag amplicons were isolated from 14 candidates (Figure 477 478 **S21C**). Disomic chromoductants stably maintained both *chrIX* and *svnIX* after restreaking onto 479 fresh lysine-deficient plates. Genotyping validation of the four disomic chromoductants suggested 480 that yWZ612 and yWZ613 had additionally acquired the HIS4 gene on chromosome III from the donor strain (the recipient strain was  $his4\Delta 0$ ), suggesting chromosome III was co-transferred along 481 with synIX (Figure S21D). 482

483

484 We next sought to destabilize chrIX in the disomic chromoductants to generate haploid strains with a single copy of synIX. We grew the disomic chromoductants (yWZ610, yWZ611) in galactose-485 containing media to activate transcription of the GAL promoter, destabilizing CEN9 and 486 consequently leading to chrIX loss during mitosis.<sup>26</sup> PCRTagging analysis of the resulting 487 488 chromoductants (yWZ618, yWZ619) showed that only synIX was present (Figure S21E). Lastly, we subjected a series of yeast strains we generated during chromoduction to a PFGE assay. As 489 490 expected, 1) synIX was slightly smaller than chrIX in both BY4741 and the recipient strain, 2) both 491 synIX and chrIX were present in disomic chromoductants, 3) chromosome III was co-transferred 492 with synIX in yWZ612 and yWZ613, and 4) chrIX was lost in chromoductants following galactose 493 induction and counterselection against chrIX (yWZ618, yWZ619) (Figure 7B). However, we were surprised to observe that two of our chromoductants, yWZ610 and yWZ611, contained an extra 494 copy of chromosome *II* in addition to copies of both *chrIX* and *synIX* (Figure 7B). Whole genome 495 sequencing of these chromoductants confirmed these PFGE results (Figure 7C). However, further 496 497 examination of SNPs differing between the donor and recipient chromosome *II* revealed that

498 yWZ610 contained one donor *II* and one recipient *II*, while yWZ611 contained two copies of
499 recipient *II* (Figure 7D). This finding suggested that in yWZ611, even though chromosome *II* was
500 not co-transferred with *synIX*, the yeast still endoreduplicated chromosome *II*, perhaps to balance
501 the gene dosage perturbation caused by *synIX* chromoduction. Interestingly, during the inactivation
502 of *chrIX* in the yWZ610 and yWZ611 strains, chromosome *II* was simultaneously lost as well
503 (Figure 7B, yWZ618 and yWZ619).

#### 504

505 To investigate whether the co-transfer and simultaneous loss of chromosome IX and II was 506 coincidental, we next isolated another 11 disomic chromoductants from a repeated chromoduction 507 experiment. Interestingly, whole genome sequencing revealed that 10 of the 11 disomic 508 chromoductants gained an extra copy of chromosome *II* (Figure 7E). Following counterselection 509 against chrIX, four disomic chromoductants lost chromosome II, while six chromoductants 510 remained disomic (Figure 7F). To rule out the possibility that this phenomenon was specific to our 511 disomic synIX donor strain, we performed chromosome substitution using one alternative disomic synIX donor strain and three monosomic synIX donor strains. As in earlier experiments, 512 chromosome II was frequently co-transferred with synIX, and the extra copy of chromosome II 513 514 often disappeared after culturing strains in galactose (Table S8). In conclusion, we successfully 515 transferred synIX from different donor strains to a kar1-1 recipient strain via chromoduction. More in-depth diagnosis of the chromoductants revealed that chromosome II was preferentially and 516 consistently co-transferred with synIX, and then was frequently lost during counterselection against 517 518 chrIX.

519



520 521

522 **Figure 7.** Chromosome substitution of *synIX*. (**A**) Chromosome substitution occurs in two steps, chromoduction, in 523 which *synIX* is transferred to a recipient strain, followed by counterselection against wild-type *chrIX* to remove the 524 native version of the chromosome and restore normal chromosome copy number in the chromoductant strain. 525 Schematic diagram of substitution steps. Black bars, wild-type chromosomes. Blue bars, synIX. Recipient strain 526 harbors ade2 mutation, and thus forms pink colonies, for easily distinguishing chromoductants. (B) Pulsed-field gel 527 electrophoresis of donor strain, recipient strain, disomic chromoductants, and chromoductants after counterselection 528 against chrIX. BY4741 serves as control. ChrIX and synIX chromosomes indicated with black arrows. Green triangles 529 indicate synIX in disomic chromoductants. Red triangles indicate the extra chromosome II or chromosome III in 530 disomic chromoductants. (C) Whole genome sequencing read depth plots of disomic chromoductants (vWZ610, 531 vWZ612) and chrIX counterselected (i.e., fully substituted) chromoductants (vWZ618). X axis: Sixteen yeast 532 chromosomes, not drawn to scale. Y axis: Relative read depth, calculated as number of reads at each position divided 533 by average read depth across the sixteen yeast chromosomes. (D) Plot of SNPs on chromosome II in chromoductant 534 strains before and after counterselection against chrlX. Recipient strain SNP highlighted in orange, donor strain SNP in 535 blue. Strains: vWZ610: chromoductant strain with one copy of chromosome // from the donor strain and one from the 536 recipient strain (mixed SNP population); yWZ611: chromoductant strain with two copies of chromosome II from the 537 recipient strain (endoreduplicated //): vWZ618; strain derived from vWZ610 after counterselection against chr/X, with 538 one copy of II harboring recipient strain SNP; yWZ619: strain derived from yWZ611 after counterselection against 539 chrlX, with one copy of I/ harboring recipient strain SNP. (E) Relative depth plot for chromosome substitution strains 540 prior to counterselection against chrIX. X axis: Sixteen yeast chromosomes, not drawn to scale. Y axis: Relative read 541 depth. (F) Relative depth plot for chromosome substitution strains after counterselection against chrIX. X axis: Sixteen 542 yeast chromosomes, not drawn to scale. Y axis: Relative read depth.

543

#### 544 **Discussion**

545

546 We successfully conclude the synthesis, debugging, and characterization of synIX. SynIX's

547 developmental history reflects advances in DNA synthesis technology over more than a decade.

548 When the project first began, each segment of DNA was painstakingly assembled by

549 undergraduates from oligonucleotides into minichunks, which in turn were integrated iteratively into

550 the yeast genome to overwrite native chromosome *IX*. While such methods enabled successful

551 incorporation of much of the synthetic sequence of interest, drawbacks included the inadvertent

retention of wild-type patches and introduction of new point mutations. Newer approaches,

553 including direct minichunk synthesis and megachunk cloning and sequencing prior to integration,

554 greatly improved assembly efficiency and reduced construction error rates. Ultimately, we

555 produced a yeast strain with a fully synthetic version of chromosome *IX* that displays near wild-type 556 fitness.

557

558 While debugging synIX, we identified a fitness defect related to deletion of DNA sequence 559 upstream of the EST3 gene containing an aspartyl tRNA gene, a Ty1 LTR, and an intervening 560 single copy DNA sequence upstream of the EST3 gene. Strains with this deletion showed a reduction in telomere length as well as reduced Est3p levels and many RNA species mapping 561 562 upstream of EST3's normal transcription start site. Reinserting all or part of the deleted sequence was sufficient to restore normal strain growth fitness. However, only a subset of the edited strains, 563 specifically those containing the tRNA gene, showed a reduction in transcriptional activity 564 upstream of EST3 and a restoration of normal or near-normal Est3p levels and normal telomere 565 566 length. We hypothesize that the tRNA sequence normally provides a buffer against production of 567 aberrant transcripts upstream of EST3, and thus, the transcripts that are produced in tRNA-568 containing strains more effectively give rise to functional, full-length Est3p. Indeed, there is considerable evidence that tRNAs suppress transcription from nearby RNA polymerase II 569 promoters in a process referred to as tRNA-mediated gene silencing.<sup>33,35,36</sup> Consistent with this 570 571 idea, nanopore RNA sequencing data suggests that a subset of transcripts include one or more

AUGs upstream of EST3, and thus presumably fail to undergo proper translation; strains with an 572 abundance of such transcripts also appear to have fewer full-length, translatable transcripts and 573 lower Est3p levels. Other changes made upstream of EST3, such as reinsertion of the gap 574 575 sequence between the tRNA and the Tv1 LTR or Tv1 LTR itself, seemingly trigger small improvements in Est3p levels vs. those found in strains with the original synIX design. These 576 577 strains have intermediate length telomeres. longer than those seen in the original svnIX design but 578 shorter than those seen in native yeast, and maintain stable fitness after repeated passaging. We 579 speculate that the changes introduced to these strains may be sufficient to maintain Est3p 580 production and consequent telomere length above a threshold level needed for normal strain

- 581 fitness, despite not fully returning either to native levels.
- 582

583 Interestingly, a similar occurrence of abnormal transcriptional patterns resulting from tRNA deletion and *loxPsym* site insertion was previously observed at the *HIS2* locus in *synVI*.<sup>7</sup> In this case, the 584 585 authors suggested that a "cryptic start site" upstream of HIS2, generated by the altered sequence in this region, might yield a subset of transcripts incapable of producing functional His2 protein, or 586 that this transcription might impair normal HIS2 promoter activity. It would be useful to examine 587 588 whether common principles underlie multiple Sc2.0 bugs following tRNA deletion.

589

590 Finally, we demonstrated a chromosome substitution technique that combines selective

591 chromoduction, or movement of a desired chromosome between yeast strains, with subsequent 592 destabilization and loss of the corresponding native chromosome. During this process of

593 chromosome substitution, we found that certain chromosomes appeared to co-transfer with synIX

at higher frequencies than did others. Specifically, we saw multiple instances of transfer of 594 595 chromosome *II* or chromosome *III* alongside *svnIX*. These extra chromosomes often disappeared

following the destabilization of *chrIX*. Prior work suggested that smaller chromosomes might be 596

597 more likely to transfer during chromoduction than larger ones, and that co-transfer of multiple

chromosomes commonly occurs.<sup>23</sup> However, as chromosome *II* is not a small chromosome, it 598

- 599 seems unlikely that chromosome size fully explains the specific co-transfer of chromosome *II*.
- 600

601 It is likely that an extra copy of chromosome *II* plays a role in balancing a gene dosage imbalance arising from the transient extra copy of chromosome IX introduced by chromoduction. In one of the 602 603 chromoductant strains, analysis of SNPs from the recipient and donor strains' copies of 604 chromosome *II* revealed that both copies of the chromosome came from the recipient strain. The 605 presence of an endoreduplicated chromosome *II* is consistent with the hypothesis that there is 606 selective pressure for a higher dosage of chromosome *II* fitness in the presence of two copies of chromosome IX, and that the co-migration of chromosomes II and synIX in the other 607

chromoductant strains may not have occurred by chance. Moreover, multiple donor strains, 608

- 609 including strains disomic and monosomic for synIX, showed co-transfer of chromosome II
- alongside synIX. It may be valuable to ascertain the mechanisms that dictate which chromosomes 610

tend to co-transfer between strains during chromoduction. 611

612

613 Ultimately, chromosome substitution has numerous applications, including transfer and

614 consolidation of chromosomes from other wild-type yeast strains, engineered strains, and non-

615 laboratory strains. Existing work profiling collections of hundreds or thousands of yeast strains has

- 616 shed light on yeast diversity and evolution;<sup>37-39</sup> chromosome substitution will facilitate further
- 617 probing of how changes on individual chromosomes or combinations of chromosomes affect strain
- 618 fitness and function. Similarly, this technology could assist in analyzing why genomic changes in
- 619 SCRaMbLEd yeast strains and other engineered strains produce particular phenotypic results.<sup>1-4</sup>
- 620 More immediately, chromosome substitution is already in use to bring the initial phase of Sc2.0 to a
- 621 close via the consolidation of all of the synthetic chromosomes into a single strain, marking a key
- 622 eukaryotic genome engineering milestone.<sup>21</sup>

### 623 STAR Methods

624

### 625 Key resources table

REAGENT or RESOURCE	SOURCE	IDENTIFIER
Antibodies		
Digoxigenin Recombinant Rabbit Monoclonal Antibody (9H27L19)	Invitrogen	Thermo Fisher Scientific Cat# 700772, RRID:AB_2532342
Monoclonal ANTI-FLAG® M2 antibody produced in mouse	Sigma-Aldrich	Sigma-Aldrich Cat# F1804, RRID:AB_262044
Rabbit anti-Histone H3 antibody	Abcam	Abcam Cat# ab1791, RRID:AB_302613
IRDye 800CW Goat anti- Mouse IgG1-Specific	LI-COR Biosciences	LI-COR Biosciences Cat# 926-32350, RRID:AB_2782997
IRDye 680RD Goat anti- Rabbit IgG	LI-COR Biosciences	LI-COR Biosciences Cat# 926-68071, RRID:AB_10956166
Bacterial and viral strains		
EPI300 <i>E. coli</i> strain	Lucigen	EC300150
TOP10 <i>E. coli</i> strain	Boeke lab collection	
Chemicals, peptides, and r	ecombinant protei	ins
Zymolyase 20T	US Biological	Z1000
Zymolyase 100T	US Biological	Z1004
EDTA-free protease inhibitor	Roche	11873580001
Lithium acetate dihydrate	Sigma-Aldrich	L6883
Polyethylene glycol	Sigma-Aldrich	81188
Herring sperm DNA	Promega	D1816
Potassium acetate	Fisher	BP364
Zinc acetate dihydrate	Sigma-Aldrich	Z0625
(S)-(+)-Camptothecin	Sigma-Aldrich	C9911
D-Sorbitol	Sigma-Aldrich	S1876
6-Azauracil	Sigma-Aldrich	A1757
Hydroxyurea	Sigma-Aldrich	H8627
Methyl methanesulfonate	Sigma-Aldrich	129925
Methyl 1-(butylcarbamoyl)- 2-benzimidazolecarbamate (Benomyl)	Sigma-Aldrich	381586
Cycloheximide	Sigma-Aldrich	01810
L-Canavanine sulfate	Sigma-Aldrich	C9758
Hydrogen peroxide	Millipore	88597
ULTRAhyb <sup>™</sup> Ultrasensitive Hybridization Buffer	Invitrogen	AM8669
SSPE Buffer 20× Concentrate	Sigma-Aldrich	S2015
Intercept (TBS) Blocking Buffer	LI-COR Biosciences	927-60001
TWEEN 20	Sigma-Aldrich	P1379
Triton <sup>™</sup> X-100	Sigma-Aldrich	X100
2-Mercaptoethanol	Sigma-Aldrich	M6250
GoTaq Green Master Mix	Promega	M7123
RNase A	Thermo Fisher Scientific	EN0531
Notl-HF	NEB	R1389L
BbsI-HF	NEB	R3539L
Afel	NEB	R0652L
Xhol	NEB	R0146L

Dynabeds <sup>™</sup> M-270 Epox, Dynabeds <sup>™</sup> M-270 Epox, Dynabeds <sup>™</sup> M-270 Epox, Dynabeds <sup>™</sup> NORode Transcriptse Deads Critical commercial assays FureLink <sup>™</sup> HPure Plasmid Midprep Kit Zypo Plasmid Miniprep Kit Scientific Oubit dsDNA HS Assay Kit Inetron Fisher Scientific Oubit dsDNA HS Assay Kit Deads Zymo-Spin I Columns Zymo-Spin I Columns Zymo-Zymo-Zymo-Zymo Zymo-Zymo Zymo-Zymo Zymo Zymo Zymo Zymo Z	Amino-terminal FLAG- BAP™ Fusion Protein	Sigma-Aldrich	P7582					
Dynabeads <sup>1M</sup> Oligo(dT) <sub>E5</sub> Jinvitrogen 18090010 Agencourt RNAClean XP Beckman Coulter A63987 Addiger Commercial assays PreLink <sup>TM</sup> HFUre Plasmid Invitrogen 463987 Critical commercial assays PreLink <sup>TM</sup> HFUre Plasmid Invitrogen 463987 Coulter 27300 Coulter 27300 Coulter 27300 Coulter 27300 Coult dsDNA HS Assay Kit Scientific 032854 Coult dsDNA HS Assay Kit Thermo Fisher Scientific 032854 Coult dsDNA HS Assay Kit Thermo Fisher Scientific 032850 NEBNext Ultra II FS DNA NEB E7805 Zymo Research 1003 Pall 60207 Minitor Part 100 Coulter 2700 Coulter 2700 Coulter 2700 Coulter 2700 Coult dsDNA BR Assay Kit Thermo Fisher Scientific 032850 Coult dsDNA BR Assay Kit Thermo Fisher Scientific 032850 Coulter 2700 Coulter 2700 Coulter 2700 Coulter 2700 Coulter 2700 Coulter 2700 Coulter 2700 Coulter 2700 Coult dsDNA BR Assay Kit Coult dsDNA BR Assay Kit Coult RNA BR Assay Kit ColaGEN 74106 Coult RNA BR Assay Kit ColaGEN 180745 Coult RNA BR Assay Kit ColaGEN 28050 Coloci ColaGEN 180745 ColaGEN 28050 Coloci ColaGEN 180745 ColaGEN	Dynabeads™ M-270 Epoxy	Invitrogen	14302D					
SuperScript M V Reverse invitrogen       18090010         Agencourt RNAClean XP beckman could invitrogen       A63967         Contler       A63967         Contler       Norgen       K2100-15         Midprep Kit       Zymo Research       D4037         Fundlink <sup>TM</sup> HiPure Plasmid       Invitrogen       K2100-15         Gubit dsDNA HS Assay Kit       Thermo Fisher Scientific       Q32854         Gubit dsDNA BR Assay Kit       Thermo Fisher Scientific       Q32850         Symo-Spin I Columns       Zymo Research       C1003         PALL 00207 Biodyne B       High Sensitivity and Low Background Nyton Transfer       Pall Corporation         Background Nyton Transfer       Pall Corporation       60207         RNeasy Mini Kit       OIAGEN       74106         RNeasy Mini Kit       OIAGEN       74106         RNeasy Mini Kit       OIAGEN       34217         NextSeq 500/550 High Output Kit V2 5 (75 Cycles)       Illumina       20024907         Output Kit v2 5 (75 Cycles)       Illumina       20024907         NastSer 2000 S00 High Output Kit V2 5 (75 Cycles)       Illumina       20024907         Oxford       Nanopore       Scientific       Q32852         Direct RNA Sequencing Kit       Thermo Fisher Scientific       Q	Dynabeads™ Oligo(dT) <sub>25</sub>	Invitrogen	61005					
Transcriptase     Introduction       Agencourt RNAClean XP     Beckman Coller       Agencourt RNAClean XP     Beckman Coller       Critical commercial assays     Invitrogen       Cubit dsDNA HS Assay Kit     Thermo Fisher Scientific     27300       Qubit dsDNA HS Assay Kit     Thermo Fisher Scientific     232854       Qubit dsDNA HS Assay Kit     Thermo Fisher Scientific     Q32850       Zymo-Spin I Columns     Zymo Research     C1003       PALL 6020 Flootyne     Pall Corporation     60207       Beakground Nyion Transfer Beakground Nyion Transfer     Pall Corporation     60207       RNeagy Mini Kit     QIAGEN     74106       RNeagy Mini Kit     QIAGEN     180745       Qubit RNA BR Assay Kit     QIAGEN     180745       RNaagy Stranded Total QIASeq Stranded Total QIASeq Stranded Total QIAGEN     QIAGEN     334217       Veast Kit     Color     20024906       Nety St 150     Illumina     20024907       Output Kit 25 (150     Coford     32862       Qiaseq Standed Total Qiater Nti     Agilent     5067.1511       Quark RNA Sasay Kit     Conford     20024906       Next Seg 500/550 High Output Kit 25 (150     QiAGEN     32862       Qiater Nti     Agilent     5067.1511       Quark RNA Sasay Kit     Coford </td <td>SuperScript™ IV Reverse</td> <td>Invitrogen</td> <td>18090010</td>	SuperScript™ IV Reverse	Invitrogen	18090010					
Agencourt RNAClean XP         Beckman         A63987           Coulter         Coulter         A63987           VerLink <sup>TM</sup> HiPure Plasmid         Invitrogen         K2100-15           Zyppy Plasmid Miniprep Kit         Zymo Research         D4037           Isolation Kit         Orgen         Z3300           Ubit dsDNA HS Assay Kit         Thermo Fisher Scientlic         Q32854           Qubit dsDNA BR Assay Kit         Thermo Fisher Scientlic         Q32850           Zymo-Spin I Columns         Zymo Research         C1003           PALL 60207 Biodyne B         High Sensitivity and Low Background Nylon Transfer Membrane, 0.45 µm Pores         Pall Corporation         60207           Size, 30 on W X 3 m L Roll         QIAGEN         74106         Columna           Qubit RNA BR Assay Kit         Thermo Fisher Scientlific         Q10210         Columna           Qubit RNA BR Assay Kit         Thermo Fisher Scientlific         Q10210         Columna           Qubit RNA BR Assay Kit         UAGEN         180745         QualceN           Rive Assay Kit         UAGEN         334217         NextSeq 500/550 High           Output Kit v2.5 (75 Cycles)         Illumina         20024906         QualceN           NextSeq 500/550 High         QualceN         30567-1511	Transcriptase	пиподен	1000010					
beads         Coller         Process           Critical commercial assays         Invitrogen         K2100-15           Zyppy Plasmid Miniprep Kit         Zymo Research         D4037           Fungi/Yeast Genomic DNA Isolation Kit         Norgen         27300           Qubit dsDNA HS Assay Kit         Thermo Fisher Scientific         Q32854           Qubit dsDNA BR Assay Kit         Thermo Fisher Scientific         Q32850           Zymo-Spin I Columns         Zymo Research         C1003           PALL 60207 Biodyne B         High Sensitivity and Low Background Nyion Transfer Membrane, 0.45 µm Pore         60207           Resay Min Kit         QIAGEN         74106            Qubit RNA BR Assay Kit         ClaGEN         74106            Qubit RNA BR Assay Kit         QIAGEN         74106            Qubit RNA BR Assay Kit         QIAGEN         74106            Qubit RNA BR Assay Kit         QIAGEN         180745            Qubit RNA BR Assay Kit         QIAGEN         334217            VestSet 500/550 High Output Kit V2.5 (75 Cycles)         Illumina         20024906            NetSeq 500/550 High Output Kit V2.5 (75 Cycles)         Illumina         20024907	Agencourt RNAClean XP	Beckman	A63987					
Critical commercial assays           PurclInK <sup>®</sup> HiPure Plasmid Midgrep Kit         Invitrogen         K2100-15           Zyppy Plasmid Miniprep Kit         Zymo Research         D4037           FunglyVeast Genomic DNA Isolation Kit         Norgen         27300           Qubit dsDNA HS Assay Kit         Thermo Fisher Scientific         Q32854           Qubit dsDNA BR Assay Kit         Thermo Fisher Scientific         Q32850           Xymo-Spin I Columns         Zymo Research         C1003           PALL 60207 Biodyne B High Sensitivity and Low Background Nyion Transfer Membrane, 0.45 µm Pore         Pall Corporation         60207           Rembrane, 0.45 µm Pore         QIAGEN         74106            Qubit RNA BR Assay Kit         Thermo Fisher Scientific         Q10210            Qubit RNA BR Assay Kit         QIAGEN         74106            Renbarge Stranded Total RNA Library Kit         QIAGEN         74106            Qubit RNA BR Assay Kit         QIAGEN         334217            Veasts of 500/550 High         QIUMINIA         Q0024906            Output Kit V2 5 (750 Cycles)         Illumina         20024907            Qalent RNA R0000 Nano Kit         Agilent RNA 6000 Nano Kit         Galent RNA 6000 Nano Kit	beads	Coulter						
PureLink"HiPure Plasmid Midiprep KitInvitrogenK2100-15Zyppy Plasmid Minjprep KitZymo ResearchD4037FunglYceast Genomic DNA Isolation KitNorgen27300Qubit dsDNA HS Assay KitThermo Fisher ScientificQ32850Qubit dsDNA BR Assay KitThermo Fisher ScientificQ32850Zymo ResearchC1003PALL 60207PALL 60207 Biodyne BFisher ScientificG0207Alta GozonaPall Corporation60207Beakground Nykon Transfer Membrane, 0.45 µm ProfPall CorporationRNeagy Mini KitQIAGEN74106Qubit RNA BR Assay KitThermo Fisher ScientificQ10210Qubit RNA BR Assay KitCIAGEN180745Qubit RNA BR Assay KitCIAGEN334217NextSeq 500/550 High Output Kit V2.5 (150Illumina20024906NextSeq 500/550 High Output Kit V2.5 (150Illumina20024907Oydelon KitLucigenMPY03100Agilent RNA 6000 Nano Kit AgilentAgilent5067-1511Qubit RNA HS Assay KitThermo Fisher ScientificQ32852Oxford MasterPure Yeast RNA Purification KitNanopore ScientificQ32852Direct RNA Sequencing KitNanopore TechnologiesSGK-RNA002Flow Cell Priming KitNanopore TechnologiesSGK-RNA002Flow Cell Priming KitNanopore TechnologiesSGK-RNA002Flow Cell Priming KitNanopore TechnologiesFLO-MIN106DMinION Flow Cell (R9.4.1)Nanopor	Critical commercial assays	1						
Midprep Kit         Zymo Research         D4037           Fung/Yeast Genomic DNA Isolation Kit         Norgen         27300           Qubit dsDNA HS Assay Kit         Thermo Fisher Scientific         Q32854           Qubit dsDNA BR Assay Kit         Thermo Fisher Scientific         Q32850           NEBNext Ultra II FS DNA Library Kit         NEB         E7805           Zymo-Spin I Columns         Zymo Research         C1003           PALL 60207 Biodyne B High Sensitivity and Low Background Nylon Transfer Membrane, 0.45 µm Pore Scientific         Q10210           Scientific         Q10210         60207           Qubit RNA BR Assay Kit         ClaGEN         74106           RNeasy Mini Kit         QIAGEN         74106           Qubit RNA BR Assay Kit         Thermo Fisher Scientific         Q10210           QIAseq FastSelect - RNA Yeast Kit         QIAGEN         334217           Veast Kit         20024906         WHY03100           Output Kit V2.5 (75 Cycles)         Illumina         20024907           Oyder Kit V2.5 (1500         Illumina         20024907           Oyder Kit V2.5 (1500         Illumina         20024907           Oyder Kit V2.5 (1500         Illumina         20024907           Oydord         Agilent Thermo Fisher Scientific <t< td=""><td>PureLink<sup>™</sup> HiPure Plasmid</td><td>Invitrogen</td><td colspan="5">K2100-15</td></t<>	PureLink <sup>™</sup> HiPure Plasmid	Invitrogen	K2100-15					
Zyppy Plasmid Winipep Xii     Zymin Research     D4037       FungiYeast Genomic DNA Isolation Kit     Norgen     27300       Qubit dsDNA HS Assay Kit     Thermo Fisher Scientific     Q32854       Qubit dsDNA BR Assay Kit     Thermo Fisher Scientific     Q32850       NEBNext Ultra II FS DNA Library Kit     NEB     E7805       Zymo-Spin I Columns     Zymo Research     C1003       PALL 60207 Biodyne B     Pall Corporation     60207       Membrane, O.45 µm Pore Size, 30 cm W x 3 m L Roll     Pall Corporation     60207       Qubit RNA BR Assay Kit     Thermo Fisher Scientific     Q10210       QlAseq Stranded Total RNA Library Kit     QIAGEN     180745       QlAseq Stranded Total RNA Library Kit     QIAGEN     334217       VestSeq 500/550 High Output Kit v2.5 (150 Cycles)     Illumina     20024906       NetSeq 500/550 High Output Kit v2.5 (150 Cycles)     Illumina     20024907       RasterPure Yeast RNA Purification Kit     Lucigen     MPY03100       Qubit RNA HS Assay Kit     Thermo Fisher Scientific     Q32852       Oxford     Nanopore Scientific     Q32862       Direct RNA 6000 Nano Kit     Agilent     5067-1511       Qubit RNA HS Assay Kit     Thermo Fisher Scientific     Q32862       Oxford     Nanopore Scientific     SQK-RNA002       Flow Cell Primi	Midiprep Kit	Zurran Dagagarah	D4027					
Pullight Teals Cellfolm DMA     Norgen     27300       Qubit dsDNA HS Assay Kit     Thermo Fisher Scientific     Q32854       Qubit dsDNA BR Assay Kit     Thermo Fisher Scientific     Q32850       NEBNext Ultra II FS DNA Library Kit     NEB     E7805       Zymo-Spin I Columns     Zymo Research     C1003       PALL 60207 Biodyne B     High Sensitivity and Low     60207       Membrane. 0.45 µm Pore     Pall Corporation     60207       Membrane. 0.45 µm Pore     Pall Corporation     60207       Membrane. 0.45 µm Pore     Pall Corporation     60207       Qubit RNA BR Assay Kit     Thermo Fisher Scientific     Q10210       QiAseg Stranded Total     QIAGEN     180745       QiAseg FastSelect -rRNA Yeast Kit     QIAGEN     334217       NextSeq 500/550 High Output Kit v2.5 (150 cycles)     Illumina     20024906       Netty V2.5 (150 cycles)     Illumina     20024907       Cycles)     MasterPure Yeast RNA Lucigen     MPY03100       Agilert RNA 6000 Nano Kit     Agilert     S067-1511       Qubrit RNA HS Assay Kit     Thermo Fisher Scientific     Q32852       Oxford     Nanopore Technologies     SQK-RNA002       Performant Kit     Agilert     Socramolgen       Quirt RNA HS Assay Kit     Thermo Fisher Scientific     Q32852	Zyppy Plasmid Miniprep Kit	Zymo Research	D4037					
Isolation Nut       Thermo Fisher Scientific       Q32854         Qubit dsDNA BR Assay Kit       Thermo Fisher Scientific       Q32850         DEBNext Utra II FS DNA Library Kit       NEB       E7805         Zymo-Spin I Columns       Zymo Research       C1003         PALL 60207 Biodyne B       60207         Background Nylon Transfer Membrane, 0.45 µm Pore Size, 30 cm W x 3 m L Roll       Pall Corporation         Qubit RNA BR Assay Kit       QIAGEN       74106         Qubit RNA BR Assay Kit       QIAGEN       74106         Qubit RNA BR Assay Kit       QIAGEN       74106         Qubit RNA BR Assay Kit       QIAGEN       180745         Qubit RNA BR Assay Kit       QIAGEN       180745         Qubit RNA BR Assay Kit       QIAGEN       134217         NextSeq 500/550 High Output Kit v2.5 (75 cycles)       Illumina       20024906         NextSeq 500/550 High Output Kit v2.5 (150       Illumina       20024907         Cycles)       MPY03100       MPY03100         Agilent RNA 6000 Nano Kit       Agilent       5067-1511         Agilent RNA 6000 Nano Kit       Agilent       5067-1511         Agilent RNA HS Assay Kit       Thermo Fisher Scientific       032852         Oxford       Nanopore       SQAK-RNA002     <	Fuligit reast Genomic DNA	Norgen	27300					
Qubit dsDNA HS Assay KitColumn in the ScientificQ32854Qubit dsDNA BR Assay KitThermo Fisher ScientificQ32850NEBNext Ultra II FS DNA Library KitNEBE7805Zymo-Spin I ColumnsZymo ResearchC1003PALL 60207 Biodyne B High Sensitivity and Low Background Nylon Transfer Membrane, 0.45 µm Pore ScientificPall CorporationRNeasy Mint KitQIAGEN74106Qubit RNA BR Assay KitCIAGEN74106Qubit RNA BR Assay KitQIAGEN180745Qubarg Stranded Total RNA Library KitQIAGEN180745Qubarg FastSelect -rRNA Veast KitQIAGEN334217NextSeq 500/550 High Output Kit v2.5 (150 Cycles)Illumina20024906Nubrary Kit Quiput Kit v2.5 (150 Cycles)Illumina20024907Output Kit v2.5 (150 Cycles)Illumina20024907Qubit RNA Bodonano Kit Agilent RNA 6000 Nano Kit ScientificQ32852Q32852Oxford Direct RNA Sequencing KitNanopore TechnologiesSQK-RNA002Flow Cell Priming KitOxford Nanopore TechnologiesSQK-RNA002Flow Cell Priming KitOxford Nanopore TechnologiesSQK-RNA002Fisher ScientificOxford Nanopore TechnologiesSQK-RNA002Fisher Sicentifier Such Tither Sicentifier SicentifierN/AAll other strains used in this study are listed in Table S1N/ANinON Flow Cell (R9.4.1)NiAN/A		Thermo Fisher						
Qubit dsDNA BR Assay Kit         Thermo Fisher Scientific         Q32850           NEBNext Ultra II FS DNA Library Kit         NEB         E7805           Zymo-Spin I Columns         Zymo Research         C1003           PALL 60207 Biodyne B High Sensitivity and Low Background Nylon Transfer Membrane, 0.45 µm Pore Size, 30 cm W x 3 m L Roll         Pall Corporation         60207           RNeasy Mini Kit         QIAGEN         74106         01210           Qubit RNA BR Assay Kit Qiage Stranded Total RNA Library Kit         QIAGEN         180745           QIAGEN Science         QIAGEN         334217           NextSeq 500/550 High Output Kit v2.5 (75 Cycles)         Illumina         20024906           NextSeq 500/550 High Output Kit v2.5 (150 Cycles)         Illumina         20024907           Qycles)         Master/Fure Yeast RNA Purification Kit         Aglient         5067-1511           Qubit RNA HS Assay Kit         Thermo Fisher Scientific         Q32852         Oxford           Qubit RNA HS Assay Kit         Thermo Fisher Scientific         Q32852         Oxford           Master/Fure Yeast RNA Purification Kit         Aglient         5067-1511         Thermo Fisher Scientific         Q32852           Direct RNA Sequencing Kit         Nanopore Technologies         EXP-FLP001         Technologies           Minopore Tech	Qubit dsDNA HS Assay Kit	Scientific	Q32854					
Qubit dsDNA BR Assay KitScientificQ32850NEBNext Ultra II FS DNA Library KitNEBE7805Zymo-Spin I ColumnsZymo ResearchC1003PALL 60207 Biodyne B High Sensitivity and Low Background Nylon Transfer Membrane, 0.45 µm PorePall CorporationRNeasy Mini KitQIAGEN74106Qubit RNA BR Assay KitClAGEN74106Qubit RNA BR Assay KitQIAGEN180745Qubit RNA BR Assay KitQIAGEN180745QIAseq Stranded TotalQIAGEN334217Yeast KitQIAGEN334217Vesat KitUlumina20024906NextSeq 500/550 High Output Kit V2.5 (75 Cycles)Illumina20024907VestSeq 500/550 High Output Kit V2.5 (75 Cycles)Illumina20024907MasterPure Yeast RNA Purification KitLucigenMPY03100Aglient RNA 6000 Nano Kit Aglient RNA Sequencing KitAglient5067-1511Qubit RNA HS Assay KitThermo Fisher ScientificQ32852Direct RNA Sequencing KitAnopore TechnologiesSQK-RNA002Flow Cell Priming KitOxford Nanopore TechnologiesSQK-RNA002Flow Cell Priming KitJef Boeke's LaboratoryJef Boeke's LaboratoryN/AAll other strains used in this study are listed in Table S1N/AN/AAll other strains used in this Study are listed in Table S1N/AN/A		Thermo Fisher						
NEBNext Ultra II FS DNA Library Kit       NEB       E7805         Zymo-Spin I Columns       Zymo Research       C1003         PALL 60207 Biodyne B High Sensitivity and Low Background Nylon Transfer Membrane, 0.45 µm Pore Size, 30 cm W x 3 m L Roll       Pall Corporation       60207         RNeasy Mini Kit       QIAGEN       74106       74106         Qubit RNA BR Assay Kit       Thermo Fisher Scientific       Q10210       0         QIAseq Stranded Total RNA Library Kit       QIAGEN       180745       0         QIAseq Solo(550 High Output Kit v2.5 (75 Cycles)       Illumina       20024906       0         NextSeq 500/550 High Output Kit v2.5 (150       Illumina       20024907       0         Yedset Kit       Lucigen       MPY03100       0         Agilent RNA 6000 Nano Kit       Agilent       5067-1511       0         Qubit RNA HS Assay Kit       Thermo Fisher Scientific       Q32852       Q32852         Oxford       Nanopore Technologies       SQK-RNA002       0         Flow Cell Priming Kit       Oxford Nanopore Technologies       SQK-RNA002       0         Sactoaromyces cerevisiae: BY4741       Jef Boeke's laboratory       N/A       N/A         All other strains used in this study are listed in Table S1       N/A       N/A	Qubit dsDNA BR Assay Kit	Scientific	Q32850					
Library KitNEBE / 805Zymo-Spin I ColumnsZymo ResearchC1003PALL 60207 Biodyne BHigh Sensitivity and Low60207Background Nylon Transfer Membrane, 0.45 µm Pore Size, 30 cm W x 3 m L RollPall Corporation60207RNeasy Mini KitOIAGEN74106Qubit RNA BR Assay KitThermo Fisher ScientificQ10210QlAseq Stranded Total RNA Library KitQIAGEN180745QlAseq Strasselect -rRNA Yeast KitQIAGEN334217NextSeq 500/550 High Output Kit v2.5 (75 Cycles)Illumina20024906NextSeq 500/550 High Output Kit v2.5 (75 Cycles)Illumina20024907CyclesUrigenMPY03100Purification KitLucigenMPY03100Agilent RNA 6000 Nano Kit Purification KitAgilent5067-1511Qubit RNA HS Assay KitThermo Fisher ScientificQ32852Direct RNA Sequencing KitNanopore TechnologiesSQK-RNA002Flow Cell Priming KitOxford Nanopore TechnologiesSQK-RNA002Flow Cell Priming KitJef Boeke's Laborapet StrainsJef Boeke's Laborapet StrainsSaccharomyces cerevisiae: BY4711Jef Boeke's LaborapetN/AAll other strains used in this study are listed in Table S1N/AN/AN/AN/A	NEBNext Ultra II FS DNA		F7005					
Zymo-Spin L ColumnsZymo ResearchC1003PALL 60207 Biodyne B High Sensitivity and Low Background Nyton Transfer Membrane, 0.45 µm Pore60207Rembrane, 0.45 µm Pore Size, 30 cm W x 3 m L RollQIAGEN74106RNeasy Mini KitQIAGEN74106Qubit RNA BR Assay KitThermo Fisher ScientificQ10210QIAseq Stranded Total Weats KitQIAGEN180745QIAseq Stranded Total Weats KitQIAGEN334217NextSeq 500/550 High Output Kit v2.5 (15 Cycles)Illumina20024906NextSeq 500/550 High Output Kit v2.5 (150Illumina20024907Output Kit v2.5 (150 Uutput Kit v2.5 (150Illumina20024907Output Kit v2.5 (150 Uutput Kit v2.5 (150Illumina20024907Output Kit v2.5 (150 Uutput Kit v2.5 (150 Cycles)MPY03100 VO24907Direct RNA Sequencing Kit Nanopore TechnologiesSQK-RNA002 VAford Nanopore TechnologiesSQK-RNA002 VAfordFlow Cell Priming Kit Saccharomyces cerevisiae: BY4741Jef Boeke's IaboratoryN/AAll other strains used in this study are listed in Table S1N/AN/AOliconculectionedN/A<	Library Kit	NEB	E7805					
PALL 60207 Biodyne B       Bigh Sensitivity and Low       Background Nylon Transfer         Background Nylon Transfer       Pall Corporation       60207         RNeasy Mini Kit       QIAGEN       74106         Qubit RNA BR Assay Kit       Thermo Fisher       Colored Colore	Zymo-Spin I Columns	Zymo Research	C1003					
High Sensitivity and Low Background Nylon Transfer Membrane, 0.45 µm Pore Size, 30 cm W x 3 m L RollPall Corporation60207RNeasy Mini KitOLAGEN74106Qubit RNA BR Assay KitThermo Fisher ScientificQ10210QIAseq Stranded Total RNA Library KitQIAGEN180745QIAseq FastSelect -rRNA Yeast KitQIAGEN334217NextSeq 500/550 High Output Kit v2.5 (75 Cycles)Illumina20024906NextSeq 500/550 High Output Kit v2.5 (150Illumina20024907Qualer RNA 6000 Nano Kit Agilent RNA 6000 Nano Kit Agilent RNA 6000 Nano Kit Nanopore TechnologiesMPY03100Oxford Nanopore TechnologiesQ32852Flow Cell Priming KitNanopore TechnologiesSQK-RNA002Flow Cell Priming KitNanopore TechnologiesEXP-FLP001Flow Cell R9.4.1)Manopore TechnologiesFLO-MIN106DFlow Cell R9.4.1)Manopore TechnologiesFLO-MIN106DFlow Cell R9.4.1)Manopore TechnologiesFLO-MIN106DTechnologiesManopore TechnologiesKYrdrMinION Flow Cell (R9.4.1)MAN/AMil other strains used in this study are listed in Table S1N/AN/AN/AN/A	PALL 60207 Biodyne B	-						
Background Nylon Transfer Membrane, 0.45 µm Pore Size, 30 cm W x 3 m L RollPall Corporation60207RNeasy Mini KitQIAGEN74106Qubit RNA BR Assay KitThermo Fisher ScientificQ10210QIAseq Stranded Total RNA Library KitQIAGEN180745QIAseq FastSelect -rRNA Yeast Kit V 2.5 (75 Cycles)QIAGEN334217NextSeq 500/550 High Output Kit v2.5 (75 Cycles)Illumina20024906Output Kit v2.5 (150 Cycles)Illumina20024907Output Kit v2.5 (150 Cycles)Illumina20024907Output Kit v2.5 (150 Cycles)Illumina20024907Output Kit NA HS Assay KitThermo Fisher ScientificQ32652Output Kit NA HS Assay KitThermo Fisher ScientificQ32652Oxford TechnologiesSQK-RNA002Pierct RNA Sequencing KitNanopore TechnologiesSQK-RNA002Flow Cell Priming KitOxford Nanopore TechnologiesSQK-RNA002Flow Cell R9.4.1)Oxford Nanopore TechnologiesNanopore TechnologiesSaccharomyces cerevisie: Stydraft 1Jef Boeke's IaboratoryN/AAll other strains used in this study are listed in Table S1M/AN/A	High Sensitivity and Low							
Membrane, 0.45 µm Pore Size, 30 cm W x 3 m L Roll       View of the second	Background Nylon Transfer	Pall Corporation	60207					
Size, 30 cm W x 3 m L Roll       QIAGEN         Rneasy Mini Kit       QIAGEN         Qubit RNA BR Assay Kit       Thermo Fisher Scientific       Q10210         QIAseq Stranded Total RNA Library Kit       QIAGEN       180745         QIAseq FastSelect -rRNA Yeast Kit       QIAGEN       334217         NextSeq 500/550 High Output Kit v2.5 (15 Cycles)       Illumina       20024906         NextSeq 500/550 High Output Kit v2.5 (150       Illumina       20024907         Output Kit v2.5 (150       Illumina       20024907         Cycles)       Lucigen       MPY03100         Purification Kit       Agilent       5067-1511         Qubit RNA HS Assay Kit       Scientific       Q32852         Direct RNA Sequencing Kit       Agilent       5067-1511         Direct RNA Sequencing Kit       Nanopore Technologies       SQK-RNA002         Flow Cell Priming Kit       Oxford Nanopore       SQK-RNA002         MinION Flow Cell (R9.4.1)       Oxford Nanopore       FLO-MIN106D         Technologies       Scharomyces carevisiae: BY4741       Jef Boeke's laboratory         Sydraf       N/A       N/A         All other strains used in this study are listed in Table S1       N/A	Membrane, 0.45 µm Pore							
RNeasy Mini Kit       QIAGEN       74106         Qubit RNA BR Assay Kit       Thermo Fisher Scientific       Q10210         QIAseq Stranded Total RNA Library Kit       QIAGEN       180745         QIAseq FastSelect -rRNA Yeast Kit       QIAGEN       334217         NextSeq 500/550 High Output Kit v2.5 (75 Cycles)       Illumina       20024906         NextSeq 500/550 High Output Kit v2.5 (150 Cycles)       Illumina       20024907         NextSeq 500/550 High Output Kit v2.5 (150       Illumina       20024907         Cycles)       MasterPure Yeast RNA Purification Kit       Lucigen       MPY03100         Agilent RNA 6000 Nano Kit       Agilent       5067-1511         Qubit RNA HS Assay Kit       Thermo Fisher Scientific       Q32852         Direct RNA Sequencing Kit       Nanopore Technologies       SQK-RNA002         Flow Cell Priming Kit       Oxford Nanopore Technologies       SQK-RNA002         MinION Flow Cell (R9.4.1)       Oxford Nanopore Technologies       SQK-RNA002         Experimental models: Organisms/strains       Saccharomyces cerevisiae: Jef Boeke's BY4741       Jaboratory         All other strains used in this study are listed in Table S1       N/A       N/A	Size, 30 cm W x 3 m L Roll							
Qubit RNA BR Assay KitIhermo Fisher ScientificQ10210QIAseq Stranded Total RNA Library KitQIAGEN180745QIAseq FastSelect -rRNA Yeast KitQIAGEN334217NextSeq 500/550 High Output Kit v2.5 (150Illumina20024906Output Kit v2.5 (150Illumina20024907Cycles)LucigenMPY03100MasterPure Yeast RNA Purification KitLucigenMPY03100Qubit RNA HS Assay KitThermo Fisher ScientificQ32852Output Kit v2.5 (150Nanopore TechnologiesQ32852Direct RNA Sequencing KitNanopore TechnologiesSQK-RNA002Flow Cell Priming KitNanopore TechnologiesSQK-RNA002Flow Cell R9.4.1)Nanopore TechnologiesSQK-RNA002Flow Cell R9.4.1)Nanopore TechnologiesN/AMinION Flow Cell (R9.4.1)N/AN/AAll other strains used in this study are listed in Table S1N/AAll other strains used in this study are listed in Table S1N/A	RNeasy Mini Kit		74106					
QIAseq Stranded Total RNA Library KitQIAGEN180745QIAseq FastSelect -rRNA Yeast KitQIAGEN334217NextSeq 500/550 High Output Kit v2.5 (75 Cycles)Illumina20024906NextSeq 500/550 High Output Kit v2.5 (150 Cycles)Illumina20024907MasterPure Yeast RNA Purification KitLucigenMPY03100Agilent RNA 6000 Nano Kit Agilent RNA 6000 Nano KitAgilent5067-1511Qubit RNA HS Assay KitThermo Fisher ScientificQ32852Direct RNA Sequencing KitNanopore TechnologiesQXford NanoporeFlow Cell Priming KitOxford Nanopore TechnologiesEXP-FLP001MinION Flow Cell (R9.4.1)Nanopore TechnologiesFLO-MIN106DExperimental models: Organisms/strainsSaccharomyces cerevisiae: IaboratoryN/AAll other strains used in this study are listed in Table S1N/AN/A	Qubit RNA BR Assay Kit	Thermo Fisher Scientific	Q10210					
RNA Library Kit     Carlocitic     From the formation of	QIAseq Stranded Total		180745					
QIAseq FastSelect -rRNA Yeast Kit     QIAGEN     334217       NextSeq 500/550 High Output Kit v2.5 (75 Cycles)     Illumina     20024906       NextSeq 500/550 High Output Kit v2.5 (150     Illumina     20024907       Oycles)     20024907     Cycles)     20024907       MasterPure Yeast RNA Purfication Kit     Lucigen     MPY03100       Agilent RNA 6000 Nano Kit     Agilent     5067-1511       Qubit RNA HS Assay Kit     Thermo Fisher Scientific     Q32852       Direct RNA Sequencing Kit     Nanopore Technologies     SQK-RNA002       Flow Cell Priming Kit     Nanopore Technologies     EXP-FLP001       MinION Flow Cell (R9.4.1)     Nanopore Technologies     FLO-MIN106D       System     Floeke's Iaboratory     N/A       All other strains used in this study are listed in Table S1     N/A     N/A	RNA Library Kit							
Yeast Kit       Illumina       20024906         NextSeq 500/550 High Output Kit v2.5 (75 Cycles)       Illumina       20024907         NextSeq 500/550 High Output Kit v2.5 (150       Illumina       20024907         MasterPure Yeast RNA Purification Kit       Lucigen       MPY03100         Agilent RNA 6000 Nano Kit       Agilent       5067-1511         Qubit RNA HS Assay Kit       Thermo Fisher Scientific       Q32852         Direct RNA Sequencing Kit       Nanopore Technologies       SQK-RNA002         Flow Cell Priming Kit       Oxford       SQK-RNA002         MinION Flow Cell (R9.4.1)       Oxford       KXP-FLP001         Experimental models: Organisms/strains       Saccharomyces cerevisiae: Jef Boeke's laboratory       N/A         All other strains used in this study are listed in Table S1       N/A       N/A	QIAseq FastSelect -rRNA	QIAGEN	334217					
NextSeq 500/550 High Output Kit v2.5 (75 Cycles)Illumina20024906NextSeq 500/550 High Output Kit v2.5 (150 Cycles)Illumina20024907MasterPure Yeast RNA Purification KitLucigenMPY03100Agilent RNA 6000 Nano KitAgilent5067-1511Qubit RNA HS Assay KitThermo Fisher ScientificQ32852Direct RNA Sequencing KitNanopore TechnologiesSQK-RNA002Flow Cell Priming KitOxford Nanopore TechnologiesEXP-FLP001MinION Flow Cell (R9.4.1)Oxford Nanopore TechnologiesFLO-MIN106DExperimental models: Organisms/strains Saccharomyces cerevisiae:Jef Boeke's laboratoryN/AAll other strains used in this study are listed in Table S1N/AN/A	Yeast Kit							
Output Kit V2.5 (75 Cycles)       20024907         NextSeq 500/550 High Output Kit V2.5 (150       Illumina       20024907         Cycles)       MasterPure Yeast RNA Purification Kit       Lucigen       MPY03100         Agilent RNA 6000 Nano Kit       Agilent       5067-1511         Qubit RNA HS Assay Kit       Thermo Fisher Scientific       Q32852         Direct RNA Sequencing Kit       Nanopore Technologies       SQK-RNA002         Flow Cell Priming Kit       Nanopore Technologies       EXP-FLP001         MinION Flow Cell (R9.4.1)       Oxford Nanopore Technologies       FLO-MIN106D         Saccharomyces cerevisiae: BY4741       Jef Beke's Iaboratory       N/A         All other strains used in this study are listed in Table S1       N/A       N/A	Output Kit v2 5 (75 Cycles)	Illumina	20024906					
Nexted Situ 22.5 (150Illumina20024907Output Kit v2.5 (150Illumina20024907MasterPure Yeast RNA Purification KitLucigenMPY03100Agilent RNA 6000 Nano KitAgilent5067-1511Qubit RNA HS Assay KitThermo Fisher ScientificQ32852Direct RNA Sequencing KitNanopore TechnologiesSQK-RNA002Flow Cell Priming KitOxford Nanopore TechnologiesEXP-FLP001MinION Flow Cell (R9.4.1)Oxford Nanopore TechnologiesFLO-MIN106DExperimental models: Organisms/strainsJef Boeke's laboratoryN/AAll other strains used in this study are listed in Table S1N/AN/A	NovtSog 500/550 High							
Outpartie       Visit (Visit)       Infinitial       20024007         MasterPure Yeast RNA       Lucigen       MPY03100         Purification Kit       Agilent       5067-1511         Qubit RNA HS Assay Kit       Thermo Fisher Scientific       Q32852         Direct RNA Sequencing Kit       Nanopore Technologies       QXford         Flow Cell Priming Kit       Nanopore Technologies       EXP-FLP001         MinION Flow Cell (R9.4.1)       Oxford Nanopore       FLO-MIN106D         Experimental models: Organisms/strains       Saccharomyces cerevisiae: BY4741       Jef Boeke's laboratory       N/A         All other strains used in this study are listed in Table S1       N/A       N/A	Output Kit $v2.5$ (150	Illumina	20024907					
JosolJunitication SiteLucigenMPY03100MasterPure Yeast RNA Purification KitLucigenMPY03100Agilent RNA 6000 Nano KitAgilent5067-1511Qubit RNA HS Assay KitThermo Fisher ScientificQ32852Direct RNA Sequencing KitNanopore TechnologiesSQK-RNA002Flow Cell Priming KitOxford Nanopore TechnologiesSQK-RNA002Flow Cell Priming KitOxford Nanopore TechnologiesEXP-FLP001MinION Flow Cell (R9.4.1)Oxford Nanopore TechnologiesFLO-MIN106DExperimental models: Organisms/strains Saccharomyces cerevisiae: BY4741Jef Boeke's laboratoryN/AAll other strains used in this study are listed in Table S1N/AN/A	Cycles)	marmina						
Autor Nation KitLucigenMPY03100Agilent RNA 6000 Nano KitAgilent5067-1511Qubit RNA HS Assay KitThermo Fisher ScientificQ32852Direct RNA Sequencing KitOxford Nanopore TechnologiesSQK-RNA002Flow Cell Priming KitOxford Nanopore TechnologiesSQK-RNA002MinION Flow Cell (R9.4.1)Oxford Nanopore TechnologiesEXP-FLP001Experimental models: Organisms/strainsSaccharomyces cerevisiae: IaboratoryJef Boeke's IaboratoryAll other strains used in this study are listed in Table S1N/AN/A	MasterPure Yeast RNA							
Agilent RNA 6000 Nano Kit       Agilent       5067-1511         Qubit RNA HS Assay Kit       Thermo Fisher Scientific       Q32852         Direct RNA Sequencing Kit       Nanopore Technologies       SQK-RNA002         Flow Cell Priming Kit       Oxford Nanopore Technologies       SQK-RNA002         MinION Flow Cell (R9.4.1)       Oxford Nanopore Technologies       EXP-FLP001         Experimental models: Organisms/strains       Oxford Nanopore       FLO-MIN106D         Saccharomyces cerevisiae: BY4741       Jef Boeke's laboratory       N/A         All other strains used in this study are listed in Table S1       N/A       N/A	Purification Kit	Lucigen	MPY03100					
Qubit RNA HS Assay Kit       Thermo Fisher Scientific       Q32852         Direct RNA Sequencing Kit       Oxford Nanopore Technologies       SQK-RNA002         Flow Cell Priming Kit       Oxford Nanopore Technologies       EXP-FLP001         MinION Flow Cell (R9.4.1)       Oxford Nanopore Technologies       FLO-MIN106D         Experimental models: Organisms/strains       Saccharomyces cerevisiae: BY4741       Jef Boeke's laboratory         All other strains used in this study are listed in Table S1       N/A       N/A	Agilent RNA 6000 Nano Kit	Agilent	5067-1511					
Qubit RINA RS Assay Kit       Scientific       Q32852         Direct RNA Sequencing Kit       Nanopore Technologies       SQK-RNA002         Flow Cell Priming Kit       Nanopore Technologies       SQK-RNA002         MinION Flow Cell (R9.4.1)       Nanopore Technologies       EXP-FLP001         MinION Flow Cell (R9.4.1)       Nanopore Technologies       FLO-MIN106D         Saccharomyces cerevisiae: BY4741       Jef Boeke's laboratory       N/A         All other strains used in this study are listed in Table S1       N/A       N/A		Thermo Fisher	020050					
Direct RNA Sequencing Kit     Oxford Nanopore Technologies     SQK-RNA002       Flow Cell Priming Kit     Oxford Nanopore Technologies     EXP-FLP001       MinION Flow Cell (R9.4.1)     Oxford Nanopore Technologies     FLO-MIN106D       Experimental models: Organisms/strains     Saccharomyces cerevisiae: Iaboratory     Jef Boeke's Iaboratory       All other strains used in this study are listed in Table S1     N/A     N/A	QUDIT KNA HS Assay Kit	Scientific	Q32852					
Direct RNA Sequencing Kit       Nanopore Technologies       SQK-RNA002         Flow Cell Priming Kit       Oxford       EXP-FLP001         Technologies       Oxford       EXP-FLP001         MinION Flow Cell (R9.4.1)       Oxford       FLO-MIN106D         Experimental models: Organisms/strains       Saccharomyces cerevisiae: BY4741       Jef Boeke's laboratory       N/A         All other strains used in this study are listed in Table S1       N/A       N/A		Oxford						
Technologies       Technologies         Flow Cell Priming Kit       Oxford       EXP-FLP001         Technologies       Oxford       EXP-FLP001         MinION Flow Cell (R9.4.1)       Oxford       FLO-MIN106D         Experimental models: Organisms/strains       Saccharomyces cerevisiae:       Jef Boeke's         BY4741       Jaboratory       N/A         All other strains used in this study are listed in Table S1       N/A       N/A	Direct RNA Sequencing Kit	Nanopore	SQK-RNA002					
Flow Cell Priming Kit       Oxford       EXP-FLP001         Technologies       Oxford       Perfection         MinION Flow Cell (R9.4.1)       Oxford       FLO-MIN106D         Experimental models: Organisms/strains       Saccharomyces cerevisiae:       Jef Boeke's         BY4741       Jaboratory       N/A         All other strains used in this study are listed in Table S1       N/A       N/A		Technologies						
Flow Cell Priming Kit     Nanopore Technologies     EXP-FLP001       MinION Flow Cell (R9.4.1)     Oxford Nanopore Technologies     FLO-MIN106D       Experimental models: Organisms/strains     FLO-MIN106D       Saccharomyces cerevisiae: BY4741     Jef Boeke's laboratory     N/A       All other strains used in this study are listed in Table S1     N/A     N/A		Oxford						
Technologies       Technologies         MinION Flow Cell (R9.4.1)       Oxford         Nanopore       FLO-MIN106D         Technologies       Technologies         Experimental models: Organisms/strains       Saccharomyces cerevisiae:         Saccharomyces cerevisiae:       Jef Boeke's         BY4741       Iaboratory         All other strains used in this study are listed in Table S1       N/A         Oligonucleotides       N/A	Flow Cell Priming Kit	Nanopore	EXP-FLP001					
MinION Flow Cell (R9.4.1)       Oxford       FLO-MIN106D         Experimental models: Organisms/strains       FLO-MIN106D         Saccharomyces cerevisiae:       Jef Boeke's       N/A         BY4741       Iaboratory       N/A         All other strains used in this study are listed in Table S1       N/A       N/A		Technologies						
MINION Flow Cell (R9.4.1)       Nanopore Technologies       FLO-MIN106D         Experimental models: Organisms/strains       Saccharomyces cerevisiae:       Jef Boeke's laboratory       N/A         BY4741       Jef Boeke's laboratory       N/A       N/A         All other strains used in this study are listed in Table S1       N/A       N/A		Oxford						
Experimental models: Organisms/strains         Saccharomyces cerevisiae:       Jef Boeke's         BY4741       laboratory         All other strains used in this       N/A         Study are listed in Table S1       N/A	MINION Flow Cell (R9.4.1)	Nanopore	FLO-MIN106D					
Experimental models: Organisms/strains         Saccharomyces cerevisiae:       Jef Boeke's laboratory         BY4741       Iaboratory         All other strains used in this study are listed in Table S1       N/A         N/A       N/A	Experimental models: Orac							
Saccharomyces cerevisiae.     Ser Doeke's       BY4741     Iaboratory       All other strains used in this study are listed in Table S1     N/A       Oligonucleotides	Experimental models: Urganisms/strains							
All other strains used in this study are listed in Table S1 N/A N/A Oligonucleotides	Saccharomyces cerevisiae:	Jel Bueke S	N/A					
study are listed in Table S1 N/A N/A	All other strains used in this	10001 at 01 y						
Oligonucleotides	study are listed in Table S1	N/A	N/A					
~	Oligonucleotides							

PCRTag primers are listed in Supplementary File 1	N/A	N/A				
Recombinant DNA						
See Tables S2-S6 for plasmids used in this study	N/A	N/A				
Software and algorithms						
GeneDesign	40					
R	R Core Team	https://www.R-project.org/				
Rstudio v1.3.1093	Rstudio	https://www.rstudio.com				
Integrated Genomics Viewer (IGV) v2.7.2	Broad Institute <sup>41</sup>	http://software.broadinstitute.org/software/igv/				
Image Studio Lite	LI-COR Biosciences	https://www.licor.com/bio/image-studio/				
GraphPad Prism v9.4.0 for macOS	GraphPad Software	http://www.graphpad.com				
Trimmomatic v0.39	42	http://www.usadellab.org/cms/index.php?page=trimmomatic				
FastQC v0.11.4	43	http://www.bioinformatics.babraham.ac.uk/projects/fastgc/				
Bowtie 2 v2.2.9	44	http://bowtie-bio.sourceforge.net/bowtie2/index.shtml				
ggplot2 v3.3.5	45	https://ggplot2.tidyverse.org/				
Seqtk v1.3	N/A	https://github.com/lh3/seqtk/				
Kallisto v0.46.0	46	https://pachterlab.github.io/kallisto/				
STAR v2.5.2a	47	https://github.com/alexdobin/STAR/				
SAMtools v1.9	48	http://samtools.sourceforge.net/				
BEDtools v2.26.0	49	https://github.com/arq5x/bedtools2/				
sleuth v0.30.0	50	https://pachterlab.github.io/sleuth/				
Other						
Resource website for Sc2.0	N/A	https://syntheticyeast.github.io/				

626 627

#### 628 **RESOURCES AVAILABILITY**

629

#### 630 Lead contact

631

Further information should be directed to and will be fulfilled by the lead contact Jef D. Boeke(jef.boeke@nyulangone.org).

## 634635 Materials availability

636

637 All unique/stable reagents generated in this study are available from the lead contact with a638 completed Materials Transfer Agreement.

# 639640 Data and code availability

- 641642 Data: All data reported in this paper will be shared by the lead contact upon request.
- 643
- 644 Code: This work did not generate any code.
- 645

646 Any additional information required to reanalyze the data reported in this paper is available from 647 the lead contact upon request.

- 648
- 649
- 650

### 651 Experimental Methods

652

### 653 SynIX design

654 Design sequences for *synIX* were developed using BioStudio, according to the same guidelines
 655 used for designing the other synthetic chromosomes for Sc2.0. Designer sequences

- 656 yeast\_chr09\_3\_54 and yeast\_chr09\_3\_55 were used as references for subsequent steps in the
- project. Version yeast\_chr09\_3\_56 was developed after completion of yeast\_chr09\_9\_02, and was
- 658 used as a reference for final strain verification.
- 659

### 660 Yeast media

- 661 Yeast strains were cultured in rich medium (YPD) or synthetic complete (SC) medium with 662 appropriate components "dropped out" e.g., SC–Ura lacks uracil. Yeast transformations were
- 663 performed using lithium acetate protocols.<sup>51</sup> YPG plates used for growth assays contained 3%
- 664 glycerol in place of dextrose as a carbon source. For expanded spot assays, the following media
- 665 types were used: pH 4 and pH 8: pH of 2X YEP + dextrose adjusted using HCl and NaOH,
- respectively before adding agar; camptothecin (Sigma-Aldrich, C9911): 0.1 μg/mL, 0.5 μg/mL, or
  1.0 μg/mL in YPD; sorbitol (Sigma-Aldrich, S1876): 0.5 M, 1.0 M, 1.5 M, or 2.0 M in YPD; 6-
- azauracil (Sigma-Aldrich, A1757): 100 μg/mL in SC medium; hydroxyurea (Sigma-Aldrich, H8627):
- 669 0.2M in YPD; MMS: methyl methanesulfonate (Sigma-Aldrich, 129925), 0.05% in YPD; benomyl
- 670 (Sigma-Aldrich, 381586): 15 μg/mL in YPD; cycloheximide (Sigma-Aldrich, 01810): 10 μg/mL in
- 671 YPD liquid medium for 2 h followed by plating to YPD; H<sub>2</sub>O<sub>2</sub> (Millipore, 88597): 1 mM in YPD liquid
- 672 medium for 2 h followed by plating to YPD. For chromoduction experiments, SC–Lys–Arg or SC–
- 673 Leu–Arg plates were used with cycloheximide (10  $\mu$ g/mL) and canavanine (12  $\mu$ g/mL) added. For
- 674 chromosome destabilization, strains were grown in YPD with 2% galactose and 1% raffinose 675 added. Sporulation medium was prepared using a 50x base consisting of 50 g potassium acetate
- 676 and 0.25 g zinc acetate dihydrate in 100 mL H2O. Final 1x sporulation media was prepared from 2
- mL of 50x sporulation medium base plus 300  $\mu$ L of 10% yeast extract, 200  $\mu$ M uracil, 2 mM leucine
- 678 300  $\mu$ M histidine, and H<sub>2</sub>O to 100 mL.
- 679

### 680 Minichunk assembly

681 Minichunks were originally assembled according to the following process, similar to the one 682 previously described for *synIII*:<sup>5</sup> First, approximately 70-mer oligonucleotides with 20-mer overlaps

- 683 were designed using GeneDesign software<sup>40</sup> and assembled by members of the Build A Genome
- class at Johns Hopkins University using polymerase cycling assembly (PCA) to produce 750 bp
- 685 building blocks. Building blocks were verified by sequencing. Overlapping building blocks with
- approximately 40 bp overlaps were then assembled into 2-4 kb minichunks by direct homologous
   recombination in yeast. Minichunks were recovered from yeast into *E. coli*, and the recovered
- 688 minichunk plasmids were verified by sequencing.
- 689 Later minichunks were ordered directly from DNA synthesis vendors with flanking *Not*l or *Bbs*l
- 690 sites, verified by sequencing, excised from the vectors in which they were built via restriction
- digestion or, for those containing internal *Not*I or *Bbs*I sites, amplified by PCR, and integrated
- 692 directly into the *synIX* strain or used for megachunk assembly.
- 693

#### 694 Megachunk assembly and verification

Unlike the subsequent chromosomes, synIX was initially designed in separate "left and right" 695 halves, which were subsequently merged. During the design process, *synIXL* was originally 696 697 divided into megachunks A-F and V-Z. Following the integration of megachunks A-F as described in the "Minichunk integration by SwAP-In" section, the minichunks comprising megachunks V-Z in 698 699 the original synIX design were reapportioned into three new megachunks, designated G. H. and I. Each megachunk plasmid was designed with flanking Notl sites for eventual excision of the 700 701 assembled megachunk, 13-16 minichunks, a terminal selectable marker (LEU2 or URA3, for use 702 with SwAP-In), and homology arms between the first and last minichunks and the plasmid 703 backbone. Overlaps between segments ranged from approximately 100-800 bp. Minichunks were 704 released from their backbone plasmids or amplified by PCR, and were co-transformed with linkers 705 and megachunk backbone plasmid into wild-type yeast (BY4741). Following plasmid assembly in yeast by homologous recombination and two days of growth on appropriate selective medium 706 707 (SC-Leu or SC-Ura) at 30°C, single colonies were isolated and assessed by PCR to detect the presence or absence of each desired assembly junction. Colonies containing all desired junctions 708 709 were recovered from yeast by phenol extraction and isopropanol precipitation, and subsequently electroporated into EPI300 E. coli cells (Lucigen, EC300150). Megachunk plasmids were extracted 710 from *E. coli* using the PureLink<sup>™</sup> HiPure Plasmid Midiprep Kit (Invitrogen, K2100-15). 711 712 To verify the sequences of the recovered megachunk plasmids, DNA sequencing libraries were 713 prepared using 300 ng of plasmid DNA as input for the NEBNext Ultra II FS DNA Library Kit (NEB. 714 E7805). DNA libraries were sequenced via Illumina NextSeq 500 using 36 bp pair-end reads. 715 Following sequencing, reads were first processed using Trimmomatic v0.39<sup>42</sup> and FastQC v0.11.4<sup>43</sup> aligned to an appropriate reference sequence using bowtie2 v2.2.9.<sup>44</sup> and visualized 716 using IGV v2.7.2.41 717

718

### 719 Minichunk integration by SwAP-In

720 Minichunks were released from their backbones by restriction enzyme digestion, and the 721 minichunks corresponding to one megachunk in the synIX design sequence (approximately 11-16 minichunks) were cotransformed into the partially completed synIX strain. Integration began at the 722 left telomere of chromosome IX and proceeded in a stepwise fashion towards the centromere. At 723 724 each step, the rightmost minichunk integrated contained either a URA3 or a LEU2 auxotrophic marker. These markers were alternated at each step according to the principles of Switching 725 Auxotrophies Progressively for Integration (SwAP-In),<sup>2,6</sup> and colonies containing the new marker 726 and lacking the previous strain's marker were identified by replica plating. Strains were assessed 727 728 by PCRTag analysis after each round of integration.

729

### 730 Megachunk integration by SwAP-In

Sequence-verified megachunks were released from their assembly backbones via *Not*l digestion.
 Approximately 1 µg of digest product was transformed into the semisynthetic *synIX* strain, and

- 733 cells were grown for two days at 30°C on appropriate SC dropout plates to select for the desired
- 734 integration product. As with minichunk integration by SwAP-In, alternating *LEU2* and *URA3*
- 735 auxotrophic markers were used each round, enabling selection of colonies containing the desired
- new integration marker and a lack of the previous round's integration marker, as determined by
- replica plating. Colonies were assessed by PCRTag analysis, and colonies containing the desired

PCRTags of interest were used for subsequent megachunk integration rounds. After integration of all three megachunks introduced by this method, yeast strain sequences were verified by whole genome sequencing, as described in the "Whole genome sequencing" section. To remove the *URA3* marker introduced during megachunk I integration, a PCR product with homology arms flanking *URA3* was transformed into the *synIX* strain, and colonies lacking the *URA3* marker were obtained following growth on 5-FOA medium.

744

### 745 **PCRTag analysis**

746 Cells from a single yeast colony were scraped from a plate using a pipette tip and resuspended in 747 30 µL of 20 mM sodium hydroxide in a 96-well PCR plate. The PCR plate was sealed and placed in a thermal cycler using the following boiling cycle: 3 cycles of 98°C for 3 m and 4°C for 1 m. 1 µL 748 of boiled product and 0.5 µM each primer were used in a 5 µL reaction using the GoTaq Green 749 750 PCR system (Promega, M7123). An acoustic liquid handler (Labcyte, Echo 550) transferred primers and DNA template. Samples were run through the following protocol in a 384-well thermal 751 752 cycler: 95°C for 5 m, 30 cycles of 95°C for 30 s, 55°C for 90 s, and 72°C for 60 s, and a final 753 extension step of 72°C for 7 m. Samples were visualized following electrophoresis on a 1% 754 agarose in 1X Tris-Taurine-EDTA gel containing ethidium bromide using the ChemiDoc XRS+ 755 System (Bio-Rad).

756

### 757 Growth assays

758 A single yeast colony for each strain was inoculated into 5 mL YPD and incubated at 30°C 759 overnight with rotation. For assays monitoring strain growth across multiple passages, overnight cultures were diluted 1:1000 in YPD and incubated at 30°C overnight with rotation until reaching 760 the desired passage number. Overnight yeast culture was diluted to A<sub>600</sub>=0.15 into 6 mL of fresh 761 762 YPD liquid medium and grown at 30°C with rotation until A<sub>600</sub> reached 0.5-0.6. Cultures were 763 serially diluted in tenfold increments in water, with a target OD of approximately 0.1 in 100  $\mu$ L total volume for the first row of the wild-type control strain used to determine the starting culture volume 764 in water for each sample. 5 µL (YPD and YPG plates) or 10 µL (SC plates) of each dilution was 765 spotted to each plate. Plates were grown at appropriate temperatures (room temperature of 766 767 approximately 22°C on YPD, 30°C on YPD, YPG, and SC, and 37°C on YPD and YPG) and photographed daily. For additional media types (Figures S18 and S19), plates were prepared as 768 described in the "Yeast media" section. Strains were serially diluted in tenfold increments of water 769 770 as described above, and 5 µL of each dilution was spotted to each plate. Plates were grown at 771 30°C and photographed daily.

772

### 773 Plate reader growth assays in liquid culture

Three single yeast colonies per strain of interest were resuspended in 1 mL YPD in a 2 mL deep
well plate and incubated for 24 h at 30°C with shaking at 800 rpm. Saturated cultures were diluted
1:100 in 200 μL YPD in a glass bottom 96-well plate with lid (Greiner Bio-One, 655892) and loaded
into a Cytation 5 cell imaging multimode reader (Agilent). The following program was run using the
Cytation 5 Gen5 software: 32 h at 37°C with continuous shaking, with A<sub>600</sub> measured every 10 m.
Data analysis and visualization were conducted using Microsoft Excel and GraphPad Prism 9.

#### 781 Genome editing using CRISPR-Cas9

782 A two-plasmid system with one guide RNA (gRNA) plasmid and one Cas9 plasmid was used for most editing.<sup>52</sup> The Cas9 coding sequence was assembled together with a *TEF1* promoter and 783 784 CYC1 terminator in a pRS415 plasmid backbone (with LEU2 marker). This plasmid was pretransformed into yeast strains prior to the introduction of the gRNA plasmid and donor DNA. Guide 785 786 RNA target sequences were selected by identifying 20 bp upstream of the desired protospacer adjacent motif (PAM) sequence (NGG), gRNAs were cloned into a 2-micron plasmid (pRS426, with 787 788 URA3 marker) under the control of a SNR52 promoter using Gibson assembly. A modified version 789 of the gRNA plasmid was used for editing two loci simultaneously. In this case, a second gRNA 790 was cloned into the plasmid under control of a RPR1 promoter and RPR1 terminator. gRNA 791 sequences were confirmed by Sanger sequencing (Genewiz). To edit the synIX strain, 792 approximately 50 ng of gRNA plasmid and 200-300 fmol (for minichunk patch replacement) or 1 793 pmol (for point mutation fixing) of donor DNA were transformed into yeast cells containing the Cas9 794 plasmid. Cells were grown on dropout medium to select for the presence of both the Cas9 and 795 gRNA plasmids, and editing was confirmed using either PCRTag analysis or PCR followed by 796 Sanger sequencing (Genewiz).

797

### 798 Sporulation and tetrad dissection experiments for bug mapping

The synIX strain yLHM1192 was crossed to a wild-type yeast strain (BY4742 with a URA3-pGAL-

- 800 *CEN9* cassette integrated near the centromere, also known as yLHM0539) to generate a
- 801 heterozygous diploid (yLHM1233) with one copy of wild-type chromosome *IX* and one copy of
- *synIX*. To prepare strains for sporulation, a single colony of each strain was inoculated into 5 mL
   YPD and incubated at 30°C overnight with rotation. Overnight cultures were diluted to an OD of ~1
- in YPD and grown to an OD of ~4, washed five times with water, and resuspended in 2 mL of 1X
- sporulation medium. Strains were incubated at room temperature for 7-10 days with rotation, and
- monitored for the presence of tetrads. For tetrad dissection, 100  $\mu$ L of these resuspended yeast
- cells in sporulation medium were washed and incubated with 25  $\mu$ L of 0.5 mg/mL zymolyase in 1M sorbitol for 8 m. 200  $\mu$ L of 1M sorbitol was added to the cells, and 10  $\mu$ L of the resulting mixture
- 809 was added to a YPD plate. Tetrads were separated and picked using a dissection microscope
- 810 (Singer Instruments) Spores were grown for 2-3 days on YPD until visible colonies emerged.
- 811 These strains were used for growth assays as described in the "Growth assays" section to classify
- 812 each spore as healthy (similar fitness to the wild-type parent) or sick (similar fitness to the *synIX*
- parent). Strains from eighteen four-spore tetrads (72 strains) were prepared for whole genome
- sequencing, as described in the "Whole genome sequencing" section below.
- 815

### 816 Plasmid cloning and transformation for *synlX* bug mapping

817 Wild-type genes of interest plus approximately 500 bp of upstream flanking sequence and 200 bp 818 of downstream flanking sequence were amplified from BY4741 by PCR. Genes were cloned into

- plasmid pRS413 linearized with *Afel* by Gibson assembly and transformed into TOP10 cells.
- Colonies were grown up in 2 mL LB + 75  $\mu$ g/mL carbenicillin and miniprepped using the Zyppy
- 821 Plasmid Miniprep Kit (Zvmo Research, D4037). Plasmid insert sequences were verified by Sanger
- sequencing. Plasmids were subsequently transformed into the *synIX* strain and grown for two days
- 823 on SC–His plates. Yeast colonies containing each plasmid of interest were isolated, grown in SC–

- His liquid medium, and used for growth assays (as described in the "Growth assays" section) to assess their relative fitness compared to the starting *synIX* strain.
- 825 826

### 827 Selective destabilization and loss of one copy of *synIX* from disomic strain

828 A cassette containing a galactose-inducible centromere and a URA3 marker (URA3-pGAL-CEN9)

- 829 with homology arms matching *synIX* was amplified by PCR and transformed into the *synIX* strain.
- After two days of growth at 30°C on SC–Ura plates, single colonies were purified and assessed by
- 331 junction PCR for the insertion of the *URA3-pGAL-CEN9* fragment. Strains were then grown in YEP + 2% galactose + 1% raffinose for two days and diluted to 5-FOA plates to obtain single colonies
- + 2% galactose + 1% raffinose for two days and diluted to 5-FOA plates to obtain single colonies
   lacking URA3. Strains were single-colony purified and prepared for whole genome sequencing to
- 834 assess chromosome copy number.
- 835

#### 836 Whole genome sequencing

A single yeast colony was inoculated into 5 mL of YPD liquid medium at 30°C overnight with

- rotation. Overnight cultured yeast cells were harvested by centrifugation at 3000 x g for 3 m, and
- 839 yeast genomic DNA was extracted using the Fungi/Yeast Genomic DNA Isolation Kit (Norgen,
- 27300). Genomic DNA concentration was determined using the Qubit dsDNA HS Assay Kit
- (Thermo Fisher Scientific, Q32854). DNA sequencing libraries were prepared using 300 ng of
- yeast genomic DNA as input for the NEBNext Ultra II FS DNA Library Kit (NEB, E7805). DNA
- libraries were sequenced via Illumina NextSeq 500 using either 75 bp or 36 bp pair-end reads.
- Following sequencing, reads were first processed using Trimmomatic v0.39<sup>42</sup> and FastQC
- v0.11.4,<sup>43</sup> aligned to an appropriate reference genome using bowtie2 v2.2.9,<sup>44</sup> and visualized using
- IGV v2.7.2.<sup>41</sup> Chromosome coverage plots were prepared in Rstudio version 1.3.1093 using the
   ggplot2 v3.3.5 package.<sup>45</sup>
- 847 848

### 849 Southern blot analysis

A single yeast colony was inoculated into 10 mL of YPD medium at 30°C overnight with rotation.
 Overnight cultured yeast cells were harvested by centrifugation at 3000 x g for 3 m, and yeast
 genomic DNA was extracted using the MasterPure<sup>™</sup> Yeast DNA Purification Kit (Lucigen,

- 853 MPY80200) according to the manufacturer's instructions, with an RNase A treatment step. DNA
- 854 was quantified using the Qubit dsDNA BR Assay Kit (Thermo Fisher Scientific, Q32850). 2 μg of
- DNA were digested with *Xhol* at 37°C overnight. Digests were electrophoresed on a 1% Tris Borate-EDTA gel in 0.5X Tris-Borate-EDTA running buffer for 10 h at 50V. Following
- electrophoresis, the gel was washed in depurination solution (0.25 M HCl) for 15 m, denaturation
- solution (0.5 M NaOH, 1.5 M NaCl) for 60 m, and neutralization solution (0.5 M Tris, 1.5 M NaCl,
- pH 7.5) for 30 m. Chromosomes were transferred by capillarity to a nylon membrane (Pall
   Corporation, 60207) using 10X SSC buffer (20X SSC: 1.5 M NaCl, 0.3 M sodium citrate, pH 7) for
- 861 18 h. After transfer, the membrane was washed in 2X SSC buffer for 5 m and baked at 80°C for 30
- m. Prehybridization was performed for 60 m in ULTRAhyb<sup>™</sup> Ultrasensitive Hybridization Buffer
   (Invitrogen, AM8669) at 55°C with rotation. For hybridization, probes were prepared by labeling a
- 864 ~360 bp PCR fragment of universal telomere cap sequence from plasmid pJS160 using the PCR
- 865 DIG Probe Synthesis Kit (Roche, 11636090910) with oligos 5'-
- 866 GCTATACGAAGTTATTAGGGTAGTGTG-3' and 5'- CTGCAGGTCGACTCTAGAGGATC-3',
- according to the manufacturer's instructions. The following thermocycler conditions were used: 1
- 868 cycle of 95°C for 2 m, 30 cycles of 95°C for 30 s, 60°C for 30 s, and 72°C for 40 s, and a final
- 869 elongation step of 72°C for 7 m. Probes were purified using Zymo-Spin I Columns (Zymo

870 Research, C1003) and denatured at 95°C for 5 m. Hybridization was performed at 55°C overnight with rotation in 10 mL ULTRAhyb<sup>™</sup> Ultrasensitive Hybridization Buffer (Invitrogen, AM8669) using 871 872 500 ng of probe per experiment. The blot was washed twice with 2X SSPE buffer (20X SSPE buffer: 0.02 M EDTA and 2.98 M NaCl in 0.2 M phosphate buffer, pH 7.4; Sigma-Aldrich, S2015) 873 874 for 5 m at room temperature, twice with 2X SSPE + 1% SDS for 30 m at 55°C, and twice with 0.1X 875 SSPE for 15 m at 55°C. Blot was incubated overnight in 1:1 Intercept (TBS) Blocking Buffer (LI-COR Biosciences, 927-60001): 1X Tris-Buffered Saline (TBS) with 0.1% TWEEN 20 (Sigma-876 877 Aldrich, P1379) (1X TBST) and 1% SDS at room temperature. The next day, the membrane was 878 washed with 1X TBST for 5 m and incubated for 1 hour at room temperature with rabbit anti-DIG 879 antibody (working concentration 1:2500 in 1:1 Intercept (TBS) Blocking Buffer: 1X TBST; 880 Invitrogen, 700772). The primary antibody solution was washed out three times with 1X TBST for 15 m at room temperature. The membrane was incubated in goat anti-rabbit IgG secondary 881 solution (LI-COR Biosciences, 926-68071, used at 1:10,000 in 1:1 Intercept (TBS) Blocking Buffer 882 883 (LI-COR Biosciences, 927-60001): 1X TBST with 1% SDS) for 1.5 h. The blot was washed three 884 times with 1X TBST for 15 m at room temperature plus one time with 1X TBS for 30 m at room 885 temperature. A LI-COR Odyssey Instrument was used to develop the blot images.<sup>53</sup>

886

#### 887 RNA extraction, sequencing, and analysis for short read sequencing

A single yeast colony was inoculated into 5 mL of YPD liquid medium at 30°C overnight with 888 rotation. Overnight yeast culture was diluted to A<sub>600</sub>=0.15-0.2 into 6 mL of fresh YPD liquid medium 889 and grown at 30°C or 37°C with rotation until A<sub>600</sub> reached 0.8-1.0. Yeast cells were harvested from 890 4 mL of culture, and total RNA was extracted using the RNeasy Mini Kit (QIAGEN, 74106) 891 according to the manufacturer's instructions. RNA was guantified using the Qubit BR RNA Assay 892 Kit (Thermo Fisher Scientific, Q10210). 1 µg total RNA was used as input for the QIAseg Stranded 893 Total RNA Library Kit (QIAGEN, 180745) plus QIAseg FastSelect -rRNA Yeast Kit (QIAGEN, 894 895 334217). RNA libraries were sequenced via Illumina NextSeq 500 using 75 bp paired-end reads. 896 For applications requiring similar numbers of reads per sample, downsampling was performed using segtk v1.3. Samples were first processed using Trimmomatic v0.3942 and FastQC v0.11.4.43 897 Processed reads were aligned to the S288C transcriptome and a custom synIX transcriptome (in 898 which chromosome IX transcripts were replaced with their synthetic versions) using Kallisto 899 v0.46.0<sup>46</sup> and STAR v2.5.2a.<sup>47</sup> STAR-aligned results were further processed using SAMtools 900 v1.9<sup>48</sup> and BEDtools v2.26.0,<sup>49</sup> and visualized in IGV v2.7.2.<sup>41</sup> Kallisto pseudoalignment data was 901 further processed in Rstudio version 1.3.1093 using the sleuth v0.30.0<sup>50</sup> and gaplot2 v3.3.5<sup>45</sup> 902 packages. For volcano plots, we calculated effect scores (approximately log2 fold change values), 903 904 tested results for significance using Wald's test, and corrected for multiple testing with the false discovery rate adjusted p-value using the Benjamini-Hochberg method. The following genes 905 deleted from the synIX design were removed from the plots in Figure 4D and Figure 4E: YIL173W 906 (VTH1), YIL171W (part of HTX12 pseudogene), YIL170W (part of HTX12 pseudogene), YIL169C 907 (CSS1), YIL167W (SDL1, blocked reading frame), YIL082W, YIL060W, YIL059C, YIL058W, 908 909 YILWTy3-1, YILCdelta3. The following dubious ORFs, retrotransposable elements, and 910 subtelomeric genes were also removed: YBL100C, YGR296C-B, YHL050W-A, YHR219C-A, YIL020C-A, YJR029W, YLL066C, YML133W-B, YMR046C, YNL339W-B, YNL339C, YOR277C, 911 YPL283W-B, YPR158W-B, YPR204C-A, YBLWdelta4, YCLWdelta3. Unfiltered versions of the 912 913 volcano plots are depicted in Figure S20. 914

915

### 916 **RNA extraction, sequencing, and analysis for nanopore direct RNA sequencing**

917 Methods recently described by us were used to perform nanopore sequencing to analyze RNA isoforms.<sup>54</sup> Total RNA was extracted from 50 mL flash-frozen cell pellets grown to mid-log (OD<sub>600</sub> 918 919 ~0.65-0.85) using the MasterPure™ Yeast RNA Purification Kit (Lucigen) including a Dnase I 920 treatment step. RNA (diluted 1:10) quality and concentration were measured by Agilent 2100 921 Bioanalyzer with the Agilent RNA 6000 Nano Kit and Qubit<sup>™</sup> RNA High Sensitivity Kit (Thermo 922 Fisher Scientific), respectively. Poly(A) mRNA was enriched from 50 µg total RNA on 132 µL Dvnabeads oligo(dT)<sub>25</sub> beads. The Direct RNA Sequencing Kit (SQK-RNA002, Oxford Nanopore 923 924 Technologies) was used to generate libraries from 500 ng poly(A) RNA. An optional reverse transcription was performed at 50°C for 50 m using SuperScript™ IV Reverse Transcriptase 925 (Invitrogen) in between the ligation of the RTA and RMX adaptors. Following reverse transcription 926 927 the RNA:cDNA was cleaned up with 1.8 volumes of Agencourt RNAClean XP beads and washed 928 with 70% ethanol. Following RMX ligation only one volume of beads was used in the clean-up, 929 and WSB (SQK-RNA002) was used in the wash steps. Direct RNA libraries (typically 150-200 ng) 930 were loaded onto primed (EXP-FLP001) MinION flow cells (FLO-MIN106D, R9 version) in RRB 931 buffer and run on the GridION with MinKNOW 5.2.2 for up to 72 h. Nanopore long reads were 932 base-called, trimmed of adapter sequences, and filtered for guality, retaining only those with the 933 best alignment scores for multi-mapping reads, as previously described.<sup>54</sup> For TSS distributions, 934 long-reads that overlapped at least 75% of the annotated gene were used.

935

### 936 Construction of FLAG-tagged EST3 strains

A DNA fragment containing a six-glycine linker plus a 3X FLAG tag (DYKDHDG-DYKDHDIDYKDDDDK) was designed for integration at the 3' end of *EST3*, in-frame with the rest of the
protein and inserted just before the *EST3* stop codon. The segment containing the tag was flanked
by 70-100 bp homology arms matching the yeast genome sequence adjacent to the desired
insertion site. Genome editing using CRISPR-Cas9 in yeast (as describe above) was used to cut
the genome near the end of the *EST3* coding sequence and insert the tag of interest. Colony PCR
and Sanger sequencing (Genewiz) were used to confirm in-frame insertion of the tag.

944

### 945 Immunoblot analysis

946 A single yeast colony was inoculated into 5 mL of YPD liquid medium at 37°C overnight with 947 rotation. Overnight culture was diluted to A<sub>600</sub>=0.15-0.2 into 50 mL of fresh YPD liquid medium and 948 grown at 37°C with rotation until A<sub>600</sub> reached 0.8-1.0. Yeast were harvested by centrifugation, and pellets were frozen overnight at -80C. Pellets were subsequently resuspended in 500 µL of lysis 949 buffer (20 mM HEPES pH 7.4, 1% Triton<sup>™</sup> X-100 (Sigma-Aldrich, X100), 2 mM MqCl<sub>2</sub>, 500 μM 950 NaCl, and 1X protease inhibitor (Roche)). Cell suspension was added to MP Biomedicals FastPrep 951 Lysing Matrix C and subjected to the following shaking program at 4°C using a FastPrep-24 5G 952 953 (MP Biomedicals): 7 cycles of 15 s shaking at 6.0 m/second + 30 s rest between cycles. Samples were centrifuged at 13000 x g for 10 m at 4°C, and 30 µL of supernatant was transferred to a new 954 955 tube and mixed with 10 µL of 4X LDS loading buffer (Invitrogen, NP0007) (pre-immunoprecipitation sample). The remaining supernatant was transferred to pre-equilibrated FLAG M2 (Sigma F1804)-956 957 conjugated M-270 Dynabeads (Invitrogen, 14302D) for immunoprecipitation.<sup>55</sup> Samples with beads

959 at 700 rpm in 50 µL 1.1X LDS loading buffer (Invitrogen, NP0007) in lysis buffer for 10 m. 2-960 mercaptoethanol (1.43 M final) was added, and samples were heated to 95°C for 5 m and loaded 961 onto an SDS-PAGE gel for immunoblotting. Mouse anti-FLAG M2 antibody (Sigma F1804, used at 1:5000) and rabbit <sup>54</sup> anti-histone H3 antibody (Abcam, 1791, used at 1:5000) were used for 962 963 blotting proteins of interest overnight at 4°C. Goat anti-mouse IgG1 (LI-COR Biosciences, 926-32350, used at 1:20,000) and goat anti-rabbit IgG (LI-COR Biosciences, 926-68071, used at 964 1:20,000) were used for blotting FLAG and H3 primary antibodies, respectively, during a one-hour 965 incubation at room temperature. A LI-COR Odyssey Instrument was used to develop the blot 966 images.53 967

968

### 969 Pulsed-field gel electrophoresis

A single yeast colony was inoculated into 7 mL of YPD liquid medium and grown at 30°C for two 970 days with rotation. Approximately 40 mg of yeast cells were harvested by centrifugation. To 971 prepare plugs, 16 µL of zymolyase solution (25 mg/mL zymolyase 20T in 10 mM KH<sub>2</sub>PO<sub>2</sub> pH 7.5) 972 973 and 360 µL of 0.5% melted low melting point agarose (Bio-Rad, 1620137) in 100 mM EDTA pH 7.5 974 were added to the cell pellet and mixed by pipetting. 90 µL of this suspension was added per well 975 of a BioRad plug mold, and the plugs were cooled for 30 m at 4°C. Solidified plugs were then 976 transferred into 15 mL tubes containing 1 mL of digestion buffer (0.5 M EDTA, 10 mM Tris, pH 7.5), 977 mixed by inversion, and incubated at 37°C with rotation overnight. The next day, 400 µL of 978 proteinase K solution (5% sarcosyl and 5 mg/mL proteinase K in 0.5 M EDTA, pH 7.5) were added 979 to the tube. Samples were incubated at 50°C for five h, and then washed once with H<sub>2</sub>O and three times with TE (2 mM Tris, 1 mM EDTA, pH 8.0). Plugs were stored at 4°C until readv for gel 980 loading. Half of one plug per sample was embedded per well of a gel consisting of 1% low melting 981 point agarose in 0.5X Tris-Borate-EDTA buffer. Samples were separated in running buffer of 0.5X 982 Tris-Borate-EDTA by clamped homogeneous electric field gel electrophoresis using the CHEF 983 Mapper XA Pulsed-Field Electrophoresis System (Bio-Rad), as previously described.<sup>56</sup> The 984 985 following program was used: auto-algorithm with size range of 200 kb to 2.5 Mb, temperature of 14°C, voltage of 6 V/cm, included angle of 120°. After electrophoresis, gels were stained using 5 986 µg/mL ethidium bromide in water for 30 m, destained in water for 30 m, and then imaged using the 987 988 ChemiDoc XRS+ System (Bio-Rad).

989

#### 990 Chromosome substitution and conditional centromere destabilization

For initial experiments, disomic synIX strain yLHM0387 was used as the donor and yWZ601 as the 991 992 recipient for kar1-1-mediated chromoduction. For follow-up experiments, disomic strain vLHM0604 and monosomic strains yLHM0721, yLHM1591, and yLHM1601 were used as donors and yWZ602 993 994 as the recipient for *kar1-1*-mediated chromoduction. yWZ601 and yWZ602 were derived from DMy044<sup>24</sup> following deletion of LYS12 from wild-type chromosome IX and insertion of a URA3-995 *pGAL* promoter cassette upstream of *CEN*9.<sup>26</sup> Strains were grown independently on YPD plates, 996 997 and subsequently mated on YPD plates via replica plating. After 14 h of growth, these plates were replica plated to SC-Lys-Arg plates with 12 µg/mL canavanine and 10 µg/mL cycloheximide 998 999 added. Plates were grown until individual colonies emerged; if a patched formed instead, patches were replica plated to the same selective medium and grown for additional time until individual 1000 colonies emerged. Single colonies were restreaked to fresh selective plates and analyzed by 1001

- 1002 PCRTagging, pulsed-field gel electrophoresis, and whole genome sequencing. To lose wild-type 1003 chromosome *IX*, strains were grown in YPD overnight, diluted 1:1000 in YEP + 2% galactose
- 1004 medium for 24 h, and plated on 5-FOA medium to select for strains that had lost wild-type
- 1005 chromosome IX. Single colonies were genotyped by PCRTag analysis and further characterized by
- 1006 pulsed-field gel electrophoresis and whole genome sequencing. For chromosome substitution to
- 1007 move *synIX* from yLHM1601 to recipient strain yLHM2231 (derived from BY4742), a LEU2 marker 1008 was integrated into *synIX* near the end of megachunk H, in minichunk 101, generating strain
- 1009 yLHM1754. The same protocol as listed above was followed using donor strain yLHM1754 and
- 1010 recipient strain yLHM2231, except after mating and 14 h of growth, plates were replica plated to
- 1011 SC-Leu-Arg (rather than SC-Lys-Arg) plates with 12  $\mu$ g/mL canavanine and 10  $\mu$ g/mL
- 1012 cycloheximide added.

### 1013 Acknowledgments

1014

1015 We thank Meghan O'Keefe for collecting the BAG student authors' information. We thank the NSF 1016 for grants MCB-1026068, MCB-1443299, MCB-1616111 and MCB-1921641 to JDB and MCB-1017 1445545 to JSB for supporting this work.

1018

### 1019 **Competing Interests**

Jef Boeke is a Founder and Director of CDI Labs, Inc., a Founder of and consultant to
Neochromosome, Inc, a Founder, SAB member of and consultant to ReOpen Diagnostics, LLC
and serves or served on the Scientific Advisory Board of the following: Logomix, Inc., Sangamo,
Inc., Modern Meadow, Inc., Rome Therapeutics, Inc., Sample6, Inc., Tessera Therapeutics, Inc.
and the Wyss Institute. Joel Bader is a founder of Neochromosome, Inc., consultant to Opentrons
Labworks, Inc. and advisor to Reflexion Pharmaceuticals, Inc.

1026

### 1027 Author contributions

J.D.B., S.C., and W.Z. conceptualized the study. W.Z., S.C., L.H.M., and J.D.B. designed the 1028 experiments, J.D.B. and J.S.B. designed svnIX. The Build A Genome class members constructed 1029 1030 building blocks. J.S., E.C., H.M., and K.Z led Build A Genome class research. V.S., N.A., and S.R. 1031 constructed minichunks and incorporated them for megachunks A to F. L.H.M. designed, 1032 constructed, and delivered megachunks G through I and performed error correction, bug mapping, 1033 bug correction, and functional studies of synIX. L.H.M. and W.Z. developed and evaluated 1034 chromosome substitution method. L.H.M., A.L.H., J.C., V.F., E.L., and G.S. performed sequencing 1035 and data analyses. L.H.M., W.Z., and J.D.B. wrote the manuscript. L.A.M., H.J., J.L., M.S.T., and W.R.B. provided advice. S.C., J.D.B., J.S.B., and L.M.S. provided lab supervision. All authors 1036 1037 reviewed and edited the manuscript.

1038

#### 1039 **References**

- 1040
- 1041 1. Dymond, J., and Boeke, J. (2012). The Saccharomyces cerevisiae SCRaMbLE system and genome minimization. Bioeng Bugs *3*, 168-171. 10.4161/bbug.19543.
- 10432.Richardson, S.M., Mitchell, L.A., Stracquadanio, G., Yang, K., Dymond, J.S., DiCarlo, J.E.,1044Lee, D., Huang, C.L., Chandrasegaran, S., Cai, Y., et al. (2017). Design of a synthetic yeast1045genome. Science 355, 1040-1044. 10.1126/science.aaf4557.
- 10463.Shen, M.J., Wu, Y., Yang, K., Li, Y., Xu, H., Zhang, H., Li, B.Z., Li, X., Xiao, W.H., Zhou, X.,1047et al. (2018). Heterozygous diploid and interspecies SCRaMbLEing. Nat Commun 9, 1934.104810.1038/s41467-018-04157-0.
- Shen, Y., Stracquadanio, G., Wang, Y., Yang, K., Mitchell, L.A., Xue, Y., Cai, Y., Chen, T.,
   Dymond, J.S., Kang, K., et al. (2016). SCRaMbLE generates designed combinatorial
   stochastic diversity in synthetic chromosomes. Genome Res 26, 36-49.
   10.1101/gr.193433.115.
- Annaluru, N., Muller, H., Mitchell, L.A., Ramalingam, S., Stracquadanio, G., Richardson,
   S.M., Dymond, J.S., Kuang, Z., Scheifele, L.Z., Cooper, E.M., et al. (2014). Total synthesis
   of a functional designer eukaryotic chromosome. Science *344*, 55-58.
   10.1126/science.1249252.
- Dymond, J.S., Richardson, S.M., Coombes, C.E., Babatz, T., Muller, H., Annaluru, N.,
   Blake, W.J., Schwerzmann, J.W., Dai, J., Lindstrom, D.L., et al. (2011). Synthetic
   chromosome arms function in yeast and generate phenotypic diversity by design. Nature
   477, 471-476. 10.1038/nature10403.
- Mitchell, L.A., Wang, A., Stracquadanio, G., Kuang, Z., Wang, X., Yang, K., Richardson, S.,
   Martin, J.A., Zhao, Y., Walker, R., et al. (2017). Synthesis, debugging, and effects of
   synthetic chromosome consolidation: synVI and beyond. Science 355.
   10.1126/science.aaf4831.
- Shen, Y., Wang, Y., Chen, T., Gao, F., Gong, J., Abramczyk, D., Walker, R., Zhao, H.,
   Chen, S., Liu, W., et al. (2017). Deep functional analysis of synII, a 770-kilobase synthetic
   yeast chromosome. Science 355. 10.1126/science.aaf4791.
- Wu, Y., Li, B.Z., Zhao, M., Mitchell, L.A., Xie, Z.X., Lin, Q.H., Wang, X., Xiao, W.H., Wang, Y., Zhou, X., et al. (2017). Bug mapping and fitness testing of chemically synthesized chromosome X. Science 355. 10.1126/science.aaf4706.
- 1071 10. Xie, Z.X., Li, B.Z., Mitchell, L.A., Wu, Y., Qi, X., Jin, Z., Jia, B., Wang, X., Zeng, B.X., Liu,
  1072 H.M., et al. (2017). "Perfect" designer chromosome V and behavior of a ring derivative.
  1073 Science 355. 10.1126/science.aaf4704.
- 1074 11. Zhang, W., Zhao, G., Luo, Z., Lin, Y., Wang, L., Guo, Y., Wang, A., Jiang, S., Jiang, Q.,
  1075 Gong, J., et al. (2017). Engineering the ribosomal DNA in a megabase synthetic
  1076 chromosome. Science *355*. 10.1126/science.aaf3981.
- 1077 12. Zhang, W., Lazar-Stefanita, L., Yamashita, H., Shen, M.J., Mitchell, L.A., Kurasawa, H., 1078 Haase, M.A.B., Sun, X., Jiang, Q., Lauer, S.L., et al. (2022). Manipulating the 3D

- 1079
   organization of the largest synthetic yeast chromosome. 2022.2004.2009.487066.

   1080
   10.1101/2022.04.09.487066 %J bioRxiv.
- 13. Blount, B.A., Lu, X., Driessen, M.R.M., Jovicevic, D., Sanchez, M.I., Ciurkot, K., Zhao, Y.,
  Lauer, S., McKiernan, R.M., Gowers, G.-O.F., et al. (2022). Synthetic yeast chromosome XI
  design enables extrachromosomal circular DNA formation on demand.
  2022.2007.2015.500197. 10.1101/2022.07.15.500197 %J bioRxiv.
- Williams, T.C., Kroukamp, H., Xu, X., Wightman, E.L.I., Llorente, B., Borneman, A.R., Carpenter, A.C., Van Wyk, N., Espinosa, M.I., Daniel, E.L., et al. (2022). Laboratory evolution and polyploid SCRaMbLE reveal genomic plasticity to synthetic chromosome defects and rearrangements. 2022.2007.2022.501046. 10.1101/2022.07.22.501046 %J bioRxiv.
- Schindler, D., Walker, R.S.K., Jiang, S., Brooks, A.N., Wang, Y., Müller, C.A., Cockram, C., Luo, Y., García, A., Schraivogel, D., et al. (2022). Design, Construction, and Functional Characterization of a tRNA Neochromosome in Yeast. 2022.2010.2003.510608.
   10.1101/2022.10.03.510608 %J bioRxiv.
- 16. Shen, Y., Gao, F., Wang, Y., Wang, Y., Zheng, J., Gong, J., Zhang, J., Luo, Z., Schindler,
  D., Deng, Y., et al. (2022). Dissecting aneuploidy phenotypes by constructing Sc2.0
  chromosome VII and SCRaMbLEing synthetic disomic yeast. 2022.2009.2001.506252.
  10.1101/2022.09.01.506252 %J bioRxiv.
- 1098
  17. Luo, J., Mercy, G., Vale-Silva, L.A., Sun, X., Agmon, N., Zhang, W., Yang, K.,
  1099
  100 Stracquadanio, G., Thierry, A., Ahn, J.Y., et al. (2018). Synthetic chromosome fusion:
  effects on genome structure and function. 381137. 10.1101/381137 %J bioRxiv.
- 1101
  18. Lendvay, T.S., Morris, D.K., Sah, J., Balasubramanian, B., and Lundblad, V. (1996).
  Senescence mutants of Saccharomyces cerevisiae with a defect in telomere replication identify three additional EST genes. Genetics *144*, 1399-1412.
  10.1093/genetics/144.4.1399.
- 110519.Morris, D.K., and Lundblad, V. (1997). Programmed translational frameshifting in a gene1106required for yeast telomere replication. Curr Biol 7, 969-976. 10.1016/s0960-11079822(06)00416-7.
- 110820.Hughes, T.R., Evans, S.K., Weilbaecher, R.G., and Lundblad, V. (2000). The Est3 protein is1109a subunit of yeast telomerase. Curr Biol *10*, 809-812. 10.1016/s0960-9822(00)00562-5.
- 21. Zhao, Y., Coelho, C., Hughes, A.L., Lazar-Stefanita, L., Yang, S., Brooks, A.N., Walker,
  R.S.K., Zhang, W., Lauer, S., Hernandez, C., et al. (2022). Debugging and consolidating
  multiple synthetic chromosomes reveals combinatorial genetic interactions.
  2022.2004.2011.486913. 10.1101/2022.04.11.486913 %J bioRxiv.
- Conde, J., and Fink, G.R. (1976). A mutant of Saccharomyces cerevisiae defective for
  nuclear fusion. Proceedings of the National Academy of Sciences of the United States of
  America 73, 3651-3655. 10.1073/pnas.73.10.3651.
- 1117 23. Dutcher, S.K. (1981). Internuclear transfer of genetic information in kar1-1/KAR1
  1118 heterokaryons in Saccharomyces cerevisiae. Mol Cell Biol *1*, 245-253.
  1119 10.1128/mcb.1.3.245-253.1981.

- Ji, H., Moore, D.P., Blomberg, M.A., Braiterman, L.T., Voytas, D.F., Natsoulis, G., and
  Boeke, J.D. (1993). Hotspots for unselected Ty1 transposition events on yeast chromosome
  III are near tRNA genes and LTR sequences. Cell 73, 1007-1018. 10.1016/00928674(93)90278-x.
- Substantiation 25.
  Guo, Z., Yin, H., Ma, L., Li, J., Ma, J., Wu, Y., and Yuan, Y. (2022). Direct Transfer and Consolidation of Synthetic Yeast Chromosomes by Abortive Mating and Chromosome Elimination. ACS Synth Biol. 10.1021/acssynbio.2c00174.
- Reid, R.J., Sunjevaric, I., Voth, W.P., Ciccone, S., Du, W., Olsen, A.E., Stillman, D.J., and Rothstein, R. (2008). Chromosome-scale genetic mapping using a set of 16 conditionally stable Saccharomyces cerevisiae chromosomes. Genetics *180*, 1799-1808.
  10.1534/genetics.108.087999.
- Johnson, D.R., Knoll, L.J., Levin, D.E., and Gordon, J.I. (1994). Saccharomyces cerevisiae
  contains four fatty acid activation (FAA) genes: an assessment of their role in regulating
  protein N-myristoylation and cellular lipid metabolism. J Cell Biol *127*, 751-762.
  10.1083/jcb.127.3.751.
- 1135 28. Knoll, L.J., Johnson, D.R., and Gordon, J.I. (1994). Biochemical studies of three
  1136 Saccharomyces cerevisiae acyl-CoA synthetases, Faa1p, Faa2p, and Faa3p. J Biol Chem
  1137 269, 16348-16356.
- Yoshikawa, K., Tanaka, T., Ida, Y., Furusawa, C., Hirasawa, T., and Shimizu, H. (2011).
  Comprehensive phenotypic analysis of single-gene deletion and overexpression strains of
  Saccharomyces cerevisiae. Yeast 28, 349-361. 10.1002/yea.1843.
- 114130.Devine, S.E., and Boeke, J.D. (1996). Integration of the yeast retrotransposon Ty1 is1142targeted to regions upstream of genes transcribed by RNA polymerase III. Genes Dev 10,1143620-633. 10.1101/gad.10.5.620.
- 114431.Geiduschek, E.P., and Tocchini-Valentini, G.P. (1988). Transcription by RNA polymerase III.1145Annu Rev Biochem 57, 873-914. 10.1146/annurev.bi.57.070188.004301.
- 114632.Newman, A.J., Ogden, R.C., and Abelson, J. (1983). tRNA gene transcription in yeast:1147effects of specified base substitutions in the intragenic promoter. Cell 35, 117-125.114810.1016/0092-8674(83)90214-3.
- 1149 33. Hull, M.W., Erickson, J., Johnston, M., and Engelke, D.R. (1994). tRNA genes as
   1150 transcriptional repressor elements. Mol Cell Biol *14*, 1266-1277. 10.1128/mcb.14.2.1266 1151 1277.1994.
- 115234.Fourel, G., Revardel, E., Koering, C.E., and Gilson, E. (1999). Cohabitation of insulators and1153silencing elements in yeast subtelomeric regions. EMBO J 18, 2522-2537.115410.1093/emboj/18.9.2522.
- 115535.Bolton, E.C., and Boeke, J.D. (2003). Transcriptional interactions between yeast tRNA1156genes, flanking genes and Ty elements: a genomic point of view. Genome Res 13, 254-263.115710.1101/gr.612203.
- 115836.Donze, D. (2012). Extra-transcriptional functions of RNA Polymerase III complexes: TFIIIC1159as a potential global chromatin bookmark. Gene 493, 169-175. 10.1016/j.gene.2011.09.018.

- 37. Peter, J., De Chiara, M., Friedrich, A., Yue, J.X., Pflieger, D., Bergstrom, A., Sigwalt, A.,
  Barre, B., Freel, K., Llored, A., et al. (2018). Genome evolution across 1,011
  Saccharomyces cerevisiae isolates. Nature *556*, 339-344. 10.1038/s41586-018-0030-5.
- 1163 38. Liti, G., Carter, D.M., Moses, A.M., Warringer, J., Parts, L., James, S.A., Davey, R.P.,
  1164 Roberts, I.N., Burt, A., Koufopanou, V., et al. (2009). Population genomics of domestic and
  1165 wild yeasts. Nature 458, 337-341. 10.1038/nature07743.
- Strope, P.K., Skelly, D.A., Kozmin, S.G., Mahadevan, G., Stone, E.A., Magwene, P.M.,
  Dietrich, F.S., and McCusker, J.H. (2015). The 100-genomes strains, an S. cerevisiae
  resource that illuminates its natural phenotypic and genotypic variation and emergence as
  an opportunistic pathogen. Genome Res 25, 762-774. 10.1101/gr.185538.114.
- 117040.Richardson, S.M., Wheelan, S.J., Yarrington, R.M., and Boeke, J.D. (2006). GeneDesign:1171rapid, automated design of multikilobase synthetic genes. Genome Res 16, 550-556.117210.1101/gr.4431306.
- 1173 41. Robinson, J.T., Thorvaldsdottir, H., Winckler, W., Guttman, M., Lander, E.S., Getz, G., and
  1174 Mesirov, J.P. (2011). Integrative genomics viewer. Nat Biotechnol 29, 24-26.
  10.1038/nbt.1754.
- 117642.Bolger, A.M., Lohse, M., and Usadel, B. (2014). Trimmomatic: a flexible trimmer for Illumina1177sequence data. Bioinformatics *30*, 2114-2120. 10.1093/bioinformatics/btu170.
- 117843.Andrews, S. (2010). FastQC: a quality control tool for high throughput sequence data.1179http://www.bioinformatics.babraham.ac.uk/projects/fastqc/.
- 118044.Langmead, B., and Salzberg, S.L. (2012). Fast gapped-read alignment with Bowtie 2. Nat1181Methods 9, 357-359. 10.1038/nmeth.1923.
- Wickham, H. (2016). ggplot2: Elegant Graphics for Data Analysis (Springer-Verlag New York).
- 118446.Bray, N.L., Pimentel, H., Melsted, P., and Pachter, L. (2016). Near-optimal probabilistic1185RNA-seq quantification. Nat Biotechnol *34*, 525-527. 10.1038/nbt.3519.
- 47. Dobin, A., Davis, C.A., Schlesinger, F., Drenkow, J., Zaleski, C., Jha, S., Batut, P.,
  Chaisson, M., and Gingeras, T.R. (2013). STAR: ultrafast universal RNA-seq aligner.
  Bioinformatics 29, 15-21. 10.1093/bioinformatics/bts635.
- 48. Li, H., Handsaker, B., Wysoker, A., Fennell, T., Ruan, J., Homer, N., Marth, G., Abecasis,
  G., Durbin, R., and Genome Project Data Processing, S. (2009). The Sequence
  Alignment/Map format and SAMtools. Bioinformatics 25, 2078-2079.
  10.1093/bioinformatics/btp352.
- 119349.Quinlan, A.R., and Hall, I.M. (2010). BEDTools: a flexible suite of utilities for comparing1194genomic features. Bioinformatics 26, 841-842. 10.1093/bioinformatics/btq033.
- 119550.Pimentel, H., Bray, N.L., Puente, S., Melsted, P., and Pachter, L. (2017). Differential1196analysis of RNA-seq incorporating quantification uncertainty. Nat Methods 14, 687-690.119710.1038/nmeth.4324.

- Brachmann, C.B., Davies, A., Cost, G.J., Caputo, E., Li, J., Hieter, P., and Boeke, J.D.
  (1998). Designer deletion strains derived from Saccharomyces cerevisiae S288C: a useful set of strains and plasmids for PCR-mediated gene disruption and other applications. Yeast 14, 115-132. 10.1002/(SICI)1097-0061(19980130)14:2<115::AID-YEA204>3.0.CO;2-2.
- 120252.DiCarlo, J.E., Norville, J.E., Mali, P., Rios, X., Aach, J., and Church, G.M. (2013). Genome1203engineering in Saccharomyces cerevisiae using CRISPR-Cas systems. Nucleic acids1204research 41, 4336-4343. 10.1093/nar/gkt135.
- 1205 53. Litovchick, L. (2020). Immunoblotting. Cold Spring Harb Protoc 2020, 098392.
   10.1101/pdb.top098392.
- 1207 54. Brooks, A.N., Hughes, A.L., Clauder-Munster, S., Mitchell, L.A., Boeke, J.D., and Steinmetz,
  1208 L.M. (2022). Transcriptional neighborhoods regulate transcript isoform lengths and
  1209 expression levels. Science 375, 1000-1005. 10.1126/science.abg0162.
- 1210 55. Domanski, M., Molloy, K., Jiang, H., Chait, B.T., Rout, M.P., Jensen, T.H., and LaCava, J.
  1211 (2012). Improved methodology for the affinity isolation of human protein complexes
  1212 expressed at near endogenous levels. Biotechniques 0, 1-6. 10.2144/000113864.
- 1213 56. Luo, J., Sun, X., Cormack, B.P., and Boeke, J.D. (2018). Karyotype engineering by
  1214 chromosome fusion leads to reproductive isolation in yeast. Nature 560, 392-396.
  10.1038/s41586-018-0374-x.

1216