



Impurity analysis of synthetic diamond for electronics and quantum physics

E. Loire, François Jomard, M.A. Pinault-Thaury

► To cite this version:

E. Loire, François Jomard, M.A. Pinault-Thaury. Impurity analysis of synthetic diamond for electronics and quantum physics. 24th International Conference on Secondary Ion Mass Spectrometry (SIMS-24), Sep 2024, La Rochelle, France. hal-04749557

HAL Id: hal-04749557

<https://hal.science/hal-04749557v1>

Submitted on 23 Oct 2024

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.





L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

Impurity analysis of synthetic diamond for electronics and quantum physics

Estelle LOIRE, François JOMARD, Marie-Amandine PINAULT-THAURY

Groupe d'Etude de la Matière Condensée (GEMaC-UMR8635),
CNRS, Université de Versailles St-Quentin-En-Yvelines (UVSQ), Université Paris-Saclay
45 av des Etats-Unis - 78035 Versailles Cedex - France


Diamond is the ultimate ultra-wide band gap semiconductor.
Its synthesis is mainly based on its impurity control.



SIMS analysis are performed with an CAMECA IMS7f dedicated to materials (metals, oxides, semiconductors...).




Our diamond conventional conditions



Element	Interferences	Mass resolution
¹ H	None	400
¹¹ B	None	400
¹³ C	¹² C ¹ H	3000
¹⁴ N ¹² C	¹³ C ₂ , ¹² C ¹³ C ¹ H	7200
²⁸ Si	¹² C ¹⁶ O, ¹² C ² ¹ H ₄	1600
³¹ P	¹² C ¹⁹ F, ¹² C ¹ H ₃ ¹⁶ O	1600
⁷⁵ As	¹² C ₆ ¹ H ₃ ¹² C ₅ ¹⁴ N ¹ H	3000


With more than 20 years of diamond synthesis and physical studies in GEMaC, the SIMS analysis is well controlled over a large range of elements with the use of dedicated implanted standards and lab samples.

High mass resolution (HMR) for ¹³C, ²⁸Si, ³¹P and ⁷⁵As




Determination of the relative sensitivity factors (RSF)

HMR & Raster size method for ¹⁴N



Element	¹ H/ ¹² C	¹¹ B/ ¹² C	¹⁴ N/ ¹² C/ ¹² C	²⁸ Si/ ¹² C	³¹ P/ ¹² C	⁷⁵ As/ ¹² C
RSF (at/cm ³) ± 10%	2.8×10^{23}	5.2×10^{24}	2.8×10^{22}	2.7×10^{23}	2.5×10^{23}	4.3×10^{23}

Low mass resolution (LMR) for ¹H and ¹¹B



Our usual detection limits (DL) by taking care of the residual vacuum of the chamber

	DL (at/cm ³)
[H]	$1.5 \pm 0.5 \times 10^{18}$
[B]	$1.5 \pm 0.5 \times 10^{16}$
[N]	$5.0 \pm 1.0 \times 10^{16}$
[Si]	$2.5 \pm 1.0 \times 10^{16}$
[P]	$4.5 \pm 1.5 \times 10^{14}$
[As]	$1.5 \pm 0.5 \times 10^{16}$

Contact: marie-amandine.pinault-thaury@uvsq.fr & estelle.loire2@uvsq.fr