

INSTITUT COCHIN ENTITY DMP - PIME Core Facility

Alain Schmitt, Thomas Guilbert

▶ To cite this version:

Alain Schmitt, Thomas Guilbert. INSTITUT COCHIN ENTITY DMP - PIME Core Facility. Institut Cochin, Inserm U1016, CNRS UMR 8104, Université Paris Cité, F-75014 Paris, France. 2024. hal-04846443

HAL Id: hal-04846443 https://hal.science/hal-04846443v1

Submitted on 18 Dec 2024

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Distributed under a Creative Commons Attribution 4.0 International License

INSTITUT COCHIN ENTITY DMP - PIME

Plan de gestion de données créé à l'aide de DMP OPIDoR, basé sur le modèle "Science Europe: structured template for research entities" fourni par DMP OPIDoR.

Plan Details

Plan title	INSTITUT COCHIN	ENTITY DMP - PIME
Deliverable		
Version	First version	
Plan purpose/scope	The aim of this DM research products p This DMP is part of DMPs of the 9 Insti PROTEOM'IC, PIV, and HISTIM) as we transfer, storage, c data produced by th The purpose of this researchers wishing wishing to call on th	P is to describe as precisely as possible the various produced by the services offered by the core facility. a wider framework, that of making available all the tut Cochin core facilities (MOUSET'IC, PIME, METABOL'IC, CYBIO, IMAG'IC, GENOM'IC II as the DMP of the CID database dedicated to the onsultation and sometimes processing of most of the ne Institut Cochin's core facilities. DMP entity is to provide drafting assistance for g to draw up a DMP for their research projects and ne expertise of our core facilities.
Fields of science and technology (from OECD classification)	1.6 Biological scien	ces
Language	eng	
Creation date	2024-07-18	
Last modification date	2024-12-18	
Identifier type	DOI	
License	Name URL	Creative Commons Attribution 4.0 International http://spdx.org/licenses/CC-BY-4.0.json

Structure Details

Entity's name	Institut Cochin Core Facility - PIME
Acronym	PIME
Identifier	https://institutcochin.fr/plateformes/pime
Identifier type	URL
Description	La plateforme de microscopie électronique est spécialisée dans l'analyse ultrastructurale des molécules biologiques, des cellules et des tissus, depuis la recherche fondamentale à la recherche appliquée. Nos domaines d'expertise couvrent la description et la caractérisation des tissus normaux et pathologiques, et l'identification et la localisation de protéines. Nos techniques sont complémentaires de la microscopie optique traditionnelle, nous permettant un pouvoir séparateur de l'ordre du nanomètre. Ce qui est particulièrement important en microscopie électronique, nous réalisons d'office l'aide à l'analyse, à l'interprétation et à la présentation des résultats et des images.
Creation date	2003-02-02

Research outputs :

- Resin Block (Physical object)
 EM Images (Image)

Contributors

Name	Affiliation	Roles
Guilbert Thomas - 0000-0001- 5069-0730	Institut Cochin - 051sk4035	• DMP manager
Schmitt Alain - 0000-0002- 8011-2277	Institut Cochin - 051sk4035	 Contact Person (ResinBlock, Images) Entity's Manager (Images, Images, ResinBlock, Images, Images, ResinBlock, Images, Images, Images, ResinBlock, Images, ResinBlock, ResinBlock, ResinBlock, ResinBlock, ResinBlock, ResinBlock) Entity's manager

INSTITUT COCHIN ENTITY DMP - PIME

Data description and collection or re-use of existing data

Resin Block		
Research output descrip	otion	
Name	Resin Block	
Description	Resin block generated from biological sample. Each samble is embedded in electron microscopy designed resins.	
Туре	Physical object	
Workpackage		
Issued Date	2024-07-22	
Persistent identifier	"project number"-"sample number"	
Identifier type	Identifiant local	
May contain personal data?	No	
May contain sensible data?	No	
May take ethical issues into account?	No	
Will existing data be reused?		
Justification	Yes, for futur research.	
How new data will be collected or produced?		

Name of the method	Embedding
Description	Embedding Protocol : Secondary fixation was performed with 1% osmium tetroxide in phosphate buffer saline for 1 hr at room temperature, then washed 3 times in H2O and further dehydrated with increasing concentrations of ethanol (70%, 90%, and 100% twice). After dehydration, samples were infiltrated with sequential 3:1, 1:1, and 1:3
	solutions of resin (Embed 812 kit; Electron Microscopy Sciences, Hatfield, PA) and 100% ethanol for 1 h separately at room temperature.
	Finally, the samples were incubated with propylene oxide for 2mn and and placed in the embedding molds with resin. The embedding resin blocks were polymerize in an oven at 60°C for 24h.
Data Nature	Experimental Data
Equipments, technical platforms used	• PIME : https://cat.opidor.fr/index.php/PIME

EM Images

Research output description

Name	EM Images
Description	Acquisition of region of interest. Analysis of the aspects of the samples in light of the problems and research of the researchers. Output consist in information linked to scientific questions
Туре	Image
Workpackage	
Keywords (free-text)	
Language	eng
Issued Date	
Persistent identifier	"project number"-"sample number"-"Unique number"
Identifier type	Identifiant local
May contain personal data?	No
May contain sensible data?	No
May take ethical issues into account?	Yes

Will existing data be reused?

Justification

Data can be analyse again at the light of new scientific questions.

How new data will be collected or produced?	
Name of the method	EM acquisition
Description	Acquisition software (digital micrograph) on the electron microscope.
Data Nature	Experimental data
Equipments, technical platforms used	• PIME : https://cat.opidor.fr/index.php/PIME

Documentation and data quality

Resin Block		
What metadata and doc data?	umentation (for example way of organising data) will accompagny the	
Description	Resine Block are referenced in the ElabOne-Pime.	
Metadata language code	fra	
Documentation software	ElabFTW	
What methods will be us	sed to ensure their scientific quality?	
Description	Visual inspection.	
EM Images		
What metadata and documentation (for example way of organising data) will accompagny the data?		
Description	Coding in the laboratory notebook.	
What methods will be used to ensure their scientific quality?		

Legal and ethical requirements, codes of conduct

Resin Block		
How will other legal issues, such as intellectual property rights and ownership, be managed? What legislation is applicable?		
Description	Researchers deposing the samples assure that all legal issues are managed by accepting the rules. For private compagny, following established contract (destruction or redonner les blocs), following NDA.	
What ethical issues and	codes of conduct are there, and how will they be taken into account?	
Description	following INSERM, CNRS and UPCité rules following APHP rules for human samples (anonymizasation, legal consent, etc)	
EM Images How will other legal issues, such as intellectual property rights and ownership, be managed? What legislation is applicable?		
Description	Researchers deposings the samples assure that all legal issues are managed by accepting the rules. For private compagny, following established contract (destruction or redonner les blocs) following NDA.	
What ethical issues and codes of conduct are there, and how will they be taken into account?		
Description	following INSERM, CNRS and UPCité rules following APHP rules for human samples (anonymization, legal consent, etc)	

Data processing and analysis

Resin Block	
How and with what reso	ources will the data be processed / analyzed?
Description	Production of sample section used for the next ME experiment. Ultrathin sections (90 nm) were cut from the block surface using an ultramicrotome (Reichert Ultracut S ultramicrotome, Leica Microsystems) then stained with uranyl acetate and lead citrate (LFG, France) and observed by TEM (JEOL 1011, JEOL, Japan)) at 80 kV of acceleration voltage. Acquisitions were performed by an Orius 1000 CCD Camera (GATAN, USA) using Digital Micrograph software (GATAN, USA).
Related references	• Ultramicrotomy of biological objects: From the beginning to the present : https://doi.org/10.1016/1047-8477(90)90074-M
Equipments, technical platforms	• PIME : https://cat.opidor.fr/index.php/PIME
EM Images	
How and with what reso	ources will the data be processed / analyzed?
Description	Data are visually analysed upon acquisition. In Electron microscopy, images provide information on size, shape, morphology and ultrastructure of all the elements of the biological sample, at a nanometric scale. Analysis is mostly qualitative, albeit some quantitative analysis could be performed with a suffisant sample size. In this case, the commonly used image analysis software equipped with an electron microscopy module can be used.
Equipments, technical platforms	• PIME : https://cat.opidor.fr/index.php/PIME

Storage and backup during the research process

Resin Block

How will data be stored and backed up during the research?

Storage needs	Dedicated container in the lab.
Measures taken for data security	Securised software (ELAB-FTW) referencing the block (Anonymous coding).
EM Images	
How will data be stored	and backed up during the research?
Storage needs	Hard drives connected to the acquisition system or CID.
Estimated volume of data	80
Unit	MB
Equipments, technical platforms	• CID : https://cat.opidor.fr/index.php/CID
Measures taken for data security	Securised software (ELAB-FTW) referencing the block (Anonymous coding).

Data sharing and long-term preservation

Resin Block		
How will data be shared	1?	
Modalities of sharing	Physical objects, so physical sharing.	
Reusability	until the block is physically exhausted.	
Data repository/catalogs	• CID : https://cat.opidor.fr/index.php/CID ()	
How will data be long-term preserved? Which data?		
Justification	Could always be useful for futur resarch. On demand or until the project is terminated.	
Start date		
End date		
Final dispositions		

EM Images	
How will data be shared?	
Modalities of sharing	By file transfering or CID.
Reusability	Future research. We recommend the use of open science sharing for future data sharing.
Data repository/catalogs	• CID : https://cat.opidor.fr/index.php/CID ()
Justification	On demand or until the project is terminated.
Estimated volume of data	80
Unit	MB
Start date	
End date	
Final dispositions	