

COSMO-RS as a cutting-edge technology for designing and pre-screening of Natural Deep Eutectic Solvents for selective extraction: Spirulina as a case study

Mervé Yagmur, Alia Chevalley, Bénédicte Montigny, Johan Jacquemin, Leslie

Boudesocque-Delaye

▶ To cite this version:

Mervé Yagmur, Alia Chevalley, Bénédicte Montigny, Johan Jacquemin, Leslie Boudesocque-Delaye. COSMO-RS as a cutting-edge technology for designing and pre-screening of Natural Deep Eutectic Solvents for selective extraction: Spirulina as a case study. International Congress on Natural Products Research (ICNPR 2024), Jul 2024, Cracow, Poland. hal-04951821

HAL Id: hal-04951821 https://hal.science/hal-04951821v1

Submitted on 17 Feb 2025


HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



COSMO-RS as a cutting-edge technology for designing and pre-screening of Natural Deep Eutectic Solvents for selective extraction: Spirulina as a case study

Mervé Yagmur^{a,*}, Alia Chevalley^b, Bénédicte Montigny^c, Johan Jacquemin^d and Leslie Boudesocque-Delaye^a

^a EA 7502 SIMBA, Pharmacy Faculty, University of de Tours, 31 Avenue Monge, 37200 Tours, France; ^b Aqua Eco Culture, 7 Rue d'Armor Maroué, 22400 Lamballe, France; ^c EA 6299 PCM2E, University of Tours, Parc de Grandmont, 37200 Tours, France; ^d MSN Department, Mohammed VI Polytechnic University, Lot 660, Hay Moulay Rachid, Ben Guerir 43150, Morocco



Green extraction 01 approach

Consider not only the biomass but also the extractive agents used at each step of the overall process

used as a probe - In silico screening of potential NaDES for selectively extracting targeted high-value compounds - Experimental validation - Eutectic window

Conventional NaDES conception 03 Experimental screening of different ratios:


a time- and solvent-consuming process

of each conformer

the solid-liquid equilibria (SLE)


rehydrated flakes






The dried biomass was rehydrated in each tube prior to running the extraction process. This was done to collect reproducible data.

Simulation for extractive solvent selection




Selected NaDES according to targeted compounds



Eutectic point determination


The SLE determination allowed estimation of each eutectic ratios (dark grey).




C8 was selected for each NaDES combination based on its high FFA extraction capacity.¹ Co-solvents were chosen based on the complementarity of their sigma profiles with C8. This strategy was also used during the *in silico* screening of NaDES formulation to efficiently extract target compounds.

Experimental results

Evaluation of the in silico NaDES selection



For each NaDES, two additional molar ratios were also investigated to experimentally evaluate the effect of the extractant formulation on the process selectivity toward targeted compounds (light grey).



Effect of the extractant formulation

■ SFA ■ PUFA ■ Total Pigments mg of target compounds/g of dry matter

These diagrams represent the 3 groups of target compounds for the different mixtures at the eutectic ratios.

 \checkmark The experimental results are consistent with the simulation hypothesis.

Conclusions and perspectives

A driver to extract nonpolar compounds from Spirulina:

Screening done according to **COSMO-RS to reach high selectivity** • SFA: C8C10 likely most preferable

- **PUFA**: C8M; C8phen & C8PDiol •
- **Pigments**: C8phen & C8PDiol due to charge delocalization

Experimentally:

- C8C10 has a higher extraction potential for SFA
- Both C8phen and C8PDiol could be used to extract nonpolar pigments and PUFA
- C8M is interesting for extracting PUFA with a low pigments content

| Total FFA C8C10 | C8M | C8phen | C8PDiol | Total | C8C10 | C8M | C8phen | C8PDiol |
|-----------------|-----|--------|---------|---------------------------------|-------|-----|-----------|---------|
| SFA PUFA | | | | pigments ■ Chlorophylls ■ Ca | | | rotenoids | |

The extractant formulation must be fine-tuned to reach a selective extraction of a specific compound. Eutectic compositions are often less effective due to strong interactions between the solvent molecules.

Eutectic effect

- Each eutectic point was determined by computing the SLE using COSMO-RS
- Eutectic composition is not always the best to be used when formulating a NaDES
- The extraction capability and selectivity can be optimized by using the sigma profiles of solvents along with their ratios to be used as extractant based on calculated SLE

Perspectives

This work will help us to go further selective in extraction of some compound to minimize the amount of the non-targeted compound categories.

