

Electrophysiological evidence for audio-visuo-lingual speech integration

Coriandre Emmanuel Vilain, Avril Treille, Marc Sato

▶ To cite this version:

Coriandre Emmanuel Vilain, Avril Treille, Marc Sato. Electrophysiological evidence for audio-visuo-lingual speech integration. IMRF 2015 - 16th International Multisensory Research Forum, Jun 2015, Pise, Italy. hal-01297678

HAL Id: hal-01297678 https://hal.science/hal-01297678

Submitted on 5 Apr 2016

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



ELECTROPHYSIOLOGICAL EVIDENCE FOR AUDIO-VISUO-LINGUAL SPEECH INTEGRATION

Coriandre Vilain¹, Avril Treille¹ & Marc Sato²

¹ GIPSA-lab, Département Parole & Cognition, CNRS & Grenoble Université, Grenoble, France ²Laboratoire Parole & Langage, CNRS & Aix Marseille Université, France

Introduction

Audio-visual speech perception is a special case of multisensory processing that interfaces with the linguistic system. One important issue is whether cross-modal interactions only depend on well-known auditory and visuo-facial modalities or, rather, might also be triggered by other sensory sources less common in speech communication. The present EEG study aimed at investigating cross-modal interactions not only between auditory, visuo-facial and audio-visuo-facial syllables but also between auditory, visuo-lingual and audio-visuo-lingual syllables.

✤ <u>Participants</u>

18 healthy adults, right-handed native French speakers.

Stimuli

/pa/, /ta/ and /ka/ syllables were individually recorded by one male and one female speakers in a sound-proof room. Two /pa/, /ta/ and /ka/ tokens were selected per speaker .



Data acquisition

- EEG data were continuously recorded from 64 scalp electrodes (international 10–20 system) using the Biosemi ActiveTwo AD-box EEG system operating at a sampling rate of 256 Hz.
- Two additional electrodes served as reference [CMS] & [DRL]
- One other external reference electrode was at the top of the nose. The electrooculogram controlling for horizontal (HEOG) and vertical (VEOG) eye movements were recorded using electrodes at the outer canthus of each eye as well as above and below the right eye. Before the experiment, the impedance of each electrode was adjusted to get low offset voltage and stable DC.

✤ Analysis

Synchronous recordings of auditory, visual and ultrasound lingual signals were acquired by a Terason T3000 ultrasound system Simultaneous recording : left : tongue movements with ultrasound system; (UltraSpeech system, Hueber et al. 2008).

right : Lips movements with camera

• 60 stimuli were created consisting of 12 /pa/, /ta/ and /ka/ syllables related to 5 conditions: an auditory condition (A), 2 visual (VL, VT) and 2 audio-visual (AVL, AVT) conditions related to either lip or tongue movements of a speaker.

Method

- 1) Before the experiment, short explanations on the ultrasound system and tongue movements required for the production of **/pa/**, **/ta/** and **/ka**/ syllables were given.
- 2) Behavioral session: a three-alternative forced-choice identification task, with participants instructed to categorize as quickly as possible each perceived syllable with their right hand.
- 3) EEG session: a three-alternative forced-choice identification task, with participants instructed to categorize each perceived syllable with their right hand, after an audio "beep".

Behavioral analyses (1):

	Behavioral experiment	EEG experiment
%	ANOVA : Environment (Noise/No Noise), Articulator (Lips/Tongue), Modality (A, AV, V)	ANOVA : Environment (Noise/No Noise), Articulator (Lips/Tongue), Modality (A, AV, V)
RT's	ANOVA : Environment (Noise/No Noise), Articulator (Lips/Tongue), Modality (A, AV, V)	

EEG analyses on fronto-central electrodes (F3/F4/C3/C4/Fz/Cz):

Pre-	- Re-referenced off-line to the nose	- Filtering : 1-30 Hz
processing	- Epochs : 1000ms (baseline from -500 to -400ms	- Rejection : ±60 µV
	to the acoustic synaple onset	

→ 1) N1 and P2 latency and amplitude analysis by ANOVA with Environment (Noise/No noise), Articulator (Lips/Tongue) and Modality (A, AV) as intra-subject variable (main effects and interactions)

2) EEG – Noise vs No Noise

→ 2) N1 and P2 latency and amplitude were also tested for integration, i.e. an ANOVA with Environment (Noise/No noise), Articulator (Lips/Tongue), Signal type (Bimodal (AV)/Sum (A+V)) as intra-subject variable (main effects and interactions)



1) Behavioral results

% Correct Response :

- Lips > tongue (p<.005) - AV > A = V (p<.0001)

- Lips: A=V < AV **Tongue:** AV=A > V (p<.05)



reflect the impact of listener's knowledge of speech production on speech perception.

Part of this research was supported by a grant from the European Research Council (FP7/2007-2013 Grant Agreement no. 339152, "Speech Unit(e)s") Correspondence: avril.treille@gipsa-lab.inpg.fr