



**HAL**  
open science

## Rheological monitoring for tau protein aggregation process

Pierre Didier, Florence Razan, Emmanuel Caplain, Christian Delamarche,  
Stephane Serfaty, Pascal Larzabal

► **To cite this version:**

Pierre Didier, Florence Razan, Emmanuel Caplain, Christian Delamarche, Stephane Serfaty, et al..  
Rheological monitoring for tau protein aggregation process. Biosensors 2016, May 2016, Göteborg,  
Sweden. hal-01393739

**HAL Id: hal-01393739**

**<https://hal.science/hal-01393739>**

Submitted on 25 Jan 2023

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

# Rheological monitoring for tau protein aggregation process

Pierre Didier<sup>\*1</sup>, Florence Razan<sup>1</sup>, Emmanuel Caplain<sup>1</sup>, Magalie Michiel<sup>1</sup>, Christian Delamarche<sup>2</sup>, Stéphane Serfaty<sup>1</sup>, Pascal Larzabal<sup>1</sup>


1. Laboratoire SATIE-UMR CNRS 8029-ENS Cachan, 61 avenue du Président Wilson – 94235 Cachan, France.

2. Laboratoire IGDR-UMR CNRS 6290 - Campus de Beaulieu, bât. 13, 35042 Rennes, France.

\* [pierre.didier@ens-rennes.fr](mailto:pierre.didier@ens-rennes.fr)

## Introduction


In the human body there are many proteins for which the polymerization can induce pathology. Diseases result from abnormalities in the structure and conformation of proteins, called conformational disease. Tau protein aggregation (Tubulin-Associated Unit) is associated with many neurological disorders grouped under the term Tauopathies. It was originally described as the "Tau factor" means an element capable of inducing polymerization of tubulin into microtubules. Neurofibrillary tangles are amassed from filaments having a straight or paired helical morphology, they are composed of the microtubule-associated protein tau under hyperphosphorylation form. Death of polymer-laden neurons is evidenced by the presence of insoluble tangle remnants to induce degeneration, they are implicated in the necrotic process of synapses. They will aggregate and probably cause cell death. The aggregation's mechanism is poorly understood and difficult to study.



In this work, a microsystem is developed for tau protein detection based on real time monitoring micro-rheological parameters by high frequency ultrasonic waves.

## Measurement principle



Tau protein polymerization is followed thanks to the thickness shear mode resonator developed in this study. The acoustic wave transducer has a resonance frequency due to its intrinsic properties (materials type, quartz dimensions, crystal cut) and mechanical properties, which will vary with biological species deposition on its surface. Complex impedance is measured by a vectorial network analyzer after introducing piezoelectric transducer in a thermostated cell during the experiment.




## Transducer characterization: Electrical response

During a deposition of chemical or biological species to the quartz surface, its resonance frequency and its electrical properties are changed and characterized by a network analyzer.

Biological species deposition is characterized by complex impedance  $Z_M$ , which depend on the viscoelastic parameters of tested solution.



## Electrical model associated with the transducer



## Extraction of viscoelastic parameters

Elastic modulus ( $G'$ ) and Viscous modulus ( $G''$ ) are extracted from  $Z_{load}$ :

$$Z_{Load} \approx i Z_0 \frac{(2n+1)\pi}{4K_{eff}^2} \frac{Z_M}{Z_{m\_eff}} \approx K \left( \sqrt{\rho_M (G'(\omega) + iG''(\omega))} + i\omega \rho_M h_{eff} \right) = R + iX + iK\omega \rho_M h_{eff}$$

Rheological parameters ( $R, X$ ) or ( $G', G''$ ) give information on the evolution of the material


$G' \propto \frac{R^2 - X^2}{\rho_M}$  elastic modulus

$G'' \propto \frac{R \times X}{\rho_M}$  viscous modulus


## Viscoelastic parameters $G'$ and $G''$ for a pure water and solution of water/glycerol:

Measurements revolves around three phases:


- **First phase** consisted in **system initialization**, with the determination of the sensor electrical parameters and quality factor.
- **Second phase** consisted in deposition of a purely Newtonian liquid (DI water) on the quartz surface to **calibrate experiment**.
- **Third phase** is the **experimental phase**. the Newtonian liquid is removed and biological solution is deposited on the quartz substrate measurement for the record in line of complex impedance.



Complex impedance measure




Viscoelastic parameters  $G'$  and  $G''$  extracted with electrical model



## Viscoelastic parameters $G'$ and $G''$ for tau protein polymerization


Protein sequence : VQIVYK





Conformational states of protein depending on peptide concentration

The choice of this peptide is justified by its presence in the tau protein in position 306-311 and is sufficient for the self-assembly of tau



Biological model of the polymerization of the tau protein. AFM micrographies of protein at different concentrations. (a) 30  $\mu$ M, (b) 120  $\mu$ M, (c) 240  $\mu$ M and (d) 480  $\mu$ M.



-> Different concentrations of protein are deposited on the quartz. Experimental temperature: 22°C (in vitro conditions), 37°C (in vivo conditions)



Increase of  $G''$  parameter: from core to fibers, viscosity increases -> Information on the concentration of the solutions tested.



Variations of  $G'$  parameter: follow proteins polymerizations and detection of peptides conformations changes are possible.

## Conclusions

Development of an instrumentation system allowing:

- Detection of different concentrations of proteins in static conditions
- Detection of different conformational states with the plasticity differentiation
- Differentiation between protein concentrations of tau ranges from 30-480  $\mu$ M by measuring their viscoelastic properties with confirmation of the biological model.

