



## Molecules, Radicals and Ions produced in an N<sub>2</sub>-H<sub>2</sub> CCP RF plasma

Nathalie Carrasco, David Dubois, Audrey Chatain, Ludovic Vettier, Guy Cernogora

### ► To cite this version:

Nathalie Carrasco, David Dubois, Audrey Chatain, Ludovic Vettier, Guy Cernogora. Molecules, Radicals and Ions produced in an N<sub>2</sub>-H<sub>2</sub> CCP RF plasma. ICPIG 2017, International Conference on Phenomena in Ionized Gases, Jul 2017, Lisbon, Portugal. hal-01486530

**HAL Id: hal-01486530**

**<https://hal.science/hal-01486530>**

Submitted on 10 Mar 2017

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

# Molecules Radicals and Ions produced in a N<sub>2</sub>-H<sub>2</sub> CCP RF

N. Carrasco<sup>1</sup>, D. Dubois<sup>1</sup>, A. Chatain<sup>1</sup>, L. Vettier<sup>1</sup>, G. Cernogora<sup>1</sup>


<sup>1</sup> LATMOS, Université Versailles St-Quentin, CNRS, 11 blvd d'Alembert, 78280 Guyancourt, France

CCP RF discharges are well known to be sources of dusty plasmas. These plasmas are used to simulate the formation of organic solid particles in planetary atmospheres, as Titan with a N<sub>2</sub> CH<sub>4</sub> mixture. As a first step for understanding these plasmas, we study here the formation of molecules, radicals and positive ions in a N<sub>2</sub> H<sub>2</sub> CCP RF plasma. mixture. Radicals and positive ions are measured by in situ mass spectrometry. Neutrals are accumulated in a cold trap downstream the plasma. These molecules are measured, after warming the trap by mass spectrometry and IR absorption spectroscopy. . When mass spectrometry gives relative values of species abundances, IR absorption gives absolute values of the most abundant molecules. A focus is done on NH<sub>3</sub>, this molecule being produced as well in the discharge as by catalytic effect on the metallic wall of the discharge.

## 1. Introduction

A Capacitively Coupled Plasma in N<sub>2</sub>-CH<sub>4</sub> mixture is used for the formation of dust to simulate the formation of solid aerosols in Titan's atmosphere. Details of this experimental device are described in details in [1]. The dissociation of CH<sub>4</sub> produces H<sub>2</sub> molecules. In order to understand the complex chemistry occurring in the N<sub>2</sub> CH<sub>4</sub> mixture, study of the N<sub>2</sub> H<sub>2</sub> mixture is developed as a first step.

## 2. Experimental device.



Experimental device

The CCP discharge is confined in a metallic cylindrical box of 13.6 cm in diameter. Two slots, diametrically opposed are managed in the cylindrical box. This box is placed into a 30 cm in diameter and 40cm height stainless steel reactor fitted with two KBr windows diametrically opposed. Gas mixture is injected continuously in the plasma and pumped with a rotary pump. The amount of H<sub>2</sub> in N<sub>2</sub> is tuned from 1 to 5%. Pressure in the plasma discharge is maintained at 1 mbar. Between the discharge and the pump, a trap cooled by liquid nitrogen condenses molecules produced in the discharge.

## 3. Measurements.

The neutral molecules, radicals and positive ions are measured using a mass spectrometer EQP Hidden placed in front of one of the slot.

After few hours of plasma run, the gas injection is turned off, valves upstream and downstream the cold trap are closed and the trap is slowly warmed up to the room temperature. Then the condensed gases are reinjected in the reactor.

The molecules densities are measured using a Nicolet 6700 Fischer FTIR spectrometer through the 40 cm of in diameter of the reactor [2]. Focus is done on the NH<sub>3</sub> molecule measured in the 967 cm<sup>-1</sup> band. In our experimental conditions, the NH<sub>3</sub> density is on the order of 10<sup>12</sup> cm<sup>-3</sup>.

## 4- Perspectives

These results obtained are now to be compared with the modelling of the plasma in our experimental conditions. The major point is the relative contribution of volume reactions versus catalytic ones for the formation of NH<sub>3</sub>.

## 5. References

- [1] Alcouffe, G. , M. Cavarroc, G. Cernogora, F. Ouni, A. Jolly, L. Boufendi and C. Szopa (2010), Plasma Sources Sci. Technol. 19(1): 015008.
- [2] Dubois D., Carrasco N., Petrucciani M., Tigrine S., Vettier L *Neutral Chemistry in Titan's Ionospheric Simulated Conditions*, Dubois D., et al., DPS-EPSC 2016, Oct 2016, Pasadena, USA