



HAL
open science

Expression of laccase IIIb from the white-rot fungus *Trametes versicolor* in the yeast *Yarrowia lipolytica* for environmental applications

Claude Jolival, Catherine Madzak, Agathe Brault, Eliane Caminade, Christian Malosse, Christian Mougin

► To cite this version:

Claude Jolival, Catherine Madzak, Agathe Brault, Eliane Caminade, Christian Malosse, et al.. Expression of laccase IIIb from the white-rot fungus *Trametes versicolor* in the yeast *Yarrowia lipolytica* for environmental applications. 2nd European Meeting in Oxizymes, Jun 2004, Naples, Italy. 2004. hal-01600125

HAL Id: hal-01600125

<https://hal.science/hal-01600125>

Submitted on 4 Jun 2020

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Distributed under a Creative Commons Attribution - ShareAlike 4.0 International License

INTRODUCTION

- Fungal laccases (EC 1.10.3.2) are involved in lignin transformation, morphogenesis, pathogenesis and fungal virulence, as well as in the oxidation of numerous xenobiotics [1-3]
- A wide range of their substrates are of environmental interest, and are recognized pollutants of waters and soils
- Laccases are potential useful tools to improve the chemical quality of polluted media through bioremediation [1,3]

Our objective is the design of laccases with:





- a high efficiency in catalysis (high redox potential)
- a wider variety of substrates
- a pH of activity optimized for natural media


Our approach is based on the knowledge of the cavity enclosing a weak reducing substrate of the enzyme, the arylamine 2,5-xylydine, and on the replacement of neighboring interacting residues

RESULTS: Expression of LacIIIb by *Yarrowia lipolytica*


- The gene encoding LacIIIb has been cloned in *Escherichia coli* using the pGEM®-T plasmid
- It has been sequenced and comprises 1563 bases (sequence accession number AF414109), including a peptide signal
- It was then expressed in the yeast *Yarrowia lipolytica* using the pINA1269 vector (Patent INRA/INA-PG), thus allowing the production of an active form of the recombinant enzyme (transformation of the substrate ABTS to a green oxidized form on solid culture media and oxidation of guaïcol after electrophoretic separation). After tryptic digestion of the transformed *Y. lipolytica* extract, the mass spectrum contains the most significative ions detected and assigned to *T. versicolor* laccase peptides.



Laccase activity detection on agar plates assay. Numbers refer to different transformant clones of *Y. lipolytica* Po1g strain having received p69TVLB plasmid.



Electrophoresis analysis of recombinant laccase IIIb. [A] Samples were boiled before loading on the gel. Lane 1 - Laccase IIIb from *T. versicolor*. Lane 2 - Concentrated crude cell free extract from a transformant (clone 4) of *Yarrowia lipolytica* expressing laccase IIIb. [B] Samples were not boiled before loading on the gel. Lanes 3 and 6 - Laccase IIIb from *T. versicolor*. Lanes 4 and 7 - Concentrated crude cell free extract from a transformant (clone 4) of *Yarrowia lipolytica* expressing laccase IIIb. Lane 5 - Molecular weight markers. Lanes 1 to 5 were stained with Coomassie bleu (R250). Laccase activity was detected in lanes 6 and 7 by coloration with guaïcol.




MALDI-TOF spectra [A] Laccase IIIb from *T. versicolor*. [B] Analysis of the band corresponding to the activity detected in electrophoresis pattern of crude supernatant from a p69TVLB-derived transformant (clone 4) of *Yarrowia lipolytica* expressing laccase IIIb.

APPLICATION: Towards the engineering of fungal laccases

1- T1 site and neighboring 2,5-xylydine ligand

2 - Sequence alignment of laccases from fungi and higher plants



- Two charged or polar residue interact with the amino group of the xylydine:
- His (H458) coordinates the copper that functions as the primary electron acceptor
 - Asp (D206) possesses acido-basic properties that could be responsible for the pH-dependent activity towards substrates

3- Strategy of directed mutagenesis

- Three mutants have been obtained comprising:
- a residue with another negative chain: D206E
 - a residue with a polar uncharged chain: D206N
 - a residue with an uncharged chain: D206A

-Their catalytic properties are currently studied.

- Biological roles of laccases from fungi and higher plants appear to be opposite: fungal laccases are involved in ligninolysis whereas enzymes from plants catalyze the initial polymerization of monolignols during lignification

-Fungal laccases have a negatively charged residue at the 206 position (Asp, D or Glu, E), whilst those of higher plants have an asparagine (Asn, N)

Fungi: *T.versicolor* = *Trametes versicolor*; *C. cinereus* = *Coprinus cinereus*; *P. cinnabarinus* = *Pycnoporus cinnabarinus*;
M. albomyces = *Melanocarpus albomyces*
Higher plants: *A. thaliana* = *Arabidopsis thaliana*; *R. vernificera* = *Rhus vernificera*; *A. pseudoplatus* = *Acer pseudoplatus*

References
[1] C. Jolivald, S. Brenon, E. Caminade, C. Mougin, M. Ponté (2000) Immobilization of laccase from *Trametes versicolor* on a modified PVDF microfiltration membrane: characterization of the grafted support and application in removing a phenylurea pesticide in wastewater. *J. Membrane Sci.* 180:103-113.
[2] C. Mougin, F.-D. Boyer, E. Caminade and R. Rama (2000). Cleavage of the diketonitrile derivative of the herbicide isoxallatolite by extracellular fungal oxidases. *J. Agric. Food Chem.* 48:4529-4534.
[3] C. Mougin, C. Jolivald, P. Briozzo, C. Madzak (2003) Fungal laccases: from structure-activity studies to environmental applications. *Environ. Chem. Lett.* 1:145-148.
[4] T. Bertrand, C. Jolivald, P. Briozzo, E. Caminade, N. Joly, C. Madzak and C. Mougin (2002) Crystal structure of a four-copper laccase complexed with an arylamine: insights into substrate recognition and correlation with kinetics. *Biochemistry* 41:7325-7333