



HAL
open science

High Speed Quantitative 3D Blood Flow Imaging by Multi illumination Holographic Microscopy

Alexey Brodoline, Nitin Rawat, Dario Donnarumma, Daniel Alexandre,
Michel Gross

► **To cite this version:**

Alexey Brodoline, Nitin Rawat, Dario Donnarumma, Daniel Alexandre, Michel Gross. High Speed Quantitative 3D Blood Flow Imaging by Multi illumination Holographic Microscopy. Journées scientifiques du LabEx NUMEV, Nov 2018, Montpellier, France. hal-01927083


HAL Id: hal-01927083

<https://hal.science/hal-01927083>

Submitted on 21 Nov 2018

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



High-Speed Quantitative 3D Blood Flow Imaging by Multi illumination Holographic Microscopy

Alexey Brodoline, Nitin Rawat, Dario Donnarumma, Daniel Alexandre, Michel Gross
 Laboratoire Charles Coulomb (L2C), UMR5221, CNRS-Université de Montpellier, Montpellier, France

Abstract:

A holographic microscopy method using two illumination beams with a 0.5 aperture objective is proposed to image blood microcirculation in zebrafish larvae. Recent achievements in 3D imaging of blood flow are presented.

OCIS codes: (090.1995) Digital holography; (170.0180) Microscopy; (170.1470) Blood or tissue constituent monitoring; (290.5850) Scattering, particles.

Background

Zebrafish
Danio rerio



Model organism:
 Developmental biology, genetics, neuroscience, biomedicine and behavioral studies.


Blood flow monitoring:

- Assess angiogenesis, tumor vascularization
- Early detection of cardiovascular diseases, age-related macular degeneration

Zebrafish used for:

- A model in drug screening
- Effect of the active substances on the cardiovascular functions
- Zebrafish heart as a model for human heart

Conventional imaging techniques




The techniques commonly used are invasive

Goals

Expected characteristics of the measurement technique

- Fast 3D Imaging**
- Position of RBCs in volume
 - 3D information in one camera frame
- Non-invasive operation**
- No contrast agents (label-free)
 - No genetic engineering
- Multimodality, information extraction:**
- Integration with conventional microscopes
 - Amplitude / phase information
 - Quantitative measurement of blood flow


Dual-illumination Holographic Microscopy Setup



The setup is an upright transmission microscope adapted for digital holography in off-axis configuration. Typical camera frequencies are 100/200 Hz.


A double illumination of the sample is performed to better localize the position of the red blood cells (RBCs) in z.

Benefit of double illumination illustrated with an element of a USAF target. Reconstruction in different planes.




water + agar + tricaine → cover slip → 2 to 5 days zebrafish embryos


Contrast enhancement by difference of images:



2D Holographic Reconstruction of Zebrafish Blood Flow




Fourier Transform of the hologram obtained with a 20x/0.5 microscope objective with a phase ring. The double illumination is visible inside the pupil of the MO.



+1 order cropped in the Fourier space. The two directions of illumination are numerically selected and represented in colors.


TF -1
 TF -1
 TF :
 Fourier Transform




Reconstruction of the field in the plane of the sample. To visualize only the moving red blood cells (RBCs) a difference of images is realized.

$$H(x, y) = \sum_{k=0}^{N-1} I_k(x, y) e^{i \frac{2\pi k}{N}}$$


Reconstructions in different planes on either sides of the sample. Averaged over 128 frames




3D Reconstruction with cleaning algorithm



3D cube of data containing the information on the light scattered around the sample. A cleaning algorithm extracts the positions of the RBCs.



The positions of the RBCs correspond to the points of maximum correlation between the two illuminations.



Reconstructed 3D images of the vascular system of a 5-days zebrafish. (a, b) Averaged in time positions of the RBCs give the shape of the vessels. (c) Positions of the RBCs in one camera frame.