

Complete Genome Sequences of 11 Type Species from the Thermococcus Genus of Hyperthermophilic and Piezophilic Archaea

Phil M. Oger

▶ To cite this version:

Phil M. Oger. Complete Genome Sequences of 11 Type Species from the Thermococcus Genus of Hyperthermophilic and Piezophilic Archaea. Genome Announcements, 2018, 6 (8). hal-02001367

HAL Id: hal-02001367 https://hal.science/hal-02001367v1

Submitted on 22 Nov 2023

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

Complete Genome Sequence of 11 type-species from the *Thermococcus* genus of hyperthermophilic and piezophilic Archaea

Philippe M. Oger

Université de Lyon, INSA Lyon, CNRS UMR5240, 11 Avenue Jean Capelle, F-69621 Villeurbanne, France

We report the genome sequences of the type strains of the species *Thermococcus barossii*, *T. celer*, *T. chitonophagus*, *T. gorgonarius*, *T. pacificus*, *T. peptonophilus*, *T. profondus*, *T. radiotolerans*, *T. siculi*,
 T. thioreducens and the prototype of a possible type strain of a novel *Thermococcus* species, strain
 P6.

5

Thermococcus celer were isolated from surface hydrothermal vent systems in the Mediterranean (1). *Thermococcus* sp. Strain P6 was isolated in 1997 by Pr. M. Thomm (University of Kiel, Germany) in 1997
from an hydrothermal site off Palaeochori Bay (Milos, Aegean Sea, Greece). *T. thioreducens* was isolated
from the Rainbow hydrothermal vent system on the Mid-Atlantic ridge, 2600m below the ocean surface
(2).All other species were isolated from hydrothermal vents systems in the Pacific ocean (3-9). All species
exhibit very similar growth characteristics, with optimal conditions for growth ranging from 83 and 88°C
(T), pH 5.8 and pH 7 and salinity between 2 and 3%.

13 DNA was extracted from cells grown in TRM medium to late exponential phase. Whole-genome shotgun sequencing was carried out by a combination of IonTorrent PGM (800 bp fragment library, 318 14 Chip, HiQ chemistry) and Illumina HiSeq run (500bp library, 100bp paired reads), which generated at 15 least 583k reads with a mean 389 bp length (IonTorrent) and 13 million paired 100 bp reads after QC 16 check and trimming (10-11) for each species. The IonTorrent data were assembled de novo with Newbler 17 2.8 (12). The paired-end Illumina data was assembled de novo with velvet (13). The pooled IonTorrent 18 19 and Illumina data was assembled de novo with MIRA (14) and Newbler. The hybrid data assembly yielded betwen 1 and 6 contigs which were connected by PCR. 20

21 Chromosome sizes range from 1.52Mb in Thermococcus sp. strain P6 to 2.07Mb in Thermococcus thioreducens. which is similar to known genomes from this family. Only two species, T. peptonophilus 22 and T. profundus, harbor an extrachromosomal element. The T. profundus plasmid belongs to the 23 24 pTBMP1 family of plasmid described in the piezophilic archaeon T. barophilus (15). To the exception of T. chitonophagus, the GC% ranges between 51.7 and 56% Similar to previously sequenced 25 Thermococcales, all genomes harbor a single copy of the 16S-23S operon, two copies of 5S operons and 26 27 45 to 48 tRNA genes, large CRISPR loci associated with cas genes (cas1, cas2, cas3, cas4, cas5, cas6, 28 *csm*²), indicating multiple contacts with mobile elements of viral or plasmidic origin. Sinteny analyzes confirm the plasticity of *Thermococcales* genomes. The gene content is consistent with an anaerobic and 29 heterotrophic lifestyle. The genome of T. barossi encodes the capacity for hydrogenogenic growth on 30 carbon monoxide, as identified previously in T. onnurienus and T. barophilus. 31

The genome sequence of *T. peptonophilus* (GC%, core genome phylogeny) confirms previous reports that it may have been misplaced into the *Thermococcus* genus, and that it should be placed in the *Pyrococcus* genus instead. Strain P6 is most closely related to *Thermococcus celer* strain Vu 13 (1), as indicated by *in silico* DNA-DNA hybridization of their genomes (predicted value = 28.90% using the genome-to-genome distance calculator GGDC2.0 (16). This genomic distance suggests that strain P6 defines a novel species inside the *Thermococcus* genus.

38 Nucleotide sequence accession number. These whole Genome Shotgun project has been deposited in
39 GenBank under the accession numbers CP014750, CP014751, CP014854, CP014855, CP014862,
40 CP014863, CP015101, CP015102, CP015103, CP015104, CP015105, CP015106, and CP015193.

41

42 ACKNOWLEDGEMENTS

43 This work was supported in part by the Agence Nationale de la Recherche (ANR-10-BLAN-44 Living deep) and a grant from the Fond Recherche of the ENS Lyon to PO. Ion Torrent 45 sequencing was performed in part at the sequencing platform of the Institut de Génomique46 Fonctionelle de Lyon by Dr. Sandrine Hughes and Benjamin Gillet.

47

48 REFERENCES

Zillig W, Holz I, Janekovic D, Schäfer W, Reiter WD. 1983. The archaebacterium *Thermococcus celer* represents, a novel genus within the thermophilic branch of the archaebacteria. Syst. Appl.
 Microbiol. 4:88-94.

Pikuta EV, Marsic D, Itoh T, Bej AK, Tang J, Whitman WB, Ng JD, Garriott OK, Hoover RB.
 2007. *Thermococcus thioreducens* sp. nov., a novel hyperthermophilic, obligately sulfur-reducing
 archaeon from a deep-sea hydrothermal vent. Int. J. Syst. Evol. Microbiol. 57:1612-8.

Grote R, Li L, Tamaoka J, Kato C, Horikoshi K, Antranikian G. 1999. *Thermococcus siculi* sp.
nov., a novel hyperthermophilic archaeon isolated from a deep-sea hydrothermal vent at the Mid-Okinawa
Trough. Extremophiles 3:55-62.

Huber R, Stöhr J, Hohenhaus S, Rachel R, Burggraf S, Jannasch HW, Stetter KO. 1995.
 Thermococcus chitonophagus sp. nov., a novel, chitin-degrading, hyperthermophilic archaeum from a
 deep-sea hydrothermal vent environment. Arch. Microbiol. 164:255-264.

Jolivet E, Corre E, L'Haridon S, Forterre P, Prieur D. 2004. *Thermococcus marinus* sp. nov. and *Thermococcus radiotolerans* sp. nov., two hyperthermophilic archaea from deep-sea hydrothermal vents
that resist ionizing radiation. Extremophiles 8:219-27.

64 6. Kobayashi T, Kwak YS, Akiba T, Kudo T, Horikoshi K. 1994. *Thermococcus profundus* sp. nov.,
65 A new hyperthermophilic archaeon isolated from a deep-sea hydrothermal vent. Syst. Appl. Microbiol.
66 17:232-236.

67 7. Miroshnichenko ML, Gongadze GM, Rainey FA, Kostyukova AS, Lysenko AM, Chernych NA,

68 Bonch-Osmolovskaya EA. 1998. *Thermococcus gorgonarius* sp. nov. and *Thermococcus pacificus* sp.

69 nov.: heterotrophic extremely thermophilic archaea from New Zealand submarine hot vents. Int. J. Syst.70 Bacteriol. 48:23-29.

Duffaud GD, d'Hennezel OB, Peek AS, Reysenbach A-L, Kelly RM. 1998. Isolation and
 characterization of *Thermococcus barossii*, sp. nov., a hyperthermophilic archaeon isolated from a
 hydrothermal vent flange formation. Syst. Appl. Microbiol. 21:40-49.

9. Gonzalez JM, Kato C, Horikoshi K. 1995. *Thermococcus peptonophilus* sp. nov., a fast-growing,
extremely thermophilic archaebacterium isolated from deep-sea hydrothermal vents. Arch. Microbiol.
164:159-164.

77 10. FastQC, A quality control tool for high throughput sequence data.

78 (http://www.bioinformatics.babraham.ac.uk/projects/)

79 11. Bolger AM, Lohse M, Usadel B. 2014. Trimmomatic: a flexible trimmer for Illumina sequence
80 data. Bioinformatics 30 (15):2114-2120.

81 12. Newbler, a software package for de novo DNA sequence assembly

82 (http://www.454.com/products/analysis-software/)

Bruijn graphs. Genome Res. 18:821-829.

85 14. Chevreux B., Wetter, T. and Suhai, S. Genome sequence assembly using trace signals and
86 additional sequence information. Computer Science and Biology: Proceedings of the German Conference
87 on Bioinformatics (GCB) 1999. 99, pp. 45-56.

88 15. Vannier P, Marteinsson V, Fridjonsson O, Oger P, Jebbar M. 2011. Complete genome sequence
89 of the hyperthermophilic piezophilic, heterotrophic and carboxydotrophic archaeon *Thermococcus*90 *barophilus* MP. J. Bacteriol. 193:1481–1482.

91 16. Meier-Kolthoff JP, Klenk H-P, Göker M. 2014. Taxonomic use of DNA G+C content and DNA–
92 DNA hybridization in the genomic age. Int. J. Syst. Evol. Microbiol. 64:352–356.