

Protein interactions in stirred yoghurt

M Renan, Fanny Guyomarc'H, Véronique Arnoult-Delest, D Paquet, Gérard Brulé, Marie-Hélène Famelart

▶ To cite this version:

M Renan, Fanny Guyomarc'H, Véronique Arnoult-Delest, D Paquet, Gérard Brulé, et al.. Protein interactions in stirred yoghurt. Food colloids 2008: creating structure, delivering functionality, Apr 2008, Le Mans, France. 2008. hal-02215957

HAL Id: hal-02215957

https://hal.science/hal-02215957

Submitted on 31 Jul 2019

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

Protein interactions in stirred yoghurt

Renan M.a, Guyomarc'h F.a, Arnoult-Delest V.b, Paquet D.b, Brulé G.a, Famelart M.H.a

^a Inra, Agrocampus Rennes, UMR 1253 Science et Technologie du Lait et de l'Œuf, 65 rue de St-Brieuc, 35042 Rennes Cedex, France

^b Danone Research, RD 128, 91767 Palaiseau Cedex, France

(*) Tel. +33.(0)2.23.48.53.43, Fax. +33.(0)2.23.48.53.50, @: marie-helene.famelart@rennes.inra.fr

Introduction

Stirred yoghurt is obtained after pumping, cooling and packaging of set yoghurt at pH ≈ 4.6. The viscoelastic properties of the gel decrease on shearing, but reincrease on storage at 4°C. This phenomenon is known as "rebodying". The nature of bonds that are responsible for this rebodying was investigated.

Materials & Methods

Reconstituted heated milk (90°C-10 min) was inoculated with non ropy yoghurt cultures. The gel was stirred through a mesh of 350-400 µm-holes and then in a home food processor (Renan et al., 2008), and characterised by low amplitude dynamic oscillation for 24 h (short-term rebodying) and viscosity (long-term rebodying) during 28 d of storage at 4°C

- ① Set yoghurts were stirred at pH 4.4 to 5.0 and stored at 4 20°C for 24 h, with or without NaN3 (using a Box-Benken design);
- ② A constant set yoghurt was stirred in different physicochemical conditions applied on stirring

Results & Discussion





Fig.1: pH and G' of set acid gels


Set yoghurt stirred at lower pH values had a higher G' value (Fig.1). This led to a higher short rebodying (Fig.2) → electrostatic interactions involved.

No effect of the storage temperature of the stirred yoghurt during 24 h on short-

term rebodying without over-acidification (Fig.2) \rightarrow hydrophobic interactions were not involved in short-term rebodying.

Over-acidification took part in rebodying (Fig.2).





pH Fig.2: Short-term rebodying of stirred yoghurt as a function of the pH on stirring and the storage temperature of stirred yoghurt during 24 h (as the increase in G' during 12 h




Fig.3. G' of stirred gels just after stirring and 12 h after stirring as a function of the pH in stirred yoghurt just after stirring and 24 h after stirring (the pH hardly decreased between 12 and 24 h)

With NaN₃ addition, the pH of the stirred gel increased by 0.2 unit on stirring.

Changes in G' during short-term rebodying very complex were (Fig.3), depending on the pH on stirring but also on its evolution during storage.

CaCl₂ (CA+) or Na-Citrate (CA-) addition on stirring

CaCl₂ addition on stirring reduced the short-term rebodying (Fig.4). same ionic strength increase with NaCl had no effect.

At 20 h:

 $tan\delta_{CA^{-}} > tan\delta_{ctrl} > tan\delta_{CA^{+}}$

- CA+ gel
- More solid-- control gel
- like gels - CA⁻ gel
- Very few significant effects on long-term rebodying:
- + NaCl → No long-term rebodying
- + Na-citrate stirred gel had the highest viscosity increase in 28 d (Fig.4).




Fig.4: Short-term rebodying: + CaCl₂ (+67 mmol.kg⁻¹) + Na-citrate 100 mmol.kg⁻¹




Fig.5: Long-term rebodying: +CaCl. (+67mmol.kg-1)

Estimations of calcium partition (Mekmene et al., 2008) give a higher content in Pserine and carboxylic acid bound calcium with CaCl2 addition and a lower one with Na-citrate addition.

Conclusions

Concerning the rebodying process, electrostatic interactions were suggested as seen from the effect of the pH on stirring on short term rebodying and the fact that NaCl addition prevented the long-term rebodying. Probably, hydrophobic interactions were of minor importance. Rebodying was partially explained by over-acidification. Calcium addition plays a negative specific role on the rebodying of stirred yoghurt. Calcium bridges between proteins likely formed, that could prevent protein rearrangements in the stirred gel.

Food Colloids 2008-Creating Structure, Delivering Functionality - Le Mans, 6-9 April









