



HAL
open science

Physico-chemical changes during acidification by a yoghurt starter of heat-treated milk

Marie-Hélène Famelart, A Laligant, D Paquet, Y Lemèner, Gérard Brulé

► **To cite this version:**

Marie-Hélène Famelart, A Laligant, D Paquet, Y Lemèner, Gérard Brulé. Physico-chemical changes during acidification by a yoghurt starter of heat-treated milk. 26^{ème} congrès international de laiterie, Sep 2002, Paris, France. 2002. hal-02244976


HAL Id: hal-02244976

<https://hal.science/hal-02244976>

Submitted on 1 Aug 2019

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



INRA

Institut National de la Recherche Agronomique

Physico-chemical changes during acidification by a yoghurt starter of heat-treated milk

Famelart M.H.^a, Laligant A.^a, Paquet D.^b, Y. Leménager^c, Brulé G.^d

^aINRA-LRTL; 65, rue de Saint Brieuc, 35042 Rennes cedex, France

^bDanone Vitapole Recherche, 5 Av. Galilée, 92350 Le Plessis-Robinson, France

^cRhodia Food, Z.I. de Buxières, 86220 Dangé Saint Romain, France

^dENSAR, 65, rue de St Brieuc, 35042 Rennes cedex, France.

Understanding the built up of acid gels by changing the fermentation temperature, as a low gelation rate can result in a higher firmness.

Gelation was followed on heat-treated milk

✓ by static measurements (at a final equilibrium)

✓ by dynamic measurements (in a continuous on-line mode)

Methods and materials

Reconstituted skim milk at 11% solid content (Ingredia, Arras, France)

Heat treatment at 90°C for 10 min

Acidification with *Streptococcus thermophilus* and *Lactobacillus delbrueckii subsp. bulgaricus* at 42 and 30°C or, when specified, with glucono-delta-lactone (GDL) at 30°C for 17 h

Static methods

pH drop stopped by addition of 2.05 g.kg⁻¹ sodium azide and stirring
Fractionation of milk into a permeate and a supernatant at 42 or 30°C
Permeate: on Centriflo membrane (CF25)
Supernatant: at 150000 g for 70 min


Soluble calcium (Ca) on permeate by atomic absorption spectroscopy [1]
Soluble casein by RP-HPLC [2]

Dynamic methods

Sinusoidal oscillation with coaxial cylinder at 2 % strain at 1 Hz (4 experiments/ temperature)
H-NMR relaxation in a low field spectrometer (Bruker, 20MHz, 0.47 T) fitted in a 1-exponential (42°C) or a 2-exponential decay (30°C) (2 experiments/ temperature)

Front-face fluorescence (FFF) of intrinsic tryptophane (TRP) at $\lambda_{exc}=290$ nm and principal component analysis (PCA) (4 experiments/ temperature)

Results and discussion







• Fermentation \Rightarrow shift of soluble Ca toward low pH values

• Ca solubilisation: 30°C = 42°C

• Heat treatment \Rightarrow shift of soluble casein toward higher pH values

• Casein solubilisation: 30°C > 42°C



¹H-RMN

- Fermentation \Rightarrow shift toward acid pH related to the shift of soluble Ca
- Sensitive to the mode of acidification
- Heterogeneity of protons : 42°C, pH 5, mono-exponential \Rightarrow a short + a long T2; 30°C, pH 4.5 \Rightarrow a break in the T2 curve
- Higher magnitude of T2 at 30°C
- Higher rate of T2 increase at 42°C

FFF + PCA : information about the environment of TRP residues during acidification

- PC1 was sensitive to fermentation time and PC2 to fermentation temperature
- 3 steps:
 - from initial pH values to the gel point : TRP environment became more hydrophobic
 - from the gel point to pH 4.5 (30°C) or 4.7-4.9 (42°C): TRP environment became less hydrophobic
 - final firmness increase: TRP environment more hydrophobic

Conclusion

- Shift in Ca and T2 with dynamic measurements due to the limited diffusion of protons
- Proteins were not in the same state during fermentation than after long time equilibrium as for acidification with GDL

[1] Le Graet et al., 1993
[3] Singh et al., 1996

[2] Laligant et al., 2002
[4] Law, 1996
[5] Mariette et al., 1998
[6] Famelart et al., 1997