

Self-assemblies of azacitidine prodrug: An innovative therapy against myelodysplastic syndromes

Milad Baroud, Elise Lepeltier, Sylvain Thépot, Yolla El Makhour, Olivier

Duval

► To cite this version:

Milad Baroud, Elise Lepeltier, Sylvain Thépot, Yolla El Makhour, Olivier Duval. Self-assemblies of azacitidine prodrug: An innovative therapy against myelodysplastic syndromes. SFNano -2019, Dec 2019, Dijon, France. hal-02415685

HAL Id: hal-02415685 https://hal.science/hal-02415685

Submitted on 17 Dec 2019

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.





Institut national de la santé et de la recherche médicale







Self-assemblies of azacitidine prodrug: An innovative therapy against myelodysplastic syndromes


<u>Baroud M¹</u>, Lepeltier E^{1*}, Thepot S², El-Makhoor Y³, Duval O¹ ¹Translational Micro and Nanomedecines (MINT, INSERM 1066 / CNRS 6021), Université d'Angers, Angers, France ²University Hospital of Angers, Hematology, 49933 Angers, France ³Laboratory of Cancer Biology and Molecular Immunology, Lebanese University, Faculty of Sciences, Hadath-Beirut, Lebanon

*corresponding author : elise.lepeltier@univ-angers.fr

ntroduction

5-azacitidine, a cytidine analogue and a hypomethylating agent, is one of the main drugs being used for the treatment of myelodysplastic syndromes ^[1]. However after administration, it exhibits several limitations including restricted diffusion and cellular internalization due to its hydrophilicity, and rapid enzymatic degradation. Our objective is to improve the drug diffusion and to protect it from metabolic degradation via the formulation of an amphiphilic prodrug and its self-assembly into a nanoparticle. Following the conjugation of the fatty acid to the azacitidine, the obtained prodrug was solubilized in acetone and mixed with water at different ratios to obtain self-assemblies by nanoprecipitation, thus protecting the active molecule from enzymatic degradation. This prodrug should be cleaved by cathepsin B^[2], an enzyme overexpressed in cancerous cells^[3], thus increasing the specificity of the drug.

Prodrug synthesis




The multi-step reaction for prodrug synthesis: starting with the protection of the alcohol groups with TBDMSCI to avoid secondary conjugations during the next step of conjugation to the chosen fatty acid. Finally, the obtained product is deprotected using TBAF to give the prodrug (azacitidine-fatty acid conjugate).

To determine size and

All of the products were identified with either ¹H NMR spectroscopy, mass spectrometry or both.

Self-assembly and characterization



Further Experiments :

2.2


Determining size, shape and zeta potential:


- To determine size and Critical Aggregation Concentration (CAC) zeta potential
- Dynamic Light Scattering (DLS)
- Nanoparticle Tracking Analysis (NTA) concentration
- Transmission electron cryo-microscopy To determine shape (cryo-TEM)

Next step: in vitro studies

The in vitro studies will be performed on several acute myeloid leukemia cell lines with varying sensitivity to azacitidine:


- HI-60
 - HEL
- Molm 13 • SKM1
- Azacitidine resistant SKM1






Studies to be performed:

Nanoprecipitation: the prodrug is solubilized in an organic solvent (acetone) then mixed with water at different ratios, leading to the formation of the self assemblies. Following that, the solvent is evaporated to obtain the self-assembled nanoparticles suspended in water. The formation of the self-assemblies is characterized by the observation of an opalescence in the solution.



- Cell viability studies (MTT, LDH)
- Cellular uptake studies
- Cell cycle study (FACS)
- Methylation quantification (LINE -1, whole DNA methylome analysis)



The Global DNA Methylation – LINE-1

is an ELISA-based assay that is specially formulated for the detection of methylation in Long Interspersed Nucleotide Element 1 (LINE-1) repeat elements of human genomic DNA. LINE-1 methylation serves as a surrogate for global DNA methylation levels and can be utilized to compare relative changes in 5methylcytosine (5-mC) levels across Add developing and [°] different sample types stop solution

Conclusion

An azacitidine prodrug was achieved following a multi-step reaction to conjugate azacitidine to a fatty acid. The obtained prodrug will undergo a process of self-assembly to produce nanoparticles, that will be fully characterized. Finally, in vitro studies will be performed to determine the biological efficacy of these nanoparticles.

References:

[1] Issa, J.-P. J. The myelodysplastic syndrome as a prototypical epigenetic disease. *Blood* 121, 3811–3817 (2013)

[2] Dheer, D., Nicolas, J. & Shankar, R. Cathepsin-sensitive nanoscale drug delivery systems for cancer therapy and other diseases. Advanced Drug Delivery Reviews (2019)

[3] Alatrash, G. et al. Cathepsin G is broadly expressed in acute myeloid leukemia and is an effective immunotherapeutic target. Leukemia 31, 234 (2016)