

NMR metabolomics for understanding the role of alkaline phosphatase in neurotransmission and inflammation processes

Stéphane Balayssac, Thomas Cruz, Marie Gleizes, Gary Lutin, Juana Como, Anne Briolay, Lionel G Nowak, Myriam Malet-Martino, David Magne, Caroline Fonta, et al.

▶ To cite this version:

Stéphane Balayssac, Thomas Cruz, Marie Gleizes, Gary Lutin, Juana Como, et al.. NMR metabolomics for understanding the role of alkaline phosphatase in neurotransmission and inflammation processes. European RFMF Metabomeeting, Jan 2020, Toulouse, France. hal-02875612

HAL Id: hal-02875612

https://hal.science/hal-02875612

Submitted on 19 Aug 2022

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.















NMR METABOLOMICS FOR UNDERSTANDING THE ROLE OF ALKALINE PHOSPHATASE IN NEUROTRANSMISSION AND INFLAMMATION PROCESSES

Stéphane BALAYSSAC¹, Thomas CRUZ¹, Marie GLEIZES², Gary LUTIN¹, Juana COMO³, Anne BRIOLAY³, Lionel G NOWAK², Myriam MALET-MARTINO¹, David MAGNE³, Caroline FONTA², Véronique GILARD¹

¹Groupe de RMN Biomédicale, Laboratoire SPCMIB, Université Paul Sabatier, Toulouse, France

²Centre de recherche Cerveau et Cognition (CerCo), CNRS-Université de Toulouse, France

³Institut de Chimie et de Biochimie Moléculaires et Supramoléculaires (ICBMS), CNRS-Université Lyon 1, Villeurbanne, France

Introduction


There are several genes for alkaline phosphatase (AP). One of them codes for TNAP (Tissue Non-specific Alkaline Phosphatase). It is a ubiquitous ectoenzyme located in different tissues in the body.

Mutations of the TNAP gene are responsible for hypophosphatasia, a rare disease characterized by defective bone mineralization caused by pyrophosphate accumulation. In the severe form, hypophosphatasia is associated with epileptic seizures. This neurological disorder could be explained by a decrease in GABA level due to lack of pyridoxal 5'phosphate (PLP) dephosphorylation.

In this study, we aimed in identifying metabolic pathways associated with TNAP activity in brain and in exploring its role in the inflammatory process. We used untargeted ¹H NMR-based metabolomics on tissue extracts from two original models: (i) Akp2 mice, which is the rodent model of hypophosphatasia by invalidation of the TNAP gene, and (ii) a model of postprandial inflammation based on food intake.

Experimental part

EXTRACTION PROCEDURE



NMR


The dried aqueous extract was redissolved in deuterated borate buffer at pH 10 and a solution of sodium 2,2,3,3-tetradeutero-3trimethylsilylpropionate (TSP) was added as an internal chemical shift and quantification reference.

¹H NMR experiments were recorded at 298K on a Bruker Avance 400 spectrometer equipped with a 5 mm BBI probe. A simple 1D NMR sequence (pulse-acquisition-relaxation delay) including water presaturation was used. Acquisition parameters were 3.7 s and 2.0 s for acquisition time and relaxation delay respectively, with a spectral width of 12 ppm, a flip angle of 30° and 768 scans.

STATISTICAL ANALYSES

The bin area method based on intelligent bucketing (KnowItAll®, BIO-RAD) was used for segmenting NMR spectra. After normalization, univariate (boxplot, p-value) and multivariate analyses (PLS-DA (Partial Least Squares-Discriminant Analysis)) were applied to NMR data using the free R and SIMCA P+ (Umetrics®) softwares. Each statistical model was validated thanks to the Q² (predictive value of the model) and CV-ANOVA (ANalysis Of VAriance of Cross-Validated predictive residuals) values.

Typical ¹H NMR spectrum




¹H NMR spectrum 4.4- 2. 9 ppm region of an aqueous extract of the brain from an Akp2 (+/+) mouse

Akp2 mice model


Metabolic changes

- 39 metabolites identified in whole brain extract.
- 8 metabolites significantly impacted by TNAP gene invalidation (cystathionine, adenosine, GABA, methionine, histidine, 3methylhistidine, N-acetylaspartate and N-acetyl-aspartylglutamate).
- Akp2 (+/-) / (+/+) mice have very **similar metabolomes**.
- Very large changes in Knockout mice for GABA, cystathionine and adenosine concentration compared to Akp2 (+/+) and Akp2 (+/-).
- GABA and adenosine levels were decreased in Akp2 (-/-) mice compared to others.
- Cystathionine levels were increased in Akp2 (-/-) mice compared to others.




Comparison of ¹H NMR spectra (zoom around 2.9 ppm) of aqueous brain extracts from Akp2+/+ (blue), Akp2+/- (green), and Akp2-/- (red) mice.

Cluster analysis revealed separation of Akp2 -/- vs Akp2 +/- and +/+ mice




Unsupervised statistical analysis (PCA) indicated that the metabolome of Akp2-/- mice was distinct from that of Akp2 +/- and +/+ mice.


3 metabolites strongly impacted by TNAP gene invalidation



Concentration decreased by 40 %




Concentration increased by 450 %



Concentration decreased by 80%

Brain extract analysis



- discriminating main metabolites were found.
- Levels of **methionine**, lactate and succinate increased in the postprandial model.

Postprandial model (Preliminary results)

PLSDA models were validated for brain ($R^2Y =$ 0.89; $Q^2 = 0.65$) and liver extract ($R^2Y = 0.49$; $Q^2 = 0.89$).

Metabolites selected according were coefficient plots (confidence intervals derived from jack-knifing ($\alpha = 0.95$) < coefficient values), VIP values (>1.5), and p-values of Wilcoxon test.

Liver extract analysis

- 10 main discriminating metabolites were identified: glucose, lactate, GSH, alanine, betaine, GPC, PC, choline, hypotaurine and succinate.
- Levels of identified discriminating metabolites except choline increase in the postprandial model.

Conclusion

- Pyridoxal phosphate (PLP) dephosphorylation by TNAP is required for cystathionine and GABA production; moreover TNAP ectonucleotidase function produces adenosine.
- → The metabolomic approach showed that TNAP deficiency reduces the level of GABA and adenosine and increases cystathionine level in the brain.
- In parallel with its regulatory role in neuronal synaptic transmission, TNAP could be involved in inflammation.
- \rightarrow The metabolic changes observed in this study suggest that TNAP is directly involved in the production of anti-inflammatory agents (adenosine) or indirectly from precursor agents (methionine and cystathionine, precursor of H_2S).

