

Cnidarian cell cryopreservation: A powerful tool for cultivating and functional assays

Clara Fricano, Eric Röttinger, Paola Furla, Stephanie Barnay-Verdier

▶ To cite this version:

Clara Fricano, Eric Röttinger, Paola Furla, Stephanie Barnay-Verdier. Cnidarian cell cryopreservation: A powerful tool for cultivating and functional assays. Cells, In press. hal-03021301

HAL Id: hal-03021301 https://hal.science/hal-03021301

Submitted on 25 Nov 2020

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.





1 Article

Cnidarian cell cryopreservation: A powerful tool for cultivation and functional assays

4 Clara Fricano¹, Eric Röttinger¹, Paola Furla^{1,\$} and Stéphanie Barnay-Verdier ^{1,2,\$,*}

- ¹ Université Co[^]te d'Azur, CNRS, INSERM, Institute for Research on Cancer and Aging (IRCAN), 28 avenue
 de Valombrose, F-06107 Nice, France
- 7 ² Sorbonne Université, UFR 927, 4 place Jussieu, F-75252 Paris, France
- 8 * Correspondence: Corresponding author: <u>stephanie.barnay-verdier@upmc.fr</u>; Tel.: +33 4 93 37 77 39
- 9 ^{\$} These authors share the seniorship position
- 10 16-digit ORCID:
- 11 Clara Fricano: 0000-0001-8496-5915
- 12 Eric Röttinger: 0000-0002-2938-6774
- 13 Paola Furla: 0000-0001-9899-942X
- 14 Stéphanie Barnay-Verdier: 0000-0002-4615-3059

15 Running Title: Cryopreservation of Cnidarian Cell Cultures

16 Abstract: Cnidarian primary cell cultures have a strong potential to become a universal tool to 17 assess stress-response mechanisms at the cellular level. However, primary cell cultures are 18 time-consuming regarding their establishment and maintenance. Cryopreservation is a commonly 19 used approach to provide stable cell stocks for experiments, but it is yet to be established for 20 Cnidarian cell cultures. The aim of this study was therefore to design a cryopreservation protocol 21 for primary cell cultures of the Cnidarian Anemonia viridis, using dimethyl sulfoxide (DMSO) as a 22 cryoprotectant, enriched or not with foetal bovine serum (FBS). We determined that DMSO 5% 23 with 25% FBS was an efficient cryosolution, resulting in 70% of post-thaw cell survival. The success 24 of this protocol was first confirmed by a constant post-thaw survival independently of the cell 25 culture age (up to 45 days old) and the storage period (up to 87 days). Finally, cryopreserved cells 26 displayed a long-term recovery with a maintenance of the primary cell culture parameters and 27 cellular functions: formation of cell aggregates, high viability and constant cell growth, and 28 unchanged intrinsic resistance to hyperthermal stress. These results will further bring new 29 opportunities for the scientific community interested in molecular, cellular and biochemical aspects 30 of cnidarian biology.

Keywords: primary cell culture; sea anemone; *Anemonia viridis*; DMSO; marine invertebrate;
 post-thaw recovery

33 1. Introduction

34 In vitro cell cultures are important tools for research in many fields, including development, 35 virology, cancer research, toxicity testing, biotechnology, biomedicine as well as for environmental 36 research [1-3]. Mammalian cells lines are well established and commonly used since decades, 37 followed by other vertebrates (e.g. zebrafish; for review Vallone et al. 2007 [4]) and insect cell lines 38 (for review Lynn, 2001 [5]). Despite much effort devoted since the 1970s (see for reviews Rinkevich, 39 2005 [6], 2011 [7] and Cai and Zhang, 2014 [8]) marine invertebrate cell cultures are not as advanced. 40 While marine invertebrate cell lines (i.e. permanently established cell cultures) are yet to be 41 available, recent reports on the establishment of primary cell cultures are encouraging with 42 maintenance and/or growing from a couple of days to weeks [9–14].

43 Once protocols for reproducibly initiating primary cell cultures are established, the next 44 important obstacle to overcome is the development of preservation procedures in order to outreach 45 primary cell cultures limitations, notably their limited lifespan. Indeed, such a preservation tool will 46 reduce the frequency of primary cell culture establishment and will minimize wastage of a valuable 47 resource at reseedings by creating cell stocks for as long as a primary cell culture is healthy. Among 48 preservation procedures, cryopreservation is considered to be the optimal long-term storage method 49 for maintaining a variety of biological materials, including cell cultures, in a state of metabolic arrest 49 for maintaining a variety of biological materials.

for considerable periods of time [16].
To date, for marine invertebrates, spermatozoa, oocytes, embryos and different larval stages

52 have been successfully cryopreserved mostly from Mollusk and Echinoderm species, but also from 53 Arthropod and Cnidarian species (see for reviews Odintsova and Boroda, 2012 [17] and Paredes, 54 2015 [18]). Other biomaterials have been studied in a cryopreservation context, such as coral 55 fragments [19] and primmorphs from sponge cells [20]. The tolerance to various cryoprotectants of 56 tissue balls from corals was also investigated by Feuillassier et al. (2015) [21]. However, considering 57 the limited advancements in marine invertebrate cell cultures, cryopreservation of marine 58 invertebrate dissociated cells is seldom reported [18]. If they are, they largely focus on Mollusks 59 [22–30], with some other studies conducted on Echinoderms (e.g. dissociated cells from sea urchins 60 larvae [24,28]), on sponge cells [31–34] and on coral cells dissociated from embryos or larvae [35,36]. 61 Besides, none of these studies, except those on sponge cells, maintained cryopreserved cells in 62 cultures for more than a few days nor use cryopreserved cells for subsequent experiments.

63 One of the major factors that determines the success of a cryopreservation protocol is the type of 64 cryoprotecting agents (CPAs) used [16]. CPAs prevent damage to the cells from changes in osmotic 65 pressure and intracellular ice crystal formation. Among the various CPAs, dimethyl sulfoxide 66 (DMSO), a penetrating CPA, is the most common and widely used cryoprotectant to maintain frozen 67 cell lines. The precise mechanism by which it protects cells remains unclear; it has been suggested 68 that DMSO depresses the freezing point of cryosolutions [37], and that it can modulate the water 69 network hydrating the membrane hence reducing the stress induced by the volume changes of 70 water during freeze-thaw [38]. Penetrating CPAs could however induce some cytotoxicity due to the 71 disruption of intracellular signalling which results in cell death [39-41]. For marine invertebrate 72 studies, Paredes (2015) [18] reported in her review that DMSO was the most effective CPA for 70% of 73 the published work on germ cells, embryos and larvae compared to others CPAs such as glycerol. 74 This trend is also found for marine invertebrate dissociated cells [22,35,21,33]. In addition, in marine 75 invertebrate studies, DMSO was found as an efficient CPA on its own [22,29,35] but more frequently 76 in combination with other CPAs or with proteins, vitamins or sugar cocktails [18,24, 26,27,33]. 77 Indeed, the preservative capacity of DMSO was long known to be increased when serum, such as 78 foetal bovine serum (containing cocktail of proteins), is added to the cryosolution [42,43].

79 We have previously reported the establishment of primary cell cultures of a soft-body 80 cnidarian, the temperate sea anemone Anemonia viridis [9]. The established cell culture protocol 81 resulted in the maintenance of primary cell cultures with gastrodermal signature [15]. These cell 82 cultures were successfully used to assess the cellular response (e.g. viability) to environmental stress 83 [15] thus creating new perspectives for further fundamental, environmental and biotechnological 84 questions. An efficient cryopreservation procedure would therefore be an essential and powerful 85 tool for facilitating research in deciphering molecular mechanisms and cellular events in cnidarian 86 cells.

The aim of this study was therefore to design a cryopreservation protocol for primary gastrodermal *A. viridis* cell cultures in order to ensure a high post-thaw cell survival, preserving long-term recovery: cell viability, cell growth and physiological responses. All these advances will participate to raise the cnidarian cell cultures as a model system for marine invertebrate research perspectives.

92 2. Material and Methods

93 2.1. Biological Material

Five individuals of *Anemonia viridis* (Forskal 1775) were collected (prefectural authorization n°107 ; 02/28/2019) from 'Plage des ondes', Antibes, France, (43°33'17"N, 7°07'17.7"E), and maintained in a closed-circuit aquarium with artificial seawater (ASW) at 36-38 ‰ with Prodibio Expert Reef Salt, at 18.0 \pm 0.5 °C with weekly water changes. A LED bar (450 nm – Deckey LED aquarium) provided light at a constant saturating irradiance of 100 µmol m⁻²s⁻¹ (measured using a special sensor QSL-100, Biospherical Instruments Inc.) on a 12h:12h (light:dark) photoperiod. Sea anemones were fed once a week with oysters.

101 2.2. Primary cell cultures

102 From each A. viridis individual, an independent primary cell culture was obtained and 103 maintained as described in Ventura et al. (2018) [15]. Briefly, cell dissociation was performed 104 enzymatically with 0.15% collagenase type I (Sigma-Aldrich). Cells were cultured at 20.0 \pm 0.5 °C and 105 in the dark, in an optimized culture medium (CM) consisted of : 20% GMIM (Gibco), 5% foetal 106 bovine serum (FBS; PAA/GE Healthcare), 1% kanamycin (100 µg/ml, Sigma-Aldrich), 1% 107 amphotericin B (2.5 µg/ml; Interchim), 1% antibiotic antimycotic solution (Sigma-Aldrich), 1% L-108 glutamate (Sigma-Aldrich) and 71% of filtered ASW. The CM was adapted in respect to the 109 Mediterranean Seawater characteristics (i.e. salinity 40 ppt and pH 8.1). From day 3, culture medium 110 was replaced weekly and cells were seeded at 250,000 cells/ml in 12 well-plates.

111 2.3. Cryopreservation protocol

As cryoprotectant, DMSO (Sigma-Aldrich) was tested at two concentrations in the final CPA solution: 5% or 10% (following Munroe et al., 2018 [33]). DMSO was dissolved in the CM or in the CM enriched with foetal bovine serum (FBS) at 25% final. Control conditions without DMSO were also tested using CM enriched or not with FBS (i.e. 'CM' or 'CM + 25% FBS').

From day 17 after dissociation, the primary cell cultures were established with reliable cellular parameters [15]. By consequence, the cultivated cells were cryopreserved at different time points, from day 17 to 45 after cell dissociation. Each cryopreserved material contained 2 million cells that were placed in a cryotube containing 1 ml of the tested solution. Cryotubes were directly placed in a -80°C freezer (Ultra-Low Temperature VIP series, SANYO) and kept there for 8 to 87 days.

For thawing, cryotubes were removed from the -80°C freezer after the defined period and immediately transferred for 1-2 min into a water bath, pre-warmed at 20°C.

For seeding the cryopreserved cells, the cryotubes were centrifuged for 5 min at 1500 rpm. The supernatant was then removed, the cell pellets resuspended in the cell culture medium and seeded at 250,000 cells/ml in 12 well-plates [15].

126 2.4. Cell survival, cell viability, cell growth rate and cell size assessment

127 Cell survival was measured right after thawing cryopreserved cells, before reseeding. It was 128 determined as the percentage of viable cells relative to the 2 million cells initially cryopreserved. To 129 assess the number of viable cells, a sub-sample (100 μ l) of cryopreserved cells was harvested after 130 the thawing phase. Cell viability was assessed by evaluating the membrane integrity thanks to the 131 Evans blue method. Therefore, viable cells (unstained) and dead cells (stained) were identified and 132 counted. on a Neubauer improved haemocytometer (Sigma-Aldrich) using an optic microscope 133 (Zeiss Axio Imager Z1).

134 Cell viability was measured every week to monitor the cell culture health state overtime. A 135 sub-sample (100μ L) of cultivated cells was harvested weekly and using Evans blue method, viable 136 cells (unstained) and dead cells (stained) were identified and counted. The cell viability was defined 137 as the percentage of viable cells relative to total cells (i.e. viable and dead cells). In addition, two 138 complementary methods for cell viability assessment, i.e. overall enzymatic activity using the

- fluorescein diacetate (FDA) staining combined with a non-vital dye (Hoechst) and cell metabolic activity with 2-(4, 5-dimethyl-2-thiazolyl)-3, 5-diphenyl-2H tetrazolium bromide (MTT) assay, were also conducted (see details in Supplementary Material and Methods).
- 142 Cell growth rate was also assessed every week using the previous viable cells counts with 143 Evans blue method. The following formula was then used to calculate the 7-day averaged daily 144 growth rate:
- 145 Cell daily growth rate = $\frac{(\text{Viable cells}_{d+7} \text{Viable cells}_d)/\text{Viable cells}_d}{(d = day)}$, (d = day)

146 Cell growth rate and cell viability were monitored for each cell culture before and after 147 cryopreservation. The monitoring of these factors for cryopreserved cells was done by considering 148 that the age of the cells at the thawing time is the same age they were at the freezing time.

- Before and after cryopreservation, during the cell counts under optic microscope (objective x20), cells were measured, and cell sizes were scored with Zeiss microscope software Zen 2 (blue edition).
- 152 2.5. Hyperthermal stress experiment

153 In order to assess the maintenance of cryopreserved cells functionality, the response of 154 cryopreserved cells to a controlled stress experiment was investigated. The cultivated A. viridis cells 155 response to hyperthermal stress was assessed following the protocol published by Ventura et al. 156 (2018) [15]. Hyperthermal stress was induced in 12-well plates exposed to two different 157 temperatures: 20°C (control) and 28°C (hyperthermal condition), for 7 days. This experiment was 158 conducted either with non-cryopreserved or cryopreserved cells, 7 days after thawing. At least four 159 independent experiments were conducted from four primary cell cultures. For each assay, we 160 analyzed 3 wells as technical replicates.

161 2.6. Statistical analyses

162 All statistical analyses were conducted using the R v.3.6.0 software[44]. In order to assess the 163 effect of the cryosolutions on cell survival, to compare global viability and growth between 164 non-cryopreserved and cryopreserved cells, and to compare viability and growth rate values after 165 hyperthermal stress, either one-way ANOVA analyses were performed when parametric analyses 166 were possible (under normality and variance equality assumptions), or Kruskal-wallis when 167 non-parametric analyses were required. These analyses were followed, if necessary, by the 168 appropriate post-hoc, i.e. Tukey for ANOVA analyses, and Dunn for Kruskal-wallis analyses. Then, 169 to investigate the effect of the storage duration and the cell culture age on the cell survival, 170 correlation tests with linear regression model were conducted. Finally, repeated measures ANOVA 171 were conducted to compare cell viability and growth through time of non-cryopreserved and 172 cryopreserved cultures.

173 3. Results and Discussion

174 3.1. Success in set up of cryopreservation protocol on cell survival

175 The efficiency of the cryopreservation solutions was first evaluated with the percentage of cells 176 that survive a cryopreservation period of 10 ± 1 days and the subsequent thawing process. When 177 cells were cryopreserved in the culture medium (CM) or in CM with DMSO 5 or 10% the mean 178 percentage of cells that survived at -80°C, was around 7%. There were no significant differences 179 between these conditions (ANOVA; p>0.05) (Figure 1). However, when cells were cryopreserved in 180 the two conditions containing an FBS supplementation, we observed a significant higher percentage 181 of cell survival compared to non-enriched medium conditions (p<0.01 and p<0.0001 respectively for 182 the CM+25% FBS and for the DMSO 5% in the CM+25% FBS). Compared to FBS enriched CM alone 183 (± 45% survival rate), adding 5% of DMSO to the FBS enriched CM significantly enhanced the 184 survival rate of the cells to 67% (p<0.05; Figure 1). Thus, the latter constituted the optimal 185 cryosolution among those tested for cryopreserving cnidarian cells.





187 Figure 1. Percentage of cell survival following 10 ± 1 days of cryopreservation and thawing. For each 188 tested cryopreservation solution, the number of viable cells were compared to the 2 million cells 189 initially cryopreserved. Mean values and standard errors are represented, with n≥3 biological 190 replicates per condition. The ANOVA revealed significant differences between data ($p=1.14.10^{-10}$) and 191 the results from the Tukey post-hoc analysis are represented with letters: $a \neq b$ (p<0.01), $a \neq c$ 192 $(p < 0.0001), b \neq c (p < 0.05).$

193 As it was observed in most of marine invertebrate studies [18,24,26,27,33], DMSO was found to 194 be an efficient CPA for A. viridis cultivated cells when it combined with serum supplementation. The 195 reason could be that carbohydrates, lipids and proteins present in serums act as membrane 196 stabilizers, therefore they may help preventing membrane damage during the freezing process 197 [45-48].

198 The survival rate of the designed cryopreservation protocol for A. viridis cultivated cells with 199 the optimal cryosolution is comparable to those determined for dissociated cells [26,33] or other 200 biomaterials [18] from marine invertebrates, as well as the ones from vertebrate in vitro cells [49].

201 Since the optimal cryopreservation solution, among those tested in this study, was found to be 202 the DMSO 5% in FBS enriched culture medium, this solution was reused for different 203 cryopreservation durations in order to determine the influence of the cryopreservation time on cell 204 survival after thawing. The representative results for a A. viridis cell culture shown in Figure 2 205 demonstrated that there was no influence of the cryopreservation duration on the cell survival. 206 Indeed, from 8 to 87 days at -80°C, the percentage of cells that survived cryopreservation and 207 thawing, did not vary significantly (linear regression model; coefficient not statistically significantly 208 different from 0, p>0.05). The statistical analyses done on all available data (at least 3 biological 209 replicates) confirmed that there were no significant differences between the cryopreservation 210





212Figure 2. Percentage of cell survival following different cryopreservation durations and thawing. For213each tested cryopreservation duration, the number of viable cells were compared to the 2 million214initially cryopreserved cells. These representative results were obtained from one primary cell215culture cryopreserved at 24 days after initial seeding in 7 cryotubes, each cryotube being analyzed at216the end of a given cryopreservation duration (from 8 to 87 days). No significant differences between217times were found (linear regression model; r²=0.037; p=0.443).

Cells could therefore be cryopreserved for almost 3 months without any impact on the post-thaw survival compared to short-term cryopreservation, further validating the efficiency of the designed protocol. This is a major and significant progress for marine invertebrate cell cryopreservation. In fact, although few studies evaluated the cryopreservation efficiency after a storage of several weeks [22,24,25,28] and a maximum of 12 months [30], the majority only cryopreserved cells for a few hours to a few days [23,26,27,35,36,29,33].

Cryopreservation of vertebrate cell lines is well-advanced, and protocols allow to keep cryopreserved cells for years [16]. Therefore, additional experiments are required to assess if longer cryopreservation durations are possible and/or to determine the maximal duration using the protocol developed in this study. If a duration limitation were to be found, a long-term storage in liquid nitrogen, or -150°C freezers following the initial -80°C freezing should be envisioned.

229 In addition, we also assessed the influence of the age (time after initial seeding) of the primary 230 cell culture on the survival at thawing. Representative results of one primary cell culture, 231 cryopreserved at five different ages (from day 17 to day 45, Figure 3) revealed no significant 232 differences (linear regression model; coefficient not statistically significantly different from 0, 233 p>0.05). The statistical analyses done on all available data (3 biological replicates) confirmed that 234 there was no influence of the culture age on the cell survival (p>0.05; data not shown). One primary 235 cell culture can therefore be cryopreserved at each reseeding, as long as the cellular parameters 236 previously defined for "healthy" A. viridis primary cell culture are maintained, i.e. cell aggregates 237 formation, high viability and constant growth rate [15]. Being able to do so, considerable amounts of 238 cell stocks for each cell culture could be created.







since its initial seeding (from day 17 to day 45). At each age tested, one cryotube was analyzed. No significant differences between ages were found (linear regression model; $r^2=0.15$; p=0.512).

244 3.2. Absence of cryopreservation impact on cell recovery and functional parameters

245 To assess long-term cell recovery, we measured weekly, at each reseeding, the cell viability of 246 all primary cell cultures cryopreserved. Viability of cryopreserved cells was monitored for a period 247 going from 6 to 12 weeks after thawing and was compared to the viability of the corresponding 248 non-cryopreserved cell culture. Results of cell viability monitoring for 6 weeks after thawing for one 249 primary cell culture cryopreserved at day 24 for 8 or 36 days are presented in Figure 4. The data 250 show that cryopreserved cells were stably viable after thawing through time (>90% viability) with a 251 cell viability equivalent to that of the corresponding non-cryopreserved cell culture (repeated 252 measures ANOVA; p>0.05). Repeated measures ANOVA conducted on all biological replicates, 253 confirmed no differences in the cell viability over time between non-cryopreserved and 254 cryopreserved cells nor with the cryopreservation storage period (p>0.05; see Figure S1 and S2). The 255 mean of viability, over time, was maintained at 95 ± 1.9 % and at 96 ± 0.99 %, respectively in the 256 different non-cryopreserved and cryopreserved cell cultures monitored (ANOVA; p>0.05; see Figure 257 S3a). In addition, data obtained with the two complementary cell viability assays performed on 258 non-cryopreserved and cryopreserved cells at different time points of the kinetics confirmed all 259 these results, i.e. no differences in cell viability between non-cryopreserved and cryopreserved cells 260 and a maintenance over time of the cell viability (ANOVA ; p>0.05; see Figure S4).





262 Figure 4. Over time cell viability of a representative cell culture, comparing cryopreserved and 263 non-cryopreserved cells. Non-cryopreserved culture in grey dotted line, the same culture 264 cryopreserved at 24 days since initial seeding and thawed after 8 days of storage in grey solid line 265 and after 36 days in black solid line. The age of the cryopreserved cells at thawing is considered the 266 same as at the freezing time. Mean values of three technical replicates are shown with standard error 267 bars (although not visible because smaller than the data point symbols). Repeated measures ANOVA 268 revealed no significant differences in cell viability at each time between the non-cryopreserved 269 culture and the cryopreserved ones (p=0.582).

As a first functional parameter, the long-term cell growth was assessed. Indeed, we considered the resumption of cell cycle after cryopreservation as an essential functional parameter to explore. Non-cryopreserved and cryopreserved cells displayed a similar daily growth rate over time and independently of the cryopreservation duration (8 or 36 days) as shown by a representative cell culture in Figure 5 (repeated measures ANOVA; p>0.05).

The same result was obtained for all cell cultures tested in this study (see Figure S1 and S2). Interestingly, the cryopreserved cells stored for 79 days (Figure S2d) displayed a low initial growth rate one week after thawing and reseeding, suggesting that cryopreserved cells may need a longer time (between one and two weeks after thawing) to fully recover after such cryopreservation storage duration. Furthermore, the mean of the daily growth rate was maintained at 1.78 ± 0.39 and at 1.73 ± 0.33 , respectively in the different non-cryopreserved cell cultures and cryopreserved cell cultures monitored, with no significant differences between these two conditions (ANOVA; *p*>0.05; see Figure S3b).



Figure 5. Over time daily growth rate of a representative cell culture, comparing cryopreserved and non-cryopreserved cells. Non-cryopreserved cell culture in grey dotted line, the same culture cryopreserved at 24 days and thawed after either 8 days of storage in grey solid line or after 36 days in black solid line. Mean values of three technical replicates and standard error bars are shown. Repeated measures ANOVA revealed no significant differences in cell growth trend between the non-cryopreserved culture and the cryopreserved ones (*p*=0.722).

291 Long-term cell viability and growth rate monitoring of cryopreserved vs non-cryopreserved 292 cell cultures corroborated the healthy state of cryopreserved cells. Therefore, these analyses were 293 completed with weekly microscope observations in order to assess the cell culture behavior. In the 294 Figure 6, the comparison of a 31-day old non cryopreserved culture with the corresponding 295 cryopreserved culture (considering the thawing age is equal to the freezing age) showed that 296 cryopreserved cells, as the non-cryopreserved cells, form adherent cell aggregates, which is the 297 characteristic architecture of A. viridis primary cell cultures[9,15]. Besides, cryopreserved cells in 298 culture presented the same mean size (5.2 \pm 0.49 μ m) that non-cryopreserved cells (5.17 \pm 0.54 μ m; 299 ANOVA; p=0.957), suggesting no volume change after thawing (see Figure S5).



300 301

302

303

304

284

Figure 6. Observation of aggregates of *A. viridis* gastrodermal cells in culture before and after cryopreservation; **a**) Cell culture at day 31 since its establishment, and **b**) the same culture cryopreserved at day 24 for 79 days, 7 days after thawing. Phase contrast microscopy (objective x20), scale bars = 10 μm.

Therefore, cryopreserved cells maintained through time identical viability, growth and shape to the corresponding non-cryopreserved cell culture, and this independently of the cryopreservation duration. These results indicate that the cryopreservation designed protocol stored cells in a healthy state allowing them to fully recover after thawing and to behave like their origin culture. This monitoring constitutes an essential part in validating the storage protocol for further use of cryopreserved cells and represents a major strength of this study. Indeed, the reseeding of thawed cells was rarely done on previous marine invertebrate studies, and cells are usually maintained only for a few hours to a few days [22–24,28,25,27], with a maximum of 15 days for bivalve cells in Dessai (2018) [30] and around 40 days for sponge cells [32,34].

314 As a second parameter of the cell functionality after cryopreservation, we investigated the 315 response of cryopreserved cells to a controlled stress experiment. In Cnidarians, and more 316 particularly in our research model, A. viridis, hyperthermal stress is well known to induce oxidative 317 damages and cell death, i.e. apoptosis [50,51]. Using A. viridis primary cell cultures, we previously 318 reported that hyperthermia (+8°C) didn't induce any oxidative damage or impact on survival but 319 provoked a drastic decrease of cell growth [15]. Thus, in this study, we compared the response of 320 cryopreserved and non-cryopreserved cells submitted to the same hyperthermal stress, in terms of 321 viability and growth. The results show that hyperthermal stress didn't change significantly the cell 322 viability neither for non-cryopreserved cells nor for cryopreserved cells (ANOVA; p>0.05) (Figure 323 7a). Moreover, cell growth rate was drastically decreased by around 80% after 7 days at 28°C 324 compared to 20°C condition (Kruskal-wallis; p<0.001), without any significant differences between 325 non-cryopreserved and cryopreserved cells (Kruskal-wallis; *p*>0.05) (Figure 7b). These results are in 326 line with data from Ventura et al. (2018) [15] and show that cryopreserved cells displayed identical 327 resistance to non-cryopreserved cells, strongly corroborating the non-alteration of cell functionality. 328 Therefore, not only A. viridis cells in vitro can be successfully cryopreserved and reseeded, but they 329 can be reliably used for further experiments, like it is done for mammalian cells. Although functional 330 analyses are sometimes conducted for marine invertebrate cells, through metabolic and enzymatic 331 activities [23,26,24,28,25,30], only some sponge cells studies conducted experiments using 332 cryopreserved cells [31,34]. However, this is a fundamental assessment in order to validate the 333 cryopreservation protocol for creating reliable models.



Figure 7. Comparison of cryopreserved and non-cryopreserved cells in response to hyperthermia. Assessment of cell viability (**a**) and growth (**b**) of *A. viridis* cells in response to a hyperthermal stress of +8°C (black bars) for 7 days, for non-cryopreserved cultures (striped bars) and for cryopreserved cultures (filled bars). Cell viability and growth are expressed relative to control condition (20°C – grey bars). Mean values with standard error bars are shown, biological replicates n≥4. The asterisks represent the significant differences between control and stress conditions (*** Kruskall-Wallis: p<0.001).

342 4. Conclusions

334

343 In this study we succeeded to design an easy and rapid cryopreservation procedure for 344 Anemonia viridis primary cell cultures. The established protocol enabled us to obtain high cell 345 survival after thawing and a full long-term recovery of the cell culture behavior. The development of 346 cryopreservation in cnidarian primary cell cultures enables us to preserve stable cell stocks available 347 shortly after thawing for experimental procedures and sharing with the scientific community. This 348 new tool will be an important asset to raise A. viridis primary cell cultures as a powerful model for 349 studying and understanding the cnidarian properties (i.e. symbiosis lifestyle, response to stress, 350 aging), difficult to study in most cnidarian models at the molecular and cellular levels.

351 **Supplementary Materials**



353 Supplementary Material and Methods

354 Complementary cell viability assays

355 FDA (Fluorescein Diacetate, 4 µg/mL; Sigma-Aldrich) and Hoechst 33342 (5 µg/mL; Sigma-Aldrich) were 356 added and incubated with cells during 15 minutes at 20°C in the dark. Viable cells (fluorescent in green) 357 and dead cells (fluorescent in blue) were identified and counted on a Neubauer improved 358 haemocytometer (Sigma-Aldrich) using a fluorescence microscope (Zeiss Axio Imager Z1). The cell 359 viability was defined as the percentage of viable cells relative to total cells (i.e. viable and dead cells).

360 MTT assay was performed following manufacturer instructions with slight modifications. Briefly, prior the 361 assay 60 000 cells were seeded in triplicate in 96-well plate in 100 µL of culture medium and incubated for 24h. 362 20 µL of 5mg/mL MTT solution (Sigma-Aldrich) is then added to each well and incubated 5h at 20°C in the 363 dark. Then, the supernatant is removed, and the yielded formazan was dissolved in the suitable detergent 364 (isopropanol) for 15 minutes. Subsequently, the plates' light absorption (OD) is read at wavelength 590 nm on 365 spectrofluorometer (SAFAS, Monaco). The cell viability was expressed as follow: Viability % = (cryopreserved 366 cells OD /non-cryopreserved cells OD) ×100





Supplementary Figures













380 Figure S2. Over time viability (a, b) and daily growth rate (c, d) of two primary cell cultures, comparing 381 non-cryopreserved cells (NC) and cells cryopreserved for a long storage period (C). The cryopreservation 382 storage period is mentioned in the legend, it goes from 25 days ('25d') to 79 days ('79d'). Each graph per 383 variable measured represent one primary cell culture and its corresponding cryopreserved cell culture(s). 384 Time points are mentioned as "Week after freezing (NC) or thawing (C) time" i.e. the monitoring 385 presented here for the non-cryopreserved cultures begins after the freezing time, and the one for 386 cryopreserved cells is done after thawing, thus considering the age at that time is the same age as at the 387 freezing time.



Figure S3. Box plots of over time viability (**a**) and daily growth rate (**b**) of all cell cultures monitored (non-cryopreserved in gray and cryopreserved in blue). Time points are identical to the ones present in Fig. S1 and S2 and are here going from 1 to 7 weeks, which represent the time points after freezing (for non-cryopreserved cells) or after thawing (for cryopreserved cells) where at least 3 biological replicates were monitored.

395Figure S4. Comparison of cell viability between 38-day old non-cryopreserved cell cultures ('NC39638d'), and cryopreserved cell cultures, either 38-day old ('C 38d') or 80-day old ('C 80d'). Cell397viability was measured with (a) FDA/Hoechst staining and with (b) MTT assay. Mean values with398standard error bars are shown (n≥3). ANOVA analyses revealed no significant differences in the399viability values between non-cryopreserved and cryopreserved cells, and over time between400cryopreserved cells (p=0.45 for panel a; p=0.378 for panel b).

407 Figure S5. Observation, on Neubauer improved haemocytometer under optic microscope (objective x20), of the same *A. viridis* gastrodermal cell culture (a) before and (b) after cryopreservation and stained with Evans blue (i.e. the dead cells stained in blue) (scale bar = 10 μm).

410

- 412 the sea anemones. Authors thank also Brigitte Poderini, Thamilla Zamoum and Maxence Burtin for aquaria 413 maintenance and taking care of the animals. Authors are grateful to Pauline Cotinat for her help in the cell 414 culture establishment.
- 415 **Funding:** This work was supported by the SATT Sud-Est (project number: 1082-SA-18-CNRS).
- 416 Author Contributions: C.F. did all the investigations (experimental work) and the statistical analyses. P.F. and

417 S.B-V. designed and supervised the research. C.F., P.F. and S.B-V. wrote the manuscript. Finally, E.R., P.F. and

- 418 S.B-V. reviewed the manuscript.
- 419 **Conflicts of Interest:** On behalf of all authors, the corresponding author states that there is no conflict of 420 interest.
- 421 References
- Nandi, S. Animal Cell Culture and Virology; New India Publishing, 2009; ISBN
 978-93-80235-05-9.
- 424 2. Maramorosch, K.; Sato, G.H. *Advances in Cell Culture*; Elsevier, 2014; ISBN
 425 978-1-4832-1506-8.
- 426 3. Hall, B.K.; Moody, S.A. Cells in Evolutionary Biology: Translating Genotypes into
 427 Phenotypes Past, Present, Future; CRC Press, 2018; ISBN 978-1-351-65202-5.
- 428 4. Vallone, D.; Santoriello, C.; Gondi, S.B.; Foulkes, N.S. Basic Protocols for Zebrafish
 429 Cell Lines. In *Circadian Rhythms: Methods and Protocols*; Rosato, E., Ed.; Methods in
 430 Molecular Biology[™]; Humana Press: Totowa, NJ, 2007; pp. 429–441 ISBN
 431 978-1-59745-257-1.
- 432 5. Lynn, D.E. Novel techniques to establish new insect cell lines. *In Vitro Cellular & Developmental Biology-Animal* 2001, *37*, 319–321.
- 434 6. Rinkevich, B. Marine Invertebrate Cell Cultures: New Millennium Trends. *Mar*435 *Biotechnol* 2005, 7, 429–439, doi:10.1007/s10126-004-0108-y.
- 436 7. Rinkevich, B. Cell Cultures from Marine Invertebrates: New Insights for Capturing
 437 Endless Stemness. *Mar Biotechnol* 2011, *13*, 345–354,
 438 Liste 4007/ 10100, 040, 00514.0
- 438 doi:10.1007/s10126-010-9354-3.
- 439 8. Cai, X.; Zhang, Y. Marine invertebrate cell culture: a decade of development. *J*440 *Oceanogr* 2014, *70*, 405–414, doi:10.1007/s10872-014-0242-8.
- Barnay-Verdier, S.; Dall'osso, D.; Joli, N.; Olivré, J.; Priouzeau, F.; Zamoum, T.; Merle,
 P.-L.; Furla, P. Establishment of primary cell culture from the temperate symbiotic
 cnidarian, Anemonia viridis. *Cytotechnology* 2013, *65*, 697–704,
 doi:10.1007/s10616-013-9566-2.
- 445 10. Huete-Stauffer, C.; Valisano, L.; Gaino, E.; Vezzulli, L.; Cerrano, C. Development of
- 446 long-term primary cell aggregates from Mediterranean octocorals. In Vitro
- 447 *Cell.Dev.Biol.-Animal* **2015**, *51*, 815–826, doi:10.1007/s11626-015-9896-9.
- Lecointe, A.; Cohen, S.; Gèze, M.; Djediat, C.; Meibom, A.; Domart-Coulon, I.
 Scleractinian coral cell proliferation is reduced in primary culture of suspended
 multicellular aggregates compared to polyps. *Cytotechnology* 2013, *65*, 705–724,
 doi:10.1007/s10616-013-9562-6.
- 452 12. Mercurio, S.; Di Benedetto, C.; Sugni, M.; Candia Carnevali, M.D. Primary cell cultures
 453 from sea urchin ovaries: a new experimental tool. *In Vitro Cell.Dev.Biol.-Animal* 2014,
 454 50, 139–145, doi:10.1007/s11626-013-9686-1.

455	13.	Rabinowitz, C.; Moiseeva, E.; Rinkevich, B. In vitro cultures of ectodermal monolayers
456		from the model sea anemone Nematostella vectensis. Cell Tissue Res 2016, 366,
457		693–705, doi:10.1007/s00441-016-2495-6.
458	14.	Vandepas, L.E.; Warren, K.J.; Amemiya, C.T.; Browne, W.E. Establishing and
459		maintaining primary cell cultures derived from the ctenophore Mnemiopsis leidyi.
460		Journal of Experimental Biology 2017 , 220, 1197–1201, doi:10.1242/jeb.152371.
461	15.	Ventura, P.; Toullec, G.; Fricano, C.; Chapron, L.; Meunier, V.; Röttinger, E.; Furla, P.;
462		Barnay-Verdier, S. Cnidarian Primary Cell Culture as a Tool to Investigate the Effect of
463		Thermal Stress at Cellular Level. Mar. Biotechnol. 2018, 20, 144–154,
464		doi:10.1007/s10126-017-9791-3.
465	16.	Cryopreservation and Freeze-Drying Protocols; Wolkers, W.F., Oldenhof, H., Eds.;
466		Methods in Molecular Biology; Springer: New York, 2015; Vol. 1257; ISBN
467		978-1-4939-2192-8.
468	17.	Odintsova, N.A.; Boroda, A.V. Cryopreservation of the cells and larvae of marine
469		organisms. Russ J Mar Biol 2012, 38, 101–111, doi:10.1134/S1063074012020083.
470	18.	Paredes, E. Exploring the evolution of marine invertebrate cryopreservation –
471		Landmarks, state of the art and future lines of research. Cryobiology 2015, 71,
472		198–209, doi:10.1016/j.cryobiol.2015.08.011.
473	19.	Hagedorn, M.; Farrell, A.; Carter, V.L. Cryobiology of coral fragments. Cryobiology
474		2013 , <i>66</i> , 17–23, doi:10.1016/j.cryobiol.2012.10.003.
475	20.	Mussino, F.; Pozzolini, M.; Valisano, L.; Cerrano, C.; Benatti, U.; Giovine, M.
476		Primmorphs Cryopreservation: A New Method for Long-Time Storage of Sponge Cells.
477		Marine biotechnology (New York, N.Y.) 2012, 15, doi:10.1007/s10126-012-9490-z.
478	21.	Feuillassier, L.; Masanet, P.; Romans, P.; Barthélémy, D.; Engelmann, F. Towards a
479		vitrification-based cryopreservation protocol for the coral Pocillopora damicornis L.:
480		Tolerance of tissue balls to 4.5 M cryoprotectant solutions. Cryobiology 2015, 71,
481		224–235, doi:10.1016/j.cryobiol.2015.07.004.
482	22.	Le Marrec-Croq, F.; Fritayre, P.; Chesné, C.; Guillouzo, A.; Dorange, G.
483		Cryopreservation of Pecten maximus Heart Cells. Cryobiology 1998, 37, 200–206,
484		doi:10.1006/cryo.1998.2113.
485	23.	Cheng, T.C.; La Peyre, J.F.; Buchanan, J.T.; Tiersch, T.R.; Cooper, R.K.
486		Cryopreservation of heart cells from the eastern oyster. In Vitro Cell.Dev.BiolAnimal
487		2001 , 37, 237–243, doi:10.1007/BF02577536.
488	24.	Odintsova, N.; Kiselev, K.; Sanina, N.; Kostetsky, E. Cryopreservation of primary cell
489		cultures of marine invertebrates. CryoLetters 2001, 22, 299-310.
490	25.	Poncet, JM.; Serpentini, A.; Boucaud-Camou, E.; Lebel, JM. Cryopreservation of
491		mantle dissociated cells from Haliotis tuberculata (Gastropoda) and postthawed
492		primary cell cultures. Cryobiology 2002, 44, 38–45,
493		doi:10.1016/S0011-2240(02)00001-9.
494	26.	Hanquet-Dufour, A.C.; Kellner, K.; Heude, C.; Naimi, A.; Mathieu, M.; Poncet, J.M.
495		Cryopreservation of Crassostrea gigas vesicular cells: Viability and metabolic activity.
496		Cryobiology 2006, 53, 28–36, doi:10.1016/j.cryobiol.2006.03.008.

497 27. Kostetsky, E.Y.; Boroda, A.V.; Odintsova, N.A. Changes in the lipid composition of 498 mussel (Mytilus trossulus) embryo cells during cryopreservation. BIOPHYSICS 2008, 499 53, 299-303, doi:10.1134/S000635090804012X. 500 28. Odintsova, N.A.; Boroda, A.V.; Velansky, P.V.; Kostetsky, E.Ya. The fatty acid profile 501 changes in marine invertebrate larval cells during cryopreservation. Cryobiology 2009, 502 59, 335–343, doi:10.1016/j.cryobiol.2009.09.006. 503 29. Daugavet, M.A.; Blinova, M.I. Culture of mussel (Mytiuls edulis L.) mantle cells. Cell 504 Tiss. Biol. 2015, 9, 233–243, doi:10.1134/S1990519X15030037. 505 30. Dessai, S.N. Cryopreservation of cultured mantle cells of Paphia malabarica for 506 perennial availability. Cryobiology 2018, 82, 93-98, 507 doi:10.1016/j.cryobiol.2018.04.002. 508 Pomponi, S.A. Development of techniques for in vitro production of bioactive natural 31. 509 products from marine sponges. Invertebrate Cell Culture : Novel Directions and 510 Biotechnology Applications 1997, 231–237. 511 32. Pomponi, S.A.; Willoughby, R. Development of sponge cell cultures for biomedical 512 application. In Aquatic invertebrate cell culture; Springer: London, 2000; pp. 323-336. 513 33. Munroe, S.; Martens, D.E.; Sipkema, D.; Pomponi, S.A. Comparison of 514 Cryopreservation Techniques For Cells of the Marine Sponge Dysidea Etheria 515 Available online: 516 https://www.ingentaconnect.com/content/cryo/cryo/2018/00000039/00000004/art000 517 07 (accessed on Dec 19, 2019). 518 34. Conkling, M.; Hesp, K.; Munroe, S.; Sandoval, K.; Martens, D.E.; Sipkema, D.; Wijffels, 519 R.H.; Pomponi, S.A. Breakthrough in Marine Invertebrate Cell Culture: Sponge Cells 520 Divide Rapidly in Improved Nutrient Medium. Sci Rep 2019, 9, 1-10, 521 doi:10.1038/s41598-019-53643-y. 522 35. Hagedorn, M.; Carter, V.; Martorana, K.; Paresa, M.K.; Acker, J.; Baums, I.B.; 523 Borneman, E.; Brittsan, M.; Byers, M.; Henley, M.; et al. Preserving and Using 524 Germplasm and Dissociated Embryonic Cells for Conserving Caribbean and Pacific 525 Coral. PLOS ONE 2012, 7, e33354, doi:10.1371/journal.pone.0033354. 526 Hagedorn, M.; van Oppen, M.J.H.; Carter, V.; Henley, M.; Abrego, D.; Puill-Stephan, 36. 527 E.; Negri, A.; Heyward, A.; MacFarlane, D.; Spindler, R. First frozen repository for the 528 Great Barrier Reef coral created. Cryobiology 2012, 65, 157–158, 529 doi:10.1016/j.cryobiol.2012.05.008. 530 37. Stolzing, A.; Naaldijk, Y.; Fedorova, V.; Sethe, S. Hydroxyethylstarch in 531 cryopreservation – Mechanisms, benefits and problems. Transfusion and Apheresis 532 Science 2012, 46, 137-147, doi:10.1016/j.transci.2012.01.007. 533 38. Cheng, C.-Y.; Song, J.; Pas, J.; Meijer, L.H.H.; Han, S. DMSO Induces Dehydration 534 near Lipid Membrane Surfaces. Biophysical Journal 2015, 109, 330-339, 535 doi:10.1016/j.bpj.2015.06.011. 536 39. Fahy, G.M. The relevance of cryoprotectant "toxicity" to cryobiology. Cryobiology **1986**, 537 23, 1-13, doi:10.1016/0011-2240(86)90013-1.

- Song, Y.C.; Khirabadi, B.S.; Lightfoot, F.; Brockbank, K.G.M.; Taylor, M.J. Vitreous
 cryopreservation maintains the function of vascular grafts. *Nature Biotechnology* 2000, *18*, 296–299, doi:10.1038/73737.
- 541 41. Elmoazzen, H.Y.; Poovadan, A.; Law, G.K.; Elliott, J.A.W.; McGann, L.E.; Jomha, N.M.
 542 Dimethyl sulfoxide toxicity kinetics in intact articular cartilage. *Cell Tissue Banking*543 2006, 8, 125, doi:10.1007/s10561-006-9023-y.
- Schaefer, V.W.; Dicke, K.A. Preservation of haemopoietic stem cells, transplantation
 potential and CFU-c activity of frozen marrow tested in mice, monkeys and man. In *Cryopreservation of Normal and Neoplastic Cells*; Weiner, R.S., Oldham, R.K.,
- 547 Schwartzenberg, L., Eds.; Paris : INSERM, 1973; pp. 63–69.
- 43. Grilli, G.; Porcellini, A.; Lucarelli, G. Role of serum on cryopreservation and subsequent
 viability of mouse bone marrow hemopoietic stem cells. *Cryobiology* **1980**, *17*,
 550 516–520, doi:10.1016/0011-2240(80)90063-2.
- 44. R Core Team *R: A language and environment for statistical computing.*; R Foundation
 for Statistical Computing: Vienna, Austria, 2020;
- 45. Anchordoguy, T.; Carpenter, J.F.; Loomis, S.H.; Crowe, J.H. Mechanisms of
 interaction of amino acids with phospholipid bilayers during freezing. *Biochimica et Biophysica Acta (BBA) Biomembranes* **1988**, *946*, 299–306,
- 556 doi:10.1016/0005-2736(88)90405-1.
- 46. Cabrita, E.; Anel, L.; Herraéz, M.P. Effect of external cryoprotectants as membrane
 stabilizers on cryopreserved rainbow trout sperm. *Theriogenology* 2001, *56*, 623–635,
 doi:10.1016/S0093-691X(01)00594-5.
- 47. He, S.; Woods III, L.C. Effects of dimethyl sulfoxide and glycine on cryopreservation
 induced damage of plasma membranes and mitochondria to striped bass (Morone
 saxatilis) sperm. *Cryobiology* 2004, *48*, 254–262, doi:10.1016/j.cryobiol.2004.01.009.
- 48. Marco-Jiménez, F.; Garzón, D.L.; Peñaranda, D.S.; Pérez, L.; Viudes-de-Castro, M.P.;
 Vicente, J.S.; Jover, M.; Asturiano, J.F. Cryopreservation of European eel (Anguilla
 anguilla) spermatozoa: Effect of dilution ratio, foetal bovine serum supplementation,
 and cryoprotectants. *Cryobiology* **2006**, *53*, 51–57, doi:10.1016/j.cryobiol.2006.03.011.
- 567 49. D. Benson, J. Modeling and Optimization of Cryopreservation. In *Cryopreservation and*568 *Freeze-Drying Protocols*; Wolkers, W.F., Oldenhof, H., Eds.; Methods in Molecular
 569 Biology; Springer: New York, NY, 2015; pp. 83–120 ISBN 978-1-4939-2193-5.
- 570 50. Richier, S.; Rodriguez-Lanetty, M.; Schnitzler, C.E.; Weis, V.M. Response of the 571 symbiotic cnidarian Anthopleura elegantissima transcriptome to temperature and UV 572 increase. *Comparative Biochemistry and Physiology Part D: Genomics and* 573 *Proteomics* 3, 283–289.
- 574 51. Richier, S.; Sabourault, C.; Courtiade, J.; Zucchini, N.; Allemand, D.; Furla, P.
 575 Oxidative stress and apoptotic events during thermal stress in the symbiotic sea
 576 anemone, Anemonia viridis. *FEBS J.* 2006, 273, 4186–4198,
- 577 doi:10.1111/j.1742-4658.2006.05414.x.
- 578 52. Moya, A.; Ganot, P.; Furla, P.; Sabourault, C. The transcriptomic response to thermal 579 stress is immediate, transient and potentiated by ultraviolet radiation in the sea

- 580 anemone *Anemonia viridis*: MOLECULAR BASIS OF CNIDARIAN BLEACHING.
- 581 *Molecular Ecology* **2012**, *21*, 1158–1174, doi:10.1111/j.1365-294X.2012.05458.x.