

In vitro evaluation of polymeric nanoparticles with a fluorine core for drug delivery triggered by focused ultrasound

L. Somaglino^{a,b}, L. Mousnier^c, A. Giron^a, W. Urbach^{d,a}, N. Tsapis^c, N. Taulier^{a,*}

^a Sorbonne Université, CNRS, INSERM, Laboratoire d'Imagerie Biomédicale, LIB, F-75006 Paris, France

^b IFREMER, La Seyne-sur-Mer, France

^c Université Paris-Saclay, CNRS, Institut Galien Paris Sud, 92296, Châtenay-Malabry, France

^d Laboratoire de Physique de l'École Normale Supérieure, ENS, Université PSL, CNRS, Sorbonne Université, Université de Paris, F-75005 Paris, France.

Figure S1. Snapshot of the magnitudes (in dB) of the Fast-Fourier Transform of the acoustic signal amplitude recorded by the focused hydrophone at acoustic peak-to-peak pressures of 0.5 MPa (in orange) and at 7.8 MPa (in blue). The last curve exhibits a broad band signal with harmonics that is characteristic of a strong inertial cavitation. No difference was observed between samples containing or not nanoparticles or liposomes.

