

# In-Cell Discovery of RNA–Protein Interactions

## Xavier Sabaté-Cadenas, Alena Shkumatava

### ▶ To cite this version:

Xavier Sabaté-Cadenas, Alena Shkumatava. In-Cell Discovery of RNA–Protein Interactions. Trends in Biochemical Sciences, 2020, 45 (3), pp.272 - 273. 10.1016/j.tibs.2019.12.003 . hal-03489898

## HAL Id: hal-03489898 https://hal.science/hal-03489898

Submitted on 21 Jul 2022

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Distributed under a Creative Commons Attribution - NonCommercial 4.0 International License

## **In-cell discovery of RNA-protein interactions**

Xavier Sabaté-Cadenas<sup>1</sup> and Alena Shkumatava<sup>1\*</sup>

<sup>1</sup>Institut Curie, PSL Research University, INSERM U934, CNRS UMR 3215, F-75005, Paris, France

\*Correspondence: alena.shkumatava@curie.fr

#### **TWO-SECTION SUMMARY**

Numerous methods have been developed to define proteins bound to an RNA of interest [1-6]. Existing techniques commonly rely on probe-based RNA capture and identification of copurified proteins by mass spectrometry, which often present limitations associated with crosslinking, RNA enrichment, and protein detection. The incPRINT (in cell protein-RNA interaction) method overcomes these challenges by using stably expressed recombinant luciferase tethered to an MS2-tagged RNA of interest through MS2 coat protein (MS2CP) [7]. If an interaction between a test RNA and a test protein occurs in the cell, it is measured by luminescence upon cell lysis and immunoprecipitation of the tagged protein.

To identify RNA-bound proteomes, high-throughput immunoprecipitation of thousands of tagged proteins [8-10] is carried out, followed by luminescence-based detection of their interactions with an individual RNA. The approach is suitable for discovery of proteins interacting with full-length transcripts or with specific RNA regions.

#### **ADVANTAGES:**

- Does not rely on direct RNA pulldowns, thus circumventing a step often plagued by low efficiency.
- Avoids crosslinking and the associated inefficiency and biases.
- Not limited by low expression levels of test RNA or proteins.
- Can identify proteomes of specific RNA regions and provide RNA-region-specific resolution of protein binding.
- Can identify transient RNA-protein interactions.
- One cell line is used for all applications.

- The protein library is not limited in size and can be expanded or reduced depending on experimental needs.
- Flexibility of throughput ranging from high to low allows usage of incPRINT as an orthogonal method for validation of individual RNA-protein interactions.
- Quantitative luciferase enables structure-function analyses and analyses of mutated or disordered RNA-protein interactions.

### **CHALLENGES:**

- incPRINT is a primary binding assay that does not consider the developmental context or timing of RNA-protein interactions.
- A subset of proteins may interact with MS2 stem loops but not with the test RNA.
- Protein tags may alter protein expression or binding.
- Some interactions may be indirect; additional validation assays are required to establish direct RNA binding.
- Some RNA-protein interactions may require additional factors that are not expressed in the cell line used for the assay.
- Ectopic expression may affect sub-cellular localization of test RNAs and proteins.

## LITERATURE

1. Chu, C. et al. (2015) Systematic discovery of Xist RNA binding proteins. Cell 161 (2), 404-16.

2. McHugh, C.A. et al. (2015) The Xist lncRNA interacts directly with SHARP to silence transcription through HDAC3. Nature 521 (7551), 232-6.

3. Ramanathan, M. et al. (2018) RNA-protein interaction detection in living cells. Nat Methods 15 (3), 207-212.

4. Lee, H.Y. et al. (2013) RNA-protein analysis using a conditional CRISPR nuclease. Proc Natl Acad Sci U S A 110 (14), 5416-21.

5. Zielinski, J. et al. (2006) In vivo identification of ribonucleoprotein-RNA interactions. Proc Natl Acad Sci U S A 103 (5), 1557-62.


6. Slobodin, B. and Gerst, J.E. (2010) A novel mRNA affinity purification technique for the identification of interacting proteins and transcripts in ribonucleoprotein complexes. RNA 16 (11), 2277-90.


7. Graindorge, A. et al. (2019) In-cell identification and measurement of RNA-protein interactions. Nat Commun 10 (1), 5317.

8. Baltz, A.G. et al. (2012) The mRNA-bound proteome and its global occupancy profile on protein-coding transcripts. Mol Cell 46 (5), 674-90.


9. Castello, A. et al. (2012) Insights into RNA biology from an atlas of mammalian mRNAbinding proteins. Cell 149 (6), 1393-406.

10. Taipale, M. et al. (2012) Quantitative analysis of HSP90-client interactions reveals principles of substrate recognition. Cell 150 (5), 987-1001.





~1500 RNA-binding proteins




~1300 transcription factors



~170 chromatin-associated proteins





 $\rightarrow$ 

Data analysis and normalization

Interaction No interaction

