



HAL
open science

Sonoporation

Spiros Kotopoulos, Anthony Delalande, Chantal Pichon, Michiel Postema

► **To cite this version:**

Spiros Kotopoulos, Anthony Delalande, Chantal Pichon, Michiel Postema. Sonoporation: Using ultrasound for targeted drug delivery. 11th US-Japan Symposium on Drug Delivery Systems, Dec 2011, Lahaina, Maui, Hawaii, United States. hal-03576189

HAL Id: hal-03576189

<https://hal.science/hal-03576189>

Submitted on 15 Feb 2022

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

Sonoporation

Using ultrasound for targeted drug delivery

Spiros Kotopoulos^{1,2}, Anthony Delalande^{3,2}, Chantal Pichon⁴

Michiel Postema²

¹ National Center for Ultrasound in Gastroenterology, Bergen, Norway

² Department of Physics and Technology, University of Bergen, Norway

³ Institute of Medicine, Hematology Section, University of Bergen, Bergen, Norway

⁴ Centre de Biophysique Moléculaire, CNRS affiliated to the University of Orléans, France

Introduction

Sonoporation is the transient permeabilisation and resealing of a cell membrane with the help of ultrasound and/or an ultrasound contrast agent.

Ultrasound contrast agents consist of elastic shelled gas-filled microbubbles. These are used to enhance the echogenicity of blood in ultrasound diagnostics. In a specific ultrasound field these microbubbles oscillate. This oscillation creates cell membrane permeabilisation allowing increased uptake of drugs.

As ultrasound can be focused, the cell membrane stimulation can be targeted, for example malignant tumors.

Our work focused on understanding the biophysical mechanisms behind sonoporation in order to improve its efficiency.




Fig. 1 - Microbubble oscillation in an ultrasound field.

Methods

We sonicated HeLa cancer cells in the presence of fluorescent labeled Micro-Marker. An MI of 0.2 at 6.6MHz was used, well within the clinical diagnostic range. We recorded the behaviour of the microbubbles using a high-speed camera and a confocal microscope (Fig.2). After sonication a z-stack was performed in order to precisely locate the microbubbles.




Fig. 2 - Schematic illustration showing configuration for simultaneous high-speed photography and fluorescent imaging.

Results

Within several seconds of sonication, microbubbles were observed to be attracted to cells and subsequently penetrate the cell membrane (Fig. 3). The z-stack proved that the bubble had entered the cell and subsequently dissolved (Fig. 4).




Fig. 3 - Microbubble entering cell recorded using a high-speed camera (top) and schematic representation (bottom). Each frame corresponds to 37 x 37 (μm)².




Fig. 4 - Location of fluorescence in two HeLa cells before and after sonoporation. After sonication fluorescence was located in the center of the cell.




Fig. 5 - Schematic representation of 5 of the phenomena involved in sonoporation.

Conclusions

It has been demonstrated that lipid-shelled microbubbles can be forced to enter cells at a low MI. Hence, if a therapeutic load is added to the bubble, ultrasound guided delivery could be facilitated at diagnostic settings. In addition, these results may have implications for the safety regulations on the use of ultrasound contrast agents for diagnostic imaging.