



HAL
open science

Phosphatidylinositol 4-phosphate: a key determinant of plasma membrane identity and function in plants

Vedrana Marković, Yvon Jaillais

► **To cite this version:**

Vedrana Marković, Yvon Jaillais. Phosphatidylinositol 4-phosphate: a key determinant of plasma membrane identity and function in plants. *New Phytologist*, inPress, 235 (3), pp.867-874. 10.1111/nph.18258 . hal-03681487

HAL Id: hal-03681487

<https://hal.science/hal-03681487v1>

Submitted on 30 May 2022

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

1 Phosphatidylinositol 4-phosphate: a key determinant of plasma membrane identity and 2 function in plants

3
4 Vedrana Marković¹ and Yvon Jaillais^{1*}

5
6 Vedrana Marković <https://orcid.org/0000-0002-0310-5102>

7 Yvon Jaillais <https://orcid.org/0000-0003-4923-883X>

8
9
10 ¹Laboratoire Reproduction et Développement des Plantes, Univ Lyon, ENS de Lyon, CNRS,
11 INRAE, F-69342, Lyon, France

12
13 * Correspondance : yvon.jaillais@ens-lyon.fr

14
15 **Summary:** Phosphatidylinositol 4-phosphate (PI4P) is an anionic phospholipid, which has
16 been described as a master regulator of the Golgi apparatus in eukaryotic cells. However, recent
17 evidence suggests that PI4P mainly accumulates at the plasma membrane in all plant cells
18 analyzed so far. Besides, many functions that are typically attributed to phosphatidylinositol
19 4,5-bisphosphate [PI(4,5)P₂] in animal and yeast cells are also supported by PI4P in plants. For
20 example, PI4P is the key anionic lipid that powers the strong electrostatic properties of the
21 plasma membrane. PI4P is also required for the establishment of stable membrane contacts
22 between the endoplasmic reticulum (ER) and the plasma membrane, for exocytosis and to
23 support signaling pathways. Thus, we propose that PI4P has a prominent role in specifying the
24 identity of the plasma membrane and in supporting some of its key functions and should be
25 considered a hallmark lipid of this compartment.

26
27 **Key words:** plant cell biology, membrane lipid, PtdIns4P, phosphoinositides, trafficking,
28 endocytosis, nanodomains, membrane contact sites.

29 30 31 **Introduction.**

32 Within the endomembrane system of eukaryotic cells, anionic lipids are low abundant lipids
33 but they are key landmarks that contribute to the identity of membrane compartments (Noack
34 & Jaillais, 2017, 2020). Anionic phospholipids have a negatively charged head group that each
35 can recruit a suite of effector proteins. Phosphatidylinositol 4-phosphate (PI4P) is a derivative
36 of phosphatidylinositol, which is phosphorylated on the fourth position of its inositol head
37 group. Early studies in animal and yeast cells established PI4P as a master regulator of the Golgi
38 apparatus (De Matteis *et al.*, 2005). Indeed, many PI4P binding proteins localize to the Golgi.
39 Such as, for example, the oxysterol-binding protein (OSBP), four-phosphate-adaptor protein 1
40 (FAPP1), ceramide transfer (CERT), Golgi-localized gamma-ear-containing ADP-ribosylation
41 factor-binding proteins (GGAs) or Golgi-phosphoprotein3 (GOLPH3) (Santiago-Tirado &
42 Bretscher, 2011). In addition, PI4P is a key trafficking regulator of the Golgi and *trans*-Golgi
43 Network (TGN) thanks to its interplay with small GTPases from the Rab and Arf families
44 (Noack & Jaillais, 2017).

45
46 Although PI4P mainly accumulates in the Golgi in animal cells, its function extends beyond the
47 Golgi (De Matteis *et al.*, 2013). New generation of PI4P biosensors together with PI4P
48 immunolocalization established that PI4P also accumulates at the plasma membrane, in late
49 endosomes and in autophagosomes (De Matteis *et al.*, 2013). In plants, the situation is even
50 more extreme since PI4P accumulates as a major pool at the plasma membrane and is present

51 as a minor pool in the TGN (Noack & Jaillais, 2020). In this review, we present how PI4P
52 localization was studied in plants and what are the evidences for a massive accumulation of
53 PI4P at the cell surface by comparison to its intracellular localization. We describe the known
54 mechanisms that regulate this PI4P balance between the plasma membrane and the TGN.
55 Finally, we summarize some of the key functions of PI4P at the plasma membrane, including
56 many cellular pathways that are traditionally associated with PI(4,5)P₂ function in animal cells.
57

58 **An inverted balance of PI4P between the TGN and the plasma membrane in plants.**

59 The plekstrin homology (PH) domain of FAPP1 and OSBP interacts with PI4P and when fused
60 to a fluorescent protein, they almost exclusively localize to the Golgi in both animal and yeast
61 cells (Levine & Munro, 1998, 2002; Simon *et al.*, 2014). By contrast, expression of the same
62 domains in plant cells, including tabaco BY-2 cells and all plant cell types analyzed in
63 transgenic *Arabidopsis thaliana*, leads to a drastically different localization (Vermeer *et al.*,
64 2009; Simon *et al.*, 2014). Indeed, PH^{FAPP1} is targeted to the TGN but also to a greater extent
65 to the plasma membrane (Simon *et al.*, 2014). PI4P sensor with higher avidity, consisting of
66 two (2xPH) or three (3xPH) repeats of the FAPP1 PH domain, are even more associated with
67 the plasma membrane than 1xPH^{FAPP1}, which hints that the plasma membrane accumulates
68 more PI4P than the TGN membrane (Simon *et al.*, 2014, 2016). This hypothesis is supported
69 by additional observations. The PH domain of FAPP1 does not only bind to PI4P but also to
70 the small GTPase ARF1, which localizes at the TGN (Levine & Munro, 2002; He *et al.*, 2011).
71 However, the two binding sites are distinct, so it is possible to mutate the ARF1 binding site
72 while leaving the PI4P interacting domain fully functional (He *et al.*, 2011). Such mutant
73 strictly localizes to the plasma membrane in plants (Simon *et al.*, 2016; Ito *et al.*, 2021). By
74 contrast, mutants still able to interact to ARF1 but not PI4P localize to the TGN and not the
75 plasma membrane (Simon *et al.*, 2016). Thus, PI4P binding is dispensable for PH^{FAPP1} to
76 associate with the TGN, while it is necessary and sufficient for plasma membrane targeting.
77 Moreover, the P4M domain from the SiDM protein of *Legionella pneumophila*, which
78 specifically binds to PI4P but not ARF1 (Hammond *et al.*, 2014), is strictly localized at the
79 plasma membrane in all cell types analyzed (i.e., various root tissues, hypocotyl, leaf epidermis,
80 shoot apical meristem) in *Arabidopsis* and in *Nicotiana benthamiana* leaves (Simon *et al.*,
81 2016).
82

83 So far, only two endogenous *Arabidopsis* protein, PLANT U-BOX13 (PUB13) and
84 VASCULAR NETWORK DEFECTIVE3/ FORKED 1 (VAN3/FKD1) have been shown to
85 bind PI4P and to localize at the TGN (**Table 1**) (Koizumi *et al.*, 2005; Naramoto *et al.*, 2009;
86 Antignani *et al.*, 2015). Of note, both proteins also require interaction with TGN-associated
87 proteins for localization, confirming that PI4P binding alone may not be sufficient for TGN
88 targeting (Table 1). In addition, recent reports found these proteins to also localize at the plasma
89 membrane (Zhou *et al.*, 2018; Adamowski & Friml, 2021), again suggesting that PI4P binding
90 favors plasma membrane targeting in plants.
91

92 **Regulation of the PI4P balance between the plasma membrane and the TGN**

93 PI4K α 1 localize at the plasma membrane and is essential for both gametophyte and sporophyte
94 development, which is consistent with a prominent functional role of PI4P at the cell surface
95 (**Figure 1a**) (Noack *et al.*, 2022). PI4K β 1 and PI4K β 2 localize to the TGN (Preuss *et al.*, 2006).
96 These two kinases are redundant, with each single mutant presenting no phenotype, while the
97 double mutant shows root hair and cytokinesis defects, as well as innate immunity-related
98 phenotypes (Preuss *et al.*, 2006; Šašek *et al.*, 2014; Antignani *et al.*, 2015; Lin *et al.*, 2019). In
99 addition, electron tomography showed that the TGN of the *pi4k β 1 β 2* double mutant has an
100 aberrant morphology (Kang *et al.*, 2011). Together, the presence of these kinases at the TGN

101 and the phenotype of the double mutant argue in favor of PI4P being synthesized and present
102 in this compartment. Note that PI4K β 1 also localizes at the cell plate in dividing cells (Lin *et al.*,
103 2019), and that this compartment strongly accumulates PI4P sensors (Simon *et al.*, 2016).

104
105 The existence of a minor pool of PI4P at the TGN was also revealed thanks to the use of a
106 specific tool designed to deplete PI4P at the plasma membrane. This tool, coined MAP-SAC1,
107 consists of a plasma membrane targeting sequence (myristoylation and palmitoylation
108 sequence, MAP), fused with a fluorescent protein and the isolated catalytic domain of the yeast
109 PI4P phosphatase SAC1 (Simon *et al.*, 2016; Platre *et al.*, 2018; Doumane *et al.*, 2022). MAP-
110 SAC1 is specifically decreasing the amount of PI4P at the plasma membrane, but is not
111 impacting the pool of PI4P at the TGN. In this condition, PH^{FAPP1} and P4M^{SidM} relocate from
112 the plasma membrane to the TGNs, demonstrating that the two PI4P pools are in competition
113 to attract PI4P-binding proteins (Simon *et al.*, 2016; Platre *et al.*, 2018). Thus, it is likely that
114 two pools of PI4P exist in plant cells, a major one at the cell surface and a minor one in the
115 TGN.

116
117 This hypothesis is also supported by the fact that the PI4P balance between the TGN and the
118 plasma membrane is actively regulated by at least two pathways (Figure 1). The first one was
119 described in root hairs. In this cell type, the PI4P phosphatase SAC7, also called ROOT HAIR
120 DEFECTIVE4 (RHD4), localizes to the TGN (Thole *et al.*, 2008). In the *sac7* mutant, PI4P
121 sensors accumulate in intracellular compartments, likely corresponding to TGNs. This
122 phenotype is even stronger in *sac6sac7* and *sac7sac8* double mutants, SAC6 and SAC8 being
123 closely related to the SAC7 protein (Dubois & Jaillais, 2020; Song *et al.*, 2021). This balance
124 is also maintained by the activity of a phosphoinositide-specific phospholipase C (PI-PLC) at
125 the TGN (Ito *et al.*, 2021). Chemical inhibition of PI-PLCs increases the amount of PI4P sensors
126 in the TGN. In addition, inhibition of very long chain fatty acid synthesis, which mainly targets
127 sphingolipids, decreases the amount of PI-PLC2/7 at the TGN and thereby increases the amount
128 of PI4P in this compartment (Ito *et al.*, 2021). These results were confirmed by direct lipid
129 measurements using mass spectrometry, thereby confirming data obtained with translocation
130 PI4P sensors. Altogether, it appears that PI4P phosphatases and PI4P phospholipases cooperate
131 to actively maintain a low concentration of PI4P at the TGN membrane (Boutté & Jaillais,
132 2020).

133 134 **PI4P is essential for many cellular functions at the cell surface**

135 One possible explanation for the continuous and active removal of PI4P at the TGN in plants,
136 is to target PI4P-binding proteins to the plasma membrane. Indeed, PI4P not only accumulates
137 at the cell surface, but it is absolutely crucial for many cellular processes that are known to rely
138 on PI(4,5)P₂ in other eukaryotic systems.

139
140 **Plasma membrane electrostatics.** The plasma membrane is highly electronegative in yeast,
141 animal and plant cells, but these organisms use different anionic lipids to power the electrostatic
142 field of this compartment (Platre & Jaillais, 2017). In animals, PI(4,5)P₂ has a central role but
143 PI4P and PI(3,4,5)P₃ are redundantly involved in the regulation of membrane surface charges.
144 In yeast, phosphatidylserine is the main lipid that drives plasma membrane electrostatics. In
145 plants, the MAP-SAC1 assay, as well as pharmacological inhibition, suggest that PI4P is the
146 main driver of plasma membrane electrostatic (Simon *et al.*, 2016). Phosphatidic acid (PA),
147 phosphatidylserine and PI(4,5)P₂ contribute to the surface charges of the plant plasma
148 membrane but to a lesser extent (Barbosa *et al.*, 2016; Platre *et al.*, 2018; Synek *et al.*, 2021).

149

150 **ER/PM contact sites.** PI4P is also required for the localization of tethering proteins involved in
151 the establishment of membrane contact sites (Noack & Jaillais, 2020; Kang *et al.*, 2022) (**Table**
152 **1**). SYNAPTOTAGMIN1 (SYT1) is a tethering protein that is embedded in the ER membrane
153 via transmembrane domains and which contacts the plasma membrane through lipid-binding
154 C2 domains (Pérez-Sancho *et al.*, 2015) (**Figure 2a**). SYT1 is homologous to mammalian
155 extended-synaptotagmin (E-SYT), which require PI(4,5)P₂ for association with the plasma
156 membrane (Giordano *et al.*, 2013). By contrast, SYT1 localize in the bulk ER when coexpressed
157 with MAP-SAC1, suggesting that PI4P is critical for SYT1 plasma membrane association
158 (Ruiz-Lopez *et al.*, 2021). PI4P is also required for the localization of MULTIPLE C2
159 DOMAIN AND TRANSMEMBRANE PROTEINs (MCTP) (Brault *et al.*, 2019). MCTPs are
160 tethering the ER and the plasma membrane within plasmodesmata, which are membrane-lined
161 pores that connect plant cells (Kang *et al.*, 2022). Upon pharmacological inhibition of PI4P,
162 MCTP4 does not accumulate in plasmodesmata anymore (Brault *et al.*, 2019).

163
164 **Exocytosis and endocytosis.** The exocyst is an octameric protein complex that tethers transport
165 vesicles to the target membrane and represents one of the fundamental regulators of exocytosis
166 in eukaryotes (Žárský *et al.*, 2013). In yeast and animal cells, the exocyst recruitment to the
167 plasma membrane depends on the direct interaction of Exo70 subunit with PI(4,5)P₂ (Liu *et al.*,
168 2007). In *Arabidopsis*, the exocyst is targeted to the plasma membrane via the EXO70A1
169 subunit (Synek *et al.*, 2021). The plasma membrane localization of EXO70A1 was completely
170 lost upon the pharmacological depletion of PI4P, and diminished upon inhibition of PA
171 synthesis, while mutants with low levels of PI(4,5)P₂ had a minor impact on exocyst membrane
172 association (**Table 1**). Thus, unlike in yeasts and animals, PI4P seem to be more important than
173 PI(4,5)P₂ in exocyst membrane targeting in *Arabidopsis* root cells (**Figure 2b**).

174
175 A similar situation may also exist for clathrin-mediated endocytosis. Indeed, PI(4,5)P₂ is
176 absolutely required for endocytosis in mammalian cells and depletion of this lipid from the
177 plasma membrane completely and acutely inhibits endocytosis (He *et al.*, 2017). In plants,
178 PI(4,5)P₂ is also involved in endocytosis but depletion of PI(4,5)P₂ from the cell surface does
179 not fully block endocytosis (Doumane *et al.*, 2021). It is thus possible that other anionic lipids
180 may replace or compensate, at least partially, PI(4,5)P₂ during endocytosis. For example, the
181 μ-homology domain (μHD) of the TML protein interacts with PI4P *in vitro*, clusters PI4P
182 molecules *in silico* and depends on PI4P for plasma membrane localization (Yperman *et al.*,
183 2021) (**Table 1**). TML is part of the T-PLATE complex, which is essential for clathrin-mediated
184 endocytosis in plants (Gadeyne *et al.*, 2014). Several plant dynamin-related proteins have also
185 been shown to interact with PI4P *in vitro* (Kim *et al.*, 2001; Lam *et al.*, 2002; Lee *et al.*, 2002;
186 Bednarek & Backues, 2010; Minami *et al.*, 2015), however the exact roles of PI4P in dynamin
187 function at the plasma membrane remain to be fully evaluated in planta (**Table 1**).

188
189 **Signaling.** A number of proteins involved in cell signaling also rely on PI4P for plasma
190 membrane association (**Table 1**). This is the case for example of AGC kinases such a PINOID
191 or D6PKs and the Rho-GTPase of Plant6 (ROP6) (Barbosa *et al.*, 2016; Simon *et al.*, 2016;
192 Platre *et al.*, 2019). In addition, NON-PHOTOTROPIC HYPCOTYL3 (NPH3), a substrate of
193 the AGC kinase and blue photoreceptor Phot1, also localizes to the plasma membrane in a PI4P-
194 dependent manner (Reuter *et al.*, 2021) (**Figure 2c**). NPH3 association with the plasma
195 membrane is dynamically regulated by lights and is dependent on NPH3 phosphorylation and
196 14-3-3 association (Reuter *et al.*, 2021). BRI1 KINASE INHIBITOR1 (BKI1) and
197 MEMBRANE ASSOCIATED KINASE REGULATOR2 (MAKR2), which are receptor kinase
198 inhibitors, also localize to the plasma membrane in a PI4P-dependent manner (Jaillais *et al.*,
199 2011; Simon *et al.*, 2016). Like for NPH3, the association between BKI1/MAKR2 and the cell

200 surface is actively regulated by phosphorylation (Jaillais *et al.*, 2011; Marquès-Bueno *et al.*,
201 2021) (**Figure 2c**).

202 Moreover, using MAP-SAC1, PI4P was shown to be required for the plasma membrane
203 localization and function of ACTIVATED DISEASE RESISTANCE1 (ADR1) and
204 RESISTANCE TO P. SYRINGAE PV MACULICOLA 1 (RPM1). ADR1 and RPM1 are
205 immune receptors (Saile *et al.*, 2021). Interestingly, the effect is specific for PI4P, since
206 depletion of PI(4,5)P₂ using MAP-OCRL, a synthetic plasma membrane-targeted 5-phosphatase,
207 did not impact the localization of ADR1 or RPM1 (Doumane *et al.*, 2021; Saile *et al.*, 2021;
208 Doumane *et al.*, 2022). Thus, PI4P is not only essential for development but also for signaling
209 events critical to immune responses.

210
211 **Nanodomain formation.** The plasma membrane is not uniform and instead is made of a myriad
212 of domains that have specific lipid and protein compositions (Jaillais & Ott, 2020). REMORINs
213 (REM) are well characterized proteins that localize in nanodomains in a sterol dependent
214 manner. Their C-terminal anchor (REM-CA) binds to PI4P in vitro (Gronnier *et al.*, 2017)
215 (Table 1). In addition, the targeting of StREM1.3 is dependent on PI4P at the plasma membrane,
216 as demonstrated by the MAP-SAC1 assay. Thus, PI4P is not only required for the localization
217 and function of many proteins at the plasma membrane, but might also be important for their
218 targeting at the nanoscale (Jaillais & Ott, 2020) (**Figure 2d**). In addition, PI4K α 1 and its protein
219 partners are not uniformly localized at the plasma membrane, but are also present in
220 nanodomains (Noack *et al.*, 2022). However, how PI4K α 1-containing nanodomains relate to
221 the function of PI4P in organizing plasma membrane nanodomains remains to be fully
222 established.

223 224 **Conclusion**

225 PI4P has an unusual accumulation pattern in plants by comparison with animal and yeast cells.
226 It is seldomly found in the TGN, thanks to an active and local removal by phosphatases and
227 phospholipases. By contrast, most PI4P concentrate at the plasma membrane, where it serves
228 essential cellular functions, in term of structural organization, intracellular trafficking and
229 signaling. However, these functions are just beginning to emerge and it is likely that many more
230 PI4P-dependent phenomena will be discovered in the near future. Another challenge will be to
231 decipher whether the balance of PI4P between the TGN and the plasma membrane is
232 dynamically regulated by signaling pathways or during development. The exact function of
233 PI4P in the Golgi/TGN also remains poorly understood. Given the specificities of the plant
234 TGN (Kang *et al.*, 2022), it is likely that PI4P may have unique functions in this compartment.
235

236 **Acknowledgements**

237 We are grateful to the SiCE group (RDP, Lyon, France) for comments and discussions. YJ has
238 received funding from the European Research Council (ERC) under the European Union's
239 Horizon 2020 research and innovation program (Grant Agreement No 101001097).

240
241 The authors declare no competing interests.
242

Protein	Localization	Mode of binding	Lipid binding method	Function of the protein	References
PUB13	TGN, plasma membrane#	<i>PI4P binding via ARM domain and co-incident binding to RabA4B</i>	<i>In vitro</i> : lipid overlay assay, liposomes	signaling	Antignani <i>et al.</i> , 2015 Zhou <i>et al.</i> , 2018
VAN3/FKD1	TGN or plasma membrane*	<i>PI4P binding via PH domain and co-incident binding to CVP2 and VAB</i>	<i>In vitro</i> : lipid overlay assay, stimulation of enzymatic activity	membrane trafficking	Koizumi <i>et al.</i> , 2005 Naramoto <i>et al.</i> , 2009
DRP2A/ADL6	TGN, plasma membrane	<i>Anionic lipids binding via PH domain</i>	<i>In vitro</i> : lipid overlay assay, liposomes	membrane trafficking	Lam <i>et al.</i> , 2002
AGD1	Plasma membrane	<i>PI4P and other monophosphorylated phosphoinositides binding via PH domain</i>	<i>In vitro</i> : lipid overlay assay	membrane trafficking	Yoo <i>et al.</i> , 2012
SYT1	Plasma membrane	<i>Anionic lipids binding via C2 domain</i>	<i>In vivo</i> : MAP-SAC1, PAO <i>In vitro</i> : lipid overlay assay, liposome	membrane contact site	Ruiz-Lopez <i>et al.</i> , 2021 Perez-Sancho <i>et al.</i> , 2015
MCTP	Plasma membrane	<i>PI4P binding via C2 domain</i>	<i>In vivo</i> : PAO <i>In silico</i> : molecular dynamics	membrane contact site	Brault <i>et al.</i> , 2019
D6PK	Plasma membrane	<i>Anionic lipids binding via polybasic stretch</i>	<i>In vitro</i> : lipid overlay assay, liposomes <i>In vivo</i> : PAO	signaling	Barbosa <i>et al.</i> , 2016
EXO70A1	Plasma membrane	<i>Anionic lipids via polybasic surface + co-incident binding?</i>	<i>In vitro</i> : lipid overlay assay, liposomes <i>In vivo</i> : PAO <i>In silico</i> : molecular dynamics	membrane trafficking	Synek <i>et al.</i> , 2021
PID	Plasma membrane	<i>Anionic lipids via polybasic stretch in kinase insertion domain</i>	<i>In vivo</i> : PAO <i>In vitro</i> : lipid overlay assay	signaling	Simon <i>et al.</i> , 2016
ROP6	Plasma membrane	<i>Anionic lipids via polybasic C-terminal tail + Geranylgeranyl lipid anchor</i>	<i>In vitro</i> : lipid overlay assay <i>In vivo</i> : PAO	signaling	Platre <i>et al.</i> , 2019
NPH3	Plasma membrane	<i>Anionic lipids binding via a C-terminal amphipathic helix with a polybasic face</i>	<i>In vivo</i> : MAP-SAC1 <i>In vitro</i> : lipid overlay assay, liposomes	signaling	Reuter <i>et al.</i> , 2021
BK11	Plasma membrane	<i>Anionic lipids binding via polybasic stretch</i>	<i>In vitro</i> : lipid overlay assay <i>In vivo</i> : BK11 localization in WT and cho1Δ yeast	signaling	Jaillais <i>et al.</i> , 2011 Simon <i>et al.</i> , 2016,
MAKR2	Plasma membrane	<i>Anionic lipids binding via polybasic stretch</i>	<i>In vitro</i> : lipid overlay assay <i>In vivo</i> : PAO; MAKR2 localization in WT and cho1Δ yeast	signaling	Simon <i>et al.</i> , 2016 Marquès-Bueno <i>et al.</i> , 2021
MAKR1, MAKR3, MAKR4	Plasma membrane	<i>Anionic lipids binding via polybasic stretch</i>	<i>In vitro</i> : lipid overlay assay	signaling	Simon <i>et al.</i> , 2016
ADR1	Plasma membrane	<i>Anionic lipids binding via CC_R domain</i>	<i>In vivo</i> : MAP-SAC1 <i>In vitro</i> : lipid overlay assay	signaling	Saile <i>et al.</i> , 2021
RPM1	Plasma membrane	<i>Anionic lipids binding via CC domain</i>	<i>In vivo</i> : MAP-SAC1 <i>In vitro</i> : lipid overlay assay	signaling	Saile <i>et al.</i> , 2021
REM1.3	Plasma membrane	<i>Phosphoinositide binding via REM-C anchor</i>	<i>In vivo</i> : MAP-SAC1 <i>In vitro</i> : liposomes <i>In silico</i> : molecular dynamics	nanodomain formation	Gronnier <i>et al.</i> , 2017
TML	Plasma membrane	<i>PI4P binding and clustering via μHD domain</i>	<i>In vivo</i> : PAO <i>In vitro</i> : lipid overlay assay <i>In silico</i> : molecular dynamics	membrane trafficking	Yperman <i>et al.</i> , 2021
POL/PLL	Plasma membrane	<i>PI4P and PS and targeting to the plasma membrane via myristoylation anchoring</i>	<i>In vitro</i> : lipid overlay assay, regulation of enzymatic activity	signaling	Gagne <i>et al.</i> , 2010

PATL2	Plasma membrane	Anionic lipid binding via SEC14 domain	In vitro: lipid overlay assay	membrane trafficking	Suzuki <i>et al.</i> , 2016
PATL3	Plasma membrane	co-incident binding to EXO70A1	In vitro: lipid overlay assay	membrane trafficking	Wu <i>et al.</i> , 2017
DRP1A	Plasma membrane	unknown	In vitro: liposomes	membrane trafficking	Backues and Bednarek, 2010
DRP1E	Plasma membrane	unknown	In vitro: lipid overlay assay	membrane trafficking	Minami <i>et al.</i> , 2015

243
244
245
246
247
248
249
250
251
252
253
254
255
256
257
258
259
260
261
262
263
264
265
266
267
268

Table 1. Summary of Arabidopsis proteins known to rely on PI4P for localization or to interact with PI4P in vitro. PLANT U-BOX13 (PUB13), VASCULAR NETWORK DEFECTIVE3/ FORKED 1 (VAN3/FKD1), DYNAMIN-RELATED PROTEIN 2A/ ARABIDOPSIS DYNAMIN-LIKE PROTEIN 6 (DRP2A/ADL6), ARF-GAP DOMAIN 1 (AGD1), SYNAPTOTAGMIN 1 (SYT1), MULTIPLE C2 DOMAIN AND TRANSMEMBRANE REGION PROTEIN (MCTP), D6 PROTEIN KINASE (D6PK), EXOCYST SUBUNIT EXO70 FAMILY PROTEIN A1 (EXO70A1), PINOID (PID), RHO-RELATED PROTEIN FROM PLANTS 6 (ROP6), NON-PHOTOTROPIC HYPOCOTYL 3 (NPH3), BRI1 KINASE INHIBITOR 1 (BKI1), MEMBRANE-ASSOCIATED KINASE REGULATOR 1-4 (MAKR1-4), ADR1 ACTIVATED DISEASE RESISTANCE 1 (ADR1), RESISTANCE TO PSEUDOMONAS SYRINGAE PV MACULICOLA 1 (RPM1), REMORIN 1.3 (REM1.3), TPLATE COMPLEX MUNISCIN-LIKE (TML), POLTERGEIST/POLTERGEIST LIKE 1 (POL/PLL1), PATELLIN 1-3 (PATL1-3), DYNAMIN-RELATED PROTEIN 1A (DRP1A), DYNAMIN-RELATED PROTEIN 1E (DRP1E), ARMADILLO DOMAIN (ARM), PLECKSTRIN HOMOLOGY DOMAIN (PH), SECRETION 14 DOMAIN (SEC14), COILED-COIL DOMAIN (CC), HELO/RPW8-LIKE COILED-COIL DOMAIN (CC_R), REMORIN C-TERMINAL ANCHOR (REM-CA), M-HOMOLOGY DOMAIN (μ HD).

PUB13 was found at the TGN in Antignani *et al.*, 2015 but Zhou *et al.*, (2018) also reported a plasma membrane localization.

* VAN3 was shown to localize in the TGN by Naramoto *et al.*, 2009 but was found at the plasma membrane in a recent preprint (Adamowski and Friml, 2021).

269 **References:**

- 270 **Adamowski M, Friml J. 2021.** Clathrin-mediated endocytosis-independent function of
271 VAN3 ARF GTPase Activating Protein at the plasma membrane. *bioRxiv*
272 11.17.468981; doi: <https://doi.org/10.1101/2021.11.17.468981>
- 273 **Antignani V, Klocko AL, Bak G, Chandrasekaran SD, Dunivin T, Nielsen E. 2015.**
274 Recruitment of PLANT U-BOX13 and the PI4K β 1/ β 2 Phosphatidylinositol-4 Kinases by the
275 Small GTPase RabA4B Plays Important Roles during Salicylic Acid-Mediated Plant Defense
276 Signaling in Arabidopsis. *The Plant Cell* **27**:243-61.
- 277 **Barbosa ICR, Shikata H, Zourelidou M, Heilmann M, Heilmann I, Schwechheimer C.**
278 **2016.** Phospholipid composition and a polybasic motif determine D6 PROTEIN KINASE
279 polar association with the plasma membrane and tropic responses. *Development* **143**: 4687-
280 4700.
- 281 **Bednarek SY, Backues SK. 2010.** Plant dynamin-related protein families DRP1 and DRP2
282 in plant development. *Biochemical Society transactions* **38**: 797–806.
- 283 **Boutté Y, Jaillais Y. 2020.** Metabolic Cellular Communications: Feedback Mechanisms
284 between Membrane Lipid Homeostasis and Plant Development. *Developmental Cell* **54**:171-
285 182.
- 286 **Braut ML, Petit JD, Immel F, Nicolas WJ, Glavier M, Brocard L, Gaston A, Fouché M,**
287 **Hawkins TJ, Crowet J-M, et al. 2019.** Multiple C2 domains and transmembrane region
288 proteins (MCTPs) tether membranes at plasmodesmata. *EMBO reports* **20**: e47182.
- 289 **De Matteis MA, Di Campli A, Godi A. 2005.** The role of the phosphoinositides at the Golgi
290 complex. *Biochimica et Biophysica Acta (BBA) - Molecular Cell Research* **1744**: 396–405.
- 291 **De Matteis MA, Wilson C, D'Angelo G. 2013.** Phosphatidylinositol-4-phosphate: The Golgi
292 and beyond. *BioEssays* **35**: 612–622.
- 293 **Doumane M, Caillaud M-C, Jaillais Y. 2022.** Experimental manipulation of
294 phosphoinositide lipids: from cells to organisms. *Trends in Cell Biology* **32**: 445-461.
- 295 **Doumane M, Lebecq A, Colin L, Fangain A, Stevens FD, Bareille J, Hamant O,**
296 **Belkhadir Y, Munnik T, Jaillais Y, et al. 2021.** Inducible depletion of PI(4,5)P₂ by the
297 synthetic iDePP system in Arabidopsis. *Nature Plants* **7**: 585-597.
- 298 **Dubois GA, Jaillais Y. 2020.** Anionic phospholipid gradients: an uncharacterized frontier of
299 the plant endomembrane network. *Plant Physiology* **185**: 577-592.
- 300 **Gadeyne A, Sánchez-Rodríguez C, Vanneste S, Di Rubbo S, Zauber H, Vanneste K, Van**
301 **Leene J, De Winne N, Eeckhout D, Persiau G, et al. 2014.** The TPLATE Adaptor Complex
302 Drives Clathrin-Mediated Endocytosis in Plants. *Cell* **156**: 691–704.
- 303 **Gagne JM, Clark SE. 2010.** The Arabidopsis stem cell factor POLTERGEIST is membrane
304 localized and phospholipid stimulated. *The Plant Cell* **22**: 729–743.
- 305 **Giordano F, Saheki Y, Idevall-Hagren O, Colombo SF, Pirruccello M, Milosevic I,**
306 **Gracheva EO, Bagriantsev SN, Borgese N, Pietro De Camilli. 2013.** PI(4,5)P₂-Dependent
307 and Ca²⁺-Regulated ER-PM Interactions Mediated by the Extended Synaptotagmins. *Cell*
308 **153**: 1494–1509.
- 309 **Gronnier J, Crowet J-M, Habenstein B, Nasir MN, Bayle V, Hosity E, Platre MP,**
310 **Gouguet P, Raffaele S, Martinez D, et al. 2017.** Structural basis for plant plasma membrane
311 protein dynamics and organization into functional nanodomains. *eLife* **6**: e26404.
- 312 **Hammond GRV, Machner MP, Balla T. 2014.** A novel probe for phosphatidylinositol 4-
313 phosphate reveals multiple pools beyond the Golgi. *The Journal of Cell Biology* **183**: 1061.
- 314 **He K, Marsland R, Upadhyayula S, Song E, Dang S, Capraro BR, Wang W, Skillern W,**
315 **Gaudin R, Ma M, et al. 2017.** Dynamics of phosphoinositide conversion in clathrin-mediated
316 endocytic traffic. *Nature*: **552**: 410-414.
- 317 **He J, Scott JL, Heroux A, Roy S, Lenoir M, Overduin M, Stahelin RV, Kutateladze TG.**

318 **2011.** Molecular basis of phosphatidylinositol 4-phosphate and ARF1 GTPase recognition by
319 the FAPP1 pleckstrin homology (PH) domain. *The Journal of biological chemistry* **286**:
320 18650–18657.

321 **Ito Y, Esnay N, Platre MP, Wattelet-Boyer V, Noack LC, Fougère L, Menzel W,**
322 **Claverol S, Fouillen L, Moreau P, et al. 2021.** Sphingolipids mediate polar sorting of PIN2
323 through phosphoinositide consumption at the trans-Golgi network. *Nature Communications*
324 **12**(1): 4267.

325 **Jaillais Y, Hothorn M, Belkhadir Y, Dabi T, Nimchuk ZL, Meyerowitz EM, Chory J.**
326 **2011.** Tyrosine phosphorylation controls brassinosteroid receptor activation by triggering
327 membrane release of its kinase inhibitor. *Genes and development* **25**: 232–237.

328 **Jaillais Y, Ott T. 2020.** The Nanoscale Organization of the Plasma Membrane and Its
329 Importance in Signaling: A Proteolipid Perspective. *Plant Physiology* **182**: 1682–1696.

330 **Kang B-H, Anderson CT, Arimura S-I, Bayer E, Bezanilla M, Botella MA, Brandizzi F,**
331 **Burch-Smith TM, Chapman KD, Dünser K, et al. 2022.** A glossary of plant cell structures:
332 Current insights and future questions. *The Plant Cell* **34**: 10-52.

333 **Kang B-H, Nielsen E, Preuss ML, Mastronarde D, Staehelin LA. 2011.** Electron
334 Tomography of RabA4b- and PI-4K β 1-Labeled Trans Golgi Network Compartments in
335 Arabidopsis. *Traffic* **12**: 313–329.

336 **Kim Y-W, Park D-S, Park S-C, Kim SH, Cheong G-W, Hwang I. 2001.** Arabidopsis
337 Dynamin-Like 2 That Binds Specifically to Phosphatidylinositol 4-Phosphate Assembles into
338 a High-Molecular Weight Complex in Vivo and in Vitro. *Plant Physiology* **127**: 1243–1255.

339 **Koizumi K, Naramoto S, Sawa S, Yahara N, Ueda T, Nakano A, Sugiyama M, Fukuda**
340 **H. 2005.** VAN3 ARF-GAP-mediated vesicle transport is involved in leaf vascular network
341 formation. *Development* **132**: 1699–1711.

342 **Lam BC-H, Sage TL, Bianchi F, Blumwald E. 2002.** Regulation of ADL6 activity by its
343 associated molecular network. *The Plant Journal* **31**: 565–576.

344 **Lee SH, Jin JB, Song J, Min MK, Park DS, Kim Y-W, Hwang I. 2002.** The Intermolecular
345 Interaction between the PH Domain and the C-terminal Domain of Arabidopsis Dynamin-like
346 6 Determines Lipid Binding Specificity. *Journal of Biological Chemistry* **277**: 31842–31849.

347 **Levine TP, Munro S. 1998.** The pleckstrin homology domain of oxysterol-binding protein
348 recognises a determinant specific to Golgi membranes. *Current Biology* **8**: 729–739.

349 **Levine TP, Munro S. 2002.** Targeting of Golgi-specific pleckstrin homology domains
350 involves both PtdIns 4-kinase-dependent and -independent components. *Current Biology* **12**:
351 695–704.

352 **Lin F, Krishnamoorthy P, Schubert V, Hause G, Heilmann M, Heilmann I. 2019.** A dual
353 role for cell plate-associated PI4K β in endocytosis and phragmoplast dynamics during plant
354 somatic cytokinesis. *The EMBO Journal* **38**: 543.

355 **Liu J, Zuo X, Yue P, Guo W. 2007.** Phosphatidylinositol 4,5-Bisphosphate Mediates the
356 Targeting of the Exocyst to the Plasma Membrane for Exocytosis in Mammalian Cells.
357 *Molecular Biology of the Cell* **18**: 4483–4492.

358 **Marquès-Bueno MM, Armengot L, Noack LC, Bareille J, Rodriguez L, Platre MP,**
359 **Bayle V, Liu M, Opdenacker D, Vanneste S, et al. 2020.** Auxin-Regulated Reversible
360 Inhibition of TMK1 Signaling by MAKR2 Modulates the Dynamics of Root Gravitropism.
361 *Current Biology* **31**: 228–237.

362 **Minami A, Tominaga Y, Furuto A, Kondo M, Kawamura Y, Uemura M. 2015.**
363 Arabidopsis dynamin-related protein 1E in sphingolipid-enriched plasma membrane domains
364 is associated with the development of freezing tolerance. *The Plant Journal* **83**: 501–514.

365 **Naramoto S, Sawa S, Koizumi K, Uemura T, Ueda T, Friml J, Nakano A, Fukuda H.**
366 **2009.** Phosphoinositide-dependent regulation of VAN3 ARF-GAP localization and activity
367 essential for vascular tissue continuity in plants. *Development* **136**: 1529–1538.

368 **Noack LC, Bayle V, Armengot L, Rozier F, Mamode-Cassim A, Stevens FD, Caillaud**
369 **M-C, Munnik T, Mongrand S, Pleskot R, et al. 2022.** A nanodomain-anchored scaffolding
370 complex is required for the function and localization of phosphatidylinositol 4-kinase alpha in
371 plants. *The Plant Cell* **34**: 302–332.

372 **Noack LC, Jaillais Y. 2017.** Precision targeting by phosphoinositides: how PIs direct
373 endomembrane trafficking in plants. *Current Opinion in Plant Biology* **40**: 22–33.

374 **Noack LC, Jaillais Y. 2020.** Functions of Anionic Lipids in Plants. *Annual Review of Plant*
375 *Biology* **71**: 71–102.

376 **Pérez-Sancho J, Vanneste S, Lee E, McFarlane HE, Esteban del Valle A, Valpuesta V,**
377 **Friml J, Botella MA, Rosado A. 2015.** The Arabidopsis Synaptotagmin1 Is Enriched in
378 Endoplasmic Reticulum-Plasma Membrane Contact Sites and Confers Cellular Resistance to
379 Mechanical Stresses. *Plant Physiology* **168**: 132–143.

380 **Platre MP, Bayle V, Armengot L, Bareille J, Marquès-Bueno MDM, Creff A, Maneta-**
381 **Peyret L, Fiche J-B, Nöllmann M, Miège C, et al. 2019.** Developmental control of plant
382 Rho GTPase nano-organization by the lipid phosphatidylserine. *Science* **364**: 57–62.

383 **Platre MP, Jaillais Y. 2017.** Anionic lipids and the maintenance of membrane electrostatics
384 in eukaryotes. *Plant Signaling & Behavior* **12**: e1282022.

385 **Platre MP, Noack LC, Doumane M, Bayle V, Simon MLA, Maneta-Peyret L, Fouillen L,**
386 **Stanislas T, Armengot L, Pejchar P, et al. 2018.** A Combinatorial Lipid Code Shapes the
387 Electrostatic Landscape of Plant Endomembranes. *Developmental Cell*: **45**: 465–480.

388 **Preuss ML, Schmitz AJ, Thole JM, Bonner HKS, Otegui MS, Nielsen E. 2006.** A role for
389 the RabA4b effector protein PI-4K 1 in polarized expansion of root hair cells in *Arabidopsis*
390 *thaliana*. *The Journal of Cell Biology* **172**: 991–998.

391 **Reuter L, Schmidt T, Manishankar P, Throm C, Keicher J, Bock A, Droste-Borel I,**
392 **Oecking C. 2021.** Light-triggered and phosphorylation-dependent 14-3-3 association with
393 NON-PHOTOTROPIC HYPOCOTYL 3 is required for hypocotyl phototropism. *Nature*
394 *Communications* **12**: 6128.

395 **Ruiz-Lopez N, Pérez-Sancho J, Esteban del Valle A, Haslam RP, Vanneste S, Catalá R,**
396 **Perea-Resa C, Van Damme D, García-Hernández S, Albert A, et al. 2021.**
397 Synaptotagmins at the endoplasmic reticulum-plasma membrane contact sites maintain
398 diacylglycerol homeostasis during abiotic stress. *The Plant Cell* **33**: 2431–2453.

399 **Saile SC, Ackermann FM, Sunil S, Keicher J, Bayless A, Bonardi V, Wan L, Doumane**
400 **M, Stöbbe E, Jaillais Y, et al. 2021.** Arabidopsis ADR1 helper NLR immune receptors
401 localize and function at the plasma membrane in a phospholipid dependent manner. *The New*
402 *phytologist* **232**(6): 2440–2456.

403 **Santiago-Tirado FH, Bretscher A. 2011.** Membrane-trafficking sorting hubs: cooperation
404 between PI4P and small GTPases at the trans-Golgi network. *Trends in Cell Biology* **21**: 515–
405 525.

406 **Šašek V, Janda M, Delage E, Puyaubert J, Guivarc’h A, López Maseda E, Dobrev PI,**
407 **Caius J, Bóka K, Valentová O, et al. 2014.** Constitutive salicylic acid accumulation in
408 pi4kIIIβ1β2 Arabidopsis plants stunts rosette but not root growth. *The New phytologist* **203**:
409 805–816.

410 **Simon MLA, Platre MP, Assil S, van Wijk R, Chen WY, Chory J, Dreux M, Munnik T,**
411 **Jaillais Y. 2014.** A multi-colour/multi-affinity marker set to visualize phosphoinositide
412 dynamics in Arabidopsis. *The Plant Journal* **77**: 322–337.

413 **Simon MLA, Platre MP, Marquès-Bueno MM, Armengot L, Stanislas T, Bayle V,**
414 **Caillaud M-C, Jaillais Y. 2016.** A PtdIns(4)P-driven electrostatic field controls cell
415 membrane identity and signalling in plants. *Nature Plants* **2**: 16089.

416 **Song L, Wang Y, Guo Z, Lam SM, Shui G, Cheng Y. 2021.** NCP2/RHD4/SAC7, SAC6
417 and SAC8 phosphoinositide phosphatases are required for PtdIns4P and PtdIns(4,5)P₂

418 homeostasis and Arabidopsis development. *The New phytologist* **231**: 713–725.

419 **Suzuki T, Matsushima C, Nishimura S, Higashiyama T, Sasabe M, Machida Y. 2016.**

420 Identification of Phosphoinositide-Binding Protein PATELLIN2 as a Substrate of Arabidopsis

421 MPK4 MAP Kinase during Septum Formation in Cytokinesis. *Plant and Cell Physiology* **57**:

422 1744–1755.

423 **Synek L, Pleskot R, Sekereš J, Serrano N, Vukašinović N, Ortmannová J, Klejchová M,**

424 **Pejchar P, Batystová K, Gutkowska M, et al. 2021.** Plasma membrane phospholipid

425 signature recruits the plant exocyst complex via the EXO70A1 subunit. *Proceedings of the*

426 *National Academy of Sciences* **118**: e2105287118.

427 **Thole JM, Vermeer JEM, Zhang Y, Gadella TWJ, Nielsen E. 2008.** ROOT HAIR

428 DEFECTIVE4 Encodes a Phosphatidylinositol-4-Phosphate Phosphatase Required for Proper

429 Root Hair Development in *Arabidopsis thaliana*. *The Plant Cell* **20**: 381–395.

430 **Vermeer JEM, Thole JM, Goedhart J, Nielsen E, Munnik T, Gadella TWJ. 2009.**

431 Imaging phosphatidylinositol 4-phosphate dynamics in living plant cells. *The Plant journal*

432 **57**: 356–372.

433 **Wu C, Tan L, van Hooren M, Tan X, Liu F, Li Y, Zhao Y, Li B, Rui Q, Munnik T, et al.**

434 **2017.** Arabidopsis EXO70A1 recruits Patellin3 to the cell membrane independent of its role as

435 an exocyst subunit. *Journal of integrative plant biology* **59**: 851–865.

436 **Yoo C-M, Quan L, Cannon AE, Wen J, Blancaflor EB. 2011.** AGD1, a class 1 ARF-GAP,

437 acts in common signaling pathways with phosphoinositide metabolism and the actin

438 cytoskeleton in controlling Arabidopsis root hair polarity. *The Plant Journal* **69**: 1064–1076.

439 **Yperman K, Wang J, Eeckhout D, Winkler J, Vu LD, Vandorpe M, Gronès P, Mylle E,**

440 **Kraus M, Merceron R, et al. 2021.** Molecular architecture of the endocytic TPLATE

441 complex. *Science Advances* **7**: eabe7999.

442 **Žárský V, Kulich I, Fendrych M, Pečenková T. 2013.** Exocyst complexes multiple

443 functions in plant cells secretory pathways. *Current Opinion in Plant Biology* **16**: 726–733.

444 **Zhou J, Liu D, Wang P, Ma X, Lin W, Chen S, Mishev K, Lu D, Kumar R, Vanhoutte I,**

445 **et al. 2018.** Regulation of Arabidopsis brassinosteroid receptor BRI1 endocytosis and

446 degradation by plant U-box PUB12/PUB13-mediated ubiquitination. *Proceedings of the*

447 *National Academy of Sciences of the United States of America* **115**: E1906–E1915.



448 **Figure legends:**

449


450 **Figure 1. Enzymes involved in balancing PI4P between the plasma membrane and the**
451 **TGN.** In *Arabidopsis thaliana*, PI4P is produced at the plasma membrane by the PI4K α 1
452 complex (shown with a higher magnification in a) and in the TGN by the redundant PI4K β 1
453 and PI4K β 2 (shown with a higher magnification in b). The concentration of PI4P in the TGN
454 is maintained at a low level thanks to the continuous action of PI-PLC2/7 and the 4-
455 phosphatases SAC7. PI, phosphatidylinositol; PI4P, phosphatidylinositol 4-phosphate; DAG,
456 diacylglycerol; IP₂, inositol (1,4)-biphosphate; PI-PLC, phosphoinositide-specific
457 phospholipase C; SAC7, SAC-DOMAIN CONTAINING7; TGN, trans-Golgi Network; EE,
458 early endosome; ALPS motif, amphipathic lipid packing sensor motif; HYC, HYCCIN-
459 CONTAINING; EFOP, EFR3-OF-PLANTS; NPG, NO POLLEN GERMINATION; NPGR,
460 NPG-RELATED; PI4K, phosphatidylinositol 4-kinase.

461


462 **Figure 2: Examples of key PI4P function at the plasma membrane.** Examples illustrating
463 some of the key functions of PI4P at the plant plasma membrane: in membrane contact site
464 establishment (a), exocytosis (b), signaling (c) and the formation of nanodomains (d). Note that
465 most of these cellular functions are mainly driven by PI(4,5)P₂ in animals, while they are
466 mediated, at least in part, by PI4P in plants. The examples shown in this figure were investigated
467 in *Arabidopsis thaliana* (a-c) and *Nicotiana benthamiana* leaf (a and d). PI4P,
468 phosphatidylinositol 4-phosphate; PI(4,5)P₂, phosphatidylinositol 4,5-bisphosphate; PA,
469 phosphatidic acid; ER, endoplasmic reticulum; SYT, SYNAPTOTAGMIN; TM,
470 transmembrane region; SMP, Synaptotagmin-like mitochondrial-lipid binding domain; C2, C2
471 domain; StREM1.3, *Solanum tuberosum* REMORIN1.3; MAKR2, MEMBRANE
472 ASSOCIATED KINASE REGULATOR2; TMK1, RECEPTOR-LIKE TRANSMEMBRANE
473 KINASE1; NPH3, NON-PHOTOTROPIC HYPOCOTYL3.




a) Membrane contact site



b) Exocytosis



c) Plasma membrane/cytoplasm switch signaling



d) Nanodomain formation

